

McsUsbNet.dll Version 5.1.15

Multi Channel Systems MCS GmbH
Aspenhaustrasse 21
72770 Reutlingen
Germany
Fon +49-71 21-90 92 5 - 0
Fax +49-71 21-90 92 5 -11
info@multichannelsystems.com
www.multichannelsystems.com

Generated by Doxygen 1.8.18

1 M	IcsUsbNet.dll for MCS USB devices	1
	1.1 Introduction	1
	1.2 System requirements	1
	1.3 Connecting to an MCS device	1
2 D	Device Classes	2
	2.1 The MCS FluidControl Device	2
	2.1.1 Introduction	2
	2.1.2 Access to the FluidControl device	2
	2.2 MCS-USB-Sw2to64 device	3
3 F	unction Classes	3
4 D	Pata ACQuisition (DACQ) Devices	4
5 T	he MCS Robo Device	5
	5.1 Introduction	5
6 S	TG200x & STG400x STimulus Generator	5
	6.1 Introduction	5
	6.2 Download mode	6
	6.2.1 Memory Layout and Trigger Setup	6
	6.3 Streaming mode	8
	6.3.1 Memory Layout and Trigger Setup	9
7 N	lamespace Index	11
	7.1 Namespace List	11
8 H	lierarchical Index	11
	8.1 Class Hierarchy	11
9 C	Class Index	16
	9.1 Class List	16
10 I	Namespace Documentation	22
	10.1 Mcs Namespace Reference	22
	10.2 Mcs::Usb Namespace Reference	22
	10.2.1 Enumeration Type Documentation	27
	10.2.2 Function Documentation	27
11 (Class Documentation	29
	11.1 CW2100_FunctionNet::AudioChannelsNet Struct Reference	29
	11.1.1 Member Data Documentation	29
	11.2 BatteryState Class Reference	29
	11.2.1 Property Documentation	30
	11.3 BesselFilterHighPassNet Class Reference	30

11.3.1 Constructor & Destructor Documentation	30
11.4 BesselFilterLowPassNet Class Reference	31
11.4.1 Constructor & Destructor Documentation	31
11.5 ButterworthFilterHighPassNet Class Reference	31
11.5.1 Constructor & Destructor Documentation	32
11.6 ButterworthFilterLowPassNet Class Reference	32
11.6.1 Constructor & Destructor Documentation	32
11.7 CChannelTestDeviceNet Class Reference	33
11.7.1 Constructor & Destructor Documentation	33
11.7.2 Member Function Documentation	33
11.8 CCMOSMea_FunctionNet Class Reference	34
11.8.1 Constructor & Destructor Documentation	36
11.8.2 Member Function Documentation	36
11.9 CCMOSMeaDeviceNet Class Reference	45
11.9.1 Constructor & Destructor Documentation	46
11.9.2 Member Function Documentation	46
11.9.3 Property Documentation	47
11.10 CCreateFilterNet Class Reference	48
11.10.1 Constructor & Destructor Documentation	48
11.10.2 Member Function Documentation	49
11.10.3 Property Documentation	49
11.11 CDacCalibrationFunctionNet Class Reference	50
11.11.1 Detailed Description	50
11.11.2 Constructor & Destructor Documentation	50
11.11.3 Member Function Documentation	51
11.12 CDacqGroupChannelGenericSelectionNet Class Reference	52
11.12.1 Constructor & Destructor Documentation	52
11.13 CDacqGroupChannelSelectionNet Class Reference	52
11.13.1 Constructor & Destructor Documentation	53
11.14 CDacqGroupChannelSelectionTemplateNet< DacqGroupChannelEnumTemplateNet, Dacq⇔	
	53
2.	54
	54
	56
	56
11.16 CDeviceGroupChannelInfoMEA2100_256Net Class Reference	56
	57
	57
11.17.1 Constructor & Destructor Documentation	57
•	57
11.18.1 Constructor & Destructor Documentation	58
11.19 CDeviceGroupChannelInfoTemplateNet< DacqGroupChannelEnumTemplateNet > Class Template Reference	58

	11.19.1 Constructor & Destructor Documentation	58
	11.19.2 Member Data Documentation	58
11.20	CDeviceGroupChannelInfoW2100Net Class Reference	59
	11.20.1 Constructor & Destructor Documentation	59
11.21	CDigOutStimulatorFunctionNet Class Reference	59
	11.21.1 Detailed Description	60
	11.21.2 Constructor & Destructor Documentation	60
	11.21.3 Member Function Documentation	61
11.22	CEncapsulatorDeviceNet Class Reference	64
	11.22.1 Detailed Description	64
	11.22.2 Constructor & Destructor Documentation	64
	11.22.3 Member Function Documentation	64
11.23	CExternDTesterDeviceNet Class Reference	64
	11.23.1 Detailed Description	65
	11.23.2 Constructor & Destructor Documentation	65
	11.23.3 Member Function Documentation	65
11.24	CFilterCoefficientsNet Class Reference	66
	11.24.1 Constructor & Destructor Documentation	67
	11.24.2 Member Function Documentation	68
	11.24.3 Property Documentation	68
11.25	CFilterConfigurationNet Class Reference	68
	11.25.1 Constructor & Destructor Documentation	69
	11.25.2 Member Function Documentation	69
11.26	CFilterConfigurationRegisterNet Class Reference	70
	11.26.1 Constructor & Destructor Documentation	71
	11.26.2 Member Function Documentation	71
11.27	CFilterPropertyNet Class Reference	72
	11.27.1 Constructor & Destructor Documentation	73
	11.27.2 Member Function Documentation	73
	11.27.3 Property Documentation	73
11.28	CFluidControlDeviceNet Class Reference	74
	11.28.1 Detailed Description	75
	11.28.2 Constructor & Destructor Documentation	75
	11.28.3 Member Function Documentation	75
	11.28.4 Property Documentation	80
11.29	CFYIDeviceNet Class Reference	80
	11.29.1 Detailed Description	80
	11.29.2 Constructor & Destructor Documentation	80
	11.29.3 Property Documentation	81
11.30	CGenericDevelopDeviceNet Class Reference	81
	11.30.1 Detailed Description	88
	11.30.2 Constructor & Destructor Documentation	88

11.30.3 Member Function Documentation	88
11.31 CGilsonDeviceNet Class Reference	98
11.31.1 Detailed Description	99
11.31.2 Constructor & Destructor Documentation	99
11.31.3 Member Function Documentation	99
11.31.4 Member Data Documentation	00
11.32 CGrapheneFunctionNet Class Reference	00
11.32.1 Detailed Description	02
11.32.2 Constructor & Destructor Documentation	02
11.32.3 Member Function Documentation	02
11.33 CHiClampDeviceNet Class Reference	09
11.33.1 Detailed Description	10
11.33.2 Constructor & Destructor Documentation	10
11.33.3 Property Documentation	10
11.34 CHLADacqNet Class Reference	10
11.34.1 Constructor & Destructor Documentation	11
11.35 CHLADeviceNet Class Reference	11
11.35.1 Detailed Description	11
11.35.2 Constructor & Destructor Documentation	11
11.35.3 Property Documentation	12
11.36 CMcsUsbDacqNet::CHWInfo Class Reference	12
11.36.1 Detailed Description	12
11.36.2 Constructor & Destructor Documentation	13
11.36.3 Member Function Documentation	13
11.37 CIntanMea_FunctionNet Class Reference	15
11.37.1 Constructor & Destructor Documentation	15
11.37.2 Member Function Documentation	16
11.38 CInterfaceboard2FunctionNet Class Reference	17
11.38.1 Detailed Description	17
11.38.2 Constructor & Destructor Documentation	18
11.38.3 Member Function Documentation	18
11.39 CInterfaceboardFunctionNet Class Reference	19
11.39.1 Detailed Description	19
11.39.2 Constructor & Destructor Documentation	19
11.39.3 Member Function Documentation	20
11.40 CLIH3DeviceNet Class Reference	21
11.40.1 Detailed Description	23
11.40.2 Constructor & Destructor Documentation	23
11.40.3 Member Function Documentation	23
11.40.4 Property Documentation	30
11.41 CMcsBus_AxisParametersNet Class Reference	30
11.41.1 Constructor & Destructor Documentation	31

11.41.2 Member Function Documentation
11.42 CMcsBus_ExtensionNet Class Reference
11.42.1 Constructor & Destructor Documentation
11.42.2 Member Function Documentation
11.43 CMcsBus_FYIExtensionNet Class Reference
11.43.1 Constructor & Destructor Documentation
11.43.2 Member Function Documentation
11.44 CMcsBus_MotorControlNet Class Reference
11.44.1 Constructor & Destructor Documentation
11.44.2 Member Function Documentation
11.45 CMcsBus_SensorNet Class Reference
11.45.1 Constructor & Destructor Documentation
11.45.2 Member Function Documentation
11.46 CMcsBus_TempSensorNet Class Reference
11.46.1 Constructor & Destructor Documentation
11.46.2 Member Function Documentation
11.47 CMcsBus_VoltageModeNet Class Reference
11.47.1 Constructor & Destructor Documentation
11.47.2 Member Function Documentation
11.48 CMcsBusNet Class Reference
11.48.1 Constructor & Destructor Documentation
11.48.2 Member Function Documentation
11.49 CMcsUsbDacqNet Class Reference
11.49.1 Detailed Description
11.49.2 Constructor & Destructor Documentation
11.49.3 Member Function Documentation
11.49.4 Member Data Documentation
11.49.5 Property Documentation
11.49.6 Event Documentation
11.50 CMcsUsbDeviceStatePushFunctionNet Class Reference
11.50.1 Constructor & Destructor Documentation
11.50.2 Member Function Documentation
11.50.3 Event Documentation
11.51 CMcsUsbDeviceStatePushNet Class Reference
11.51.1 Constructor & Destructor Documentation
11.51.2 Member Function Documentation
11.51.3 Event Documentation
11.52 CMcsUsbFactoryNet Class Reference
11.52.1 Constructor & Destructor Documentation
11.52.2 Member Function Documentation
11.52.3 Member Data Documentation
11.53 CMcsUsbFunctionNet Class Reference

11.53.1 Constructor & Destructor Documentation
11.53.2 Member Function Documentation
11.53.3 Member Data Documentation
11.54 CMcsUsbFunctionPointerContainer Class Reference
11.55 CMcsUsbListEntryNet Class Reference
11.55.1 Detailed Description
11.55.2 Constructor & Destructor Documentation
11.55.3 Member Function Documentation
11.55.4 Property Documentation
11.56 CMcsUsbListNet Class Reference
11.56.1 Detailed Description
11.56.2 Constructor & Destructor Documentation
11.56.3 Member Function Documentation
11.56.4 Property Documentation
11.56.5 Event Documentation
11.57 CMcsUsbNet Class Reference
11.57.1 Detailed Description
11.57.2 Constructor & Destructor Documentation
11.57.3 Member Function Documentation
11.57.4 Member Data Documentation
11.57.5 Property Documentation
11.58 CMcsUsbPointerContainer Class Reference
11.59 CMEA2100_256DacqGroupChannelSelectionNet Class Reference
11.59.1 Constructor & Destructor Documentation
11.60 CMEA2100x256FunctionNet Class Reference
11.60.1 Detailed Description
11.60.2 Constructor & Destructor Documentation
11.60.3 Member Function Documentation
11.61 CMeaAudioFunctionNet Class Reference
11.61.1 Constructor & Destructor Documentation
11.61.2 Member Function Documentation
11.62 CMeaCleanDeviceNet Class Reference
11.62.1 Detailed Description
11.62.2 Constructor & Destructor Documentation
11.62.3 Member Function Documentation
11.63 CMeaCoatDeviceNet Class Reference
11.63.1 Detailed Description
11.63.2 Constructor & Destructor Documentation
11.63.3 Member Function Documentation
11.64 CMeaDeviceNet Class Reference
11.64.1 Detailed Description
11.64.2 Constructor & Destructor Documentation

11.64.3 Member Function Documentation
11.64.4 Property Documentation
11.65 CMeaDigitalDataFunctionNet Class Reference
11.65.1 Constructor & Destructor Documentation
11.65.2 Member Function Documentation
11.66 CMeaFeedbackFunctionNet Class Reference
11.66.1 Constructor & Destructor Documentation
11.66.2 Member Function Documentation
11.67 CMealmpedanceDeviceNet Class Reference
11.67.1 Constructor & Destructor Documentation
11.67.2 Member Function Documentation
11.68 CMeasureTableDeviceNet Class Reference
11.68.1 Detailed Description
11.68.2 Constructor & Destructor Documentation
11.68.3 Property Documentation
11.69 CMeaSwitchDeviceNet Class Reference
11.69.1 Detailed Description
11.69.2 Constructor & Destructor Documentation
11.69.3 Member Function Documentation
11.70 CMeaUSBDeviceNet Class Reference
11.70.1 Detailed Description
11.70.2 Constructor & Destructor Documentation
11.71 CMeFunctionNet Class Reference
11.71.1 Detailed Description
11.71.2 Constructor & Destructor Documentation
11.71.3 Member Function Documentation
11.72 CMultiBatteryChargerDeviceNet Class Reference
11.72.1 Detailed Description
11.72.2 Constructor & Destructor Documentation
11.72.3 Member Function Documentation
11.73 CMultiwellCallbackFunctionNet Class Reference
11.73.1 Detailed Description
11.73.2 Constructor & Destructor Documentation
11.73.3 Member Function Documentation
11.73.4 Event Documentation
11.74 CMultiwellDeviceNet Class Reference
11.74.1 Detailed Description
11.74.2 Constructor & Destructor Documentation
11.74.3 Member Function Documentation
11.75 CMultiwellOptoStimFunctionNet Class Reference
11.75.1 Detailed Description
11.75.2 Constructor & Destructor Documentation

11.75.3 Member Function Documentation	15
11.76 CNF_GenDeviceNet Class Reference	19
11.76.1 Constructor & Destructor Documentation	19
11.76.2 Member Function Documentation	19
11.77 COctoPotDeviceNet Class Reference	20
11.77.1 Constructor & Destructor Documentation	20
11.77.2 Member Function Documentation	21
11.78 COkuvisionStimulatorDeviceNet Class Reference	24
11.78.1 Constructor & Destructor Documentation	25
11.78.2 Member Function Documentation	25
11.79 CPatchServerDeviceNet Class Reference	27
11.79.1 Detailed Description	28
11.79.2 Constructor & Destructor Documentation	28
11.79.3 Property Documentation	28
11.80 CPathIdentDeviceNet Class Reference	28
11.80.1 Constructor & Destructor Documentation	29
11.80.2 Member Function Documentation	29
11.81 CPedoterDeviceNet Class Reference	29
11.81.1 Detailed Description	30
11.81.2 Constructor & Destructor Documentation	30
11.81.3 Member Function Documentation	30
11.82 CPeristalticPumpDeviceNet Class Reference	31
11.82.1 Detailed Description	31
11.82.2 Constructor & Destructor Documentation	31
11.82.3 Property Documentation	32
11.83 CPgaDeviceNet Class Reference	32
11.83.1 Constructor & Destructor Documentation	32
11.83.2 Member Function Documentation	33
11.84 CPositionIIDeviceNet Class Reference	34
11.84.1 Detailed Description	36
11.84.2 Constructor & Destructor Documentation	36
11.84.3 Member Function Documentation	36
11.84.4 Property Documentation	43
11.85 CPositionImpDeviceNet Class Reference	43
11.85.1 Detailed Description	1 3
11.85.2 Constructor & Destructor Documentation	14
11.85.3 Member Function Documentation	14
11.86 CPPCDeviceNet Class Reference	16
11.86.1 Constructor & Destructor Documentation	16
11.86.2 Property Documentation	17
11.87 CPPCFunctionNet Class Reference	17
11.87.1 Datailed Description	12

11.87.2 Constructor & Destructor Documentation	8
11.87.3 Member Function Documentation	١9
11.88 CPPS_DeviceNet Class Reference	6
11.88.1 Constructor & Destructor Documentation	6
11.88.2 Property Documentation	6
11.89 CPPS_FunctionNet Class Reference	57
11.89.1 Constructor & Destructor Documentation	57
11.89.2 Member Function Documentation	58
11.90 CPPSDeviceNet Class Reference	31
11.90.1 Detailed Description	31
11.90.2 Constructor & Destructor Documentation	32
11.91 CProgramPressureCurveNet Class Reference	32
11.91.1 Detailed Description	32
11.91.2 Constructor & Destructor Documentation	32
11.91.3 Member Function Documentation	3
11.92 CPulseGeneratorFunctionNet Class Reference	3
11.92.1 Detailed Description	64
11.92.2 Constructor & Destructor Documentation	34
11.92.3 Member Function Documentation	35
11.93 CRadioControledDevicesNet Class Reference	6
11.93.1 Constructor & Destructor Documentation	57
11.93.2 Member Function Documentation	57
11.94 CCMOSMeaDeviceNet::CRegionOfInterestRect Class Reference	8
11.94.1 Constructor & Destructor Documentation	8
11.94.2 Member Function Documentation	8
11.94.3 Member Data Documentation	
11.95 CRetinaLedDeviceNet Class Reference	9
11.95.1 Constructor & Destructor Documentation	'0
11.95.2 Member Function Documentation	'0
11.96 CRFFunctionNet Class Reference	
11.96.1 Detailed Description	'2
11.96.2 Constructor & Destructor Documentation	'2
11.96.3 Member Function Documentation	'2
11.97 CRobo_FYIProgram_FunctionNet Class Reference	'6
11.97.1 Constructor & Destructor Documentation	'6
11.97.2 Member Function Documentation	'6
11.98 CRobo_FYITemp_FunctionNet Class Reference	7
11.98.1 Constructor & Destructor Documentation	'8
11.98.2 Member Function Documentation	
11.99 CRoboDacqNet Class Reference	'9
11.99.1 Constructor & Destructor Documentation	32
11.99.2 Member Function Documentation	32

11.100 CRoboDeviceNet Class Reference
11.100.1 Detailed Description
11.100.2 Constructor & Destructor Documentation
11.100.3 Member Function Documentation
11.100.4 Member Data Documentation
11.100.5 Property Documentation
11.100.6 Event Documentation
11.101 CRoboFluidDeviceNet Class Reference
11.101.1 Constructor & Destructor Documentation
11.101.2 Member Function Documentation
11.101.3 Member Data Documentation
11.101.4 Property Documentation
11.102 CRobolnjectDeviceNet Class Reference
11.102.1 Detailed Description
11.102.2 Constructor & Destructor Documentation
11.103 CRoboocyte2DeviceNet Class Reference
11.103.1 Detailed Description
11.103.2 Constructor & Destructor Documentation
11.103.3 Member Function Documentation
11.104 CRoboStatorDeviceNet Class Reference
11.104.1 Constructor & Destructor Documentation
11.104.2 Member Function Documentation
11.104.3 Property Documentation
11.105 CSafeISDeviceNet Class Reference
11.105.1 Detailed Description
11.105.2 Constructor & Destructor Documentation
11.105.3 Member Function Documentation
11.105.4 Property Documentation
11.106 CSCUDacqGroupChannelSelectionNet Class Reference
11.106.1 Constructor & Destructor Documentation
11.107 CSCUFunctionNet Class Reference
11.107.1 Detailed Description
11.107.2 Constructor & Destructor Documentation
11.107.3 Member Function Documentation
11.107.4 Event Documentation
11.108 CSerialPortNet Class Reference
11.108.1 Constructor & Destructor Documentation
11.108.2 Member Function Documentation
11.109 CStg200xBasicNet Class Reference
11.109.1 Detailed Description
11.109.2 Constructor & Destructor Documentation
11.109.3 Member Function Documentation

11.110 CStg200xDownloadBasicNet Class Reference
11.110.1 Detailed Description
11.110.2 Member Function Documentation
11.110.3 Property Documentation
11.111 CStg200xDownloadNet Class Reference
11.111.1 Detailed Description
11.111.2 Constructor & Destructor Documentation
11.111.3 Member Function Documentation
11.111.4 Event Documentation
11.112 CStimulusFunctionNet Class Reference
11.112.1 Constructor & Destructor Documentation
11.112.2 Member Function Documentation
11.112.3 Event Documentation
11.113 CSw2to64DeviceNet Class Reference
11.113.1 Detailed Description
11.113.2 Constructor & Destructor Documentation
11.113.3 Member Function Documentation
11.114 CTcxDeviceNet Class Reference
11.114.1 Detailed Description
11.114.2 Constructor & Destructor Documentation
11.114.3 Member Function Documentation
11.115 CTEERFunctionNet Class Reference
11.115.1 Detailed Description
11.115.2 Constructor & Destructor Documentation
11.115.3 Member Function Documentation
11.116 CTEERMachineDeviceNet Class Reference
11.116.1 Constructor & Destructor Documentation
11.116.2 Property Documentation
11.117 CUsbDeviceConfigurationFunctionNet Class Reference
11.117.1 Detailed Description
11.117.2 Constructor & Destructor Documentation
11.117.3 Member Function Documentation
11.118 CUsbExceptionNet Class Reference
11.118.1 Detailed Description
11.118.2 Constructor & Destructor Documentation
11.118.3 Property Documentation
11.119 CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet Class Reference
11.119.1 Constructor & Destructor Documentation
11.119.2 Member Data Documentation
11.120 CW2100_FunctionNet Class Reference
11.120.1 Constructor & Destructor Documentation
11.120.2 Member Function Documentation

11.120.3 Property Documentation	34
11.121 CW2100_StimulatorFunctionNet Class Reference	34
11.121.1 Constructor & Destructor Documentation	36
11.121.2 Member Function Documentation	36
11.121.3 Member Data Documentation	41
11.121.4 Event Documentation	42
11.122 CW2100DacqGroupChannelSelectionNet Class Reference	42
11.122.1 Constructor & Destructor Documentation	42
11.123 CWarnerUssingDeviceNet Class Reference	43
11.123.1 Detailed Description	43
11.123.2 Constructor & Destructor Documentation	43
11.123.3 Property Documentation	44
11.124 CWarnerUssingFunctionNet Class Reference	44
11.124.1 Detailed Description	46
11.124.2 Constructor & Destructor Documentation	46
11.124.3 Member Function Documentation	46
11.125 CWarnerValveControllerDeviceNet Class Reference	60
11.125.1 Detailed Description	64
11.125.2 Constructor & Destructor Documentation	64
11.125.3 Member Function Documentation	65
11.125.4 Event Documentation	82
11.126 CWarnerValveControllerDeviceTesterFunctionNet Class Reference	85
11.126.1 Detailed Description	86
11.126.2 Constructor & Destructor Documentation	86
11.126.3 Member Function Documentation	86
11.127 CWClassicFunctionNet Class Reference	88
11.127.1 Constructor & Destructor Documentation	89
11.127.2 Member Function Documentation	89
11.128 CWirelessBaseFunctionNet Class Reference	93
11.128.1 Constructor & Destructor Documentation	
	93
11.128.2 Member Function Documentation	93 93
11.129 DeviceIdNet Struct Reference	93
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55	93 93
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55	93 93 94 94
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 55	93 93 94 94
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 55 11.129.3 Member Function Documentation 55 11.129.4 Member Data Documentation 55	93 93 94 94
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 55 11.129.3 Member Function Documentation 55 11.129.4 Member Data Documentation 55 11.130 DigitalSource 55 11.130 DigitalSource 55	93 93 94 94 95
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 55 11.129.3 Member Function Documentation 55 11.129.4 Member Data Documentation 55 11.130 DigitalSource digitalsourceenum > Class Template Reference 55 11.130.1 Constructor & Destructor Documentation 55	93 93 94 94 95 95
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 56 11.129.3 Member Function Documentation 56 11.129.4 Member Data Documentation 56 11.130 DigitalSource digitalsourceenum Class Template Reference 56 11.130.1 Constructor & Destructor Documentation 56 11.130.2 Member Function Documentation 56	93 93 94 94 95 95
11.129 DeviceIdNet Struct Reference 55 11.129.1 Detailed Description 55 11.129.2 Constructor & Destructor Documentation 56 11.129.3 Member Function Documentation 56 11.129.4 Member Data Documentation 57 11.130 DigitalSource digitalsourceenum > Class Template Reference 56 11.130.1 Constructor & Destructor Documentation 57 11.130.2 Member Function Documentation 57 11.130.3 Property Documentation 57	93 93 94 94 95 95 96

11.131.2 Member Function Documentation	597
11.131.3 Property Documentation	597
11.132 DriverVersionNet Class Reference	598
11.132.1 Detailed Description	599
11.132.2 Constructor & Destructor Documentation	599
11.132.3 Member Function Documentation	599
11.133 FirmwareDestinationNames Class Reference	603
11.133.1 Member Data Documentation	603
11.134 HeadStageIDType Class Reference	606
11.134.1 Member Enumeration Documentation	607
11.134.2 Constructor & Destructor Documentation	607
11.134.3 Member Function Documentation	607
11.134.4 Property Documentation	608
11.135 HeadstageIDTypeObject Class Reference	609
11.135.1 Constructor & Destructor Documentation	609
11.135.2 Member Function Documentation	610
11.135.3 Member Data Documentation	610
11.135.4 Property Documentation	610
11.136 HeadStageIDTypeState Class Reference	610
11.136.1 Property Documentation	611
11.137 mkfilterNet Class Reference	611
11.137.1 Member Function Documentation	612
11.138 CRoboDeviceNet::RoboMainLowLevelCommands Class Reference	614
11.138.1 Member Function Documentation	615
11.139 CRoboStatorDeviceNet::RoboMainStatorLowLevelCommands Class Reference	621
11.139.1 Member Function Documentation	621
11.140 CFilterCoefficientsNet::s_FilterAttributesNet Struct Reference	621
11.140.1 Constructor & Destructor Documentation	622
11.140.2 Member Function Documentation	622
11.140.3 Member Data Documentation	622
11.141 CMeaAudioFunctionNet::s_setaudionet Struct Reference	623
11.141.1 Member Data Documentation	623
11.142 CStimulusFunctionNet::SidebandData Class Reference	623
11.142.1 Constructor & Destructor Documentation	623
11.142.2 Property Documentation	624
11.143 StgStatusNet Class Reference	624
11.143.1 Member Function Documentation	624
11.143.2 Member Data Documentation	625
11.144 CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData Class Reference	625
11.144.1 Constructor & Destructor Documentation	625
11.144.2 Property Documentation	626
11 1/5 uchSatunPackat + Class Reference	626

Index	628
11.146.1 Member Data Documentation	
11.146 W2100_StimulusParametersNet Struct Refer	nce
11.145.1 Member Data Documentation	

1 McsUsbNet.dll for MCS USB devices

1.1 Introduction

This DLL provides the .NET interface to MCS devices

The most important options are accessing our stimulator and data acquisition devices:

- STG200x & STG400x STimulus Generator
- Data ACQuisition (DACQ) Devices

See here for a list of our other devices: Device Classes.

And here for a list of function classes addressing groups of features that might be shared between different devices: Function Classes.

1.2 System requirements

The DLL can be used with any .NET compatible language.

The DLL needs the .NET Framework 4.7.2.

It requires the Microsoft Visual C++ Redistributable for Visual Studio 2019 to be installed.

It also requires the **USB driver** to be installed.

The simplest way to achieve this is to install the latest **Multi Channel Experimenter** setup (will install 64bit redistributable).

All examples assume that the Mcs.Usb namespace is loaded:

```
using namespace Mcs.Usb;
```

Include the file McsUsbNet.dll into the references of your project.

1.3 Connecting to an MCS device

A connection to a DAQ device is established by Mcs.Usb.CMcsUsbNet.Connect. When this function is called without argument, the first DAQ device found on the USB bus is used:

```
CMcsUsbNet device = new CMcsUsbNet();
device.Connect();
```

When more than one DAQ device of the specific type is connected, you can use the Mcs.Usb.CMcsUsbListNet class to get a list of available devices:

```
CMcsUsbListNet usblist = new CMcsUsbListNet(DeviceEnumNet.MCS_DEVICE_USB);
var entry = usblist.GetUsbListEntry((uint)0);
CMcsUsbNet device = new CMcsUsbNet();
device.Connect(entry);
```

After you are finished with the device, you can disconnect the device object from the device by:

```
device.Disonnect();
```

2 Device Classes 2

2 Device Classes

- For FluidControl device see MCS FluidControl
- For SW2TO64 device see MCS-USB-Sw2to64
- For TCx device see Mcs.Usb.CTcxDeviceNet

2.1 The MCS FluidControl Device

2.1.1 Introduction

The FluidControl Device can control up to 24 valves. The nominal voltage is 24V.

8 TTL level digital output ports are available and 8 TTL inputs can be read in.

The device has 8 ADC inputs with a rage from 0V to 3.3V.

2.1.2 Access to the FluidControl device

For connecting to a FluidControl device see Connecting to an MCS device.*

```
CFluidControlDevice* m_dacq;
m_fluidcontrol = new CFluidControlDevice;
status = m_fluidcontrol->Connect();
```

The valves are controlled with the CFluidControlDevice::SetValve call. The argument given is a bit pattern of all valves which should be open.

The digital outputs can be controlled with the CFluidControlDevice::SetDigout call. Again, a bit pattern of all digital output pins which should be set to a logic high level is given as an argument.

The current state of the valves and the digital outputs can be read back with the CFluidControlDevice::GetValve and CFluidControlDevice::GetDigout

The command to read an ADC-Channel is CFluidControlDevice::GetAdc. Here the channelnummer which should be read in is given as an argument and the return value is the current Adc level.

The state of the digital inputs is read with the CFluidControlDevice::GetDigin call. Here the return value is the bit pattern of the digital inputs.

The connection to the device is closed with the CFluidControlDevice::Disconnect call.

2.2 MCS-USB-Sw2to64 device

The class Mcs.Usb.CSw2to64DeviceNet controls the setting of the switches in the MCS-USB-Sw2to64 device.

First construct an object of the class:

```
CSw2to64DeviceNet device = new CSw2to64DeviceNet();
```

For connecting to an MCS-USB-Sw2to64 device see Connecting to an MCS device.

To get the number of channels the device handles:

```
int number = device.GetNumber()
```

Set all channel switches at once:

```
byte 2 = 1;
byte[] pattern = new byte[number];
for(int i = 0; i < number; i++)
{
   pattern[i] = z; // pattern you want to switch this channel to
}
device.SetChannels(pattern);</pre>
```

Get all channel switches at once:

```
byte[] pattern = device.GetChannels();
```

Set one channel switch:

```
ushort index = 10;
byte pattern = 1;
device.SetChannel(index, pattern)
```

Get one channel switch:

```
ushort index = 10;
byte pattern = device.GetChannel(index);
```

3 Function Classes

- Mcs.Usb.CCMOSMea_FunctionNet
- Mcs.Usb.CDacCalibrationFunctionNet
- Mcs.Usb.CDigOutStimulatorFunctionNet
- Mcs.Usb.CIntanMea_FunctionNet
- Mcs.Usb.CInterfaceboardFunctionNet
- Mcs.Usb.CMcsBus_MotorControlNet
- Mcs.Usb.CMcsBus_VoltageModeNet
- Mcs.Usb.CMcsBus_AxisParametersNet
- Mcs.Usb.CMcsBus_SensorNet
- Mcs.Usb.CMcsBus_TempSensorNet
- Mcs.Usb.CMcsBus_ExtensionNet
- Mcs.Usb.CMcsBus_FYIExtensionNet
- Mcs.Usb.CMcsUsbDeviceStatePushFunctionNet
- Mcs.Usb.CMEA2100x256FunctionNet
- Mcs.Usb.CMeaAudioFunctionNet
- Mcs.Usb.CMeaDigitalDataFunctionNet

- Mcs.Usb.CMeaFeedbackFunctionNet
- Mcs.Usb.CMeFunctionNet
- · Mcs.Usb.CMultiwellCallbackFunctionNet
- Mcs.Usb.CMultiwellOptoStimFunctionNet
- Mcs.Usb.CPPCFunctionNet
- Mcs.Usb.CPPS_FunctionNet
- · Mcs.Usb.CPPS FunctionNet
- Mcs.Usb.CPulseGeneratorFunctionNet
- · Mcs.Usb.CRFFunctionNet
- Mcs.Usb.CRobo_FYITemp_FunctionNet
- Mcs.Usb.CRobo_FYIProgram_FunctionNet
- Mcs.Usb.CRobo_FYITemp_FunctionNet
- Mcs.Usb.CRobo_FYIProgram_FunctionNet
- Mcs.Usb.CSCUFunctionNet
- Mcs.Usb.CStimulusFunctionNet
- Mcs.Usb.CTEERFunctionNet
- Mcs.Usb.CW2100_FunctionNet
- Mcs.Usb.CW2100_StimulatorFunctionNet
- Mcs.Usb.CWarnerUssingFunctionNet
- Mcs.Usb.CWarnerValveControllerDeviceTesterFunctionNet
- Mcs.Usb.CWClassicFunctionNet
- Mcs.Usb.CWirelessBaseFunctionNet

4 Data ACQuisition (DACQ) Devices

There are different device types of (MEA) data acquisition (DACQ) devices. All of them are supported by this class.

This library does **not** support the writing of the MCD (MC_Rack), MSRD (Multi Channel Experimenter) or HDF5 file format!

The class Mcs.Usb.CMeaDeviceNet is the base class for DACQ devices.

The base class Mcs.Usb.CMeaDeviceNet constructs actually the underlying classes for USB-MEA devices (Mcs.Usb.CMeaUSBDeviceNet).

```
.
CMeaDeviceNet device = new CMeaDeviceNet (McsBusTypeEnumNet.MCS_USB_BUS, OnChannelData, OnError);
```

For connecting to a DACQ device see Connecting to an MCS device.

Get the number of available analog hardware channels and set the number of channels to the maximum.

```
int hwchannels;
device.HWInfo().GetNumberOfHWADCChannels(out hwchannels);
device.SetNumberOfChannels(hwchannels);
int samplingrate = 1000;
device.SetSamplerate(samplingrate, 1, 0);
device.EnableDigitalIn(true, 0);
```

5 The MCS Robo Device 5

Get the layout to know how the data look like that you receive

```
int ana, digi, che, tim, block;
device.GetChannelLayout(out ana, out digi, out che, out tim, out block);
```

For the Mcs.Usb.OnChannelData callback function you have to provide a definition of the channels you want to receive.

```
bool[] selChannels = new bool[block];
for (int i = 0; i < block; i++)
{
    selChannels[i] = true; // With true channel i is selected
    // selChannels[i] = false; // With false the channel i is deselected
}
channelblocksize = samplingrate / 10;
// queue size and threshold should be selected carefully
device.SetSelectedChannels(selChannels, 10 * channelblocksize, channelblocksize);</pre>
```

The Mcs.Usb.OnChannelData callback function gets a callback for each channelblock that is defined. In this example a callback for each channel.

```
void OnChannelData(CMcsUsbDacqNet d, int cbHandle, int numSamples)
{
   int size_ret;
   ushort[] channeldata = device.ChannelBlock_ReadFramesUI16(CbHandle, numSamples, out size_ret);
}
void OnError(String msg, int info)
{
   MessageBox.Show("Mea Device Error: " + msg);
}
```

see MEA Recording in the Examples directory.

5 The MCS Robo Device

5.1 Introduction

Up to now two MCS devices exist that base on the Robo platform.

- The MCS Roboinject device is controlled by the Mcs.Usb.CRobolnjectDeviceNet class.
- The MCS Roboocyte2 device is controlled by the Mcs.Usb.CRoboocyte2DeviceNet class.

Both classes are derived from Mcs.Usb.CRoboDeviceNet

6 STG200x & STG400x STimulus Generator

6.1 Introduction

The STG200x & STG400x Series Stimulus Generators have two distinct modes of operation, the Download mode and the Streaming mode.

6.2 Download mode 6

6.2 Download mode

The Download mode is the "classic" mode of operation, as used by the MC Stimulus software. In this mode, one or multiple waveforms are defined in PC memory and downloaded to the STG. The waveforms are stored in STG device onboard memory and can be sent to the analog and sync outputs once or multiple times. The STG can operate independently from the PC (without computer connection) after the download. Output is triggered either by the front panel start/stop button, the digital trigger inputs or under software control.

In the Download mode, there are up to eight independent triggers available (depending on the device). The user can assign each of the analog outputs and sync (digital) outputs to any of the triggers.

The analog output waveform is stored sample by sample in the STG memory. To reduce memory usage, this data can be compressed: whenever a given output value is to be held for more than one sample period, it has only to be given once. The user can define the number of sample periods for that a pattern should remain active. Compression is done for each channel independently of the others, thus the algorithm to compress the data is very easy to implement.

A new feature of the Download mode is the segmentation of the STG memory. The onboard memory can be devided into up to 100 segments. Each segment can hold its own waveform pattern. Under software control, the user can switch between the defined segments within milliseconds. Another option is to use the four trigger inputs to select between four predefined segments. This option is accessible from the MC_Stimulus Software as the "Multi-File mode", and can start each of up to four defined waveforms within microseconds. This feature allows a predefinied flexible response (feedback) to recorded data.

Mcs.Usb.CStg200xDownloadNet is the class for using the STG in download mode.

6.2.1 Memory Layout and Trigger Setup

The class to be used for the Download mode is Mcs.Usb.CStg200xDownloadNet, which is derived from Mcs.Usb.CStg200xBasicNet. You can add a poll handler delegate (Mcs.Usb.OnStg200xPollStatus) to the constructor Mcs.Usb.CStg200xDownloadNet.

For connecting to an STG see Connecting to an MCS device.

To use the Download mode, the memory layout of the STG200x can be set up, if the default is not sufficient. The total amount of memory available in the STG is obtained by the Mcs.Usb.CStg200xDownloadNet.GetTotalMemory call. With Mcs.Usb.CStg200xDownloadNet.SendSegmentDefine the segment sizes are assigned.

```
uint32_t memory = device.GetTotalMemory(); // obtain total memory available
uint[] segmentmemory = new uint[2]; // each segments has half of total memory
segmentmemory[0] = memory / 2;
segmentmemory[1] = memory / 2;
device.SendSegmentDefine(segmentmemory);// setup the STG
```

Next, for each segment, one has to assign the amount of memory to be used for each channel and sync output. This is done by Mcs.Usb.CStg200xDownloadBasicNet.SetCapacity. Its arguments contain a list of memory sizes, with one entry per channel and one entry per sync output. Again, the total memory assigned to the channels and sync outputs must not exceed the memory assigned to the segment.

6.2 Download mode 7

}

Before the STG can start, the trigger has to be configured. This is done by the Mcs.Usb.CStg200xDownloadNet.SetupTrigger call. Its arguments are a list of channelmaps, syncoutmaps and repeats, one for each of the four available triggers. channelmap is a bitmap, each bit representing one of the available channels. To assign channel 1 and syncout 1 to trigger 1 and channel 3 to trigger 2 use:

```
uint32_t TriggerInputs = device.GetNumberOfTriggerInputs();
uint[] channelmap = new uint[TriggerInputs];
uint[] syncoutmap = new uint[TriggerInputs];
uint[] repeat = new uint[TriggerInputs];
for (int i = 0; i < TriggerInputs; i++)
{
    channelmap[i] = 0;
    syncoutmap[i] = 0;
    repeat[i] = 0;
}
// Trigger 0
    channelmap[0] = 1; // Channel 1
    syncoutmap[0] = 1; // Syncout 1
    repeat[0] = 0; // forever
// Trigger 1
    channelmap[1] = 4; // Channel 3
    device.SetupTrigger(channelmap, syncoutmap, repeat);</pre>
```

For the STG400x series you have to set the output mode of the channels. Mcs.Usb.CStg200xDownloadNet.SetVoltageMode interprets the values as voltages. Mcs.Usb.CStg200xDownloadNet.SetCurrentMode as currents.

```
// Only meaningfull for STG400x
device.SetVoltageMode();
```

For each segment, data can be sent to each of the defined channels and sync outputs using the Mcs.Usb.CStg200xDownloadNet.SendChannelData and Mcs.Usb.CStg200xDownloadNet.SendSyncData calls. channeldata and syncdata are a list of analog and digital samples as a list of two byte values (unsigned short). Multiple calls to Mcs.Usb.CStg200xDownloadNet.SendSyncData to the same channel append data to that channel.

If the Multi-File mode of the STG is enabled using the Mcs.Usb.CStg200xDownloadNet.EnableMultiFileMode call, the four trigger inputs are used to switch between four segments. A hardware trigger signal (TTL) on trigger input 1 selects the first segment and starts all pulses in this segment. Thus with the Multi-File mode, one can predefine four stimulus patterns and switch between them without a connection to the PC.

The STG200x series has an analog resolution of 13 bits, thus the analog data contains the information in bits 0 to 12 of each sample. Bits 13 to 15 have to be 0.

```
int DACResolution = device GetDACResolution();
// Data for Channel 0
    device.ClearChannelData(0);
    double factor = 0.1;
    const int 1 = 1000;
   ushort[] pData = new ushort[1];
        Uint64_t[] tData = new Uint64_t[1];
        for (int i = 0; i < 1; i++)
             // calculate Sin-Wave
            double sin = factor * (Math.Pow(2, DACResolution - 1) - 1.0) *
                Math.Sin(2.0 * (double)i * Math.PI / (double)1);
             // calculate sign
            pData[i] = sin >= 0 ? (ushort)sin : (ushort)((int)Math.Abs(sin) +
                 (int)Math.Pow(2, DACResolution - 1));
            tData[i] = (Uint64_t)20; // duration in \mus
        device. SendChannelData(0, pData, tData);
// Data for Channel 3
    device.ClearChannelData(2);
    double factor = 0.1;
    const int 1 = 700;
    // without compression
    ushort[] pData = new ushort[1];
    Uint64_t[] tData = new Uint64_t[1];
    for (int i = 0; i < 1; i++)
        // calculate Sin-Wave
        double sin = factor * (Math.Pow(2, DACResolution - 1) - 1.0) *
           Math.Sin(2.0 * (double)i * Math.PI / (double)1);
```

6.3 Streaming mode 8

Start the trigger by pushing the front button or by software

```
// Start Trigger 1 and 2
device.SendStart(1 + 2); // Trigger 1 und 2
```

see the StgDownloadExampleNet in the example directory.

6.3 Streaming mode

The other mode of operation is the Streaming mode. Here the analog output is sent to the STG device in "real time". The PC has to be connected to the STG all the time. The data that is sent to the analog output is downloaded from the PC to the STG on the fly.

The Streaming mode is useful for applications where flexible feedback is needed as well for applications where very long waveforms which are not repeated (such as white noise) are used.

The Streaming mode works by use of two ring buffers which hold data. One is in PC memory and managed by the DLL, and one is in on-board STG memory. Data is transferred from PC memory to the STG via the USB bus in time slices of one millisecond.

The user can define both the size of the ring buffer in DLL memory and in the STG memory. Once the Streaming mode is started, the STG request data from the PC. The data rate from PC to STG is variable and controlled by the STG. The STG request data from the PC at a rate to keep its internal ringbuffer at about half full.

It is the responsibility of the user to keep the ring buffer in the memory of the PC filled, so the DLL can supply sufficient data to the STG. To do so, the Windows DLL allows to define a "callback" function which is called whenever new data is needed, or more precise, as soon as the ring buffer in the memory of the PC falls below the user defined threshold.

Small buffers have the advantage of a low latency between data generation in the callback funtion and its output as a analog signal from the STG. However for low latency to work, the user-written callback function has to be fast and to produce a steady flow of data.

In the Streaming mode, all triggers are available as well. Each of the eight analog and sync outputs can be assigned to one of the triggers.

The output rate is user defined with a maximum of 50 kHz

Mcs.Usb.CStg200xStreamingNet is the class for using the STG in streaming mode.

6.3 Streaming mode 9

6.3.1 Memory Layout and Trigger Setup

With the constructor for Mcs.Usb.CStg200xStreamingNet.CStg200xStreamingNet, the name of the callback function for the data handler is provided. The data handler function is called automatically, whenever the STG needs new data. This data is first written to a ring buffer in the memory of the PC. The size for this ring buffer is defined as first argument in the constructor. The user provided delegate gets the trigger number which needs new data as argument

```
CStg200xStreamingNet device = new CStg200xStreamingNet(10000, dataHandler, errorHandler);
```

The callback funtion, which is defined in the constructor, is called whenever the STG needs new data for a trigger, or more precise, whenever the ring buffer in PC memory falls below the defined threshold.

The user can query the amount of space available for queuing by use of the Mcs.Usb.CStg200xStreamingNet. ← GetDataQueueSpace call. Its return value is the number of samples that can be send to the STG.

User code is required to fill an array analog and sync out data, sample by sample for up to the maximum number of samples as obtained by Mcs.Usb.CStg200xStreamingNet.GetDataQueueSpace or Mcs.Usb.CStg200xStreaming

Net.GetSyncoutQueueSpace.

The values for the analog outputs are 16 bits signed integers. The lower bits are trunctated according to the resolution of the STG. This behaviour is different to the behaviour in download mode.

Note: Compression as described in the download mode can NOT be used for the streaming mode.

The new data is sent to the STG by using the Mcs.Usb.CStg200xStreamingNet.EnqueueData call.

```
void dataHandler(uint32_t trigger)
    double factor = 1;
if (trigger == 0) // Callback for Trigger 1
          {// Handle Channel 1
               uint32_t channel = 0;
                    uint32 t space = device.GetDataOueueSpace(channel);
                    if (space < 1000)
                    short[] data = new short[1000];
                    for (int i = 0; i < 1000; i++)
                         // Calc Sin-Wave (16 bits) lower bits will be removed according resolution double sin = factor \star (Math.Pow(2, 16 - 1) - 1.0) \star Math.Sin(2.0 \star (double)i \star Math.PI / (double)1000);
                         data[i] = (short)sin;
                    uint32_t engueued = device.EngueueData(channel, data);
          {// Handle Channel 3
               uint32_t channel = 2;
               for (; ; )
                    uint32_t space = device.GetDataQueueSpace(channel);
                    if (space < 700)</pre>
                         break:
                    short[] data = new short[700];
                    for (int i = 0; i < 700; i++)
                         // Calc Sin-Wave (16 bits) lower bits will be removed according resolution double sin = factor * (Math.Pow(2, 16 - 1) - 1.0) * Math.Sin(2.0 * (double)i * Math.PI / (double)700);
                         data[i] = (short)sin;
                    uint32_t enqueued = device.EnqueueData(channel, data);
          {// Handle Syncout 1
               uint32_t channel = 0;
               for (; ; )
                    uint32_t space = device.GetSyncoutQueueSpace(channel);
                    if (space < 1000)</pre>
                         break:
                    ushort[] data = new ushort[1000];
                    for (int i = 0; i < 1000; i++)</pre>
```

6.3 Streaming mode 10

For connecting to an STG device see Connecting to an MCS device.

With enabling or disabling the continuous mode it can be selected how the STG handles an "out of data" situation.

When Mcs.Usb.CStg200xStreamingNet.EnableContinousMode is used, the STG does not stop when it runs out of data, but it keeps running and sends a zero voltage to its outputs.

When Mcs.Usb.CStg200xStreamingNet.DisableContinousMode is used, the STG stops when it runs out of data. It has to be retriggered to resume the output.

```
device.EnableContinousMode();
```

Mcs.Usb.CStg200xStreamingNet.SetOutputRate is used to set the sampling rate. device.SetOutputRate(50000);

To use the Streaming mode, the memory layout of the STG has to be set up. To total amount of memory available in the STG is obtained by the Mcs.Usb.CStg200xStreamingNet.GetTotalMemory call.

This memory can be assigned to four ring buffers (one per trigger) which buffer the data received from the PC via USB cable. This is done with the CStg200xStreaming::SetCapacity call. The total amount of memory must not exceed the total memory size as obtained by Mcs.Usb.CStg200xStreamingNet.GetTotalMemory.

This internal ring buffer is crucial for proper operation of the Streaming mode. The size of the ring buffer determines the latency of the Streaming mode. The firmware of the STG requests data from the PC in order to keep the ring buffer about half full. Thus the average latency is:

```
latency = (ringbuffersize in bytes/4) / output rate
```

If the ring buffer size is too big, the latency of the STG might be too long. If the ring buffer size is too low, an overflow or underflow of data in the STG ringbuffer might occur, resulting in data jumps of the output signals or the "out of data" situation described erlier.

The following example divides the total memory equally amoung the four triggers:

Before the STG can start, the trigger has to be configured. This is done by the Mcs.Usb.CStg200xStreaming
Net.SetupTrigger call. Its arguments are a list of channelmaps, syncoutmaps, digoutmap, autostart and callback
_threshold, with one entry for each of the available triggers. channelmap is a bitmap, each bit representing one of the available channels. To assign channel 1 and 3 and syncout 1 to trigger 1 use:

```
uint32_t ntrigger = device.GetNumberOfTriggerInputs();  // obtain number of triggers in this STG
uint[] channelmap = new uint[ntrigger];
uint[] syncoutmap = new uint[ntrigger];
uint[] digoutmap = new uint[ntrigger];
uint[] autostart = new uint[ntrigger];
uint[] callback_threshold = new uint[ntrigger];
for (int i = 0; i < ntrigger; i++)
{
    channelmap[i] = 0;</pre>
```

7 Namespace Index 11

```
syncoutmap[i] = 0;
digoutmap[i] = 0;
autostart[i] = 0;
callback_threshold[i] = 0;
}
channelmap[0] = 0x1 + 0x4; // Channel 1 und Channel 3 to Trigger 1
syncoutmap[0] = 0x1; // Syncout 1 to Trigger 1
autostart[0] = 1;
callback_threshold[0] = 50; // 50% of buffer size
device.SetupTrigger(channelmap, syncoutmap, digoutmap, autostart, callback_threshold);
device.StartLoop();
System.Threading.Thread.Sleep(1000); // Give StartLoop some time
```

Start Trigger by pushing the front button or by Software

device.SendStart(1);

see the StgStreamingExampleNet in the example directory.

7 Namespace Index

7.1 Namespace List

Here is a list of all namespaces with brief descriptions:

Mcs	22
Mcs::Usb	22

8 Hierarchical Index

8.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

CW2100_FunctionNet::AudioChannelsNet	
BatteryState	29
CCreateFilterNet	48
BesselFilterHighPassNet	30
BesselFilterLowPassNet	31
ButterworthFilterHighPassNet	31
ButterworthFilterLowPassNet	32
${\tt CDeviceGroupChannelInfoTemplateNet} < {\tt DacqGroupChannelEnumTemplateNet} >$	58
${\bf CDeviceGroupChannelInfoTemplateNet} < {\bf DacqGroupChannelEnumNet} >$	58
CDeviceGroupChannelInfoNet	57
${\tt CDeviceGroupChannelInfoTemplateNet} < {\tt int} >$	58
CDeviceGroupChannelInfoGenericNet	56

8.1 Class Hierarchy 12

${\tt CDeviceGroupChannelInfoTemplateNet} < {\tt MEA2100_256DacqGroupChannelEnumNet} >$	58
CDeviceGroupChannelInfoMEA2100_256Net	56
${\tt CDeviceGroupChannelInfoTemplateNet} < {\tt SCUDacqGroupChannelEnumNet} >$	58
CDeviceGroupChannelInfoSCUNet	57
${\tt CDeviceGroupChannelInfoTemplateNet} < {\tt W2100DacqGroupChannelEnumNet} >$	58
CDeviceGroupChannelInfoW2100Net	59
CFilterCoefficientsNet	66
CFilterPropertyNet	72
CMcsUsbDacqNet::CHWInfo	112
CMcsUsbFunctionNet	232
$\label{lem:continuous} \textbf{CDacqGroupChannelEnumNet}, \textbf{DacqGroup} \leftarrow \\ \textbf{ChannelEnum, CDeviceGroupChannelInfoNet} >$	53
CDacqGroupChannelSelectionNet	52
${\tt CDacqGroupChannelSelectionTemplateNet} {< } {\tt int, int, CDeviceGroupChannelInfoGenericNet} {< }$	53
CDacqGroupChannelGenericSelectionNet	52
CDacqGroupChannelSelectionTemplateNet< MEA2100_256DacqGroupChannelEnumNet, M← EA2100_256DacqGroupChannelEnum, CDeviceGroupChannelInfoMEA2100_256Net >	53
CMEA2100_256DacqGroupChannelSelectionNet	263
${\bf CDacqGroupChannelSelectionTemplateNet} < {\bf SCUDacqGroupChannelEnumNet}, {\bf SCUDacq} \leftarrow {\bf GroupChannelEnum, CDeviceGroupChannelInfoSCUNet} >$	53
CSCUDacqGroupChannelSelectionNet	422
$\label{lem:continuous} \textbf{CDacqGroupChannelSelectionTemplateNet} < & \textbf{W2100DacqGroupChannelEnumNet}, & \textbf{W2100} \leftarrow \textbf{DacqGroupChannelEnum}, & \textbf{CDeviceGroupChannelInfoW2100Net} > \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ $	53
CW2100DacqGroupChannelSelectionNet	542
CCMOSMea_FunctionNet	34
CDacCalibrationFunctionNet	50
$\textbf{CDacqGroupChannelSelectionTemplateNet} < \textbf{DacqGroupChannelEnumTemplateNet}, \textbf{Dacq} \leftarrow \textbf{GroupChannelEnumTemplate}, \textbf{CDeviceGroupChannelInfoTemplateNet} >$	53
CDigOutStimulatorFunctionNet	59
CFilterConfigurationNet	68
CFilterConfigurationRegisterNet	70
CGrapheneFunctionNet	100
CIntanMea_FunctionNet	115
CInterfaceboardFunctionNet	119

CInterfaceboard2FunctionNet	117
CMcsBus_AxisParametersNet	130
CMcsBus_ExtensionNet	132
CMcsBus_FYIExtensionNet	133
CMcsBus_MotorControlNet	135
CMcsBus_SensorNet	152
CMcsBus_TempSensorNet	162
CMcsBus_VoltageModeNet	164
CMcsBusNet	169
CMcsUsbDeviceStatePushFunctionNet	221
CMultiwellCallbackFunctionNet	306
CSCUFunctionNet	423
CMEA2100x256FunctionNet	264
CMeaAudioFunctionNet	265
CMeaDigitalDataFunctionNet	286
CMeaFeedbackFunctionNet	288
CMeFunctionNet	297
CMultiwellOptoStimFunctionNet	314
CPPCFunctionNet	347
CPPS_FunctionNet	357
CProgramPressureCurveNet	362
CPulseGeneratorFunctionNet	363
CRFFunctionNet	371
CRobo_FYIProgram_FunctionNet	376
CRobo_FYITemp_FunctionNet	377
CStimulusFunctionNet	485
CTEERFunctionNet	512
CUsbDeviceConfigurationFunctionNet	523
CW2100_StimulatorFunctionNet	534
CWarnerUssingFunctionNet	544
CWarnerValveControllerDeviceTesterFunctionNet	585
CWirelessBaseFunctionNet	593

8.1	Class	Hiera	rchv
-----	-------	-------	------

CW2100_FunctionNet	527
CWClassicFunctionNet	588
CMcsUsbFunctionPointerContainer	233
CMcsUsbListEntryNet	233
CMcsUsbListNet	239
CMcsUsbNet	242
CExternDTesterDeviceNet	64
CFluidControlDeviceNet	74
CGenericDevelopDeviceNet	81
CGilsonDeviceNet	98
CMcsUsbDacqNet	173
CMeaDeviceNet	278
CMeaUSBDeviceNet	296
CCMOSMeaDeviceNet	45
CHLADacqNet	110
CLIH3DeviceNet	121
CMultiwellDeviceNet	308
CWarnerUssingDeviceNet	543
COctoPotDeviceNet	320
CRoboDacqNet	379
CMcsUsbDeviceStatePushNet	222
CWarnerValveControllerDeviceNet	560
CMcsUsbFactoryNet	223
CMeaCleanDeviceNet	268
CMeaCoatDeviceNet	272
CMealmpedanceDeviceNet	292
CMeaSwitchDeviceNet	294
CChannelTestDeviceNet	33
CMultiBatteryChargerDeviceNet	299
CNF_GenDeviceNet	319
COkuvisionStimulatorDeviceNet	324
CPathIdentDeviceNet	328

	CPedoterDeviceNet	329
	CPeristalticPumpDeviceNet	331
	CPgaDeviceNet	332
	CPositionIIDeviceNet	334
	CPositionImpDeviceNet	343
	CPPCDeviceNet	346
	CPPS_DeviceNet	356
	CRadioControledDevicesNet	366
	CRetinaLedDeviceNet	369
	CRoboDeviceNet	393
	CEncapsulatorDeviceNet	64
	CFYIDeviceNet	80
	CHiClampDeviceNet	109
	CHLADeviceNet	111
	CMeasureTableDeviceNet	293
	CPatchServerDeviceNet	327
	CPPSDeviceNet	361
	CRoboStatorDeviceNet	413
	CRobolnjectDeviceNet	411
	CRoboocyte2DeviceNet	411
	CTEERMachineDeviceNet	522
	CRoboFluidDeviceNet	407
	CSafeISDeviceNet	419
	CSerialPortNet	438
	CStg200xBasicNet	439
	CStg200xDownloadBasicNet	469
	CStg200xDownloadNet	478
	CSw2to64DeviceNet	496
	CTcxDeviceNet	499
CN	IcsUsbPointerContainer	263
CC	MOSMeaDeviceNet::CRegionOfInterestRect	368
CN	IcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet	526

9 Class Index 16

	DeviceIdNet	593
	DigitalSource< digitalsourceenum >	595
	DigitalSourceGeneral	596
	DriverVersionNet Exception	598
	CUsbExceptionNet	525
	FirmwareDestinationNames	603
	HeadstageIDTypeObject	609
	HeadStageIDTypeState IComparable	610
	HeadStageIDType	606
	mkfilterNet	611
	CRoboDeviceNet::RoboMainLowLevelCommands	614
	CRoboStatorDeviceNet::RoboMainStatorLowLevelCommands	621
	CFilterCoefficientsNet::s_FilterAttributesNet	621
	CMeaAudioFunctionNet::s_setaudionet	623
	CStimulusFunctionNet::SidebandData	623
	StgStatusNet stgstreaming	624
	CStg200xBasicNet	439
	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData	625
	usbSetupPacket_t	626
	W2100_StimulusParametersNet	627
9	Class Index	
9.	1 Class List	
Нє	ere are the classes, structs, unions and interfaces with brief descriptions:	
	CW2100_FunctionNet::AudioChannelsNet	29
	BatteryState	29
	BesselFilterHighPassNet	30
	BesselFilterLowPassNet	31
	ButterworthFilterHighPassNet	31

9

ButterworthFilterLowPassNet	32
CChannelTestDeviceNet	33
CCMOSMea_FunctionNet	34
CCMOSMeaDeviceNet	45
CCreateFilterNet	48
CDacCalibrationFunctionNet 50	
CDacqGroupChannelGenericSelectionNet	52
CDacqGroupChannelSelectionNet	52
${\bf CDacqGroupChannelSelectionTemplateNet} < {\bf DacqGroupChannelEnumTemplateNet}, {\bf$	annelEnumTemplate
CDeviceGroupChannelInfoGenericNet	56
CDeviceGroupChannelInfoMEA2100_256Net	56
CDeviceGroupChannelInfoNet	57
CDeviceGroupChannelInfoSCUNet	57
${\bf CDeviceGroupChannelInfoTemplateNet} < {\bf DacqGroupChannelEnumTemplateNet} >$	58
CDeviceGroupChannelInfoW2100Net	59
CDigOutStimulatorFunctionNet CDigOutStimulatorFunctionNet is the class of the DigOut stimulator function class	59
CEncapsulatorDeviceNet CEncapsulatorDeviceNet is the to control the MCS HiClamp device	64
CExternDTesterDeviceNet CExternDTesterDeviceNet is the class to access the ExternD Tester (Handheld Device Tester D)	64
CFilterCoefficientsNet	66
CFilterConfigurationNet	68
CFilterConfigurationRegisterNet	70
CFilterPropertyNet	72
CFluidControlDeviceNet CFluidControlDeviceNet is the class to control MCS FluidControl (FCB and FCX) device	74
CFYIDeviceNet CFYIDeviceNet is the class to control the MCS FYI device	80
CGenericDevelopDeviceNet CGenericDevelopDeviceNet is the class to use during development of a new device	81
CGilsonDeviceNet CGilsonDeviceNet is the class to control a Gilson device	98
CGrapheneFunctionNet CGrapheneFunctionNet is the class to control the TEER device	100

CHiClampDeviceNet	100
CHiClampDeviceNet is the to control the MCS HiClamp device	109
CHLADacqNet	110
CHLADeviceNet CHLADeviceNet is the to control the MCS HLA device	111
CMcsUsbDacqNet::CHWInfo Class to provide hardware information about the device	112
CIntanMea_FunctionNet	115
CInterfaceboard2FunctionNet CInterfaceboard2FunctionNet is the class to control the Interfaceboard	117
CInterfaceboardFunctionNet CInterfaceboardFunctionNet is the class to control the Interfaceboard	119
CLIH3DeviceNet CLIH3DeviceNet is the class to access the HEKA LIH3 device	121
CMcsBus_AxisParametersNet	130
CMcsBus_ExtensionNet	132
CMcsBus_FYIExtensionNet	133
CMcsBus_MotorControlNet	135
CMcsBus_SensorNet	152
CMcsBus_TempSensorNet	162
CMcsBus_VoltageModeNet	164
CMcsBusNet	169
CMcsUsbDacqNet Base class for data acquisition devices	173
CMcsUsbDeviceStatePushFunctionNet	221
CMcsUsbDeviceStatePushNet	222
CMcsUsbFactoryNet	223
CMcsUsbFunctionNet	232
CMcsUsbFunctionPointerContainer	233
CMcsUsbListEntryNet McsUsbListEntryNet identifies a connected device	233
CMcsUsbListNet Class to handle a list of connected MCS USB devices	239
CMcsUsbNet Base class to handle MCS USB devices. All device classes are derived from this class. Functionality that is provided by all MCS devices is handled by this class	242
CMcsUsbPointerContainer	263

CMEA2100_256DacqGroupChannelSelectionNet	263
CMEA2100x256FunctionNet CMEA2100x256FunctionNet is the class to control the MEA2100-256 device needs #include "← Stg200xNet.h" to resolve documentation reference	264
CMeaAudioFunctionNet	265
CMeaCleanDeviceNet CMeaCleanDeviceNet is the class to access the MEA Clean device	268
CMeaCoatDeviceNet CMeaCoatDeviceNet is the class to access the MEA Coat device	272
CMeaDeviceNet Base class for MEA data acquisition devices	278
CMeaDigitalDataFunctionNet	286
CMeaFeedbackFunctionNet	288
CMealmpedanceDeviceNet	292
CMeasureTableDeviceNet CMeasureTableDeviceNet is the to control the MCS HLA device	293
CMeaSwitchDeviceNet The class to control the USB-MEA-Switch	294
CMeaUSBDeviceNet Class for data acquisition via ME and MEA USB amplifiers	296
CMeFunctionNet 297	
CMultiBatteryChargerDeviceNet CMultiBatteryChargerDeviceNet is the class to access the MBC-08 device	299
CMultiwellCallbackFunctionNet CMultiwellCallbackFunctionNet is the class to access the Multiwell-Mini-Stimulator	306
CMultiwellDeviceNet CMultiwellDeviceNet is the class to access the Multiwell device	308
CMultiwellOptoStimFunctionNet CMultiwellOptoStimFunctionNet is the class to access the optical properties of the Multiwell Optostim device	314
CNF_GenDeviceNet	319
COctoPotDeviceNet	320
COkuvisionStimulatorDeviceNet	324
CPatchServerDeviceNet CPatchServerDeviceNet is the class to control the MCS PatchServer device	327
CPathIdentDeviceNet	328
CPedoterDeviceNet 329	

CPeristalticPumpDeviceNet CPeristalticPumpDeviceNet is the class to control a Persistaltic Pump	331
CPgaDeviceNet	332
CPositionIIDeviceNet CPositionIIDeviceNet is the class to control PositionIII devices	334
CPositionImpDeviceNet CPositionImpDeviceNet is the class to access the Position/Imp devices	343
CPPCDeviceNet	346
CPPCFunctionNet CPPCFunctionNet is the class to access the PPC (high precision Patch Peristalic patch Pump	347
CPPS_DeviceNet	356
CPPS_FunctionNet	357
CPPSDeviceNet CPPS4plus1DeviceNet is the to control the MCS HLA device	361
CProgramPressureCurveNet CProgramPressureCurveNet is the class to program pressure curves	362
CPulseGeneratorFunctionNet CPulseGeneratorFunctionNet is the class to control the pulse generator for video tracking	363
CRadioControledDevicesNet	366
CCMOSMeaDeviceNet::CRegionOfInterestRect	368
CRetinaLedDeviceNet	369
CRFFunctionNet CRFFunctionNet is the class to control RF devices	371
CRobo_FYIProgram_FunctionNet	376
CRobo_FYITemp_FunctionNet	377
CRoboDacqNet	379
CRoboDeviceNet CRoboDeviceNet is the base class for all Robo platform based devices	393
CRoboFluidDeviceNet	407
CRobolnjectDeviceNet CRobolnjectDeviceNet is the to control the MCS Robolnject device	411
CRoboocyte2DeviceNet CRoboocyte2DeviceNet is the class to control the MCS Roboocyte2 device	411
CRoboStatorDeviceNet	413
CSafeISDeviceNet 419	
CSCUDacqGroupChannelSelectionNet	422

CSCUFunctionNet CSCUFunctionNet is the class to control the SCU device	423
CSerialPortNet	438
	430
CStg200xBasicNet Base class for the Stg200x	439
CStg200xDownloadBasicNet CStg200xDownloadBasicNet is the base class to control the download mode of the MCS STG device	469
CStg200xDownloadNet Main class for the STG download mode This class implements the STG download mode interface.	478
CStimulusFunctionNet	485
CSw2to64DeviceNet The class to control the MCS-USB-Sw2to64 device	496
CTcxDeviceNet Class to control a Temperature Controller (TCX)	499
CTEERFunctionNet CTEERFunctionNet is the class to control the TEER device	512
CTEERMachineDeviceNet	522
CUsbDeviceConfigurationFunctionNet CUsbDeviceConfigurationFunctionNet is the class to configure the USB firmware	523
CUsbExceptionNet Exception class that is thrown in case of an USB error	525
CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet	526
CW2100_FunctionNet	527
CW2100_StimulatorFunctionNet	534
CW2100DacqGroupChannelSelectionNet	542
CWarnerUssingDeviceNet CWarnerUssingDeviceNet is the class to control the Ussing device	543
CWarnerUssingFunctionNet CWarnerUssingFunctionNet is the class to control the Ussing device	544
CWarnerValveControllerDeviceNet CWarnerValveControllerDeviceNet is the class to access the Warner Valve Controller	560
CWarnerValveControllerDeviceTesterFunctionNet CWarnerValveControllerDeviceTesterFunctionNet is the class to access the functions for the Warner Valve Controller Device Tester	585
CWClassicFunctionNet	588
CWirelessBaseFunctionNet	593
DeviceIdNet Device Id	593

DigitalSource< digitalsourceenum >	595
DigitalSourceGeneral	596
DriverVersionNet Class gives firmware versions of the device's firmware destinations	598
FirmwareDestinationNames	603
HeadStageIDType	606
HeadstageIDTypeObject	609
HeadStageIDTypeState	610
mkfilterNet	611
CRoboDeviceNet::RoboMainLowLevelCommands	614
CRoboStatorDeviceNet::RoboMainStatorLowLevelCommands	621
CFilterCoefficientsNet::s_FilterAttributesNet	621
CMeaAudioFunctionNet::s_setaudionet	623
CStimulusFunctionNet::SidebandData	623
StgStatusNet	624
CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData	625
usbSetupPacket_t	626
W2100_StimulusParametersNet	627

10 Namespace Documentation

10.1 Mcs Namespace Reference

Namespaces

• Usb

10.2 Mcs::Usb Namespace Reference

Classes

- class BatteryState
- class BesselFilterHighPassNet
- class BesselFilterLowPassNet
- class ButterworthFilterHighPassNet
- class ButterworthFilterLowPassNet
- class CChannelTestDeviceNet
- class CCMOSMea_FunctionNet
- class CCMOSMeaDeviceNet

- class CCreateFilterNet
- class CDacCalibrationFunctionNet
- class CDacqGroupChannelGenericSelectionNet
- class CDacqGroupChannelSelectionNet
- class CDacqGroupChannelSelectionTemplateNet
- class CDeviceGroupChannelInfoGenericNet
- class CDeviceGroupChannelInfoMEA2100_256Net
- · class CDeviceGroupChannelInfoNet
- class CDeviceGroupChannelInfoSCUNet
- class CDeviceGroupChannelInfoTemplateNet
- class CDeviceGroupChannelInfoW2100Net
- · class CDigOutStimulatorFunctionNet

CDigOutStimulatorFunctionNet is the class of the DigOut stimulator function class.

· class CEncapsulatorDeviceNet

CEncapsulatorDeviceNet is the to control the MCS HiClamp device

class CExternDTesterDeviceNet

CExternDTesterDeviceNet is the class to access the ExternD Tester (Handheld Device Tester D)

- · class CFilterCoefficientsNet
- · class CFilterConfigurationNet
- · class CFilterConfigurationRegisterNet
- class CFilterPropertyNet
- class CFluidControlDeviceNet

CFluidControlDeviceNet is the class to control MCS FluidControl (FCB and FCX) device.

class CFYIDeviceNet

CFYIDeviceNet is the class to control the MCS FYI device

class CGenericDevelopDeviceNet

CGenericDevelopDeviceNet is the class to use during development of a new device.

class CGilsonDeviceNet

CGilsonDeviceNet is the class to control a Gilson device.

class CGrapheneFunctionNet

CGrapheneFunctionNet is the class to control the TEER device

class CHiClampDeviceNet

CHiClampDeviceNet is the to control the MCS HiClamp device

- · class CHLADacqNet
- class CHLADeviceNet

CHLADeviceNet is the to control the MCS HLA device

- · class CIntanMea FunctionNet
- · class CInterfaceboard2FunctionNet

CInterfaceboard2FunctionNet is the class to control the Interfaceboard

• class CInterfaceboardFunctionNet

CInterfaceboardFunctionNet is the class to control the Interfaceboard

class CLIH3DeviceNet

CLIH3DeviceNet is the class to access the HEKA LIH3 device.

- class CMcsBus_AxisParametersNet
- · class CMcsBus_ExtensionNet
- class CMcsBus_FYIExtensionNet
- class CMcsBus_MotorControlNet
- class CMcsBus_SensorNet
- class CMcsBus_TempSensorNet
- class CMcsBus_VoltageModeNet
- class CMcsBusNet
- class CMcsUsbDacqNet

Base class for data acquisition devices.

- class CMcsUsbDeviceStatePushFunctionNet
- · class CMcsUsbDeviceStatePushNet
- · class CMcsUsbFactoryNet
- class CMcsUsbFunctionNet
- class CMcsUsbFunctionPointerContainer
- class CMcsUsbListEntryNet

McsUsbListEntryNet identifies a connected device.

class CMcsUsbListNet

Class to handle a list of connected MCS USB devices.

class CMcsUsbNet

Base class to handle MCS USB devices. All device classes are derived from this class. Functionality that is provided by all MCS devices is handled by this class.

- class CMcsUsbPointerContainer
- class CMEA2100_256DacqGroupChannelSelectionNet
- class CMEA2100x256FunctionNet

CMEA2100x256FunctionNet is the class to control the MEA2100-256 device needs #include "Stg200xNet.h" to resolve documentation reference

- class CMeaAudioFunctionNet
- class CMeaCleanDeviceNet

CMeaCleanDeviceNet is the class to access the MEA Clean device.

class CMeaCoatDeviceNet

CMeaCoatDeviceNet is the class to access the MEA Coat device.

class CMeaDeviceNet

Base class for MEA data acquisition devices.

- class CMeaDigitalDataFunctionNet
- class CMeaFeedbackFunctionNet
- class CMealmpedanceDeviceNet
- class CMeasureTableDeviceNet

CMeasureTableDeviceNet is the to control the MCS HLA device

class CMeaSwitchDeviceNet

The class to control the USB-MEA-Switch.

class CMeaUSBDeviceNet

Class for data acquisition via ME and MEA USB amplifiers

- · class CMeFunctionNet
- · class CMultiBatteryChargerDeviceNet

CMultiBatteryChargerDeviceNet is the class to access the MBC-08 device.

class CMultiwellCallbackFunctionNet

CMultiwellCallbackFunctionNet is the class to access the Multiwell-Mini-Stimulator

· class CMultiwellDeviceNet

CMultiwellDeviceNet is the class to access the Multiwell device.

class CMultiwellOptoStimFunctionNet

CMultiwellOptoStimFunctionNet is the class to access the optical properties of the Multiwell Optostim device

- · class CNF GenDeviceNet
- · class COctoPotDeviceNet
- class COkuvisionStimulatorDeviceNet
- class CPatchServerDeviceNet

CPatchServerDeviceNet is the class to control the MCS PatchServer device

- class CPathIdentDeviceNet
- · class CPedoterDeviceNet
- class CPeristalticPumpDeviceNet

CPeristalticPumpDeviceNet is the class to control a Persistaltic Pump.

- · class CPgaDeviceNet
- · class CPositionIIDeviceNet

CPositionIIDeviceNet is the class to control PositionII devices

class CPositionImpDeviceNet

CPositionImpDeviceNet is the class to access the Position/Imp devices

- class CPPCDeviceNet
- · class CPPCFunctionNet

CPPCFunctionNet is the class to access the PPC (high precision Patch Peristalic patch Pump

- class CPPS DeviceNet
- class CPPS FunctionNet
- · class CPPSDeviceNet

CPPS4plus1DeviceNet is the to control the MCS HLA device

· class CProgramPressureCurveNet

CProgramPressureCurveNet is the class to program pressure curves

class CPulseGeneratorFunctionNet

CPulseGeneratorFunctionNet is the class to control the pulse generator for video tracking

- · class CRadioControledDevicesNet
- · class CRetinaLedDeviceNet
- class CRFFunctionNet

CRFFunctionNet is the class to control RF devices

- class CRobo_FYIProgram_FunctionNet
- class CRobo_FYITemp_FunctionNet
- class CRoboDacqNet
- · class CRoboDeviceNet

CRoboDeviceNet is the base class for all Robo platform based devices

- · class CRoboFluidDeviceNet
- · class CRobolnjectDeviceNet

CRobolnjectDeviceNet is the to control the MCS Robolnject device

class CRoboocyte2DeviceNet

CRoboocyte2DeviceNet is the class to control the MCS Roboocyte2 device

- · class CRoboStatorDeviceNet
- · class CSafeISDeviceNet
- class CSCUDacqGroupChannelSelectionNet
- · class CSCUFunctionNet

CSCUFunctionNet is the class to control the SCU device

- · class CSerialPortNet
- class CStg200xBasicNet

Base class for the Stg200x.

class CStg200xDownloadBasicNet

CStg200xDownloadBasicNet is the base class to control the download mode of the MCS STG device.

• class CStg200xDownloadNet

Main class for the STG download mode This class implements the STG download mode interface.

- class CStimulusFunctionNet
- · class CSw2to64DeviceNet

The class to control the MCS-USB-Sw2to64 device.

class CTcxDeviceNet

Class to control a Temperature Controller (TCX)

class CTEERFunctionNet

CTEERFunctionNet is the class to control the TEER device

- · class CTEERMachineDeviceNet
- class CUsbDeviceConfigurationFunctionNet

CUsbDeviceConfigurationFunctionNet is the class to configure the USB firmware

class CUsbExceptionNet

Exception class that is thrown in case of an USB error.

- class CW2100 FunctionNet
- class CW2100 StimulatorFunctionNet
- class CW2100DacqGroupChannelSelectionNet
- class CWarnerUssingDeviceNet

CWarnerUssingDeviceNet is the class to control the Ussing device

class CWarnerUssingFunctionNet

CWarnerUssingFunctionNet is the class to control the Ussing device

class CWarnerValveControllerDeviceNet

CWarnerValveControllerDeviceNet is the class to access the Warner Valve Controller

class CWarnerValveControllerDeviceTesterFunctionNet

CWarnerValveControllerDeviceTesterFunctionNet is the class to access the functions for the Warner Valve Controller Device Tester

- · class CWClassicFunctionNet
- class CWirelessBaseFunctionNet
- struct DeviceIdNet

Device Id.

- class DigitalSource
- · class DigitalSourceGeneral
- · class DriverVersionNet

Class gives firmware versions of the device's firmware destinations.

- class FirmwareDestinationNames
- · class HeadStageIDType
- · class HeadstageIDTypeObject
- class HeadStageIDTypeState
- class mkfilterNet
- class StgStatusNet
- class usbSetupPacket_t
- struct W2100_StimulusParametersNet

Enumerations

```
    enum enCMosMeaChipType {
        unknown = 0,
        nMos16LV = 1,
        nMos32LV = 3,
        nMos36LN = 6,
        nMos64LN = 7 }
    enum EnSTG200x_STATUS {
        OK,
        NOT_CONNECTED,
        DEVICE_NOT_FOUND }
```

Functions

- public delegate void OnMcsUsbDeviceState (usbSetupPacket t[∧] request)
- private delegate void OnMcsUsbDeviceStateCallback (IntPtr pThis, uint32_t size, IntPtr buffer)
- public delegate void OnUpdateFirmwareStatusChange (String[^])
- public delegate void OnUpdateFirmwareProgress (int)
- public delegate void OnDeviceArrivalRemoval (CMcsUsbListEntryNet^ entry)

Delegate to show a device arrival or removal.

- public delegate void OnStgPollStatus (unsigned int status, StgStatusNet[^] stgStatusNet, array< int >[^] index_list)
- public delegate void OnMwPollStatus (unsigned int CurrentTemp, unsigned int PlateState, unsigned int SwitchState)
- public delegate void RoboStatusEventDelegate (array< unsigned char >^ buffer)
- public delegate void OnStg200xDataHandler (uint32_t trigger)
- public delegate void OnStg200xErrorHandler ()
- public delegate void OnChannelData (CMcsUsbDacqNet^ dacq, int CbHandle, int numFrames)
- public delegate void OnError (String[^] msg, int action)

10.2.1 Enumeration Type Documentation

10.2.1.1 enCMosMeaChipType enum enCMosMeaChipType [strong]

Enumerator

unknown	
nMos16LV	
nMos32LV	
nMos36LN	
nMos64LN	

10.2.1.2 EnSTG200x_STATUS enum EnSTG200x_STATUS [strong]

Enumerator

OK	
NOT_CONNECTED	
DEVICE_NOT_FOUND	

10.2.2 Function Documentation

```
10.2.2.2 OnDeviceArrivalRemoval() public delegate void Mcs::Usb::OnDeviceArrivalRemoval (
             CMcsUsbListEntryNet^ entry )
Delegate to show a device arrival or removal.
10.2.2.3 OnError() public delegate void Mcs::Usb::OnError (
             String^{\wedge} msg,
             int action )
10.2.2.4 OnMcsUsbDeviceState() public delegate void OnMcsUsbDeviceState (
             usbSetupPacket_t^ request )
10.2.2.5 OnMcsUsbDeviceStateCallback() private delegate void OnMcsUsbDeviceStateCallback (
             IntPtr pThis,
             uint32_t size,
             IntPtr buffer )
10.2.2.6 OnMwPollStatus() public delegate void Mcs::Usb::OnMwPollStatus (
             unsigned int CurrentTemp,
             unsigned int PlateState,
             unsigned int SwitchState )
10.2.2.7 OnStg200xDataHandler() public delegate void Mcs::Usb::OnStg200xDataHandler (
             uint32_t trigger )
10.2.2.8 OnStg200xErrorHandler() public delegate void Mcs::Usb::OnStg200xErrorHandler ( )
10.2.2.9 OnStgPollStatus() public delegate void Mcs::Usb::OnStgPollStatus (
             unsigned int status,
             StgStatusNet^ stgStatusNet,
             array < int >^  index_list )
```

11 Class Documentation 29

```
10.2.2.10 OnUpdateFirmwareProgress() public delegate void Mcs::Usb::OnUpdateFirmwareProgress
             int )
10.2.2.11 OnUpdateFirmwareStatusChange() public delegate void Mcs::Usb::OnUpdateFirmware←
StatusChange (
            String^ )
10.2.2.12 RoboStatusEventDelegate() public delegate void Mcs::Usb::RoboStatusEventDelegate (
             array< unsigned char >^{\wedge} buffer )
11
     Class Documentation
11.1 CW2100_FunctionNet::AudioChannelsNet Struct Reference
Public Attributes

    W2100DacqGroupChannelEnumNet dacqgroup

    int channel

   · int amplification
11.1.1 Member Data Documentation
11.1.1.1 amplification int amplification
11.1.1.2 channel int channel
11.1.1.3 dacqgroup W2100DacqGroupChannelEnumNet dacqgroup
11.2 BatteryState Class Reference
Properties
   • unsigned int Charge [get]
   • unsigned int Voltage [get]
   • System::String ChargeString [get]
   • System::String^ ChargeRegionString [get]
   • System::String^ VoltageString [get]
```

11.2.1 Property Documentation

11.2.1.1 Charge unsigned int Charge [get]

11.2.1.2 ChargeRegionString System:: String^ ChargeRegionString [get]


11.2.1.3 ChargeString System:: String ChargeString [get]

11.2.1.4 Voltage unsigned int Voltage [get]

11.2.1.5 VoltageString System:: String^ VoltageString [get]

11.3 BesselFilterHighPassNet Class Reference

Inheritance diagram for BesselFilterHighPassNet:



Public Member Functions

• BesselFilterHighPassNet (int numCoefSets, int order, double sampleRate, double cutoffFrequency, double scale)

Additional Inherited Members


11.3.1 Constructor & Destructor Documentation

11.3.1.1 BesselFilterHighPassNet() BesselFilterHighPassNet (

```
int numCoefSets,
int order,
double sampleRate,
double cutoffFrequency,
double scale )
```

11.4 BesselFilterLowPassNet Class Reference

Inheritance diagram for BesselFilterLowPassNet:



Public Member Functions

 BesselFilterLowPassNet (int numCoefSets, int order, double sampleRate, double cutoffFrequency, double scale)

Additional Inherited Members


11.4.1 Constructor & Destructor Documentation

11.4.1.1 BesselFilterLowPassNet() BesselFilterLowPassNet (

```
int numCoefSets,
int order,
double sampleRate,
double cutoffFrequency,
double scale )
```

11.5 ButterworthFilterHighPassNet Class Reference

Inheritance diagram for ButterworthFilterHighPassNet:



Public Member Functions

ButterworthFilterHighPassNet (int numCoefSets, int order, double sampleRate, double cutoffFrequency, double scale)

Additional Inherited Members


11.5.1 Constructor & Destructor Documentation

11.5.1.1 ButterworthFilterHighPassNet() ButterworthFilterHighPassNet (

```
int numCoefSets,
int order,
double sampleRate,
double cutoffFrequency,
double scale )
```

11.6 ButterworthFilterLowPassNet Class Reference

Inheritance diagram for ButterworthFilterLowPassNet:



Public Member Functions

ButterworthFilterLowPassNet (int numCoefSets, int order, double sampleRate, double cutoffFrequency, double scale)

Additional Inherited Members


11.6.1 Constructor & Destructor Documentation

11.6.1.1 ButterworthFilterLowPassNet() ButterworthFilterLowPassNet (

```
int numCoefSets,
int order,
double sampleRate,
double cutoffFrequency,
double scale )
```

11.7 CChannelTestDeviceNet Class Reference

Inheritance diagram for CChannelTestDeviceNet:



Public Member Functions

- CChannelTestDeviceNet ()
- ∼CChannelTestDeviceNet ()
- void SetWaveform (unsigned int Waveform)
- void SetAmplitude (unsigned int Amplitude)
- void SetFrequency (unsigned int Frequency)
- · void SetAttenuation (unsigned int Attenuation)

Additional Inherited Members

11.7.1 Constructor & Destructor Documentation


```
11.7.1.1 CChannelTestDeviceNet() CChannelTestDeviceNet ( )
```

```
11.7.1.2 \sim CChannelTestDeviceNet() \sim CChannelTestDeviceNet ( )
```

11.7.2 Member Function Documentation

11.8 CCMOSMea FunctionNet Class Reference

Inheritance diagram for CCMOSMea_FunctionNet:



Public Member Functions

- CCMOSMea_FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] cMOSMea_←
 FunctionPointerContainer)
- CCMOSMea_FunctionNet (CMcsUsbNet^ mcsusb)
- void SetADCInputOffset (int32_t offset)
- int32_t GetADCInputOffset ()
- void SetSourceDrain (int32_t voltage)
- int32 t GetSourceDrain ()
- void SetSourceGate (int32_t voltage)
- int32_t GetSourceGate ()
- void SetSourceBulk (int32_t voltage)
- int32 t GetSourceBulk ()
- void SetGate (int32_t voltage)
- int32 t GetGate ()
- void SetBath (int32_t voltage)
- int32_t GetBath ()
- int32 t GetGNDI ()
- int32_t GetVDDI ()
- int32_t GetVDD3I ()
- void UpdateTransistorVoltages ()
- bool AreTransistorVoltagesSet ()
- void PowerChip (bool on)
- bool IsChipPowered ()
- enCMosMeaChipType DetectChipType ()
- void SetGateToVOP ()
- void SetGateFloating ()
- bool IsGateFloating ()
- void VOPSTimerSetResetTimes (uint32_t ResetTime, uint32_t IntervalTime)
- void VOPSTimerSetResetTimes (uint32_t ResetTime, uint32_t IntervalTime, uint32_t HPFilterResetTime)
- void SetBathMode (CMOSMeaBathModeEnumNet Mode)
- CMOSMeaBathModeEnumNet GetBathMode ()
- void SetNeurochipMemoryData (uint16_t MemAddress, uint32_t MemData)

- void SetNeurochipMemoryData (uint16_t MemAddress, array< uint32_t >^ MemData)
- uint32_t GetNeurochipMemoryData (uint16_t MemAddress)
- array< uint32 t > ^ GetNeurochipMemoryData (uint16 t MemAddress, uint32 t RegestLength)
- uint32 t GetNeurochipMemorySize ()
- uint32_t GetMaxNumOfColumns (uint32_t Samplerate)
- void SetStimulusSites (List< int16_t >^ SwitchPosition)
- List< int16 t > ^ GetStimulusSites ()
- void ClearSTGOutput (uint32 t Channel)
- uint32 t GetNumberOfSupportedGroups ()
- uint32 t GetNumberOfSupportedGroups (uint32 t virtualDevice)
- DacqGroupChannelEnumNet GetGroupID (uint32 t Index)
- DacqGroupChannelEnumNet GetGroupID (uint32_t Index, uint32_t virtualDevice)
- uint32_t GetGroupNumberOfChannels (DacqGroupChannelEnumNet GroupID)
- uint32_t GetGroupNumberOfChannels (DacqGroupChannelEnumNet GroupID, uint32_t virtualDevice)
- DacqMeaGroupTypeEnumNet GetGroupType (DacqGroupChannelEnumNet GroupID)
- void EnableChannelsInGroup (DacqGroupChannelEnumNet GroupID, List< bool >^ EnabledChannelsBit← Map)
- void EnableChannelsInGroup (DacqGroupChannelEnumNet GroupID, List< bool >[^] EnabledChannelsBit←
 Map, uint32 t virtualDevice)
- List< bool > ^ GetEnabledChannelsInGroup (DacqGroupChannelEnumNet GroupID)
- List< bool > ^ GetEnabledChannelsInGroup (DacqGroupChannelEnumNet GroupID, uint32_t virtualDevice)
- SampleSizeNet GetGroupSampleSize (DacqGroupChannelEnumNet GroupID)
- SampleSizeNet GetGroupSampleSize (DacqGroupChannelEnumNet GroupID, uint32_t virtualDevice)
- uint32 t GetGroupResolutionPerDigit (DacqGroupChannelEnumNet GroupID)
- uint32 t GetGroupResolutionPerDigit (DacqGroupChannelEnumNet GroupID, uint32 t virtualDevice)
- CMOSMeaValueUnitEnumNet GetGroupUnit (DacqGroupChannelEnumNet GroupID)
- CMOSMeaValueUnitEnumNet GetGroupUnit (DacqGroupChannelEnumNet GroupID, uint32 t virtualDevice)
- int32_t GetGroupDCOffset (DacqGroupChannelEnumNet GroupID)
- int32_t GetGroupDCOffset (DacqGroupChannelEnumNet GroupID, uint32_t virtualDevice)
- int32_t GetGroupADCBits (DacqGroupChannelEnumNet GroupID)
- int32 t GetGroupADCBits (DacqGroupChannelEnumNet GroupID, uint32 t virtualDevice)
- uint32_t GetGroupChannelBitmaskBySelect (DacqGroupChannelEnumNet GroupID, uint32_t Channel← Number)
- uint32_t GetGroupChannelBitmaskBySelect (DacqGroupChannelEnumNet GroupID, uint32_t Channel
 — Number, uint32_t virtualDevice)
- CMOSMeaInterfaceADCEnumNet GetGroupChannelBitmaskInterfaceADC (uint32_t ChannelNumber)
- CMOSMeaInterfaceADCEnumNet GetGroupChannelBitmaskInterfaceADC (uint32_t ChannelNumber, uint32_t virtualDevice)
- CMOSMealFDigChannelEnumNet GetGroupChannelBitmaskIFDigChannels (uint32 t ChannelNumber)
- CMOSMealFDigChannelEnumNet GetGroupChannelBitmasklFDigChannels (uint32_t ChannelNumber, uint32_t virtualDevice)
- CMOSMeaHeadstage1NCBathCurrentEnumNet GetGroupChannelBitmaskHS1NCBathCurrent (uint32_t ChannelNumber)
- CMOSMeaHeadstage1NCBathCurrentEnumNet GetGroupChannelBitmaskHS1NCBathCurrent (uint32_t ChannelNumber, uint32_t virtualDevice)
- CMOSMeaHeadstage1NCCol2CurrentEnumNet GetGroupChannelBitmaskHS1NCCol2Current (uint32_t ChannelNumber)
- CMOSMeaHeadstage1NCCol2CurrentEnumNet GetGroupChannelBitmaskHS1NCCol2Current (uint32_t ChannelNumber, uint32_t virtualDevice)
- CMOSMeaHeadstage1NChipTempEnumNet GetGroupChannelBitmaskHS1NChipTemp (uint32_t Channel
 — Number)
- CMOSMeaHeadstage1NChipTempEnumNet GetGroupChannelBitmaskHS1NChipTemp (uint32_t Channel
 — Number, uint32_t virtualDevice)

- CMOSMeaSTG1DACSignalEnumNet GetGroupChannelBitmaskSTG1DACSignal (uint32_t Channel↔ Number)
- CMOSMeaSTG1DACSignalEnumNet GetGroupChannelBitmaskSTG1DACSignal (uint32_t Channel
 — Number, uint32_t virtualDevice)
- CMOSMeaHS1SidebandEnumNet GetGroupChannelBitmaskHS1Sidebands (uint32_t ChannelNumber)
- CMOSMeaHS1SidebandEnumNet GetGroupChannelBitmaskHS1Sidebands (uint32_t ChannelNumber, uint32_t virtualDevice)
- CMOSMeaHS1TriggerStatusEnumNet GetGroupChannelBitmaskHS1TriggerStatus (uint32_t Channel
 — Number)
- CMOSMeaHS1TriggerStatusEnumNet GetGroupChannelBitmaskHS1TriggerStatus (uint32_t Channel Number, uint32_t virtualDevice)
- CMOSMeaPacketFrameContextGroupEnumNet GetGroupChannelBitmaskPacketFrameContext (uint32_← t ChannelNumber)
- CMOSMeaPacketFrameContextGroupEnumNet GetGroupChannelBitmaskPacketFrameContext (uint32_← t ChannelNumber, uint32_t virtualDevice)

Additional Inherited Members

11.8.1 Constructor & Destructor Documentation

```
11.8.2.2 ClearSTGOutput() void ClearSTGOutput ( uint32_t Channel)
```

11.8.2.3 DetectChipType() enCMosMeaChipType DetectChipType ()

```
11.8.2.4 EnableChannelsInGroup() [1/2] void EnableChannelsInGroup (
              DacqGroupChannelEnumNet GroupID,
              List<br/>< bool >^ EnabledChannelsBitMap )
\textbf{11.8.2.5} \quad \textbf{EnableChannelsInGroup() [2/2]} \quad \texttt{void EnableChannelsInGroup ()}
              DacqGroupChannelEnumNet GroupID,
              List < bool >^{\wedge} EnabledChannelsBitMap,
              uint32_t virtualDevice )
11.8.2.6 GetADCInputOffset() int32_t GetADCInputOffset ( )
11.8.2.7 GetBath() int32_t GetBath ()
11.8.2.8 GetBathMode() CMOSMeaBathModeEnumNet GetBathMode ( )
11.8.2.9 GetEnabledChannelsInGroup() [1/2] List<br/>bool> ^ GetEnabledChannelsInGroup (
              DacqGroupChannelEnumNet GroupID )
11.8.2.10 GetEnabledChannelsInGroup() [2/2] List<bool> ^ GetEnabledChannelsInGroup (
              DacqGroupChannelEnumNet GroupID,
              uint32_t virtualDevice )
11.8.2.11 GetGate() int32_t GetGate ()
11.8.2.12 GetGNDI() int32_t GetGNDI ()
11.8.2.13 GetGroupADCBits() [1/2] int32_t GetGroupADCBits (
              {\tt DacqGroupChannelEnumNet} \  \, \textit{GroupID} \  \, )
```

```
11.8.2.14 GetGroupADCBits() [2/2] int32_t GetGroupADCBits (
            DacqGroupChannelEnumNet GroupID,
            uint32_t virtualDevice )
11.8.2.15 GetGroupChannelBitmaskBySelect() [1/2] uint32_t GetGroupChannelBitmaskBySelect (
            DacqGroupChannelEnumNet GroupID,
            uint32_t ChannelNumber )
11.8.2.16 GetGroupChannelBitmaskBySelect() [2/2] uint32_t GetGroupChannelBitmaskBySelect (
            DacqGroupChannelEnumNet GroupID,
            uint32_t ChannelNumber,
            uint32_t virtualDevice )
11.8.2.17 GetGroupChannelBitmaskHS1NCBathCurrent() [1/2] CMOSMeaHeadstage1NCBathCurrentEnum
Net GetGroupChannelBitmaskHS1NCBathCurrent (
            uint32_t ChannelNumber )
11.8.2.18 GetGroupChannelBitmaskHS1NCBathCurrent() [2/2] CMOSMeaHeadstagelNCBathCurrentEnum↔
Net GetGroupChannelBitmaskHS1NCBathCurrent (
            uint32_t ChannelNumber,
            uint32_t virtualDevice )
11.8.2.19 GetGroupChannelBitmaskHS1NCCol2Current() [1/2] CMOSMeaHeadstagelNCCol2CurrentEnum←
Net GetGroupChannelBitmaskHS1NCCol2Current (
            uint32_t ChannelNumber )
11.8.2.20 GetGroupChannelBitmaskHS1NCCol2Current() [2/2] CMOSMeaHeadstage1NCCol2CurrentEnum←
Net GetGroupChannelBitmaskHS1NCCol2Current (
            uint32_t ChannelNumber,
            uint32_t virtualDevice )
11.8.2.21 GetGroupChannelBitmaskHS1NChipTemp() [1/2] CMOSMeaHeadstage1NChipTempEnumNet Get↔
GroupChannelBitmaskHS1NChipTemp (
            uint32_t ChannelNumber )
```

```
11.8.2.22 GetGroupChannelBitmaskHS1NChipTemp() [2/2] CMOSMeaHeadstagelNChipTempEnumNet Get←
GroupChannelBitmaskHS1NChipTemp (
                                      uint32_t ChannelNumber,
                                       uint32_t virtualDevice )
11.8.2.23 GetGroupChannelBitmaskHS1Sidebands() [1/2] CMOSMeaHS1SidebandEnumNet GetGroup↔
ChannelBitmaskHS1Sidebands (
                                       uint32_t ChannelNumber )
\textbf{11.8.2.24} \quad \textbf{GetGroupChannelBitmaskHS1Sidebands() [2/2]} \quad \texttt{CMOSMeaHS1SidebandEnumNet} \quad \texttt{GetGroup} \leftarrow \texttt{CMOSMeaHS1SidebandEnumNet} \quad \texttt{CMOSMeaHS1Si
ChannelBitmaskHS1Sidebands (
                                       uint32_t ChannelNumber,
                                       uint32_t virtualDevice )
11.8.2.25 GetGroupChannelBitmaskHS1TriggerStatus() [1/2] CMOSMeaHS1TriggerStatusEnumNet Get↔
GroupChannelBitmaskHS1TriggerStatus (
                                       uint32_t ChannelNumber )
11.8.2.26 GetGroupChannelBitmaskHS1TriggerStatus() [2/2] CMOSMeaHS1TriggerStatusEnumNet Get↔
GroupChannelBitmaskHS1TriggerStatus (
                                       uint32_t ChannelNumber,
                                       uint32_t virtualDevice )
11.8.2.27 GetGroupChannelBitmaskIFDigChannels() [1/2] CMOSMeaIFDigChannelEnumNet GetGroup↔
ChannelBitmaskIFDigChannels (
                                     uint32_t ChannelNumber )
11.8.2.28 GetGroupChannelBitmaskIFDigChannels() [2/2] CMOSMeaIFDigChannelEnumNet GetGroup↔
ChannelBitmaskIFDigChannels (
                                       uint32_t ChannelNumber,
                                       uint32_t virtualDevice )
11.8.2.29 GetGroupChannelBitmaskInterfaceADC() [1/2] CMOSMeaInterfaceADCEnumNet GetGroup↔
ChannelBitmaskInterfaceADC (
                                       uint32_t ChannelNumber )
```

```
11.8.2.30 GetGroupChannelBitmaskInterfaceADC() [2/2] CMOSMeaInterfaceADCEnumNet GetGroup←
ChannelBitmaskInterfaceADC (
             uint32_t ChannelNumber,
             uint32_t virtualDevice )
11.8.2.31 GetGroupChannelBitmaskPacketFrameContext() [1/2] CMOSMeaPacketFrameContextGroup↔
EnumNet GetGroupChannelBitmaskPacketFrameContext (
              uint32_t ChannelNumber )
\textbf{11.8.2.32} \quad \textbf{GetGroupChannelBitmaskPacketFrameContext()} \  \  \texttt{[2/2]} \quad \texttt{CMOSMeaPacketFrameContextGroup} \leftarrow \\
{\tt EnumNet~GetGroupChannelBitmaskPacketFrameContext~(}
              uint32_t ChannelNumber,
              uint32_t virtualDevice )
11.8.2.33 GetGroupChannelBitmaskSTG1DACSignal() [1/2] CMOSMeaSTG1DACSignalEnumNet GetGroup←
ChannelBitmaskSTG1DACSignal (
              uint32_t ChannelNumber )
11.8.2.34 GetGroupChannelBitmaskSTG1DACSignal() [2/2] CMOSMeaSTG1DACSignalEnumNet GetGroup↔
ChannelBitmaskSTG1DACSignal (
             uint32_t ChannelNumber,
              uint32_t virtualDevice )
11.8.2.35 GetGroupDCOffset() [1/2] int32_t GetGroupDCOffset (
              DacqGroupChannelEnumNet GroupID )
\textbf{11.8.2.36} \quad \textbf{GetGroupDCOffset() [2/2]} \quad \texttt{int32\_t GetGroupDCOffset ()}
              DacqGroupChannelEnumNet GroupID,
              uint32_t virtualDevice )
11.8.2.37 GetGroupID() [1/2] DacqGroupChannelEnumNet GetGroupID (
              uint32_t Index )
```

```
11.8.2.38 GetGroupID() [2/2] DacqGroupChannelEnumNet GetGroupID (
             uint32_t Index,
             uint32_t virtualDevice )
11.8.2.39 GetGroupNumberOfChannels() [1/2] uint32_t GetGroupNumberOfChannels (
             DacqGroupChannelEnumNet GroupID )
11.8.2.40 GetGroupNumberOfChannels() [2/2] uint32_t GetGroupNumberOfChannels (
             DacqGroupChannelEnumNet GroupID,
             uint32_t virtualDevice )
11.8.2.41 GetGroupResolutionPerDigit() [1/2] uint32_t GetGroupResolutionPerDigit (
             DacqGroupChannelEnumNet GroupID )
11.8.2.42 GetGroupResolutionPerDigit() [2/2] uint32_t GetGroupResolutionPerDigit (
             DacqGroupChannelEnumNet GroupID,
             uint32_t virtualDevice )
11.8.2.43 GetGroupSampleSize() [1/2] SampleSizeNet GetGroupSampleSize (
             DacqGroupChannelEnumNet GroupID )
11.8.2.44 GetGroupSampleSize() [2/2] SampleSizeNet GetGroupSampleSize (
             DacqGroupChannelEnumNet GroupID,
             uint32_t virtualDevice )
11.8.2.45 GetGroupType() [1/2] DacqMeaGroupTypeEnumNet GetGroupType (
             {\tt DacqGroupChannelEnumNet}\ \textit{GroupID}\ )
11.8.2.46 GetGroupType() [2/2] DacqMeaGroupTypeEnumNet GetGroupType (
             DacqGroupChannelEnumNet GroupID,
             uint32_t virtualDevice )
```

```
11.8.2.47 GetGroupUnit() [1/2] CMOSMeaValueUnitEnumNet GetGroupUnit (
             DacqGroupChannelEnumNet GroupID )
11.8.2.48 GetGroupUnit() [2/2] CMOSMeaValueUnitEnumNet GetGroupUnit (
             DacqGroupChannelEnumNet GroupID,
             uint32_t virtualDevice )
11.8.2.49 GetMaxNumOfColumns() uint32_t GetMaxNumOfColumns (
             uint32_t Samplerate )
\textbf{11.8.2.50} \quad \textbf{GetNeurochipMemoryData() [1/2]} \quad \texttt{uint} \ \texttt{32\_t GetNeurochipMemoryData} \ \ \textbf{(}
             uint16_t MemAddress )
11.8.2.51 GetNeurochipMemoryData() [2/2] array<uint32_t> ^ GetNeurochipMemoryData (
             uint16_t MemAddress,
             uint32_t ReqestLength )
11.8.2.52 GetNeurochipMemorySize() uint32_t GetNeurochipMemorySize ( )
11.8.2.53 GetNumberOfSupportedGroups() [1/2] uint32_t GetNumberOfSupportedGroups ( )
11.8.2.54 GetNumberOfSupportedGroups() [2/2] uint32_t GetNumberOfSupportedGroups (
             uint32_t virtualDevice )
11.8.2.55 GetSourceBulk() int32_t GetSourceBulk ( )
11.8.2.56 GetSourceDrain() int32_t GetSourceDrain ( )
```

```
11.8.2.57 GetSourceGate() int32_t GetSourceGate ( )
11.8.2.58 GetStimulusSites() List<int16_t> ^{\land} GetStimulusSites ( )
11.8.2.59 GetVDD3I() int32_t GetVDD3I ()
11.8.2.60 GetVDDI() int32_t GetVDDI ()
11.8.2.61 IsChipPowered() bool IsChipPowered ( )
\textbf{11.8.2.62} \quad \textbf{IsGateFloating()} \quad \texttt{bool IsGateFloating ()} \\
11.8.2.63 PowerChip() void PowerChip (
             bool on )
11.8.2.64 SetADCInputOffset() void SetADCInputOffset (
              int32_t offset )
11.8.2.65 SetBath() void SetBath (
              int32_t voltage )
11.8.2.66 SetBathMode() void SetBathMode (
              CMOSMeaBathModeEnumNet Mode )
```

```
11.8.2.67 SetGate() void SetGate (
             int32_t voltage )
11.8.2.68 SetGateFloating() void SetGateFloating ( )
11.8.2.69 SetGateToVOP() void SetGateToVOP ()
11.8.2.70 SetNeurochipMemoryData() [1/2] void SetNeurochipMemoryData (
             uint16_t MemAddress,
              array < uint32_t >^{\land} MemData)
11.8.2.71 SetNeurochipMemoryData() [2/2] void SetNeurochipMemoryData (
             uint16_t MemAddress,
             uint32_t MemData )
11.8.2.72 SetSourceBulk() void SetSourceBulk (
             int32_t voltage )
11.8.2.73 SetSourceDrain() void SetSourceDrain (
             int32_t voltage )
11.8.2.74 SetSourceGate() void SetSourceGate (
             int32_t voltage )
11.8.2.75 SetStimulusSites() void SetStimulusSites (
             List< int16_t >^{\land} SwitchPosition )
\textbf{11.8.2.76} \quad \textbf{UpdateTransistorVoltages()} \quad \texttt{void UpdateTransistorVoltages ()}
```


```
11.8.2.77 VOPSTimerSetResetTimes() [1/2] void VOPSTimerSetResetTimes ( uint32_t ResetTime, uint32_t IntervalTime )
```

11.8.2.78 VOPSTimerSetResetTimes() [2/2] void VOPSTimerSetResetTimes (uint32_t ResetTime,

uint32_t IntervalTime,
uint32_t HPFilterResetTime)

11.9 CCMOSMeaDeviceNet Class Reference

Inheritance diagram for CCMOSMeaDeviceNet:



Classes

class CRegionOfInterestRect

Public Member Functions

- CCMOSMeaDeviceNet (void)
- ∼CCMOSMeaDeviceNet ()
- void SetBaseSamplerate (int BaseSamplerate)
- int GetBaseSamplerate ()
- array< int > ^ GetAvailableBaseSamplerates ()
- void SetRegionOfInterests (System::Collections::Generic::Dictionary< int, CRegionOfInterestRect^>^ rois)
- void UpdateChannelBlock (int queuesize, int threshold, int channels_in_block)
- System::Collections::Generic::Dictionary< int, array< array< int16_t >^> ^ GetCMOSDataDictionary (int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
- System::Collections::Generic::Dictionary< int, array< uint16_t >^> ^ GetChannelDataUI16 (DacqGroup← ChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
- System::Collections::Generic::Dictionary< int, array< int16_t >^> ^ GetChannelDatal16 (DacqGroup← ChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
- System::Collections::Generic::Dictionary< int, array< uint32_t >^> ^ GetChannelDataUl32 (DacqGroup← ChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames ret)
- System::Collections::Generic::Dictionary< int, array< int32_t >^> ^ GetChannelDatal32 (DacqGroup← ChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Properties

```
• CCMOSMea_FunctionNet^ CMosMea [get]
```

```
• CStimulusFunctionNet^ Stimulus [get]
```

Additional Inherited Members

```
11.9.1 Constructor & Destructor Documentation
```

```
11.9.1.1 CCMOSMeaDeviceNet() CCMOSMeaDeviceNet (
            void )
11.9.1.2 ~CCMOSMeaDeviceNet() ~CCMOSMeaDeviceNet ()
11.9.2 Member Function Documentation
11.9.2.1 GetAvailableBaseSamplerates() array<int> ^ GetAvailableBaseSamplerates ()
11.9.2.2 GetBaseSamplerate() int GetBaseSamplerate ( )
11.9.2.3 GetChannelDatal16() System::Collections::Generic::Dictionary<int, array<int16_t>^> ^
GetChannelDataI16 (
            DacqGroupChannelEnumNet group,
            int frames,
            [System::Runtime::InteropServices::Out] int % frames_ret )
11.9.2.4 GetChannelDatal32() System::Collections::Generic::Dictionary<int, array<int32_t>^> ^
GetChannelDataI32 (
```

DacqGroupChannelEnumNet group,


[System::Runtime::InteropServices::Out] int % frames_ret)

int frames,

```
11.9.2.5 GetChannelDataUl16() System::Collections::Generic::Dictionary<int, array<uint16_t>^>
^ GetChannelDataUI16 (
             DacqGroupChannelEnumNet group,
             int frames,
             [System::Runtime::InteropServices::Out] int % frames_ret )
11.9.2.6 GetChannelDataUl32() System::Collections::Generic::Dictionary<int, array<uint32_t>^>
^ GetChannelDataUI32 (
             DacqGroupChannelEnumNet group,
             int frames,
             [System::Runtime::InteropServices::Out] int % frames_ret )
11.9.2.7 GetCMOSDataDictionary() System::Collections::Generic::Dictionary<int, array<array<int16↔
_{t}^{\wedge} \(^{\} \) \(^{\} \) GetCMOSDataDictionary (
             int frames,
             [System::Runtime::InteropServices::Out] int % frames_ret )
11.9.2.8 SetBaseSamplerate() void SetBaseSamplerate (
             int BaseSamplerate )
11.9.2.9 SetRegionOfInterests() void SetRegionOfInterests (
             {\tt System::Collections::Generic::Dictionary<~int,~CRegionOfInterestRect^>^~\it rois~)}
11.9.2.10 UpdateChannelBlock() void UpdateChannelBlock (
             int queuesize,
             int threshold,
             int channels_in_block )
11.9.3 Property Documentation
11.9.3.1 CMosMea CCMOSMea_FunctionNet^ CMosMea [get]
11.9.3.2 Stimulus CStimulusFunctionNet^ Stimulus [get]
```

11.10 CCreateFilterNet Class Reference

Inheritance diagram for CCreateFilterNet:



Public Member Functions

- · CCreateFilterNet (int numCoefSets, int order, double sampleRate, double cutoffFrequency, double scale)
- ∼CCreateFilterNet ()
- CFilterCoefficientsNet ^ GetBiQuad (int index)
- array< CFilterCoefficientsNet[^]> [^] GetBiQuads ()

Static Public Member Functions

- static int FindFilter (array< CFilterCoefficientsNet^>^ coef, array< CCreateFilterNet^>^ param)
- static int FindFilter (array< array< uint64_t >^> coef, array< CCreateFilterNet^> param, CFilterCoefficientsNet::s_FilterAttributesNet^ FiltAttr, bool DoMCSLegacyCompare)

Protected Member Functions

• CCreateFilterNet (int numCoefSets, CCreateFilter *pCreateFilter)

Properties

- int NumCoefSets [get]
- int Order [get]
- double SampleRate [get]
- double CutoffFrequency [get]
- double Scale [get]

11.10.1 Constructor & Destructor Documentation

11.10.1.1 CCreateFilterNet() [1/2] CCreateFilterNet (

```
int numCoefSets,
int order,
double sampleRate,
double cutoffFrequency,
double scale )
```

```
11.10.1.2 ~CCreateFilterNet() ~CCreateFilterNet ()
11.10.1.3 CCreateFilterNet() [2/2] CCreateFilterNet (
             int numCoefSets,
             CCreateFilter * pCreateFilter ) [protected]
11.10.2 Member Function Documentation
11.10.2.1 FindFilter() [1/2] static int FindFilter (
            array< uint64_t >^{\wedge}>^{\wedge} coef,
             array< CCreateFilterNet^>^ param,
            CFilterCoefficientsNet::s_FilterAttributesNet^ FiltAttr,
            bool DoMCSLegacyCompare ) [static]
11.10.2.2 FindFilter() [2/2] static int FindFilter (
            array< CFilterCoefficientsNet^>^ coef,
             array< CCreateFilterNet^>^ param ) [static]
11.10.2.3 GetBiQuad() CFilterCoefficientsNet ^ GetBiQuad (
            int index )
11.10.2.4 GetBiQuads() array<CFilterCoefficientsNet^> ^ GetBiQuads ()
11.10.3 Property Documentation
11.10.3.1 CutoffFrequency double CutoffFrequency [get]
11.10.3.2 NumCoefSets int NumCoefSets [get]
```


11.10.3.3 Order int Order [get]

11.10.3.4 SampleRate double SampleRate [get]

11.10.3.5 Scale double Scale [get]

11.11 CDacCalibrationFunctionNet Class Reference

Inheritance diagram for CDacCalibrationFunctionNet:



Public Member Functions

CDacCalibrationFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pDac
 — CalibrationFunctionPointerContainer)

Initializes a new instance of the CDacCalibrationFunctionNet class.

- CDacCalibrationFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CDacCalibrationFunctionNet ()
- !CDacCalibrationFunctionNet ()
- void SetDacOffset (uint16_t dacChannel, uint32_t offset)

Sets the offset of a DAC channel.

• uint32_t GetDacOffset (uint16_t dacChannel)

Gets the offset of a DAC channel.

void BurnDacOffset (uint16_t dacChannel)

Writes the offset of a DAC channel to permanent memory.

Additional Inherited Members

11.11.1 Detailed Description

11.11.2 Constructor & Destructor Documentation

```
11.11.2.1 CDacCalibrationFunctionNet() [1/2] CDacCalibrationFunctionNet (

CMcsUsbNet^ mcsusb,

CMcsUsbFunctionPointerContainer^ pDacCalibrationFunctionPointerContainer)
```

Initializes a new instance of the CDacCalibrationFunctionNet class.

```
11.11.2.2 CDacCalibrationFunctionNet() [2/2] CDacCalibrationFunctionNet (
CMcsUsbNet^ mcsusb)
```

```
11.11.2.3 ~CDacCalibrationFunctionNet() virtual ~CDacCalibrationFunctionNet () [virtual]
```

```
11.11.2.4 "!CDacCalibrationFunctionNet() !CDacCalibrationFunctionNet ( )
```

11.11.3 Member Function Documentation

```
11.11.3.1 BurnDacOffset() void BurnDacOffset (
uint16_t dacChannel)
```

Writes the offset of a DAC channel to permanent memory.

Parameters

```
dacChannel The DAC channel number.
```

Gets the offset of a DAC channel.

Parameters

dacChannel	The DAC channel number.

Returns

The offset in digits.


Sets the offset of a DAC channel.

Parameters

dacChannel	The DAC channel number.
offset	The offset in digits.

11.12 CDacqGroupChannelGenericSelectionNet Class Reference

Inheritance diagram for CDacqGroupChannelGenericSelectionNet:



Public Member Functions

CDacqGroupChannelGenericSelectionNet (CMcsUsbNet[^] mcsusb)


Additional Inherited Members

11.12.1 Constructor & Destructor Documentation

```
11.12.1.1 CDacqGroupChannelGenericSelectionNet() CDacqGroupChannelGenericSelectionNet ( CMcsUsbNet^ mcsusb )
```

11.13 CDacqGroupChannelSelectionNet Class Reference

Inheritance diagram for CDacqGroupChannelSelectionNet:



Public Member Functions

CDacqGroupChannelSelectionNet (CMcsUsbNet[^] mcsusb)

Additional Inherited Members

11.13.1 Constructor & Destructor Documentation

11.13.1.1 CDacqGroupChannelSelectionNet() CDacqGroupChannelSelectionNet (CMcsUsbNet^ mcsusb)

11.14 CDacqGroupChannelSelectionTemplateNet < DacqGroupChannelEnumTemplateNet, DacqGroupChannelEnumTemplate, CDeviceGroupChannelInfoTemplateNet > Class Template Reference

Inheritance diagram for CDacqGroupChannelSelectionTemplateNet< DacqGroupChannelEnumTemplateNet, DacqGroupChannelEnumTemplate, CDeviceGroupChannelInfoTemplateNet >:



Public Member Functions

- CDacqGroupChannelSelectionTemplateNet (CMcsUsbNet^{\(\Lambda\)} mcsusb)
- uint32_t GetNumberOfSupportedGroups ()
- uint32_t GetNumberOfSupportedGroups (uint32_t virtualDevice)
- DacqGroupChannelEnumTemplateNet GetGroupID (uint32_t Index)
- DacqGroupChannelEnumTemplateNet GetGroupID (uint32_t Index, uint32_t virtualDevice)
- uint32_t GetGroupNumberOfChannels (DacqGroupChannelEnumTemplateNet GroupID)
- DacqMeaGroupTypeEnumNet GetGroupType (DacqGroupChannelEnumTemplateNet GroupID)
- DacqMeaGroupTypeEnumNet GetGroupType (DacqGroupChannelEnumTemplateNet GroupID, uint32_

 t virtualDevice)
- void EnableChannelsInGroup (DacqGroupChannelEnumTemplateNet GroupID, List< bool >[^] Enabled←
 ChannelsBitMap)
- void EnableChannelsInGroup (DacqGroupChannelEnumTemplateNet GroupID, List< bool >[^] Enabled←
 ChannelsBitMap, uint32_t virtualDevice)
- $\bullet \ \ \, \mathsf{List} \! < \mathsf{bool} > \land \ \, \mathsf{GetEnabledChannelsInGroup} \, (\mathsf{DacqGroupChannelEnumTemplateNet} \, \mathsf{GroupID})$
- List< bool > [^] GetEnabledChannelsInGroup (DacqGroupChannelEnumTemplateNet GroupID, uint32_← t virtualDevice)
- SampleSizeNet GetGroupSampleSize (DacqGroupChannelEnumTemplateNet GroupID)
- SampleSizeNet GetGroupSampleSize (DacqGroupChannelEnumTemplateNet GroupID, uint32_t virtual ← Device)
- List< CDeviceGroupChannelInfoTemplateNet^> ^ GetDeviceGroupChannelInfos ()
- $\bullet \ \, \mathsf{List} \! < \mathsf{CDeviceGroupChannelInfoTemplateNet}^{\wedge} \! > \ \, ^{\wedge} \ \, \mathsf{GetDeviceGroupChannelInfoS} \ \, (\mathsf{uint32_t} \ \mathsf{virtualDevice}) \\$

Additional Inherited Members

11.14.1 Constructor & Destructor Documentation

```
11.14.1.1 CDacqGroupChannelSelectionTemplateNet() CDacqGroupChannelSelectionTemplateNet (
CMcsUsbNet^ mcsusb )
```

11.14.2 Member Function Documentation

```
11.14.2.1 EnableChannelsInGroup() [1/2] void EnableChannelsInGroup (

DacqGroupChannelEnumTemplateNet GroupID,

List< bool >^ EnabledChannelsBitMap )
```

```
11.14.2.2 EnableChannelsInGroup() [2/2] void EnableChannelsInGroup (

DacqGroupChannelEnumTemplateNet GroupID,

List< bool >^ EnabledChannelsBitMap,

uint32_t virtualDevice )
```


```
11.14.2.5 GetEnabledChannelsInGroup() [1/2] List<bool> ^ GetEnabledChannelsInGroup (
DacqGroupChannelEnumTemplateNet GroupID )
```

```
\textbf{11.14.2.7} \quad \textbf{GetGroupID()} \; \texttt{[1/2]} \quad \texttt{DacqGroupChannelEnumTemplateNet} \; \texttt{GetGroupID} \; \; \texttt{(}
               uint32_t Index )
11.14.2.8 GetGroupID() [2/2] DacqGroupChannelEnumTemplateNet GetGroupID (
               uint32_t Index,
               uint32_t virtualDevice )
\textbf{11.14.2.9} \quad \textbf{GetGroupNumberOfChannels() [1/2]} \quad \texttt{uint32\_t GetGroupNumberOfChannels ()}
               DacqGroupChannelEnumTemplateNet GroupID )
11.14.2.10 GetGroupNumberOfChannels() [2/2] uint32_t GetGroupNumberOfChannels (
               DacqGroupChannelEnumTemplateNet GroupID,
               uint32_t virtualDevice )
11.14.2.11 GetGroupSampleSize() [1/2] SampleSizeNet GetGroupSampleSize (
               {\tt DacqGroupChannelEnumTemplateNet} \ \ {\tt GroupID} \ )
11.14.2.12 GetGroupSampleSize() [2/2] SampleSizeNet GetGroupSampleSize (
               DacqGroupChannelEnumTemplateNet GroupID,
               uint32_t virtualDevice )
\textbf{11.14.2.13} \quad \textbf{GetGroupType()} \; \texttt{[1/2]} \quad \texttt{DacqMeaGroupTypeEnumNet} \; \texttt{GetGroupType} \; \; \texttt{(}
               DacqGroupChannelEnumTemplateNet GroupID )
11.14.2.14 GetGroupType() [2/2] DacqMeaGroupTypeEnumNet GetGroupType (
               DacqGroupChannelEnumTemplateNet GroupID,
               uint32_t virtualDevice )
11.14.2.15 GetNumberOfSupportedGroups() [1/2] uint32_t GetNumberOfSupportedGroups ( )
```

11.14.2.16 **GetNumberOfSupportedGroups()** [2/2] uint32_t GetNumberOfSupportedGroups (uint32_t virtualDevice)

11.15 CDeviceGroupChannelInfoGenericNet Class Reference

Inheritance diagram for CDeviceGroupChannelInfoGenericNet:



Public Member Functions


• CDeviceGroupChannelInfoGenericNet (int id, int channels, DacqMeaGroupTypeEnumNet type)

Additional Inherited Members

11.15.1 Constructor & Destructor Documentation

11.16 CDeviceGroupChannelInfoMEA2100_256Net Class Reference

 $Inheritance\ diagram\ for\ CDevice Group Channel Info MEA 2100_256 Net:$



Public Member Functions


 CDeviceGroupChannelInfoMEA2100_256Net (MEA2100_256DacqGroupChannelEnumNet id, int channels, DacqMeaGroupTypeEnumNet type)

Additional Inherited Members

11.16.1 Constructor & Destructor Documentation

11.17 CDeviceGroupChannelInfoNet Class Reference

Inheritance diagram for CDeviceGroupChannelInfoNet:



Public Member Functions


CDeviceGroupChannelInfoNet (DacqGroupChannelEnumNet id, int channels, DacqMeaGroupTypeEnumNet type)

Additional Inherited Members

11.17.1 Constructor & Destructor Documentation

11.18 CDeviceGroupChannelInfoSCUNet Class Reference

Inheritance diagram for CDeviceGroupChannelInfoSCUNet:



Public Member Functions

CDeviceGroupChannelInfoSCUNet (SCUDacqGroupChannelEnumNet id, int channels, DacqMeaGroup
 —
 ТуреEnumNet type)

Additional Inherited Members

11.18.1 Constructor & Destructor Documentation

11.19 CDeviceGroupChannelInfoTemplateNet < DacqGroupChannelEnumTemplateNet > Class Template Reference

Public Member Functions

 CDeviceGroupChannelInfoTemplateNet (DacqGroupChannelEnumTemplateNet id, int channels, DacqMea← GroupTypeEnumNet type)

Public Attributes

- DacqGroupChannelEnumTemplateNet GroupID
- int NumberOfChannels
- DacqMeaGroupTypeEnumNet GroupType

11.19.1 Constructor & Destructor Documentation

11.19.2 Member Data Documentation


11.19.2.1 GroupID DacqGroupChannelEnumTemplateNet GroupID

11.19.2.2 GroupType DacqMeaGroupTypeEnumNet GroupType

11.19.2.3 NumberOfChannels int NumberOfChannels

11.20 CDeviceGroupChannelInfoW2100Net Class Reference

Inheritance diagram for CDeviceGroupChannelInfoW2100Net:



Public Member Functions

 CDeviceGroupChannelInfoW2100Net (W2100DacqGroupChannelEnumNet id, int channels, DacqMea← GroupTypeEnumNet type)


Additional Inherited Members

11.20.1 Constructor & Destructor Documentation

11.21 CDigOutStimulatorFunctionNet Class Reference

CDigOutStimulatorFunctionNet is the class of the DigOut stimulator function class.

Inheritance diagram for CDigOutStimulatorFunctionNet:



Public Member Functions

CDigOutStimulatorFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pDigOut
 — StimulatorFunctionPointerContainer)

Initializes a new instance of the CDigOutStimulatorFunctionNet class.

- CDigOutStimulatorFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CDigOutStimulatorFunctionNet ()
- !CDigOutStimulatorFunctionNet ()
- void ClearChannel (int32_t NrChannel)

clear stimulation pattern

CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData ^ PrepareChannelData (array< int32_t >^ Amplitude, array< uint64_t >^ Duration)

prepares the channel data for the device and unrolles the data for the GUI

void SendChannelData (int32_t NrChannel, CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData[^] device_data_and_unrolled)

send or append stimulation pattern

• int32_t GetNumberOfChannels ()

get the number of channels available on the device

void SetGlobalRepeat (int32_t NrChannel, bool value)

set repeat whole stimulation pattern

bool GetGlobalRepeat (int32_t NrChannel)

get repeat whole stimulation pattern

- void SetStartTriggerSlope (int32_t NrChannel, DigitalStimulatorTriggerSlopeEnumNet Condition) sets start condition of digital out stimulator
- DigitalStimulatorTriggerSlopeEnumNet GetStartTriggerSlope (int32_t NrChannel)

queries start condition of digital out stimulator

- void SetStopTriggerSlope (int32_t NrChannel, DigitalStimulatorTriggerSlopeEnumNet Condition)
 sets stop condition of digital out stimulator
- DigitalStimulatorTriggerSlopeEnumNet GetStopTriggerSlope (int32_t NrChannel)
 queries stop condition of digital out stimulator

Additional Inherited Members

11.21.1 Detailed Description

CDigOutStimulatorFunctionNet is the class of the DigOut stimulator function class.

11.21.2 Constructor & Destructor Documentation

Initializes a new instance of the CDigOutStimulatorFunctionNet class.

```
11.21.2.2 CDigOutStimulatorFunctionNet() [2/2] CDigOutStimulatorFunctionNet (
             CMcsUsbNet^ mcsusb )
11.21.2.3 ~CDigOutStimulatorFunctionNet() virtual ~CDigOutStimulatorFunctionNet () [virtual]
11.21.2.4 "!CDigOutStimulatorFunctionNet() !CDigOutStimulatorFunctionNet ( )
11.21.3 Member Function Documentation
11.21.3.1 ClearChannel() void ClearChannel (
             int32_t NrChannel )
clear stimulation pattern
Parameters
 NrChannel
             the channel to clear
11.21.3.2 GetGlobalRepeat() bool GetGlobalRepeat (
             int32_t NrChannel )
get repeat whole stimulation pattern
Parameters
 NrChannel
             channel number
Returns
    current value
11.21.3.3 GetNumberOfChannels() int32_t GetNumberOfChannels ( )
get the number of channels available on the device
Returns
```

the number of channels

11.21.3.4 GetStartTriggerSlope() DigitalStimulatorTriggerSlopeEnumNet GetStartTriggerSlope (int32_t NrChannel)

queries start condition of digital out stimulator

Parameters

```
NrChannel channel number
```

Returns

start condition (rising or falling edge)

```
11.21.3.5 GetStopTriggerSlope() DigitalStimulatorTriggerSlopeEnumNet GetStopTriggerSlope ( int32_t NrChannel )
```

queries stop condition of digital out stimulator

Parameters

NrChannel channel number

Returns

stop condition (rising or falling edge)

prepares the channel data for the device and unrolles the data for the GUI

 $array < uint64_t >^{\land} Duration)$

Parameters

Amplitude	array of amplitudes
Duration	array of durations

Returns

send or append stimulation pattern

Parameters

NrChannel	the channel to send data to
device_data_and_unrolled	internal, use value obtained from PrepareChannelData

set repeat whole stimulation pattern

Parameters

NrChannel	channel number
value	new value

sets start condition of digital out stimulator

Parameters

NrChannel	channel number
Condition	start condition (rising or falling edge)


sets stop condition of digital out stimulator

NrChannel	channel number
Condition	stop condition (rising or falling edge)

11.22 CEncapsulatorDeviceNet Class Reference

CEncapsulatorDeviceNet is the to control the MCS HiClamp device

Inheritance diagram for CEncapsulatorDeviceNet:



Public Member Functions

- CEncapsulatorDeviceNet (void)
- CRoboFluidDeviceNet ^ GetRoboFluidDevice ()

Additional Inherited Members

11.22.1 Detailed Description

CEncapsulatorDeviceNet is the to control the MCS HiClamp device

11.22.2 Constructor & Destructor Documentation


```
11.22.2.1 CEncapsulatorDeviceNet() CEncapsulatorDeviceNet (
void )
```

11.22.3 Member Function Documentation

11.22.3.1 GetRoboFluidDevice() CRoboFluidDeviceNet ^ GetRoboFluidDevice ()

11.23 CExternDTesterDeviceNet Class Reference

CExternDTesterDeviceNet is the class to access the ExternD Tester (Handheld Device Tester D)
Inheritance diagram for CExternDTesterDeviceNet:



Public Member Functions

CExternDTesterDeviceNet ()

Initializes a new instance of the CExternDTesterDeviceNet class.

- virtual ~CExternDTesterDeviceNet ()
- !CExternDTesterDeviceNet ()
- array< uint8_t > ^ Read (int configString_Length)

Reads the config string from the device.

String ^ Read2 ()

Reads the config string from the device.

void Write (array< uint8_t >^ configString)

Reads the config string from the device.

void Write2 (String[^] configString)

Reads the config string from the device.

Additional Inherited Members

11.23.1 Detailed Description

CExternDTesterDeviceNet is the class to access the ExternD Tester (Handheld Device Tester D)

11.23.2 Constructor & Destructor Documentation

```
11.23.2.1 CExternDTesterDeviceNet() CExternDTesterDeviceNet ()
```

Initializes a new instance of the CExternDTesterDeviceNet class.

```
11.23.2.2 ~CExternDTesterDeviceNet() virtual ~CExternDTesterDeviceNet () [virtual]
```

```
11.23.2.3 "!CExternDTesterDeviceNet() !CExternDTesterDeviceNet ()
```

11.23.3 Member Function Documentation

Reads the config string from the device.

configString_Length	The maximal length of configString.
---------------------	-------------------------------------

Returns

The config string.

11.23.3.2 Read2() String $^{\land}$ Read2 ()

Reads the config string from the device.

Returns

The config string.

Reads the config string from the device.

Parameters

configString	The config string.

11.23.3.4 Write2() void Write2 ($String^{\wedge} configString$)

Reads the config string from the device.

Parameters

configString	The config string.

11.24 CFilterCoefficientsNet Class Reference

Classes

• struct s_FilterAttributesNet

Public Member Functions

- CFilterCoefficientsNet ()
- CFilterCoefficientsNet (double b0, double b1, double b2, double a1, double a2)
- CFilterCoefficientsNet (double b0, double b1, double a1)
- CFilterCoefficientsNet (array< double >[^] b, array< double >[^] a)
- ∼CFilterCoefficientsNet ()
- bool IsEqual (array< uint64_t >^ coefficients, s_FilterAttributesNet^ FiltAttr)
- bool IsEqual (array< uint64_t >^ coefficients, s_FilterAttributesNet^ FiltAttr, bool DoMCSLegacyCompare)
- uint64_t GetUintB (int index, s_FilterAttributesNet[^] FiltAttr)
- uint64 t GetUintA (int index, s FilterAttributesNet[^] FiltAttr)

Properties

```
• array< double >^ A [get]
```

• array< double >^ B [get]

11.24.1 Constructor & Destructor Documentation

```
11.24.1.1 CFilterCoefficientsNet() [1/4] CFilterCoefficientsNet ( )
```

```
11.24.1.2 CFilterCoefficientsNet() [2/4] CFilterCoefficientsNet (
```

```
double b0,
double b1,
double b2,
double a1,
double a2)
```

```
11.24.1.3 CFilterCoefficientsNet() [3/4] CFilterCoefficientsNet (
```

```
double b0,
double b1,
double a1 )
```

11.24.1.4 CFilterCoefficientsNet() [4/4] CFilterCoefficientsNet (

```
array< double >^{\wedge} b, array< double >^{\wedge} a)
```

11.24.1.5 \sim CFilterCoefficientsNet() \sim CFilterCoefficientsNet ()

11.24.2 Member Function Documentation

```
11.24.2.1 GetUintA() uint64_t GetUintA (
    int index,
    s_FilterAttributesNet^ FiltAttr )

11.24.2.2 GetUintB() uint64_t GetUintB (
    int index,
    s_FilterAttributesNet^ FiltAttr )

11.24.2.3 IsEqual() [1/2] bool IsEqual (
    array< uint64_t >^ coefficients,
    s_FilterAttributesNet^ FiltAttr )

11.24.2.4 IsEqual() [2/2] bool IsEqual (
    array< uint64_t >^ coefficients,
    s_FilterAttributesNet^ FiltAttr,
    bool DoMCSLegacyCompare )


11.24.3 Property Documentation
```

11.25 CFilterConfigurationNet Class Reference

Inheritance diagram for CFilterConfigurationNet:

11.24.3.1 A array< double>^ A [get]

11.24.3.2 B array< double>^ B [get]



Public Member Functions

- CFilterConfigurationNet (CMcsUsbNet[^] mcsusb)
- void SetFilterParameter (DacqGroupChannelEnumNet GroupID, uint32_t FilterNumber, CFilterCoefficientsNet[^] Coefficients, CFilterPropertyNet[^] FilterProp)
- void SetFilterParameter (DacqGroupChannelEnumNet GroupID, uint32 t FilterNumber, CFilterCoefficientsNet[^] CoefficientsSet1, CFilterCoefficientsNet[^] CoefficientsSet2, CFilterPropertyNet[^] FilterProp)
- void SetFilterParameterPermanent (DacqGroupChannelEnumNet GroupID, uint32 t FilterNumber)
- void EraseFilterParameterPermanent (DacqGroupChannelEnumNet GroupID, uint32_t FilterNumber)
- void SetHighpassFilterEnable (bool enable)
- bool GetHighpassFilterEnable ()
- void ResetHighpassFilter ()
- EnumNet index)
- CFilterCoefficientsNet::s_FilterAttributesNet ^ GetFilterAttributes (DacqGroupChannelEnumNet GroupID, uint32 t FilterNumber)

Additional Inherited Members

11.25.1 Constructor & Destructor Documentation

```
11.25.1.1 CFilterConfigurationNet() CFilterConfigurationNet (
             CMcsUsbNet^ mcsusb )
```

11.25.2 Member Function Documentation

```
11.25.2.1 EraseFilterParameterPermanent() void EraseFilterParameterPermanent (
             DacqGroupChannelEnumNet GroupID,
             uint32_t FilterNumber )
```

```
11.25.2.2 GetFilterAttributes() [1/2] CFilterCoefficientsNet::s_FilterAttributesNet ^ GetFilter↔
Attributes (
             DacqGroupChannelEnumNet GroupID,
```

```
uint32_t FilterNumber )
```

```
11.25.2.3 GetFilterAttributes() [2/2] uint32_t GetFilterAttributes (
             DacqGroupChannelEnumNet GroupID,
```


```
uint32_t FilterNumber,
FilterAttributeEnumNet index )
```

```
11.25.2.4 GetHighpassFilterEnable() bool GetHighpassFilterEnable ( )
11.25.2.5 ResetHighpassFilter() void ResetHighpassFilter ( )
11.25.2.6 SetFilterParameter() [1/2] void SetFilterParameter (
             DacqGroupChannelEnumNet GroupID,
             uint32_t FilterNumber,
             CFilterCoefficientsNet<sup>∧</sup> Coefficients,
             CFilterPropertyNet^ FilterProp )
11.25.2.7 SetFilterParameter() [2/2] void SetFilterParameter (
             DacqGroupChannelEnumNet GroupID,
             uint32_t FilterNumber,
             CFilterCoefficientsNet^ CoefficientsSet1,
             CFilterCoefficientsNet<sup>∧</sup> CoefficientsSet2,
             CFilterPropertyNet<sup>∧</sup> FilterProp )
11.25.2.8 SetFilterParameterPermanent() void SetFilterParameterPermanent (
             DacqGroupChannelEnumNet GroupID,
             uint32_t FilterNumber )
11.25.2.9 SetHighpassFilterEnable() void SetHighpassFilterEnable (
```

11.26 CFilterConfigurationRegisterNet Class Reference

Inheritance diagram for CFilterConfigurationRegisterNet:

bool enable)



Public Member Functions

- CFilterConfigurationRegisterNet (CMcsUsbNet[∧] mcsusb)
- void SetFilterParameter (uint32_t FilterCoefRegBase, CFilterCoefficientsNet[^] Coefficients, uint32_t Filter←
 InfoRegBase, CFilterPropertyNet[^] FilterProp)
- void SetFilterParameter (uint32_t FilterCoefSet1RegBase, CFilterCoefficientsNet[^] CoefficientsSet1, uint32_t FilterCoefSet2RegBase, CFilterCoefficientsNet[^] CoefficientsSet2, uint32_t FilterInfoRegBase, CFilterPropertyNet[^] FilterProp)
- void SetFilterParameterPermanent (uint32_t FilterCoefRegBase, uint32_t FilterCoefDmaReg, uint32_←
 t FilterInfoRegBase, uint32_t FilterInfoDmaReg, uint32_t EEPROMBize)
- void SetFilterParameterPermanent (uint32_t FilterCoefSet1RegBase, uint32_t FilterCoefSet1DmaReg, uint32_t FilterCoefSet2RegBase, uint32_t FilterCoefSet2DmaReg, uint32_t FilterInfoRegBase, uint32_t FilterInfoDmaReg, uint32_t EEPROMBase, uint32_t EEPROMSize)
- void EraseFilterParameterPermanent (uint32_t FilterCoefDmaReg, uint32_t FilterInfoDmaReg, uint32_t E←
 EPROMBase, uint32_t EEPROMSize)
- void EraseFilterParameterPermanent (uint32_t FilterCoefSet1DmaReg, uint32_t FilterCoefSet2DmaReg, uint32_t FilterInfoDmaReg, uint32_t EEPROMBase, uint32_t EEPROMSize)

Additional Inherited Members

11.26.1 Constructor & Destructor Documentation

11.26.2 Member Function Documentation

uint32_t EEPROMSize)

```
11.26.2.3 SetFilterParameter() [1/2] void SetFilterParameter (
             uint32_t FilterCoefRegBase,
             CFilterCoefficientsNet^ Coefficients,
             uint32_t FilterInfoRegBase,
             CFilterPropertyNet^ FilterProp )
11.26.2.4 SetFilterParameter() [2/2] void SetFilterParameter (
             uint32_t FilterCoefSet1RegBase,
             CFilterCoefficientsNet^ CoefficientsSet1,
             uint32_t FilterCoefSet2RegBase,
             CFilterCoefficientsNet<sup>∧</sup> CoefficientsSet2,
             uint32_t FilterInfoRegBase,
             CFilterPropertyNet^ FilterProp )
11.26.2.5 SetFilterParameterPermanent() [1/2] void SetFilterParameterPermanent (
             uint32_t FilterCoefRegBase,
             uint32_t FilterCoefDmaReg,
             uint32_t FilterInfoRegBase,
             uint32_t FilterInfoDmaReg,
             uint32_t EEPROMBase,
             uint32_t EEPROMSize )
11.26.2.6 SetFilterParameterPermanent() [2/2] void SetFilterParameterPermanent (
             uint32_t FilterCoefSet1RegBase,
             uint32_t FilterCoefSet1DmaReg,
             uint32_t FilterCoefSet2RegBase,
             uint32_t FilterCoefSet2DmaReg,
             uint32_t FilterInfoRegBase,
             uint32_t FilterInfoDmaReg,
             uint32_t EEPROMBase,
             uint32_t EEPROMSize )
```

11.27 CFilterPropertyNet Class Reference

Public Member Functions

- CFilterPropertyNet (uint32_t CornerFrequenzymHz, uint32_t Order, FilterBandEnumNet FilterBand, Filter FamilyEnumNet FilterTypeEnumNet FilterType, bool Active)
- ∼CFilterPropertyNet ()
- virtual System::String ^ ToString () override

Properties

```
    uint32_t CornerFrequencymHz [get]
    double CornerFrequency [get]
    uint32_t Order [get]
    FilterBandEnumNet FilterBand [get]
    FilterFamilyEnumNet FilterFamily [get]
```

FilterTypeEnumNet FilterType [get]

• bool FilterActive [get]

11.27.1 Constructor & Destructor Documentation

```
11.27.1.1 CFilterPropertyNet() CFilterPropertyNet (
            uint32_t CornerFrequenzymHz,
            uint32_t Order,
            FilterBandEnumNet FilterBand,
            FilterFamilyEnumNet FilterFamily,
            FilterTypeEnumNet FilterType,
             bool Active )
11.27.1.2 ~CFilterPropertyNet() ~CFilterPropertyNet ()
11.27.2 Member Function Documentation
11.27.2.1 ToString() virtual System::String ^ ToString ( ) [override], [virtual]
11.27.3 Property Documentation
11.27.3.1 CornerFrequency double CornerFrequency [get]
11.27.3.2 CornerFrequencymHz uint32_t CornerFrequencymHz [get]
11.27.3.3 FilterActive bool FilterActive [get]
11.27.3.4 FilterBand FilterBandEnumNet FilterBand [get]
11.27.3.5 FilterFamily FilterFamilyEnumNet FilterFamily [get]
```


11.27.3.6 FilterType FilterTypeEnumNet FilterType [get]

11.27.3.7 Order uint32_t Order [get]

11.28 CFluidControlDeviceNet Class Reference

CFluidControlDeviceNet is the class to control MCS FluidControl (FCB and FCX) device.

Inheritance diagram for CFluidControlDeviceNet:



Public Member Functions

CFluidControlDeviceNet ()

Initialize a new instance of the CFluidControlDeviceNet class.

∼CFluidControlDeviceNet ()

Default destructor.

void SetValve (unsigned int value)

Open or Close valves.

void SetSingleValve (unsigned short valve, unsigned short onoff)

Opens or Closes a valve.

• void SetDigout (unsigned int value)

Define the pattern on the Digital Output.

void SetPWM (unsigned int channel, unsigned int value)

Sets the duty cycle of the PWM output.

• void CalibrateThermocouple (unsigned int channel)

Calibrates the ADC which is used for the Thermocouple. For the calibration, Short circuit the Thermocouple and use this function to correct a possible offset of the ADC which measures the thermocouple.

• void SetThermocoupleNanovoltPerKelvin (unsigned int channel, unsigned int value)

Sets the proportinal constant for the Thermocouple.

• unsigned int GetValve ()

Gets the state of the valves.

• unsigned short GetSingleValve (unsigned short valve)

Gets the state of a valve.

• unsigned int GetDigout ()

Gets the state of the digital output.

• unsigned int GetPWM (unsigned int channel)

Gets the state of the PWM output.

unsigned int GetAdc (unsigned int channel)

Reads an ADC Value.

• unsigned int GetDigin ()

Reads the digital input.

• int GetThermocoupleTemperature (unsigned int channel)

Reads the temperature from Thermocouple. The functions gives the temperature difference between both Thermocouple junctions. To get the absolute temperature, add the reference temperature.

• int GetReferenceTemperature (unsigned int channel)

Reads the reference temperature for the Thermocouple.

unsigned int GetThermocoupleCalibration (unsigned int channel)

Gets the calibration constant for the Thermocouple ADC.

unsigned int GetThermocoupleNanovoltPerKelvin (unsigned int channel)

Reads the proportional constant for the Thermocouple.

Properties

CMcsBus_VoltageModeNet[^] McsBus_VoltageMode [get]

Additional Inherited Members

11.28.1 Detailed Description

CFluidControlDeviceNet is the class to control MCS FluidControl (FCB and FCX) device.

11.28.2 Constructor & Destructor Documentation

```
11.28.2.1 CFluidControlDeviceNet() CFluidControlDeviceNet ()
```

Initialize a new instance of the CFluidControlDeviceNet class.

```
11.28.2.2 ~CFluidControlDeviceNet() ~CFluidControlDeviceNet ()
```

Default destructor.

11.28.3 Member Function Documentation

```
11.28.3.1 CalibrateThermocouple() void CalibrateThermocouple ( unsigned int channel )
```

Calibrates the ADC which is used for the Thermocouple. For the calibration, Short circuit the Thermocouple and use this function to correct a possible offset of the ADC which measures the thermocouple.

channel	Thermocouple channel number.
---------	------------------------------

11.28.3.2 GetAdc() unsigned int GetAdc (unsigned int channel)

Reads an ADC Value.

Parameters

channel The ADC channel number to query	/.
---	----

Returns

The current ADC value.

11.28.3.3 GetDigin() unsigned int GetDigin ()

Reads the digital input.

Returns

The bit pattern of the state of the digital inputs.

11.28.3.4 GetDigout() unsigned int GetDigout ()

Gets the state of the digital output.

Returns

The current state of the digital outputs as a bit pattern.

```
11.28.3.5 GetPWM() unsigned int GetPWM ( unsigned int channel)
```

Gets the state of the PWM output.

Returns

The current state of the PWM outputs duty cycle in permille.

```
11.28.3.6 GetReferenceTemperature() int GetReferenceTemperature ( unsigned int channel )
```

Reads the reference temperature for the Thermocouple.

channel	Thermocouple channel number.
---------	------------------------------

Returns

The temperature from the Thermocouple in 1/100 °C.

11.28.3.7 **GetSingleValve()** unsigned short GetSingleValve (unsigned short *valve*)

Gets the state of a valve.

Parameters

valve	number of valve
vaive	number of valve

Returns

state of the valve

11.28.3.8 GetThermocoupleCalibration() unsigned int GetThermocoupleCalibration (unsigned int *channel*)

Gets the calibration constant for the Thermocouple ADC.

Parameters

channel	Thermocouple channel number.

Returns

The calibration constant for the Thermocouple ADC.

11.28.3.9 **GetThermocoupleNanovoltPerKelvin()** unsigned int GetThermocoupleNanovoltPerKelvin (unsigned int *channel*)

Reads the proportional constant for the Thermocouple.

	channel	Thermocouple channel number.
--	---------	------------------------------

Returns

The proportional constant in Nanovolt per Kelvin.

11.28.3.10 **GetThermocoupleTemperature()** int GetThermocoupleTemperature (unsigned int *channel*)

Reads the temperature from Thermocouple. The functions gives the temperature difference between both Thermocouple junctions. To get the absolute temperature, add the reference temperature.

Parameters

channel	Thermocouple channel number.
---------	------------------------------

Returns

The temperature difference between both Thermocouple junctions in 1/100 ℃.

11.28.3.11 GetValve() unsigned int GetValve ()

Gets the state of the valves.

Returns

The current state of the valves as a bit pattern.

```
11.28.3.12 SetDigout() void SetDigout (
unsigned int value)
```

Define the pattern on the Digital Output.

Parameters

```
value bit pattern on the digital output.
```

```
11.28.3.13 SetPWM() void SetPWM (
unsigned int channel,
unsigned int value)
```

Sets the duty cycle of the PWM output.

channel	PWM channel number.
value	duty cycle of the PWM output in permille.

```
11.28.3.14 SetSingleValve() void SetSingleValve (
unsigned short valve,
unsigned short onoff)
```

Opens or Closes a valve.

Parameters

valve number of valve to be change

Parameters

onoff open or close the valve.

```
11.28.3.15 SetThermocoupleNanovoltPerKelvin() void SetThermocoupleNanovoltPerKelvin (
unsigned int channel,
unsigned int value)
```

Sets the proportinal constant for the Thermocouple.

Parameters

channel	Thermocouple channel number.
value	proportinal constant for the Thermocouple in Nanovolt per Kelvin.

```
11.28.3.16 SetValve() void SetValve (
unsigned int value)
```

Open or Close valves.

value	bit pattern of valves which should be open.


11.28.4 Property Documentation

11.28.4.1 McsBus_VoltageMode CMcsBus_VoltageModeNet^ McsBus_VoltageMode [get]

11.29 CFYIDeviceNet Class Reference

CFYIDeviceNet is the class to control the MCS FYI device

Inheritance diagram for CFYIDeviceNet:



Public Member Functions

• CFYIDeviceNet (void)

Properties

- CRobo_FYITemp_FunctionNet^ FYITemp [get]
- CRobo_FYIProgram_FunctionNet^ FYIProgram [get]
- CMcsBus_SensorNet^ Sensor [get]

Additional Inherited Members

11.29.1 Detailed Description

CFYIDeviceNet is the class to control the MCS FYI device

11.29.2 Constructor & Destructor Documentation

```
11.29.2.1 CFYIDeviceNet() CFYIDeviceNet (
```

11.29.3 Property Documentation

```
11.29.3.1 FYIProgram CRobo_FYIProgram_FunctionNet^ FYIProgram [get]
```


```
11.29.3.2 FYITemp CRobo_FYITemp_FunctionNet^ FYITemp [get]
```

```
11.29.3.3 Sensor CMcsBus_SensorNet^ Sensor [get]
```

11.30 CGenericDevelopDeviceNet Class Reference

CGenericDevelopDeviceNet is the class to use during development of a new device.

Inheritance diagram for CGenericDevelopDeviceNet:



Public Member Functions

- CGenericDevelopDeviceNet (void)
 - Initialize a new instance of the CGenericDevelopDeviceNet class.
- ~CGenericDevelopDeviceNet (void)
- void SetValue (uint16_t value, uint16_t index)

Sets .

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

• template<typename C > void SetBuffer (uint16_t value, uint16_t index, array< C > buffer)

void SetUByteBuffer (uint16_t value, uint16_t index, array< unsigned char >^ buffer)
 Sends an array of type unsigned char to the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

Parameters

buffer	The buffer to send.

• void SetByteBuffer (uint16_t value, uint16_t Index, array< char $>^{\wedge}$ buffer) Sends an array of type char to the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

Index The index of the request.

Parameters

ſ	buffer	The buffer to send.

• void SetUShortBuffer (uint16_t value, uint16_t index, array< unsigned short >^ buffer)

Sends an array of type unsigned short to the device.

Parameters

value The value of the request.	
---------------------------------	--

index The index of the request.	_
---------------------------------	---

void SetShortBuffer (uint16_t value, uint16_t index, array< short >^ buffer)
 Sends an array of type short to the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

Parameters

ouffer The buffer to send.

• void SetUIntBuffer (uint16_t value, uint16_t index, array< unsigned int >^ buffer)

Sends an array of unsigned int to the device.

Parameters

value	The value of the request.

Parameters

	index	The index of the request.
- 1		

Parameters

buffer	The buffer to send.
Dunci	The build to send.

• void SetIntBuffer (uint16_t value, uint16_t index, array< int $>^{\land}$ buffer) Sends an array of type int to the device.

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

Parameters

- template < typename C > array < C > $^{\land}$ GetBuffer (uint16_t value, uint16_t index, int size)
- array< unsigned char > ^ GetUByteBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type unsigned char from the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index The index of the reques	st.
-------------------------------	-----


Parameters

size The size of the	array.
----------------------	--------

Returns

The array of data from the device.

array< char > ^ GetByteBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type char from the device.



index	The index of the request.
-------	---------------------------

Parameters

Returns

The array of data from the device.

array< unsigned short > ^ GetUShortBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type unsigned short from the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index The index of the requ	uest.
-----------------------------	-------

Parameters

size The	size of the array.
----------	--------------------

Returns

The array of data from the device.

array< short > ^ GetShortBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type short from the device.

Parameters

value The value of the request.	_
---------------------------------	---

st.
st.

size	The size of the array.
------	------------------------

Returns

The array of data from the device.

array< unsigned int > ^ GetUIntBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type unsigned int from the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

Parameters

Returns

The array of data from the device.

array< int > ^ GetIntBuffer (uint16_t value, uint16_t index, int size)
 Gets an array of type int from the device.

Parameters

value The value	of the request.
-----------------	-----------------

Parameters

index	The index of the request.
-------	---------------------------

Returns

The array of data from the device.

• IntPtr OpenPipe (uint8_t endpointAddress)

Open a Pipe to an USB Endpoint.

Parameters

Returns

A handle to the endpoint.

• void ClosePipe (IntPtr pipeHandle)

Close a Pipe to an USB Endpoint.

Parameters

pipeHandle	The endpoint handle.
------------	----------------------

• void ResetPipe (IntPtr pipeHandle)

Reset a Pipe to an USB Endpoint.

Parameters

pipeHandle	The endpoint handle.
------------	----------------------

• template<typename C > array< C > $^{\wedge}$ ReadPipe (IntPtr pipeHandle, uint32_t size)

Read data from an USB Endpoint.

Parameters

pipeHandle	The endpoint handle.
------------	----------------------

Parameters

size	Number of elements to read.
------	-----------------------------

Returns

An array of data read.

• template<typename C > void WritePipe (IntPtr pipeHandle, array< C $>^{\wedge}$ buffer)

Write data to an USB Endpoint.

pipeHandle	The endpoint handle.
pipenariole	тпе епаропп папате.

buffer	An array of data to write.
--------	----------------------------

Additional Inherited Members

11.30.1 Detailed Description

CGenericDevelopDeviceNet is the class to use during development of a new device.

11.30.2 Constructor & Destructor Documentation

```
11.30.2.1 CGenericDevelopDeviceNet() CGenericDevelopDeviceNet (
void )
```

Initialize a new instance of the CGenericDevelopDeviceNet class.

```
11.30.2.2 ~CGenericDevelopDeviceNet() ~CGenericDevelopDeviceNet ( void )
```

11.30.3 Member Function Documentation

Close a Pipe to an USB Endpoint.

pipeHandle	The endpoint handle.

Gets an array of type char from the device.

Parameters

value The value of the request.

Parameters

	index	The index of the request.
--	-------	---------------------------

Parameters

```
size The size of the array.
```

Returns

The array of data from the device.

Gets an array of type int from the device.

value	The value of the request.
-------	---------------------------

index	The index of the request.
-------	---------------------------

Parameters

size	The size of the array.
------	------------------------

Returns

The array of data from the device.

Gets an array of type short from the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index The index of the reques	st.
-------------------------------	-----

size	The size of the array.
0,20	The oize of the array.

Returns

The array of data from the device.

Gets an array of type unsigned char from the device.

Parameters

value The value of the request.

Parameters

	index	The index of the request.
--	-------	---------------------------

Parameters

```
size The size of the array.
```

Returns

The array of data from the device.

Gets an array of type unsigned int from the device.

value	The value of the request.
-------	---------------------------

index	The index of the request.
-------	---------------------------

Parameters

Returns

The array of data from the device.

Gets an array of type unsigned short from the device.

Parameters

	value	The value of the request.
--	-------	---------------------------

Parameters

index The index of the reques	t.
-------------------------------	----

size	The size of the array.

Returns

The array of data from the device.

Open a Pipe to an USB Endpoint.

Parameters

Returns

A handle to the endpoint.

Read data from an USB Endpoint.

Parameters

pipeHandle	The endpoint handle.
piporiariaro	i i i o o i aponiti i anaio.

Parameters

size	Number of elements to read.	
0,20	i tambor or oronnome to road.	

Returns

An array of data read.

Reset a Pipe to an USB Endpoint.

Parameters

pipeHandle The endpoint handle.

Sends an array of type char to the device.

Parameters

	value	The value of the request.
--	-------	---------------------------

Parameters

	Index	The index of the request.
--	-------	---------------------------

Parameters

buffer The buffer to send.

Sends an array of type int to the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index The index of the req	uest.
----------------------------	-------

Parameters

buffer	The buffer to send.
Dunci	i i i c banci to scria.

Sends an array of type short to the device.

Parameters

lue of the request.	value
---------------------	-------

	index	The index of the request.
--	-------	---------------------------

Parameters

Sends an array of type unsigned char to the device.

Parameters

value The value of the request	t.
--------------------------------	----

Parameters

index	The index of the request.
-------	---------------------------

Parameters

```
buffer The buffer to send.
```

Sends an array of unsigned int to the device.

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

Parameters

buffer The buffer to send.

Sends an array of type unsigned short to the device.

Parameters

value	The value of the request.
-------	---------------------------

Parameters

index	The index of the request.
-------	---------------------------

buffer The buffer to	send.
----------------------	-------

Sets.

Parameters

est.

Parameters

index The index of the re

```
11.30.3.20 WritePipe() void WritePipe (

IntPtr pipeHandle,

array < C >^ buffer)
```

Write data to an USB Endpoint.

Parameters


pipeHandle	The endpoint handle.
------------	----------------------

Parameters

11.31 CGilsonDeviceNet Class Reference

CGilsonDeviceNet is the class to control a Gilson device.

Inheritance diagram for CGilsonDeviceNet:



Public Member Functions

- CGilsonDeviceNet (void)
 - Initialize a new instance of the CGilsonDeviceNet class.
- ∼CGilsonDeviceNet (void)
- void ConnectSlave (byte ID)
- void SendImmediate (wchar_t command)
- String \(^\) SendImmediateGetResponse (wchar_t command)
- void SendBuffered (String[^] command)
- String ^ GetLastAnswer ()

Protected Attributes

• CGilsonDevice * m pGilsonDevice

Additional Inherited Members

11.31.1 Detailed Description

CGilsonDeviceNet is the class to control a Gilson device.

11.31.2 Constructor & Destructor Documentation

```
11.31.2.1 CGilsonDeviceNet() CGilsonDeviceNet (
void )
```

Initialize a new instance of the CGilsonDeviceNet class.

```
11.31.2.2 ~CGilsonDeviceNet() ~CGilsonDeviceNet (
```


11.31.3 Member Function Documentation

```
11.31.3.1 ConnectSlave() void ConnectSlave (
             byte ID )
11.31.3.2 GetLastAnswer() String ^ GetLastAnswer ( )
11.31.3.3 SendBuffered() void SendBuffered (
             String^{\wedge} command)
11.31.3.4 SendImmediate() void SendImmediate (
             wchar_t command )
11.31.3.5 SendImmediateGetResponse() String ^ SendImmediateGetResponse (
             wchar_t command )
11.31.4 Member Data Documentation
11.31.4.1 m_pGilsonDevice CGilsonDevice* m_pGilsonDevice [protected]
```

11.32 CGrapheneFunctionNet Class Reference

CGrapheneFunctionNet is the class to control the TEER device

Inheritance diagram for CGrapheneFunctionNet:



Public Member Functions

CGrapheneFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pGraphene←
 FunctionPointerContainer)

Initializes a new instance of the CGrapheneFunctionNet class.

- CGrapheneFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ∼CGrapheneFunctionNet ()
- !CGrapheneFunctionNet ()
- void GetVdVsDAC ([System::Runtime::InteropServices::Out]int16_t% Vd, [System::Runtime::Interop←
 Services::Out]int16_t% Vs)

Gets Vd and Vs

 void GetVdVsDAC (uint32_t Headstage, [System::Runtime::InteropServices::Out]int16_t% Vd, [System::← Runtime::InteropServices::Out]int16_t% Vs)

Gets Vd and Vs

void SetVdVsDAC (int16 t Vd, int16 t Vs)

Sets Vd and Vs

void SetVdVsDAC (uint32_t Headstage, int16_t Vd, int16_t Vs)

Sets Vd and VS

· bool GetVoltageReached ()

Gets the reached voltage

bool GetVoltageReached (uint32_t Headstage)

Gets the reached voltage

int32_t GetVoltageRange ()

Gets the voltage range

int32_t GetVoltageRange (uint32_t Headstage)

Gets the voltage range

void SetVoltageRange (int32_t range)

Sets the voltage range

void SetVoltageRange (uint32_t Headstage, int32_t range)

Sets the voltage range

• int32_t GetVoltageResolution ()

Gets the voltage resolution

int32_t GetVoltageResolution (uint32_t Headstage)

Gets the voltage resolution

void SetVoltageResolution (int32_t resolution)

Sets the voltage resolution

void SetVoltageResolution (uint32_t Headstage, int32_t resolution)

Sets the voltage resolution

void GetDACOffset ([System::Runtime::InteropServices::Out]int16_t% offset_vd, [System::Runtime::
 — InteropServices::Out]int16_t% offset_vs)

Gets the DAC offset

 void GetDACOffset (uint32_t Headstage, [System::Runtime::InteropServices::Out]int16_t% offset_vd, [System::Runtime::InteropServices::Out]int16_t% offset_vs)

Gets the DAC offset

void SetDACOffset (int16_t offset_vd, int16_t offset_vs)

Sets the DAC offset

void SetDACOffset (uint32_t Headstage, int16_t offset_vd, int16_t offset_vs)

Set the DAC offset

void GetVdVs ([System::Runtime::InteropServices::Out]int32_t% Vd, [System::Runtime::InteropServices::
 — Out]int32_t% Vs)

Gets Vd and Vs

```
102
11.32 CGrapheneFunctionNet Class Reference

    void GetVdVs (uint32_t Headstage, [System::Runtime::InteropServices::Out]int32_t% Vd, [System::

      Runtime::InteropServices::Out]int32_t% Vs)
          Gets Vd and Vs

    void SetVdVs (int32_t Vd, int32_t Vs)

          Sets Vd and Vs

    void SetVdVs (uint32_t Headstage, int32_t Vd, int32_t Vs)

         Sets Vd and Vs

    void SetVgs (int32_t Vgs)

          Sets Vgs

    void SetVgs (uint32 t Headstage, int32 t Vgs)

    void SetVds (int32_t Vds)

         Sets Vds

    void SetVds (uint32_t Headstage, int32_t Vds)

         Sets Vds
Additional Inherited Members
11.32.1 Detailed Description
CGrapheneFunctionNet is the class to control the TEER device
11.32.2 Constructor & Destructor Documentation
11.32.2.1 CGrapheneFunctionNet() [1/2] CGrapheneFunctionNet (
               CMcsUsbNet^ mcsusb,
               {\tt CMcsUsbFunctionPointerContainer}^{\wedge} \ p{\tt GrapheneFunctionPointerContainer} \ )
Initializes a new instance of the CGrapheneFunctionNet class.
```

```
11.32.2.2 CGrapheneFunctionNet() [2/2] CGrapheneFunctionNet (
            CMcsUsbNet^ mcsusb )
11.32.2.3 ~CGrapheneFunctionNet() virtual ~CGrapheneFunctionNet ( ) [virtual]
11.32.2.4 "!CGrapheneFunctionNet() !CGrapheneFunctionNet ( )
11.32.3 Member Function Documentation
```

```
11.32.3.1 GetDACOffset() [1/2] void GetDACOffset (
             [System::Runtime::InteropServices::Out] int16_t% offset_vd,
             [System::Runtime::InteropServices::Out] int16_t% offset_vs )
```

Gets the DAC offset

Parameters

offset_vd	Vd offset
offset_vs	Vs offset

$\textbf{11.32.3.2} \quad \textbf{GetDACOffset() [2/2]} \quad \texttt{void GetDACOffset ()}$

```
uint32_t Headstage,
[System::Runtime::InteropServices::Out] int16_t% offset_vd,
[System::Runtime::InteropServices::Out] int16_t% offset_vs )
```

Gets the DAC offset

Parameters

Headstage	The headstage to query.
offset_vd	Vd offset
offset_vs	Vs offset

11.32.3.3 GetVdVs() [1/2] void GetVdVs (

```
[System::Runtime::InteropServices::Out] int32_t% Vd, [System::Runtime::InteropServices::Out] int32_t% Vs)
```

Gets Vd and Vs

Parameters

Vd	Vd
Vs	Vs

11.32.3.4 GetVdVs() [2/2] void GetVdVs (

```
uint32_t Headstage,
[System::Runtime::InteropServices::Out] int32_t% Vd,
[System::Runtime::InteropServices::Out] int32_t% Vs )
```

Gets Vd and Vs

Headstage	The headstage to query.
Vd	Vd
Vs	Vs

11.32.3.5 GetVdVsDAC() [1/2] void GetVdVsDAC (

```
[System::Runtime::InteropServices::Out] int16_t% Vd, [System::Runtime::InteropServices::Out] int16_t% Vs)
```

Gets Vd and Vs

Parameters

Vd	Vd
Vs	Vs

11.32.3.6 GetVdVsDAC() [2/2] void GetVdVsDAC (

```
uint32_t Headstage,
[System::Runtime::InteropServices::Out] int16_t% Vd,
[System::Runtime::InteropServices::Out] int16_t% Vs )
```

Gets Vd and Vs

Parameters

Headstage	The headstage to query.
Vd	Vd
Vs	Vs

11.32.3.7 GetVoltageRange() [1/2] int32_t GetVoltageRange ()

Gets the voltage range

Returns

The voltage range

11.32.3.8 GetVoltageRange() [2/2] int32_t GetVoltageRange (

uint32_t Headstage)

Gets the voltage range

Headstage	The headstage to query.

Returns

The voltage range

11.32.3.9 GetVoltageReached() [1/2] bool GetVoltageReached ()

Gets the reached voltage

Returns

the reached voltage

11.32.3.10 GetVoltageReached() [2/2] bool GetVoltageReached (uint32_t Headstage)

Gets the reached voltage

Parameters

Headstage	The headstage to query.
-----------	-------------------------

Returns

The reached voltage

11.32.3.11 GetVoltageResolution() [1/2] int32_t GetVoltageResolution ()

Gets the voltage resolution

Returns

The voltage resolution

```
11.32.3.12 GetVoltageResolution() [2/2] int32_t GetVoltageResolution ( uint32_t Headstage )
```

Gets the voltage resolution

Headstage	The headstage to query.
-----------	-------------------------

Returns

The voltage resolution

Sets the DAC offset

Parameters

offset_vd	Vd
offset_vs	Vs

Set the DAC offset

Parameters

Headstage	The headstage to query.
offset_vd	Vd
offset_vs	Vs

```
11.32.3.15 SetVds() [1/2] void SetVds ( int32_t Vds )
```

Sets Vds

Parameters

```
Vds Vds
```

Sets Vds

Parameters

Headstage	The headstage to query.
Vds	Vds

Sets Vd and Vs

Parameters

Vd	Vd
Vs	Vs

Sets Vd and Vs

Parameters

Headstage	The headstage to query.
Vd	Vd
Vs	Vs

Sets Vd and Vs

Vd	Vd
Vs	Vs

Sets Vd and VS

Parameters

Headstage	The headstage to query.
Vd	Vd
Vs	Vs

```
11.32.3.21 SetVgs() [1/2] void SetVgs ( int32_t Vgs )
```

Sets Vgs

Parameters

Sets Vgs

Parameters

Headstage	The headstage to query.
Vgs	Vgs

Sets the voltage range

range	The voltage range

Sets the voltage range

Parameters

Headstage	The headstage to query.
range	The voltage range

Sets the voltage resolution

Parameters

resolution	The voltage resolution
------------	------------------------

11.32.3.26 SetVoltageResolution() [2/2] void SetVoltageResolution (uint32_t Headstage, int32_t resolution)

Sets the voltage resolution


Parameters

Headstage	The headstage to query.
resolution	The voltage resolution

11.33 CHiClampDeviceNet Class Reference

CHiClampDeviceNet is the to control the MCS HiClamp device

Inheritance diagram for CHiClampDeviceNet:



Public Member Functions

CHiClampDeviceNet (void)

Properties

• CRoboDacqNet^ RoboDacq [get]

Additional Inherited Members

11.33.1 Detailed Description

CHiClampDeviceNet is the to control the MCS HiClamp device

11.33.2 Constructor & Destructor Documentation


```
11.33.2.1 CHiClampDeviceNet() CHiClampDeviceNet ( void )
```

11.33.3 Property Documentation

```
11.33.3.1 RoboDacq CRoboDacqNet^ RoboDacq [get]
```

11.34 CHLADacqNet Class Reference

Inheritance diagram for CHLADacqNet:



Public Member Functions

CHLADacqNet (void)

Additional Inherited Members


11.34.1 Constructor & Destructor Documentation

```
11.34.1.1 CHLADacqNet() CHLADacqNet (
void )
```

11.35 CHLADeviceNet Class Reference

CHLADeviceNet is the to control the MCS HLA device

Inheritance diagram for CHLADeviceNet:



Public Member Functions

CHLADeviceNet (void)

Properties

- CHLADacqNet^ HLADacq [get]
- CSerialPortNet^ SerialPort [get]

Additional Inherited Members

11.35.1 Detailed Description

CHLADeviceNet is the to control the MCS HLA device

11.35.2 Constructor & Destructor Documentation

```
11.35.2.1 CHLADeviceNet() CHLADeviceNet (
```

11.35.3 Property Documentation

```
11.35.3.1 HLADacq CHLADacqNet^ HLADacq [get]
```

```
11.35.3.2 SerialPort CSerialPortNet^ SerialPort [get]
```

11.36 CMcsUsbDacqNet::CHWInfo Class Reference

Class to provide hardware information about the device.

Classes

· class CVoltageRangeInfoNet

Public Member Functions

- CHWInfo (CMcsUsbDacqNet[^] device)

Get the number of analog channels the device supports.

Get the number of digital channels the device supports.

· virtual bool IsDigitalChannelDedicated ()

Query if the digital channel replaces an analog channel when enabled (e.g. on MC_Card) or adds a channel link on USB devices.

- virtual uint32_t GetAvailableSampleRates ([System::Runtime::InteropServices::Out]System::Collections::
 Generic::List< int32_t >^% sampleRates)
- virtual System::Collections::Generic::List< int32_t > ^ GetAvailableVoltageRangesInMicroVolt (int milliGain)
 Gets a List of voltage ranges the device supports.
- virtual System::Collections::Generic::List< CVoltageRangeInfoNet[^]> [^] GetAvailableVoltageRangesInMicroVoltAndStringsInM (int milliGain)

Gets a List of voltage ranges the device supports.

11.36.1 Detailed Description

Class to provide hardware information about the device.

11.36.2 Constructor & Destructor Documentation

```
11.36.2.1 CHWInfo() CHWInfo (

CMcsUsbDacqNet^ device)
```

11.36.3 Member Function Documentation

Gets a List of voltage ranges the device supports.

The List is scaled by the gain factor given as argument to this function. Use "1000" as scale factor for backwards compatibility. To get a list of voltage ranges for the headstage, obtain the gain of the headstage with the Get← Gain() call and use the result in the milliGain parameter. To get a list of voltage ranges for the analog inputs of the interfaceboard, obtain the gain of the analog inputs with the GetAnalogGain() call and use the result in the milliGain parameter.

Parameters

milliGain The gain factor (in milliGain) used to scale the list of voltage ranges.

Returns

List of voltage ranges in µV.

Gets a List of voltage ranges the device supports.

The List is scaled by the gain factor given as argument to this function. Use "1000" as scale factor for backwards compatibility. Each list entry contains the voltage range as an integer and as a string. To get a list of voltage ranges for the headstage, obtain the gain of the headstage with the GetGain() call and use the result in the milliGain parameter. To get a list of voltage ranges for the analog inputs of the interfaceboard, obtain the gain of the analog inputs with the GetAnalogGain() call and use the result in the milliGain parameter.

Parameters

milliGain	The gain factor (in milliGain) used to scale the list of voltage ranges.
-----------	--

Returns

List of voltage ranges in µV.

11.36.3.4 **GetNumberOfHWADCChannels()** virtual uint32_t GetNumberOfHWADCChannels (
[System::Runtime::InteropServices::Out] int% numberOfChannels) [virtual]

Get the number of analog channels the device supports.

Parameters

numberOfChannels	Number of analog channels the device supports.
------------------	--

Returns

Error Status. 0 on success.

11.36.3.5 **GetNumberOfHWDigitalChannels()** virtual uint32_t GetNumberOfHWDigitalChannels ([System::Runtime::InteropServices::Out] int% numberOfChannels) [virtual]

Get the number of digital channels the device supports.

Parameters

numberOfChannels	Number of digital channels the device supports.

Returns

Error Status. 0 on success.

11.36.3.6 IsDigitalChannelDedicated() virtual bool IsDigitalChannelDedicated () [virtual]


Query if the digital channel replaces an analog channel when enabled (e.g. on MC_Card) or adds a channel link on USB devices.

Returns

false when the digital channel replaces an analog channel when enabled, true when the digital channels is appended to the analog channels when enabled.

11.37 CIntanMea_FunctionNet Class Reference

Inheritance diagram for CIntanMea FunctionNet:



Public Member Functions

- CIntanMea_FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] intalMea_Function
 —
 PointerContainer)
- CIntanMea_FunctionNet (CMcsUsbNet[^] mcsusb)
- int GetUpperFrequencyByIndex (unsigned short index)
- int GetLowerFrequencyByIndex (unsigned short index)
- int64_t GetDSPHighPassByIndex (unsigned short index)
- int GetIntanRegister (unsigned short chip, unsigned short registernumber)
- int GetImpedanceResult (unsigned short channel)
- void SetBandwidthByIndex (int upper_index, int lower_index)
- void SetDSPHighPassByIndex (int index)
- void AmplifierSettle ()
- · void SetIntanRegister (unsigned short register number, int value)
- · void SetDiagnosticMode (unsigned char onoff)
- void BeginImpedanceCheck (array< int >[∧] config_values)

Additional Inherited Members


11.37.1 Constructor & Destructor Documentation

```
11.37.1.2 CIntanMea_FunctionNet() [2/2] CIntanMea_FunctionNet (
             CMcsUsbNet^ mcsusb )
11.37.2 Member Function Documentation
11.37.2.1 AmplifierSettle() void AmplifierSettle ( )
11.37.2.2 BeginImpedanceCheck() void BeginImpedanceCheck (
             array < int >^ config_values )
11.37.2.3 GetDSPHighPassByIndex() int64_t GetDSPHighPassByIndex (
             unsigned short index)
11.37.2.4 GetImpedanceResult() int GetImpedanceResult (
             unsigned short channel )
11.37.2.5 GetIntanRegister() int GetIntanRegister (
             unsigned short chip,
             unsigned short registernumber )
\textbf{11.37.2.6} \quad \textbf{GetLowerFrequencyByIndex()} \quad \texttt{int GetLowerFrequencyByIndex} \quad \textbf{(}
             unsigned short index)
11.37.2.7 GetUpperFrequencyByIndex() int GetUpperFrequencyByIndex (
             unsigned short index)
11.37.2.8 SetBandwidthByIndex() void SetBandwidthByIndex (
             int upper_index,
             int lower_index )
```

11.38 CInterfaceboard2FunctionNet Class Reference

CInterfaceboard2FunctionNet is the class to control the Interfaceboard

Inheritance diagram for CInterfaceboard2FunctionNet:



Public Member Functions

 CInterfaceboard2FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pInterfaceboard2← FunctionPointerContainer)

Initializes a new instance of the CInterfaceboard2FunctionNet class.

- CInterfaceboard2FunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CInterfaceboard2FunctionNet ()
- !CInterfaceboard2FunctionNet ()
- void SetloVoltage (IoVoltageEnumNet ioVoltage)

Sets the I/O Voltage level for the IFB2 digital and AUX ports, default is 3.3V.

• IoVoltageEnumNet GetIoVoltage ()

Gets the I/O Voltage level for the IFB2 digital and AUX ports, default is 3.3V.

Additional Inherited Members

11.38.1 Detailed Description

CInterfaceboard2FunctionNet is the class to control the Interfaceboard

11.38.2 Constructor & Destructor Documentation

```
11.38.2.1 CInterfaceboard2FunctionNet() [1/2] CInterfaceboard2FunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pInterfaceboard2FunctionPointerContainer)
```

Initializes a new instance of the CInterfaceboard2FunctionNet class.

```
11.38.2.2 CInterfaceboard2FunctionNet() [2/2] CInterfaceboard2FunctionNet ( CMcsUsbNet^{\land} mcsusb )
```

```
11.38.2.3 ~CInterfaceboard2FunctionNet() virtual ~CInterfaceboard2FunctionNet ( ) [virtual]
```

```
11.38.2.4 "!CInterfaceboard2FunctionNet() !CInterfaceboard2FunctionNet ()
```

11.38.3 Member Function Documentation

```
11.38.3.1 GetloVoltage() IoVoltageEnumNet GetIoVoltage ( )
```

Gets the I/O Voltage level for the IFB2 digital and AUX ports, default is 3.3V.

Returns

Enum for the IO Voltage (3.3V or 5.0V).

Sets the I/O Voltage level for the IFB2 digital and AUX ports, default is 3.3V.


Parameters

ioVoltage Enum for the I/O Voltage (3.3V or 5.0V).

11.39 CInterfaceboardFunctionNet Class Reference

CInterfaceboardFunctionNet is the class to control the Interfaceboard

Inheritance diagram for CInterfaceboardFunctionNet:



Public Member Functions

CInterfaceboardFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pInterfaceboard
 —
 FunctionPointerContainer)

Initializes a new instance of the CInterfaceboardFunctionNet class.

- CInterfaceboardFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CInterfaceboardFunctionNet ()
- !CInterfaceboardFunctionNet ()
- void SetCardinalDacqSamplerate (uint32_t samplerate)

Sets the fundamental/cardinal data aquisition samplerate of the Interfaceboard, default is 50 kHz

uint32_t GetCardinalDacqSamplerate ()

Gets the fundamental/cardinal data aquisition samplerate of the Interfaceboard, default is 50 kHz

void SetCardinalStgOutputrate (uint32_t outputrate)

Sets the fundamental/cardinal STG output rate of the Interfaceboard, default is 50 kHz

• uint32 t GetCardinalStgOutputrate ()

Gets the fundamental/cardinal STG output rate of the Interfaceboard, default is 50 kHz

Additional Inherited Members

11.39.1 Detailed Description

CInterfaceboardFunctionNet is the class to control the Interfaceboard

11.39.2 Constructor & Destructor Documentation

```
11.39.2.1 CInterfaceboardFunctionNet() [1/2] CInterfaceboardFunctionNet (

CMcsUsbNet^ mcsusb,

CMcsUsbFunctionPointerContainer^ pInterfaceboardFunctionPointerContainer)
```

Initializes a new instance of the CInterfaceboardFunctionNet class.

```
11.39.2.2 CInterfaceboardFunctionNet() [2/2] CInterfaceboardFunctionNet (
               CMcsUsbNet^ mcsusb )
11.39.2.3 ~CInterfaceboardFunctionNet() virtual ~CInterfaceboardFunctionNet () [virtual]
11.39.2.4 "!CInterfaceboardFunctionNet() !CInterfaceboardFunctionNet ( )
11.39.3 Member Function Documentation
11.39.3.1 GetCardinalDacqSamplerate() uint32_t GetCardinalDacqSamplerate ( )
Gets the fundamental/cardinal data aquisition samplerate of the Interfaceboard, default is 50 kHz
Returns
     The samplerate in Hz.
\textbf{11.39.3.2} \quad \textbf{GetCardinalStgOutputrate()} \quad \texttt{uint32\_t GetCardinalStgOutputrate ()} \\
Gets the fundamental/cardinal STG output rate of the Interfaceboard, default is 50 kHz
Returns
     The output rate in Hz.
\textbf{11.39.3.3} \quad \textbf{SetCardinalDacqSamplerate()} \quad \texttt{void SetCardinalDacqSamplerate ()}
               uint32_t samplerate )
Sets the fundamental/cardinal data aquisition samplerate of the Interfaceboard, default is 50 kHz
Parameters
               The samplerate in Hz.
 samplerate
```

11.39.3.4 SetCardinalStgOutputrate() void SetCardinalStgOutputrate (uint32_t outputrate)

Sets the fundamental/cardinal STG output rate of the Interfaceboard, default is 50 kHz


Parameters

outputrate The output rate in Hz.

11.40 CLIH3DeviceNet Class Reference

CLIH3DeviceNet is the class to access the HEKA LIH3 device.

Inheritance diagram for CLIH3DeviceNet:



Public Member Functions

· CLIH3DeviceNet ()

Initializes a new instance of the CLIH3DeviceNet class.

- virtual ~CLIH3DeviceNet ()
- !CLIH3DeviceNet ()
- void DummyCommand (uint32_t dummyParameter)

Dummy command to show how to use the DLL.

void SetEEpromPage (uint32_t EEpromStartAddress, array< int8_t >^ EEpromData, LIH30_EPC10_Bus
 _EnumNet epc10bus)

Writes into EEprom on the EPC10 EEPROM

• array< int8_t > ^ GetEEpromPage (uint32_t EEpromStartAddress, int EEpromData_Length, LIH30_EP \leftarrow C10_Bus_EnumNet epc10bus)

Reads the requested amount of EEprom byte from the EPC10 EEPROM

void SetSampleInterval (uint32_t SampleInterval)

Sets the Sample Interval for the DACQ and Stimulation

• uint32 t GetSampleInterval ()

Gets the Sample Interval for the DACQ and Stimulation

void SetAdcOffset (LIH30_ADC_Channel_EnumNet AdcChannel, int32_t Offset)

Sets the ADC offset of the DACQ for a single channel

int32_t GetAdcOffset (LIH30_ADC_Channel_EnumNet AdcChannel)

Gets the ADC offset of the DACQ for a single channel

void SetAdcOffsetPermanent (LIH30_ADC_Channel_EnumNet AdcChannel)

Writes the ADC offset of the DACQ for a single channel to permanent EEProm memory

void ErasePermanentAdcOffset (LIH30_ADC_Channel_EnumNet AdcChannel)

Delets the ADC offset of the DACQ for a single channel in permanent EEProm memory

uint32_t ReadClipping (LIH30_EPC10_Bus_EnumNet epc10bus)

Gets the clipping information

void SetDigOutState (uint16 t DigOutState)

Writes to the LIH30 digital output

uint16_t GetDigInState ()

Reads from the LIH30 digital input

• void SendCommand (LIH30_EPC10_Bus_EnumNet epc10bus, uint16_t Command)

Send command to the EPC10

uint16 t GetDacqRunStatus ()

Gets the data acquisition running status

void SetDacUseIdleValue (uint32 t DacChannel, bool UseIdle)

Sets if the DAC Idle value is used after stimulation

bool GetDacUseIdleValue (uint32_t DacChannel)

Gets if the DAC Idle value is used after stimulation

void SetDacIdleValue (uint32_t DacChannel, int32_t IdleValue)

Sets the DAC Idle value

int32 t GetDacIdleValue (uint32 t DacChannel)

Gets the DAC Idle value

void EnableUserTrigger (bool enable)

Enables the User Trigger

bool IsUserTriggerEnabled ()

Is the User Trigger enabled

void SetDacOffset (LIH30_DAC_Channel_EnumNet DacChannel, int32_t Offset)

Sets the offset of a DAC channel.

• int32_t GetDacOffset (LIH30_DAC_Channel_EnumNet DacChannel)

Gets the offset of a DAC channel.

void SetDacOffsetPermanent (LIH30_DAC_Channel_EnumNet DacChannel)

Writes the DAC offset of the STG for a single channel to permanent EEProm memory

void ErasePermanentDacOffset (LIH30_DAC_Channel_EnumNet DacChannel)

Delets the DAC offset of the STG for a single channel in permanent EEProm memory

void SetAudioOutDacParameter (uint32_t Frequency, uint32_t Amplification)

Sets the parameter of the audio DAC output.

 void GetAudioOutDacParameter ([System::Runtime::InteropServices::Out]uint32_t% Frequency, [System::← Runtime::InteropServices::Out]uint32_t% Amplification)

Gets the parameter of the audio DAC output.

String ^ ReadUARTData ()

Reads the config string from the device.

void WriteUARTData (String[^] commandString)

Write the command string to the device.

Properties

CStimulusFunctionNet[^] StimulusFunction [get]

Additional Inherited Members

11.40.1 Detailed Description

CLIH3DeviceNet is the class to access the HEKA LIH3 device.

11.40.2 Constructor & Destructor Documentation

```
11.40.2.1 CLIH3DeviceNet() CLIH3DeviceNet ()
```

Initializes a new instance of the CLIH3DeviceNet class.

```
11.40.2.2 ~CLIH3DeviceNet() virtual ~CLIH3DeviceNet () [virtual]
```

```
11.40.2.3 "!CLIH3DeviceNet() !CLIH3DeviceNet ()
```

11.40.3 Member Function Documentation

```
11.40.3.1 DummyCommand() void DummyCommand (
uint32_t dummyParameter)
```

Dummy command to show how to use the DLL.

Parameters

dummyParameter | parameter to send to the device

```
11.40.3.2 EnableUserTrigger() void EnableUserTrigger ( bool enable )
```

Enables the User Trigger

Parameters

enable Enable

Delets the ADC offset of the DACQ for a single channel in permanent EEProm memory

Parameters

AdcChannel	The ADC channel
riaconanici	THE ADD CHAINCE

Delets the DAC offset of the STG for a single channel in permanent EEProm memory

Parameters

DacChannel The DAC channel

Gets the ADC offset of the DACQ for a single channel

Parameters

AdcChannel	The ADC channel

Returns

The offset for the given channel number

Gets the parameter of the audio DAC output.

Frequency	Frequency(1 - 25000 Hz)
Amplification	Amplification(0 - 0xFFFF)

Gets the DAC Idle value

Parameters

```
DacChannel The DAC channel
```

Returns

The idle value

Gets the offset of a DAC channel.

Parameters

```
DacChannel The DAC channel
```

Returns

The offset for the given channel number

```
11.40.3.9 GetDacqRunStatus() uint16_t GetDacqRunStatus ( )
```

Gets the data acquisition running status

Returns

The status (1: running / 0: stopped)

```
11.40.3.10 GetDacUseldleValue() bool GetDacUseIdleValue ( uint32_t DacChannel )
```

Gets if the DAC Idle value is used after stimulation

Parameters

DacChannel	The DAC channel
------------	-----------------

Returns

Use idle value

11.40.3.11 GetDigInState() uint16_t GetDigInState ()

Reads from the LIH30 digital input

Returns

The bit mask defining the digital input state

Reads the requested amount of EEprom byte from the EPC10 EEPROM

Parameters

EEpromStartAddress	start address of memory area to read from
EEpromData_Length	The maximal length of EEpromData.
epc10bus	The EPC10 bus

Returns

pointer to internal memory for the requested amount of data

11.40.3.13 GetSampleInterval() uint32_t GetSampleInterval ()

Gets the Sample Interval for the DACQ and Stimulation

Returns

Sample Interval configured on the device

Is the User Trigger enabled

Returns

Enabled

Gets the clipping information

Parameters

```
epc10bus The EPC10 bus
```

Returns

The clipping value

11.40.3.16 ReadUARTData() String $^{\wedge}$ ReadUARTData ()

Reads the config string from the device.

Returns

The config string.

Send command to the EPC10

epc10bus	The EPC10 bus
Command	The command

Sets the ADC offset of the DACQ for a single channel

Parameters

AdcChannel	The ADC channel
Offset	The offset for the given channel number

Writes the ADC offset of the DACQ for a single channel to permanent EEProm memory

Parameters

AdcChannel The ADC channe

Sets the parameter of the audio DAC output.

Parameters

Frequency	Frequency(1 - 25000 Hz)
Amplification	Amplification(0 - 0xFFFF)

Sets the DAC Idle value

DacChannel	The DAC channel
IdleValue	The idle value

Sets the offset of a DAC channel.

Parameters

DacChannel	The DAC channel
Offset	The offset for the given channel number

Writes the DAC offset of the STG for a single channel to permanent EEProm memory

Parameters

DacChannel The DAC channel

Sets if the DAC Idle value is used after stimulation

Parameters

DacChannel	The DAC channel
Useldle	Use idle value

```
11.40.3.25 SetDigOutState() void SetDigOutState ( uint16_t DigOutState )
```

Writes to the LIH30 digital output

Parameters

DigOutState	The bit mask defining the digital output state

11.40.3.26 SetEEpromPage() void SetEEpromPage (

```
uint32_t EEpromStartAddress,
array< int8_t >^ EEpromData,
LIH30_EPC10_Bus_EnumNet epc10bus )
```

Writes into EEprom on the EPC10 EEPROM

Parameters

EEpromStartAddress	start address of memory area to write to
EEpromData	pointer to internal memory for the supported amount of data
epc10bus	The EPC10 bus

11.40.3.27 SetSampleInterval() void SetSampleInterval (uint32_t SampleInterval)

Sets the Sample Interval for the DACQ and Stimulation

Parameters

11.40.3.28 WriteUARTData() void WriteUARTData ($String^{\wedge} commandString$)

Write the command string to the device.

Parameters


commandString	The config string.
---------------	--------------------

11.40.4 Property Documentation

11.40.4.1 StimulusFunction CStimulusFunctionNet^ StimulusFunction [get]

11.41 CMcsBus_AxisParametersNet Class Reference

Inheritance diagram for CMcsBus_AxisParametersNet:



Public Member Functions

- CMcsBus AxisParametersNet (CMcsUsbNet[^] device)
- ~CMcsBus_AxisParametersNet (void)
- void SetAxisParametersEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short index, unsigned int parameter)
- void SetAxisParametersEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short index, int parameter)
- unsigned int GetAxisParametersUnsignedEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short index)
- int GetAxisParametersSignedEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short index)

Additional Inherited Members

11.41.1 Constructor & Destructor Documentation

11.41.2 Member Function Documentation

```
11.41.2.1 GetAxisParametersSignedEeprom() int GetAxisParametersSignedEeprom (
          unsigned char busaumber,
          unsigned char busaddress,
          unsigned char axis,
          unsigned short index )
```

```
11.41.2.2 GetAxisParametersUnsignedEeprom() unsigned int GetAxisParametersUnsignedEeprom ( unsigned char busnumber,
```

```
unsigned char busaddress,
unsigned char axis,
unsigned char axis,
unsigned short index)
```

$\textbf{11.41.2.3} \quad \textbf{SetAxisParametersEeprom()} \; \texttt{[1/2]} \quad \texttt{void SetAxisParametersEeprom ()}$


```
unsigned char busnumber,
unsigned char busaddress,
unsigned char axis,
unsigned short index,
int parameter)
```

11.41.2.4 SetAxisParametersEeprom() [2/2] void SetAxisParametersEeprom (

```
unsigned char busnumber,
unsigned char busaddress,
unsigned char axis,
unsigned short index,
unsigned int parameter)
```

11.42 CMcsBus_ExtensionNet Class Reference

Inheritance diagram for CMcsBus_ExtensionNet:



Public Member Functions

- CMcsBus_ExtensionNet (CMcsUsbNet[^] device)
- ~CMcsBus_ExtensionNet (void)
- void SetLEDSwitch (unsigned char busnumber, unsigned char busaddress, unsigned short LEDSwitch)
- unsigned short GetLEDSwitch (unsigned char busnumber, unsigned char busaddress)

Additional Inherited Members

11.42.1 Constructor & Destructor Documentation

```
11.42.1.2 ~CMcsBus_ExtensionNet() ~CMcsBus_ExtensionNet (
```


11.42.2 Member Function Documentation

```
11.42.2.1 GetLEDSwitch() unsigned short GetLEDSwitch (
    unsigned char busnumber,
    unsigned char busaddress)
```

```
11.42.2.2 SetLEDSwitch() void SetLEDSwitch (
unsigned char busnumber,
unsigned char busaddress,
unsigned short LEDSwitch )
```

11.43 CMcsBus_FYIExtensionNet Class Reference

Inheritance diagram for CMcsBus_FYIExtensionNet:



Public Member Functions

- CMcsBus FYIExtensionNet (CMcsUsbNet[^] device)
- ~CMcsBus_FYIExtensionNet (void)
- void SetValves (unsigned char busnumber, unsigned char busaddress, unsigned int states)
- unsigned int GetValves (unsigned char busnumber, unsigned char busaddress)
- void SetDIO (unsigned char busnumber, unsigned char busaddress, unsigned short io)
- unsigned short GetDIO (unsigned char busnumber, unsigned char busaddress)
- void SetSingleHeater (unsigned char busnumber, unsigned char busaddress, short index, unsigned short power)
- · unsigned short GetSingleHeater (unsigned char busnumber, unsigned char busaddress, short index)

Additional Inherited Members

11.43.1 Constructor & Destructor Documentation


```
11.43.1.1 CMcsBus_FYIExtensionNet() CMcsBus_FYIExtensionNet (
CMcsUsbNet^ device )
```

```
11.43.1.2 ~CMcsBus_FYIExtensionNet() ~CMcsBus_FYIExtensionNet (
             void )
11.43.2 Member Function Documentation
11.43.2.1 GetDIO() unsigned short GetDIO (
             unsigned char busnumber,
             unsigned char busaddress )
11.43.2.2 GetSingleHeater() unsigned short GetSingleHeater (
             unsigned char busnumber,
             unsigned char busaddress,
             short index )
11.43.2.3 GetValves() unsigned int GetValves (
             unsigned char busnumber,
             unsigned char busaddress )
11.43.2.4 SetDIO() void SetDIO (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short io )
11.43.2.5 SetSingleHeater() void SetSingleHeater (
             unsigned char busnumber,
             unsigned char busaddress,
             short index,
             unsigned short power )
11.43.2.6 SetValves() void SetValves (
             unsigned char busnumber,
             unsigned char busaddress,
```

unsigned int *states*)

11.44 CMcsBus MotorControlNet Class Reference

Inheritance diagram for CMcsBus_MotorControlNet:



Public Member Functions

- CMcsBus MotorControlNet (CMcsUsbNet[^] device)
- ~CMcsBus_MotorControlNet (void)
- void SetMCScalingFactorEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, int factor)
- int GetMCScalingFactorEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCScalingFactor (unsigned char busnumber, unsigned char busaddress, unsigned char axis, int factor)
- int GetMCScalingFactor (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxSpeedEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short speed)
- unsigned short GetMCMaxSpeedEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxSpeed (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short speed)
- unsigned short GetMCMaxSpeed (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxTravelEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, int travel)
- int GetMCMaxTravelEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxTravel (unsigned char busnumber, unsigned char busaddress, unsigned char axis, int travel)
- int GetMCMaxTravel (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCMaxCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCMaxCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCRegulatorGainEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short gain)
- short GetMCRegulatorGainEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCRegulatorGain (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short gain)
- short GetMCRegulatorGain (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxAccelerationEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short acceleration)
- unsigned short GetMCMaxAccelerationEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)

- void SetMCMaxAcceleration (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short acceleration)
- unsigned short GetMCMaxAcceleration (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCStandbyCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short percent)
- short GetMCStandbyCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCStandbyCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short percent)
- short GetMCStandbyCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCStandbyTimeEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short t)
- short GetMCStandbyTimeEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCStandbyTime (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short t)
- short GetMCStandbyTime (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCBreakCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCBreakCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCBreakCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCBreakCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCConfigEeprom (unsigned char busnumber, unsigned char busnumber, unsigned
- unsigned short GetMCConfigEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCConfig (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short config)
- unsigned short GetMCConfig (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCSpeedEeprom (unsigned char busnumber, unsigned c
- unsigned short GetMCSpeedEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCSpeed (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short speed)
- short GetMCSpeed (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCAccelerationEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short acceleration)
- unsigned short GetMCAccelerationEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCAcceleration (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short acceleration)
- unsigned short GetMCAcceleration (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCReferenceCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCReferenceCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCReferenceCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCReferenceCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrentModeEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, RoboCurrentModeEnumNet mode)
- RoboCurrentModeEnumNet GetMCCurrentModeEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)

- void SetMCCurrentMode (unsigned char busnumber, unsigned char busaddress, unsigned char axis, Robo

 CurrentModeEnumNet mode)
- RoboCurrentModeEnumNetGetMCCurrentMode (unsigned char busnumber, unsigned char busn
- void SetMCAxisRevisionEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short revision)
- unsigned short GetMCAxisRevisionEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCSpeedUnitEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis, int32 t speedunit)
- int32 t GetMCSpeedUnitEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCOutputOnOff (unsigned char busnumber, unsigned char busnumber, under busn
- bool GetMCOutputOnOff (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCSpeedShortCommand (unsigned char busnumber, unsigned char busnumber, unsi
- short GetMCSpeedShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCAccelerationShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short acceleration)
- unsigned short GetMCAccelerationShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrentShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis, short current)
- short GetMCCurrentShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCMaxTravelShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis. int travel)
- int GetMCMaxTravelShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrentPosition (unsigned char busnumber, unsigned char busnumber, under busnumber, under busnumber, under busnumber, under busnumber, under busnumb
- int GetMCCurrentPosition (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCNewPosition (unsigned char busnumber, unsigned char busnumber, under busn
- int GetMCNewPosition (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- short GetMCCurrentSpeed (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void StartMCMovement (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCRotation (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned char onoff)
- unsigned short GetMCMovement (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCReference (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned char switch_enable, unsigned char switch_polarity)
- unsigned char GetMCReference (unsigned char busnumber, unsigned char busaddress, unsigned char axis, [System::Runtime::InteropServices::Out]unsigned char% switch_port)
- void StopMCMovement (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetMCCurrentModeShortCommand (unsigned char busnumber, unsigned char busaddress, unsigned char axis, RoboCurrentModeEnumNet mode)
- RoboCurrentModeEnumNet GetMCCurrentModeShortCommand (unsigned char busnumber, unsigned char busnumber, unsigned char axis)
- unsigned short GetMCPhase (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- unsigned short GetMCPhaseOffset (unsigned char busnumber, unsigned char busaddress, unsigned char axis)
- void SetSubChannel (unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short subchannel)
- unsigned short GetSubChannel (unsigned char busnumber, unsigned char busaddress, unsigned char axis)

Additional Inherited Members

11.44.1 Constructor & Destructor Documentation

```
11.44.1.1 CMcsBus_MotorControlNet() CMcsBus_MotorControlNet (
             CMcsUsbNet^ device )
11.44.1.2 ~CMcsBus_MotorControlNet() ~CMcsBus_MotorControlNet (
             void )
11.44.2 Member Function Documentation
11.44.2.1 GetMCAcceleration() unsigned short GetMCAcceleration (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.2 GetMCAccelerationEeprom() unsigned short GetMCAccelerationEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.3 GetMCAccelerationShortCommand() unsigned short GetMCAccelerationShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.4 GetMCAxisRevisionEeprom() unsigned short GetMCAxisRevisionEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
```

```
11.44.2.5 GetMCBreakCurrent() short GetMCBreakCurrent (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
\textbf{11.44.2.6} \quad \textbf{GetMCBreakCurrentEeprom()} \quad \texttt{short GetMCBreakCurrentEeprom ()}
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
\textbf{11.44.2.7} \quad \textbf{GetMCConfig()} \quad \texttt{unsigned short GetMCConfig ()}
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.8 GetMCConfigEeprom() unsigned short GetMCConfigEeprom (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.9 GetMCCurrent() short GetMCCurrent (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.10 GetMCCurrentEeprom() short GetMCCurrentEeprom (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.11 GetMCCurrentMode() RoboCurrentModeEnumNet GetMCCurrentMode (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
```

```
11.44.2.12 GetMCCurrentModeEeprom() RoboCurrentModeEnumNet GetMCCurrentModeEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.13 GetMCCurrentModeShortCommand() RoboCurrentModeEnumNet GetMCCurrentModeShort←
Command (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.14 GetMCCurrentPosition() int GetMCCurrentPosition (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.15 GetMCCurrentShortCommand() short GetMCCurrentShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.16 GetMCCurrentSpeed() short GetMCCurrentSpeed (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.17 GetMCMaxAcceleration() unsigned short GetMCMaxAcceleration (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.18 GetMCMaxAccelerationEeprom() unsigned short GetMCMaxAccelerationEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
```

```
11.44.2.19 GetMCMaxCurrent() short GetMCMaxCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
\textbf{11.44.2.20} \quad \textbf{GetMCMaxCurrentEeprom()} \quad \texttt{short GetMCMaxCurrentEeprom ()}
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.21 GetMCMaxSpeed() unsigned short GetMCMaxSpeed (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.22 GetMCMaxSpeedEeprom() unsigned short GetMCMaxSpeedEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.23 GetMCMaxTravel() int GetMCMaxTravel (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.24 GetMCMaxTravelEeprom() int GetMCMaxTravelEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.25 GetMCMaxTravelShortCommand() int GetMCMaxTravelShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
```

```
{\bf 11.44.2.26} \quad {\bf GetMCMovement()} \quad {\tt unsigned \ short \ GetMCMovement \ (}
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
{\bf 11.44.2.27} \quad {\bf GetMCNewPosition()} \quad {\tt int GetMCNewPosition} \quad (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.28 GetMCOutputOnOff() bool GetMCOutputOnOff (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.29 GetMCPhase() unsigned short GetMCPhase (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.30 GetMCPhaseOffset() unsigned short GetMCPhaseOffset (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
11.44.2.31 GetMCReference() unsigned char GetMCReference (
              unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis,
              [System::Runtime::InteropServices::Out] unsigned char% switch_port )
11.44.2.32 GetMCReferenceCurrent() short GetMCReferenceCurrent (
             unsigned char busnumber,
              unsigned char busaddress,
              unsigned char axis )
```

```
\textbf{11.44.2.33} \quad \textbf{GetMCReferenceCurrentEeprom()} \quad \texttt{short GetMCReferenceCurrentEeprom} \quad \textbf{(}
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.34 GetMCRegulatorGain() short GetMCRegulatorGain (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.35 GetMCRegulatorGainEeprom() short GetMCRegulatorGainEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.36 GetMCScalingFactor() int GetMCScalingFactor (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.37 GetMCScalingFactorEeprom() int GetMCScalingFactorEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.38 GetMCSpeed() short GetMCSpeed (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.39 GetMCSpeedEeprom() unsigned short GetMCSpeedEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
```

```
11.44.2.40 GetMCSpeedShortCommand() short GetMCSpeedShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.41 GetMCSpeedUnitEeprom() int32_t GetMCSpeedUnitEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.42 GetMCStandbyCurrent() short GetMCStandbyCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.43 GetMCStandbyCurrentEeprom() short GetMCStandbyCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.44 GetMCStandbyTime() short GetMCStandbyTime (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.45 GetMCStandbyTimeEeprom() short GetMCStandbyTimeEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
11.44.2.46 GetSubChannel() unsigned short GetSubChannel (
            unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis )
```

```
11.44.2.47 SetMCAcceleration() void SetMCAcceleration (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short acceleration )
\textbf{11.44.2.48} \quad \textbf{SetMCAccelerationEeprom()} \quad \texttt{void SetMCAccelerationEeprom} \quad \textbf{(}
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short acceleration )
11.44.2.49 SetMCAccelerationShortCommand() void SetMCAccelerationShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short acceleration )
11.44.2.50 SetMCAxisRevisionEeprom() void SetMCAxisRevisionEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short revision )
11.44.2.51 SetMCBreakCurrent() void SetMCBreakCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.52 SetMCBreakCurrentEeprom() void SetMCBreakCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
```

```
11.44.2.53 SetMCConfig() void SetMCConfig (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short config )
\textbf{11.44.2.54} \quad \textbf{SetMCConfigEeprom()} \quad \texttt{void SetMCConfigEeprom ()}
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short config )
11.44.2.55 SetMCCurrent() void SetMCCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.56 SetMCCurrentEeprom() void SetMCCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.57 SetMCCurrentMode() void SetMCCurrentMode (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             RoboCurrentModeEnumNet mode )
11.44.2.58 SetMCCurrentModeEeprom() void SetMCCurrentModeEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             RoboCurrentModeEnumNet mode )
```

```
11.44.2.59 SetMCCurrentModeShortCommand() void SetMCCurrentModeShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             {\tt RoboCurrentModeEnumNet}\ \textit{mode}\ )
11.44.2.60 SetMCCurrentPosition() void SetMCCurrentPosition (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int position )
11.44.2.61 SetMCCurrentShortCommand() void SetMCCurrentShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.62 SetMCMaxAcceleration() void SetMCMaxAcceleration (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short acceleration )
11.44.2.63 SetMCMaxAccelerationEeprom() void SetMCMaxAccelerationEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short acceleration )
11.44.2.64 SetMCMaxCurrent() void SetMCMaxCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
```

```
11.44.2.65 SetMCMaxCurrentEeprom() void SetMCMaxCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.66 SetMCMaxSpeed() void SetMCMaxSpeed (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short speed )
11.44.2.67 SetMCMaxSpeedEeprom() void SetMCMaxSpeedEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short speed )
11.44.2.68 SetMCMaxTravel() void SetMCMaxTravel (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int travel )
11.44.2.69 SetMCMaxTravelEeprom() void SetMCMaxTravelEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int travel )
11.44.2.70 SetMCMaxTravelShortCommand() void SetMCMaxTravelShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int travel )
```


```
11.44.2.71 SetMCNewPosition() void SetMCNewPosition (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int position )
11.44.2.72 SetMCOutputOnOff() void SetMCOutputOnOff (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             bool OnOff_status )
11.44.2.73 SetMCReference() void SetMCReference (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned char switch_enable,
             unsigned char switch_polarity )
11.44.2.74 SetMCReferenceCurrent() void SetMCReferenceCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.75 SetMCReferenceCurrentEeprom() void SetMCReferenceCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short current )
11.44.2.76 SetMCRegulatorGain() void SetMCRegulatorGain (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short gain )
```

```
11.44.2.77 SetMCRegulatorGainEeprom() void SetMCRegulatorGainEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short gain )
11.44.2.78 SetMCRotation() void SetMCRotation (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned char onoff )
11.44.2.79 SetMCScalingFactor() void SetMCScalingFactor (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int factor )
11.44.2.80 SetMCScalingFactorEeprom() void SetMCScalingFactorEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int factor )
11.44.2.81 SetMCSpeed() void SetMCSpeed (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short speed )
11.44.2.82 SetMCSpeedEeprom() void SetMCSpeedEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short speed )
```

```
11.44.2.83 SetMCSpeedShortCommand() void SetMCSpeedShortCommand (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short speed )
11.44.2.84 SetMCSpeedUnitEeprom() void SetMCSpeedUnitEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int32_t speedunit )
11.44.2.85 SetMCStandbyCurrent() void SetMCStandbyCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short percent )
11.44.2.86 SetMCStandbyCurrentEeprom() void SetMCStandbyCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short percent )
11.44.2.87 SetMCStandbyTime() void SetMCStandbyTime (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short t )
11.44.2.88 SetMCStandbyTimeEeprom() void SetMCStandbyTimeEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short t )
```

11.45 CMcsBus SensorNet Class Reference

Inheritance diagram for CMcsBus SensorNet:



Public Member Functions

- CMcsBus SensorNet (CMcsUsbNet[^] device)
- ∼CMcsBus_SensorNet (void)
- void SetMinimalThreshold (unsigned char busnumber, unsigned char busaddress, unsigned short threshold)
- unsigned short GetMinimalThreshold (unsigned char busnumber, unsigned char busaddress)
- void SetDetectionThreshold (unsigned char busnumber, unsigned char busaddress, unsigned short threshold)
- unsigned short GetDetectionThreshold (unsigned char busnumber, unsigned char busaddress)
- void SetLatency (unsigned char busnumber, unsigned char busaddress, unsigned short latency)
- unsigned short GetLatency (unsigned char busnumber, unsigned char busaddress)
- unsigned short GetBubbleStatus (unsigned char busnumber, unsigned char busaddress)
- unsigned short GetLatencyCounter (unsigned char busnumber, unsigned char busaddress)
- unsigned short GetDetectorValue (unsigned char busnumber, unsigned char busaddress)
- array< int > ^ GetPressure (unsigned char busnumber, unsigned char busaddress, int n)
- int GetPressure (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetRegulatorOnOff (unsigned char busnumber, unsigned char busaddress, unsigned short index, unsigned char onoff)
- unsigned char GetRegulatorOnOff (unsigned char busnumber, unsigned char busaddress, unsigned short index)

- void SetSollPressure (unsigned char busnumber, unsigned char busaddress, unsigned short index, int pressure)
- int GetSollPressure (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetRegulatorFactor (unsigned char busnumber, unsigned char busaddress, unsigned short index, int factor)
- int GetRegulatorFactor (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetPressureOffset (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- array< unsigned short > ^ GetPressureOffset (unsigned char busnumber, unsigned char busaddress)
- int GetPressureOffset (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- unsigned int GetRegulatorStatus (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetRotatePump (unsigned char busnumber, unsigned char busaddress, unsigned short index, short speed)
- short GetRotatePump (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetMovePump (unsigned char busnumber, unsigned char busaddress, unsigned short index, unsigned short speed, int position)
- void SetRegulationTimeouts (unsigned char busnumber, unsigned char busaddress, unsigned short Max
 SpeedWait, unsigned short MaxSignChange)
- array< int > ^ Get4ADC (unsigned char busnumber, unsigned char busaddress)
- array< int > ^ Get4ADCAverage (unsigned char busnumber, unsigned char busaddress)
- void Set4DAC (unsigned char busnumber, unsigned char busaddress, array< unsigned short >^ dac)
- array< unsigned short > ^ Get4DAC (unsigned char busnumber, unsigned char busaddress)
- void Set4ADCMode (unsigned char busnumber, unsigned char busaddress, PatchServAdcModeEnumNet mode)
- PatchServAdcModeEnumNet Get4ADCMode (unsigned char busnumber, unsigned char busaddress)
- void Set4ADCCatchampAverageShift (unsigned char busnumber, unsigned char busaddress, unsigned int shift)
- unsigned int Get4ADCCatchampAverageShift (unsigned char busnumber, unsigned char busaddress)
- array< unsigned short > ^ Get2AnalogInput (unsigned char busnumber, unsigned char busaddress)
- unsigned short Get2DigitalInput (unsigned char busnumber, unsigned char busaddress)
- array< unsigned short > ^ GetADCs (unsigned char busnumber, unsigned char busaddress, int n)
- array< unsigned short > ^ GetADCsLoop (unsigned char busnumber, unsigned char busaddress, int n)
- void SetPiezoState (unsigned char busnumber, unsigned char busaddress, int state)
- void GetPiezoState (unsigned char busnumber, unsigned char busaddress, [System::Runtime::Interop←
 Services::Out]int% state, [System::Runtime::InteropServices::Out]int% reason)
- void SetDACs (unsigned char busnumber, unsigned char busaddress, unsigned short index, array< unsigned short >^ dac_times_voltages)
- array< unsigned short > ^ GetDACs (unsigned char busnumber, unsigned char busaddress, unsigned short index)
- void SetSamplePeriode (unsigned char busnumber, unsigned char busaddress, unsigned short periode)
- unsigned short GetSamplePeriode (unsigned char busnumber, unsigned char busaddress)
- void StartSync (unsigned char busnumber, unsigned char busaddress)
- unsigned short GetSyncState (unsigned char busnumber, unsigned char busaddress)
- void CatchAmpSetDacAmplitude (unsigned char busnumber, unsigned char busaddress, unsigned short dacAmplitude)
- unsigned short CatchAmpGetDacAmplitude (unsigned char busnumber, unsigned char busaddress)
- void CatchAmpSetDacOffset (unsigned char busnumber, unsigned char busaddress, short dacOffset)
- short CatchAmpGetDacOffset (unsigned char busnumber, unsigned char busaddress)
- int CatchAmpGetAdcMean (unsigned char busnumber, unsigned char busaddress)
- int CatchAmpGetAdcValue (unsigned char busnumber, unsigned char busaddress)

- int CatchAmpGetAdcValueH (unsigned char busnumber, unsigned char busaddress)
- int CatchAmpGetAdcValueL (unsigned char busnumber, unsigned char busaddress)
- void CatchAmpSetPwmEnable (unsigned char busnumber, unsigned char busaddress, bool pwmEnable)
- bool CatchAmpGetPwmEnable (unsigned char busnumber, unsigned char busaddress)
- void CatchAmpSetDacEnable (unsigned char busnumber, unsigned char busaddress, bool dacEnable)
- bool CatchAmpGetDacEnable (unsigned char busnumber, unsigned char busaddress)
- int TactSwitchGetState (unsigned char busnumber, unsigned char busaddress)
- · void TactSwitchSetDisplay (unsigned char busnumber, unsigned char busaddress, int Melody)

Additional Inherited Members

11.45.1 Constructor & Destructor Documentation

11.45.2 Member Function Documentation

```
11.45.2.1 CatchAmpGetAdcMean() int CatchAmpGetAdcMean (
unsigned char busnumber,
unsigned char busaddress)
```

```
11.45.2.2 CatchAmpGetAdcValue() int CatchAmpGetAdcValue ( unsigned char busnumber,
```

unsigned char busaddress)

unsigned char busaddress)

```
11.45.2.3 CatchAmpGetAdcValueH() int CatchAmpGetAdcValueH ( unsigned char busnumber,
```

```
11.45.2.4 CatchAmpGetAdcValueL() int CatchAmpGetAdcValueL (
             unsigned char busnumber,
             unsigned char busaddress )
\textbf{11.45.2.5} \quad \textbf{CatchAmpGetDacAmplitude()} \quad \texttt{unsigned short CatchAmpGetDacAmplitude ()}
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.6 CatchAmpGetDacEnable() bool CatchAmpGetDacEnable (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.7 CatchAmpGetDacOffset() short CatchAmpGetDacOffset (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.8 CatchAmpGetPwmEnable() bool CatchAmpGetPwmEnable (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.9 CatchAmpSetDacAmplitude() void CatchAmpSetDacAmplitude (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short dacAmplitude )
11.45.2.10 CatchAmpSetDacEnable() void CatchAmpSetDacEnable (
             unsigned char busnumber,
             unsigned char busaddress,
             bool dacEnable )
11.45.2.11 CatchAmpSetDacOffset() void CatchAmpSetDacOffset (
             unsigned char busnumber,
             unsigned char busaddress,
             short dacOffset )
```

```
11.45.2.12 CatchAmpSetPwmEnable() void CatchAmpSetPwmEnable (
             unsigned char busnumber,
             unsigned char busaddress,
             bool pwmEnable )
11.45.2.13 Get2AnalogInput() array<unsigned short> ^ Get2AnalogInput (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.14 Get2DigitalInput() unsigned short Get2DigitalInput (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.15 Get4ADC() array<int> ^ Get4ADC (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.16 Get4ADCAverage() array<int> ^ Get4ADCAverage (
             unsigned char busnumber,
             unsigned char busaddress )
\textbf{11.45.2.17} \quad \textbf{Get4ADCCatchampAverageShift()} \quad \textbf{unsigned int Get4ADCCatchampAverageShift ()}
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.18 Get4ADCMode() PatchServAdcModeEnumNet Get4ADCMode (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.19 Get4DAC() array<unsigned short> ^{\land} Get4DAC (
             unsigned char busnumber,
             unsigned char busaddress )
```

```
11.45.2.20 GetADCs() array<unsigned short> ^{\wedge} GetADCs (
             unsigned char busnumber,
             unsigned char busaddress,
             int n)
11.45.2.21 GetADCsLoop() array<unsigned short> ^{\wedge} GetADCsLoop (
             unsigned char busnumber,
             unsigned char busaddress,
             int n)
11.45.2.22 GetBubbleStatus() unsigned short GetBubbleStatus (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.23 GetDACs() array<unsigned short> ^ GetDACs (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.24 GetDetectionThreshold() unsigned short GetDetectionThreshold (
             unsigned char busnumber,
             unsigned char busaddress )
{\bf 11.45.2.25} \quad {\bf GetDetectorValue()} \quad {\tt unsigned short GetDetectorValue ()}
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.26 GetLatency() unsigned short GetLatency (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.27 GetLatencyCounter() unsigned short GetLatencyCounter (
             unsigned char busnumber,
             unsigned char busaddress )
```

```
11.45.2.28 GetMinimalThreshold() unsigned short GetMinimalThreshold (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.29 GetMovePump() void GetMovePump (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             [System::Runtime::InteropServices::Out] unsigned short% speed,
             [System::Runtime::InteropServices::Out] int% position )
11.45.2.30 GetPiezoState() void GetPiezoState (
             unsigned char busnumber,
             unsigned char busaddress,
             [System::Runtime::InteropServices::Out] int% state,
             [System::Runtime::InteropServices::Out] int% reason )
11.45.2.31 GetPressure() [1/2] array<int> ^ GetPressure (
             unsigned char busnumber,
             unsigned char busaddress,
             int n)
11.45.2.32 GetPressure() [2/2] int GetPressure (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index )
11.45.2.33 GetPressureOffset() [1/2] array<unsigned short> ^ GetPressureOffset (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.34 GetPressureOffset() [2/2] int GetPressureOffset (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index )
```

```
11.45.2.35 GetRegulationTimeouts() void GetRegulationTimeouts (
             unsigned char busnumber,
             unsigned char busaddress,
             [System::Runtime::InteropServices::Out] unsigned short% MaxSpeedWait,
             [System::Runtime::InteropServices::Out] unsigned short% MaxSignChange )
11.45.2.36 GetRegulatorFactor() int GetRegulatorFactor (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.37 GetRegulatorOnOff() unsigned char GetRegulatorOnOff (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.38 GetRegulatorStatus() unsigned int GetRegulatorStatus (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.39 GetRotatePump() short GetRotatePump (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.40 GetSamplePeriode() unsigned short GetSamplePeriode (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.41 GetSollPressure() int GetSollPressure (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
```


```
11.45.2.42 GetSyncState() unsigned short GetSyncState (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.43 Set4ADCCatchampAverageShift() void Set4ADCCatchampAverageShift (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned int shift )
11.45.2.44 Set4ADCMode() void Set4ADCMode (
             unsigned char busnumber,
             unsigned char busaddress,
             PatchServAdcModeEnumNet mode )
11.45.2.45 Set4DAC() void Set4DAC (
             unsigned char busnumber,
             unsigned char busaddress,
             array< unsigned short >^{\wedge} dac )
11.45.2.46 SetDACs() void SetDACs (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             array< unsigned short >^{\wedge} dac\_times\_voltages )
11.45.2.47 SetDetectionThreshold() void SetDetectionThreshold (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short threshold )
11.45.2.48 SetLatency() void SetLatency (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short latency )
```

```
11.45.2.49 SetMinimalThreshold() void SetMinimalThreshold (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short threshold )
11.45.2.50 SetMovePump() void SetMovePump (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             unsigned short speed,
             int position )
11.45.2.51 SetPiezoState() void SetPiezoState (
             unsigned char busnumber,
             unsigned char busaddress,
             int state )
11.45.2.52 SetPressureOffset() void SetPressureOffset (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index)
11.45.2.53 SetRegulationTimeouts() void SetRegulationTimeouts (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short MaxSpeedWait,
             unsigned short MaxSignChange )
11.45.2.54 SetRegulatorFactor() void SetRegulatorFactor (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             int factor )
11.45.2.55 SetRegulatorOnOff() void SetRegulatorOnOff (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             unsigned char onoff )
```

```
11.45.2.56 SetRotatePump() void SetRotatePump (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             short speed )
11.45.2.57 SetSamplePeriode() void SetSamplePeriode (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short periode )
11.45.2.58 SetSollPressure() void SetSollPressure (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short index,
             int pressure )
11.45.2.59 StartSync() void StartSync (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.60 TactSwitchGetState() int TactSwitchGetState (
             unsigned char busnumber,
             unsigned char busaddress )
11.45.2.61 TactSwitchSetDisplay() void TactSwitchSetDisplay (
             unsigned char busnumber,
             unsigned char busaddress,
             int Melody )
```

11.46 CMcsBus_TempSensorNet Class Reference

Inheritance diagram for CMcsBus_TempSensorNet:



Public Member Functions

- CMcsBus_TempSensorNet (CMcsUsbNet[^] device)
- ~CMcsBus_TempSensorNet (void)
- short GetTemperatur (unsigned char busnumber, unsigned char busaddress)
- short GetTemperatur (unsigned char busnumber, unsigned char busaddress, short index)
- void SetNanoVoltsPerKelvin (unsigned char busnumber, unsigned char busaddress, int nanovoltsperkelvin)
- int GetNanoVoltsPerKelvin (unsigned char busnumber, unsigned char busaddress)
- short GetThermoVoltage (unsigned char busnumber, unsigned char busaddress, short index)
- short GetThermoTemp (unsigned char busnumber, unsigned char busaddress, short index)
- void SetThermoOffset (unsigned char busnumber, unsigned char busaddress, short index, short offset)
- short GetThermoOffset (unsigned char busnumber, unsigned char busaddress, short index)

Additional Inherited Members


11.46.1 Constructor & Destructor Documentation

```
11.46.2.3 GetTemperatur() [2/2] short GetTemperatur (
unsigned char busnumber,
unsigned char busaddress,
short index)
```

```
11.46.2.4 GetThermoOffset() short GetThermoOffset (
             unsigned char busnumber,
             unsigned char busaddress,
             short index )
11.46.2.5 GetThermoTemp() short GetThermoTemp (
             unsigned char busnumber,
             unsigned char busaddress,
             short index )
11.46.2.6 GetThermoVoltage() short GetThermoVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             short index )
11.46.2.7 SetNanoVoltsPerKelvin() void SetNanoVoltsPerKelvin (
             unsigned char busnumber,
             unsigned char busaddress,
             int nanovoltsperkelvin )
11.46.2.8 SetThermoOffset() void SetThermoOffset (
             unsigned char busnumber,
             unsigned char busaddress,
             short index,
             short offset )
```

11.47 CMcsBus_VoltageModeNet Class Reference

Inheritance diagram for CMcsBus_VoltageModeNet:



Public Member Functions

- CMcsBus_VoltageModeNet (CMcsUsbNet[^] device)
- ~CMcsBus_VoltageModeNet (void)
- void SetVMMaxPositiveCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short current)
- short GetVMMaxPositiveCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxPositiveCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short current)
- short GetVMMaxPositiveCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxNegativeCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short current)
- short GetVMMaxNegativeCurrentEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxNegativeCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short current)
- short GetVMMaxNegativeCurrent (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxPositiveVoltageEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short voltage)
- short GetVMMaxPositiveVoltageEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxPositiveVoltage (unsigned char busnumber, unsigned char
- short GetVMMaxPositiveVoltage (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxNegativeVoltageEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short voltage)
- short GetVMMaxNegativeVoltageEeprom (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMMaxNegativeVoltage (unsigned char busnumber, unsigned char busaddress, unsigned char channel, short voltage)
- short GetVMMaxNegativeVoltage (unsigned char busnumber, unsigned char busaddress, unsigned char channel)
- void SetVMOutputOnOff (unsigned char busnumber, unsigned char busaddress, unsigned char channel, unsigned short status)
- unsigned short GetVMOutputOnOff (unsigned char busnumber, unsigned char busnumber, under busnumber, under busnumber, under busnumber, under busnumber, under busnumber, under busnumber
- void SetVMVoltage (unsigned char busnumber, unsigned char busnumber, under busnum
- short GetVMVoltage (unsigned char busnumber, unsigned char busaddress, unsigned char channel)

Additional Inherited Members

11.47.1 Constructor & Destructor Documentation


```
11.47.1.2 ~ CMcsBus_VoltageModeNet() ~ CMcsBus_VoltageModeNet (
             void )
11.47.2 Member Function Documentation
\textbf{11.47.2.1} \quad \textbf{GetVMMaxNegativeCurrent()} \quad \texttt{short GetVMMaxNegativeCurrent ()}
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.2 GetVMMaxNegativeCurrentEeprom() short GetVMMaxNegativeCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.3 GetVMMaxNegativeVoltage() short GetVMMaxNegativeVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.4 GetVMMaxNegativeVoltageEeprom() short GetVMMaxNegativeVoltageEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.5 GetVMMaxPositiveCurrent() short GetVMMaxPositiveCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.6 GetVMMaxPositiveCurrentEeprom() short GetVMMaxPositiveCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
```


unsigned char channel)

```
11.47.2.7 GetVMMaxPositiveVoltage() short GetVMMaxPositiveVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.8 GetVMMaxPositiveVoltageEeprom() short GetVMMaxPositiveVoltageEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.9 GetVMOutputOnOff() unsigned short GetVMOutputOnOff (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.10 GetVMVoltage() short GetVMVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel )
11.47.2.11 SetVMMaxNegativeCurrent() void SetVMMaxNegativeCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short current )
11.47.2.12 SetVMMaxNegativeCurrentEeprom() void SetVMMaxNegativeCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short current )
11.47.2.13 SetVMMaxNegativeVoltage() void SetVMMaxNegativeVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short voltage )
```

```
11.47.2.14 SetVMMaxNegativeVoltageEeprom() void SetVMMaxNegativeVoltageEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short voltage )
11.47.2.15 SetVMMaxPositiveCurrent() void SetVMMaxPositiveCurrent (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short current )
11.47.2.16 SetVMMaxPositiveCurrentEeprom() void SetVMMaxPositiveCurrentEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short current )
11.47.2.17 SetVMMaxPositiveVoltage() void SetVMMaxPositiveVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short voltage )
11.47.2.18 SetVMMaxPositiveVoltageEeprom() void SetVMMaxPositiveVoltageEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short voltage )
11.47.2.19 SetVMOutputOnOff() void SetVMOutputOnOff (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             unsigned short status )
11.47.2.20 SetVMVoltage() void SetVMVoltage (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char channel,
             short voltage )
```

11.48 CMcsBusNet Class Reference

Inheritance diagram for CMcsBusNet:



Public Member Functions

- CMcsBusNet (CMcsUsbNet[^] device)
- virtual ∼CMcsBusNet (void)
- void SetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned short value)
- void SetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, short value)
- void SetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, unsigned int value)
- void SetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, int value)
- void GetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, [System::Runtime::InteropServices::Out]unsigned short% value)
- void GetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, [System::Runtime::InteropServices::Out]short% value)
- void GetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, [System::Runtime::InteropServices::Out]unsigned int% value)
- void GetCommand (unsigned char command, unsigned char busnumber, unsigned char busaddress, unsigned char axis, [System::Runtime::InteropServices::Out]int% value)
- void SetBusAddressEeprom (unsigned char busnumber, unsigned char busaddress, unsigned short newaddress)
- unsigned short GetBusAddressEeprom (unsigned char busnumber, unsigned char busaddress)
- void SetBusAddress (unsigned char busnumber, unsigned char busaddress, unsigned short newaddress)
- unsigned short GetBusAddress (unsigned char busnumber, unsigned char busaddress)
- void CMcsBusNet::SetModeEeprom (unsigned char busnumber, unsigned char busaddress, unsigned short mode)
- unsigned short CMcsBusNet::GetModeEeprom (unsigned char busnumber, unsigned char busaddress)
- void CMcsBusNet::SetMode (unsigned char busnumber, unsigned char busaddress, unsigned short mode)
- unsigned short CMcsBusNet::GetMode (unsigned char busnumber, unsigned char busaddress)
- void SetHWRevisionEeprom (unsigned char busnumber, unsigned char busaddress, unsigned short revision)
- unsigned short GetHWRevisionEeprom (unsigned char busnumber, unsigned char busaddress)

Additional Inherited Members

11.48.1 Constructor & Destructor Documentation

```
11.48.1.1 CMcsBusNet() CMcsBusNet (
CMcsUsbNet^ device)
```

```
11.48.1.2 \sim CMcsBusNet() virtual \sim CMcsBusNet (
             void ) [virtual]
11.48.2 Member Function Documentation
11.48.2.1 CMcsBusNet::GetMode() unsigned short CMcsBusNet::GetMode (
             unsigned char busnumber,
             unsigned char busaddress )
11.48.2.2 CMcsBusNet::GetModeEeprom() unsigned short CMcsBusNet::GetModeEeprom (
             unsigned char busnumber,
             unsigned char busaddress )
11.48.2.3 CMcsBusNet::SetMode() void CMcsBusNet::SetMode (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short mode )
11.48.2.4 CMcsBusNet::SetModeEeprom() void CMcsBusNet::SetModeEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short \mathit{mode} )
11.48.2.5 GetBusAddress() unsigned short GetBusAddress (
             unsigned char busnumber,
             unsigned char busaddress )
```

 ${\bf 11.48.2.6} \quad {\bf GetBusAddressEeprom()} \quad {\tt unsigned \ short \ GetBusAddressEeprom \ ()}$

unsigned char busnumber,
unsigned char busaddress)


```
11.48.2.7 GetCommand() [1/4] void GetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             [System::Runtime::InteropServices::Out] int% value )
11.48.2.8 GetCommand() [2/4] void GetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             [System::Runtime::InteropServices::Out] short% value )
11.48.2.9 GetCommand() [3/4] void GetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             [{\tt System::Runtime::InteropServices::Out}] \  \, {\tt unsigned int} \% \  \, {\tt value} \ )
11.48.2.10 GetCommand() [4/4] void GetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             [System::Runtime::InteropServices::Out] unsigned short% value )
11.48.2.11 GetHWRevisionEeprom() unsigned short GetHWRevisionEeprom (
             unsigned char busnumber,
             unsigned char busaddress )
11.48.2.12 SetBusAddress() void SetBusAddress (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short newaddress )
```

```
11.48.2.13 SetBusAddressEeprom() void SetBusAddressEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short newaddress )
11.48.2.14 SetCommand() [1/4] void SetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             int value )
11.48.2.15 SetCommand() [2/4] void SetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             short value )
11.48.2.16 SetCommand() [3/4] void SetCommand (
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned int value )
\textbf{11.48.2.17} \quad \textbf{SetCommand()} \; \textbf{[4/4]} \quad \text{void SetCommand (}
             unsigned char command,
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned char axis,
             unsigned short value )
11.48.2.18 SetHWRevisionEeprom() void SetHWRevisionEeprom (
             unsigned char busnumber,
             unsigned char busaddress,
             unsigned short revision )
```

11.49 CMcsUsbDacqNet Class Reference

Base class for data acquisition devices.

Inheritance diagram for CMcsUsbDacqNet:



Classes

· class CHWInfo

Class to provide hardware information about the device.

Public Member Functions

- CMcsUsbDacqNet ()
- ∼CMcsUsbDacqNet ()
- virtual uint32 t GetVoltageRangeIndex (unsigned int virtualDevice)
- virtual void SetVoltageRangeByIndex (int32_t voltageRangeIndex, unsigned int virtualDevice)

Sets the voltage range on devices which support multiple voltage ranges.

virtual void SetVoltageRangeInMicroVolt (int32_t voltageRange, unsigned int virtualDevice)

Sets the voltage range on devices which support multiple voltage ranges.

virtual int32_t GetVoltageRangeInMicroVolt (unsigned int virtualDevice)

Gets the currently selected voltage range on devices which support multiple voltage ranges.

virtual int32_t GetVoltageRangeInMilliVolt ()

Gets the currently selected voltage range on devices which support multiple voltage ranges.

virtual void SetDataMode (DataModeEnumNet dataMode, unsigned int virtualDevice)

Sets the data mode, can be 16, 24 or 32bit, all signed or unsigned on the MEA2100 device.

virtual DataModeEnumNet GetDataMode (unsigned int virtualDevice)

Gets the data mode, can be 16, 24 or 32bit, all signed or unsigned on the MEA2100 device.

void SetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, DigitalSourceEnumNet source, int bitnumber_offset)

Sets the function/source of an digital output bit.

 void SetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, W2100DigitalSourceEnumNet source, int bitnumber_offset)

Sets the function/source of an digital output bit.

 void SetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, SCUDigitalSourceEnumNet source, int bitnumber offset)

Sets the function/source of an digital output bit.

void SetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, MEA2100_256DigitalSource
 — EnumNet source, int bitnumber_offset)

Sets the function/source of an digital output bit.

template<typename digitalsourceenum >
 void SetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, DigitalSource< digital-sourceenum >^ source, int bitnumber_offset)

Sets the function/source of an digital output bit.

void GetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, [System::Runtime::Interop
 — Services::Out]DigitalSourceEnumNet% source, [System::Runtime::InteropServices::Out]int% bitnumber
 — offset)

Gets the function/source of an digital output bit.

void GetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, [System::Runtime::
 — InteropServices::Out]W2100DigitalSourceEnumNet% source, [System::Runtime::InteropServices::Out]int% bitnumber offset)

Gets the function/source of an digital output bit.

void GetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, [System::Runtime::
 InteropServices::Out]SCUDigitalSourceEnumNet% source, [System::Runtime::InteropServices::Out]int% bitnumber offset)

Gets the function/source of an digital output bit.

Gets the function/source of an digital output bit.

• template<typename digitalsourceenum >

void GetDigitalSource (DigitalTargetEnumNet digitaltarget, int32_t NrChannel, [System::Runtime::Interop ← Services::Out]DigitalSource < digitalsourceenum > % source, [System::Runtime::InteropServices::Out]int% bitnumber_offset)

Gets the function/source of an digital output bit.

virtual AdapterTypeEnumNet GetAdapterType ()

Gets the adapter which is connected to the MEA2100 device.

virtual MeaLayoutEnumNet GetMeaLayout ()

Gets the MEA layout which is connected to the MEA2100 device.

virtual uint32_t GetAdcDataFormat (uint32_t virtualDevice)

Gets the ADC data format, 16 means 16 bits, 24 means 24 bits, 32 means 32 bits.

- virtual uint32_t GetResolutionPerDigit (uint32_t virtualDevice, DacqGroupChannelEnumNet group, [System::Runtime::InteropServices::Out] int% res, [System::Runtime::InteropServices::Out] int% resUnit)
- virtual uint32_t GetHardwareMinRange (uint32_t virtualDevice, DacqGroupChannelEnumNet group, [System::Runtime::InteropServices::Out] int% r, [System::Runtime::InteropServices::Out] int% rUnit)
- virtual uint32_t GetHardwareMaxRange (uint32_t virtualDevice, DacqGroupChannelEnumNet group, [System::Runtime::InteropServices::Out] int% r, [System::Runtime::InteropServices::Out] int% rUnit)
- virtual uint32_t GetDataFormat (uint32_t virtualDevice, DacqGroupChannelEnumNet group, [System::
 Runtime::InteropServices::Out] int% numberOfBits)
- virtual uint32_t GetNumberOfDataBits (uint32_t virtualDevice, DacqGroupChannelEnumNet group, [System::Runtime::InteropServices::Out] int% numberOfBits)

Get the real number of data bits.

• virtual void SetSamplerate (int32_t rate, unsigned int oversample, unsigned int virtualDevice)

Sets the sampling frequency of the device.

virtual int32 t GetSamplerate (unsigned int virtualDevice)

Gets the sampling frequency of the device.

virtual uint32 t GetMaxSamplingFrequency (int virtualDevice)

Gets the maximal sampling frequency of the device.

virtual uint32 t GetMinSamplingFrequencyStepsize ()

Gets the minimal sampling frequency step size increment value of the device.

virtual int32_t GetChannelsInBlock (unsigned int virtualDevice)

Get the number of 16 bit datawords which will be collected per sample frame, use after the device is configured.

- virtual void SendStartDacq ()

Start sampling.

virtual void SendStartDacq (int VirtualDacqMap)

Start sampling.

virtual void SendStartStgAndDacq (uint32_t trigger_map, int VirtualDacqMap)

Start sampling together with the STG.

virtual void SendStopDacq ()

Stop sampling.

virtual void SendStopDacq (int VirtualDacqMap)

Stop sampling.

Parameters

VirtualDacqMap

virtual void SendStopStgAndDacq (uint32_t trigger_map, int VirtualDacqMap)

Stop sampling together with the STG.

virtual void SendStopStgAndDacqWithOptions (uint32_t trigger_map, int VirtualDacqMap, int options)

Stop sampling together with the STG and options.

• virtual void StartLoop ()

Start the data acquisition thread.

virtual void StartLoop (int32_t timeout)

Start the data acquisition thread.

virtual void StartLoop (int32_t timeout, int32_t numSubmittedUsbBuffers, int32_t numUsbBuffers, int32_← t packetsInUrb)

Start the data acquisition thread.

virtual void StartLoop (int32_t timeout, int32_t numSubmittedUsbBuffers, int32_t numUsbBuffers, int32_t packetsInUrb, uint32_t virtualDevice)

Start the data acquisition thread.

- virtual void StopLoop ()
- virtual void ClearBuffers ()
- virtual void StartDacq ()

Start the data acquisition thread and sampling.

virtual void StartDacq (int32 t timeout)

Start the data acquisition thread and sampling.

virtual void StartDacq (int32_t timeout, int32_t numSubmittedUsbBuffers, int32_t numUsbBuffers, int32_← t packetsInUrb)

Start the data acquisition thread and sampling.

 virtual void StartDacq (int32_t timeout, int32_t numSubmittedUsbBuffers, int32_t numUsbBuffers, int32_← t packetsInUrb, uint32_t virtualDevice)

Start the data acquisition thread and sampling.

virtual void StopDacq ()

Stop the data acquisition thread and sampling.

virtual void StopDacq (uint32 t virtualDevice)

Stop the data acquisition thread and sampling.

- virtual uint32_t SetPoti (uint32_t channel, uint32_t value, bool write_nvram)
- virtual uint32_t GetPoti (uint32_t channel, [System::Runtime::InteropServices::Out]uint32_t% value)
- virtual CFilterPropertyNet ^ GetFilterProperty (DacqGroupChannelEnumNet GroupID, unsigned int index)
- virtual array< CFilterPropertyNet[^]> [^] CMcsUsbDacqNet::GetFilterProperties (DacqGroupChannelEnum← Net GroupID)
- int GetChannelDataFillSize ()
- virtual void SetSelectedChannels (int nChannels, int queuesize, int threshold, SampleSizeNet samplesize, int ChannelsInBlock)

Create a FIFO queue per channel. Each channel will have its own FIFO and Callback function.

- virtual void SetSelectedChannels (int nChannels, int queuesize, int threshold, SampleSizeNet sampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual void SetSelectedChannels (array< bool >[^] selectedChannels, int queuesize, int threshold, Sample
 SizeNet samplesize, int ChannelsInBlock)

Create a FIFO queue per channel. Each channel will have its own FIFO and Callback function.

- virtual void SetSelectedChannels (array < bool > ^ selectedChannels, int queuesize, int threshold, Sample ← SizeNet sampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual void SetSelectedData (int nChannels, int queuesize, int threshold, SampleSizeNet samplesize, int ChannelsInBlock)

Create a common FIFO queue for all channels. Use handle = 0 in the ChannelBlock_ReadFrames... functions.

- virtual void SetSelectedData (int nChannels, int queuesize, int threshold, SampleSizeNet sampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual void SetSelectedData (array < bool >^ selectedChannels, int queuesize, int threshold, SampleSizeNet samplesize, int ChannelsInBlock)

Create a common FIFO queue for all channels. Use handle = 0 in the ChannelBlock_ReadFrames... functions.

- virtual void SetSelectedData (array < bool >^ selectedChannels, int queuesize, int threshold, SampleSizeNet sampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual int AddSelectedChannelsQueue (int nByteOffset, int nChannelOffset, int nChannels, int queuesize, int threshold, SampleSizeNet samplesize)

Adds a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_Read← FramesDict... with handle = 0 to read the data.

- virtual int AddSelectedChannelsQueue (int nByteOffset, int nChannelOffset, int nChannels, int queuesize, int threshold, SampleSizeNet samplesize, SampleDstSizeNet sampleDstSize)
- virtual int AddSelectedChannelsQueue (int nByteOffset, int nChannelOffset, array< bool >[∧] selected ← Channels, int queuesize, int threshold, SampleSizeNet samplesize)

Adds a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_Read← FramesDict... with handle = 0 to read the data.

- virtual int AddSelectedChannelsQueue (int nByteOffset, int nChannelOffset, array< bool >[^] selected←
 Channels, int queuesize, int threshold, SampleSizeNet sampleDstSizeNet sampleDstSizeNet
- virtual void SetSelectedChannelsQueue (int nChannels, int queuesize, int threshold, SampleSizeNet sample-size, int ChannelsInBlock)

Create a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_← ReadFramesDict... with handle = 0 to read the data.

- virtual void SetSelectedChannelsQueue (int nChannels, int queuesize, int threshold, SampleSizeNet samplesize, size, SampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual void SetSelectedChannelsQueue (array< bool >^ selectedChannels, int queuesize, int threshold, SampleSizeNet samplesize, int ChannelsInBlock)

Create a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_← ReadFramesDict... with handle = 0 to read the data.

- virtual void SetSelectedChannelsQueue (array< bool >^ selectedChannels, int queuesize, int threshold, SampleSizeNet sampleDstSizeNet sampleDstSizeNet sampleDstSize, int ChannelsInBlock)
- virtual uint32_t ChannelBlock_AvailFrames (int handle)

Get the number of sample frames already available in the FIFO.

- virtual uint32_t ChannelBlock_AvailFrames (int handle, int queue)

Read data from a FIFO queue in uint16_t data format

• virtual uint32_t ChannelBlock_ReadFramesUI16 (int handle, array< uint16_t >^ buffer, int frames_pos, int frames, [System::Runtime::InteropServices::Out]int % frames ret)

Read data from a FIFO queue in uint16 t data format

virtual array< int16_t > ^ ChannelBlock_ReadFramesI16 (int handle, int frames, [System::Runtime::←
 InteropServices::Out]int % frames ret)

Read data from a FIFO queue in int16_t data format

• virtual uint32_t ChannelBlock_ReadFramesI16 (int handle, array< int16_t >^ buffer, int frames_pos, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue in int16_t data format

Read data from a FIFO queue in uint32_t data format

• virtual uint32_t ChannelBlock_ReadFramesUl32 (int handle, array< uint32_t >^ buffer, int frames_pos, int frames, [System::Runtime::InteropServices::Out]int % frames ret)

Read data from a FIFO queue in uint32_t data format

virtual array< int32_t > ^ ChannelBlock_ReadFramesl32 (int handle, int frames, [System::Runtime::←
 InteropServices::Out]int % frames ret)

Read data from a FIFO queue in uint32_t data format

• virtual uint32_t ChannelBlock_ReadFramesI32 (int handle, array< int32_t >^ buffer, int frames_pos, int frames, [System::Runtime::InteropServices::Out]int % frames ret)

Read data from a FIFO queue in uint32_t data format

• virtual array< array< uint16_t >^> ^ ChannelBlock_ReadAsFrameArrayUI16 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

• virtual array< array< uint16_t >^> ^ ChannelBlock_ReadAsFrameArrayUI16 (int handle, int queue, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

virtual array < array < int16_t >^ > ^ ChannelBlock_ReadAsFrameArrayI16 (int handle, int frames, [System ← ::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

virtual array< array< int16_t >^> ^ ChannelBlock_ReadAsFrameArrayI16 (int handle, int queue, int frames,
 [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

• virtual array< array< uint32_t $>^{^{\wedge}}>$ ChannelBlock_ReadAsFrameArrayUI32 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

virtual array< array< uint32_t >^> ^ ChannelBlock_ReadAsFrameArrayUl32 (int handle, int queue, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

virtual array < array < int32_t >^> ^ ChannelBlock_ReadAsFrameArrayl32 (int handle, int frames, [System ← ::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

• virtual array< array< int32_t >^> ^ ChannelBlock_ReadAsFrameArrayl32 (int handle, int queue, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue as array of uint16_t data frame arrays

virtual System::Collections::Generic::Dictionary< int, array< uint16_t >^> ^ ChannelBlock_ReadFramesDictUI16 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames ret)

Read data from a FIFO queue in uint16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

virtual System::Collections::Generic::Dictionary< int, array< int16_t >^> ^ ChannelBlock_ReadFramesDictI16
 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)

Read data from a FIFO queue in int16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

- virtual System::Collections::Generic::Dictionary< int, array< uint32_t >^> ^ ChannelBlock_ReadFramesDictUI32 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames ret)
 - Read data from a FIFO queue in uint32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- virtual System::Collections::Generic::Dictionary< int, array< int32_t >^> ^ ChannelBlock_ReadFramesDictl32 (int handle, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
 - Read data from a FIFO queue in int32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- virtual System::Collections::Generic::Dictionary< int, array< uint16_t >^> ^ GetGroupChannelDataUI16 (DacqGroupChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
 - Read data from a FIFO queue in uint16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- virtual System::Collections::Generic::Dictionary< int, array< int16_t >^> ^ GetGroupChannelDatal16 (DacqGroupChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames ret)
 - Read data from a FIFO queue in int16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- virtual System::Collections::Generic::Dictionary< int, array< uint32_t >^> ^ GetGroupChannelDataUl32 (DacqGroupChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames_ret)
 - Read data from a FIFO queue in uint32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- virtual System::Collections::Generic::Dictionary< int, array< int32_t >^> ^ GetGroupChannelDatal32 (DacqGroupChannelEnumNet group, int frames, [System::Runtime::InteropServices::Out]int % frames ret)
 - Read data from a FIFO queue in int32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number
- void SetupGroupDacqQueue (int queuesize, int threshold)
- void SetupGroupDacqQueue (int queuesize, int threshold, unsigned int virtualDevice)
- CHWInfo ^ HWInfo ()

Static Public Attributes

- static const int Error_Callback_Queue_Full = 0x100
- static const int Error_Callback_Aquisition_Stopped = 0x200
- static const int Error_Callback_Packet_Error = 1
- static const int Error_Callback_RingQueue_Full = 3
- static const int Error_Callback_Frames_Lost = 4
- static const int Error_Callback_Data_lost = 5

Properties

• virtual int Samplerate [get, set]

The sampling frequency of the device in Hz.

Events

- OnChannelData^ ChannelDataEvent
- OnError[^] ErrorEvent

Additional Inherited Members

11.49.1 Detailed Description

Base class for data acquisition devices.

11.49.2 Constructor & Destructor Documentation

11.49.2.1 CMcsUsbDacqNet() CMcsUsbDacqNet ()

```
11.49.2.2 ~CMcsUsbDacqNet() ~CMcsUsbDacqNet ()
```

11.49.3 Member Function Documentation

11.49.3.1 AddSelectedChannelsQueue() [1/4] virtual int AddSelectedChannelsQueue (

```
int nByteOffset,
int nChannelOffset,
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize ) [virtual]
```

Adds a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_Read← FramesDict... with handle = 0 to read the data.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

Parameters

ByteOffset Number of bytes to start with.

Parameters

nChannelOffset	Number of channel to start with (counted in samplesize bytes).	

selectedChannels	List of channels to be collected in the FIFO.

queuesize	Size of sample frames the FIFO can hold.
-----------	--

Parameters

Parameters

samplesize size of the datawords, either 16 or 32b
--

Returns

The handle to the Queue.

$\textbf{11.49.3.2} \quad \textbf{AddSelectedChannelsQueue() [2/4]} \quad \text{virtual int AddSelectedChannelsQueue (} \\$

```
int nByteOffset,
int nChannelOffset,
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize ) [virtual]
```

11.49.3.3 AddSelectedChannelsQueue() [3/4] virtual int AddSelectedChannelsQueue (

```
int nByteOffset,
int nChannelOffset,
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize ) [virtual]
```

Adds a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_Read ← FramesDict... with handle = 0 to read the data.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

nByteOffset	Number of bytes to start with.
-------------	--------------------------------

Parameters

nChannelOffset Number of channel to start with (counted in samplesize

Parameters

nChannels	Number of channels to be collected in the FIFO.
-----------	---

Parameters

Parameters

threshold Number of sample frames the FIFO must acquire before the callback function is called
--

Parameters

	camplaciza	size of the datawords, either 16 or 32bit.	
ı	Samplesize	Size of the datawords, either to or szbit	

Returns

The handle to the Queue.

```
11.49.3.4 AddSelectedChannelsQueue() [4/4] virtual int AddSelectedChannelsQueue (
```

```
int nByteOffset,
int nChannelOffset,
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize ) [virtual]
```

```
11.49.3.5 ChannelBlock_AvailFrames() [1/2] virtual uint32_t ChannelBlock_AvailFrames ( int handle ) [virtual]
```

Get the number of sample frames already available in the FIFO.

Parameters

handle Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.

Returns

Number of sample frames available in the FIFO.

```
11.49.3.6 ChannelBlock_AvailFrames() [2/2] virtual uint32_t ChannelBlock_AvailFrames (
    int handle,
    int queue ) [virtual]
```

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

handle Handle of the FIFO queue. Zero when the SetSelectedData call was used.

frames

Parameters

frames_ret Number of sample frames which were read, might be
--

Returns

Array of int16_t frame arrays.

11.49.3.8 ChannelBlock_ReadAsFrameArrayl16() [2/2] virtual array<array<int16_t>^> ^ Channel \leftarrow

```
Block_ReadAsFrameArrayI16 (
    int handle,
    int queue,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

Parameters

queue	Number of the sub queue.
frames	Number of sample frames to read.

frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Array of int16_t frame arrays.

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

handle	Handle of the FIFO queue. Zero when the SetSelectedData call was used.
--------	--

Parameters

```
frames Number of sample frames to read.
```

Parameters

```
frames_ret | Number of sample frames which were read, might be smaller than frames.
```

Returns

Array of int32_t frame arrays.

[System::Runtime::InteropServices::Out] int % frames_ret) [virtual]

Read data from a FIFO queue as array of uint16_t data frame arrays

handle	Handle of the FIFO queue. Zero when the SetSelectedData call was used.
queue	Number of the sub queue.
frames	Number of sample frames to read.

Parameters

Returns

Array of int32_t frame arrays.

11.49.3.11 ChannelBlock_ReadAsFrameArrayUl16() [1/2] virtual array<array<uint16_t>^> ^

```
ChannelBlock_ReadAsFrameArrayUI16 (
    int handle,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

handle	Handle of the FIFO gueue. Zero when the SetSelectedData call was used.
--------	--

Parameters

Г	-	
П	framaa rat	Number of comple frames which were read might be smaller than frames
1	mannes rec	Number of sample frames which were read, might be smaller than frames.
П		,

Returns

Array of uint16_t frame arrays.

11.49.3.12 ChannelBlock_ReadAsFrameArrayUl16() [2/2] virtual array<array<uint16_t> $^{^{\wedge}}$ > $^{^{\wedge}}$

```
ChannelBlock_ReadAsFrameArrayUI16 (
    int handle,
    int queue,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

Parameters

queue	Number of the sub queue.
frames	Number of sample frames to read.

Parameters

frames_ret	Number of sample frames which were read, might be smaller than frames.
------------	--

Returns

Array of uint16_t frame arrays.

11.49.3.13 ChannelBlock_ReadAsFrameArrayUl32() [1/2] virtual array<array<uint32_t>^> ^

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

frames	Number of sample frames to read.
--------	----------------------------------

Parameters

trames ret Number of sample frames which were read, might be smaller than fra	frames ret	Number of sample frames which were read, might be smaller than frames.
---	------------	--

Returns

Array of uint32_t frame arrays.

11.49.3.14 ChannelBlock_ReadAsFrameArrayUl32() [2/2] virtual array<array<uint32_t>^> ^

```
ChannelBlock_ReadAsFrameArrayUI32 (
    int handle,
    int queue,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue as array of uint16_t data frame arrays

Parameters

FIFO queue. Zero when the SetSelectedData call was used.
--

queue	Number of the sub queue.
frames	Number of sample frames to read.

frames_ret	Number of sample frames which were read, might be smaller than frames.
------------	--

Returns

Array of uint32_t frame arrays.

Read data from a FIFO queue in int16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

	handle	Handle of the FIFO queue. Zero when the SetSelectedChannelsQueue call was used.	
--	--------	---	--

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of int16_t arrays and hardware channel as key.

Read data from a FIFO queue in int32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

handle	Handle of the FIFO queue. Zero when the SetSelectedChannelsQueue call was used.
--------	---

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of int32_t arrays and hardware channel as key.

Read data from a FIFO queue in uint16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

queue. Zero when the SetSelectedChannelsQueue call was used.	handle Handle of the FIFO queue. Zer
--	--------------------------------------

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of uint16_t arrays and hardware channel as key.

Read data from a FIFO queue in uint32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

handle	Handle of the FIFO queue. Zero when the SetSelectedChannelsQueue call was used.
--------	---

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of uint32_t arrays and hardware channel as key.

```
11.49.3.19 ChannelBlock_ReadFramesI16() [1/2] virtual uint32_t ChannelBlock_ReadFramesI16 (
    int handle,
    array< int16_t >^ buffer,
    int frames_pos,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in int16_t data format

Parameters

handle	Handle of the FIFO gueue. Either zero when the SetSelectedData call was used or the channel number.
--------	---

buffer	Buffer to put the data from the device in.
frames_pos	Position in buffer where to put the data.
frames	Number of sample frames to read.

frames_ret | Number of sample frames which were read, might be smaller than frames.

Returns

Error Status. 0 on success.

Read data from a FIFO queue in int16_t data format

Parameters

handle Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.

Parameters

frames Number of sample frames to read.

Parameters

frames_ret | Number of sample frames which were read, might be smaller than frames.

```
11.49.3.21 ChannelBlock_ReadFramesI32() [1/2] virtual uint32_t ChannelBlock_ReadFramesI32 (
    int handle,
    array< int32_t >^ buffer,
    int frames_pos,
```

```
int frames,
[System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in uint32_t data format

Parameters

handle	Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.
--------	---

Parameters

buffer	Buffer to put the data from the device in.
frames_pos	Position in buffer where to put the data.
frames	Number of sample frames to read.

Parameters

Returns

Error Status. 0 on success.

```
11.49.3.22 ChannelBlock_ReadFramesl32() [2/2] virtual array<int32_t> ^ ChannelBlock_Read←

FramesI32 (

int handle,
int frames,

[System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in uint32_t data format

Parameters

handle Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.

frames	Number of sample frames to read.
--------	----------------------------------

Parameters

frames ret	Number of sample frames which were read, might be smaller than frames.

```
11.49.3.23 ChannelBlock_ReadFramesUl16() [1/2] virtual uint32_t ChannelBlock_ReadFramesUI16 (
    int handle,
    array< uint16_t >^ buffer,
    int frames_pos,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in uint16_t data format

Parameters

handle	Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel	
	number.	
buffer	Buffer to put the data from the device in.	
frames_pos	S Position in buffer where to put the data.	
frames	Number of sample frames to read.	

Parameters

fram	es_ret	Number of sample frames which were read, might be smaller than frames.
------	--------	--

Returns

Error Status. 0 on success.

```
11.49.3.24 ChannelBlock_ReadFramesUl16() [2/2] virtual array<uint16_t> ^ ChannelBlock_Read←
FramesUI16 (
    int handle,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in uint16_t data format

	handle	Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.	
--	--------	---	--

Parameters

frames	Number of sample frames to read.
--------	----------------------------------

Parameters

trames ret Number of sample frames which were read, might be smaller than fra	frames ret	Number of sample frames which were read, might be smaller than frames.
---	------------	--

Returns

Array of data from the device.

```
11.49.3.25 ChannelBlock_ReadFramesUl32() [1/2] virtual uint32_t ChannelBlock_ReadFramesUl32 (
    int handle,
    array< uint32_t >^ buffer,
    int frames_pos,
    int frames,
    [System::Runtime::InteropServices::Out] int % frames_ret ) [virtual]
```

Read data from a FIFO queue in uint32_t data format

Parameters

```
handle Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.
```

buffer	Buffer to put the data from the device in.
frames_pos	Position in buffer where to put the data.
frames	Number of sample frames to read.

frames_ret Number of sample frames which were read, might be smaller than frames.

Returns

Error Status. 0 on success.

Read data from a FIFO queue in uint32_t data format

Parameters

handle | Handle of the FIFO queue. Either zero when the SetSelectedData call was used or the channel number.

Parameters

frames Number of sample frames to read.

Parameters

frames_ret Number of sample frames which were read, might be smaller than frames.

11.49.3.27 ClearBuffers() virtual void ClearBuffers () [virtual]

```
11.49.3.28 CMcsUsbDacqNet::GetFilterProperties() virtual array<CFilterPropertyNet^> ^ CMcs↔ UsbDacqNet::GetFilterProperties (

DacqGroupChannelEnumNet GroupID ) [virtual]
```

```
11.49.3.29 GetAdapterType() virtual AdapterTypeEnumNet GetAdapterType ( ) [virtual]
```

Gets the adapter which is connected to the MEA2100 device.

Returns

AdapterTypeEnumNet which enumerates the possible adapters.

```
11.49.3.30 GetAdcDataFormat() virtual uint32_t GetAdcDataFormat ( uint32_t virtualDevice ) [virtual]
```

Gets the ADC data format, 16 means 16 bits, 24 means 24 bits, 32 means 32 bits.

Returns

The data format in bits.

11.49.3.33 GetChannelDataFillSize() int GetChannelDataFillSize ()

```
11.49.3.35 GetChannelsInBlock() virtual int32_t GetChannelsInBlock ( unsigned int virtualDevice ) [virtual]
```

Get the number of 16 bit datawords which will be collected per sample frame, use after the device is configured.

Returns

Number of 16 bit datawords per sample frame.

```
11.49.3.37 GetDataMode() virtual DataModeEnumNet GetDataMode (
unsigned int virtualDevice ) [virtual]
```

Gets the data mode, can be 16, 24 or 32bit, all signed or unsigned on the MEA2100 device.

Parameters

```
virtualDevice Virtual device to use.
```

Returns

DataModeEnumNet which enumerates the possible data modes.

Gets the function/source of an digital output bit.

This is the templated generic implementation.

digitaltarget	The digital target to query.
NrChannel	The channel/bit of target to query.
source	The source/function assignd to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Gets the function/source of an digital output bit.

This overload is for the MEA2100 device.

Parameters

digitaltarget	The digital target to query.
NrChannel	The channel/bit of target to query.
source	The source/function assignd to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Gets the function/source of an digital output bit.

This overload is for the MEA2100-256 device.

digitaltarget	The digital target to query.
NrChannel	The channel/bit of target to query.
source	The source/function assignd to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

```
11.49.3.41 GetDigitalSource() [4/5] void GetDigitalSource (
DigitalTargetEnumNet digitaltarget,
```

```
int32_t NrChannel,
[System::Runtime::InteropServices::Out] SCUDigitalSourceEnumNet% source,
[System::Runtime::InteropServices::Out] int% bitnumber_offset )
```

Gets the function/source of an digital output bit.

This overload is for the SCU device.

Parameters

digitaltarget	The digital target to query.
NrChannel	The channel/bit of target to query.
source	The source/function assignd to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Gets the function/source of an digital output bit.

This overload is for the W2100 device.

Parameters

digitaltarget	The digital target to query.
NrChannel	The channel/bit of target to query.
source	The source/function assignd to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

```
11.49.3.43 GetFilterProperty() virtual CFilterPropertyNet ^ GetFilterProperty (
DacqGroupChannelEnumNet GroupID,
unsigned int index ) [virtual]
```

Read data from a FIFO queue in int16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

group	Group selector supported by the device.
-------	---

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of int16_t arrays and hardware channel as key.

Read data from a FIFO queue in int32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

group G	Group selector supported by the device.
---------	---

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of int32_t arrays and hardware channel as key.

Read data from a FIFO queue in uint16_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

Parameters

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of uint16_t arrays and hardware channel as key.

Read data from a FIFO queue in uint32_t data format, that contains subqueues, each populates an entry in the dictionary by hardware channel number

Parameters

ſ	group	Group selector supported by the device.

frames	Number of sample frames to read.
frames_ret	Number of sample frames which were read, might be smaller than frames.

Returns

Dictonary of uint32_t arrays and hardware channel as key.

```
11.49.3.50 GetMaxSamplingFrequency() virtual uint32_t GetMaxSamplingFrequency ( int virtualDevice ) [virtual]
```

Gets the maximal sampling frequency of the device.

Returns

Sampling frequency in Hz.

```
11.49.3.51 GetMeaLayout() virtual MeaLayoutEnumNet GetMeaLayout ( ) [virtual]
```

Gets the MEA layout which is connected to the MEA2100 device.

Returns

MeaLayoutEnumNet which enumerates the MEA types.

```
11.49.3.52 GetMinSamplingFrequencyStepsize() virtual uint32_t GetMinSamplingFrequencyStepsize ( ) [virtual]
```

Gets the minimal sampling frequency step size increment value of the device.

Returns

Sampling frequency step size in Hz.

Get the real number of data bits.

This value may be different from the value returned by GetDataFormat, e.g. in MC_Card the data are shifted 2 bits so the real number is 14 while the data format is 16 bits

Gets the sampling frequency of the device.

Returns

Sampling frequency in Hz.

```
11.49.3.57 GetVoltageRangeIndex() virtual uint32_t GetVoltageRangeIndex ( unsigned int virtualDevice ) [virtual]
```

Gets the currently selected voltage range on devices which support multiple voltage ranges.

Returns

The Voltage Range in uV.

```
11.49.3.59 GetVoltageRangeInMilliVolt() virtual int32_t GetVoltageRangeInMilliVolt () [virtual]
```

Gets the currently selected voltage range on devices which support multiple voltage ranges.

Returns

The rounded Voltage Range in mV.

```
11.49.3.60 HWInfo() CHWInfo ^ HWInfo ()
```

```
11.49.3.61 SendStartDacq() [1/2] virtual void SendStartDacq ( ) [virtual]
```

Start sampling.

```
11.49.3.62 SendStartDacq() [2/2] virtual void SendStartDacq ( int VirtualDacqMap ) [virtual]
```

Start sampling.

Parameters

VirtualDacqMap

Start sampling together with the STG.

trigger_map
VirtualDacqMap

11.49.3.64 SendStopDacq() [1/2] virtual void SendStopDacq () [virtual]

Stop sampling.

```
11.49.3.65 SendStopDacq() [2/2] virtual void SendStopDacq ( int VirtualDacqMap ) [virtual]
```

Stop sampling.

Parameters

VirtualDacqMap

Stop sampling together with the STG.

Parameters

trigger_map

Stop sampling together with the STG and options.

trigger_map

Parameters

options

Parameters

VirtualDacqMap

```
11.49.3.68 SetDataMode() virtual void SetDataMode (

DataModeEnumNet dataMode,

unsigned int virtualDevice ) [virtual]
```

Sets the data mode, can be 16, 24 or 32bit, all signed or unsigned on the MEA2100 device.

Parameters

dataMode	DataModeEnumNet enumerates the possible data modes.
virtualDevice	Virtual device to use.

Sets the function/source of an digital output bit.

This is the templated generic implementation.

digitaltarget	The digital target to change.
NrChannel	The channel/bit of target to change.
source	The source/function to assign to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Sets the function/source of an digital output bit.

This overload is for the MEA2100 device.

Parameters

digitaltarget	The digital target to change.
NrChannel	The channel/bit of target to change.
source	The source/function to assign to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Sets the function/source of an digital output bit.

This overload is for the MEA2100-256 device.

Parameters

digitaltarget	The digital target to change.
NrChannel	The channel/bit of target to change.
source	The source/function to assign to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Sets the function/source of an digital output bit.

This overload is for the SCU device.

digitaltarget	The digital target to change.
NrChannel	The channel/bit of target to change.
source	The source/function to assign to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Sets the function/source of an digital output bit.

This overload is for the W2100 device.

Parameters

digitaltarget	The digital target to change.
NrChannel	The channel/bit of target to change.
source	The source/function to assign to the digital target.
bitnumber_offset	An offset / bit number with the source/function.

Sets the sampling frequency of the device.

Parameters

```
rate Sampling frequency in Hz.
```

11.49.3.76 SetSelectedChannels() [1/4] virtual void SetSelectedChannels (

```
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock ) [virtual]
```

Create a FIFO queue per channel. Each channel will have its own FIFO and Callback function.

When using a 32bit sample size, the number obtained from GetChannelsInBlock must be devided by 2 to be used here, since GetChannelsInBlock returns the number of 16 bit datapoints per sample frame, while this functions uses the number of sample frames in its own data format.

Parameters

selectedChannels	List of channels to be collected in the FIFO.
------------------	---

Parameters

queuesize	Size of sample frames the FIFO can hold.
-----------	--

Parameters

Parameters

Parameters

```
ChannelsInBlock value obtained from GetChannelsInBlock.
```

11.49.3.77 SetSelectedChannels() [2/4] virtual void SetSelectedChannels (

```
array<br/><br/> bool >^{\wedge} selectedChannels,<br/>int queuesize,<br/>int threshold,<br/>SampleSizeNet samplesize,
```

```
SampleDstSizeNet sampleDstSize,
int ChannelsInBlock ) [virtual]
```

11.49.3.78 SetSelectedChannels() [3/4] virtual void SetSelectedChannels (

```
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock ) [virtual]
```

Create a FIFO queue per channel. Each channel will have its own FIFO and Callback function.

When using a 32bit sample size, the number obtained from GetChannelsInBlock must be devided by 2 to be used here, since GetChannelsInBlock returns the number of 16 bit datapoints per sample frame, while this functions uses the number of sample frames in its own data format.

Parameters

nChannels 1	Number of channels to be collected in the FIFO.
-------------	---

Parameters

Parameters

Parameters

samplesize size of the		size of the datawords, either 16 or 32bit.
	ChannelsInBlock	value obtained from GetChannelsInBlock.

11.49.3.79 SetSelectedChannels() [4/4] virtual void SetSelectedChannels (

```
int nChannels,
int queuesize,
int threshold,
```

```
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize,
int ChannelsInBlock ) [virtual]
```

11.49.3.80 SetSelectedChannelsQueue() [1/4] virtual void SetSelectedChannelsQueue (

```
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock) [virtual]
```

Create a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_← ReadFramesDict... with handle = 0 to read the data.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

Parameters

selectedChannels	List of channels to be collected in the FIFO.
------------------	---

Parameters

Parameters

threshold	Number of sample frames the FIFO must acquire before the callback function is called.

Parameters

ChannelsInBlock	value obtained from GetChannelsInBlock.
-----------------	---

11.49.3.81 SetSelectedChannelsQueue() [2/4] virtual void SetSelectedChannelsQueue (

```
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize,
int ChannelsInBlock ) [virtual]
```

11.49.3.82 SetSelectedChannelsQueue() [3/4] virtual void SetSelectedChannelsQueue (

```
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock ) [virtual]
```

Create a common FIFO queue for all channels. Data in callback will be a list per channel. Use ChannelBlock_← ReadFramesDict... with handle = 0 to read the data.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

Parameters

nChannole	Number of channels to be collected in the FIFO.
HUHAHHUUS	Number of charmers to be confected in the FIFO.

Parameters

queuesize Size of sample frames the FIFO	can hold.
--	-----------

Parameters

threshold	Number of sample frames the FIFO must acquire before the callback function is called.
-----------	---

aamplaaiza	size of the datawords.	oithar 16 ar 20hit
Samuesize	Size of the datawords.	enner to or azon.

ChannelsInBlock	value obtained from GetChannelsInBlock.
-----------------	---

11.49.3.83 SetSelectedChannelsQueue() [4/4] virtual void SetSelectedChannelsQueue (

```
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize,
int ChannelsInBlock ) [virtual]
```

11.49.3.84 SetSelectedData() [1/4] virtual void SetSelectedData (

```
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock ) [virtual]
```

Create a common FIFO queue for all channels. Use handle = 0 in the ChannelBlock_ReadFrames... functions.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

Parameters

selectedChannels List of channe	Is to be collected in the FIFO.
---------------------------------	---------------------------------

Parameters

	queuesize	Size of sample frames the FIFO can hold.
--	-----------	--

لما ممام مسطة	Number of sample frames the FIFO must acquire before the callback function is called
Inresnoia	Number of sample trames the FIFO must acquire before the caliback function is called.

samplesize	size of the datawords, either 16 or 32bit.
ChannelsInBlock	value obtained from GetChannelsInBlock.

11.49.3.85 SetSelectedData() [2/4] virtual void SetSelectedData (

```
array< bool >^ selectedChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
SampleDstSizeNet sampleDstSize,
int ChannelsInBlock ) [virtual]
```

11.49.3.86 SetSelectedData() [3/4] virtual void SetSelectedData (

```
int nChannels,
int queuesize,
int threshold,
SampleSizeNet samplesize,
int ChannelsInBlock ) [virtual]
```

Create a common FIFO queue for all channels. Use handle = 0 in the ChannelBlock_ReadFrames... functions.

When using 32 bit data format, ChannelsInBlock is still the number of 16 bit channels per frame, as obtained from GetChannelsInBlock, while nChannels is the number of 32 bit channels to be read from the device. So when all channels from a device are read in 32 bit data format nChannels = ChannelsInBlock/2

Parameters

	nChannels	Number of channels to be collected in the FIFO.
--	-----------	---

Parameters

queuesize Size of sample frames the FIFO can hold.
--

threshold	Number of sample frames the FIFO must acquire before the callback function is called.
unicanola	I Number of sample frames the Fit O must acquire before the camback function is called.

samplesize	size of the datawords, either 16 or 32bit.
ChannelsInBlock	value obtained from GetChannelsInBlock.

```
11.49.3.87 SetSelectedData() [4/4] virtual void SetSelectedData (
              int nChannels,
              int queuesize,
              int threshold,
              SampleSizeNet samplesize,
              SampleDstSizeNet sampleDstSize,
              int ChannelsInBlock ) [virtual]
\textbf{11.49.3.88} \quad \textbf{SetupGroupDacqQueue() [1/2]} \quad \texttt{void SetupGroupDacqQueue} \quad \textbf{(}
              int queuesize,
              int threshold )
11.49.3.89 SetupGroupDacqQueue() [2/2] void SetupGroupDacqQueue (
              int queuesize,
              int threshold,
              unsigned int virtualDevice )
11.49.3.90 SetVoltageRangeByIndex() virtual void SetVoltageRangeByIndex (
              int32_t voltageRangeIndex,
```

Sets the voltage range on devices which support multiple voltage ranges.

unsigned int virtualDevice) [virtual]

Parameters

```
voltageRangeIndex Voltage Range to use as index, smaller values are larger voltage ranges.
```

Sets the voltage range on devices which support multiple voltage ranges.

This replaces SetVoltageRange, where the value of the range was in mV!

```
11.49.3.92 StartDacq() [1/4] virtual void StartDacq ( ) [virtual]
```

Start the data acquisition thread and sampling.

```
11.49.3.93 StartDacq() [2/4] virtual void StartDacq ( int32_t timeout ) [virtual]
```

Start the data acquisition thread and sampling.

Parameters

```
timeout Timeout in ms.
```

Start the data acquisition thread and sampling.

Parameters

```
timeout Timeout in ms.
```

Parameters

Start the data acquisition thread and sampling.

Parameters

numSubmittedUsbBuffers	Number of USB Buffers that are simultaniously	submitted.

Parameters

Parameters

numl lab Duffara	Number of LICE Buffers to use
numuspbutters	Number of USB Buffers to use.

Parameters

packetsInUrb	Packets in each URB.
packetsiiioib	rackets ill each onb.

virtualDevice	Virtual Device to start.
vii luaiDevice	VII luai Device lu siai i.

```
11.49.3.96 StartLoop() [1/4] virtual void StartLoop ( ) [virtual]
```

Start the data acquisition thread.

Start the data acquisition thread.

Parameters

```
timeout Timeout in ms.
```

Start the data acquisition thread.

Parameters

Parameters

11.15.00	N I (HODD W I
niimi ishkiitters	Number of LISB Buffers to use

packetsInUrb	Packets in each URB.
--------------	----------------------

11.49.3.99 StartLoop() [4/4] virtual void StartLoop (

```
int32_t timeout,
int32_t numSubmittedUsbBuffers,
int32_t numUsbBuffers,
int32_t packetsInUrb,
uint32_t virtualDevice ) [virtual]
```

Start the data acquisition thread.

Parameters

numSubmittedUsbBuffers Number of USB Buffers that are simultanion

Parameters

timeout	Timeout in ms.
unicour	HIIIIGOULIII IIIG.

Parameters

numl IshBuffers	Number of USB Buffers to use.

Parameters

```
11.49.3.100 StopDacq() [1/2] virtual void StopDacq ( ) [virtual]
```

Stop the data acquisition thread and sampling.

```
11.49.3.101 StopDacq() [2/2] virtual void StopDacq ( uint32_t virtualDevice ) [virtual]
```

Stop the data acquisition thread and sampling.

Parameters

virtualDevice Virtual Device to start

```
11.49.3.102 StopLoop() virtual void StopLoop ( ) [virtual]
```

11.49.4 Member Data Documentation

11.49.4.1 Error_Callback_Aquisition_Stopped const int Error_Callback_Aquisition_Stopped = 0x200 [static]

```
11.49.4.2 Error_Callback_Data_lost const int Error_Callback_Data_lost = 5 [static]
```

11.49.4.3 Error_Callback_Frames_Lost const int Error_Callback_Frames_Lost = 4 [static]

11.49.4.4 Error_Callback_Packet_Error const int Error_Callback_Packet_Error = 1 [static]

11.49.4.5 Error_Callback_Queue_Full const int Error_Callback_Queue_Full = 0x100 [static]

11.49.4.6 Error_Callback_RingQueue_Full const int Error_Callback_RingQueue_Full = 3 [static]

11.49.5 Property Documentation

11.49.5.1 Samplerate virtual int Samplerate [get], [set]

The sampling frequency of the device in Hz.


11.49.6 Event Documentation

11.49.6.1 ChannelDataEvent OnChannelData^ ChannelDataEvent

11.49.6.2 ErrorEvent OnError^ ErrorEvent

11.50 CMcsUsbDeviceStatePushFunctionNet Class Reference

Inheritance diagram for CMcsUsbDeviceStatePushFunctionNet:



Public Member Functions

• void TriggerStatus ()

Protected Member Functions

CMcsUsbDeviceStatePushFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] p
 Device)

Events

 $\bullet \ \ OnMcsUsbDeviceState ^{\wedge} \ McsUsbDeviceStateEvent \ \ [add, remove, raise]\\$

Additional Inherited Members

11.50.1 Constructor & Destructor Documentation

11.50.1.1 CMcsUsbDeviceStatePushFunctionNet() CMcsUsbDeviceStatePushFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pDevice) [protected]

11.50.2 Member Function Documentation


```
11.50.2.1 TriggerStatus() void TriggerStatus ()
```

11.50.3 Event Documentation

11.50.3.1 McsUsbDeviceStateEvent OnMcsUsbDeviceState^ McsUsbDeviceStateEvent [add], [remove], [raise]

11.51 CMcsUsbDeviceStatePushNet Class Reference

 $Inheritance\ diagram\ for\ CMcsUsbDeviceStatePushNet:$



Public Member Functions

• void TriggerStatus ()

Protected Member Functions

CMcsUsbDeviceStatePushNet (CMcsUsbPointerContainer[^] pDevice)

Events

• OnMcsUsbDeviceState^ McsUsbDeviceStateEvent [add, remove, raise]

Additional Inherited Members

11.51.1 Constructor & Destructor Documentation

```
11.51.1.1 CMcsUsbDeviceStatePushNet() CMcsUsbDeviceStatePushNet (
CMcsUsbPointerContainer^ pDevice) [protected]
```

11.51.2 Member Function Documentation


```
11.51.2.1 TriggerStatus() void TriggerStatus ()
```

11.51.3 Event Documentation

```
11.51.3.1 McsUsbDeviceStateEvent OnMcsUsbDeviceState^ McsUsbDeviceStateEvent [add], [remove], [raise]
```

11.52 CMcsUsbFactoryNet Class Reference

Inheritance diagram for CMcsUsbFactoryNet:



Public Member Functions

- CMcsUsbFactoryNet ()
- ∼CMcsUsbFactoryNet ()
- unsigned int GetNumDestinations ()
- String \(^{\text{ GetDestinationName}}\) (unsigned int index)
- String \(^{\text{GetDestinationName}}\) (CFirmwareDestinationNet dest)
- void SetDestinationSerialNumber (CFirmwareDestinationNet dest, String[^] serialnumber)
- String \(^\) GetDestinationSerialNumber (CFirmwareDestinationNet dest)
- CFirmwareDestinationNet GetDestination (unsigned int index)
- CFirmwareDestinationNet GetDestination (String[^] Key)
- unsigned int GetDestinationTargetAddress (CFirmwareDestinationNet destination)

Gets the target base address for the destination.

- uint32_t ChangeSerialNumber (String^ serial)
- bool LoadUserFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry)

Send the DSP Firmware to the MEA21 device.

- bool LoadUserFirmware (String^ FirmwareFile, CMcsUsbListEntryNet^ listEntry, uint32 t LockMask)
- bool UpdateFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestinationNet Dest, OnUpdateFirmwareStatusChange[^] deleg, OnUpdateFirmwareProgress[^] progress, bool SkipWait)

Flashes a firmware file to the device.

- bool UpdateFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestinationNet Dest, OnUpdateFirmwareStatusChange[^] deleg, OnUpdateFirmwareProgress[^] progress, bool SkipWait, unsigned int LockMask)
- bool UpdateFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestinationNet dest)

Flashes a firmware file to the device.

 bool UpdateFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestinationNet dest, bool SkipWait)

Flashes a firmware file to the device.

- bool UpdateFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestinationNet dest, bool SkipWait, uint32_t LockMask)
- bool CompareFirmware (String[^] FirmwareFile, CMcsUsbListEntryNet[^] listEntry, CFirmwareDestination
 — Net Dest, OnUpdateFirmwareStatusChange[^] deleg, OnUpdateFirmwareProgress[^] progress, String[^]
 MessagePrefix, unsigned int LockMask, [System::Runtime::InteropServices::Out] String[^]% ErrorText,
 [System::Runtime::InteropServices::Out] String[^]% Protokoll)
- uint32_t Coldstart (CFirmwareDestinationNet dest)
- int32 t GetXilinxFlashOffset (CFirmwareDestinationNet dest)
- uint32 t GetXilinxFlashReadCommand (CFirmwareDestinationNet dest)
- array< uint8 t > ^ DownloadFirmware (CFirmwareDestinationNet Dest, uint32 t Address, uint32 t length)
- bool GetUsercodeFromFlash (unsigned int FPGA, unsigned int Address, [System::Runtime::Interop
 — Services::Out] unsigned int% Usercode)
- array< unsigned char > ^ ReadBlockFromFlash (unsigned int FPGA, unsigned int Address)
- void ReadBlockFromFlash (unsigned int FPGA, unsigned int Address, array< unsigned char >^ buffer, int position)
- array< unsigned char > ^ ReadBlockFromIFBGlobalEEprom (unsigned int Address)
- array< unsigned char > ^ ReadBlockFromNVMEM (unsigned int FPGA, unsigned int Offset, unsigned int Address)

Static Public Member Functions

- static String \(^\) GetDestinationDisplayLabel (String\(^\) RawLabel, CFirmwareDestinationNet dest)
- static String ^ FindFirmwareVersionMagicInBuffer (array< unsigned char >^ buffer, int length, [System::

 Runtime::InteropServices::Out]int% position)
- static bool GetFirmwareVersionFromFile (String[^] FirmwareFile, [System::Runtime::InteropServices::Out] uint32 t% Version)

Retrives version info from a Firmware update file.

- static bool GetFirmwareVersionFromFile (String^ FirmwareFile, [System::Runtime::InteropServices::Out] uint32_t% Version, [System::Runtime::InteropServices::Out] uint32_t% Position)
- static bool GetFirmwareVersionFromHexFile (String[^] FirmwareFile, [System::Runtime::InteropServices::Out] uint32_t% Version)
- static uint32_t GetChecksumFromFX3Image (String[^] FirmwareFile)
- static uint32_t GetUSBDeviceIDFromFX3Image (String[^] FirmwareFile)
- static bool GetUsercodeFromBitFile (String[^] FirmwareFile, [System::Runtime::InteropServices::Out] unsigned int% Usercode)

Static Public Attributes

- static const uint32 t FX3MCSDataAddress = 0x40037E00
- static const uint32_t FX3MCSDataDeviceIdOffset = 0x4
- static const uint32_t FX3MCSDataVersionOffset = 0x8
- static const uint32 t FX3MCSDataIFB2ImageOffset = 0xC
- static const uint32_t FX3MCSDataIFB1ImageOffset = 0x2C

Additional Inherited Members

11.52.1 Constructor & Destructor Documentation

```
11.52.1.1 CMcsUsbFactoryNet() CMcsUsbFactoryNet ( )
```

11.52.1.2 \sim CMcsUsbFactoryNet() \sim CMcsUsbFactoryNet ()

11.52.2 Member Function Documentation

```
11.52.2.1 ChangeSerialNumber() uint32_t ChangeSerialNumber ( String^{\wedge} serial)
```

```
11.52.2.2 Coldstart() uint32_t Coldstart (
              CFirmwareDestinationNet dest)
11.52.2.3 CompareFirmware() bool CompareFirmware (
              String^ FirmwareFile,
              CMcsUsbListEntryNet^ listEntry,
              CFirmwareDestinationNet Dest,
              OnUpdateFirmwareStatusChange^ deleg,
              OnUpdateFirmwareProgress^ progress,
              String^{\wedge} MessagePrefix,
              unsigned int LockMask,
               [System::Runtime::InteropServices::Out] String^% ErrorText,
               [System::Runtime::InteropServices::Out] String^% Protokoll )
11.52.2.4 DownloadFirmware() array<uint8_t> ^ DownloadFirmware (
              CFirmwareDestinationNet Dest,
              uint32_t Address,
              uint32_t length )
11.52.2.5 FindFirmwareVersionMagicInBuffer() static String ^ FindFirmwareVersionMagicInBuffer (
              array< unsigned char >^{\wedge} buffer,
              int length,
               [System::Runtime::InteropServices::Out] int% position ) [static]
11.52.2.6 GetChecksumFromFX3Image() static uint32_t GetChecksumFromFX3Image (
              String<sup>∧</sup> FirmwareFile ) [static]
\textbf{11.52.2.7} \quad \textbf{GetDestination()} \; \texttt{[1/2]} \quad \texttt{CFirmwareDestinationNet GetDestination} \; (
              String^{\wedge} Key )
11.52.2.8 GetDestination() [2/2] CFirmwareDestinationNet GetDestination (
              unsigned int index)
\textbf{11.52.2.9} \quad \textbf{GetDestinationDisplayLabel()} \quad \texttt{static String} \; \land \; \texttt{GetDestinationDisplayLabel} \; \; \textbf{(}
              String^ RawLabel,
              CFirmwareDestinationNet dest ) [static]
```

```
11.52.2.10 GetDestinationName() [1/2] String ^ GetDestinationName (
             CFirmwareDestinationNet dest )
11.52.2.11 GetDestinationName() [2/2] String ^ GetDestinationName (
             unsigned int index )
11.52.2.12 GetDestinationSerialNumber() String ^ GetDestinationSerialNumber (
             CFirmwareDestinationNet dest)
11.52.2.13 GetDestinationTargetAddress() unsigned int GetDestinationTargetAddress (
             CFirmwareDestinationNet destination )
Gets the target base address for the destination.
Parameters
 destination
             The destination to be queried.
Returns
     The base address as a 32 bit number, only the lower 16 bit represent the address.
11.52.2.14 GetFirmwareVersionFromFile() [1/2] static bool GetFirmwareVersionFromFile (
             String^ FirmwareFile,
             [System::Runtime::InteropServices::Out] uint32_t% Version ) [static]
Retrives version info from a Firmware update file.
11.52.2.15 GetFirmwareVersionFromFile() [2/2] static bool GetFirmwareVersionFromFile (
             String Firmware File,
             [System::Runtime::InteropServices::Out] uint32_t% Version,
             [System::Runtime::InteropServices::Out] uint32_t% Position ) [static]
11.52.2.16 GetFirmwareVersionFromHexFile() static bool GetFirmwareVersionFromHexFile (
             String^ FirmwareFile,
```

[System::Runtime::InteropServices::Out] uint32_t% Version) [static]

```
11.52.2.17 GetNumDestinations() unsigned int GetNumDestinations ( )
\textbf{11.52.2.18} \quad \textbf{GetUSBDeviceIDFromFX3Image()} \quad \texttt{static uint32\_t GetUSBDeviceIDFromFX3Image ()} \\
               String FirmwareFile ) [static]
11.52.2.19 GetUsercodeFromBitFile() static bool GetUsercodeFromBitFile (
               String Firmware File,
               [{\tt System::Runtime::InteropServices::Out}] \ \ unsigned \ \ int \ \ \textit{Usercode} \ ) \ \ \ [{\tt static}]
\textbf{11.52.2.20} \quad \textbf{GetUsercodeFromFlash()} \quad \texttt{bool GetUsercodeFromFlash ()}
               unsigned int FPGA,
               unsigned int Address,
               [System::Runtime::InteropServices::Out] \ unsigned \ int \% \ \textit{Usercode} \ )
11.52.2.21 GetXilinxFlashOffset() int32_t GetXilinxFlashOffset (
               CFirmwareDestinationNet dest )
11.52.2.22 GetXilinxFlashReadCommand() uint32_t GetXilinxFlashReadCommand (
               CFirmwareDestinationNet dest )
11.52.2.23 LoadUserFirmware() [1/2] bool LoadUserFirmware (
               String^{\wedge} FirmwareFile,
               CMcsUsbListEntryNet^ listEntry )
Send the DSP Firmware to the MEA21 device.
Parameters
 FirmwareFile
                 Filename of the DSP Firmware (*.bin) file.
```

listEntry Device to use for the connection. See CMcsUsbListNet.

```
11.52.2.24 LoadUserFirmware() [2/2] bool LoadUserFirmware (
             String^ FirmwareFile,
             CMcsUsbListEntryNet^ listEntry,
             uint32_t LockMask )
11.52.2.25 ReadBlockFromFlash() [1/2] array<unsigned char> ^ ReadBlockFromFlash (
             unsigned int FPGA,
             unsigned int Address )
11.52.2.26 ReadBlockFromFlash() [2/2] void ReadBlockFromFlash (
             unsigned int FPGA,
             unsigned int Address,
             array< unsigned char >^{\wedge} buffer,
             int position )
11.52.2.27 ReadBlockFromIFBGlobalEEprom() array<unsigned char> ^ ReadBlockFromIFBGlobalE←
Eprom (
             unsigned int Address )
11.52.2.28 ReadBlockFromNVMEM() array<unsigned char> ^ ReadBlockFromNVMEM (
             unsigned int FPGA,
             unsigned int Offset,
             unsigned int Address )
11.52.2.29 SetDestinationSerialNumber() void SetDestinationSerialNumber (
             CFirmwareDestinationNet dest,
             String^ serialnumber )
11.52.2.30 UpdateFirmware() [1/5] bool UpdateFirmware (
             String^ FirmwareFile,
             CMcsUsbListEntryNet<sup>^</sup> listEntry,
             CFirmwareDestinationNet dest )
```

Flashes a firmware file to the device.

FirmwareFile	Filename of the Firmware file.	
--------------	--------------------------------	--

Parameters

```
listEntry Device to use for the connection.
```

Flashes a firmware file to the device.

Parameters

Parameters

```
listEntry Device to use for the connection.
```

bool SkipWait,
uint32_t LockMask)

```
11.52.2.33 UpdateFirmware() [4/5] bool UpdateFirmware (
```

```
String FirmwareFile,
CMcsUsbListEntryNet listEntry,
CFirmwareDestinationNet Dest,
```

```
OnUpdateFirmwareStatusChange^ deleg,
OnUpdateFirmwareProgress^ progress,
bool SkipWait )
```

Flashes a firmware file to the device.

Parameters

```
FirmwareFile Filename of the Firmware file.
```

```
11.52.2.34 UpdateFirmware() [5/5] bool UpdateFirmware (
String^ FirmwareFile,
CMcsUsbListEntryNet^ listEntry,
```

```
CMcsUsbListEntryNet^ listEntry,
CFirmwareDestinationNet Dest,
OnUpdateFirmwareStatusChange^ deleg,
OnUpdateFirmwareProgress^ progress,
bool SkipWait,
unsigned int LockMask)
```

11.52.3 Member Data Documentation

```
11.52.3.1 FX3MCSDataAddress const uint32_t FX3MCSDataAddress = 0x40037E00 [static]
```

- 11.52.3.2 FX3MCSDataDeviceIdOffset const uint32_t FX3MCSDataDeviceIdOffset = 0x4 [static]
- 11.52.3.3 FX3MCSDataIFB1ImageOffset const uint32_t FX3MCSDataIFB1ImageOffset = 0x2C [static]
- 11.52.3.4 FX3MCSDataIFB2ImageOffset const uint32_t FX3MCSDataIFB2ImageOffset = 0xC [static]
- 11.52.3.5 FX3MCSDataVersionOffset const uint32_t FX3MCSDataVersionOffset = 0x8 [static]

11.53 CMcsUsbFunctionNet Class Reference

Inheritance diagram for CMcsUsbFunctionNet:

CMcsUshFunctionNet	
	CDecqGroupChannelSelectionTemphiteNet's DecqGroupChannelEtzunNet, DecqGroupChannelEtzun, CDeviceGroupChannelInfoNet >
	CDocyGroup Channel Selection TempiatoNet visr, int., CDovice Group Channel InfeGenericNet >
	·
	CDacqGroupChannelSelectionTemplateNets MEA2100_256DacqGroupChannelEnumNet, MEA2100_256DacqGroupChannelEnum, CDeviceGroupChannelInfoMEA2100_256DacqGroupChan
	CDacqGroupChannelSelectionTemplateNet+SCUDacqGroupChannelBraunNet,SCUDacqGroupChannelBraun,CDeviceGroupChannelBnfoSCUNet+
	CDacqGroupChannelSelectionTemplateNets W2100DacqGroupChannelEzuanNet, W2100DacqGroupChannelEnam, CDeviceGroupChannelInfoW2100Net >
	CCMOSMea_FunctionNet
	CDucCalibrationFunctionNet
	CDacqGroupChannelSelectiosTemplateNet+ DacqGroupChannelEzumTemplateNet, DacqGroupChannelEnumTemplate, CDeviceGroupChannelInfoTemplateNet+
	CDigOut Stimulater FunctionNet
	CFtherConfigurationNet
	CPitarConfiguraionRegisterNet
	COmplemel'unctionNet
	ClasaMa_FaccionNet
	ClaterfaceboantPunctionNet
	CModisu, AxisParanstersNet
	CMcdau, ExtensionNet
	CMcdau_FYExtensioeNet
	CMcdlau_MossControlNet
	CMcollius, SemonNet
	CMcdBu_TempSenseeNet
	CMcsBur, VoltageModeNet
	CMcdlanNet
	Olio Urb Device di sur Probl'ancion Net
	CMLA2100c256PanctionNet
	CMeaAudieFunctionNet
	CMeaDigitaDataFunctionNet
	CMtaFeedbackFunctionNet
	CMcFunctionNet
	CMultiwellOstoStitusGuactionNet
	CPICFunctionNet
	C198_FunctionNet
	ChogualhessanCuveNet
	CPutseGeneratorFunctionNet
	CRFFunctionNet
	CRobo_FYIProgram_FunctionNet
	CRobo_FYTTemp_FunctionNet
	CSimila PanciseNet
	L
	CTE329 unctionNet
	CUshDeviceConfigurationPsuccionNet
	CW2100_StimulatorFunctionNet
	CWarnerUssingFunctionNet
	CWarnerValveControllerDeviceTesterFunctionNet
	CWireless Ruse Function Net

Public Member Functions

- CMcsUsbFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CMcsUsbFunctionNet (void)
- !CMcsUsbFunctionNet ()
- void ThrowCUsbExceptionNetOnError (uint32_t status)

Protected Member Functions

• CMcsUsbFunctionNet (CMcsUsbNet^ mcsusb, CMcsUsbFunctionPointerContainer^ mcsusbfunction)

Protected Attributes

- CMcsUsbNet ^ m_pMcsUsb
- CMcsUsbFunction * m_pMcsUsbFunction

11.53.1 Constructor & Destructor Documentation

```
11.53.1.1 CMcsUsbFunctionNet() [1/2] CMcsUsbFunctionNet (
            CMcsUsbNet^ mcsusb )
11.53.1.2 ~CMcsUsbFunctionNet() virtual ~CMcsUsbFunctionNet (
            void ) [virtual]
11.53.1.3 "!CMcsUsbFunctionNet() !CMcsUsbFunctionNet ()
11.53.1.4 CMcsUsbFunctionNet() [2/2] CMcsUsbFunctionNet (
            CMcsUsbNet^ mcsusb,
            CMcsUsbFunctionPointerContainer<sup>∧</sup> mcsusbfunction ) [protected]
11.53.2 Member Function Documentation
11.53.2.1 ThrowCUsbExceptionNetOnError() void ThrowCUsbExceptionNetOnError (
            uint32_t status )
11.53.3 Member Data Documentation
11.53.3.1 m pMcsUsb CMcsUsbNet ^ m_pMcsUsb [protected]
11.53.3.2 m pMcsUsbFunction CMcsUsbFunction* m_pMcsUsbFunction [protected]
11.54 CMcsUsbFunctionPointerContainer Class Reference
11.55 CMcsUsbListEntryNet Class Reference
```

McsUsbListEntryNet identifies a connected device.

Public Member Functions

- ∼CMcsUsbListEntryNet ()
- virtual bool Equals (Object[^] obj) override

Checks weather two CMcsUsbListEntryNet represent the same USB device.

void SetStringFormat (String ^ format)

Specify the text the CMcsUsbListEntryNet.ToString() function should return. The special code N expands to the device name and S expands to the serial number of the device.

• virtual String ^ ToString () override

Static Public Member Functions

static CMcsUsbListEntryNet [^] GetEntry ()

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

static CMcsUsbListEntryNet ^ GetEntry (DeviceEnumNet McsUsbDevice)

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

static CMcsUsbListEntryNet ^ GetEntry (DeviceEnumNet McsUsbDevice, unsigned int index)

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

• static unsigned int GetEntryCount ()

Returns the number of devices connected to the computer.

static unsigned int GetEntryCount (DeviceEnumNet McsUsbDevice)

Returns the number of devices connected to the computer.

Properties

• String Manufacturer [get]

The Manufacturer ID of the device represented by this CMcsUsbListEntryNet.

String[^] Product [get]

 $\textit{The Product ID of the device represented by this \textit{CMcsUsbListEntryNet}.}$

• String DeviceName [get]

The device name of the device represented by this CMcsUsbListEntryNet.

String[^] SerialNumber [get]

The serial number of the device represented by this CMcsUsbListEntryNet.

• String HwVersion [get]

The hardware revision of the device represented by this CMcsUsbListEntryNet.

• DeviceIdNet^ DeviceId [get]

11.55.1 Detailed Description

McsUsbListEntryNet identifies a connected device.

11.55.2 Constructor & Destructor Documentation

11.55.2.1 ~CMcsUsbListEntryNet() ~CMcsUsbListEntryNet ()

11.55.3 Member Function Documentation

```
11.55.3.1 Equals() virtual bool Equals (
Object^ obj ) [override], [virtual]
```

Checks weather two CMcsUsbListEntryNet represent the same USB device.

obj The CMcsUsbListEntryNet to compare with.

```
11.55.3.2 GetEntry() [1/3] static CMcsUsbListEntryNet ^ GetEntry ( ) [static]
```

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

Returns

A CMcsUsbListEntryNet to be used to connect to the device.

```
11.55.3.3 GetEntry() [2/3] static CMcsUsbListEntryNet ^ GetEntry (

DeviceEnumNet McsUsbDevice ) [static]
```

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

Parameters

McsUsbDevice | Specifies the type of devices to look for.

Returns

A CMcsUsbListEntryNet to be used to connect to the device.

```
11.55.3.4 GetEntry() [3/3] static CMcsUsbListEntryNet ^ GetEntry (

DeviceEnumNet McsUsbDevice,

unsigned int index ) [static]
```

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

Parameters

McsUsbDevice Specifies the type of devices to look for.

index	number of the entry to use.
-------	-----------------------------

Returns

A CMcsUsbListEntryNet to be used to connect to the device.

11.55.3.5 GetEntryCount() [1/2] static unsigned int GetEntryCount () [static]

Returns the number of devices connected to the computer.

Returns

The number of devices.

```
11.55.3.6 GetEntryCount() [2/2] static unsigned int GetEntryCount (

DeviceEnumNet McsUsbDevice ) [static]
```

Returns the number of devices connected to the computer.

Parameters

McsUsbDevice	Specifies the type of devices to look for.
--------------	--

Returns

The number of devices.

Specify the text the CMcsUsbListEntryNet.ToString() function should return. The special code N expands to the device name and S expands to the serial number of the device.

format A String containing the format template. Default is "%N (%S)".

11.55.3.8 ToString() virtual String ^ ToString () [override], [virtual]

11.55.4 Property Documentation

11.55.4.1 DeviceId DeviceIdNet^ DeviceId [get]

11.55.4.2 DeviceName String DeviceName [get]

The device name of the device represented by this CMcsUsbListEntryNet.

11.55.4.3 HwVersion String MwVersion [get]

The hardware revision of the device represented by this CMcsUsbListEntryNet.

11.55.4.4 Manufacturer String^ Manufacturer [get]

The Manufacturer ID of the device represented by this CMcsUsbListEntryNet.

11.55.4.5 **Product** String^ Product [get]

The Product ID of the device represented by this CMcsUsbListEntryNet.

11.55.4.6 SerialNumber String SerialNumber [get]

The serial number of the device represented by this CMcsUsbListEntryNet.

11.56 CMcsUsbListNet Class Reference

Class to handle a list of connected MCS USB devices.

Public Member Functions

• CMcsUsbListNet (DeviceEnumNet McsUsbDevice)

Initializes a new instance of CMcsUsbListNet class.

CMcsUsbListNet (array< DeviceIdNet[^]>[^] DeviceIdList)

Initializes a new instance of CMcsUsbListNet class.

∼CMcsUsbListNet ()

Destructor: called by Dispose()

• !CMcsUsbListNet ()

Finalizer: called by GC before collecting

void SetStringFormat (String ^ format)

Specify the text the CMcsUsbListEntryNet.ToString() function should return. The special code N expands to the device name and S expands to the serial number of the device.

• uint32 t GetNumberOfDevices ()

Gets the number of devices currently in the list.

CMcsUsbListEntryNet ^ GetUsbListEntry (unsigned int index)

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

array< CMcsUsbListEntryNet[^]> [^] GetUsbListEntries ()

Returns all entries from the list of USB Devices connected to the computer.

bool IsDeviceTypeOf (CMcsUsbListEntryNet[^] entry, DeviceEnumNet McsUsbDevice)

Properties

• uint32_t Count [get]

Gets the number of devices currently in the list.

Events

- OnDeviceArrivalRemoval[^] DeviceArrival
- OnDeviceArrivalRemoval[^] DeviceRemoval

11.56.1 Detailed Description

Class to handle a list of connected MCS USB devices.

11.56.2 Constructor & Destructor Documentation

```
11.56.2.1 CMcsUsbListNet() [1/2] CMcsUsbListNet (

DeviceEnumNet McsUsbDevice )
```

Initializes a new instance of CMcsUsbListNet class.

```
11.56.2.2 CMcsUsbListNet() [2/2] CMcsUsbListNet (

array< DeviceIdNet^>^ DeviceIdList )
```

Initializes a new instance of CMcsUsbListNet class.

```
11.56.2.3 ~CMcsUsbListNet() ~CMcsUsbListNet ()
```

Destructor: called by Dispose()

11.56.2.4 "!CMcsUsbListNet() !CMcsUsbListNet ()

Finalizer: called by GC before collecting

11.56.3 Member Function Documentation

```
11.56.3.1 GetNumberOfDevices() uint32_t GetNumberOfDevices ( )
```

Gets the number of devices currently in the list.

Returns

The number of devices currently in the list.

```
11.56.3.2 GetUsbListEntries() array<CMcsUsbListEntryNet^> ^ GetUsbListEntries ( )
```

Returns all entries from the list of USB Devices connected to the computer.

```
11.56.3.3 GetUsbListEntry() CMcsUsbListEntryNet ^ GetUsbListEntry ( unsigned int index )
```

Returns one CMcsUsbListEntryNet from the list of USB Devices connected to the computer.

Parameters

index | number of the entry to use.

```
11.56.3.5 SetStringFormat() void SetStringFormat ( String ^{\wedge} format )
```

Specify the text the CMcsUsbListEntryNet.ToString() function should return. The special code N expands to the device name and S expands to the serial number of the device.

Parameters

format A String containing the format template. Default is "%	N (%S)".
---	----------

11.56.4 Property Documentation

11.56.4.1 Count uint32_t Count [get]

Gets the number of devices currently in the list.

11.56.5 Event Documentation


11.56.5.1 DeviceArrival OnDeviceArrivalRemoval^ DeviceArrival

11.56.5.2 DeviceRemoval OnDeviceArrivalRemoval^ DeviceRemoval

11.57 CMcsUsbNet Class Reference

Base class to handle MCS USB devices. All device classes are derived from this class. Functionality that is provided by all MCS devices is handled by this class.

Inheritance diagram for CMcsUsbNet:



Public Member Functions

• CMcsUsbNet ()

Initializes a new instance of the base class to handle MCS USB devices.

CMcsUsbNet (McsBusTypeEnumNet bustype)

Initializes a new instance of the base class to handle MCS USB devices.

- virtual ∼CMcsUsbNet ()
- !CMcsUsbNet ()
- DeviceEnumNet GetDeviceEnum ()
- virtual uint32_t Connect ()

Opens a connection to the device.

virtual uint32_t Connect (unsigned int LockMask)

Opens a connection to the device.

virtual uint32_t Connect (CMcsUsbListEntryNet[^] entry)

Opens a connection to the device.

virtual uint32_t Connect (CMcsUsbListEntryNet[^] entry, unsigned int LockMask)

Opens a connection to the device.

- virtual uint32 t GetStatus ([System::Runtime::InteropServices::Out]uint32 t% iStatus)
- virtual bool IsConnected ()

Check if a device is Connected.

virtual void Disconnect ()

Disconnect from a device.

- CMcsUsbListEntryNet ^ GetUsbListEntry ()
- virtual String ^ GetSerialNumber ()

Query the Serial Number of the device.

- DriverVersionNet ^ GetVersion ()
- DriverVersionNet ^ GetVersion (CFirmwareDestinationNet dest)
- DeviceIdNet ^ GetDeviceId ()
- uint32_t GetIdent ([System::Runtime::InteropServices::Out]String^% Answer)
- void MultibootSelectImage (unsigned int sector)

Select the multiboot image specified by "sector" (range: 0..2) for IFB FPGA.

String \(^\) MultibootGetImageId (unsigned int sector)

Query the multiboot image id of the device located in specified sector (range: 0..2 / 0..9) of IFB1 / IFB2 FPGA.

uint32 t MultibootGetCypressImageId (unsigned int sector)

Query the multiboot image id of the device located in specified sector (range: 0..9 0..9) of IFB2 Cypress.

uint32_t MultibootGetSelectedImage ()

Gets sector index of selected FPGA boot image on IFB

uint32_t GetMea21UsbPort ()

Gets the USB port if an IFB that is used by this connection

HeadstageIdEnumNet GetHeadstageID (uint32_t headstage)

Gets the ID of a connected headstage.

bool GetHeadstagePresent (uint32 t headstage)

queries whether a headstage is present

bool GetHeadstageActive (uint32_t headstage)

queries whether a headstage is active

• void RescanHeadstage (uint32 t headstage)

rescans and activates a headstage

- array< BYTE > ^ GetSoftwareKey (unsigned int index)
- void SetSoftwareKey (unsigned int index, array< BYTE >^ buffer)
- void RemoveSoftwareKey (unsigned int index)
- void AddSoftwareKey (String[^] key)
- bool EmptyKey (String[^] key)
- bool ValidKey (String[^] key, [System::Runtime::InteropServices::Out]String[^]% serial_number)
- bool HasSoftwareKey (uint8_t ProgrammID, uint8_t majorversion)
- bool HasSoftwareKey (SoftwareKeyProgrammIdsNet::ProgrammIdsNet ProgrammID, uint8_t majorversion)
- String \(^{\text{GetSoftwareKeyString}}\) (uint8 t ProgrammID, uint8 t majorversion)
- String ^ GetSoftwareKeyString (SoftwareKeyProgrammIdsNet::ProgrammIdsNet ProgrammID, uint8_t majorversion)
- bool IsDeviceHighSpeedCapable ()
- bool IsDeviceHighSpeed ()
- McsUsbSpeedEnumNet GetDeviceCapableSpeed ()
- McsUsbSpeedEnumNet GetDeviceSpeed ()

Query the Connection Speed of the device.

- unsigned int TxnTestMemoryWrite (unsigned short index)
- unsigned int TxnTestMemoryReadAndCheck (unsigned short index)
- void TxnSetSerialNumber (unsigned int number)
- unsigned int TxnGetSerialNumber ()
- unsigned int ReadRegister (unsigned int reg)

- array< uint32 t > ^ ReadRegister (unsigned int reg, int length)
- unsigned int ReadRegister32 (unsigned int adr)
- unsigned int ReadRegisterTimeSlot (unsigned int reg, int TimeSlot)
- void WriteRegister (unsigned int reg, unsigned int value)
- void WriteRegisterValue (unsigned int reg, unsigned int value)
- void WriteRegister32 (unsigned int adr, unsigned int value)
- void WriteRegister (unsigned int reg, array< unsigned int >[∧] values)
- void WriteRegisterArray (unsigned int reg, array< unsigned int >^ values)
- · void WriteRegisterTimeSlot (unsigned int reg, unsigned int value, int TimeSlot)
- void WriteRegisterTimeSlot (unsigned int reg, array< unsigned int >^ values, int TimeSlot)
- bool ReadEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, [System::Runtime::
 —
 InteropServices::Out]uint32 t% DMA value)
- bool ReadEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, [System::Runtime::
 —
 InteropServices::Out]uint32_t% DMA_value, uint32_t EEPROMSize)
- bool ReadEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, [System::Runtime::
 —
 InteropServices::Out]uint32_t% DMA_value, uint32_t EEPROMSize, uint32_t EepromStartAddress)
- void WriteEepromRegisterPreconfig (uint32 t EEPROMBase, uint32 t DMA reg, uint32 t DMA value)
- void WriteEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, uint32_t DMA_value, uint32_t EEPROMSize)
- void WriteEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, uint32_t DMA_value, uint32_t EEPROMSize, uint32_t EepromStartAddress)
- void EraseEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg)
- void EraseEepromRegisterPreconfig (uint32 t EEPROMBase, uint32 t DMA reg, uint32 t EEPROMSize)
- void EraseEepromRegisterPreconfig (uint32_t EEPROMBase, uint32_t DMA_reg, uint32_t EEPROMSize, uint32_t EepromStartAddress)
- unsigned int GetLastUSBError ()
- void ThrowCUsbExceptionNetOnError (uint32 t status)
- bool GetDeviceCannotStallOutRequests ()
- String \(^\) GetHardwareRevision ()
- unsigned int GetFirmwareVersion (CFirmwareDestinationNet destination)

Gets the firmware version for the destination.

- uint8 t GetNumConfigurations ()
- uint8 t GetConfiguration ()
- void SetConfiguration (uint8 t config)
- uint32 t GetDeviceRootHubVendorID ()

Gets the Vendor ID of the USB root hub the device is connected to.

UsbVendorIdEnumNet GetDeviceRootHubVendorEnum ()

Gets the Vendor ID of the USB root hub the device is connected to.

String \(^\) GetDeviceRootHubVendorName ()

Gets the Vendor Name of the USB root hub the device is connected to.

void EnableExceptions (bool enable)

Enables or Disables Exceptions for calls to McsUsb Devices. If Exceptions are disabled, the return value of a command can be queries with the GetStatusOfLastCommand call instead.

- bool IsExceptionsEnabled ()
- uint32_t GetStatusOfLastCommand ()

Gets the status of the last call to the McsUsb Library.

- uint32 t CyclePort ()
- void AssociateToThis (CMcsUsbNet[^] device)

Static Public Member Functions

static String \(^\) GetErrorText (unsigned int Status)

Gets the error text string that belongs to a status number.

Static Public Attributes

```
    static const uint32 t Status Crc = (0xE0100001L)

    static const uint32 t Status Btstuff = (0xE0100002L)

    static const uint32 t Status DataToggleMismatch = (0xE0100003L)

    static const uint32 t Status Stall = (0xE0100004L)

    static const uint32_t Status_DevNotResponding = (0xE0100005L)

    static const uint32 t Status PidCheckFailure = (0xE0100006L)

    static const uint32_t Status_UnexpectedPid = (0xE0100007L)

• static const uint32_t Status_DataOverrun = (0xE0100008L)

    static const uint32 t Status DataUnderrun = (0xE0100009L)

    static const uint32 t Status BufferOverrun = (0xE010000CL)

    static const uint32 t Status BufferUnderrun = (0xE010000DL)

    static const uint32 t Status NotAccessed = (0xE010000FL)

    static const uint32_t Status_Fifo = (0xE0100010L)

    static const uint32 t Status EndpointHalted = (0xE0100030L)

    static const uint32 t Status NoMemory = (0xE0100100L)

    static const uint32_t Status_InvalidUrbFunction = (0xE0100200L)

    static const uint32 t Status InvalidParameter = (0xE0100300L)

    static const uint32_t Status_InvalidDeviceHandle = (0xE0100013L)

    static const uint32 t Status InvalidHandle = (0xE0100012L)

    static const uint32 t Status ErrorBusy = (0xE0100400L)

    static const uint32_t Status_RequestFailed = (0xE0100500L)

    static const uint32 t Status InvalidPipeHandle = (0xE0100600L)

    static const uint32 t Status NoBandwidth = (0xE0100700L)

    static const uint32_t Status_InternalHcError = (0xE0100800L)

    static const uint32 t Status ErrorShortTransfer = (0xE0100900L)

    static const uint32 t Status BadStartFrame = (0xE0100A00L)

    static const uint32 t Status IsochRequestFailed = (0xE0100B00L)

    static const uint32_t Status_FrameControlOwned = (0xE0100C00L)

    static const uint32_t Status_ControlNotOwned = (0xE0100D00L)

    static const uint32 t Status Canceled = (0xE0110000L)

    static const uint32 t Status Canceling = (0xE0120000L)

    static const uint32 t Status AlreadyConfigured = (0xE0110001L)

    static const uint32 t Status Unconfigured = (0xE0110002L)

    static const uint32_t Status_NoSuchDevice = (0xE01F0002L)

    static const uint32 t Status DeviceNotFound = (0xE01F0003L)

    static const uint32 t Status NotSupported = (0xE01F0005L)

    static const uint32 t Status IoPending = (0xE01F0006L)

    static const uint32 t Status IoTimeout = (0xE01F0007L)

    static const uint32_t Status_DeviceRemoved = (0xE01F0008L)

    static const uint32_t Status_PipeNotLinked = (0xE01F0009L)

    static const uint32 t Status ConnectedPipes = (0xE01F000AL)

    static const uint32_t Status_DeviceLocked = (0xE01F0010L)

    static const uint32 t Status RequestMutexTimeout = (0xE01F0020L)

    static const uint32 t Status RequestMutexFailed = (0xE01F0021L)

    static const uint32 t Status LastUsbErrorMismatch = (0xE01F0022L)

    static const uint32 t WPAError ScanningIsPending = ( (0xA0220000L) | 0x0036 )
```

Properties

virtual String[^] SerialNumber [get]

11.57.1 Detailed Description

Base class to handle MCS USB devices. All device classes are derived from this class. Functionality that is provided by all MCS devices is handled by this class.

11.57.2 Constructor & Destructor Documentation

```
11.57.2.1 CMcsUsbNet() [1/2] CMcsUsbNet ( )
```

Initializes a new instance of the base class to handle MCS USB devices.

```
11.57.2.2 CMcsUsbNet() [2/2] CMcsUsbNet (

McsBusTypeEnumNet bustype)
```

Initializes a new instance of the base class to handle MCS USB devices.

Parameters

```
bustype Type of device to use, either USB or PCI.
```

```
11.57.2.3 \sim CMcsUsbNet() virtual \sim CMcsUsbNet () [virtual]
```

```
11.57.2.4 "!CMcsUsbNet() !CMcsUsbNet ()
```

11.57.3 Member Function Documentation

```
11.57.3.1 AddSoftwareKey() void AddSoftwareKey ( String^{\land} key )
```

```
11.57.3.2 AssociateToThis() void AssociateToThis (

CMcsUsbNet^ device )
```

11.57.3.3 Connect() [1/4] virtual uint32_t Connect () [virtual]

Opens a connection to the device.

Returns

Error Status. 0 on success.

```
11.57.3.4 Connect() [2/4] virtual uint32_t Connect (
CMcsUsbListEntryNet^ entry ) [virtual]
```

Opens a connection to the device.

Parameters

entry	The Device List Entry for the device to be connected.
-------	---

Returns

Error Status. 0 on success.

```
11.57.3.5 Connect() [3/4] virtual uint32_t Connect (

CMcsUsbListEntryNet^ entry,

unsigned int LockMask ) [virtual]
```

Opens a connection to the device.

Parameters

entry	The Device List Entry for the device to be connected.
LockMask	The Lock Mask for this connection.

Returns

Error Status. 0 on success.

```
11.57.3.6 Connect() [4/4] virtual uint32_t Connect (
unsigned int LockMask) [virtual]
```

Opens a connection to the device.

Parameters

LockMask The Lock Mask for this connection	۱.
--	----

Returns

Error Status. 0 on success.

```
11.57.3.7 CyclePort() uint32_t CyclePort ()
```

```
11.57.3.8 Disconnect() virtual void Disconnect ( ) [virtual]
```

Disconnect from a device.

```
11.57.3.9 EmptyKey() bool EmptyKey ( String^{\wedge} key )
```

```
11.57.3.10 EnableExceptions() void EnableExceptions (
bool enable)
```

Enables or Disables Exceptions for calls to McsUsb Devices. If Exceptions are disabled, the return value of a command can be queries with the GetStatusOfLastCommand call instead.

Parameters

```
enable True to enable Exceptions, False to disable.
```

```
11.57.3.11 EraseEepromRegisterPreconfig() [1/3] void EraseEepromRegisterPreconfig ( uint32_t EEPROMBase, uint32_t DMA_reg )
```

```
11.57.3.12 EraseEepromRegisterPreconfig() [2/3] void EraseEepromRegisterPreconfig (
```

```
uint32_t EEPROMBase,
uint32_t DMA_reg,
uint32_t EEPROMSize )
```

```
\textbf{11.57.3.13} \quad \textbf{EraseEepromRegisterPreconfig() [3/3]} \quad \texttt{void EraseEepromRegisterPreconfig ()} \\
               uint32_t EEPROMBase,
               uint32_t DMA_reg,
               uint32_t EEPROMSize,
               uint32_t EepromStartAddress )
\textbf{11.57.3.14} \quad \textbf{GetConfiguration()} \quad \texttt{uint8\_t GetConfiguration ()}
\textbf{11.57.3.15} \quad \textbf{GetDeviceCannotStallOutRequests()} \quad \texttt{bool GetDeviceCannotStallOutRequests ()} \\
11.57.3.16 GetDeviceCapableSpeed() McsUsbSpeedEnumNet GetDeviceCapableSpeed ( )
11.57.3.17 GetDeviceEnum() DeviceEnumNet GetDeviceEnum ( )
11.57.3.18 GetDeviceId() DeviceIdNet ^ GetDeviceId ( )
11.57.3.19 GetDeviceRootHubVendorEnum() UsbVendorIdEnumNet GetDeviceRootHubVendorEnum ( )
Gets the Vendor ID of the USB root hub the device is connected to.
Returns
     An enum which enumerates the PCI Vendor ID.
11.57.3.20 GetDeviceRootHubVendorID() uint32_t GetDeviceRootHubVendorID ( )
Gets the Vendor ID of the USB root hub the device is connected to.
Returns
```

The PCI Vendor ID, 0x8086 for Intel, 0x1912 for Renesas, 0x1b21 for ASMedia.

11.57.3.21 GetDeviceRootHubVendorName() String ^ GetDeviceRootHubVendorName ()

Gets the Vendor Name of the USB root hub the device is connected to.

Returns

The PCI Vendor Name, either "Intel", "Renesas", "ASMedia" or "unknown".

11.57.3.22 GetDeviceSpeed() McsUsbSpeedEnumNet GetDeviceSpeed ()

Query the Connection Speed of the device.

Returns

0 for Low-Speed, 1 for Full-Speed, 2 for High-Speed and 3 for SuperSpeed.

Gets the error text string that belongs to a status number.

Parameters

Status	The status number you want the text for.
--------	--

Returns

The error text string that belongs to the status number.

11.57.3.24 **GetFirmwareVersion()** unsigned int GetFirmwareVersion (CFirmwareDestinationNet destination)

Gets the firmware version for the destination.

Parameters

destination	η.	The destination to be queried.

Returns

The firmware version as a 32 bit number, the upper 16 bit contain the majaor version number, the lower 16 bit the minor version number.

11.57.3.25 GetHardwareRevision() String $^{\wedge}$ GetHardwareRevision ()

```
11.57.3.26 GetHeadstageActive() bool GetHeadstageActive ( uint32_t headstage )
```

queries whether a headstage is active

Parameters

in <i>headstage</i> th	e headstage number (0 or 1)
------------------------	-----------------------------

Returns

true if the headstage is active

11.57.3.27 **GetHeadstagelD()** HeadstageIdEnumNet GetHeadstageID (uint32_t headstage)

Gets the ID of a connected headstage.

Parameters

ſ			
١	in	headstage	the headstage number (0 or 1)

Returns

enumerated Headstage ID

11.57.3.28 GetHeadstagePresent() bool GetHeadstagePresent (uint32_t headstage)

queries whether a headstage is present

Parameters

in	headstage	the headstage number (0 or 1)

Returns

true if the headstage is present

```
11.57.3.29 GetIdent() uint32_t GetIdent (
              [System::Runtime::InteropServices::Out] String^% Answer )
11.57.3.30 GetLastUSBError() unsigned int GetLastUSBError ()
11.57.3.31 GetMea21UsbPort() uint32_t GetMea21UsbPort ( )
Gets the USB port if an IFB that is used by this connection
Returns
     number of used port; range: 0..1
11.57.3.32 GetNumConfigurations() uint8_t GetNumConfigurations ( )
11.57.3.33 GetSerialNumber() virtual String ^ GetSerialNumber ( ) [virtual]
Query the Serial Number of the device.
Returns
     The Serial Number.
11.57.3.34 GetSoftwareKey() array<br/>SYTE> ^{\wedge} GetSoftwareKey (
             unsigned int index)
11.57.3.35 GetSoftwareKeyString() [1/2] String ^ GetSoftwareKeyString (
             SoftwareKeyProgrammIdsNet::ProgrammIdsNet ProgrammID,
             uint8_t majorversion )
11.57.3.36 GetSoftwareKeyString() [2/2] String ^ GetSoftwareKeyString (
             uint8_t ProgrammID,
             uint8_t majorversion )
```

```
11.57.3.37 GetStatus() virtual uint32_t GetStatus (
             [System::Runtime::InteropServices::Out] uint32_t% iStatus ) [virtual]
11.57.3.38 GetStatusOfLastCommand() uint32_t GetStatusOfLastCommand ( )
Gets the status of the last call to the McsUsb Library.
Returns
    The Error Status of the last McsUsb command. 0 on success.
11.57.3.39 GetUsbListEntry() CMcsUsbListEntryNet ^ GetUsbListEntry ( )
11.57.3.40 GetVersion() [1/2] DriverVersionNet ^ GetVersion ( )
11.57.3.41 GetVersion() [2/2] DriverVersionNet ^ GetVersion (
             CFirmwareDestinationNet dest )
11.57.3.42 HasSoftwareKey() [1/2] bool HasSoftwareKey (
             SoftwareKeyProgrammIdsNet::ProgrammIdsNet ProgrammID,
             uint8_t majorversion )
11.57.3.43 HasSoftwareKey() [2/2] bool HasSoftwareKey (
             uint8_t ProgrammID,
             uint8_t majorversion )
11.57.3.44 IsConnected() virtual bool IsConnected ( ) [virtual]
Check if a device is Connected.
```

Generated by Doxygen

true if the device is connected.

Returns

```
11.57.3.45 IsDeviceHighSpeed() bool IsDeviceHighSpeed ()
```

```
11.57.3.46 IsDeviceHighSpeedCapable() bool IsDeviceHighSpeedCapable ( )
```

```
11.57.3.48 MultibootGetCypressImageId() uint32_t MultibootGetCypressImageId ( unsigned int sector )
```

Query the multiboot image id of the device located in specified sector (range: 0..9 0..9) of IFB2 Cypress.

Returns

The magic ident code of the image.

```
11.57.3.49 MultibootGetImageId() String ^ MultibootGetImageId ( unsigned int sector )
```

Query the multiboot image id of the device located in specified sector (range: 0..2 / 0..9) of IFB1 / IFB2 FPGA.

Returns

The magic ident code of the image.

```
11.57.3.50 MultibootGetSelectedImage() uint32_t MultibootGetSelectedImage ()
```

Gets sector index of selected FPGA boot image on IFB

Returns

Sector index of image; range: 0..2

```
11.57.3.51 MultibootSelectImage() void MultibootSelectImage ( unsigned int sector )
```

Select the multiboot image specified by "sector" (range: 0..2) for IFB FPGA.

Returns

```
Throws exception on error.
11.57.3.52 ReadEepromRegisterPreconfig() [1/3] bool ReadEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             [System::Runtime::InteropServices::Out] uint32_t% DMA_value )
11.57.3.53 ReadEepromRegisterPreconfig() [2/3] bool ReadEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             [System::Runtime::InteropServices::Out] uint32_t% DMA_value,
             uint32_t EEPROMSize )
11.57.3.54 ReadEepromRegisterPreconfig() [3/3] bool ReadEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             [System::Runtime::InteropServices::Out] uint32_t% DMA_value,
             uint32_t EEPROMSize,
             uint32_t EepromStartAddress )
11.57.3.55 ReadRegister() [1/2] unsigned int ReadRegister (
             unsigned int reg )
11.57.3.56 ReadRegister() [2/2] array<uint32_t> ^ ReadRegister (
             unsigned int reg,
             int length )
11.57.3.57 ReadRegister32() unsigned int ReadRegister32 (
```

unsigned int adr)

```
11.57.3.58 ReadRegisterTimeSlot() unsigned int ReadRegisterTimeSlot (
              unsigned int reg,
              int TimeSlot )
11.57.3.59 RemoveSoftwareKey() void RemoveSoftwareKey (
             unsigned int index )
\textbf{11.57.3.60} \quad \textbf{RescanHeadstage()} \quad \texttt{void RescanHeadstage ()}
             uint32_t headstage )
rescans and activates a headstage
Parameters
       headstage
                    the headstage number (0 or 1)
   in
11.57.3.61 SetConfiguration() void SetConfiguration (
             uint8_t config )
11.57.3.62 SetSoftwareKey() void SetSoftwareKey (
             unsigned int index,
              array < BYTE >^{\wedge} buffer )
11.57.3.63 ThrowCUsbExceptionNetOnError() void ThrowCUsbExceptionNetOnError (
             uint32_t status )
11.57.3.64 TxnGetSerialNumber() unsigned int TxnGetSerialNumber ()
11.57.3.65 TxnSetSerialNumber() void TxnSetSerialNumber (
              unsigned int number )
```

```
11.57.3.66 TxnTestMemoryReadAndCheck() unsigned int TxnTestMemoryReadAndCheck (
             unsigned short index)
11.57.3.67 TxnTestMemoryWrite() unsigned int TxnTestMemoryWrite (
             unsigned short index)
11.57.3.68 ValidKey() [1/2] bool ValidKey (
             String^{\wedge} key,
             [System::Runtime::InteropServices::Out] \ String^{\ } \ serial\_number \ )
11.57.3.69 ValidKey() [2/2] bool ValidKey (
             String^{\wedge} key,
             uint8_t ProgrammID,
             uint8_t majorversion,
             [System::Runtime::InteropServices::Out] String^% serial_number )
11.57.3.70 WriteEepromRegisterPreconfig() [1/3] void WriteEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             uint32_t DMA_value )
11.57.3.71 WriteEepromRegisterPreconfig() [2/3] void WriteEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             uint32_t DMA_value,
             uint32_t EEPROMSize )
11.57.3.72 WriteEepromRegisterPreconfig() [3/3] void WriteEepromRegisterPreconfig (
             uint32_t EEPROMBase,
             uint32_t DMA_reg,
             uint32_t DMA_value,
             uint32_t EEPROMSize,
             uint32_t EepromStartAddress )
```

```
11.57.3.73 WriteRegister() [1/2] void WriteRegister (
             unsigned int reg,
             array< unsigned int >^{\wedge} values)
11.57.3.74 WriteRegister() [2/2] void WriteRegister (
             unsigned int reg,
             unsigned int value )
11.57.3.75 WriteRegister32() void WriteRegister32 (
             unsigned int adr,
             unsigned int value )
11.57.3.76 WriteRegisterArray() void WriteRegisterArray (
             unsigned int reg,
             array< unsigned int >^{\land} values)
11.57.3.77 WriteRegisterTimeSlot() [1/2] void WriteRegisterTimeSlot (
             unsigned int reg,
             array< unsigned int >^{\wedge} values,
             int TimeSlot )
11.57.3.78 WriteRegisterTimeSlot() [2/2] void WriteRegisterTimeSlot (
             unsigned int reg,
             unsigned int value,
             int TimeSlot )
11.57.3.79 WriteRegisterValue() void WriteRegisterValue (
             unsigned int reg,
             unsigned int value )
```

11.57.4 Member Data Documentation

```
11.57.4.1 Status_AlreadyConfigured const uint32_t Status_AlreadyConfigured = (0xE0110001L)
[static]
```

- 11.57.4.2 Status_BadStartFrame const uint32_t Status_BadStartFrame = (0xE0100A00L) [static]
- 11.57.4.3 Status_Btstuff const uint32_t Status_Btstuff = (0xE0100002L) [static]
- 11.57.4.4 Status_BufferOverrun const uint32_t Status_BufferOverrun = (0xE010000CL) [static]
- 11.57.4.5 Status_BufferUnderrun const uint32_t Status_BufferUnderrun = (0xE010000DL) [static]
- 11.57.4.6 Status_Canceled const uint32_t Status_Canceled = (0xE0110000L) [static]
- 11.57.4.7 Status_Canceling const uint32_t Status_Canceling = (0xE0120000L) [static]
- 11.57.4.8 Status_ConnectedPipes const uint32_t Status_ConnectedPipes = (0xE01F000AL) [static]
- 11.57.4.9 Status_ControlNotOwned const uint32_t Status_ControlNotOwned = (0xE0100D00L) [static]
- 11.57.4.10 Status_Crc const uint32_t Status_Crc = (0xE01000001L) [static]
- 11.57.4.11 Status_DataOverrun const uint32_t Status_DataOverrun = (0xE0100008L) [static]

```
11.57.4.12 Status_DataToggleMismatch const uint32_t Status_DataToggleMismatch = (0xE0100003L)
[static]
11.57.4.13 Status_DataUnderrun const uint32_t Status_DataUnderrun = (0xE0100009L) [static]
11.57.4.14 Status_DeviceLocked const uint32_t Status_DeviceLocked = (0xE01F0010L) [static]
11.57.4.15 Status_DeviceNotFound const uint32_t Status_DeviceNotFound = (0xE01F0003L) [static]
11.57.4.16 Status DeviceRemoved const uint32_t Status_DeviceRemoved = (0xE01F0008L) [static]
11.57.4.17 Status_DevNotResponding const uint32_t Status_DevNotResponding = (0xE0100005L)
[static]
11.57.4.18 Status_EndpointHalted const uint32_t Status_EndpointHalted = (0xE0100030L) [static]
11.57.4.19 Status_ErrorBusy const uint32_t Status_ErrorBusy = (0xE0100400L) [static]
```

- 11.57.4.20 Status_ErrorShortTransfer const uint32_t Status_ErrorShortTransfer = (0xE0100900L) [static]
- 11.57.4.21 Status_Fifo const uint32_t Status_Fifo = (0xE0100010L) [static]
- 11.57.4.22 Status_FrameControlOwned const uint32_t Status_FrameControlOwned = (0xE0100C00L) [static]

```
11.57.4.23 Status_InternalHcError const uint32_t Status_InternalHcError = (0xE0100800L) [static]
```

11.57.4.24 Status_InvalidDeviceHandle const uint32_t Status_InvalidDeviceHandle = (0xE0100013L) [static]

11.57.4.25 Status_InvalidHandle const uint32_t Status_InvalidHandle = (0xE0100012L) [static]

11.57.4.26 Status_InvalidParameter const uint32_t Status_InvalidParameter = (0xE0100300L) [static]

11.57.4.27 Status_InvalidPipeHandle const uint32_t Status_InvalidPipeHandle = (0xE0100600L) [static]

11.57.4.28 Status_InvalidUrbFunction const uint32_t Status_InvalidUrbFunction = (0xE0100200L) [static]

11.57.4.29 Status loPending const uint32_t Status_IoPending = (0xE01F0006L) [static]

11.57.4.30 Status_IoTimeout const uint32_t Status_IoTimeout = (0xE01F0007L) [static]

11.57.4.31 Status_IsochRequestFailed const uint32_t Status_IsochRequestFailed = (0xE0100B00L) [static]

11.57.4.32 Status_LastUsbErrorMismatch const uint32_t Status_LastUsbErrorMismatch = (0xE01← F0022L) [static]

```
11.57.4.33 Status_NoBandwidth const uint32_t Status_NoBandwidth = (0xE0100700L) [static]
11.57.4.34 Status_NoMemory const uint32_t Status_NoMemory = (0xE0100100L) [static]
11.57.4.35 Status NoSuchDevice const uint32_t Status_NoSuchDevice = (0xE01F0002L) [static]
11.57.4.36 Status_NotAccessed const uint32_t Status_NotAccessed = (0xE010000FL) [static]
11.57.4.37 Status_NotSupported const uint32_t Status_NotSupported = (0xE01F0005L) [static]
11.57.4.38 Status_PidCheckFailure const uint32_t Status_PidCheckFailure = (0xE0100006L) [static]
11.57.4.39 Status_PipeNotLinked const uint32_t Status_PipeNotLinked = (0xE01F0009L) [static]
11.57.4.40 Status_RequestFailed const uint32_t Status_RequestFailed = (0xE0100500L) [static]
11.57.4.41 Status_RequestMutexFailed const uint32_t Status_RequestMutexFailed = (0xE01F0021L)
[static]
11.57.4.42 Status_RequestMutexTimeout const uint32_t Status_RequestMutexTimeout = (0xE01↔
F0020L) [static]
```

11.57.4.43 Status_Stall const uint32_t Status_Stall = (0xE0100004L) [static]

```
11.57.4.44 Status_Unconfigured const uint32_t Status_Unconfigured = (0xE0110002L) [static]
```

```
11.57.4.45 Status_UnexpectedPid const uint32_t Status_UnexpectedPid = (0xE0100007L) [static]
```

```
11.57.4.46 WPAError_ScanninglsPending const uint32_t WPAError_ScanningIsPending = ( (0x \leftarrow A0220000L) \mid 0x0036 ) [static]
```

11.57.5 Property Documentation

```
11.57.5.1 SerialNumber virtual String^ SerialNumber [get]
```

11.58 CMcsUsbPointerContainer Class Reference

11.59 CMEA2100_256DacqGroupChannelSelectionNet Class Reference

Inheritance diagram for CMEA2100_256DacqGroupChannelSelectionNet:



Public Member Functions

CMEA2100_256DacqGroupChannelSelectionNet (CMcsUsbNet[^] mcsusb)

Additional Inherited Members

11.59.1 Constructor & Destructor Documentation


```
11.59.1.1 CMEA2100_256DacqGroupChannelSelectionNet() CMEA2100_256DacqGroupChannelSelectionNet (

CMcsUsbNet^ mcsusb )
```

11.60 CMEA2100x256FunctionNet Class Reference

CMEA2100x256FunctionNet is the class to control the MEA2100-256 device needs #include "Stg200xNet.h" to resolve documentation reference

Inheritance diagram for CMEA2100x256FunctionNet:



Public Member Functions

CMEA2100x256FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pME
 — A2100x256FunctionPointerContainer)

Initializes a new instance of the CMEA2100x256FunctionNet class.

- CMEA2100x256FunctionNet (CMcsUsbNet[∧] mcsusb)
- virtual ~CMEA2100x256FunctionNet ()
- !CMEA2100x256FunctionNet ()
- StimulationLayoutConfigurationEnumNet GetLayoutConfiguration ()

Gets the stimulation layout configuration. Can be single well, 6-well or 9-well. The number of $D \leftarrow AC$ channels available per well is Mcs::Usb::Cstg200xBasicNet::GetNumberOfAnalogChannels divided by Mcs::Usb::Cstg200xBasicNet::GetNumberOfStimulationSourcesPerElectrode.

void SetLayoutConfiguration (StimulationLayoutConfigurationEnumNet LayoutConfiguration)

Sets the stimulation layout configuration. Can be single well, 6-well or 9-well. The number of D← AC channels available per well is Mcs::Usb::CStg200xBasicNet::GetNumberOfAnalogChannels divided by Mcs::Usb::CStg200xBasicNet::GetNumberOfStimulationSourcesPerElectrode.

Additional Inherited Members

11.60.1 Detailed Description

CMEA2100x256FunctionNet is the class to control the MEA2100-256 device needs #include "Stg200xNet.h" to resolve documentation reference

11.60.2 Constructor & Destructor Documentation

```
11.60.2.1 CMEA2100x256FunctionNet() [1/2] CMEA2100x256FunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pMEA2100x256FunctionPointerContainer)
```

Initializes a new instance of the CMEA2100x256FunctionNet class.

```
11.60.2.2 CMEA2100x256FunctionNet() [2/2] CMEA2100x256FunctionNet ( CMcsUsbNet^{\wedge} mcsusb )
```

11.60.2.3 ~CMEA2100x256FunctionNet() virtual ~CMEA2100x256FunctionNet () [virtual]

11.60.2.4 "!CMEA2100x256FunctionNet() !CMEA2100x256FunctionNet ()

11.60.3 Member Function Documentation

```
11.60.3.1 GetLayoutConfiguration() StimulationLayoutConfigurationEnumNet GetLayoutConfiguration ()
```

Gets the stimulation layout configuration. Can be single well, 6-well or 9-well. The number of D← AC channels available per well is Mcs::Usb::CStg200xBasicNet::GetNumberOfAnalogChannels divided by Mcs::Usb::CStg200xBasicNet::GetNumberOfStimulationSourcesPerElectrode.

Returns

The currently active stimulation layout configuration.

```
11.60.3.2 SetLayoutConfiguration() void SetLayoutConfiguration (
StimulationLayoutConfigurationEnumNet LayoutConfiguration)
```


Sets the stimulation layout configuration. Can be single well, 6-well or 9-well. The number of D← AC channels available per well is Mcs::Usb::CStg200xBasicNet::GetNumberOfAnalogChannels divided by Mcs::Usb::CStg200xBasicNet::GetNumberOfStimulationSourcesPerElectrode.

Parameters

LawautCanfiguration	The new stimulation layout configuration.
i avoniconnonranon	Line new sumulation lavour confiduration
Layoutooningaration	ino non dimalation layout donngaration.

11.61 CMeaAudioFunctionNet Class Reference

Inheritance diagram for CMeaAudioFunctionNet:



Classes

· struct s setaudionet

Public Member Functions

- CMeaAudioFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] meaAudioFunction
 —
 PointerContainer)
- CMeaAudioFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual uint32_t GetNumberOfAudioChannels ()

Gets the number of available audio channels.

virtual uint32_t SetAudioChannels (array< s_setaudionet^>^ channels)

Sets the electrode to monitor and amplification for the audio channels.

- $\bullet \ \ virtual \ uint 32_t \ Set Audio Channels \ (array < s_set audionet ^> ^ \ channels, \ unsigned \ int \ virtual Device)$
 - Sets the electrode to monitor and amplification for the audio channels.
- virtual uint32_t GetAudioChannels ([System::Runtime::InteropServices::Out]array< s_setaudionet^>^% channels)

Gets the electrode to monitor and amplification for the audio channels.

virtual uint32_t GetAudioChannels ([System::Runtime::InteropServices::Out]array< s_setaudionet^>^% channels, unsigned int virtualDevice)

Gets the electrode to monitor and amplification for the audio channels.

Additional Inherited Members

11.61.1 Constructor & Destructor Documentation

11.61.2 Member Function Documentation

Gets the electrode to monitor and amplification for the audio channels.

Parameters

Returns

Error Status. 0 on success.

Gets the electrode to monitor and amplification for the audio channels.

Parameters

Parameters

virtualDevice	Virtual device to use.
---------------	------------------------

Returns

Error Status. 0 on success.

11.61.2.3 GetNumberOfAudioChannels() virtual uint32_t GetNumberOfAudioChannels () [virtual]

Gets the number of available audio channels.

Returns

The number of audio channels available, 0 when there are none.

Sets the electrode to monitor and amplification for the audio channels.

Parameters

channels	Struct which defines the electrode (channel) and amplification.	
----------	---	--

Returns

Error Status. 0 on success.

Sets the electrode to monitor and amplification for the audio channels.

Parameters

channels	Struct which defines the electrode (channel) and amplification.
----------	---

Parameters


Returns

Error Status. 0 on success.

11.62 CMeaCleanDeviceNet Class Reference

CMeaCleanDeviceNet is the class to access the MEA Clean device.

Inheritance diagram for CMeaCleanDeviceNet:



Public Member Functions

CMeaCleanDeviceNet ()

Initializes a new instance of the CMeaCleanDeviceNet class.

- virtual ~CMeaCleanDeviceNet ()
- !CMeaCleanDeviceNet ()
- void Start ()

Starts a MEA Clean run.

• void Stop ()

Stops a MEA Clean run.

void SetSlope (uint32_t voltageSlope)

Sets the voltage slope.

void SetCycles (uint32_t cycles)

Sets the number of cycles.

void SetMinVoltage (int32_t voltageMin)

Sets the lower voltage level.

void SetMaxVoltage (int32 t voltageMax)

Sets the upper voltage level.

• bool IsRunning ()

Gets if the MEA Clean device is running.

• uint32_t GetSlope ()

Gets the voltage slope.

uint32_t GetCycles ()

Gets the number of cycles.

• int32_t GetMinVoltage ()

Gets the lower voltage level.

int32_t GetMaxVoltage ()

Gets the upper voltage level

int32_t GetOutputVoltage ()

Gets the output voltage.

int32_t GetCycle ()

Gets the current cycle.

Additional Inherited Members

11.62.1 Detailed Description

CMeaCleanDeviceNet is the class to access the MEA Clean device.

11.62.2 Constructor & Destructor Documentation

11.62.2.1 CMeaCleanDeviceNet() CMeaCleanDeviceNet ()

Initializes a new instance of the CMeaCleanDeviceNet class.

```
11.62.2.2 ~CMeaCleanDeviceNet() virtual ~CMeaCleanDeviceNet () [virtual]
11.62.2.3 "!CMeaCleanDeviceNet() !CMeaCleanDeviceNet ()
11.62.3 Member Function Documentation
11.62.3.1 GetCycle() int32_t GetCycle ()
Gets the current cycle.
Returns
     The cycle number.
11.62.3.2 GetCycles() uint32_t GetCycles ()
Gets the number of cycles.
Returns
     The number of cycles to run for.
11.62.3.3 GetMaxVoltage() int32_t GetMaxVoltage ()
Gets the upper voltage level
Returns
     The upper voltage level in mV.
11.62.3.4 GetMinVoltage() int32_t GetMinVoltage ()
Gets the lower voltage level.
Returns
```

The lower voltage level in mV.

11.62.3.5 GetOutputVoltage() int32_t GetOutputVoltage ()

Gets the output voltage.

Returns

The output voltage in mV.

11.62.3.6 GetSlope() uint32_t GetSlope ()

Gets the voltage slope.

Returns

The voltage slope in mV/s.

11.62.3.7 IsRunning() bool IsRunning ()

Gets if the MEA Clean device is running.

Returns

"true" when a run is in progress, otherwise "false".

```
11.62.3.8 SetCycles() void SetCycles (
uint32_t cycles)
```

Sets the number of cycles.

Parameters

cycles The number of cycles to run for (0 .. 99).

Sets the upper voltage level.

Parameters

voltaneMay	The upper voltage level in mV ((-16 16 V/)
vuitautiviax	THE UDDEL VUILAGE IEVEL III III V (-1.0 1.0 v /.

Sets the lower voltage level.

Parameters

In The lower voltage level in mV (-1.6 1.6 \	Min The lower voltage level in mV (-1.6 1.6 V).	mV (-1.6 1	1.6 1.6 V)
--	---	------------	------------

```
11.62.3.11 SetSlope() void SetSlope ( uint32_t voltageSlope )
```

Sets the voltage slope.

Parameters

voltageSlope	The voltage slope in mV/s (range 0 60 V/s).

```
11.62.3.12 Start() void Start ()
```

Starts a MEA Clean run.


```
11.62.3.13 Stop() void Stop ()
```

Stops a MEA Clean run.

11.63 CMeaCoatDeviceNet Class Reference

CMeaCoatDeviceNet is the class to access the MEA Coat device.

Inheritance diagram for CMeaCoatDeviceNet:



Public Member Functions

CMeaCoatDeviceNet ()

Initializes a new instance of the CMeaCoatDeviceNet class.

- virtual ∼CMeaCoatDeviceNet ()
- !CMeaCoatDeviceNet ()
- void Start ()

Starts a MEA Coat run.

• void Stop ()

Stops a MEA Coat run.

void SetSlope (int32_t currentSlope)

Sets the current slope.

void SetDuration (uint32_t duration)

Sets the duration of a MEA Coat run.

void SetMaxCurrent (uint32_t currentMax)

Sets the limit of the current ramp (absolute value).

void SetOffsetCurrent (int32 t currentOffset)

Sets the offset of the current.

• bool IsRunning ()

Gets if the MEA Clean device is running.

int32 t GetSlope ()

Gets the current slope.

uint32_t GetDuration ()

Gets the duration of a MEA Coat run.

uint32_t GetMaxCurrent ()

Gets the limit of the current ramp (absolute value).

• int32_t GetOffsetCurrent ()

Gets the offset of the current.

int32_t GetOutputCurrent ()

Gets the output current.

• int32_t GetTimeInPlateau ()

Gets the time in the plateau.

void SetPauseDuration (uint32_t pauseDuration)

Sets the duration of the pause between MEA Coat pulses.

uint32_t GetPauseDuration ()

Gets the duration of the pause between MEA Coat pulses.

• int32_t GetTimeInPause ()

Gets the time in the pause.

void SetCycles (uint32_t cycles)

Sets the number of cycles.

uint32_t GetCycles ()

Gets the number of cycles.

• int32_t GetCurrentCycle ()

Gets the current cycle.

Additional Inherited Members

11.63.1 Detailed Description

CMeaCoatDeviceNet is the class to access the MEA Coat device.

11.63.2 Constructor & Destructor Documentation

11.63.2.1 CMeaCoatDeviceNet() CMeaCoatDeviceNet () Initializes a new instance of the CMeaCoatDeviceNet class. 11.63.2.2 ~CMeaCoatDeviceNet() virtual ~CMeaCoatDeviceNet () [virtual] 11.63.2.3 "!CMeaCoatDeviceNet() !CMeaCoatDeviceNet () 11.63.3 Member Function Documentation 11.63.3.1 GetCurrentCycle() int32_t GetCurrentCycle ()

Gets the current cycle.

Returns

The cycle number.

```
11.63.3.2 GetCycles() uint32_t GetCycles ()
```

Gets the number of cycles.

Returns

The number of cycles to run for.

11.63.3.3 GetDuration() uint32_t GetDuration ()

Gets the duration of a MEA Coat run.

Returns

The duration in ms.

11.63.3.4 GetMaxCurrent() uint32_t GetMaxCurrent ()

Gets the limit of the current ramp (absolute value).

Returns

The limit of the current ramp in pA (absolute value).

11.63.3.5 GetOffsetCurrent() int32_t GetOffsetCurrent ()

Gets the offset of the current.

Returns

The offset of the current in pA.

11.63.3.6 GetOutputCurrent() int32_t GetOutputCurrent ()

Gets the output current.

Returns

The output current in pA.

11.63.3.7 GetPauseDuration() uint32_t GetPauseDuration ()

Gets the duration of the pause between MEA Coat pulses.

Returns

The duration in ms.

11.63.3.8 GetSlope() int32_t GetSlope ()

Gets the current slope.

Returns

The current slope in pA/s.

11.63.3.9 GetTimeInPause() int32_t GetTimeInPause ()

Gets the time in the pause.

Returns

The time in the pause in ms.

11.63.3.10 GetTimeInPlateau() int32_t GetTimeInPlateau ()

Gets the time in the plateau.

Returns

The time in the plateau in ms.

11.63.3.11 IsRunning() bool IsRunning ()

Gets if the MEA Clean device is running.

Returns

"true" when a run is in progress, otherwise "false".

```
11.63.3.12 SetCycles() void SetCycles (
uint32_t cycles)
```

Sets the number of cycles.

Parameters

cycles The number of cycles to run for (0 .. 99).

```
11.63.3.13 SetDuration() void SetDuration ( uint32_t duration)
```

Sets the duration of a MEA Coat run.

n in ms (range 0 65 s).

```
11.63.3.14 SetMaxCurrent() void SetMaxCurrent ( uint32_t currentMax )
```

Sets the limit of the current ramp (absolute value).

Parameters

currentMax | The limit of the current ramp in pA (absolute value, 0 .. 18 nA).

Sets the offset of the current.

Parameters

```
currentOffset The offset of the current in pA (-10 .. 10 nA).
```

```
11.63.3.16 SetPauseDuration() void SetPauseDuration ( uint32_t pauseDuration )
```

Sets the duration of the pause between MEA Coat pulses.

Parameters

```
pauseDuration The duration in ms (range 0 .. 65 s).
```

```
11.63.3.17 SetSlope() void SetSlope ( int32_t currentSlope )
```

Sets the current slope.

Parameters

```
currentSlope The current slope in pA/s (range -65 .. 65 nA/s).
```

11.63.3.18 Start() void Start ()

Starts a MEA Coat run.


11.63.3.19 Stop() void Stop ()

Stops a MEA Coat run.

11.64 CMeaDeviceNet Class Reference

Base class for MEA data acquisition devices.

Inheritance diagram for CMeaDeviceNet:



Public Member Functions

CMeaDeviceNet (McsBusTypeEnumNet bustype)

Initializes a new instance of CMeaDeviceNet class.

- CMeaDeviceNet (McsBusTypeEnumNet bustype, OnChannelData^ channelData, OnError^ error)
 - Initializes a new instance of CMeaDeviceNet class.
- ∼CMeaDeviceNet ()
- virtual int32_t GetGain ()

Gets the amplifier gain of the device.

- int32_t GetEnumerationSpeed ()
- virtual int32 t GetAnalogGain ()

Gets the gain of the analog inputs of the device.

• virtual uint32_t EnableDigitalIn (bool enable, unsigned int virtualDevice)

Enable the digital data word in the datastream.

- virtual uint32_t EnableDigitalIn (DigitalDatastreamEnableEnumNet enable, unsigned int virtualDevice)
 - Enable digital data words in the datastream.
- virtual uint32_t EnableTimestamp (bool enable, unsigned int virtualDevice)

Enable the timestamp data word in the datastream. The timestamp is a 64 bit counter.

virtual uint32 t EnableChecksum (bool enable, unsigned int virtualDevice)

Enable the checksum data word in the datastream. The checksum is a 32 bit counter and 2x16 bit magic numbers.

virtual void SetDigitalOut (unsigned int digout_value, int pulselength)

Generate a pulse on the digital output.

• virtual uint32_t SetNumberOfChannels (int NumberOfChannels)

Sets the number of analog channels in the datastream.

• virtual uint32_t SetNumberOfChannels (int NumberOfChannels, unsigned int virtualDevice)

Sets the number of analog channels in the datastream.

virtual uint32_t SetNumberOfAnalogChannels (unsigned int NumberOfChannels_HS1, unsigned int NumberOfChannels_HS2, unsigned int NumberOfChannels_DSP, unsigned int NumberOfChannels_IF, unsigned int virtualDevice)

Sets the number of analog channels in the datastream for the MEA2100 device.

- virtual uint32_t SetTriggerPeriod (int samples, unsigned int virtualDevice)
 - Sets the maximum number of samples per trigger.
- virtual uint32_t SetTriggerMaskValue (unsigned int mask, unsigned int value, unsigned int virtualDevice)

Defines a pattern on the digital dataword which will start a trigger when found.

Properties

- CWClassicFunctionNet[^] WClassicFunctionNet [get]
- CW2100 FunctionNet[^] W2100 FunctionNet [get]
- CMeaAudioFunctionNet[^] MeaAudioFunctionNet [get]
- CMeaDigitalDataFunctionNet[^] MeaDigitalDataFunctionNet [get]
- CMeaFeedbackFunctionNet[^] MeaFeedbackFunctionNet [get]
- virtual int Gain [get]

The amplifier gain of the device. Value is gain times 1000, a value of 1000 corresponds to a gain of 1.0.

virtual int AnalogGain [get]

The gain of the analog inputs of the device. Value is gain times 1000, a value of 1000 corresponds to a gain of 1.0.

Additional Inherited Members

11.64.1 Detailed Description

Base class for MEA data acquisition devices.

There are two different device types for MEA data aquistion devices. There are the USB-MEA devices and the MC← _Card. In .NET both classes can be accessed by the contructor of the base class CMeaDeviceNet, which contructs the correct underlying C++ class for the USB-MEA device on the one hand or the MC_Card device on the other hand. Through this interface both device types USB-MEA devices and MC Card devices can be accessed

11.64.2 Constructor & Destructor Documentation

```
11.64.2.1 CMeaDeviceNet() [1/2] CMeaDeviceNet (
McsBusTypeEnumNet bustype)
```

Initializes a new instance of CMeaDeviceNet class.

Parameters

bustype Type of device to use, either USB or PCI.

11.64.2.2 CMeaDeviceNet() [2/2] CMeaDeviceNet (

```
McsBusTypeEnumNet bustype,
OnChannelData^ channelData,
OnError^ error )
```

Initializes a new instance of CMeaDeviceNet class.

Parameters

bustype Type of device to use, eithe	r USB or PCI.
--------------------------------------	---------------

Parameters

Parameters

```
error | Callback to call when an error occurred.
```

11.64.2.3 \sim CMeaDeviceNet() \sim CMeaDeviceNet ()

11.64.3 Member Function Documentation

Enable the checksum data word in the datastream. The checksum is a 32 bit counter and 2x16 bit magic numbers.

Parameters

enable	True to enable, False to disable.
virtualDevice	virtual device to use.

Returns

Error Status. 0 on success.

Enable the digital data word in the datastream.

Parameters

enable	True to enable, False to disable.
virtualDevice	virtual device to use.

Returns

Error Status. 0 on success.

Enable digital data words in the datastream.

Parameters

enable	True to enable, False to disable.
virtualDevice	virtual device to use.

Returns

Error Status. 0 on success.

Enable the timestamp data word in the datastream. The timestamp is a 64 bit counter.

enable	True to enable, False to disable.
virtualDevice	virtual device to use.

Returns

Error Status. 0 on success.

11.64.3.5 GetAnalogGain() virtual int32_t GetAnalogGain () [virtual]

Gets the gain of the analog inputs of the device.

Returns

Gain times 1000, a value of 1000 corresponds to a gain of 1.0.

11.64.3.6 GetEnumerationSpeed() int32_t GetEnumerationSpeed ()

```
11.64.3.7 GetGain() virtual int32_t GetGain ( ) [virtual]
```

Gets the amplifier gain of the device.

Returns

Gain times 1000, a value of 1000 corresponds to a gain of 1.0.

Generate a pulse on the digital output.

Parameters

digout_value	Bitmask to set on the digital out.
--------------	------------------------------------

11.64.3.9 SetNumberOfAnalogChannels() virtual uint32_t SetNumberOfAnalogChannels (

```
unsigned int NumberOfChannels_HS1,
unsigned int NumberOfChannels_HS2,
unsigned int NumberOfChannels_DSP,
unsigned int NumberOfChannels_IF,
unsigned int virtualDevice ) [virtual]
```

Sets the number of analog channels in the datastream for the MEA2100 device.

Parameters

NumberOfChannels_HS1	Number of analog channels from the Headstage 1.
----------------------	---

Parameters

NumberOfChannels_HS2

Parameters

NumberOfChannels_DSP	Number of data words from the DSP.
----------------------	------------------------------------

Parameters

NumberOfChannels↔	Number of analog channels from the Interfaceboard.
_IF	

Parameters

virtualDevice virtualDevice to use.

Returns

Error Status. 0 on success.

11.64.3.10 SetNumberOfChannels() [1/2] virtual uint32_t SetNumberOfChannels (int NumberOfChannels) [virtual]

Sets the number of analog channels in the datastream.

Parameters

NumberOfChannels	Number of analog channels.
------------------	----------------------------

Returns

Error Status. 0 on success.

11.64.3.11 SetNumberOfChannels() [2/2] virtual uint32_t SetNumberOfChannels (int NumberOfChannels, unsigned int virtualDevice) [virtual]

Sets the number of analog channels in the datastream.

Parameters

NumberOfChannels	Number of analog channels.
virtualDevice	virtual device to use.

Returns

Error Status. 0 on success.

```
11.64.3.12 SetTriggerMaskValue() virtual uint32_t SetTriggerMaskValue (
    unsigned int mask,
    unsigned int value,
    unsigned int virtualDevice ) [virtual]
```

Defines a pattern on the digital dataword which will start a trigger when found.

Parameters

mask	Bits in the digital dataword which are monitored for a match with value.
------	--

value	Pattern which must match for the trigger to start.
-------	--

Returns

Error Status. 0 on success.

Sets the maximum number of samples per trigger.

Parameters

samples	Number of samples to acquire after the trigger condition is met.

Returns

Error Status. 0 on success.

11.64.4 Property Documentation

11.64.4.1 AnalogGain virtual int AnalogGain [get]

The gain of the analog inputs of the device. Value is gain times 1000, a value of 1000 corresponds to a gain of 1.0.

```
11.64.4.2 Gain virtual int Gain [get]
```

The amplifier gain of the device. Value is gain times 1000, a value of 1000 corresponds to a gain of 1.0.

11.64.4.3 MeaAudioFunctionNet CMeaAudioFunctionNet^ MeaAudioFunctionNet [get]

11.64.4.4 MeaDigitalDataFunctionNet CMeaDigitalDataFunctionNet^ MeaDigitalDataFunctionNet [get]

```
11.64.4.5 MeaFeedbackFunctionNet CMeaFeedbackFunctionNet^ MeaFeedbackFunctionNet [get]
```


```
11.64.4.6 MeFunctionNet CMeFunctionNet^ MeFunctionNet [get]
```

```
11.64.4.7 W2100_FunctionNet CW2100_FunctionNet^ W2100_FunctionNet [get]
```

11.64.4.8 WClassicFunctionNet CWClassicFunctionNet^ WClassicFunctionNet [get]

11.65 CMeaDigitalDataFunctionNet Class Reference

Inheritance diagram for CMeaDigitalDataFunctionNet:



Public Member Functions

- CMeaDigitalDataFunctionNet (CMcsUsbNet[^] mcsusb)
- void SetDigitalData (unsigned int digital_value, unsigned int digital_value_mask)

Generate a value on the digital output.

• void SetDigitalData (unsigned int bit_number, bool value)

Generate a value on the digital output.

unsigned int GetDigitalData ()

Get the value of the digital output.

Additional Inherited Members

11.65.1 Constructor & Destructor Documentation

```
11.65.1.1 CMeaDigitalDataFunctionNet() [1/2] CMeaDigitalDataFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ meaDigitalFunctionPointerContainer)
```

11.65.1.2 CMeaDigitalDataFunctionNet() [2/2] CMeaDigitalDataFunctionNet (CMcsUsbNet^ mcsusb)

11.65.2 Member Function Documentation

11.65.2.1 GetDigitalData() unsigned int GetDigitalData ()

Get the value of the digital output.

Returns

Value on the digital data register.

Generate a value on the digital output.

Parameters

```
bit_number Bit number to change.
```

Parameters

```
value Bit value.
```


Generate a value on the digital output.


digital value	Value to get
alallal value	value to set.

Parameters

digital_value_mask	Mask for change.
--------------------	------------------

11.66 CMeaFeedbackFunctionNet Class Reference

Inheritance diagram for CMeaFeedbackFunctionNet:



Public Member Functions

- CMeaFeedbackFunctionNet (CMcsUsbNet[∧] mcsusb, CMcsUsbFunctionPointerContainer[∧] meaFeedback
 FunctionNet)
- CMeaFeedbackFunctionNet (CMcsUsbNet[∧] mcsusb)
- void FeedbackSetFeedback (unsigned char on, unsigned short digoutmask, unsigned short diginmask)
- unsigned int FeedbackGetSampleTimerCount ([System::Runtime::InteropServices::Out]unsigned int% CurrentCount, [System::Runtime::InteropServices::Out]unsigned int% LastKnownCount, [System::Runtime::InteropServices::Out]bool% On)
- void FeedbackSetDigitalMapping (unsigned short channel, unsigned short outmapping, unsigned short inmapping)
- void FeedbackSetFilterParameter (unsigned char filter, array< short >[∧] parameters)
- void FeedbackSetFilterParameter32 (unsigned char filter, array< int >^ parameters)
- void FeedbackSetIIRFilterParameter (unsigned char filter, int length, array< double >^ parameters)
- void FeedbackSetMkFilter (unsigned char filter, String[^] filtertype, double cheb_ribble, String[^] passtype, int order, double alpha1, double alpha2)
- void FeedbackSetChannelFilter (short channel, char filter)
- void FeedbackSetGlobalChannelFilter (char filter, unsigned short firstchannel, unsigned short lastchannel)
- void FeedbackSetFilterOff ()
- void FeedbackSetNumberOfSpikeDetectors (unsigned short number)
- void FeedbackSetSpikeDetectorThreshold (unsigned short position, unsigned short sourcechannel, unsigned short resultchannel, unsigned short trigger, unsigned short totzeit, int threshold1, int threshold2, short slope)
- void FeedbackSetNumberOfRateCounter (unsigned short number)
- void FeedbackSetRateCounter (unsigned short position, unsigned short sourcechannel, unsigned short resultchannel)
- void FeedbackSetNumberOfRateDetectors (unsigned short number)
- void FeedbackSetRateDetector (unsigned short position, unsigned short resultchannel, unsigned short trigger, unsigned short totzeit, unsigned short pulses, unsigned int duration1, unsigned int duration2)
- void FeedbackSetNumberOfLogics (unsigned short number)
- void FeedbackSetLogic (unsigned short position, array< unsigned short >^ sourcechannel, unsigned short resultchannel, unsigned int lookup)
- void FeedbackSetNumberOfTriggers (unsigned short number)
- void FeedbackSetTrigger (unsigned short position, unsigned short sourcechannel, unsigned short resultchannel, unsigned short trigger, unsigned short totzeit)
- void FeedbackSetAnalogSource (AnalogSourceEnumNet AnalogSource, unsigned int Channels, unsigned int Offset)

Additional Inherited Members

11.66.1 Constructor & Destructor Documentation

```
\textbf{11.66.1.1} \quad \textbf{CMeaFeedbackFunctionNet() [1/2]} \quad \texttt{CMeaFeedbackFunctionNet ()}
             CMcsUsbNet^ mcsusb,
             CMcsUsbFunctionPointerContainer<sup>∧</sup> meaFeedbackFunctionNet )
11.66.1.2 CMeaFeedbackFunctionNet() [2/2] CMeaFeedbackFunctionNet (
             CMcsUsbNet^ mcsusb )
11.66.2 Member Function Documentation
11.66.2.1 FeedbackGetSampleTimerCount() unsigned int FeedbackGetSampleTimerCount (
              [System::Runtime::InteropServices::Out] unsigned int% CurrentCount,
              [System::Runtime::InteropServices::Out] unsigned int% LastKnownCount,
              [System::Runtime::InteropServices::Out] bool% On )
11.66.2.2 FeedbackSetAnalogSource() void FeedbackSetAnalogSource (
             AnalogSourceEnumNet AnalogSource,
             unsigned int Channels,
             unsigned int Offset )
11.66.2.3 FeedbackSetChannelFilter() void FeedbackSetChannelFilter (
             short channel,
             char filter )
11.66.2.4 FeedbackSetDigitalMapping() void FeedbackSetDigitalMapping (
             unsigned short channel,
             unsigned short outmapping,
             unsigned short inmapping)
```

```
11.66.2.5 FeedbackSetFeedback() void FeedbackSetFeedback (
             unsigned char on,
             unsigned short digoutmask,
             unsigned short diginmask )
11.66.2.6 FeedbackSetFilterOff() void FeedbackSetFilterOff ()
11.66.2.7 FeedbackSetFilterParameter() void FeedbackSetFilterParameter (
             unsigned char filter,
             array < short >^{\wedge} parameters)
11.66.2.8 FeedbackSetFilterParameter32() void FeedbackSetFilterParameter32 (
             unsigned char filter,
             array< int >^{\wedge} parameters )
11.66.2.9 FeedbackSetGlobalChannelFilter() void FeedbackSetGlobalChannelFilter (
             char filter,
             unsigned short firstchannel,
             unsigned short lastchannel )
11.66.2.10 FeedbackSetIIRFilterParameter() void FeedbackSetIIRFilterParameter (
             unsigned char filter,
             int length,
             array< double >^{\land} parameters )
11.66.2.11 FeedbackSetLogic() void FeedbackSetLogic (
             unsigned short position,
             array< unsigned short >^{\wedge} sourcechannel,
             unsigned short resultchannel,
             unsigned int lookup )
```

```
11.66.2.12 FeedbackSetMkFilter() void FeedbackSetMkFilter (
               unsigned char filter,
               String^ filtertype,
               double cheb_ribble,
               String^{\wedge} passtype,
               int order,
               double alpha1,
               double alpha2 )
11.66.2.13 FeedbackSetNumberOfLogics() void FeedbackSetNumberOfLogics (
               unsigned short number )
11.66.2.14 FeedbackSetNumberOfRateCounter() void FeedbackSetNumberOfRateCounter (
               unsigned short number )
\textbf{11.66.2.15} \quad \textbf{FeedbackSetNumberOfRateDetectors()} \quad \texttt{void} \quad \texttt{FeedbackSetNumberOfRateDetectors} \quad \textbf{(}
               unsigned short number )
\textbf{11.66.2.16} \quad \textbf{FeedbackSetNumberOfSpikeDetectors()} \quad \texttt{void FeedbackSetNumberOfSpikeDetectors} \quad \textbf{(}
               unsigned short number )
11.66.2.17 FeedbackSetNumberOfTriggers() void FeedbackSetNumberOfTriggers (
               unsigned short number )
\textbf{11.66.2.18} \quad \textbf{FeedbackSetRateCounter()} \quad \texttt{void FeedbackSetRateCounter} \quad \textbf{(}
               unsigned short position,
               unsigned short sourcechannel,
               unsigned short resultchannel )
11.66.2.19 FeedbackSetRateDetector() void FeedbackSetRateDetector (
               unsigned short position,
               unsigned short resultchannel,
               unsigned short trigger,
               unsigned short totzeit,
               unsigned short pulses,
               unsigned int duration1,
               unsigned int duration2 )
```

11.66.2.20 FeedbackSetSpikeDetectorThreshold() void FeedbackSetSpikeDetectorThreshold (


```
unsigned short position,
unsigned short sourcechannel,
unsigned short resultchannel,
unsigned short trigger,
unsigned short totzeit,
int threshold1,
int threshold2,
short slope )
```

$\textbf{11.66.2.21} \quad \textbf{FeedbackSetTrigger()} \quad \texttt{void} \ \texttt{FeedbackSetTrigger} \ \ \textbf{(}$

```
unsigned short position,
unsigned short sourcechannel,
unsigned short resultchannel,
unsigned short trigger,
unsigned short totzeit)
```

11.67 CMealmpedanceDeviceNet Class Reference

Inheritance diagram for CMealmpedanceDeviceNet:



Public Member Functions

- CMealmpedanceDeviceNet ()
- ∼CMeaImpedanceDeviceNet ()
- virtual void StartMeasurement (unsigned short channel)
- virtual unsigned short GetReady ()
- virtual unsigned short GetArraySize ()
- virtual array< unsigned short > ^ GetResult ()
- unsigned short GetAdapterCode ()
- virtual unsigned int GetImpedanceTestFrequency ()
- virtual void SetImpedanceTestFrequency (unsigned int TestFrequency_Hertz)

Additional Inherited Members

11.67.1 Constructor & Destructor Documentation


11.67.1.1 CMealmpedanceDeviceNet() CMealmpedanceDeviceNet ()

```
11.67.1.2 ~ CMealmpedanceDeviceNet() ~ CMealmpedanceDeviceNet ()
11.67.2 Member Function Documentation
11.67.2.1 GetAdapterCode() unsigned short GetAdapterCode ( )
11.67.2.2 GetArraySize() virtual unsigned short GetArraySize ( ) [virtual]
11.67.2.3 GetImpedanceTestFrequency() virtual unsigned int GetImpedanceTestFrequency ( )
[virtual]
11.67.2.4 GetReady() virtual unsigned short GetReady ( ) [virtual]
11.67.2.5 GetResult() virtual array<unsigned short> ^ GetResult ( ) [virtual]
11.67.2.6 SetImpedanceTestFrequency() virtual void SetImpedanceTestFrequency (
             unsigned int TestFrequency_Hertz ) [virtual]
11.67.2.7 StartMeasurement() virtual void StartMeasurement (
             {\tt unsigned \ short \ } {\it channel \ }) \quad [{\tt virtual}]
```

11.68 CMeasureTableDeviceNet Class Reference

CMeasureTableDeviceNet is the to control the MCS HLA device

Inheritance diagram for CMeasureTableDeviceNet:



Public Member Functions

CMeasureTableDeviceNet (void)

Properties

• CMcsBus_SensorNet^ Sensor [get]

Additional Inherited Members

11.68.1 Detailed Description

CMeasureTableDeviceNet is the to control the MCS HLA device

11.68.2 Constructor & Destructor Documentation

```
11.68.2.1 CMeasureTableDeviceNet() CMeasureTableDeviceNet ( void )
```


11.68.3 Property Documentation

```
11.68.3.1 Sensor CMcsBus_SensorNet^ Sensor [get]
```

11.69 CMeaSwitchDeviceNet Class Reference

The class to control the USB-MEA-Switch.

Inheritance diagram for CMeaSwitchDeviceNet:



Public Member Functions

CMeaSwitchDeviceNet ()

Constructor.

∼CMeaSwitchDeviceNet ()

Destructor.

unsigned short GetNumber ()

Gets the number of boards in the device.

array< unsigned char > ^ GetPattern ()

Gets the pattern of the switches that are currently set in the device as char array.

array< bool > ^ GetPatternBool ()

Gets the pattern of the switches that are currently set in he device as bools.

void SetPattern (array< unsigned char >^ pattern)

Sets the pattern of switches from a char array.

void SetPatternBool (array< bool >^ pattern)

Sets the pattern of switches from a.

Additional Inherited Members

11.69.1 Detailed Description

The class to control the USB-MEA-Switch.

This class controls the settings of the USB-MEA-Switch. The box has two inputs for signals from a MEA amplifier. Each of the 64 outputs can be connected to one of the MEAs at the same channel.

11.69.2 Constructor & Destructor Documentation

```
11.69.2.1 CMeaSwitchDeviceNet() CMeaSwitchDeviceNet ( )
```

Constructor.

$\textbf{11.69.2.2} \quad \sim \textbf{CMeaSwitchDeviceNet()} \quad \sim \texttt{CMeaSwitchDeviceNet()} \quad \\$

Destructor.

11.69.3 Member Function Documentation

```
 \textbf{11.69.3.1} \quad \textbf{GetNumber()} \quad \texttt{unsigned short GetNumber ()} \\
```

Gets the number of boards in the device.

The MEA-Switch are delivered with 64 or 128 channels

```
11.69.3.2 GetPattern() array<unsigned char> ^{\wedge} GetPattern ( )
```

Gets the pattern of the switches that are currently set in the device as char array.

```
11.69.3.3 GetPatternBool() array<bool> ^ GetPatternBool ( )
```

Gets the pattern of the switches that are currently set in he device as bools.

```
11.69.3.4 SetPattern() void SetPattern ( array < unsigned char >^{\wedge} pattern)
```


Sets the pattern of switches from a char array.

Sets the pattern of switches from a.

11.70 CMeaUSBDeviceNet Class Reference

Class for data acquisition via ME and MEA USB amplifiers

Inheritance diagram for CMeaUSBDeviceNet:



Public Member Functions

- CMeaUSBDeviceNet (OnChannelData[^] channelData, OnError[^] error)
 - Initializes a new instance of CMeaDeviceNet class.
- CMeaUSBDeviceNet ()

Initializes a new instance of CMeaDeviceNet class.

• \sim CMeaUSBDeviceNet ()

Additional Inherited Members

11.70.1 Detailed Description

Class for data acquisition via ME and MEA USB amplifiers

11.70.2 Constructor & Destructor Documentation

```
11.70.2.1 CMeaUSBDeviceNet() [1/2] CMeaUSBDeviceNet (
OnChannelData^ channelData,
OnError^ error )
```

Initializes a new instance of CMeaDeviceNet class.

Parameters

channelData	Handler to call when new data is available.
-------------	---

Parameters

error Handler to call when an error occurs.


11.70.2.2 CMeaUSBDeviceNet() [2/2] CMeaUSBDeviceNet ()

Initializes a new instance of CMeaDeviceNet class.

11.70.2.3 ~CMeaUSBDeviceNet() ~CMeaUSBDeviceNet ()

11.71 CMeFunctionNet Class Reference

Inheritance diagram for CMeFunctionNet:



Public Member Functions

CMeFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] meFunctionPointer
 — Container)

Initializes a new instance of the CDacCalibrationFunctionNet class.

- CMeFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ∼CMeFunctionNet (void)
- !CMeFunctionNet (void)
- void SetTransformer (unsigned int index, bool onoff)

Additional Inherited Members

11.71.1 Detailed Description

11.71.2 Constructor & Destructor Documentation

Initializes a new instance of the CDacCalibrationFunctionNet class.

```
11.71.2.2 CMeFunctionNet() [2/2] CMeFunctionNet (
CMcsUsbNet^ mcsusb )
```

```
11.71.2.3 ~CMeFunctionNet() virtual ~CMeFunctionNet (
void ) [virtual]
```


```
11.71.2.4 "!CMeFunctionNet() !CMeFunctionNet (
```

11.71.3 Member Function Documentation

11.72 CMultiBatteryChargerDeviceNet Class Reference

CMultiBatteryChargerDeviceNet is the class to access the MBC-08 device.

Inheritance diagram for CMultiBatteryChargerDeviceNet:



Public Member Functions

CMultiBatteryChargerDeviceNet ()

Initializes a new instance of the CMultiBatteryChargerDeviceNet class.

- virtual ~CMultiBatteryChargerDeviceNet ()
- !CMultiBatteryChargerDeviceNet ()
- uint32_t GetChargeCurrent (uint32_t NrChannel)

gets the charge current; unit: mA

• uint32_t GetDischargeCurrent (uint32_t NrChannel)

gets the discharge current; unit: mA

void SetDischargeCurrentSetPoint (uint32_t NrChannel, uint32_t DischargeCurrent_mA)

sets the setpoint for the discharge current; unit: mA

uint32_t GetDischargeCurrentSetPoint (uint32_t NrChannel)

gets the setpoint for the discharge current; unit: mA

• void SetFinalDischargeVoltage (uint32_t NrChannel, uint32_t FinalDischargeVoltage mV)

sets the final discharge voltage; unit: mV

uint32_t GetFinalDischargeVoltage (uint32_t NrChannel)

gets the final discharge voltage; unit: mV

uint32_t GetDischargeCapacity (uint32_t NrChannel)

gets the discharge capacity; unit: ?Ah

uint32_t GetChargeCapacity (uint32_t NrChannel)

gets the charge capacity; unit: ?Ah

uint32_t GetBatteryVoltage (uint32_t NrChannel)

gets the battery voltage; unit: mV

uint32_t GetChannels ()

gets number of channels

void SetRatedCapacityVolatile (uint32_t NrChannel, MbcRatedCapacityEnumNet NewRatedCapacity)

sets the rated capacity (i.e. charge current) without storing it persistently

void SetChargingMode (uint32 t NrChannel, MbcChargingModeEnumNet NewOperatingMode)

sets the charging mode: StorageCharge, LowCurrentCharge and HighCurrentCharge

MbcChargingModeEnumNet GetChargingMode (uint32_t NrChannel)

gets the charging mode: StorageCharge, LowCurrentCharge and HighCurrentCharge

MbcChannelStateEnumNet GetChannelState (uint32_t NrChannel)

gets the channel state: IdleNoBattery, IdleChargeFinished, CapacityTestPreCharge, CapacityTestDischarge, StorageCharge, LowCurrentCharge, HighCurrentCharge

void CapacityTest (uint32_t NrChannel)

start capacity test on channel

· void ChannelReset (uint32 t NrChannel)

cancel charging and capacity test functions; check if battery is connected

- void SetChargingPCoefficient (uint32_t pCoefficient)
 sets the p-coefficient for charging in mA/V / nominal charging current
- uint32 t GetChargingPCoefficient ()

gets the p-coefficient for charging in mA/V / nominal charging current

void SetRatedCapacity (uint32_t NrChannel, MbcRatedCapacityEnumNet NewRatedCapacity)
 sets the rated capacity

 MbcRatedCapacityEnumNet GetRatedCapacity (uint32_t NrChannel) gets the rated capacity

Additional Inherited Members

11.72.1 Detailed Description

CMultiBatteryChargerDeviceNet is the class to access the MBC-08 device.

11.72.2 Constructor & Destructor Documentation

11.72.2.1 CMultiBatteryChargerDeviceNet() CMultiBatteryChargerDeviceNet ()

Initializes a new instance of the CMultiBatteryChargerDeviceNet class.

```
11.72.2.2 ~CMultiBatteryChargerDeviceNet() virtual ~CMultiBatteryChargerDeviceNet ( ) [virtual]
```

```
11.72.2.3 "!CMultiBatteryChargerDeviceNet() !CMultiBatteryChargerDeviceNet ( )
```

11.72.3 Member Function Documentation

start capacity test on channel

Parameters

NrChannel the channel number

```
11.72.3.2 ChannelReset() void ChannelReset ( uint32_t NrChannel)
```

cancel charging and capacity test functions; check if battery is connected

Parameters

NrChannel the channel number

11.72.3.3 GetBatteryVoltage() uint32_t GetBatteryVoltage (uint32_t NrChannel)

gets the battery voltage; unit: mV

Parameters

Returns

the battery voltage in mV

11.72.3.4 GetChannels() uint32_t GetChannels ()

gets number of channels

Returns

number of channels

11.72.3.5 GetChannelState() MbcChannelStateEnumNet GetChannelState (uint32_t NrChannel)

gets the channel state: IdleNoBattery, IdleChargeFinished, CapacityTestPreCharge, CapacityTestDischarge, StorageCharge, LowCurrentCharge, HighCurrentCharge

N 01 1	
NrChannei	the channel number

Returns

the current state

gets the charge capacity; unit: ?Ah

Parameters

NrChannel	the channel number
Michallie	the charmer number

Returns

the capacity in uAh

gets the charge current; unit: mA

Parameters

NrChannel	the channel number
i ii Onamo	tilo olialilloi lialillooi

Returns

the measured charge current in mA

11.72.3.8 GetChargingMode() MbcChargingModeEnumNet GetChargingMode (uint32_t NrChannel)

gets the charging mode: StorageCharge, LowCurrentCharge and HighCurrentCharge

Parameters

NrChannel the channel number

Returns

the charging mode

11.72.3.9 GetChargingPCoefficient() uint32_t GetChargingPCoefficient ()

gets the p-coefficient for charging in mA/V / nominal charging current

Returns

the p-coefficient

11.72.3.10 GetDischargeCapacity() uint32_t GetDischargeCapacity (uint32_t NrChannel)

gets the discharge capacity; unit: ?Ah

Parameters

NrChannel the channel number

Returns

the capacity in uAh

11.72.3.11 GetDischargeCurrent() uint32_t GetDischargeCurrent (uint32_t NrChannel)

gets the discharge current; unit: mA

Parameters

NrChannel | the channel number

Returns

the measured discharge current in mA

11.72.3.12 GetDischargeCurrentSetPoint() uint32_t GetDischargeCurrentSetPoint (uint32_t NrChannel)

gets the setpoint for the discharge current; unit: mA

Parameters

NrChannel the channel number

Returns

the discharge current in mA

```
11.72.3.13 GetFinalDischargeVoltage() uint32_t GetFinalDischargeVoltage ( uint32_t NrChannel )
```

gets the final discharge voltage; unit: mV

Parameters

Returns

the battery voltage in mV at the end of discharge

11.72.3.14 **GetRatedCapacity()** MbcRatedCapacityEnumNet GetRatedCapacity (uint32_t NrChannel)

gets the rated capacity

Parameters

NrChannel	the channel number

Returns

the capacity

```
11.72.3.15 SetChargingMode() void SetChargingMode (
```

uint32_t NrChannel,
MbcChargingModeEnumNet NewOperatingMode)

sets the charging mode: StorageCharge, LowCurrentCharge and HighCurrentCharge

NrChannel	the channel number
NewOperatingMode	the charging mode

11.72.3.16 SetChargingPCoefficient() void SetChargingPCoefficient (uint32_t pCoefficient)

sets the p-coefficient for charging in mA/V / nominal charging current

Parameters

pCoefficient	the p-coefficient
--------------	-------------------

sets the setpoint for the discharge current; unit: mA

Parameters

NrChannel	the channel number
DischargeCurrent_mA	the discharge current in mA

sets the final discharge voltage; unit: mV

Parameters

NrChannel	the channel number
FinalDischargeVoltage_mV	the battery voltage in mV at the end of discharge

11.72.3.19 SetRatedCapacity() void SetRatedCapacity (uint32_t NrChannel,

 ${\tt MbcRatedCapacityEnumNet}\ \textit{NewRatedCapacity}\)$

sets the rated capacity

NrChannel	the channel number
NewRatedCapacity	the capacity

sets the rated capacity (i.e. charge current) without storing it persistently


Parameters

NrChannel	the channel number
NewRatedCapacity	the capacity

11.73 CMultiwellCallbackFunctionNet Class Reference

CMultiwellCallbackFunctionNet is the class to access the Multiwell-Mini-Stimulator

Inheritance diagram for CMultiwellCallbackFunctionNet:



Public Member Functions

- delegate void OnGetPlateClampStateByHeadstage (uint32_t Headstage, PlateClampEnumNet plateState)
- CMultiwellCallbackFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pMultiwell←
 CallbackFunctionPointerContainer)

Initializes a new instance of the CMultiwellCallbackFunctionNet class.

- CMultiwellCallbackFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CMultiwellCallbackFunctionNet ()
- !CMultiwellCallbackFunctionNet ()
- PlateClampEnumNet GetPlateClampStateByHeadstage (uint32_t Headstage)

Gets the state of the plate

Events

Event fires when the plate state for the headstage number has changed

Additional Inherited Members

11.73.1 Detailed Description

CMultiwellCallbackFunctionNet is the class to access the Multiwell-Mini-Stimulator

11.73.2 Constructor & Destructor Documentation

```
11.73.2.1 CMultiwellCallbackFunctionNet() [1/2] CMultiwellCallbackFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pMultiwellCallbackFunctionPointerContainer)
```

Initializes a new instance of the CMultiwellCallbackFunctionNet class.

```
11.73.2.2 CMultiwellCallbackFunctionNet() [2/2] CMultiwellCallbackFunctionNet ( CMcsUsbNet^ mcsusb )
```

```
11.73.2.3 ~CMultiwellCallbackFunctionNet() virtual ~CMultiwellCallbackFunctionNet ( ) [virtual]
```

```
11.73.2.4 "!CMultiwellCallbackFunctionNet() !CMultiwellCallbackFunctionNet ( )
```

11.73.3 Member Function Documentation

11.73.3.1 GetPlateClampStateByHeadstage() PlateClampEnumNet GetPlateClampStateByHeadstage (uint32_t Headstage)

Gets the state of the plate

Parameters

Headstage	The headstage number
-----------	----------------------

Returns

The plate state

11.73.4 Event Documentation


11.73.4.1 GetPlateClampStateByHeadstageEvent OnGetPlateClampStateByHeadstage[∧] GetPlateClamp↔ StateByHeadstageEvent [add], [remove], [raise]

Event fires when the plate state for the headstage number has changed

11.74 CMultiwellDeviceNet Class Reference

CMultiwellDeviceNet is the class to access the Multiwell device.

Inheritance diagram for CMultiwellDeviceNet:



Public Member Functions

• CMultiwellDeviceNet ()

Initializes a new instance of the CMultiwellDeviceNet class.

- virtual ∼CMultiwellDeviceNet ()
- !CMultiwellDeviceNet ()
- PlateClampEnumNet GetPlateClampState ()

Gets the state of the Multiwell plate clamp.

• PlateClampEnumNet GetPlateClampState (uint32_t Headstage)

Gets the state of the plate

• void OpenPlateClamp ()

Opens the plate clamp.

void ClosePlateClamp ()

Closes the plate clamp.

void StopPlateClamp ()

Stops the plate clamp movement.

• uint32 t GetPlateClampLockState ()

Gets the state of the plate clamp lock.

void LockPlateClamp ()

Locks the plate clamp.

void UnlockPlateClamp ()

Unlocks the plate clamp.

• MultiwellPlateTypeEnumNet GetPlateType ()

Gets the plate type.

MultiwellPlateTypeEnumNet GetPlateType (uint32_t Headstage)

Gets the plate type.

void SetPlateType (MultiwellPlateTypeEnumNet plateType)

Sets the plate type.

• void SetPlateType (uint32_t Headstage, MultiwellPlateTypeEnumNet plateType)

Sets the plate type.

void SetPlateMux (uint32_t muxSelection)

Selects a one quarter of the electrodes on a high density Multiwell plate.

void SetPlateMux (uint32_t Headstage, uint32_t muxSelection)

Selects a one quarter of the electrodes on a high density Multiwell plate.

uint32 t GetPlateMux ()

Gets the selected quarter of the electrodes on a high density Multiwell plate.

uint32 t GetPlateMux (uint32 t Headstage)

Gets the selected quarter of the electrodes on a high density Multiwell plate.

bool IsPlateTypeValid ()

Checks whether the plate type is valid, meaning all pins have contact.

bool IsPlateTypeValid (uint32_t Headstage)

Checks whether the plate type is valid, meaning all pins have contact.

void SetPowerMuxPlate (uint32_t Headstage, bool powerOn)

On the Multiwell Mini device, turn Power to the MUX Plate On or Off.

bool GetPowerMuxPlate (uint32_t Headstage)

On the Multiwell Mini device, Query if Power to the MUX Plate is On or Off.

Additional Inherited Members

11.74.1 Detailed Description

CMultiwellDeviceNet is the class to access the Multiwell device.

11.74.2 Constructor & Destructor Documentation

11.74.2.1 CMultiwellDeviceNet() CMultiwellDeviceNet ()

Initializes a new instance of the CMultiwellDeviceNet class.

$\textbf{11.74.2.2} \quad \sim \textbf{CMultiwellDeviceNet()} \quad \text{virtual} \quad \sim \texttt{CMultiwellDeviceNet ()} \quad [\texttt{virtual}]$

11.74.2.3 "!CMultiwellDeviceNet() !CMultiwellDeviceNet ()

11.74.3 Member Function Documentation

11.74.3.1 ClosePlateClamp() void ClosePlateClamp ()

Closes the plate clamp.

11.74.3.2 GetPlateClampLockState() uint32_t GetPlateClampLockState ()

Gets the state of the plate clamp lock.

Returns

the state of the plate lock (unlocked/locked)

11.74.3.3 GetPlateClampState() [1/2] PlateClampEnumNet GetPlateClampState ()

Gets the state of the Multiwell plate clamp.

Returns

the state of the plate clamp (open/closed)

11.74.3.4 GetPlateClampState() [2/2] PlateClampEnumNet GetPlateClampState (uint32_t Headstage)

Gets the state of the plate

Parameters

Headstage	The headstage number

Returns

The plate state

11.74.3.5 GetPlateMux() [1/2] uint32_t GetPlateMux ()

Gets the selected quarter of the electrodes on a high density Multiwell plate.

the selected quarter

```
11.74.3.6 GetPlateMux() [2/2] uint32_t GetPlateMux (
             uint32_t Headstage )
```

Gets the selected quarter of the electrodes on a high density Multiwell plate.

Parameters

Headstage	The headstage to query.
-----------	-------------------------

Returns

the selected quarter

11.74.3.7 GetPlateType() [1/2] MultiwellPlateTypeEnumNet GetPlateType ()

Gets the plate type.

Returns

the plate type

```
11.74.3.8 GetPlateType() [2/2] MultiwellPlateTypeEnumNet GetPlateType (
             uint32_t Headstage )
```

Gets the plate type.

Parameters

Headstage The headstage to query

Returns

the plate type

```
11.74.3.9 GetPowerMuxPlate() bool GetPowerMuxPlate (
```

uint32_t Headstage)

On the Multiwell Mini device, Query if Power to the MUX Plate is On or Off.

Headstage The headstage to que

Returns

"true" Power is On, "false" Power is Off

11.74.3.10 IsPlateTypeValid() [1/2] bool IsPlateTypeValid ()

Checks whether the plate type is valid, meaning all pins have contact.

Returns

"true" when all pins have contact, otherwise "false".

11.74.3.11 **IsPlateTypeValid()** [2/2] bool IsPlateTypeValid (uint32_t Headstage)

Checks whether the plate type is valid, meaning all pins have contact.

Parameters

Headstage	The headstage to query.

Returns

"true" when all pins have contact, otherwise "false".

11.74.3.12 LockPlateClamp() void LockPlateClamp ()

Locks the plate clamp.

11.74.3.13 OpenPlateClamp() void OpenPlateClamp ()

Opens the plate clamp.

Selects a one quarter of the electrodes on a high density Multiwell plate.

Headstage	The headstage to query.
muxSelection	the selected quarter

11.74.3.15 SetPlateMux() [2/2] void SetPlateMux (uint32_t muxSelection)

Selects a one quarter of the electrodes on a high density Multiwell plate.

Parameters

muxSelection	the selected quarter
--------------	----------------------

Sets the plate type.

Parameters

```
plateType the plate type
```

Sets the plate type.

Parameters

Headstage	The headstage to query.
plateType	the plate type

On the Multiwell Mini device, turn Power to the MUX Plate On or Off.

Headstage	The headstage to query.
powerOn	"true" to turn Power On, "false" to turn Power Off

11.74.3.19 StopPlateClamp() void StopPlateClamp ()

Stops the plate clamp movement.


11.74.3.20 UnlockPlateClamp() void UnlockPlateClamp ()

Unlocks the plate clamp.

11.75 CMultiwellOptoStimFunctionNet Class Reference

CMultiwellOptoStimFunctionNet is the class to access the optical properties of the Multiwell Optostim device

Inheritance diagram for CMultiwellOptoStimFunctionNet:



Public Member Functions

- CMultiwellOptoStimFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pMultiwell←
 OptoStimFunctionPointerContainer)
 - Initializes a new instance of the CMultiwellOptoStimFunctionNet class.
- CMultiwellOptoStimFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CMultiwellOptoStimFunctionNet ()
- !CMultiwellOptoStimFunctionNet ()
- uint32_t GetWaveLengthInNanometer (uint16_t channel)
- uint32_t GetAbsMaxCurrentInMicroAmp (uint16_t channel)
- uint32_t GetMaxDurationHighCurrentInMicroSec (uint16_t channel)
- uint32_t GetMaxDutyCycleHighCurrent (uint16_t channel)
- uint32_t GetPermanentCurrentInMicroAmp (uint16_t channel)
- uint32_t GetColorRgb (uint16_t channel)
- String \(^\) GetColorStr (uint16_t channel)
- void SetWaveLengthInNanometer (uint16_t channel, uint32_t WaveLength_nm)
- void SetAbsMaxCurrentInMicroAmp (uint16_t channel, uint32_t AbsoluteMaxCurrent_uA)
- void SetMaxDurationHighCurrentInMicroSec (uint16_t channel, uint32_t AbsoluteMaxDuration_us)
- void SetMaxDutyCycleHighCurrent (uint16 t channel, uint32 t MaxDutyCycleHighCurrent)
- void SetPermanentCurrentInMicroAmp (uint16_t channel, uint32_t PermanentCurrent_uA)
- void SetColorRgb (uint16_t channel, uint32_t ColorRGB)
- void SetColorStr (uint16_t channel, String[^] ColorString)

Additional Inherited Members

11.75.1 Detailed Description

CMultiwellOptoStimFunctionNet is the class to access the optical properties of the Multiwell Optostim device

11.75.2 Constructor & Destructor Documentation

```
11.75.2.1 CMultiwellOptoStimFunctionNet() [1/2] CMultiwellOptoStimFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pMultiwellOptoStimFunctionPointerContainer )
```

Initializes a new instance of the CMultiwellOptoStimFunctionNet class.

```
11.75.2.2 CMultiwellOptoStimFunctionNet() [2/2] CMultiwellOptoStimFunctionNet (
CMcsUsbNet^ mcsusb)
```

11.75.2.3 ~CMultiwellOptoStimFunctionNet() virtual ~CMultiwellOptoStimFunctionNet () [virtual]

```
11.75.2.4 "!CMultiwellOptoStimFunctionNet() !CMultiwellOptoStimFunctionNet ( )
```

11.75.3 Member Function Documentation

```
11.75.3.1 GetAbsMaxCurrentInMicroAmp() uint32_t GetAbsMaxCurrentInMicroAmp ( uint16_t channel)
```

Parameters

```
channel the (analog) channel number
```

Returns

absolute max. current; unit: uA

channel	the (analog) channel number
---------	-----------------------------

Returns

RGB-value of LED color

Parameters

channel	the (analog) channel number
---------	-----------------------------

Returns

LED color as string

11.75.3.4 **GetMaxDurationHighCurrentInMicroSec()** uint32_t GetMaxDurationHighCurrentInMicroSec (uint16_t channel)

Parameters

channel	the (analog) channel number
---------	-----------------------------

Returns

max. duration the LED can stand the abs. max current; unit: us

11.75.3.5 GetMaxDutyCycleHighCurrent() uint32_t GetMaxDutyCycleHighCurrent (uint16_t *channel*)

Parameters

channel	the (analog) channel number

max. duty cycle at max. current; unit: 100*%

$\textbf{11.75.3.6} \quad \textbf{GetPermanentCurrentInMicroAmp()} \quad \texttt{uint32_t} \quad \texttt{GetPermanentCurrentInMicroAmp} \quad \textbf{(}$

uint16_t channel)

Parameters

channel	the (analog) channel number
---------	-----------------------------

Returns

max. current the LED can stand when always switched on; unit: uA

$\textbf{11.75.3.7} \quad \textbf{GetWaveLengthInNanometer()} \quad \texttt{uint32_t GetWaveLengthInNanometer} \quad \textbf{(}$

uint16_t channel)

Parameters

channel the (analog) channel number

Returns

wavelength of this channel's LEDs; unit: nm

11.75.3.8 SetAbsMaxCurrentInMicroAmp() void SetAbsMaxCurrentInMicroAmp (

```
uint16_t channel,
uint32_t AbsoluteMaxCurrent_uA )
```

Parameters

channel	the (analog) channel number
AbsoluteMaxCurrent_uA	absolute max. current; unit: uA

11.75.3.9 SetColorRgb() void SetColorRgb (

```
uint16_t channel,
uint32_t ColorRGB )
```

channel	the (analog) channel number
ColorRGB	RGB-value of LED color

Parameters

channel	the (analog) channel number
ColorString	LED color as string

Parameters

channel	the (analog) channel number
AbsoluteMaxDuration_us	max. duration the LED can stand the abs. max current; unit: us

Parameters

channel	the (analog) channel number
MaxDutyCycleHighCurrent	max. duty cycle at max. current; unit: 100*%

Parameters


channel	the (analog) channel number
PermanentCurrent_uA	max. current the LED can stand when always switched on; unit: uA

Parameters

channel	the (analog) channel number
WaveLength_nm	wavelength of this channel's LEDs; unit: nm

11.76 CNF_GenDeviceNet Class Reference

Inheritance diagram for CNF_GenDeviceNet:



Public Member Functions

- CNF_GenDeviceNet (void)
- ∼CNF_GenDeviceNet (void)
- void Set_Values (unsigned int frequency, unsigned int amplitude)

Additional Inherited Members

11.76.1 Constructor & Destructor Documentation

```
11.76.1.1 CNF_GenDeviceNet() CNF_GenDeviceNet (

void )

11.76.1.2 ~CNF_GenDeviceNet() ~CNF_GenDeviceNet (
```


11.76.2 Member Function Documentation

void)

```
11.76.2.1 Set_Values() void Set_Values (
          unsigned int frequency,
          unsigned int amplitude )
```

11.77 COctoPotDeviceNet Class Reference

Inheritance diagram for COctoPotDeviceNet:



Public Member Functions

- COctoPotDeviceNet (void)
- COctoPotDeviceNet (OnChannelData^ channelData, OnError^ error)
- uint32_t SetOutputRate (uint32_t rate)
- uint32_t SetBathclamp (unsigned int block, bool enable)
- uint32 t SetDacValue (int channel, int value)
- uint32_t SetDacAutoControl (unsigned int channel)
- uint32_t SetPidParameter (unsigned int channel, int const_p, int const_i, int shift_p, int shift_i)
- uint32_t SetRampParameter (unsigned int channel, int start, int min, int max, int slope, int slope2, int pause, unsigned int samples)
- uint32_t RampStart (int channelmap)
- uint32_t SetSineParameter (unsigned int channel, int amplitude)
- uint32_t SineStart (int channelmap)
- uint32_t SetPatternListEntry (unsigned int channel, unsigned int position, unsigned int duration, int value)
- uint32_t PatternListStart (int channelmap)
- uint32_t SetAdcOffset (unsigned int channel, int offset)
- uint32_t SetDacOffset (unsigned int channel, int offset)
- uint32_t ResetAdcOffset (unsigned int channel)
- uint32_t ResetDacOffset (unsigned int channel)
- uint32 t BurnAdcOffset ()
- uint32_t BurnDacOffset ()
- uint32_t GetAdcOffset (unsigned int channel, [System::Runtime::InteropServices::Out] int ^ offset)
- uint32 t GetDacOffset (unsigned int channel, [System::Runtime::InteropServices::Out] int ^ offset)
- uint32 t SetAmplificationSwitch (unsigned int channel, unsigned int state)
- uint32_t SetChannelSwitch (unsigned int channel, unsigned int state)
- uint32_t SetNumberOfChannels (unsigned int NumberOfChannels)
- uint32_t EnableDigitalIn (bool enable)
- uint32_t EnableTimestamp (bool enable)
- uint32 t EnableChecksum (bool enable)

Additional Inherited Members

11.77.1 Constructor & Destructor Documentation


```
11.77.1.1 COctoPotDeviceNet() [1/2] COctoPotDeviceNet (
             void )
11.77.1.2 COctoPotDeviceNet() [2/2] COctoPotDeviceNet (
             OnChannelData^ channelData,
             OnError<sup>∧</sup> error )
11.77.2 Member Function Documentation
11.77.2.1 BurnAdcOffset() uint32_t BurnAdcOffset ( )
11.77.2.2 BurnDacOffset() uint32_t BurnDacOffset ( )
11.77.2.3 EnableChecksum() uint32_t EnableChecksum (
             bool enable )
11.77.2.4 EnableDigitalIn() uint32_t EnableDigitalIn (
             bool enable )
11.77.2.5 EnableTimestamp() uint32_t EnableTimestamp (
             bool enable )
11.77.2.6 GetAdcOffset() uint32_t GetAdcOffset (
             unsigned int channel,
             [System::Runtime::InteropServices::Out] int ^{\land} offset )
11.77.2.7 GetDacOffset() uint32_t GetDacOffset (
             unsigned int channel,
             [System::Runtime::InteropServices::Out] int ^{\land} offset )
```

```
11.77.2.8 PatternListStart() uint32_t PatternListStart (
              int channelmap )
11.77.2.9 RampStart() uint32_t RampStart (
              int channelmap )
\textbf{11.77.2.10} \quad \textbf{ResetAdcOffset()} \quad \texttt{uint32\_t} \;\; \texttt{ResetAdcOffset} \;\; \texttt{(}
              unsigned int channel )
11.77.2.11 ResetDacOffset() uint32_t ResetDacOffset (
              unsigned int channel )
11.77.2.12 SetAdcOffset() uint32_t SetAdcOffset (
              unsigned int channel,
              int offset )
11.77.2.13 SetAmplificationSwitch() uint32_t SetAmplificationSwitch (
              unsigned int channel,
              unsigned int state )
11.77.2.14 SetBathclamp() uint32_t SetBathclamp (
              unsigned int block,
              bool enable )
11.77.2.15 SetChannelSwitch() uint32_t SetChannelSwitch (
              unsigned int channel,
              unsigned int state )
11.77.2.16 SetDacAutoControl() uint32_t SetDacAutoControl (
              unsigned int channel )
```

```
11.77.2.17 SetDacOffset() uint32_t SetDacOffset (
             unsigned int channel,
             int offset )
11.77.2.18 SetDacValue() uint32_t SetDacValue (
             int channel,
             int value )
11.77.2.19 SetNumberOfChannels() uint32_t SetNumberOfChannels (
             unsigned int NumberOfChannels )
11.77.2.20 SetOutputRate() uint32_t SetOutputRate (
             uint32_t rate )
11.77.2.21 SetPatternListEntry() uint32_t SetPatternListEntry (
             unsigned int channel,
             unsigned int position,
             unsigned int duration,
             int value )
11.77.2.22 SetPidParameter() uint32_t SetPidParameter (
             unsigned int channel,
             int const_p,
             int const_i,
             int shift_p,
             int shift_i )
11.77.2.23 SetRampParameter() uint32_t SetRampParameter (
             unsigned int channel,
             int start,
             int min,
             int max,
             int slope,
             int slope2,
             int pause,
             unsigned int samples )
```

11.78 COkuvisionStimulatorDeviceNet Class Reference

Inheritance diagram for COkuvisionStimulatorDeviceNet:



Public Member Functions

- · COkuvisionStimulatorDeviceNet (void)
- ~COkuvisionStimulatorDeviceNet (void)
- void SetPulseform (int channel, int current, int pulsewidth, int periode, int duration)
- void GetPulseform (int channel, [System::Runtime::InteropServices::Out] int% current, [System::Runtime ← ::InteropServices::Out] int% pulsewidth, [System::Runtime::InteropServices::Out] int% periode, [System::← Runtime::InteropServices::Out] int% duration)
- void SetMaxPower (int channel, int power)
- int GetMaxPower (int channel)
- void SetMaxVoltage (int channel, int voltage)
- int GetMaxVoltage (int channel)
- void SetCheckVoltage (int channel, int voltage)
- int GetCheckVoltage (int channel)
- int GetVoltage (int channel)
- · void SetDACOffset (int channel, int part, int offset)
- int GetDACOffset (int channel, int part)
- void SetRTC (uint8_t year, uint8_t month, uint8_t day, uint8_t hour, uint8_t minute, uint8_t second)
- void GetRTC ([System::Runtime::InteropServices::Out] uint8_t% year, [System::Runtime::InteropServices
 ::Out] uint8_t% month, [System::Runtime::InteropServices::Out] uint8_t% day, [System::Runtime::Interop
 Services::Out] uint8_t% hour, [System::Runtime::InteropServices::Out] uint8_t% minute, [System::Runtime
 ::InteropServices::Out] uint8_t% second)
- void SetRTC (DateTime timestamp)
- DateTime GetRTC ()
- void GetStimulatorStatus ([System::Runtime::InteropServices::Out] int% startstop, [System::Runtime::

 InteropServices::Out] int% last_error, [System::Runtime::InteropServices::Out] int% battery_status)
- void SetCurrentFactor (int channel, int factor)
- · int GetCurrentFactor (int channel)

Additional Inherited Members

11.78.1 Constructor & Destructor Documentation

```
11.78.1.1 COkuvisionStimulatorDeviceNet() COkuvisionStimulatorDeviceNet (
            void )
11.78.1.2 ~COkuvisionStimulatorDeviceNet() ~COkuvisionStimulatorDeviceNet (
             void )
11.78.2 Member Function Documentation
11.78.2.1 GetCheckVoltage() int GetCheckVoltage (
             int channel )
11.78.2.2 GetCurrentFactor() int GetCurrentFactor (
             int channel )
11.78.2.3 GetDACOffset() int GetDACOffset (
            int channel,
             int part )
11.78.2.4 GetMaxPower() int GetMaxPower (
            int channel )
11.78.2.5 GetMaxVoltage() int GetMaxVoltage (
             int channel )
```


```
\textbf{11.78.2.6} \quad \textbf{GetPulseform()} \quad \texttt{void GetPulseform ()}
              int channel,
              [System::Runtime::InteropServices::Out] int% current,
              [System::Runtime::InteropServices::Out] int% pulsewidth,
              [System::Runtime::InteropServices::Out] int% periode,
              [System::Runtime::InteropServices::Out] int% duration )
11.78.2.7 GetRTC() [1/2] DateTime GetRTC ( )
11.78.2.8 GetRTC() [2/2] void GetRTC (
              [System::Runtime::InteropServices::Out] uint8_t% year,
              [System::Runtime::InteropServices::Out] uint8_t% month,
              [System::Runtime::InteropServices::Out] uint8_t% day,
              [System::Runtime::InteropServices::Out] uint8_t% hour,
              [System::Runtime::InteropServices::Out] uint8_t% minute,
              [System::Runtime::InteropServices::Out] uint8_t% second )
11.78.2.9 GetStimulatorStatus() void GetStimulatorStatus (
              [System::Runtime::InteropServices::Out] int% startstop,
              [{\tt System::Runtime::InteropServices::Out}] \  \, {\tt int} \\ \textit{ last\_error,}
              [System::Runtime::InteropServices::Out] int% battery_status )
11.78.2.10 GetVoltage() int GetVoltage (
              int channel )
11.78.2.11 SetCheckVoltage() void SetCheckVoltage (
              int channel,
              int voltage )
11.78.2.12 SetCurrentFactor() void SetCurrentFactor (
              int channel,
              int factor )
```

```
\textbf{11.78.2.13} \quad \textbf{SetDACOffset()} \quad \texttt{void SetDACOffset ()}
              int channel,
              int part,
              int offset )
11.78.2.14 SetMaxPower() void SetMaxPower (
              int channel,
              int power )
11.78.2.15 SetMaxVoltage() void SetMaxVoltage (
              int channel,
              int voltage )
11.78.2.16 SetPulseform() void SetPulseform (
              int channel,
              int current,
              int pulsewidth,
              int periode,
              int duration )
11.78.2.17 SetRTC() [1/2] void SetRTC (
              DateTime timestamp )
11.78.2.18 SetRTC() [2/2] void SetRTC (
              uint8_t year,
              uint8_t month,
              uint8_t day,
              uint8_t hour,
              uint8_t minute,
              uint8_t second )
```

11.79 CPatchServerDeviceNet Class Reference

CPatchServerDeviceNet is the class to control the MCS PatchServer device

Inheritance diagram for CPatchServerDeviceNet:



Public Member Functions

CPatchServerDeviceNet (void)

Properties

• CMcsBus_SensorNet^ Sensor [get]

Additional Inherited Members

11.79.1 Detailed Description

CPatchServerDeviceNet is the class to control the MCS PatchServer device

11.79.2 Constructor & Destructor Documentation


```
11.79.2.1 CPatchServerDeviceNet() CPatchServerDeviceNet (
void )
```

11.79.3 Property Documentation

```
11.79.3.1 Sensor CMcsBus_SensorNet^ Sensor [get]
```

11.80 CPathIdentDeviceNet Class Reference

Inheritance diagram for CPathIdentDeviceNet:



Public Member Functions

- CPathIdentDeviceNet (void)
- ∼CPathIdentDeviceNet (void)
- void Set_Values (unsigned int frequency, unsigned int amplitude)
- void Measure ([System::Runtime::InteropServices::Out] unsigned int% phase, [System::Runtime::Interop←
 Services::Out] unsigned int% amplitude)


Additional Inherited Members

11.80.1 Constructor & Destructor Documentation

11.81 CPedoterDeviceNet Class Reference

unsigned int amplitude)

Inheritance diagram for CPedoterDeviceNet:



Public Member Functions

CPedoterDeviceNet ()

Initializes a new instance of the CPedoterDeviceNet class.

- virtual ∼CPedoterDeviceNet ()
- !CPedoterDeviceNet ()
- uint32_t GetCommand (uint16_t Argument)

Get value from the pedoter device

void SetCommand (uint16_t Argument, uint32_t pData)

Set value on the pedoter device

Additional Inherited Members

11.81.1 Detailed Description

11.81.2 Constructor & Destructor Documentation

```
11.81.2.1 CPedoterDeviceNet() CPedoterDeviceNet ()
```

Initializes a new instance of the CPedoterDeviceNet class.

```
11.81.2.2 ~ CPedoterDeviceNet() virtual ~ CPedoterDeviceNet () [virtual]
```

```
11.81.2.3 "!CPedoterDeviceNet() !CPedoterDeviceNet ()
```

11.81.3 Member Function Documentation

```
11.81.3.1 GetCommand() uint32_t GetCommand ( uint16_t Argument )
```

Get value from the pedoter device

Parameters

```
Argument argument
```

Returns

value


Set value on the pedoter device

Argument	argument
pData	value

11.82 CPeristalticPumpDeviceNet Class Reference

CPeristalticPumpDeviceNet is the class to control a Persistaltic Pump.

Inheritance diagram for CPeristalticPumpDeviceNet:



Public Member Functions

- CPeristalticPumpDeviceNet (void)

 Initialize a new instance of the CPeristalticPumpDeviceNet class.
- ~CPeristalticPumpDeviceNet (void)

Properties

CMcsBus_MotorControlNet[^] McsBus_MotorControl [get]

Additional Inherited Members

11.82.1 Detailed Description

CPeristalticPumpDeviceNet is the class to control a Persistaltic Pump.

11.82.2 Constructor & Destructor Documentation

11.82.2.1 CPeristalticPumpDeviceNet() CPeristalticPumpDeviceNet (void)

Initialize a new instance of the CPeristalticPumpDeviceNet class.


```
11.82.2.2 ~CPeristalticPumpDeviceNet() ~CPeristalticPumpDeviceNet (
```

11.82.3 Property Documentation

11.82.3.1 McsBus_MotorControl CMcsBus_MotorControlNet^ McsBus_MotorControl [get]

11.83 CPgaDeviceNet Class Reference

Inheritance diagram for CPgaDeviceNet:



Public Member Functions

- · CPgaDeviceNet ()
- ∼CPgaDeviceNet ()
- uint32_t GetNumFrequencyRanges ([System::Runtime::InteropServices::Out]int% numRanges)
- uint32_t GetFrequencyRange (int rangeIndex, [System::Runtime::InteropServices::Out]int% low, [System::

 Runtime::InteropServices::Out]int% high, [System::Runtime::InteropServices::Out]int% channels, [System
 ::Runtime::InteropServices::Out]int% gain)
- uint32_t GetNumAmplifications ([System::Runtime::InteropServices::Out]int% number)
- uint32_t GetAmplification (int index, [System::Runtime::InteropServices::Out]int% amplification, [System:: \leftarrow Runtime::InteropServices::Out]int% poti1, [System::Runtime::InteropServices::Out]int% poti2)
- uint32_t DefineNumFrequencyRanges (int rnum)
- uint32_t DefineFrequencyRange (int index, int low, int high, int channels, int gain)
- uint32_t DefineNumAmplifications (int number)
- uint32_t DefineAmplification (int index, int amplification, int poti1, int poti2)
- uint32_t SetGain (int channel, int Gain, int poti1, int poti2)
- uint32_t GetGain (int channel, [System::Runtime::InteropServices::Out]int% Gain, [System::Runtime::
 InteropServices::Out]int% poti1, [System::Runtime::InteropServices::Out]int% poti2)
- uint32 t ApplyGains ()

Additional Inherited Members

11.83.1 Constructor & Destructor Documentation

11.83.1.1 CPgaDeviceNet() CPgaDeviceNet ()

```
11.83.1.2 ~CPgaDeviceNet() ~CPgaDeviceNet ()
11.83.2 Member Function Documentation
11.83.2.1 ApplyGains() uint32_t ApplyGains ()
11.83.2.2 DefineAmplification() uint32_t DefineAmplification (
              int index,
              int amplification,
              int poti1,
              int poti2 )
11.83.2.3 DefineFrequencyRange() uint32_t DefineFrequencyRange (
              int index,
              int low,
              int high,
              int channels,
              int gain )
11.83.2.4 DefineNumAmplifications() uint32_t DefineNumAmplifications (
              int number )
\textbf{11.83.2.5} \quad \textbf{DefineNumFrequencyRanges()} \quad \texttt{uint32\_t DefineNumFrequencyRanges} \quad \textbf{(}
              int rnum )
11.83.2.6 GetAmplification() uint32_t GetAmplification (
              int index,
              [System::Runtime::InteropServices::Out] int% amplification,
              [System::Runtime::InteropServices::Out] int% poti1,
```


[System::Runtime::InteropServices::Out] int% poti2)

```
11.83.2.7 GetFrequencyRange() uint32_t GetFrequencyRange (
             int rangeIndex,
             [System::Runtime::InteropServices::Out] int% low,
             [System::Runtime::InteropServices::Out] int% high,
             [System::Runtime::InteropServices::Out] int% channels,
             [System::Runtime::InteropServices::Out] int% gain )
11.83.2.8 GetGain() uint32_t GetGain (
             int channel,
             [System::Runtime::InteropServices::Out] int% Gain,
             [System::Runtime::InteropServices::Out] int% poti1,
             [System::Runtime::InteropServices::Out] int% poti2 )
11.83.2.9 GetNumAmplifications() uint32_t GetNumAmplifications (
             [System::Runtime::InteropServices::Out] int% number )
11.83.2.10 GetNumFrequencyRanges() uint32_t GetNumFrequencyRanges (
             [System::Runtime::InteropServices::Out] int% numRanges )
11.83.2.11 SetGain() uint32_t SetGain (
             int channel,
             int Gain,
             int poti1,
             int poti2 )
```

11.84 CPositionIIDeviceNet Class Reference

CPositionIIDeviceNet is the class to control PositionII devices

Inheritance diagram for CPositionIIDeviceNet:



Public Member Functions

CPositionIIDeviceNet ()

Initializes a new instance of the CPositionIIDeviceNet class.

- virtual ∼CPositionIIDeviceNet ()
- !CPositionIIDeviceNet ()
- uint32_t GetCoilCommunication (uint16_t coil)

get if the communication to the coil is working

uint32 t GetOnOff (uint16 t coil)

get if the coil is switched on/off

void SwitchOnOff (uint16_t coil, uint32_t on)

switched the coild on of

uint32_t GetImplantState (uint16_t coil)

gets the implantat state

• uint32_t GetImplantCurrentSetpoint (uint16_t coil)

sets the implant current setpoint

void SetImplantCurrentSetpoint (uint16 t coil, uint32 t current)

gets the implant current setpoint

uint32 t GetPowerStrength (uint16 t coil)

sets the power for the trigger pulses

void SetPowerStrength (uint16_t coil, uint32_t power)

gets the power for the trigger pulses

uint32 t GetImplantResult (uint16 t coil)

gets the last result of the implant pulse trigger

void GetRTC ([System::Runtime::InteropServices::Out]uint8_t% year, [System::Runtime::InteropServices
::Out]uint8_t% month, [System::Runtime::InteropServices::Out]uint8_t% day, [System::Runtime::Interop
Services::Out]uint8_t% hour, [System::Runtime::InteropServices::Out]uint8_t% minute, [System::Runtime
::InteropServices::Out]uint8_t% second)

Get the RTC

• void SetRTC (uint8_t year, uint8_t month, uint8_t day, uint8_t hour, uint8_t minute, uint8_t second)

Set the RTC

uint32_t GetStateDebugData (uint16_t coil)

get the debug queue state

void SetStateDebugData (uint16_t coil, uint32_t state)

clears/starts/stops the debug queue for a certain coil

void GetDebugData (uint16_t coil, [System::Runtime::InteropServices::Out]uint16_t% index, [System::Runtime::InteropServices::Out]uint16_t% voltage, [System::Runtime::InteropServices::Out]uint16_t% numberofpulses, [System::Runtime::InteropServices::Out]uint16_t% mediantime)

get the oldest debug entry for a certain coil

uint32_t GetStateEventData ()

get the event queue state

void SetStateEventData (uint32_t state)

clears/starts/stops the event queue for a certain coil

get the oldest event entry

Properties

• CRFFunctionNet^ RFFunction [get]

Additional Inherited Members

11.84.1 Detailed Description

CPositionIIDeviceNet is the class to control PositionII devices

11.84.2 Constructor & Destructor Documentation

```
11.84.2.1 CPositionIIDeviceNet() CPositionIIDeviceNet ()
```

Initializes a new instance of the CPositionIIDeviceNet class.

```
11.84.2.2 ~CPositionIIDeviceNet() virtual ~CPositionIIDeviceNet () [virtual]
```

```
11.84.2.3 "!CPositionIIDeviceNet() !CPositionIIDeviceNet ( )
```

11.84.3 Member Function Documentation

```
11.84.3.1 GetCoilCommunication() uint32_t GetCoilCommunication ( uint16_t coil)
```

get if the communication to the coil is working

Parameters

coil the coil

Returns

is communicating

11.84.3.2 GetDebugData() void GetDebugData (

```
uint16_t coil,
[System::Runtime::InteropServices::Out] uint16_t% index,
[System::Runtime::InteropServices::Out] uint16_t% voltage,
[System::Runtime::InteropServices::Out] uint16_t% numberofpulses,
[System::Runtime::InteropServices::Out] uint16_t% mediantime)
```

get the oldest debug entry for a certain coil

Parameters

coil	the coil
index	the debug entry index number
voltage	the voltage applied
numberofpulses	the number of pulses detected
mediantime	the median time between pulses

11.84.3.3 GetEventData() void GetEventData (

```
[System::Runtime::InteropServices::Out] uint16_t% index,
[System::Runtime::InteropServices::Out] uint8_t% year,
[System::Runtime::InteropServices::Out] uint8_t% month,
[System::Runtime::InteropServices::Out] uint8_t% day,
[System::Runtime::InteropServices::Out] uint8_t% hour,
[System::Runtime::InteropServices::Out] uint8_t% minute,
[System::Runtime::InteropServices::Out] uint8_t% second,
[System::Runtime::InteropServices::Out] uint16_t% coil,
[System::Runtime::InteropServices::Out] uint16_t% type,
[System::Runtime::InteropServices::Out] uint16_t% value)
```

get the oldest event entry

Parameters

index	the event index number
year	the year
month	the month
day	the day
hour	the hour
minute	the minute
second	the second
coil	the coil
type	the event type
value	the even value

$\textbf{11.84.3.4} \quad \textbf{GetImplantCurrentSetpoint()} \quad \texttt{uint32_t} \quad \texttt{GetImplantCurrentSetpoint} \quad \textbf{(}$

uint16_t coil)

sets the implant current setpoint

Da	ro	m	Δi	ŀΔ	rc
га	17		ы	ш	15

coil t	he coil
--------	---------

the current

11.84.3.5 GetImplantResult() uint32_t GetImplantResult (uint16_t coil)

gets the last result of the implant pulse trigger

Parameters

```
coil the coil
```

Returns

the result

11.84.3.6 GetImplantState() uint32_t GetImplantState (uint16_t coil)

gets the implantat state

Parameters

```
coil the coil
```

Returns

the state

11.84.3.7 GetOnOff() uint32_t GetOnOff (uint16_t coil)

get if the coil is switched on/off

Parameters

coil the coil

0 = off, 1 = on

```
11.84.3.8 GetPowerStrength() uint32_t GetPowerStrength ( uint16_t coil )
```

sets the power for the trigger pulses

Parameters

coil	the coil
------	----------

Returns

the power

```
11.84.3.9 GetRTC() void GetRTC (
```

```
[System::Runtime::InteropServices::Out] uint8_t% year,
[System::Runtime::InteropServices::Out] uint8_t% month,
[System::Runtime::InteropServices::Out] uint8_t% day,
[System::Runtime::InteropServices::Out] uint8_t% hour,
[System::Runtime::InteropServices::Out] uint8_t% minute,
[System::Runtime::InteropServices::Out] uint8_t% second )
```

Get the RTC

Parameters

year	the year
month	the month
day	the day
hour	the hour
minute	the minute
second	the second

```
11.84.3.10 GetStateDebugData() uint32_t GetStateDebugData ( uint16_t coil )
```

get the debug queue state

Parameters

coil the coil

the state

```
11.84.3.11 GetStateEventData() uint32_t GetStateEventData ( )
```

get the event queue state

Returns

the state

gets the implant current setpoint

Parameters

coil	the coil
current	the current

gets the power for the trigger pulses

Parameters

coil	the coil
power	the power

Set the RTC

year	the year
month	the month
day	the day
hour	the hour
minute	the minute
second	the second

clears/starts/stops the debug queue for a certain coil

Parameters

coil	the coil
state	clear/start/stop

11.84.3.16 SetStateEventData() void SetStateEventData (uint32_t state)

clears/starts/stops the event queue for a certain coil

Parameters

state	clear/start/stop
-------	------------------

switched the coild on of

Parameters

coil	the coil
on	0 = off, 1 = on


11.84.4 Property Documentation

11.84.4.1 RFFunction CRFFunctionNet^ RFFunction [get]

11.85 CPositionImpDeviceNet Class Reference

CPositionImpDeviceNet is the class to access the Position/Imp devices

Inheritance diagram for CPositionImpDeviceNet:



Public Member Functions

• CPositionImpDeviceNet ()

Initializes a new instance of the CPositionImpDeviceNet class.

- virtual \sim CPositionImpDeviceNet ()
- !CPositionImpDeviceNet ()
- void ConnectImp (uint32_t id)

Connect to a Imp device with a certain ID

• uint32_t ConnectedImp ()

The ID of the connected Imp device

• int32_t GetRFFrequency ()

Gets currently used RF frequency

void SetRFFrequency (int32_t frequency)

Sets the current RF frequency

• uint32_t GetDeviceList (int32_t index)

Gets the device list

• void SetDeviceList (int32_t index, uint32_t id)

Sets the device list

• uint32_t GetImpId ()

Gets the ID of the impedance measure device

void SetImpId (uint32_t id)

Sets the ID of the impedance measure device

Additional Inherited Members

11.85.1 Detailed Description

CPositionImpDeviceNet is the class to access the Position/Imp devices

11.85.2 Constructor & Destructor Documentation

```
11.85.2.1 CPositionImpDeviceNet() CPositionImpDeviceNet ()
```

Initializes a new instance of the CPositionImpDeviceNet class.

```
11.85.2.2 ~CPositionImpDeviceNet() virtual ~CPositionImpDeviceNet () [virtual]
```

```
11.85.2.3 "!CPositionImpDeviceNet() !CPositionImpDeviceNet ()
```

11.85.3 Member Function Documentation

```
11.85.3.1 ConnectedImp() uint32_t ConnectedImp ()
```

The ID of the connected Imp device

Returns

The ID

```
11.85.3.2 ConnectImp() void ConnectImp ( uint32_t id )
```

Connect to a Imp device with a certain ID

Parameters

id The ID

Gets the device list

index 1	the index
---------	-----------

Returns

the ID

11.85.3.4 GetImpld() uint32_t GetImpld ()

Gets the ID of the impedance measure device

Returns

the ID

11.85.3.5 GetRFFrequency() int32_t GetRFFrequency ()

Gets currently used RF frequency

Returns

The frequency

Sets the device list

Parameters

index	the index
id	the ID

```
11.85.3.7 SetImpId() void SetImpId ( uint32_t id )
```

Sets the ID of the impedance measure device

id the ID


Sets the current RF frequency

Parameters

frequency The frequency

11.86 CPPCDeviceNet Class Reference

Inheritance diagram for CPPCDeviceNet:



Public Member Functions

• CPPCDeviceNet (void)

Properties

- CPPCFunctionNet^ PPCFunction [get]
- CMcsBusNet^ McsBus [get]
- CMcsBus_MotorControlNet^ McsBus_MotorControl [get]
- CMcsBus_SensorNet^ McsBus_Sensor [get]

Additional Inherited Members

11.86.1 Constructor & Destructor Documentation

11.86.1.1 CPPCDeviceNet() CPPCDeviceNet (void)

11.86.2 Property Documentation

```
11.86.2.1 McsBus CMcsBusNet^ McsBus [get]
```

11.86.2.2 McsBus_MotorControl CMcsBus_MotorControlNet^ McsBus_MotorControl [get]


11.86.2.3 McsBus_Sensor CMcsBus_SensorNet^ McsBus_Sensor [get]

11.86.2.4 PPCFunction CPPCFunctionNet^ PPCFunction [get]

11.87 CPPCFunctionNet Class Reference

CPPCFunctionNet is the class to access the PPC (high precision Patch Peristalic patch Pump

Inheritance diagram for CPPCFunctionNet:



Public Member Functions

CPPCFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pPPCFunctionPointer←
 Container)

Initializes a new instance of the CPPCFunctionNet class.

- CPPCFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ∼CPPCFunctionNet ()
- !CPPCFunctionNet ()
- int GetPumpSpeedUnit (uint16_t channel)

Reads the Pump Speed Unit

void SetPumpSpeedUnit (uint16_t channel, int SpeedUnit)

Writes the Pump Speed Unit

PP_Pump_Mode_Type_EnumNet GetPumpModeType (uint16_t channel)

Reads the Pump Mode Type.

• void SetPumpModeType (uint16_t channel, PP_Pump_Mode_Type_EnumNet PumpMode)

Writes the config string from the device.

void GetAnalogVoltageRange (uint16_t channel, [System::Runtime::InteropServices::Out]uint16_t% min_
 voltage, [System::Runtime::InteropServices::Out]uint16_t% max_voltage)

Reads the Analog Input Voltage Range

void SetAnalogVoltageRange (uint16_t channel, uint16_t min_voltage, uint16_t max_voltage)

Writes the Analog Input Voltage Range

void GetPressureRange (uint16_t channel, [System::Runtime::InteropServices::Out]int32_t% lower_
 pressure, [System::Runtime::InteropServices::Out]int32_t% upper_pressure)

Get the pressure range that is used between the analog voltage or the digital states

void SetPressureRange (uint16_t channel, int32_t lower_pressure, int32_t upper_pressure)

Get the pressure range that is used between the analog voltage or the digital states

uint16_t GetSupplyVoltage ()

Reads the current supply voltage in mV

uint16 t GetAnalogVoltage (uint16 t channel)

Reads the current analog voltage

uint16 t GetDigitalIn (uint16 t channel)

Reads the digital input state

• int GetValveActive (uint16 t valve)

Gets the valve active/inactive state

void SetValveActive (uint16_t valve, int valveActive)

Sets the valve active/inactive state

• void SetPressureOffset ()

Sets the pressure offset

void LoadPressure (int32_t pressure, uint32_t options)

Loads the reservoir with a pressure

void IsBusy ([System::Runtime::InteropServices::Out]int16_t% task, [System::Runtime::InteropServices::
 Out]int16_t% wait)

Is the PPC busy with a task

void FirePressurePulse (int32_t duration, int32_t nextpressure)

Fire a pressure pulse from the reservoir

int32_t MeasureReservoir ()

Measures the reservoir pressure

Additional Inherited Members

11.87.1 Detailed Description

CPPCFunctionNet is the class to access the PPC (high precision Patch Peristalic patch Pump

11.87.2 Constructor & Destructor Documentation

Initializes a new instance of the CPPCFunctionNet class.

```
11.87.2.2 CPPCFunctionNet() [2/2] CPPCFunctionNet (
CMcsUsbNet^ mcsusb )
```

11.87.2.3 ~CPPCFunctionNet() virtual ~CPPCFunctionNet () [virtual]

```
11.87.2.4 "!CPPCFunctionNet() !CPPCFunctionNet ( )
```

11.87.3 Member Function Documentation

Fire a pressure pulse from the reservoir

Parameters

duration	The pulse duration (valves open)
nextpressure	The next pressure

Reads the current analog voltage

Parameters

channel	The Channel Number
---------	--------------------

Returns

The Analog Voltage

11.87.3.3 **GetAnalogVoltageRange()** void GetAnalogVoltageRange (uint16_t *channel*,

[System::Runtime::InteropServices::Out] uint16_t% min_voltage, [System::Runtime::InteropServices::Out] uint16_t% max_voltage)

Reads the Analog Input Voltage Range

Parameters

channel	The Channel Number
min_voltage	The voltage that should be seen as the minimum voltage
max_voltage	The voltage that should be seen as the maximum voltage

Reads the digital input state

Parameters

channel The Channel Number

Returns

The Digital State

Get the pressure range that is used between the analog voltage or the digital states

Parameters

channel	The Channel Number
lower_pressure	The lower border of the pressure range
upper_pressure	The upper border of the pressure range

```
11.87.3.6 GetPumpModeType() PP_Pump_Mode_Type_EnumNet GetPumpModeType ( uint16_t channel )
```

Reads the Pump Mode Type.

Parameters

channel	The Channel Number

Returns

The Pump Mode Type.

Reads the Pump Speed Unit

Parameters

Returns

The Speed Unit

11.87.3.8 GetSupplyVoltage() uint16_t GetSupplyVoltage ()

Reads the current supply voltage in mV

Returns

The supply voltage

```
11.87.3.9 GetValveActive() int GetValveActive ( uint16_t valve )
```

Gets the valve active/inactive state

Parameters

```
valve The valve number
```

Returns

The valve state

Is the PPC busy with a task

Parameters

task	The task state
wait	The wait state

Loads the reservoir with a pressure

Parameters

pressure	The pressure
options	The options: end with 0=regulate on patch 1=regulate on reservoir

11.87.3.12 MeasureReservoir() int32_t MeasureReservoir ()

Measures the reservoir pressure

Returns

The pressure

Writes the Analog Input Voltage Range

Parameters

channel	The Channel Number
min_voltage	The voltage that should be seen as the minimum voltage
max_voltage	The voltage that should be seen as the maximum voltage

11.87.3.14 SetPressureOffset() void SetPressureOffset ()

Sets the pressure offset

Get the pressure range that is used between the analog voltage or the digital states

Parameters

channel	The Channel Number
lower_pressure	The lower border of the pressure range
upper_pressure	The upper border of the pressure range

Writes the config string from the device.

Parameters

channel	The Channel Number
PumpMode	The Pump Mode Type.

Writes the Pump Speed Unit

Parameters

channel	The Channel Number
SpeedUnit	The Speed Unit


Sets the valve active/inactive state

Parameters

valve	The valve number
valveActive	The valve state

11.88 CPPS_DeviceNet Class Reference

Inheritance diagram for CPPS_DeviceNet:



Public Member Functions

CPPS_DeviceNet (void)

Properties

- CPPS_FunctionNet^ PPS_Function [get]
- CMcsBusNet^ McsBus [get]
- CMcsBus_MotorControlNet^ McsBus_MotorControl [get]
- CMcsBus_SensorNet^ McsBus_Sensor [get]

Additional Inherited Members

11.88.1 Constructor & Destructor Documentation

```
11.88.1.1 CPPS_DeviceNet() CPPS_DeviceNet (
void )
```

11.88.2 Property Documentation

11.88.2.1 McsBus CMcsBusNet^ McsBus [get]


11.88.2.2 McsBus_MotorControl CMcsBus_MotorControlNet^ McsBus_MotorControl [get]

11.88.2.3 McsBus_Sensor CMcsBus_SensorNet^ McsBus_Sensor [get]

11.88.2.4 PPS_Function CPPS_FunctionNet^ PPS_Function [get]

11.89 CPPS_FunctionNet Class Reference

Inheritance diagram for CPPS_FunctionNet:



Public Member Functions

- CPPS_FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] cPPS_FunctionPointer←
 Container)
- CPPS FunctionNet (CMcsUsbNet[^] mcsusb)
- · void SetPumpMaxSpeed (unsigned short index, unsigned short maxspeed)
- unsigned short GetPumpMaxSpeed (unsigned short index)
- void SetPumpSpeedUnit (unsigned short index, int speedunit)
- int GetPumpSpeedUnit (unsigned short index)
- void SetPumpModeType (unsigned short index, PP_Pump_Mode_Type_EnumNet type)
- PP_Pump_Mode_Type_EnumNet GetPumpModeType (unsigned short index)
- void SetPumpCouple (unsigned int i)
- unsigned int GetPumpCouple ()
- void SetPumpEnableSpeedRatio (unsigned int enable)
- unsigned int GetPumpEnableSpeedRatio ()
- · void SetPumpManualOnOff (unsigned short index, unsigned int onoff)
- unsigned int GetPumpManualOnOff (unsigned short index)
- void SetPumpFunctionSpeeds (unsigned short index, short offspeed, short onspeed)
- void GetPumpFunctionSpeeds (unsigned short index, [System::Runtime::InteropServices::Out]short% offspeed, [System::Runtime::InteropServices::Out]short% onspeed)
- void SetPumpSpeedRatio (int ratio)
- int GetPumpSpeedRatio ()
- void SetPumpFastOnOff (unsigned short index, unsigned int onoff)
- unsigned int GetPumpFastOnOff (unsigned short index)
- void SetPumpFastSpeed (unsigned short index, short fastspeed)
- short GetPumpFastSpeed (unsigned short index)
- void SetAnalogVoltages (unsigned short index, unsigned short minvoltage, unsigned short maxvoltage)
- void GetAnalogVoltages (unsigned short index, [System::Runtime::InteropServices::Out]unsigned short% minvoltage, [System::Runtime::InteropServices::Out]unsigned short% maxvoltage)
- void SetUseBubble (unsigned short index, unsigned int usebubble)
- unsigned int GetUseBubble (unsigned short index)
- unsigned short GetSupplyVoltage ()
- unsigned short GetAnalogVoltage (unsigned short index)
- unsigned short GetDigitalIn (unsigned short index)
- · unsigned short GetBubbleState ()

Additional Inherited Members

11.89.1 Constructor & Destructor Documentation

```
11.89.1.1 CPPS_FunctionNet() [1/2] CPPS_FunctionNet (
             CMcsUsbNet^ mcsusb,
             {\tt CMcsUsbFunctionPointerContainer}^{\wedge} \ \ cPPS\_FunctionPointerContainer \ )
11.89.1.2 CPPS_FunctionNet() [2/2] CPPS_FunctionNet (
             CMcsUsbNet^ mcsusb )
11.89.2 Member Function Documentation
11.89.2.1 GetAnalogVoltage() unsigned short GetAnalogVoltage (
             unsigned short index)
11.89.2.2 GetAnalogVoltages() void GetAnalogVoltages (
             unsigned short index,
             [System::Runtime::InteropServices::Out] unsigned short% minvoltage,
             [System::Runtime::InteropServices::Out] unsigned short% maxvoltage )
11.89.2.3 GetBubbleState() unsigned short GetBubbleState ( )
11.89.2.4 GetDigitalIn() unsigned short GetDigitalIn (
             unsigned short index)
11.89.2.5 GetPumpCouple() unsigned int GetPumpCouple ()
11.89.2.6 GetPumpEnableSpeedRatio() unsigned int GetPumpEnableSpeedRatio ( )
11.89.2.7 GetPumpFastOnOff() unsigned int GetPumpFastOnOff (
             unsigned short index)
```


```
11.89.2.8 GetPumpFastSpeed() short GetPumpFastSpeed (
              unsigned short index)
11.89.2.9 GetPumpFunctionSpeeds() void GetPumpFunctionSpeeds (
              unsigned short index,
              [System::Runtime::InteropServices::Out] short% offspeed,
              [System::Runtime::InteropServices::Out] short% onspeed )
\textbf{11.89.2.10} \quad \textbf{GetPumpManualOnOff()} \quad \texttt{unsigned int GetPumpManualOnOff} \quad \textbf{(}
              unsigned short index)
11.89.2.11 GetPumpMaxSpeed() unsigned short GetPumpMaxSpeed (
              unsigned short index)
11.89.2.12 GetPumpModeType() PP_Pump_Mode_Type_EnumNet GetPumpModeType (
              unsigned short index)
11.89.2.13 GetPumpSpeedRatio() int GetPumpSpeedRatio ( )
11.89.2.14 GetPumpSpeedUnit() int GetPumpSpeedUnit (
              unsigned short index)
11.89.2.15 GetSupplyVoltage() unsigned short GetSupplyVoltage ( )
{\bf 11.89.2.16} \quad {\bf GetUseBubble()} \quad {\tt unsigned int GetUseBubble ()}
              unsigned short index)
```

```
\textbf{11.89.2.17} \quad \textbf{SetAnalogVoltages()} \quad \texttt{void SetAnalogVoltages ()}
             unsigned short index,
             unsigned short minvoltage,
             unsigned short maxvoltage )
11.89.2.18 SetPumpCouple() void SetPumpCouple (
             unsigned int i )
11.89.2.19 SetPumpEnableSpeedRatio() void SetPumpEnableSpeedRatio (
             unsigned int enable )
11.89.2.20 SetPumpFastOnOff() void SetPumpFastOnOff (
             unsigned short index,
             unsigned int onoff )
11.89.2.21 SetPumpFastSpeed() void SetPumpFastSpeed (
             unsigned short index,
             short fastspeed )
11.89.2.22 SetPumpFunctionSpeeds() void SetPumpFunctionSpeeds (
             unsigned short index,
             short offspeed,
             short onspeed )
11.89.2.23 SetPumpManualOnOff() void SetPumpManualOnOff (
             unsigned short index,
             unsigned int onoff )
11.89.2.24 SetPumpMaxSpeed() void SetPumpMaxSpeed (
             unsigned short index,
             unsigned short {\it maxspeed} )
```

11.90 CPPSDeviceNet Class Reference

CPPS4plus1DeviceNet is the to control the MCS HLA device

Inheritance diagram for CPPSDeviceNet:



Public Member Functions

• CPPSDeviceNet (void)

Additional Inherited Members

11.90.1 Detailed Description

CPPS4plus1DeviceNet is the to control the MCS HLA device


11.90.2 Constructor & Destructor Documentation

```
11.90.2.1 CPPSDeviceNet() CPPSDeviceNet (
void )
```

11.91 CProgramPressureCurveNet Class Reference

CProgramPressureCurveNet is the class to program pressure curves

Inheritance diagram for CProgramPressureCurveNet:



Public Member Functions

- CProgramPressureCurveNet (CMcsUsbNet[^] mcsusb)
 Initializes a new instance of the CPulseGeneratorFunctionNet class.
- virtual ~CProgramPressureCurveNet (void)
- !CProgramPressureCurveNet (void)
- void Program (unsigned char busnumber, unsigned char busaddress, int32_t channel, array< int32_t $>^{\land}$ pressures, array< int32_t $>^{\land}$ durations)
- void SetRepeats (unsigned char busnumber, unsigned char busaddress, int32_t channel, uint32_t repeats)
- unsigned int GetRepeats (unsigned char busnumber, unsigned char busaddress, int32 t channel)

Additional Inherited Members

11.91.1 Detailed Description

CProgramPressureCurveNet is the class to program pressure curves

11.91.2 Constructor & Destructor Documentation

```
11.91.2.1 CProgramPressureCurveNet() CProgramPressureCurveNet (
CMcsUsbNet^ mcsusb )
```

Initializes a new instance of the CPulseGeneratorFunctionNet class.


```
11.91.2.2 ~CProgramPressureCurveNet() virtual ~CProgramPressureCurveNet (
             void ) [virtual]
11.91.2.3 "!CProgramPressureCurveNet() !CProgramPressureCurveNet (
             void )
11.91.3 Member Function Documentation
11.91.3.1 GetRepeats() unsigned int GetRepeats (
             unsigned char busnumber,
             unsigned char busaddress,
             int32_t channel )
11.91.3.2 Program() void Program (
             unsigned char busnumber,
             unsigned char busaddress,
             int32_t channel,
             array< int32_t >^{\wedge} pressures,
             array< int32_t >^{\wedge} steps,
             array< int16_t >^{\land} durations )
11.91.3.3 SetRepeats() void SetRepeats (
             unsigned char busnumber,
             unsigned char busaddress,
             int32_t channel,
```

11.92 CPulseGeneratorFunctionNet Class Reference

CPulseGeneratorFunctionNet is the class to control the pulse generator for video tracking

Inheritance diagram for CPulseGeneratorFunctionNet:

uint32_t repeats)



Public Member Functions

CPulseGeneratorFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pPulse←
 GeneratorFunctionPointerContainer)

Initializes a new instance of the CPulseGeneratorFunctionNet class.

- CPulseGeneratorFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CPulseGeneratorFunctionNet ()
- !CPulseGeneratorFunctionNet ()
- int32_t GetPeriod (int32_t generator_number)

Reads the generator period

• void SetPeriod (int32_t generator_number, int32_t period_in_samples)

Writes the generator period

int32_t GetPulseLength (int32_t generator_number)

Reads the generator pulse length

void SetPulseLength (int32 t generator number, int32 t pulselength in 10us)

Writes the generator pulse length

void GetModeSelect (int32_t generator_number, [System::Runtime::InteropServices::Out]PulseGenerator
 —Mode_EnumNet% mode, [System::Runtime::InteropServices::Out]int32_t% digitalchannel)

Reads the generator mode

void SetModeSelect (int32_t generator_number, PulseGenerator_Mode_EnumNet mode, int32_t digitalchannel)

Writes the generator mode

Additional Inherited Members

11.92.1 Detailed Description

CPulseGeneratorFunctionNet is the class to control the pulse generator for video tracking

11.92.2 Constructor & Destructor Documentation

```
11.92.2.1 CPulseGeneratorFunctionNet() [1/2] CPulseGeneratorFunctionNet (

CMcsUsbNet^ mcsusb,

CMcsUsbFunctionPointerContainer^ pPulseGeneratorFunctionPointerContainer)
```

Initializes a new instance of the CPulseGeneratorFunctionNet class.

```
11.92.2.2 CPulseGeneratorFunctionNet() [2/2] CPulseGeneratorFunctionNet (
CMcsUsbNet^ mcsusb )
```

```
11.92.2.3 ~CPulseGeneratorFunctionNet() virtual ~CPulseGeneratorFunctionNet ( ) [virtual]
```

11.92.2.4 "!CPulseGeneratorFunctionNet() !CPulseGeneratorFunctionNet ()

11.92.3 Member Function Documentation

Reads the generator mode

Parameters

generator_number	The generator number
mode	The generator mode
digitalchannel	The digital in channel used as gate

Reads the generator period

Parameters

generator_number	The generator number
------------------	----------------------

Returns

The period

Reads the generator pulse length

Parameters

generator_number The generator number

Returns

The pulse length

Writes the generator mode

Parameters

generator_number	The generator number
mode	The generator mode
digitalchannel	The digital in channel used as gate

Writes the generator period

Parameters

generator_number	The generator number
period_in_samples	The period


Writes the generator pulse length

Parameters

generator_number	The generator number
pulselength_in_10us	The pulse length

11.93 CRadioControledDevicesNet Class Reference

Inheritance diagram for CRadioControledDevicesNet:



Public Member Functions

- CRadioControledDevicesNet (void)
- bool HasRadioControl ()
- array< unsigned short > ^ GetDeviceNames ()
- void ConnectDevice (unsigned short sn)
- void DisConnectDevice ()
- bool StillConnected ()
- void SetFrequency (unsigned short frequency)
- unsigned short GetFrequency ()

Protected Member Functions

• CRadioControledDevicesNet (CRadioControledDevices *pRadioControled)

Additional Inherited Members

11.93.1 Constructor & Destructor Documentation

```
11.93.1.1 CRadioControledDevicesNet() [1/2] CRadioControledDevicesNet ( void )
```

```
11.93.1.2 CRadioControledDevicesNet() [2/2] CRadioControledDevicesNet (
CRadioControledDevices * pRadioControled) [protected]
```

11.93.2 Member Function Documentation

```
11.93.2.1 ConnectDevice() void ConnectDevice ( unsigned short sn )
```

11.93.2.2 DisConnectDevice() void DisConnectDevice ()

11.94 CCMOSMeaDeviceNet::CRegionOfInterestRect Class Reference

Public Member Functions

- CRegionOfInterestRect (int left, int top, int right, int bottom)
- CRegionOfInterestRect ^ DeepCopy ()

Public Attributes

```
• int m_Left
```

- int m_Top
- int m_Right
- int m_Bottom

11.94.1 Constructor & Destructor Documentation

```
11.94.1.1 CRegionOfInterestRect() CRegionOfInterestRect (
          int left,
          int top,
          int right,
          int bottom )
```

11.94.2 Member Function Documentation

11.94.2.1 DeepCopy() CRegionOfInterestRect ^ DeepCopy ()

11.94.3 Member Data Documentation

11.94.3.1 m_Bottom int m_Bottom


11.94.3.2 m_Left int m_Left

11.94.3.3 m_Right int m_Right

11.94.3.4 m_Top int m_Top

11.95 CRetinaLedDeviceNet Class Reference

Inheritance diagram for CRetinaLedDeviceNet:



Public Member Functions

- CRetinaLedDeviceNet ()
- ∼CRetinaLedDeviceNet ()
- unsigned int SetTrigger (int enable)
- unsigned int SetLED (unsigned long long pattern)
- unsigned int SetTablepointer (int position)
- unsigned int GetTablepointer (int % position)
- unsigned int ClearTable ()
- unsigned int AddTableEntry (unsigned long long pattern)
- unsigned int AddLoopEntry (unsigned short repeats, unsigned short steps_back)
- unsigned int SetRepeat (int repeat)
- unsigned int SetLumi (int lumi)
- unsigned int SetPersistency (unsigned int persistency)

Additional Inherited Members


```
11.95.1 Constructor & Destructor Documentation
```

```
11.95.1.1 CRetinaLedDeviceNet() CRetinaLedDeviceNet ()
11.95.1.2 \sim CRetinaLedDeviceNet() \sim CRetinaLedDeviceNet ()
11.95.2 Member Function Documentation
11.95.2.1 AddLoopEntry() unsigned int AddLoopEntry (
             unsigned short repeats,
             unsigned short steps_back )
11.95.2.2 AddTableEntry() unsigned int AddTableEntry (
             unsigned long long pattern )
11.95.2.3 ClearTable() unsigned int ClearTable ( )
11.95.2.4 GetTablepointer() unsigned int GetTablepointer (
             int % position )
11.95.2.5 SetLED() unsigned int SetLED (
             unsigned long long pattern )
11.95.2.6 SetLumi() unsigned int SetLumi (
             int lumi )
```

11.96 CRFFunctionNet Class Reference

CRFFunctionNet is the class to control RF devices

Inheritance diagram for CRFFunctionNet:



Public Member Functions

CRFFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pRFFunctionPointer←
 Container)

Initializes a new instance of the CRFFunctionNet class.

- CRFFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CRFFunctionNet ()
- !CRFFunctionNet ()
- uint32_t GetBaseFrequency (CFirmwareDestinationNet destination)

gets the base advertise frequency

void SetBaseFrequency (CFirmwareDestinationNet destination, uint32_t frequency)

sets the base advertise frequency

• uint32 t GetWorkingFrequency ()

gets the working frequency

void SetWorkingFrequency (uint32_t frequency)

sets the working frequency

array< uint32_t > ^ GetAvailableDeviceListEx (int list_Length)

get a list of available devices

array< uint32_t > ^ GetAvailableDeviceList ()

```
get a list of available devices

    array< uint32_t > ^ GetAvailableStateListEx (int list_Length)

          get a list of the states of the available devices

    array< uint32_t > ^ GetAvailableStateList ()

          get a list of the states of the available devices
    • void Connect (uint32_t sn)
          connect to a RF device, use 0 to disconnect

    uint32_t GetConnectedDevice ()

          get connect RF device, 0 = no device connected
    • uint32_t GetState ()
          get connection state

    void SetTestMode (uint32_t mode)

          set test mode

    uint32_t GetTestMode ()

          gets test mode
Additional Inherited Members
11.96.1 Detailed Description
CRFFunctionNet is the class to control RF devices
11.96.2 Constructor & Destructor Documentation
11.96.2.1 CRFFunctionNet() [1/2] CRFFunctionNet (
               CMcsUsbNet<sup>∧</sup> mcsusb,
               CMcsUsbFunctionPointerContainer^{\wedge} pRFFunctionPointerContainer)
Initializes a new instance of the CRFFunctionNet class.
```

```
11.96.2.2 CRFFunctionNet() [2/2] CRFFunctionNet (
CMcsUsbNet^ mcsusb )

11.96.2.3 ~CRFFunctionNet() virtual ~CRFFunctionNet ( ) [virtual]

11.96.2.4 "!CRFFunctionNet() !CRFFunctionNet ( )

11.96.3 Member Function Documentation
```

```
Generated by Doxygen
```

11.96.3.1 Connect() void Connect (
uint32_t sn)

connect to a RF device, use 0 to disconnect

Parameters

sn the serial number

```
11.96.3.2 GetAvailableDeviceList() array<uint32_t> ^ GetAvailableDeviceList ( )
```

get a list of available devices

Returns

array of devices

```
11.96.3.3 GetAvailableDeviceListEx() array<uint32_t> ^{\land} GetAvailableDeviceListEx ( int list\_Length)
```

get a list of available devices

Parameters

list_Length	The maximal length of list.
-------------	-----------------------------

Returns

array of devices

```
11.96.3.4 GetAvailableStateList() array<uint32_t> ^ GetAvailableStateList ( )
```

get a list of the states of the available devices

Returns

array of states

```
11.96.3.5 GetAvailableStateListEx() array<uint32_t> ^{\land} GetAvailableStateListEx ( int list_Length)
```

get a list of the states of the available devices

Parameters

list_Length	The maximal length of list.
-------------	-----------------------------

Returns

array of states

11.96.3.6 GetBaseFrequency() uint32_t GetBaseFrequency (CFirmwareDestinationNet destination)

gets the base advertise frequency

Parameters

destination	the destination to query
-------------	--------------------------

Returns

the frequency

11.96.3.7 GetConnectedDevice() uint32_t GetConnectedDevice ()

get connect RF device, 0 = no device connected

Returns

the serial number

11.96.3.8 GetState() uint32_t GetState ()

get connection state

Returns

the state

```
11.96.3.9 GetTestMode() uint32_t GetTestMode ( )
```

gets test mode

Returns

the mode

11.96.3.10 GetWorkingFrequency() uint32_t GetWorkingFrequency ()

gets the working frequency

Returns

the frequency

sets the base advertise frequency

Parameters

destination	the destination to set
frequency	the frequency

```
11.96.3.12 SetTestMode() void SetTestMode ( uint32_t mode )
```

set test mode

Parameters

mode the mode

```
11.96.3.13 SetWorkingFrequency() void SetWorkingFrequency ( uint32_t frequency )
```


sets the working frequency

Parameters

frequency the frequency

11.97 CRobo_FYIProgram_FunctionNet Class Reference

Inheritance diagram for CRobo_FYIProgram_FunctionNet:



Public Member Functions

- CRobo_FYIProgram_FunctionNet (CMcsUsbNet[∧] mcsusb, CMcsUsbFunctionPointerContainer[∧] robo_FY → IProgram_FunctionPointerContainer)
- CRobo FYIProgram FunctionNet (CMcsUsbNet[∧] mcsusb)
- void SetValve1 (unsigned char index, unsigned int valve1)
- unsigned int GetValve1 (unsigned char index)
- void SetValve2 (unsigned char index, unsigned int valve2)
- unsigned int GetValve2 (unsigned char index)
- void SetLength (unsigned char index, int length)
- int GetLength (unsigned char index)
- void Start ()
- · int GetState ()

Additional Inherited Members

11.97.1 Constructor & Destructor Documentation


```
11.97.1.2 CRobo_FYIProgram_FunctionNet() [2/2] CRobo_FYIProgram_FunctionNet (
CMcsUsbNet^ mcsusb )
```

11.97.2 Member Function Documentation

```
\textbf{11.97.2.1} \quad \textbf{GetLength()} \quad \texttt{int GetLength ()}
              unsigned char index )
11.97.2.2 GetState() int GetState ()
11.97.2.3 GetValve1() unsigned int GetValve1 (
              unsigned char index)
11.97.2.4 GetValve2() unsigned int GetValve2 (
              unsigned char index)
11.97.2.5 SetLength() void SetLength (
              unsigned char index,
              int length )
11.97.2.6 SetValve1() void SetValve1 (
              unsigned char index,
              unsigned int valve1 )
11.97.2.7 SetValve2() void SetValve2 (
              unsigned char index,
              unsigned int valve2)
11.97.2.8 Start() void Start ()
```

11.98 CRobo_FYITemp_FunctionNet Class Reference

Inheritance diagram for CRobo_FYITemp_FunctionNet:



Public Member Functions

- CRobo_FYITemp_FunctionNet (CMcsUsbNet[^] mcsusb)
- void SetRegulatorOnOff (unsigned char index, int onoff)
- int GetRegulatorOnOff (unsigned char index)
- void SetSollTemp (unsigned char index, int temp)
- int GetSollTemp (unsigned char index)
- void SetPCoeff (unsigned char index, int pcoeff)
- int GetPCoeff (unsigned char index)
- void SetICoeff (unsigned char index, int icoeff)
- int GetlCoeff (unsigned char index)
- void SetMaxPower (unsigned char index, int power)
- int GetMaxPower (unsigned char index)

Additional Inherited Members

11.98.1 Constructor & Destructor Documentation

```
11.98.1.1 CRobo_FYITemp_FunctionNet() CRobo_FYITemp_FunctionNet (
CMcsUsbNet^ mcsusb )
```

11.98.2 Member Function Documentation

```
11.98.2.1 GetlCoeff() int GetICoeff (
          unsigned char index )
```

```
11.98.2.2 GetMaxPower() int GetMaxPower ( unsigned char index )
```


```
11.98.2.3 GetPCoeff() int GetPCoeff (
unsigned char index)
```

```
11.98.2.4 GetRegulatorOnOff() int GetRegulatorOnOff (
    unsigned char index )
```

```
11.98.2.5 GetSollTemp() int GetSollTemp (
             unsigned char index )
11.98.2.6 SetlCoeff() void SetICoeff (
             unsigned char index,
             int icoeff )
11.98.2.7 SetMaxPower() void SetMaxPower (
             unsigned char index,
             int power )
11.98.2.8 SetPCoeff() void SetPCoeff (
             unsigned char index,
             int pcoeff )
11.98.2.9 SetRegulatorOnOff() void SetRegulatorOnOff (
             unsigned char index,
             int onoff )
11.98.2.10 SetSollTemp() void SetSollTemp (
             unsigned char index,
             int temp )
```

11.99 CRoboDacqNet Class Reference

Inheritance diagram for CRoboDacqNet:



Public Member Functions

- CRoboDacqNet (void)
- CRoboDacqNet (CRoboDeviceNet[^] robodevice)
- void RunTable ()
- void RunTable (int timeout)
- void StopTable ()
- void StopTable (int timeout)
- void CancelTableLoop ()
- void CancelTableLoopAndStopTable ()
- void SetConfigurationBit (unsigned short bit, bool value)
- void SetConfigurationBitSupply (bool value)
- void SetConfigurationBitRelais (bool value)
- void SetConfigurationBitStream (bool value)
- void SetConfigurationBitAxc (bool value)
- void SetConfigurationBitCC_Gen (bool value)
- void SetConfigurationBitCV Gen (bool value)
- void SetConfigurationBitRC_Gen (bool value)
- void SetConfigurationBitRV Gen (bool value)
- void SetConfigurationBitBlu_Led (bool value)
- void GotGoringarationEttEld_Edd (BGG) value)
- void SetConfigurationBitRed_Led (bool value)
- void SetConfigurationBitBlu_LedToggleSlow (bool value)
- void SetConfigurationBitRed LedToggleSlow (bool value)
- void SetConfigurationBitBlu_LedToggleFast (bool value)
- void SetConfigurationBitRed_LedToggleFast (bool value)
- void SetConfigurationBitRed_LedSaturation (bool value)
- void SetSimulation (unsigned int enable)
- void SetUClamp (int uClamp)
- void SetIClamp (int iClamp)
- void SetPGain (int pGain)
- void SetlGain (int iGain)
- void SetFilter (int filter)
- void SetUVOffset (int UVOffset)
- void SetUCOffset (int UCOffset)
- void SetICOffset (int ICOffset)
- void SetCrossTalkOffset (int CrossTalk)
- void SetXGain (int xGain)
- void SetCrossTalkOptimum (int cxOptimum)
- void SetRecordingNumber (unsigned int recordingNumber)
- void ClampAmpRestart ()
- void DoRamp (int startValue, int endValue, int duration, int mode)
- unsigned int GetClampAmpSerialNumber ()
- unsigned int GetConfigurationBits ()
- bool GetConfigurationBit (unsigned short bit)
- bool GetConfigurationBitSupply ()
- bool GetConfigurationBitRelais ()
- · bool GetConfigurationBitStream ()
- bool GetConfigurationBitAxc ()
- bool GetConfigurationBitCC Gen ()
- bool GetConfigurationBitCV_Gen ()
- bool GetConfigurationBitRC Gen ()
- bool GetConfigurationBitRV Gen ()
- bool GetConfigurationBitBlu Led ()
- bool GetConfigurationBitRed_Led ()
- bool GetConfigurationBitBlu_LedToggleSlow ()

- bool GetConfigurationBitRed_LedToggleSlow ()
- bool GetConfigurationBitBlu_LedToggleFast ()
- bool GetConfigurationBitRed LedToggleFast ()
- bool GetConfigurationBitRed_LedSaturation ()
- unsigned int GetSimulation ()
- int GetUClamp ()
- int GetlClamp ()
- int GetPGain ()
- int GetlGain ()
- int GetFilter ()
- int GetUVOffset ()
- · int GetUCOffset ()
- int GetICOffset ()
- int GetCrossTalkOffset ()
- int GetXGain ()
- int GetCrossTalkOptimum ()
- unsigned int GetRecordingNumber ()
- int GetResistanceC ()
- int GetResistanceV ()
- int GetCapacityC ()
- · int GetCapacityV ()
- int GetCapacityX ()
- int GetUV ()
- int GetUC ()
- int GetIC ()
- int GetNUV MS ()
- int GetNUC MS ()
- int GetNIC MS ()
- void SetAllDigout (uint32_t value)
- uint32_t GetAllDigout ()
- · void SetCommand (unsigned char command, int value)
- int GetCommand (unsigned char command)
- void SetDigout (uint16_t index, bool enable)
- bool GetDigout (uint16 t index)
- · void TableDefBegin ()
- void TableDefEnd ()
- void Table_Wait (unsigned int tableWait)
- void SetDownsampleFactor (int index, int downsample_factor)
- void SetFilterCoeffs (int index, array< int >^ coeffs)
- · void SetNoFilterCoeffs (int index)
- int GetDownsampleFactor (int index)
- array< int > ^ GetFilterCoeffs (int index)
- void Emu_SetElectrodeResists (int emuElectrodeResist)
- void Emu SetCellResists (int emuCellResist)
- void Emu_SetCellCapacity (int emuCellCapacity)
- void Emu_SetCellPotential (int emuCellPotential)
- void Emu_SetNoise (int emuNoise)
- int Emu_GetElectrodeResists ()
- int Emu_GetCellResists ()
- int Emu_GetCellCapacity ()int Emu_GetCellPotential ()
- int Emu_GetNoise ()
- void SetDisplayText (int index, String[^] displayText)
- void SetScreen (int screen)
- void UpdateDisplay ()
- String \(^\) GetDisplayText (int index)
- int GetScreen ()
- int GetUpdateDisplay ()

Additional Inherited Members

```
11.99.1 Constructor & Destructor Documentation
```

```
11.99.1.1 CRoboDacqNet() [1/2] CRoboDacqNet (
            void )
11.99.1.2 CRoboDacqNet() [2/2] CRoboDacqNet (
             CRoboDeviceNet^ robodevice )
11.99.2 Member Function Documentation
11.99.2.1 CancelTableLoop() void CancelTableLoop ( )
11.99.2.2 CancelTableLoopAndStopTable() void CancelTableLoopAndStopTable ( )
11.99.2.3 ClampAmpRestart() void ClampAmpRestart ( )
11.99.2.4 DoRamp() void DoRamp (
            int startValue,
            int endValue,
            int duration,
            int mode )
11.99.2.5 Emu_GetCellCapacity() int Emu_GetCellCapacity ( )
11.99.2.6 Emu_GetCellPotential() int Emu_GetCellPotential ( )
```

```
11.99.2.7 Emu_GetCellResists() int Emu_GetCellResists ( )
\textbf{11.99.2.8} \quad \textbf{Emu\_GetElectrodeResists()} \quad \texttt{int} \ \texttt{Emu\_GetElectrodeResists()} \quad \texttt{int} \ \texttt{Emu\_GetElectrodeResists()}
11.99.2.9 Emu_GetNoise() int Emu_GetNoise ( )
11.99.2.10 Emu_SetCellCapacity() void Emu_SetCellCapacity (
               int emuCellCapacity )
11.99.2.11 Emu_SetCellPotential() void Emu_SetCellPotential (
               int emuCellPotential )
11.99.2.12 Emu_SetCellResists() void Emu_SetCellResists (
               int emuCellResist )
11.99.2.13 Emu_SetElectrodeResists() void Emu_SetElectrodeResists (
               int emuElectrodeResist )
11.99.2.14 Emu_SetNoise() void Emu_SetNoise (
               int emuNoise )
11.99.2.15 GetAllDigout() uint32_t GetAllDigout ( )
11.99.2.16 GetCapacityC() int GetCapacityC ( )
```

```
11.99.2.17 GetCapacityV() int GetCapacityV ( )
11.99.2.18 GetCapacityX() int GetCapacityX ( )
11.99.2.19 GetClampAmpSerialNumber() unsigned int GetClampAmpSerialNumber ( )
11.99.2.20 GetCommand() int GetCommand (
             unsigned char command )
11.99.2.21 GetConfigurationBit() bool GetConfigurationBit (
             unsigned short bit )
\textbf{11.99.2.22} \quad \textbf{GetConfigurationBitAxc()} \quad \texttt{bool GetConfigurationBitAxc ()}
11.99.2.23 GetConfigurationBitBlu_Led() bool GetConfigurationBitBlu_Led ( )
11.99.2.24 GetConfigurationBitBlu_LedToggleFast() bool GetConfigurationBitBlu_LedToggleFast ( )
11.99.2.25 GetConfigurationBitBlu_LedToggleSlow() bool GetConfigurationBitBlu_LedToggleSlow ( )
11.99.2.26 GetConfigurationBitCC_Gen() bool GetConfigurationBitCC_Gen ( )
11.99.2.27 GetConfigurationBitCV_Gen() bool GetConfigurationBitCV_Gen ( )
```

11.99.2.28	GetConfigurationBitRC_Gen() bool GetConfigurationBitRC_Gen ()
11.99.2.29	GetConfigurationBitRed_Led() bool GetConfigurationBitRed_Led ()
11.99.2.30	GetConfigurationBitRed_LedSaturation() bool GetConfigurationBitRed_LedSaturation ()
11.99.2.31	GetConfigurationBitRed_LedToggleFast() bool GetConfigurationBitRed_LedToggleFast ()
11.99.2.32	GetConfigurationBitRed_LedToggleSlow() bool GetConfigurationBitRed_LedToggleSlow (
11.99.2.33	GetConfigurationBitRelais() bool GetConfigurationBitRelais ()
11.99.2.34	GetConfigurationBitRV_Gen() bool GetConfigurationBitRV_Gen ()
11.99.2.35	GetConfigurationBits() unsigned int GetConfigurationBits ()
11.99.2.36	GetConfigurationBitStream() bool GetConfigurationBitStream ()
11.99.2.37	GetConfigurationBitSupply() bool GetConfigurationBitSupply ()
11.99.2.38	<pre>GetCrossTalkOffset() int GetCrossTalkOffset ()</pre>

```
11.99.2.39 GetCrossTalkOptimum() int GetCrossTalkOptimum ( )
\textbf{11.99.2.40} \quad \textbf{GetDigout()} \quad \texttt{bool GetDigout} \quad (
             uint16_t index )
11.99.2.41 GetDisplayText() String ^ GetDisplayText (
              int index )
11.99.2.42 GetDownsampleFactor() int GetDownsampleFactor (
              int index )
11.99.2.43 GetFilter() int GetFilter ( )
11.99.2.44 GetFilterCoeffs() array<int> ^{\wedge} GetFilterCoeffs (
             int index )
11.99.2.45 GetIC() int GetIC ( )
11.99.2.46 GetlClamp() int GetIClamp ()
11.99.2.47 GetICOffset() int GetICOffset ( )
11.99.2.48 GetlGain() int GetlGain ()
11.99.2.49 GetNIC_MS() int GetNIC_MS ()
```

```
\textbf{11.99.2.50} \quad \textbf{GetNUC\_MS()} \quad \texttt{int GetNUC\_MS ()}
11.99.2.51 GetNUV_MS() int GetNUV_MS ()
11.99.2.52 GetPGain() int GetPGain ()
11.99.2.53 GetRecordingNumber() unsigned int GetRecordingNumber ( )
11.99.2.54 GetResistanceC() int GetResistanceC ( )
11.99.2.55 GetResistanceV() int GetResistanceV ()
11.99.2.56 GetScreen() int GetScreen ()
11.99.2.57 GetSimulation() unsigned int GetSimulation ()
11.99.2.58 GetUC() int GetUC ( )
11.99.2.59 GetUClamp() int GetUClamp ()
11.99.2.60 GetUCOffset() int GetUCOffset ()
```

```
11.99.2.61 GetUpdateDisplay() int GetUpdateDisplay ( )
11.99.2.62 GetUV() int GetUV ( )
11.99.2.63 GetUVOffset() int GetUVOffset ( )
11.99.2.64 GetXGain() int GetXGain ( )
11.99.2.65 RunTable() [1/2] void RunTable ( )
11.99.2.66 RunTable() [2/2] void RunTable (
            int timeout )
11.99.2.67 SetAllDigout() void SetAllDigout (
            uint32_t value )
11.99.2.68 SetCommand() void SetCommand (
             unsigned char command,
             int value )
11.99.2.69 SetConfigurationBit() void SetConfigurationBit (
            unsigned short bit,
             bool value )
11.99.2.70 SetConfigurationBitAxc() void SetConfigurationBitAxc (
             bool value )
```

```
11.99.2.71 SetConfigurationBitBlu_Led() void SetConfigurationBitBlu_Led (
              bool value )
11.99.2.72 SetConfigurationBitBlu_LedToggleFast() void SetConfigurationBitBlu_LedToggleFast (
             bool value )
\textbf{11.99.2.73} \quad \textbf{SetConfigurationBitBlu\_LedToggleSlow()} \quad \texttt{void SetConfigurationBitBlu\_LedToggleSlow} \quad \textbf{(}
              bool value )
11.99.2.74 SetConfigurationBitCC_Gen() void SetConfigurationBitCC_Gen (
             bool value )
11.99.2.75 SetConfigurationBitCV_Gen() void SetConfigurationBitCV_Gen (
              bool value )
11.99.2.76 SetConfigurationBitRC_Gen() void SetConfigurationBitRC_Gen (
             bool value )
11.99.2.77 SetConfigurationBitRed_Led() void SetConfigurationBitRed_Led (
              bool value )
11.99.2.78 SetConfigurationBitRed LedSaturation() void SetConfigurationBitRed_LedSaturation (
              bool value )
11.99.2.79 SetConfigurationBitRed_LedToggleFast() void SetConfigurationBitRed_LedToggleFast (
             bool value )
11.99.2.80 SetConfigurationBitRed_LedToggleSlow() void SetConfigurationBitRed_LedToggleSlow (
              bool value )
```

```
\textbf{11.99.2.81} \quad \textbf{SetConfigurationBitRelais()} \quad \texttt{void SetConfigurationBitRelais} \quad \textbf{(}
                bool value )
11.99.2.82 SetConfigurationBitRV_Gen() void SetConfigurationBitRV_Gen (
                bool value )
\textbf{11.99.2.83} \quad \textbf{SetConfigurationBitStream()} \quad \texttt{void SetConfigurationBitStream} \ \ \textbf{(}
                bool value )
\textbf{11.99.2.84} \quad \textbf{SetConfigurationBitSupply()} \quad \texttt{void SetConfigurationBitSupply} \quad \textbf{(}
                bool value )
11.99.2.85 SetCrossTalkOffset() void SetCrossTalkOffset (
               int CrossTalk )
11.99.2.86 SetCrossTalkOptimum() void SetCrossTalkOptimum (
                int cxOptimum )
11.99.2.87 SetDigout() void SetDigout (
                uint16_t index,
                bool enable )
11.99.2.88 SetDisplayText() void SetDisplayText (
                int index,
                String^{\wedge} displayText)
11.99.2.89 SetDownsampleFactor() void SetDownsampleFactor (
                int index,
                int downsample_factor )
```

```
11.99.2.90 SetFilter() void SetFilter (
             int filter )
11.99.2.91 SetFilterCoeffs() void SetFilterCoeffs (
             int index,
              array< int >^{\land} coeffs )
11.99.2.92 SetlClamp() void SetIClamp (
             int iClamp )
11.99.2.93 SetICOffset() void SetICOffset (
             int ICOffset )
11.99.2.94 SetlGain() void SetIGain (
             int iGain )
11.99.2.95 SetNoFilterCoeffs() void SetNoFilterCoeffs (
             int index )
11.99.2.96 SetPGain() void SetPGain (
             int pGain )
\textbf{11.99.2.97} \quad \textbf{SetRecordingNumber()} \quad \texttt{void SetRecordingNumber ()}
              unsigned int recordingNumber )
11.99.2.98 SetScreen() void SetScreen (
            int screen )
```


```
\textbf{11.99.2.99} \quad \textbf{SetSimulation()} \quad \texttt{void SetSimulation ()}
             unsigned int enable )
11.99.2.100 SetUClamp() void SetUClamp (
             int uClamp )
11.99.2.101 SetUCOffset() void SetUCOffset (
             int UCOffset )
11.99.2.102 SetUVOffset() void SetUVOffset (
             int UVOffset )
11.99.2.103 SetXGain() void SetXGain (
             int xGain )
11.99.2.104 StopTable() [1/2] void StopTable ( )
11.99.2.105 StopTable() [2/2] void StopTable (
             int timeout )
11.99.2.106 Table_Wait() void Table_Wait (
             unsigned int tableWait )
11.99.2.107 TableDefBegin() void TableDefBegin ( )
11.99.2.108 TableDefEnd() void TableDefEnd ( )
```

11.99.2.109 UpdateDisplay() void UpdateDisplay ()

11.100 CRoboDeviceNet Class Reference

CRoboDeviceNet is the base class for all Robo platform based devices

Inheritance diagram for CRoboDeviceNet:



Classes

class RoboMainLowLevelCommands

Public Member Functions

- CRoboDeviceNet (void)
- ∼CRoboDeviceNet (void)
- void SetInMovement ()

Low level command, sets the internal state to "In Movement"

bool GetInMovement ()

Low level command, gets the internal state "In Movement"

uint32_t GetMovementError ()

Low level command, gets the error of the last movement end

- void FindReference (unsigned char busaddress, char axes)
- void FindReference (unsigned char busaddress, char axes, int timeout)

Searches the reference position of the motor

- void MoveAbs (unsigned char busaddress, char axes, int x, int y)
- void MoveAbs (unsigned char busaddress, char axes, int x, int y, int timeout)

Moves the motor to the new absolute position

- void StopMovement (unsigned char busaddress, char axes)
- void StopMovement (unsigned char busaddress, char axes, int timeout)

Stops the current movement

- void SetCurrentAndAir (unsigned char busaddress, char axes, unsigned short onoff)
- void SetCurrentAndAir (unsigned char busaddress, char axes, unsigned short onoff, int timeout)
- void CancelPoolLoop ()
- void CancelPoolLoopAndStopMovement ()
- void GetCurrentPosition (unsigned char busaddress, char axes, [System::Runtime::InteropServices::Out]int%
 x, [System::Runtime::InteropServices::Out]int% y)

Gets the current position of motors

- · void SetAirValve (unsigned int onoff)
- unsigned int GetAirValve ()
- unsigned int GetVoltageValves ()
- unsigned int GetVoltageRs485A ()
- unsigned int GetVoltageRs485B ()
- · unsigned int GetVoltageAirvalve ()
- unsigned int GetCurrentAirvalve ()
- unsigned int GetVoltage12V ()
- unsigned int GetAirpressure ()
- unsigned int GetVoltage5V ()
- unsigned int GetErrorVoltageValves ()
- unsigned int GetErrorVoltageRs485A ()
- unsigned int GetErrorVoltageRs485B ()
- unsigned int GetErrorVoltageAirvalve ()
- unsigned int GetErrorCurrentAirvalve ()
- unsigned int GetErrorVoltage12V ()
- unsigned int GetErrorAirpressure ()
- unsigned int GetErrorVoltage5V ()
- void SetVoltageValvesLimit (unsigned int lowervoltage, unsigned int uppervoltage)
- void SetVoltageRs485ALimit (unsigned int lowervoltage, unsigned int uppervoltage)
- void SetVoltageRs485BLimit (unsigned int lowervoltage, unsigned int uppervoltage)
- void SetVoltageAirvalveLimit (unsigned int lowervoltage, unsigned int uppervoltage)
- void SetCurrentAirvalveLimit (unsigned int lowercurrent, unsigned int uppercurrent)
- void SetVoltage12VLimit (unsigned int lowervoltage, unsigned int uppervoltage)
 void SetAirpressureLimit (unsigned int lowerpressure, unsigned int upperpressure)
- void SetVoltage5VLimit (unsigned int lowervoltage, unsigned int uppervoltage)
- void GetVoltageRs485ALimit ([System::Runtime::InteropServices::Out]unsigned int% lowervoltage, [System::Runtime::InteropServices::Out]unsigned int% uppervoltage)
- void GetVoltageRs485BLimit ([System::Runtime::InteropServices::Out]unsigned int% lowervoltage, [System::Runtime::InteropServices::Out]unsigned int% uppervoltage)
- void GetVoltageAirvalveLimit ([System::Runtime::InteropServices::Out]unsigned int% lowervoltage, [System::Runtime::InteropServices::Out]unsigned int% uppervoltage)
- void GetCurrentAirvalveLimit ([System::Runtime::InteropServices::Out]unsigned int% lowercurrent,
 [System::Runtime::InteropServices::Out]unsigned int% uppercurrent)
- void GetVoltage12VLimit ([System::Runtime::InteropServices::Out]unsigned int% lowervoltage, [System::
 Runtime::InteropServices::Out]unsigned int% uppervoltage)
- void GetAirpressureLimit ([System::Runtime::InteropServices::Out]unsigned int% lowerpressure, [System ← ::Runtime::InteropServices::Out]unsigned int% upperpressure)
- void GetVoltage5VLimit ([System::Runtime::InteropServices::Out]unsigned int% lowervoltage, [System::
 Runtime::InteropServices::Out]unsigned int% uppervoltage)
- void SetMinPressure (int pressure)
- int GetMinPressure ()

Static Public Attributes

```
    static const uint32 t RoboError Base = (0xA0110000L)

    static const uint32 t RoboError UnknownCommand = ( (0xA0110000L) )

static const uint32_t RoboError_Timeout = ( (0xA0110000L) | 0x0001 )

    static const uint32 t RoboError Pressure = ( (0xA0110000L) | 0x0002 )

    static const uint32 t RoboError RangeExceeded = ( (0xA0110000L) | 0x0003 )

    static const uint32_t RoboError_CommunicationTimeout = ( (0xA0110000L) | 0x0004 )

    static const uint32_t RoboError_AnotherMaster = ( (0xA0110000L) | 0x0005 )

• static const uint32 t RoboError FindReferenceMethod = ( (0xA0110000L) | 0x0006 )

    static const uint32 t RoboError NoSpeedOrAcceleration = ( (0xA0110000L) | 0x0007 )

    static const uint32 t RoboError NoEndSwitch = ( (0xA0110000L) | 0x0008 )

    static const uint32_t RoboError_CannotEscapeEndSwitch = ( (0xA0110000L) | 0x0009 )

    static const uint32 t RoboError CommandAlreadyInProgress = ( (0xA0110000L) | 0x000A )

    static const uint32_t RoboError_NoReference = ( (0xA0110000L) | 0x000B )

    static const uint32 t RoboError OverPressure = ( (0xA0110000L) | 0x000C )

    static const uint32 t RoboError Phase0OutOfRange = ( (0xA0110000L) | 0x000D )

    static const uint32_t RoboError_PeristalticTimeout = ( (0xA0110000L) | 0x000E )

    static const uint32 t RoboError GilsonTimeout = ( (0xA0110000L) | 0x000F )

    static const uint32_t RoboError_GilsonWrondID = ( (0xA0110000L) | 0x0010 )

    static const uint32 t RoboError GilsonCommandPending = ((0xA0110000L) | 0x0011)

    static const uint32 t RoboError ParameterNotAllowed = ( (0xA0110000L) | 0x0012 )

    static const uint32_t RoboError_StateChangeNotPossible = ( (0xA0110000L) | 0x0013 )

    static const uint32 t RoboError CommandNotPossible = ( (0xA0110000L) | 0x0014 )

    static const uint32_t RoboError_DacqNotReady = ( (0xA0110000L) | 0x0015 )

static const uint32_t RoboError_NoMoreData = ( (0xA0110000L) | 0x0016 )

    static const uint32 t RoboError McsBus UnknownCommand = ( (0xA0110000L) | 0x003F)

    static const uint32 t RoboError DLLMovementTimeout = ( (0xA0110000L) | 0x1001)

    static const uint32 t RoboError PollLoopCanceled = ( (0xA0110000L) | 0x1002)

    static const uint32_t RoboError_PollLoopCanceledAndStopMovement = ( (0xA0110000L) | 0x1003)

• static const byte McsBus_XY = 1
     McsBus address for the xy-plane

    static const byte McsBus ZI = 2

     McsBus address for the z and i axes
• static const byte Axis X = 0
     Axis number of x for axis argument
static const byte Axis_Y = 1
     Axis number of y for axis argument
• static const byte Axis Z = 0
     Axis number of z for axis argument

    static const byte Axis I = 1

     Axis number of i for axis argument

    static const char Axes X = 1

     Bit pattern for x axis for axes argument

    static const char Axes Y = 2

     Bit pattern for y axis for axes argument

    static const char Axes Z = 1

     Bit pattern for z axis for axes argument
static const char Axes_I = 2
     Bit pattern for i axis for axes argument
```

Properties

- CMcsBusNet^ McsBus [get]
- CMcsBus_MotorControlNet[^] McsBus_MotorControl [get]
- RoboMainLowLevelCommands RoboMainLowLevelCommand [get]

Events

RoboStatusEventDelegate[^] RoboStatusEvent

Additional Inherited Members

11.100.1 Detailed Description

CRoboDeviceNet is the base class for all Robo platform based devices

11.100.2 Constructor & Destructor Documentation

11.100.3 Member Function Documentation

```
11.100.3.1 CancelPoolLoop() void CancelPoolLoop ( )
```

11.100.3.2 CancelPoolLoopAndStopMovement() void CancelPoolLoopAndStopMovement ()

```
11.100.3.3 FindReference() [1/2] void FindReference (
          unsigned char busaddress,
          char axes )
```

```
11.100.3.4 FindReference() [2/2] void FindReference (
          unsigned char busaddress,
          char axes,
          int timeout )
```

Searches the reference position of the motor

Parameters

busaddress	Address of the McsBus
axes	Bit pattern of axes to drive
timeout	Timeout of maximal waiting for the end of the command (-1 is forever)

```
11.100.3.5 GetAirpressure() unsigned int GetAirpressure ()
11.100.3.6 GetAirpressureLimit() void GetAirpressureLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowerpressure,
              [{\tt System::Runtime::InteropServices::Out}] \ unsigned \ {\tt int} \$ \ upperpressure \ )
\textbf{11.100.3.7} \quad \textbf{GetAirValve()} \quad \texttt{unsigned int GetAirValve ()}
11.100.3.8 GetCurrentAirvalve() unsigned int GetCurrentAirvalve ()
11.100.3.9 GetCurrentAirvalveLimit() void GetCurrentAirvalveLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowercurrent,
              [System::Runtime::InteropServices::Out] unsigned int% uppercurrent )
11.100.3.10 GetCurrentPosition() void GetCurrentPosition (
              unsigned char busaddress,
              char axes,
              [System::Runtime::InteropServices::Out] int% x,
```

Gets the current position of motors

Parameters

busaddress	Address of the McsBus
axes	Bit pattern of axes to drive
X	Current position of first axis if pattern in axes is set
У	Current position of second axis if pattern in axes is set

[System::Runtime::InteropServices::Out] int% y)

11.100.3.11	GetErrorAirpressure() unsigned int GetErrorAirpressure ()
11.100.3.12	GetErrorCurrentAirvalve() unsigned int GetErrorCurrentAirvalve ()
11.100.3.13	GetErrorVoltage12V() unsigned int GetErrorVoltage12V ()
11.100.3.14	GetErrorVoltage5V() unsigned int GetErrorVoltage5V ()
11.100.3.15	GetErrorVoltageAirvalve() unsigned int GetErrorVoltageAirvalve ()
11.100.3.16	GetErrorVoltageRs485A() unsigned int GetErrorVoltageRs485A ()
11.100.3.17	<pre>GetErrorVoltageRs485B() unsigned int GetErrorVoltageRs485B ()</pre>
11.100.3.18	GetErrorVoltageValves() unsigned int GetErrorVoltageValves ()
	GetInMovement() bool GetInMovement () nmand, gets the internal state "In Movement"
11.100.3.20	GetMinPressure() int GetMinPressure ()

```
\textbf{11.100.3.21} \quad \textbf{GetMovementError()} \quad \texttt{uint32\_t} \quad \texttt{GetMovementError} \quad \textbf{( )}
Low level command, gets the error of the last movement end
11.100.3.22 GetVoltage12V() unsigned int GetVoltage12V ()
11.100.3.23 GetVoltage12VLimit() void GetVoltage12VLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowervoltage,
              [System::Runtime::InteropServices::Out] unsigned int% uppervoltage )
11.100.3.24 GetVoltage5V() unsigned int GetVoltage5V ()
11.100.3.25 GetVoltage5VLimit() void GetVoltage5VLimit (
              [{\tt System::Runtime::InteropServices::Out}] \ {\tt unsigned int \$} \ {\tt lowervoltage,}
              [System::Runtime::InteropServices::Out] unsigned int% uppervoltage )
11.100.3.26 GetVoltageAirvalve() unsigned int GetVoltageAirvalve ()
11.100.3.27 GetVoltageAirvalveLimit() void GetVoltageAirvalveLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowervoltage,
              [System::Runtime::InteropServices::Out] unsigned int% uppervoltage )
11.100.3.28 GetVoltageRs485A() unsigned int GetVoltageRs485A ( )
11.100.3.29 GetVoltageRs485ALimit() void GetVoltageRs485ALimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowervoltage,
              [System::Runtime::InteropServices::Out] unsigned int% uppervoltage )
```

```
11.100.3.30 GetVoltageRs485B() unsigned int GetVoltageRs485B ( )
11.100.3.31 GetVoltageRs485BLimit() void GetVoltageRs485BLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowervoltage,
              [{\tt System::Runtime::InteropServices::Out}] \  \, unsigned \  \, int {\tt \textit{uppervoltage}} \  \, )
11.100.3.32 GetVoltageValves() unsigned int GetVoltageValves ()
11.100.3.33 GetVoltageValvesLimit() void GetVoltageValvesLimit (
              [System::Runtime::InteropServices::Out] unsigned int% lowervoltage,
              [System::Runtime::InteropServices::Out] unsigned int% uppervoltage )
11.100.3.34 MoveAbs() [1/2] void MoveAbs (
             unsigned char busaddress,
             char axes,
             int x,
             int y)
11.100.3.35 MoveAbs() [2/2] void MoveAbs (
             unsigned char busaddress,
             char axes,
             int x,
             int y,
             int timeout )
```

Moves the motor to the new absolute position

Parameters

busaddress	Address of the McsBus
axes	Bit pattern of axes to drive
X	Position of first axis, if pattern in axes is set
У	Position of second axis if pattern in axes is set
timeout	Timeout of maximal waiting for the end of the command (-1 is forever)

11.100.3.36 SetAirpressureLimit() void SetAirpressureLimit (

```
unsigned int lowerpressure,
             unsigned int upperpressure )
11.100.3.37 SetAirValve() void SetAirValve (
             unsigned int onoff )
11.100.3.38 SetCurrentAirvalveLimit() void SetCurrentAirvalveLimit (
             unsigned int lowercurrent,
             unsigned int uppercurrent )
11.100.3.39 SetCurrentAndAir() [1/2] void SetCurrentAndAir (
             unsigned char busaddress,
             char axes,
             unsigned short onoff )
11.100.3.40 SetCurrentAndAir() [2/2] void SetCurrentAndAir (
             unsigned char busaddress,
             char axes,
             unsigned short onoff,
             int timeout )
11.100.3.41 SetInMovement() void SetInMovement ()
Low level command, sets the internal state to "In Movement"
11.100.3.42 SetMinPressure() void SetMinPressure (
             int pressure )
11.100.3.43 SetVoltage12VLimit() void SetVoltage12VLimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
```

```
11.100.3.44 SetVoltage5VLimit() void SetVoltage5VLimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
11.100.3.45 SetVoltageAirvalveLimit() void SetVoltageAirvalveLimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
11.100.3.46 SetVoltageRs485ALimit() void SetVoltageRs485ALimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
11.100.3.47 SetVoltageRs485BLimit() void SetVoltageRs485BLimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
11.100.3.48 SetVoltageValvesLimit() void SetVoltageValvesLimit (
             unsigned int lowervoltage,
             unsigned int uppervoltage )
11.100.3.49 StopMovement() [1/2] void StopMovement (
             unsigned char busaddress,
             char axes )
11.100.3.50 StopMovement() [2/2] void StopMovement (
             unsigned char busaddress,
             char axes,
             int timeout )
```

Stops the current movement

Parameters

busaddress	Address of the McsBus
axes	Bit pattern of axes to drive
timeout	Timeout of maximal waiting for the end of the command (-1 is forever)

11.100.4 Member Data Documentation

11.100.4.1 Axes_I const char Axes_I = 2 [static]

Bit pattern for i axis for axes argument

11.100.4.2 Axes_X const char Axes_X = 1 [static]

Bit pattern for x axis for axes argument

11.100.4.3 Axes_Y const char Axes_Y = 2 [static]

Bit pattern for y axis for axes argument

11.100.4.4 Axes_Z const char $Axes_Z = 1$ [static]

Bit pattern for z axis for axes argument

11.100.4.5 Axis_I const byte Axis_I = 1 [static]

Axis number of i for axis argument

11.100.4.6 Axis_X const byte $Axis_X = 0$ [static]

Axis number of x for axis argument

11.100.4.7 Axis_Y const byte Axis_Y = 1 [static]

Axis number of y for axis argument

```
11.100.4.8 Axis_Z const byte Axis_Z = 0 [static]
Axis number of z for axis argument
11.100.4.9 McsBus_XY const byte McsBus_XY = 1 [static]
McsBus address for the xy-plane
11.100.4.10 McsBus_ZI const byte McsBus_ZI = 2 [static]
McsBus address for the z and i axes
11.100.4.11 RoboError_AnotherMaster const uint32_t RoboError_AnotherMaster = ( (0xA0110000L) |
0x0005 ) [static]
11.100.4.12 RoboError_Base const uint32_t RoboError_Base = (0xA0110000L) [static]
11.100.4.13 RoboError_CannotEscapeEndSwitch const uint32_t RoboError_CannotEscapeEndSwitch =
( (0xA0110000L) | 0x0009 ) [static]
11.100.4.14 RoboError_CommandAlreadyInProgress const uint32_t RoboError_CommandAlreadyIn↔
Progress = ( (0xA0110000L) | 0x000A ) [static]
\textbf{11.100.4.15} \quad \textbf{RoboError\_CommandNotPossible} \quad \texttt{const uint32\_t RoboError\_CommandNotPossible} = (
```

11.100.4.16 RoboError_CommunicationTimeout const uint32_t RoboError_CommunicationTimeout = (

Generated by Doxygen

(0xA0110000L) | 0x0014) [static]

(0xA0110000L) | 0x0004) [static]

```
11.100.4.17 RoboError_DacqNotReady const uint32_t RoboError_DacqNotReady = ( (0xA0110000L) |
0x0015 ) [static]
11.100.4.18 RoboError_DLLMovementTimeout const uint32_t RoboError_DLLMovementTimeout = (
(0xA0110000L) | 0x1001) [static]
11.100.4.19 RoboError_FindReferenceMethod const uint32_t RoboError_FindReferenceMethod = (
(0xA0110000L) | 0x0006 ) [static]
11.100.4.20 RoboError_GilsonCommandPending const uint32_t RoboError_GilsonCommandPending =
( (0xA0110000L) | 0x0011 ) [static]
11.100.4.21 RoboError_GilsonTimeout const uint32_t RoboError_GilsonTimeout = ( (0xA0110000L) |
0x000F ) [static]
11.100.4.22 RoboError_GilsonWrondID const uint32_t RoboError_GilsonWrondID = ( (0xA0110000L)
| 0x0010 ) [static]
11.100.4.23 RoboError_McsBus_UnknownCommand const uint32_t RoboError_McsBus_Unknown←
Command = ( (0xA0110000L) | 0x003F) [static]
11.100.4.24 RoboError NoEndSwitch const uint32_t RoboError_NoEndSwitch = ( (0xA0110000L) |
0x0008 ) [static]
11.100.4.25 RoboError_NoMoreData const uint32_t RoboError_NoMoreData = ( (0xA0110000L) |
0x0016 ) [static]
11.100.4.26 RoboError_NoReference const uint32_t RoboError_NoReference = ( (0xA0110000L) |
0x000B ) [static]
```

```
11.100.4.27 RoboError_NoSpeedOrAcceleration const uint32_t RoboError_NoSpeedOrAcceleration =
  ( (0xA0110000L) | 0x0007 ) [static]
11.100.4.28 RoboError_OverPressure const uint32_t RoboError_OverPressure = ( (0xA0110000L) |
0x000C ) [static]
11.100.4.29 RoboError_ParameterNotAllowed const uint32_t RoboError_ParameterNotAllowed = (
  (0xA0110000L) | 0x0012 ) [static]
11.100.4.30 RoboError_PeristalticTimeout const uint32_t RoboError_PeristalticTimeout = ( (0x←
A0110000L) | 0x000E ) [static]
\textbf{11.100.4.31} \quad \textbf{RoboError\_Phase0OutOfRange} \quad \texttt{const uint} \\ 32\_t \\ \ \texttt{RoboError\_Phase0OutOfRange} \\ = \text{ ( } \texttt{(0x} \\ \leftarrow \texttt{(0
A0110000L) | 0x000D ) [static]
11.100.4.32 RoboError_PollLoopCanceled const uint32_t RoboError_PollLoopCanceled = ( (0x←)
A0110000L) | 0x1002) [static]
11.100.4.33 RoboError_PollLoopCanceledAndStopMovement const uint32_t RoboError_PollLoop←
CanceledAndStopMovement = ((0xA0110000L) | 0x1003) [static]
11.100.4.34 RoboError Pressure const uint32_t RoboError_Pressure = ( (0xA0110000L) | 0x0002 )
 [static]
11.100.4.35 RoboError_RangeExceeded const uint32_t RoboError_RangeExceeded = ( (0xA0110000L)
| 0x0003 ) [static]
11.100.4.36 RoboError_StateChangeNotPossible const uint32_t RoboError_StateChangeNotPossible
= ( (0xA0110000L) | 0x0013 ) [static]
```

11.100.4.37 RoboError_Timeout const uint32_t RoboError_Timeout = ((0xA0110000L) | 0x0001) [static]

11.100.4.38 RoboError_UnknownCommand const uint32_t RoboError_UnknownCommand = ($(0x \leftarrow A0110000L)$) [static]

11.100.5 Property Documentation

11.100.5.1 McsBus CMcsBusNet^ McsBus [get]

11.100.5.2 McsBus_MotorControl CMcsBus_MotorControlNet^ McsBus_MotorControl [get]


11.100.5.3 RoboMainLowLevelCommand RoboMainLowLevelCommands^ RoboMainLowLevelCommand [get]

11.100.6 Event Documentation

11.100.6.1 RoboStatusEvent RoboStatusEventDelegate^ RoboStatusEvent

11.101 CRoboFluidDeviceNet Class Reference

Inheritance diagram for CRoboFluidDeviceNet:



Public Member Functions

- CRoboFluidDeviceNet (void)
- ∼CRoboFluidDeviceNet (void)
- void SetValve (int value)

Open or Close valves.

void SetSingleValve (int valve, bool onoff)

Opens or Closes a valve.

• int GetValve ()

Query the state of the values.

bool GetSingleValve (int valve)

Query the state of a valve.

- void CloseAllValves ()
- void PumpOn (int index, short speed)
- void SetPumpSpeed (int index, short speed)
- void PumpOff (int index)
- short GetPumpSpeed (int index)
- bool IsPumpMotorOn (int index)

Protected Attributes

- CRoboFluidDevice * m_pRoboFluidDevice
- CMcsBus_MotorControlNet ^ m_pMcsBus_MotorControlNet

Properties

CMcsBus MotorControlNet[^] McsBus MotorControl [get]

Additional Inherited Members

11.101.1 Constructor & Destructor Documentation

11.101.2 Member Function Documentation

```
\textbf{11.101.2.1} \quad \textbf{CloseAllValves()} \quad \texttt{void CloseAllValves ()}
11.101.2.2 GetPumpSpeed() short GetPumpSpeed (
              int index )
11.101.2.3 GetSingleValve() bool GetSingleValve (
              int valve )
Query the state of a valve.
Parameters
 valve
        number of valve /*!
Returns
     state of the valve
11.101.2.4 GetValve() int GetValve ( )
Query the state of the values.
Returns
     the current state of the valves as a bit pattern.
11.101.2.5 IsPumpMotorOn() bool IsPumpMotorOn (
              int index )
11.101.2.6 PumpOff() void PumpOff (
              int index )
11.101.2.7 PumpOn() void PumpOn (
              int index,
              short speed )
```

Opens or Closes a valve.

Parameters

valve	number of valve to be changed /*!
onoff	open or close the valve

```
11.101.2.10 SetValve() void SetValve ( int value )
```

Open or Close valves.

Parameters

value	bit pattern of valves which should be open.
-------	---

11.101.3 Member Data Documentation

11.101.3.1 m_pMcsBus_MotorControlNet CMcsBus_MotorControlNet ^ m_pMcsBus_MotorControlNet [protected]

11.101.3.2 m_pRoboFluidDevice CRoboFluidDevice* m_pRoboFluidDevice [protected]


11.101.4 Property Documentation

11.101.4.1 McsBus_MotorControl CMcsBus_MotorControlNet^ McsBus_MotorControl [get]

11.102 CRobolnjectDeviceNet Class Reference

CRobolnjectDeviceNet is the to control the MCS Robolnject device

Inheritance diagram for CRobolnjectDeviceNet:



Public Member Functions

CRobolnjectDeviceNet (void)

Additional Inherited Members

11.102.1 Detailed Description


CRobolnjectDeviceNet is the to control the MCS Robolnject device

11.102.2 Constructor & Destructor Documentation

11.103 CRoboocyte2DeviceNet Class Reference

CRoboocyte2DeviceNet is the class to control the MCS Roboocyte2 device

Inheritance diagram for CRoboocyte2DeviceNet:



Public Member Functions

- CRoboocyte2DeviceNet (void)
- void SetAxisLED (bool onoff)
- bool GetAxisLED ()
- CRoboDacqNet ^ GetRoboDacq ()
- CRoboFluidDeviceNet ^ GetRoboFluidDevice ()
- CGilsonDeviceNet ^ GetGilsonDevice ()
- CMcsBus_ExtensionNet ^ GetMcsBus_Extension ()

Additional Inherited Members

11.103.1 Detailed Description

CRoboocyte2DeviceNet is the class to control the MCS Roboocyte2 device

11.103.2 Constructor & Destructor Documentation

```
11.103.2.1 CRoboocyte2DeviceNet() CRoboocyte2DeviceNet (
void )
```

11.103.3 Member Function Documentation

```
11.103.3.1 GetAxisLED() bool GetAxisLED ( )
```

```
11.103.3.2 GetGilsonDevice() CGilsonDeviceNet ^ GetGilsonDevice ( )
```

```
11.103.3.3 GetMcsBus_Extension() CMcsBus_ExtensionNet ^ GetMcsBus_Extension ( )
```


11.103.3.4 GetRoboDacq() CRoboDacqNet ^ GetRoboDacq ()

11.103.3.5 GetRoboFluidDevice() CRoboFluidDeviceNet ^ GetRoboFluidDevice ()

```
11.103.3.6 SetAxisLED() void SetAxisLED (
bool onoff)
```

11.104 CRoboStatorDeviceNet Class Reference

Inheritance diagram for CRoboStatorDeviceNet:



Classes

· class RoboMainStatorLowLevelCommands

Public Member Functions

- CRoboStatorDeviceNet (void)
- void FindReferenceXY ()
- void FindReferenceXY (int timeout)
- void FindReferenceZ ()
- void FindReferenceZ (int timeout)
- void FindReferencel ()
- void FindReferencel (int timeout)
- unsigned char HasRefXY ()
- unsigned char HasRefZ ()
- unsigned char HasRefl ()
- void MoveAbsXY (int x, int y)
- void MoveAbsXY (int x, int y, int timeout)
- void MoveAbsZ (int z)
- void MoveAbsZ (int z, int timeout)
- void MoveAbsI (int i)
- void MoveAbsI (int i, int timeout)
- void StopMovementXY ()
- void StopMovementXY (int timeout)
- void StopMovementZ ()
- void StopMovementZ (int timeout)
- void StopMovementI ()
- void StopMovementI (int timeout)
- void SetCurrentAndAirXY (unsigned short onoff)

- void SetCurrentAndAirXY (unsigned short onoff, int timeout)
- void GetCurrentPositionXY ([System::Runtime::InteropServices::Out]int% x, [System::Runtime::Interop←
 Services::Out]int% y)
- void GetCurrentPositionZ ([System::Runtime::InteropServices::Out]int% z)
- void GetCurrentPositionI ([System::Runtime::InteropServices::Out]int% i)
- void SetVelocityXY (int v)
- void SetVelocityZ (int v)
- void SetVelocityI (int v)
- void SetSpeedXY (int v)
- void SetSpeedZ (int v)
- void SetSpeedI (int v)
- void SetSpeedNativeXY (int v)
- void SetSpeedNativeZ (int v)
- void SetSpeedNativel (int v)
- void SetAccelerationXY (int a)
- void SetAccelerationZ (int a)
- void SetAccelerationI (int a)
- void SetAccelerationNativeXY (int a)
- void SetAccelerationNativeZ (int a)
- void SetAccelerationNativel (int a)

Properties

• RoboMainStatorLowLevelCommands RoboMainStatorLowLevelCommand [get]

Additional Inherited Members

11.104.1 Constructor & Destructor Documentation

```
11.104.1.1 CRoboStatorDeviceNet() CRoboStatorDeviceNet (
```

11.104.2 Member Function Documentation

```
11.104.2.1 FindReferencel() [1/2] void FindReferenceI ( )
```

```
11.104.2.2 FindReferencel() [2/2] void FindReferenceI ( int timeout )
```

```
11.104.2.3 FindReferenceXY() [1/2] void FindReferenceXY ( )
11.104.2.4 FindReferenceXY() [2/2] void FindReferenceXY (
             int timeout )
11.104.2.5 FindReferenceZ() [1/2] void FindReferenceZ ( )
11.104.2.6 FindReferenceZ() [2/2] void FindReferenceZ (
             int timeout )
11.104.2.7 GetCurrentPositionI() void GetCurrentPositionI (
             [System::Runtime::InteropServices::Out] int% i )
11.104.2.8 GetCurrentPositionXY() void GetCurrentPositionXY (
             [System::Runtime::InteropServices::Out] int% x,
             [System::Runtime::InteropServices::Out] int% y )
11.104.2.9 GetCurrentPositionZ() void GetCurrentPositionZ (
             [System::Runtime::InteropServices::Out] int% z )
11.104.2.10 HasRefl() unsigned char HasRefI ()
11.104.2.11 HasRefXY() unsigned char HasRefXY ( )
11.104.2.12 HasRefZ() unsigned char HasRefZ ( )
```

```
11.104.2.13 MoveAbsI() [1/2] void MoveAbsI (
             int i)
11.104.2.14 MoveAbsI() [2/2] void MoveAbsI (
             int i,
             int timeout )
11.104.2.15 MoveAbsXY() [1/2] void MoveAbsXY (
             int x_{i}
             int y )
11.104.2.16 MoveAbsXY() [2/2] void MoveAbsXY (
             int x,
              int y,
             int timeout )
11.104.2.17 MoveAbsZ() [1/2] void MoveAbsZ (
             int z )
11.104.2.18 MoveAbsZ() [2/2] void MoveAbsZ (
             int z,
             int timeout )
11.104.2.19 SetAccelerationI() void SetAccelerationI (
             int a )
\textbf{11.104.2.20} \quad \textbf{SetAccelerationNativeI()} \quad \texttt{void SetAccelerationNativeI} \quad (
             int a)
11.104.2.21 SetAccelerationNativeXY() void SetAccelerationNativeXY (
             int a)
```

```
11.104.2.22 SetAccelerationNativeZ() void SetAccelerationNativeZ (
              int a)
11.104.2.23 SetAccelerationXY() void SetAccelerationXY (
             int a)
11.104.2.24 SetAccelerationZ() void SetAccelerationZ (
             int a)
11.104.2.25 SetCurrentAndAirXY() [1/2] void SetCurrentAndAirXY (
              unsigned short onoff )
11.104.2.26 SetCurrentAndAirXY() [2/2] void SetCurrentAndAirXY (
              unsigned short onoff,
              int timeout )
11.104.2.27 SetSpeedI() void SetSpeedI (
             int v)
11.104.2.28 SetSpeedNativel() void SetSpeedNativeI (
              int v)
11.104.2.29 SetSpeedNativeXY() void SetSpeedNativeXY (
             int v )
\textbf{11.104.2.30} \quad \textbf{SetSpeedNativeZ()} \quad \texttt{void SetSpeedNativeZ} \quad \textbf{(}
              int v)
```

```
\textbf{11.104.2.31} \quad \textbf{SetSpeedXY()} \quad \texttt{void SetSpeedXY} \ \ (
             int \ v )
11.104.2.32 SetSpeedZ() void SetSpeedZ (
             int v)
11.104.2.33 SetVelocityI() void SetVelocityI (
             int v)
11.104.2.34 SetVelocityXY() void SetVelocityXY (
             int v)
11.104.2.35 SetVelocityZ() void SetVelocityZ (
             int v)
11.104.2.36 StopMovementI() [1/2] void StopMovementI ( )
11.104.2.37 StopMovementI() [2/2] void StopMovementI (
             int timeout )
11.104.2.38 StopMovementXY() [1/2] void StopMovementXY ( )
11.104.2.39 StopMovementXY() [2/2] void StopMovementXY (
             int timeout )
11.104.2.40 StopMovementZ() [1/2] void StopMovementZ ( )
```


11.104.2.41 StopMovementZ() [2/2] void StopMovementZ (int timeout)

11.104.3 Property Documentation

11.104.3.1 RoboMainStatorLowLevelCommand RoboMainStatorLowLevelCommands ↑ RoboMainStator← LowLevelCommand [get]

11.105 CSafeISDeviceNet Class Reference

Inheritance diagram for CSafeISDeviceNet:



Public Member Functions

CSafeISDeviceNet (void)

Initializes a new instance of the CSafeISDeviceNet class.

~CSafeISDeviceNet (void)

Releases unmanaged resources and performs other cleanup operations before the CSafelSDeviceNet is reclaimed by garbage collection.

void SetSwitches (unsigned short switches)

Sets the switches for all electrodes on the device. Do not use during measurement

• void SetAdcChannels (unsigned char channels)

Sets the ADC channels you want to be sampled

void SetAdcSamplePos (array< unsigned short >^ positions)

Sets the sample position of the ADC.

• void SetDacMode (unsigned char mode)

Sets the DAC mode.

void SetDacPulseform (array< short >^ pulseform)

Sets the DAC pulseform.

• void SetDacPeriode (unsigned int periode)

Sets the DAC periode.

Properties

CRoboDeviceNet[^] RoboDevice [get]

Gets the CRoboDeviceNet. Use this to control the syringe.

CFluidControlDeviceNet[^] FluidControlDevice [get]

Gets the CFluidControlDeviceNet. Use this to control the valves. Only SetSingleValve is implemented for CSafeISDeviceNet.

CMcsUsbDacqNet[^] DacqDevice [get]

Gets the CMcsUsbDacqNet. Use this to control the data aquisition.

Additional Inherited Members

11.105.1 Detailed Description

11.105.2 Constructor & Destructor Documentation

```
11.105.2.1 CSafeISDeviceNet() CSafeISDeviceNet (
```

Initializes a new instance of the CSafeISDeviceNet class.

```
11.105.2.2 ~CSafeISDeviceNet() ~CSafeISDeviceNet (
```

Releases unmanaged resources and performs other cleanup operations before the CSafeISDeviceNet is reclaimed by garbage collection.

11.105.3 Member Function Documentation

```
11.105.3.1 SetAdcChannels() void SetAdcChannels (
unsigned char channels )
```

Sets the ADC channels you want to be sampled

Parameters

channels The bitmap of the 8 channels. Set bit to 1 for the channels you want measure

Sets the sample position of the ADC.

Parameters

positions | The positions in units of 0.1μs.

```
11.105.3.3 SetDacMode() void SetDacMode ( unsigned char mode )
```

Sets the DAC mode.

Parameters

```
mode The mode: 0 = Impedance; 1 = Amperometry
```

Sets the DAC periode.

Parameters

periode The periode in units of 10μs.

Sets the DAC pulseform.

Parameters

pulseform The pulseform.

Sets the switches for all electrodes on the device. Do not use during measurement

Parameters

switches The switches: See Schematics for the meaning

11.105.4 Property Documentation

```
11.105.4.1 DacqDevice CMcsUsbDacqNet^ DacqDevice [get]
```

Gets the CMcsUsbDacqNet. Use this to control the data aquisition.

```
11.105.4.2 FluidControlDevice CFluidControlDeviceNet^ FluidControlDevice [get]
```


Gets the CFluidControlDeviceNet. Use this to control the valves. Only SetSingleValve is implemented for CSafeISDeviceNet.

```
11.105.4.3 RoboDevice CRoboDeviceNet^ RoboDevice [get]
```

Gets the CRoboDeviceNet. Use this to control the syringe.

11.106 CSCUDacqGroupChannelSelectionNet Class Reference

Inheritance diagram for CSCUDacqGroupChannelSelectionNet:



Public Member Functions

CSCUDacqGroupChannelSelectionNet (CMcsUsbNet[^] mcsusb)

Additional Inherited Members


11.106.1 Constructor & Destructor Documentation

11.106.1.1 CSCUDacqGroupChannelSelectionNet() CSCUDacqGroupChannelSelectionNet (CMcsUsbNet^ mcsusb)

11.107 CSCUFunctionNet Class Reference

CSCUFunctionNet is the class to control the SCU device

Inheritance diagram for CSCUFunctionNet:



Public Member Functions

- delegate void OnGetAvailableHeadstages (uint32 t AvailableHeadstages)
- delegate void OnlsHeadstageAvailable (uint32_t Headstage, bool available)
- CSCUFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pSCUFunctionPointer
 — Container)

Initializes a new instance of the CSCUFunctionNet class.

- CSCUFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ∼CSCUFunctionNet ()
- !CSCUFunctionNet ()
- uint32_t GetAvailableHeadstages ()

Gets a bitmap of available headstages.

bool IsInDacqLegacyMode ()

Is the SCU in legacy mode

void SetDacqLegacyMode (bool enable)

Enable the SCU legacy mode

• uint32_t GetMaxStimulusChannelsPerHeadstage ()

Gets the maximal number of stimulation channels a headstage can have.

uint32_t GetMaxNumberOfHeadstages ()

Gets the maximal number of headstages.

SCU_HeadstageIdEnumNet GetHeadstageID (uint32_t Headstage)

Gets the headstage fpga ID.

bool IsHeadstageAvailable (uint32_t Headstage)

Checks whether the given headstage is available.

void PowerHS (uint32 t Headstage, bool power)

Power the HS

bool IsHSPowered (uint32_t Headstage)

Is the HS powered

bool HasHSPowerSwitch ()

Has SCU HS power switch

String \(^\) GetHeadstageSerialNumber (uint32 t Headstage)

Gets the serial number of a given headstage.

uint32_t GetHeadstageNumberOfAnalogChannels (uint32_t Headstage)

Gets the number of analog channels for a given headstage.

• uint32_t GetHeadstageNumberOfStimulationChannels (uint32_t Headstage)

Gets the number of stimulation channels for a given headstage.

uint32_t GetHeadstageGainInPermille (uint32_t Headstage)

Gets the gain factor in permille for a given headstage.

uint32 t GetHeadstageAdcRangeInMicroVolt (uint32 t Headstage)

Gets the ADC Range in uV for a given headstage.

uint32_t GetHeadstageAdcBits (uint32_t Headstage)

Gets the Number of ADC bits for a given headstage.

uint32_t GetHeadstageDacVoltageRangeInMilliVolt (uint32_t Headstage)

Gets the DAC Voltage Range in mV for a given headstage.

uint32 t GetHeadstageDacVoltageResolutionInMicroVolt (uint32 t Headstage)

Gets the DAC Voltage Resolution in uV for a given headstage.

uint32_t GetHeadstageDacCurrentRangeInMicroAmpere (uint32_t Headstage)

Gets the DAC Current Range in uA for a given headstage.

• uint32_t GetHeadstageDacCurrentResolutionInNanoAmpere (uint32_t Headstage)

Gets the DAC Current Resolution in nA for a given headstage.

• uint32_t GetHeadstageDacBits (uint32_t Headstage)

Gets the Number of DAC bits for a given headstage.

uint32 t GetHeadstageSamplerate (uint32 t Headstage)

Gets the Samplerate of a given headstage.

bool GetHeadstagePowerStateAtStart (uint32 t Headstage)

Gets the Power Status at SCU Power on of a given headstage.

void SetHeadstagePowerStateAtStart (uint32 t Headstage, bool Powerstatus)

Sets the Power Status at SCU Power on of a given headstage.

bool HasGalvanicIsolation ()

Has galvanic isolated hardware

• bool HasAnalogOut ()

Has AnalogOut hardware

void EnableAnalogOut (bool enable)

Enables AnalogOut globally

bool IsAnalogOutEnabled ()

Is AnalogOut enabled

void SetAnalogOutDACRange (AnalogOut DAC Range EnumNet range)

Sets the analog out DAC range

AnalogOut_DAC_Range_EnumNet GetAnalogOutDACRange ()

Gets the analog out DAC range

void SetAnalogOutADCRange (uint32 t range)

Sets the analog out ADC range

uint32_t GetAnalogOutADCRange ()

Gets the analog out ADC range

void AutomaticAnalogOut (bool automatic)

Sets automatic source channel selection

bool IsAutomaticAnalogOut ()

Is Automatic source channel selection selected

void SetAnalogOutChannels (uint32_t out_channel, uint32_t source_channel)

Set the source channel number for a certain output channel

uint32_t GetAnalogOutChannels (uint32_t out_channel)

Get the connected source channel number for a certain output channel

 void SetReferenceElectrodeSwitchState (uint32_t Headstage, ReferenceElectrodeSwitchPositionEnumNet NewSwitchPos)

Sets the position of the switch for the reference electrode

• ReferenceElectrodeSwitchPositionEnumNet GetReferenceElectrodeSwitchState (uint32_t Headstage)

Gets the position of the switch for the reference electrode

void SetReferenceElectrodeMode (uint32_t Headstage, ReferenceElectrodeModeEnumNet NewValue)

Sets the mode for the reference electrode

- ReferenceElectrodeModeEnumNet GetReferenceElectrodeMode (uint32_t Headstage)
 - Gets the mode for the reference electrode
- CFilterPropertyNet ^ GetFilterProperty (SCUDacqGroupChannelEnumNet GroupID, uint32_t FilterNumber)
 Gets the filter property
- array< CFilterPropertyNet[^]> [^] GetFilterProperties (SCUDacqGroupChannelEnumNet GroupID, int filter
 — Configurations_Length)

Gets multiple filter properties

Events

- OnGetAvailableHeadstages^ GetAvailableHeadstagesEvent [add, remove, raise]

 Event fires when the bitmap of available headstages has changed
- OnlsHeadstageAvailable IsHeadstageAvailableEvent [add, remove, raise]

 Event fires when 'true' if the headstage is connected for the headstage to guery has changed

Additional Inherited Members

11.107.1 Detailed Description

CSCUFunctionNet is the class to control the SCU device

11.107.2 Constructor & Destructor Documentation

```
11.107.2.1 CSCUFunctionNet() [1/2] CSCUFunctionNet (

CMcsUsbNet^ mcsusb,

CMcsUsbFunctionPointerContainer^ pSCUFunctionPointerContainer)
```

Initializes a new instance of the CSCUFunctionNet class.

```
11.107.2.2 CSCUFunctionNet() [2/2] CSCUFunctionNet (

CMcsUsbNet^ mcsusb )

11.107.2.3 ~CSCUFunctionNet() virtual ~CSCUFunctionNet ( ) [virtual]
```

11.107.3 Member Function Documentation

11.107.2.4 "!CSCUFunctionNet() !CSCUFunctionNet ()

Sets automatic source channel selection

```
{\bf 11.107.3.2} \quad {\bf Enable Analog Out ()} \quad {\tt void Enable Analog Out} \quad (
```

bool enable)

Enables AnalogOut globally

Parameters

enable Enable

11.107.3.3 GetAnalogOutADCRange() uint32_t GetAnalogOutADCRange ()

Gets the analog out ADC range

Returns

Range

Get the connected source channel number for a certain output channel

Parameters

out_channel	Output channel number

Returns

Source channel number

11.107.3.5 GetAnalogOutDACRange() AnalogOut_DAC_Range_EnumNet GetAnalogOutDACRange ()

Gets the analog out DAC range

Returns

Range

$\textbf{11.107.3.6} \quad \textbf{GetAvailableHeadstages()} \quad \texttt{uint32_t GetAvailableHeadstages ()}$

Gets a bitmap of available headstages.

Returns

The bitmap of available headstages.

Gets multiple filter properties

Parameters

GroupID	The group ID
filterConfigurations_Length	The maximal length of filterConfigurations.

Returns

array of filter properties

Gets the filter property

Parameters

GroupID	The group ID
FilterNumber	The filter number

Returns

The filter property

```
11.107.3.9 GetHeadstageAdcBits() uint32_t GetHeadstageAdcBits ( uint32_t Headstage )
```

Gets the Number of ADC bits for a given headstage.

Headstage The headstage to query.

Returns

The number of bits the ADC has for the given headstage.

11.107.3.10 GetHeadstageAdcRangeInMicroVolt() uint32_t GetHeadstageAdcRangeInMicroVolt (uint32_t Headstage)

Gets the ADC Range in uV for a given headstage.

Parameters

Headstage The headst	age to query.
----------------------	---------------

Returns

The ADC Range in uV for the given headstage.

11.107.3.11 GetHeadstageDacBits() uint32_t GetHeadstageDacBits (uint32_t Headstage)

Gets the Number of DAC bits for a given headstage.

Parameters

Headstage	The headstage to query.

Returns

The number of bits the DAC has for the given headstage.

Gets the DAC Current Range in uA for a given headstage.

Headstage T	he headstage to query.
-------------	------------------------

Returns

The DAC Current Range in uA for the given headstage.

11.107.3.13 GetHeadstageDacCurrentResolutionInNanoAmpere() uint32_t GetHeadstageDacCurrent↔ ResolutionInNanoAmpere (

```
uint32_t Headstage )
```

Gets the DAC Current Resolution in nA for a given headstage.

Parameters

dstage The headstage to	to query.
-------------------------	-----------

Returns

The DAC Current Resolution in nA for the given headstage.

Gets the DAC Voltage Range in mV for a given headstage.

Parameters

Headstage	The headstage to query.

Returns

The DAC Voltage Range in mV for the given headstage.

Gets the DAC Voltage Resolution in uV for a given headstage.

Returns

The DAC Voltage Resolution in uV for the given headstage.

11.107.3.16 GetHeadstageGainInPermille() uint32_t GetHeadstageGainInPermille (uint32_t Headstage)

Gets the gain factor in permille for a given headstage.

Parameters

Returns

The gain factor in permille for the given headstage.

11.107.3.17 **GetHeadstageID()** SCU_HeadstageIdEnumNet GetHeadstageID (uint32_t Headstage)

Gets the headstage fpga ID.

Parameters

Headstage	The headstage to query.
	, ,

Returns

The headstage fpga ID.

Gets the number of analog channels for a given headstage.

Headstage	The headstage to query.
-----------	-------------------------

Returns

The number of analog channels the headstage has.

Gets the number of stimulation channels for a given headstage.

Parameters

dstage The headstage to	to query.
-------------------------	-----------

Returns

The number of stimulation channels the headstage has.

11.107.3.20 GetHeadstagePowerStateAtStart() bool GetHeadstagePowerStateAtStart (uint32_t Headstage)

Gets the Power Status at SCU Power on of a given headstage.

Parameters

Headstage	The headstage to query.

Returns

The Power State at startup for the given headstage: bool false \rightarrow off, bool true \rightarrow on.

```
11.107.3.21 GetHeadstageSamplerate() uint32_t GetHeadstageSamplerate ( uint32_t Headstage)
```

Gets the Samplerate of a given headstage.

Returns

The samplerate in Hz for the given headstage.

11.107.3.22 GetHeadstageSerialNumber() String ^ GetHeadstageSerialNumber (uint32_t Headstage)

Gets the serial number of a given headstage.

Parameters

Returns

The serial number of the headstage.

11.107.3.23 GetMaxNumberOfHeadstages() uint32_t GetMaxNumberOfHeadstages ()

Gets the maximal number of headstages.

Returns

The maximal number of headstages.

11.107.3.24 GetMaxStimulusChannelsPerHeadstage() uint32_t GetMaxStimulusChannelsPerHeadstage

Gets the maximal number of stimulation channels a headstage can have.

Returns

The maximal number of stimulation channels a headstage can have.

Gets the mode for the reference electrode

Da	ro	m	Δi	ŀΔ	rc
Εи	17		ы	ш	15

Headstage	The headstage number
-----------	----------------------

Returns

The mode

$\textbf{11.107.3.26} \quad \textbf{GetReferenceElectrodeSwitchState()} \quad \texttt{ReferenceElectrodeSwitchPositionEnumNet} \quad \texttt{Get} \leftarrow \textbf{Constant} \quad \texttt{Constant} \quad \texttt$

Gets the position of the switch for the reference electrode

Parameters

Headstage	The headstage number
-----------	----------------------

Returns

The switch position

11.107.3.27 HasAnalogOut() bool HasAnalogOut ()

Has AnalogOut hardware

Returns

Enabled

11.107.3.28 HasGalvanicIsolation() bool HasGalvanicIsolation ()

Has galvanic isolated hardware

Returns

Enabled

```
11.107.3.29 HasHSPowerSwitch() bool HasHSPowerSwitch ( )
Has SCU HS power switch
Returns
     Has Switch
11.107.3.30 IsAnalogOutEnabled() bool IsAnalogOutEnabled ( )
Is AnalogOut enabled
Returns
     Enabled
11.107.3.31 IsAutomaticAnalogOut() bool IsAutomaticAnalogOut ( )
Is Automatic source channel selection selected
Returns
     Automatic
11.107.3.32 | IsHeadstageAvailable() bool IsHeadstageAvailable (
             uint32_t Headstage )
Checks whether the given headstage is available.
Parameters
 Headstage
              The headstage to query.
Returns
     'true' if the headstage is connected.
```

```
Generated by Doxygen
```

Is the HS powered

11.107.3.33 IsHSPowered() bool IsHSPowered (uint32_t Headstage)

Returns

'true' if the headstage is powered.

11.107.3.34 | IsInDacqLegacyMode() | bool | IsInDacqLegacyMode ()

Is the SCU in legacy mode

Returns

Is Enabled

```
11.107.3.35 OnGetAvailableHeadstages() delegate void OnGetAvailableHeadstages ( uint32_t AvailableHeadstages )
```

11.107.3.36 OnlsHeadstageAvailable() delegate void OnlsHeadstageAvailable (uint32_t Headstage, bool available)

Power the HS

Parameters

Headstage	The headstage to query.
power	'true' if the headstage is powered.

${\bf 11.107.3.38} \quad {\bf SetAnalogOutADCRange()} \quad {\tt void SetAnalogOutADCRange ()} \quad {\tt$

uint32_t range)

Sets the analog out ADC range

range I	Range
-----------	-------

Set the source channel number for a certain output channel

Parameters

out_channel	Output channel number
source_channel	Source channel number

11.107.3.40 SetAnalogOutDACRange() void SetAnalogOutDACRange (AnalogOut_DAC_Range_EnumNet range)

Sets the analog out DAC range

Parameters

range	Range
-------	-------

11.107.3.41 SetDacqLegacyMode() void SetDacqLegacyMode (bool enable)

Enable the SCU legacy mode

Parameters

```
enable Enable
```

11.107.3.42 SetHeadstagePowerStateAtStart() void SetHeadstagePowerStateAtStart (uint32_t Headstage, bool Powerstatus)

Sets the Power Status at SCU Power on of a given headstage.

Headstage	The headstage number
Powerstatus	The Power State at startup for the given headstage: bool false -> off, bool true -> on.

11.107.3.43 SetReferenceElectrodeMode() void SetReferenceElectrodeMode (uint32_t Headstage, ReferenceElectrodeModeEnumNet NewValue)

Sets the mode for the reference electrode

Parameters

Headstage	The headstage number
NewValue	The mode

Sets the position of the switch for the reference electrode

Parameters

Headstage	The headstage number
NewSwitchPos	The switch position

11.107.4 Event Documentation

```
11.107.4.1 GetAvailableHeadstagesEvent OnGetAvailableHeadstages^ GetAvailableHeadstagesEvent [add], [remove], [raise]
```


Event fires when the bitmap of available headstages has changed

```
11.107.4.2 IsHeadstageAvailableEvent OnIsHeadstageAvailable^ IsHeadstageAvailableEvent [add], [remove], [raise]
```

Event fires when 'true' if the headstage is connected for the headstage to guery has changed

11.108 CSerialPortNet Class Reference

Inheritance diagram for CSerialPortNet:



Public Member Functions

- CSerialPortNet (void)
- void Send (array< byte >^ buffer)
- void Send (String[^] command)
- array< byte > ^ Receive (void)
- array< byte > ^ Receive (int length)
- String ^ ReceiveString (void)
- String ^ ReceiveString (int length)
- int GetBytesAvailable (void)

Additional Inherited Members

11.108.1 Constructor & Destructor Documentation

11.108.2 Member Function Documentation

```
11.108.2.1 GetBytesAvailable() int GetBytesAvailable ( void )
```


11.108.2.2 Receive() [1/2] array
byte>
$$^{\land}$$
 Receive (int length)

11.108.2.3 Receive() [2/2] array
byte>
$$^{\wedge}$$
 Receive (void)

11.109 CStg200xBasicNet Class Reference

Base class for the Stg200x.

Inheritance diagram for CStg200xBasicNet:



Public Member Functions

virtual ~CStg200xBasicNet ()

The destructor.

void SetOutputRate (uint32_t rate)

Change the output rate of the STG. Valid rates are from 1000 Hz to 50000 Hz.

• uint32_t GetOutputRate ()

Queries the output rate of the STG. Valid rates are from 1000 Hz to 50000 Hz.

void SendStart (uint32_t triggermap)

Start (Trigger) the STG. The startup delay is in the range of a few ms.

void SendStop (uint32_t triggermap)

Stop some or all triggers of the STG.

void SendStop (uint32_t triggermap, int options)

Stop some or all triggers of the STG.

void GetStgVersionInfo ([Out]String^% SwVersion, [Out]String^% HwVersion)

Queries software and hardware version.

• void GetAnalogRanges (int channel, [Out]int% URange, [Out]int% IRange)

Gets the range of the analog outputs.

void GetAnalogResolution (int channel, [Out]int% URes, [Out]int% IRes)

Gets the resolution of the analog outputs.

virtual int32 t GetDACResolution ()

Gets number of bits of the DAC resolution.

virtual int32 t GetVoltageRangeInMicroVolt (uint32 t channel)

Gets the Voltage Range of the specified channel in Microvolts.

virtual int32 t GetVoltageResolutionInMicroVolt (uint32 t channel)

Gets the Voltage Resolution of the specified channel in Microvolts.

virtual int32 t GetCurrentRangeInNanoAmp (uint32 t channel)

Gets the Current Range of the specified channel in Nanoamps.

virtual int32 t GetCurrentResolutionInNanoAmp (uint32 t channel)

Gets the Current Resolution of the specified channel in Nanoamps.

Queries Download information from the STG. If download information was stored by the use of SetStgProgramInfo, this function can be used to retrieve it.

void GetStgProgramInfo ([Out]bool% IsProgrammed, [Out]DateTime% timestamp, [Out]String^% filename, [Out]Guid% guid)

Queries Download information from the STG. If download information was stored by the use of SetStgProgramInfo, this function can be used to retrieve it.

void SetStgProgramInfo (DateTime timestamp, String[^] filename, Guid guid)

Store Download information in the STG. This function can be used to store the filename and timestamp of the last download for later query.

• uint32 t GetAvailableMemory ()

Gets the amount of memory available in the currently selected segment of the STG.

uint32_t GetTotalMemory ()

Gets the total amount of memory available on the STG (all segments).

virtual uint32_t GetNumberOfAnalogChannels ()

Gets the Number of available analog channels of the device.

virtual uint32_t GetNumberOfSyncoutChannels ()

Gets the Number of available syncout channels of the device.

virtual uint32_t GetNumberOfTriggerInputs ()

Gets the Number of trigger inputs of the device.

virtual uint32_t GetNumberOfHWDACPaths ()

Gets the Number of HW Stimulation DACs of the device.

virtual uint32 t GetNumberOfStimulationSourcesPerElectrode ()

Gets the number of stimulation sources (DACs) per electrode.

virtual void SetVoltageMode (unsigned int channel)

Sets a channel to voltage mode (STG3008-FA and STG400x only).

virtual void SetCurrentMode (unsigned int channel)

Sets a channel to current mode (STG3008-FA and STG400x only).

• virtual void SetVoltageMode ()

Sets all channels to voltage mode (STG3008-FA and STG400x only).

virtual void SetCurrentMode ()

Sets all channels to current mode (STG3008-FA and STG400x only).

virtual void SetMeasurementMode (unsigned int channel)

Sets a channel to measurement mode (STG3008-FA)

- virtual void SetFAAmplification (unsigned int amplification)
- virtual uint32 t GetFAAmplification ()
- · virtual void SetAutocalibrationDisabled (unsigned int channel, bool disable)

Sets the autocalibration configuration.

virtual bool GetAutocalibrationDisabled (unsigned int channel)

Gets the autocalibration configuration.

virtual void SetElectrodeMode (uint32 t electrode, array< ElectrodeModeEnumNet >^ mode)

Puts an electrode in either automatic or manual mode.

virtual void SetElectrodeMode (uint32 t electrode, ElectrodeModeEnumNet mode)

Puts an electrode in either automatic or manual mode.

virtual void SetElectrodeMode (uint32_t Scu_HS, uint32_t electrode, array< ElectrodeModeEnumNet >^
mode)

Puts an electrode in either automatic or manual mode.

virtual void SetElectrodeMode (uint32 t Scu HS, uint32 t electrode, ElectrodeModeEnumNet mode)

Puts an electrode in either automatic or manual mode.

virtual uint32 t GetElectrodeMode (uint32 t electrode)

Gets the mode an electrode is in.

virtual uint32_t GetElectrodeMode (uint32_t Scu_HS, uint32_t electrode)

Gets the mode an electrode is in.

Defines the DAC to use for an electrode.

virtual void SetElectrodeDacMux (uint32_t electrode, uint32_t listmodeIndex, ElectrodeDacMuxEnumNet dacMux)

Defines the DAC to use for an electrode.

virtual void SetElectrodeDacMux (uint32_t Scu_HS, uint32_t electrode, uint32_t listmodeIndex, Electrode
 —
 DacMuxEnumNet dacMux)

Defines the DAC to use for an electrode.

virtual void SetElectrodeDacMux (uint32_t Scu_HS, uint32_t electrode, uint32_t listmodeIndex, array
 ElectrodeDacMuxEnumNet >^ dacMux)

Defines the DAC to use for an electrode.

virtual ElectrodeDacMuxEnumNet GetElectrodeDacMux (uint32_t electrode, uint32_t listmodeIndex)

Gets the DAC which is used for an electrode.

virtual ElectrodeDacMuxEnumNet GetElectrodeDacMux (uint32_t Scu_HS, uint32_t electrode, uint32_
 t listmodeIndex)

Gets the DAC which is used for an electrode.

virtual void SetElectrodeEnable (uint32 t electrode, uint32 t listmodeIndex, array < bool >^ enable)

Enables or disables the stimulation switch for an electrode.

virtual void SetElectrodeEnable (uint32_t electrode, uint32_t listmodeIndex, bool enable)

Enables or disables the stimulation switch for an electrode.

virtual void SetElectrodeEnable (uint32_t Scu_HS, uint32_t electrode, uint32_t listmodeIndex, bool enable)

Enables or disables the stimulation switch for an electrode.

virtual void SetElectrodeEnable (uint32_t Scu_HS, uint32_t electrode, uint32_t listmodeIndex, array< bool
 <p>^ enable)

Enables or disables the stimulation switch for an electrode.

• virtual bool GetElectrodeEnable (uint32 t electrode, uint32 t listmodeIndex)

Gets weather an electrode is enabled or disabled for stimulation.

virtual bool GetElectrodeEnable (uint32_t Scu_HS, uint32_t electrode, uint32_t listmodeIndex)

Gets weather an electrode is enabled or disabled for stimulation.

virtual void SetExternalElectrodeEnable (uint32_t electrode, uint32_t listmodeIndex, array< bool >^ enable)

Enables or disables the stimulation switch for an external electrode.

virtual void SetExternalElectrodeEnable (uint32_t electrode, uint32_t listmodeIndex, bool enable)

Enables or disables the stimulation switch for an external electrode.

virtual bool GetExternalElectrodeEnable (uint32 t electrode, uint32 t listmodeIndex)

Gets weather an electrode is enabled or disabled for stimulation.

virtual void SetBlankingEnable (uint32 t electrode, bool enable)

Defines whether an electrode should be blanked while stimulation is in progress.

virtual void SetBlankingEnable (uint32_t electrode, array< bool >^ enable)

Defines whether an electrode should be blanked while stimulation is in progress.

• virtual void SetBlankingEnable (uint32_t Scu_HS, uint32_t electrode, bool enable)

Defines whether an electrode should be blanked while stimulation is in progress.

virtual void SetBlankingEnable (uint32_t Scu_HS, uint32_t electrode, array< bool >^ enable)

Defines whether an electrode should be blanked while stimulation is in progress.

virtual bool GetBlankingEnable (uint32 t electrode)

Gets whether an electrode should be blanked while stimulation is in progress.

virtual bool GetBlankingEnable (uint32 t Scu HS, uint32 t electrode)

Gets whether an electrode should be blanked while stimulation is in progress.

• virtual void SetEnableAmplifierProtectionSwitch (uint32 t electrode, bool enable)

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

virtual void SetEnableAmplifierProtectionSwitch (uint32_t electrode, array< bool >^ enable)

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

• virtual void SetEnableAmplifierProtectionSwitch (uint32_t Scu_HS, uint32_t electrode, bool enable)

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

 virtual void SetEnableAmplifierProtectionSwitch (uint32_t Scu_HS, uint32_t electrode, array< bool >^ enable)

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

virtual bool GetEnableAmplifierProtectionSwitch (uint32_t electrode)

Gets whether the Amplifier Protection Switch is openend while stimulation is in progress.

virtual bool GetEnableAmplifierProtectionSwitch (uint32 t Scu HS, uint32 t electrode)

Gets whether the Amplifier Protection Switch is openend while stimulation is in progress.

- virtual uint32 t GetNumberOfStimulationElectrodes ()
- template<typename digitalsourceenum >

virtual void SetTriggerSource (unsigned int triggernum, DigitalSource< digitalsourceenum >^ triggersource, int bitnum_offset)

- virtual void SetTriggerSource (unsigned int triggernum, TriggerSourceEnumNet triggersource, int bitnum_

 offset)
- virtual void SetTriggerSource (unsigned int triggernum, TriggerSourceEnumNet triggersource)
- virtual TriggerSourceEnumNet GetTriggerSource (unsigned int triggernum)
- virtual void SetListmodeIndexRange (unsigned int sideband, unsigned int startIndex, unsigned int endIndex, unsigned int mode)
- virtual void GetListmodeIndexRange (unsigned int sideband, unsigned int &startIndex, unsigned int &end
 —
 Index, unsigned int &mode)
- virtual void SetListmodeTriggerSource (unsigned int sideband, TriggerSourceEnumNet triggersource, int bitnumOffset)
- virtual void SetListmodeTriggerSource (unsigned int sideband, TriggerSourceEnumNet triggersource)
- · virtual TriggerSourceEnumNet GetListmodeTriggerSource (unsigned int sideband)
- virtual void ListModeSendStart (unsigned int sidebandMask)
- virtual void ListModeSendStop (unsigned int sidebandMask)
- virtual void SetHeadstage (unsigned int headstage)
- virtual uint32 t GetHeadstage ()
- virtual void SetDacAmplificationFactor (uint32_t DacNumber, double Factor)

Set the amplification factor for a DAC.

virtual double GetDacAmplificationFactor (uint32_t DacNumber)

Get the amplification factor for a DAC.

• virtual void SetDigoutMode (Stg200xDigoutModeEnumNet digoutMode)

Sets the operation mode of the digital outport port, can be Monitor, Manual or SyncOut

virtual Stg200xDigoutModeEnumNet GetDigoutMode()

Gets the operation mode of the digital outport port, can be Monitor, Manual or SyncOut

virtual void SetDigoutValue (uint32_t digoutValue)

Sets the Value on the digital output port when in manual mode.

virtual uint32_t GetDigoutValue ()

Gets the Value on the digital output port.

virtual uint32_t GetDiginValue ()

Gets the Value on the digital input port.

virtual void SetSyncoutMap (uint32 t channel, uint32 t syncoutMap)

Sets the mapping between external syncout outputs and internal syncout channels.

virtual uint32_t GetSyncoutMap (uint32_t channel)

Gets the mapping between external syncout outputs and internal syncout channels.

Additional Inherited Members

11.109.1 Detailed Description

Base class for the Stg200x.

From this class all STG related classes are derived: Mcs.Usb.CStg200xDownloadBasicNet Mcs.Usb.CStg200xDownloadNet for Download Mode and Mcs.Usb.CStg200xStreamingNet for Streaming Mode.

CStg200xBasicNet is the base class to control MCS STG device.

11.109.2 Constructor & Destructor Documentation

```
11.109.2.1 ~CStg200xBasicNet() virtual ~CStg200xBasicNet () [virtual]
```

The destructor.

11.109.3 Member Function Documentation

Gets the range of the analog outputs.

channel	The channel which is queried.
URange	The Voltage range in mV.
IRange	The Current range in uA.


```
[Out] int% URes,
[Out] int% IRes)
```

Gets the resolution of the analog outputs.

Parameters

channel	The channel which is queried.
URes	The Voltage resolution in mV.
IRes	The Current resolution in uA.

11.109.3.3 GetAutocalibrationDisabled() virtual bool GetAutocalibrationDisabled (unsigned int *channel*) [virtual]

Gets the autocalibration configuration.

Parameters

channel	The channel number.

Returns

true if autocalibration is disabled.

11.109.3.4 GetAvailableMemory() uint32_t GetAvailableMemory ()

Gets the amount of memory available in the currently selected segment of the STG.

Returns

The memory available in the currently selected segment in bytes.

```
11.109.3.5 GetBlankingEnable() [1/2] virtual bool GetBlankingEnable ( uint32_t electrode ) [virtual]
```

Gets whether an electrode should be blanked while stimulation is in progress.

electrode The electrode number.

Returns

true if blanking is enabled while stimulation is in progress.

Gets whether an electrode should be blanked while stimulation is in progress.

Parameters

Parameters

electrode	The electrode number.
CICCLI DUC	i ilie electione mumbei.

Returns

true if blanking is enabled while stimulation is in progress.

```
11.109.3.7 GetCurrentRangeInNanoAmp() virtual int32_t GetCurrentRangeInNanoAmp ( uint32_t channel) [virtual]
```

Gets the Current Range of the specified channel in Nanoamps.

Parameters

channel	Channel which is queried.
oriariror	Onamor Willom to quonou.

Returns

The Current Range of the specified channel in Nanoamps.

11.109.3.8 GetCurrentResolutionInNanoAmp() virtual int32_t GetCurrentResolutionInNanoAmp (uint32_t channel) [virtual]

Gets the Current Resolution of the specified channel in Nanoamps.

Parameters

channel	Channel which is queried.
---------	---------------------------

Returns

The Current Resolution of the specified channel in Nanoamps.

11.109.3.9 **GetDacAmplificationFactor()** virtual double GetDacAmplificationFactor (uint32_t DacNumber) [virtual]

Get the amplification factor for a DAC.

Parameters

DacNumber The number of the DAC.

Returns

the amplification factor for the DAC queried, range is from -1.99999 to +1.99999.

```
11.109.3.10 GetDACResolution() virtual int32_t GetDACResolution ( ) [virtual]
```

Gets number of bits of the DAC resolution.

Returns

The DAC resolution in bits.

```
11.109.3.11 GetDiginValue() virtual uint32_t GetDiginValue ( ) [virtual]
```

Gets the Value on the digital input port.

Returns

The current value on the digital inputs.

11.109.3.12 GetDigoutMode() virtual Stg200xDigoutModeEnumNet GetDigoutMode () [virtual]

Gets the operation mode of the digital outport port, can be Monitor, Manual or SyncOut

Returns

The current operation mode.

11.109.3.13 GetDigoutValue() virtual uint32_t GetDigoutValue () [virtual]

Gets the Value on the digital output port.

Returns

The current value on the digital outputs.

Gets the DAC which is used for an electrode.

Parameters

electrode	The electrode number.
listmodeIndex	The index for listmode.

Returns

The DAC in use, can be 1, 2 or 3. If the electrode is grounded 0 is returned.

Gets the DAC which is used for an electrode.

Parameters

Scu HS	The SCU headstage number.
004_770	The coo headstage nambor.

electrode	The electrode number.
listmodeIndex	The index for listmode.

Returns

The DAC in use, can be 1, 2 or 3. If the electrode is grounded 0 is returned.

```
11.109.3.16 GetElectrodeEnable() [1/2] virtual bool GetElectrodeEnable ( uint32_t electrode, uint32_t listmodeIndex ) [virtual]
```

Gets weather an electrode is enabled or disabled for stimulation.

Parameters

electrode	The electrode number.
listmodeIndex	The index for listmode.

Returns

true if the electrode is enabled, false if it is disabled.

Gets weather an electrode is enabled or disabled for stimulation.

Parameters

Scu_HS	The SCU headstage number.

Parameters

electrode	The electrode number.
listmodeIndex	The index for listmode.

Returns

true if the electrode is enabled, false if it is disabled.

```
11.109.3.18 GetElectrodeMode() [1/2] virtual uint32_t GetElectrodeMode ( uint32_t electrode ) [virtual]
```

Gets the mode an electrode is in.

Parameters

electrode	The electrode number.
CICCLIOUC	THE ELECTIONS HUTTING.

Returns

0 for automatic and 3 for manual mode.

```
11.109.3.19 GetElectrodeMode() [2/2] virtual uint32_t GetElectrodeMode ( uint32_t Scu_HS, uint32_t electrode ) [virtual]
```

Gets the mode an electrode is in.

Parameters

Scu_HS	The SCU headstage number.
--------	---------------------------

Parameters

electrode	The electrode number.
-----------	-----------------------

Returns

0 for automatic and 3 for manual mode.

Gets whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

electrode	The electrode number.
CICCIIOUC	THE ELECTIONE HUITIDEL.

Returns

true if the switch is to be opened, false if it is closed while stimulation is in progress.

Gets whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

Scu_HS	The SCU headstage number.
--------	---------------------------

Parameters

electrode	The electrode number.
-----------	-----------------------

Returns

true if the switch is to be opened, false if it is closed while stimulation is in progress.

Gets weather an electrode is enabled or disabled for stimulation.

Parameters

electrode	The electrode number.
listmodeIndex	The index for listmode.

Returns

true if the electrode is enabled, false if it is disabled.

```
11.109.3.23 GetFAAmplification() virtual uint32_t GetFAAmplification ( ) [virtual]
```

```
11.109.3.24 GetHeadstage() virtual uint32_t GetHeadstage ( ) [virtual]
```

```
\textbf{11.109.3.25} \quad \textbf{GetListmodeIndexRange()} \quad \texttt{virtual void GetListmodeIndexRange ()}
```

```
unsigned int sideband,
unsigned int & startIndex,
unsigned int & endIndex,
unsigned int & mode ) [virtual]
```

```
11.109.3.26 GetListmodeTriggerSource() virtual TriggerSourceEnumNet GetListmodeTriggerSource ( unsigned int sideband ) [virtual]
```

11.109.3.27 GetNumberOfAnalogChannels() virtual uint32_t GetNumberOfAnalogChannels () [virtual]

Gets the Number of available analog channels of the device.

Returns

The number of analog channels.

11.109.3.28 GetNumberOfHWDACPaths() virtual uint32_t GetNumberOfHWDACPaths () [virtual]

Gets the Number of HW Stimulation DACs of the device.

Returns

The number of independent HW Stimulation outputs.

```
11.109.3.29 GetNumberOfStimulationElectrodes() virtual uint32_t GetNumberOfStimulationElectrodes ( ) [virtual]
```

```
11.109.3.30 GetNumberOfStimulationSourcesPerElectrode() virtual uint32_t GetNumberOfStimulation← SourcesPerElectrode ( ) [virtual]
```

Gets the number of stimulation sources (DACs) per electrode.

Returns

The number of stimulation sources (DACs) per electrode.

```
11.109.3.31 GetNumberOfSyncoutChannels() virtual uint32_t GetNumberOfSyncoutChannels ( ) [virtual]
```

Gets the Number of available syncout channels of the device.

Returns

The number of analog channels.

11.109.3.32 GetNumberOfTriggerInputs() virtual uint32_t GetNumberOfTriggerInputs () [virtual] Gets the Number of trigger inputs of the device.

Returns

The number of trigger inputs.

```
11.109.3.33 GetOutputRate() uint32_t GetOutputRate ( )
```

Queries the output rate of the STG. Valid rates are from 1000 Hz to 50000 Hz.

Returns

Returns the current output rate in Hz.

11.109.3.34 GetStgProgramInfo() [1/2] void GetStgProgramInfo (

```
[Out] bool% IsProgrammed,
[Out] DateTime% timestamp,
[Out] String^% filename,
[Out] Guid% guid)
```

Queries Download information from the STG. If download information was stored by the use of SetStgProgramInfo, this function can be used to retrieve it.

IsProgrammed	Flag wether download information is valid.
timestamp	The timestamp of last download.
filename	The filename of the downlaoded waveform.
guid	A GUID.

Queries Download information from the STG. If download information was stored by the use of SetStgProgramInfo, this function can be used to retrieve it.

Parameters

IsProgrammed	Flag wether download information is valid.
timestamp	The timestamp of last download.
filename	The filename of the downlaoded waveform.

Queries software and hardware version.

Parameters

SwVersion	The current Software Version of the STG.
HwVersion	The Hardware Revision of the STG.

Gets the mapping between external syncout outputs and internal syncout channels.

Parameters

channel	The external syncout output channel number (zero based).

Returns

The bitmap of internal syncout channels mapped to channel.

11.109.3.38 GetTotalMemory() uint32_t GetTotalMemory ()

Gets the total amount of memory available on the STG (all segments).

Returns

The total memory available on the STG in bytes.

```
11.109.3.39 GetTriggerSource() virtual TriggerSourceEnumNet GetTriggerSource (
unsigned int triggernum) [virtual]
```

```
11.109.3.40 GetVoltageRangeInMicroVolt() virtual int32_t GetVoltageRangeInMicroVolt ( uint32_t channel) [virtual]
```

Gets the Voltage Range of the specified channel in Microvolts.

Parameters

channel	Channel which is queried.

Returns

The Voltage Range of the specified channel in Microvolts.

```
11.109.3.41 GetVoltageResolutionInMicroVolt() virtual int32_t GetVoltageResolutionInMicroVolt ( uint32_t channel ) [virtual]
```

Gets the Voltage Resolution of the specified channel in Microvolts.

Parameters

-		
ſ	channel	Channel which is queried.

Returns

The Voltage Resolution of the specified channel in Microvolts.

```
11.109.3.42 ListModeSendStart() virtual void ListModeSendStart ( unsigned int sidebandMask ) [virtual]
```

```
11.109.3.43 ListModeSendStop() virtual void ListModeSendStop ( unsigned int sidebandMask ) [virtual]
```

```
11.109.3.44 SendStart() void SendStart ( uint32_t triggermap )
```

Start (Trigger) the STG. The startup delay is in the range of a few ms.

Parameters

triggermap	A bitmap of triggers which will be started.
------------	---

```
11.109.3.45 SendStop() [1/2] void SendStop ( uint32_t triggermap )
```

Stop some or all triggers of the STG.

Parameters

trigge	rmap	A bitmap of triggers which will be stopped.
--------	------	---

Stop some or all triggers of the STG.

triggermap	A bitmap of triggers which will be stopped.
options	bitmap of options, currently only STOP_OPTION_SAVESTOP (0x80) is defined, which bypasses the stop commands when a syncout assossiated with a given sync-out has bit 1 (0x02) set. Can be used e.g. to prevent a stop while a biphasic stimulation pulse is active

11.109.3.47 SetAutocalibrationDisabled() virtual void SetAutocalibrationDisabled (unsigned int channel, bool disable) [virtual]

Sets the autocalibration configuration.

Parameters

channel	The channel number.
disable	true if autocalibration is to be disabled.

```
11.109.3.48 SetBlankingEnable() [1/4] virtual void SetBlankingEnable ( uint32_t electrode, array< bool >^ enable ) [virtual]
```

Defines whether an electrode should be blanked while stimulation is in progress.

Parameters

electrode	The electrode number.
enable	True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.

```
11.109.3.49 SetBlankingEnable() [2/4] virtual void SetBlankingEnable ( uint32_t electrode, bool enable ) [virtual]
```

Defines whether an electrode should be blanked while stimulation is in progress.

Parameters

electrode	The electrode number.
enable	True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.

Defines whether an electrode should be blanked while stimulation is in progress.

Scu_HS	The SCU headstage number.
--------	---------------------------

electrode	The electrode number.
enable	True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.

Defines whether an electrode should be blanked while stimulation is in progress.

Parameters

Parameters

electrode	The electrode number.
enable	True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.

```
11.109.3.52 SetCurrentMode() [1/2] virtual void SetCurrentMode ( ) [virtual]
```

Sets all channels to current mode (STG3008-FA and STG400x only).

```
11.109.3.53 SetCurrentMode() [2/2] virtual void SetCurrentMode (
unsigned int channel) [virtual]
```

Sets a channel to current mode (STG3008-FA and STG400x only).

channel	The channel to change.
Charine	The charmer to change.

Set the amplification factor for a DAC.

Parameters

DacNumber	The number of the DAC.
Factor	the amplification factor for that DAC, range is from -1.99999 to +1.99999.

11.109.3.55 SetDigoutMode() virtual void SetDigoutMode (Stg200xDigoutModeEnumNet digoutMode) [virtual]

Sets the operation mode of the digital outport port, can be Monitor, Manual or SyncOut

Parameters

digoutMode	The new operation mode.
------------	-------------------------

```
11.109.3.56 SetDigoutValue() virtual void SetDigoutValue ( uint32_t digoutValue ) [virtual]
```

Sets the Value on the digital output port when in manual mode.

Parameters

digoutValue	The new value on the digital outputs.
-------------	---------------------------------------

Defines the DAC to use for an electrode.

listmodeIndex	The index for listmode.	
dacMux	The DAC to use, can be ElectrodeDacMuxEnumNet.Stg1 (1), ElectrodeDacMuxEnumNet.Stg2	
	(2) or ElectrodeDacMuxEnumNet.Stg3 (3). To ground an electrode, use	
	(2) or ElectrodeDacMuxEnumNet.Stg3 (3). To ground an electrode, use ElectrodeDacMuxEnumNet.Ground (0).	

Defines the DAC to use for an electrode.

Parameters

Parameters

listmodeIndex	The index for listmode.	
dacMux	The DAC to use, can be ElectrodeDacMuxEnumNet.Stg1 (1), ElectrodeDacMuxEnumNet.Stg2 (2) or ElectrodeDacMuxEnumNet.Stg3 (3). To ground an electrode, use	
	ElectrodeDacMuxEnumNet.Ground (0).	

Defines the DAC to use for an electrode.

Scu_HS	The SCU headstage number.

electrode The electrode number.	r.
---------------------------------	----

Parameters

listmodeIndex	The index for listmode.	
dacMux	The DAC to use, can be ElectrodeDacMuxEnumNet.Stg1 (1), ElectrodeDacMuxEnumNet.Stg2	
	(2) or ElectrodeDacMuxEnumNet.Stg3 (3). To ground an electrode, use	
	ElectrodeDacMuxEnumNet.Ground (0).	

$\textbf{11.109.3.60} \quad \textbf{SetElectrodeDacMux() [4/4]} \quad \text{virtual void SetElectrodeDacMux ()}$

```
uint32_t Scu_HS,
uint32_t electrode,
uint32_t listmodeIndex,
ElectrodeDacMuxEnumNet dacMux ) [virtual]
```

Defines the DAC to use for an electrode.

Parameters

Scu_HS The SCU headstage nu

Parameters

electrode	The electrode number.
	i ilie electione iluilibei.

Parameters

listmodeIndex	The index for listmode.	
dacMux	The DAC to use, can be ElectrodeDacMuxEnumNet.Stg1 (1), ElectrodeDacMuxEnumNet.Stg2	
	(2) or ElectrodeDacMuxEnumNet.Stg3 (3). To ground an electrode, use	
	ElectrodeDacMuxEnumNet.Ground (0).	

11.109.3.61 SetElectrodeEnable() [1/4] virtual void SetElectrodeEnable (uint32_t electrode,

```
uint32_t listmodeIndex, array < bool >^  enable ) [virtual]
```

Enables or disables the stimulation switch for an electrode.

Parameters

Parameters

listmodeIndex	The index for listmode.
enable	1 to enable the electrode, 0 to disable.

Enables or disables the stimulation switch for an electrode.

Parameters

electrode	The electrode number.
-----------	-----------------------

Parameters

listmodeIndex	The index for listmode.
enable	1 to enable the electrode, 0 to disable.

Enables or disables the stimulation switch for an electrode.

Scu_HS	The SCU headstage number.
--------	---------------------------

Parameters

Parameters

listmodeIndex	The index for listmode.
enable	1 to enable the electrode, 0 to disable.

11.109.3.64 SetElectrodeEnable() [4/4] virtual void SetElectrodeEnable (

```
uint32_t Scu_HS,
uint32_t electrode,
uint32_t listmodeIndex,
bool enable ) [virtual]
```

Enables or disables the stimulation switch for an electrode.

Parameters

Parameters

listmodeIndex	The index for listmode.
enable	1 to enable the electrode, 0 to disable.

```
11.109.3.65 SetElectrodeMode() [1/4] virtual void SetElectrodeMode ( uint32_t electrode, array< ElectrodeModeEnumNet >^ mode ) [virtual]
```

Puts an electrode in either automatic or manual mode.

Parameters

electrode	The electrode number.
-----------	-----------------------

Returns

0 for automatic and 3 for manual mode.

Puts an electrode in either automatic or manual mode.

Parameters

electrode	The electrode number.
electione	I THE EIECHOUE HUITIDEL.

Parameters

mode	0 for automatic and 3 for manual mode.
moue	i o ioi automatic and 3 ioi mandai mode.

Puts an electrode in either automatic or manual mode.

Scu_HS The SCU headstage number.	Scu_HS	The SCU headstage number.
------------------------------------	--------	---------------------------

Parameters

electrode The electrode number.	
---------------------------------	--

Returns

0 for automatic and 3 for manual mode.

11.109.3.68 SetElectrodeMode() [4/4] virtual void SetElectrodeMode (

```
uint32_t Scu_HS,
uint32_t electrode,
ElectrodeModeEnumNet mode ) [virtual]
```

Puts an electrode in either automatic or manual mode.

Parameters

Scu_HS The SCU headstage number.

Parameters

electrode	The electrode number.
-----------	-----------------------

Parameters

mode 0 for automatic and 3 for manual mode.

11.109.3.69 SetEnableAmplifierProtectionSwitch() [1/4] virtual void SetEnableAmplifierProtection← Switch (

```
uint32_t electrode,
array< bool >^{\wedge} enable ) [virtual]
```

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

electrode	de The electrode number.	
enable True if the switch is to be opened, false if it is to remain closed while stimulation is in prog		

```
11.109.3.70 SetEnableAmplifierProtectionSwitch() [2/4] virtual void SetEnableAmplifierProtection←
Switch (
uint32_t electrode,
bool enable) [virtual]
```

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

electrod	The electrode number.	
enable True if the switch is to be opened, false if it is to remain closed while stimulation is in progres		

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

Scu_HS	The SCU headstage number.
--------	---------------------------

electrode	The electrode number.	
enable	True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.	

```
11.109.3.72 SetEnableAmplifierProtectionSwitch() [4/4] virtual void SetEnableAmplifierProtection← Switch (
```

```
uint32_t Scu_HS,
uint32_t electrode,
bool enable ) [virtual]
```

Defines whether the Amplifier Protection Switch is openend while stimulation is in progress.

Parameters

Parameters

	electrode	The electrode number.	
enable True if the switch is to be opened, false if it is to remain closed while stimulation is in prog		True if the switch is to be opened, false if it is to remain closed while stimulation is in progress.	

Enables or disables the stimulation switch for an external electrode.

Parameters

electrode The elect	trode number.
---------------------	---------------

Parameters

listmodeIndex	The index for listmode.
enable	1 to enable the electrode, 0 to disable.

Enables or disables the stimulation switch for an external electrode.

electrode	The electrode number.
-----------	-----------------------

Parameters

listmodeIndex	The index for listmode.	
enable	1 to enable the electrode, 0 to disable.	

```
11.109.3.75 SetFAAmplification() virtual void SetFAAmplification (
             unsigned int amplification ) [virtual]
11.109.3.76 SetHeadstage() virtual void SetHeadstage (
             unsigned int headstage ) [virtual]
11.109.3.77 SetListmodeIndexRange() virtual void SetListmodeIndexRange (
             unsigned int sideband,
             unsigned int startIndex,
             unsigned int endIndex,
             unsigned int mode ) [virtual]
11.109.3.78 SetListmodeTriggerSource() [1/2] virtual void SetListmodeTriggerSource (
             unsigned int sideband,
             {\tt TriggerSourceEnumNet}\ triggersource\ )\quad [{\tt virtual}]
11.109.3.79 SetListmodeTriggerSource() [2/2] virtual void SetListmodeTriggerSource (
             unsigned int sideband,
             TriggerSourceEnumNet triggersource,
             int bitnumOffset ) [virtual]
```

Sets a channel to measurement mode (STG3008-FA).

11.109.3.80 SetMeasurementMode() virtual void SetMeasurementMode (

unsigned int channel) [virtual]

Change the output rate of the STG. Valid rates are from 1000 Hz to 50000 Hz.

Parameters

rate The new output rate in F	łz.
-------------------------------	-----

Store Download information in the STG. This function can be used to store the filename and timestamp of the last download for later query.

Parameters

timestamp	The timestamp of last download.
filename	The filename of the downlaoded waveform.

Sets the mapping between external syncout outputs and internal syncout channels.

channel	The external syncout output channel number (zero based).
syncoutMap	A bitmap of internal syncout channels to map to channel.

```
11.109.3.84 SetTriggerSource() [1/3] virtual void SetTriggerSource ( unsigned int triggernum,
```

```
\label{eq:digitalsource} \mbox{DigitalSource} < \mbox{digitalsourceenum} >^{\wedge} triggersource, \\ \mbox{int } \mbox{bitnum\_offset} \mbox{)} \mbox{ [virtual]}
```

```
11.109.3.85 SetTriggerSource() [2/3] virtual void SetTriggerSource ( unsigned int triggernum,
```

 ${\tt TriggerSourceEnumNet}\ triggersource\)\quad [{\tt virtual}]$

```
11.109.3.86 SetTriggerSource() [3/3] virtual void SetTriggerSource (
```

```
unsigned int triggernum,
TriggerSourceEnumNet triggersource,
int bitnum_offset ) [virtual]
```

```
11.109.3.87 SetVoltageMode() [1/2] virtual void SetVoltageMode ( ) [virtual]
```

Sets all channels to voltage mode (STG3008-FA and STG400x only).

```
11.109.3.88 SetVoltageMode() [2/2] virtual void SetVoltageMode (
unsigned int channel) [virtual]
```

Sets a channel to voltage mode (STG3008-FA and STG400x only).


Parameters

channel	The channel to change.

11.110 CStg200xDownloadBasicNet Class Reference

CStg200xDownloadBasicNet is the base class to control the download mode of the MCS STG device.

Inheritance diagram for CStg200xDownloadBasicNet:



Public Member Functions

virtual void SetupTrigger (uint32_t first_trigger, array< uint32_t >^ channelmap, array< uint32_t >^ syncoutmap, array< uint32_t >^ repeat)

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

- virtual void SetupTriggerSingle (uint32_t trigger, uint32_t channelmap, uint32_t syncoutmap, uint32_t repeat)

 Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.
- void GetTrigger ([Out] array< uint32_t >^% channelmap, [Out] array< uint32_t >^% syncoutmap, [Out] array< uint32_t >^% repeat)

Queries the trigger settings for the STG. Note that all memory segments have their own trigger setting.

void GetSweepCount ([Out] array< uint32 t >^% sweeps, [Out] array< uint32 t >^% triggers)

Get the sweep and trigger count of the STG.

- The triggercount tells how many times each trigger was active and is reset to zero on download of new channel data.
- The sweepcount tells how many times each trigger was already repeated. This count is set to zero on trigger and counts up to repeat in CStg200xDownloadBasicNet::SetupTrigger.
- void ForceStatusEvent ()

Force a status event.

void ResetStatus (uint32 t triggermap)

Reset the status flag.

uint32 t GetMemoryUsageDAC (uint32 t Channel)

Queries the memory usage of the current segment and selected analog DAC channel.

uint32_t GetMemoryUsageSyncout (uint32_t Channel)

Queries the memory usage of the current segment and selected syncout channel.

virtual void ClearSyncData (uint32_t channel)

Delete a SyncOut pattern for a channel from STG memory.

virtual void SendSyncData (uint32_t channel, array< uint16_t >^ pData, array< uint64_t >^ tData)

Uploads sync output data to the STG.

Sends sync output data to a given channel on the STG. The list of datapoints will be sent to the selected sync output channel. Sync output data previously sent to the channel is overwritten.

Each datapoint is represented by an integer value and can be either 0 or 1.

The duration is given as a list of 64 bit integers. Durations are given in units of μs. The STG has a resolution of 20 μs. If your application can not handle 64 bit integers, use the STG200x_SendSyncData32() call instead.

virtual void ClearChannelData (uint32_t channel)

Delete a stimulus pattern for a channel from STG memory

virtual void SendChannelData (uint32_t channel, array< uint16_t >^ pData, array< uint64_t >^ tData)

Uploads analog data (stimulus patterns) to the STG.

Sends datapoints to a given channel on the STG. The list of datapoints will be sent to the selected channel. Data previously sent to the channel is overwritten.

Each datapoint is represented by an integer value in the range from 0 to 4095 (bit 0 to 11), its sign is taken from bit 12, 0 is for positive amplitude, and 1 for negative amplitude Bits 13 to 15 have to be zero.

The duration is given as a list of 64 bit integers. Durations are given in units of μs. The STG has a resolution of 20 μs.

virtual void EnableAutoReset ()

Enable AutoReset of the STG Status.

virtual void DisableAutoReset ()

Disable AutoReset of the STG Status.

 virtual void SetupRetriggerMode (int8_t trigger, RetriggerActionEnumNet same_trigger, RetriggerAction← EnumNet other trigger)

Define the action on triggers while the STG is running.

The STG has three options how to handle a successive trigger while a trigger is active.

- stop this trigger (default action)
- restart this trigger
- ignore the signal
- virtual void SetupRetriggerMode (RetriggerActionEnumNet same_trigger, RetriggerActionEnumNet other_
 trigger)

Define the action on triggers while the STG is running.

The STG has three options how to handle a successive trigger while a trigger is active.

- stop this trigger (default action)
- restart this trigger
- ignore the signal

Properties

• CStimulusFunctionNet[^] Stimulus [get]

Additional Inherited Members

11.110.1 Detailed Description

CStg200xDownloadBasicNet is the base class to control the download mode of the MCS STG device.

11.110.2 Member Function Documentation

```
11.110.2.1 ClearChannelData() virtual void ClearChannelData ( uint32_t channel) [virtual]
```

Delete a stimulus pattern for a channel from STG memory

Parameters

channel	Specifies the channel to clear.
---------	---------------------------------

```
11.110.2.2 ClearSyncData() virtual void ClearSyncData ( uint32_t channel ) [virtual]
```

Delete a SyncOut pattern for a channel from STG memory.

Parameters

channel Specifies the syncout channel to clear.

11.110.2.3 DisableAutoReset() virtual void DisableAutoReset () [virtual]

Disable AutoReset of the STG Status.

If autoreset is disabled, the STG status switches to FINISHED after the defined number of sweeps is finished. To switch back to the IDLE status, use CStg200xDownload::ResetStatus()

```
11.110.2.4 EnableAutoReset() virtual void EnableAutoReset ( ) [virtual]
```

Enable AutoReset of the STG Status.

This is the default on power up. If autoreset is enabled, the STG status switches to FINISHED only for one poll cycle after this, it switches to IDLE automatically.

```
11.110.2.5 ForceStatusEvent() void ForceStatusEvent ()
```

Force a status event.

Force the DLL to create a PollMessage event and to call the pPollCallback function, even if no new status information is available.

```
11.110.2.6 GetMemoryUsageDAC() uint32_t GetMemoryUsageDAC ( uint32_t Channel)
```

Queries the memory usage of the current segment and selected analog DAC channel.

The currently used memory is reported for the requested channel.

Parameters

Channel	channel for the amount of interested usage.
---------	---

Returns

Returns the usage in STG memory.

```
11.110.2.7 GetMemoryUsageSyncout() uint32_t GetMemoryUsageSyncout ( uint32_t Channel)
```

Queries the memory usage of the current segment and selected syncout channel.

The currently used memory is reported for the requested channel.

Channel	channel for the amount of interested usage.

Returns

Returns the usage in STG memory.

```
11.110.2.8 GetSweepCount() void GetSweepCount (

[Out] array< uint32_t >^% sweeps,

[Out] array< uint32_t >^% triggers )
```

Get the sweep and trigger count of the STG.

- The triggercount tells how many times each trigger was active and is reset to zero on download of new channel data.
- The sweepcount tells how many times each trigger was already repeated. This count is set to zero on trigger and counts up to repeat in CStg200xDownloadBasicNet::SetupTrigger.

Parameters

sweeps	on return contains the number of sweeps for each trigger.
triggers	on return contains the number of trigger events seen for each trigger.

```
11.110.2.9 GetTrigger() void GetTrigger (

[Out] array< uint32_t >^% channelmap,

[Out] array< uint32_t >^% syncoutmap,

[Out] array< uint32_t >^% repeat )
```

Queries the trigger settings for the STG. Note that all memory segments have their own trigger setting.

Parameters

```
channelmap For each trigger, a bitmap of channels that belong to this trigger.
```

syncoutmap	For each trigger, a bitmap of syncouts that belong to this trigger.
repeat	For each trigger, define the number of times this trigger should be repeated.

```
11.110.2.10 ResetStatus() void ResetStatus (
uint32_t triggermap )
```

Reset the status flag.

	triggermap	bitmap of trigger for which to reset the status.
--	------------	--

Uploads analog data (stimulus patterns) to the STG.

Sends datapoints to a given channel on the STG. The list of datapoints will be sent to the selected channel. Data previously sent to the channel is overwritten.

Each datapoint is represented by an integer value in the range from 0 to 4095 (bit 0 to 11), its sign is taken from bit 12, 0 is for positive amplitude, and 1 for negative amplitude Bits 13 to 15 have to be zero.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Parameters

channel	Specifies the channel to append the data to.
pData	A list of datapoints.
tData	A list of durations as int64_t. The time is given in units of μs.

Uploads sync output data to the STG.

Sends sync output data to a given channel on the STG. The list of datapoints will be sent to the selected sync output channel. Sync output data previously sent to the channel is overwritten.

Each datapoint is represented by an integer value and can be either 0 or 1.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s. If your application can not handle 64 bit integers, use the STG200x_SendSyncData32() call instead.

channel	Specifies the sync output channel to append the data to.
pData	A list of datapoints.
tData	A list of durations as int64_t. The time is given in units of μs.

Define the action on triggers while the STG is running.

The STG has three options how to handle a successive trigger while a trigger is active.

- stop this trigger (default action)
- · restart this trigger
- · ignore the signal

Parameters

trigger	The trigger to change.
same_trigger	Action for successive triggers in Normal Mode, and for triggers to the currently selected segment in Multi - File Mode.
other_trigger	Action for successive triggers in Multi-File Mode for a trigger on a segment not currently selected.Not used in Normal Mode.

```
11.110.2.14 SetupRetriggerMode() [2/2] virtual void SetupRetriggerMode (
RetriggerActionEnumNet same_trigger,
RetriggerActionEnumNet other_trigger) [virtual]
```

Define the action on triggers while the STG is running.

The STG has three options how to handle a successive trigger while a trigger is active.

- stop this trigger (default action)
- · restart this trigger
- · ignore the signal

Parameters

same_trigger	Action for successive triggers in Normal Mode, and for triggers to the currently selected segment in Multi - File Mode.
other_trigger	Action for successive triggers in Multi-File Mode for a trigger on a segment not currently selected. Not used in Normal Mode.

11.110.2.15 SetupTrigger() virtual void SetupTrigger (

```
uint32_t first_trigger,
array< uint32_t >^ channelmap,
array< uint32_t >^ syncoutmap,
array< uint32_t >^ repeat ) [virtual]
```

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

Parameters

first_trigger	The number of the first trigger to change.
---------------	--

Parameters

channelmap	For each trigger, a bitmap of channels that belong to this trigger.
------------	---

Parameters

syncoutmap	For each trigger, a bitmap of syncouts that belong to this trigger.
repeat	For each trigger, define the number of times this trigger should be repeated.

uint32_t repeat) [virtual]

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

Parameters

trigger	The trigger to change.

channelmap	A bitmap of channels that belong to this trigger.

syncoutmap	A bitmap of syncouts that belong to this trigger.
repeat	The number of times this trigger should be repeated.


11.110.3 Property Documentation

11.110.3.1 Stimulus CStimulusFunctionNet^ Stimulus [get]

11.111 CStg200xDownloadNet Class Reference

Main class for the STG download mode This class implements the STG download mode interface.

Inheritance diagram for CStg200xDownloadNet:



Public Member Functions

CStg200xDownloadNet ()

Use this constructor if you do not want to use the status callback.

CStg200xDownloadNet (OnStgPollStatus[^] pollStatus)

Use this constructor if you want to use the status callback.

- ∼CStg200xDownloadNet ()
- void PrepareAndSendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType)

Prepare and send data to a given channel on the STG. Previous data sent to that channel is erased first.

 void PrepareAndAppendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType)

Prepare and append data to a given channel on the STG.

void ClearChannel_PrepareAndSendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType, bool doClear)

Prepare and append data to a given channel on the STG.

void SendSegmentDefine (array< uint32_t >^ segment_list)

Defines the segment memory layout of the STG.

 void SendSegmentStart (uint32_t triggermap, uint32_t segment, Stg200xSegmentFlagsEnumNet segmentflags)

Switchs segment and starts trigger.

- void SendSegmentSelect (uint32_t segment, Stg200xSegmentFlagsEnumNet segmentflags)
 Switchs segment.
- void EnableMultiFileMode (uint32_t submode)

Enable the Multi-File mode of the STG.

void DisableMultiFileMode ()

Disable the Multi-File mode of the STG

- StgStatusNet ^ QueryTriggerstatus ()
- void SetOutputMap (array< uint32_t $>^{\land}$ ChannelLayout)
- int32 t GetModuleTemp (unsigned int channel)
- uint32_t GetModuleCurrent (unsigned int channel)

Events

- OnStgPollStatus Stg200xPollStatusEvent [add, remove, raise]
- OnMwPollStatus^ MwPollStatusEvent [add, remove, raise]

Additional Inherited Members

11.111.1 Detailed Description

Main class for the STG download mode This class implements the STG download mode interface.

11.111.2 Constructor & Destructor Documentation

```
11.111.2.1 CStg200xDownloadNet() [1/2] CStg200xDownloadNet ( )
```

Use this constructor if you do not want to use the status callback.

```
11.111.2.2 CStg200xDownloadNet() [2/2] CStg200xDownloadNet ( OnStgPollStatus^{\wedge} pollStatus )
```

Use this constructor if you want to use the status callback.

11.111.2.3 ~CStg200xDownloadNet() ~CStg200xDownloadNet ()

11.111.3 Member Function Documentation

$\textbf{11.111.3.1} \quad \textbf{ClearChannel_PrepareAndSendData()} \quad \texttt{void ClearChannel_PrepareAndSendData} \quad \textbf{(}$

```
uint32_t channel,
array< int32_t >^ amplitude,
array< uint64_t >^ duration,
STG_DestinationEnumNet destType,
bool doClear )
```

Prepare and append data to a given channel on the STG.

Each datapoint is represented by an signed 32bit integer value. When using voltage stimulation, the values are in multiple of 1 uV, thus the possible range is += 2000 V. When using current stimulation, the values are in multiple of 1 nA, this the possible range is += 2000 mA.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Blocks of data which should repeat can be defined by prepending such a block with an entry in the arrays where both amplitude and duration is zero. The end of such an block is marked by an entry where the duration is set to zero and the amplitude beeing set to the number of times the block should run. Blocks can be nested.

Parameters

channel T	he channel number to send data to.
-----------	------------------------------------

Parameters

amplitude	A list of amplitudes in units of μV and nA in voltage and current mode, respectively.

Parameters

duration A list of durations in units of μs.		A list of durations in units of μs.
	destType	specifies wheather the data is for syncout, current or voltage stimulation.

11.111.3.2 DisableMultiFileMode() void DisableMultiFileMode ()

Disable the Multi-File mode of the STG

Switch the STG back to normal mode. In this mode, trigger inputs are assigned to channels, not to segments.

```
11.111.3.3 EnableMultiFileMode() void EnableMultiFileMode (
             uint32_t submode )
```

Enable the Multi-File mode of the STG.

In Multi-File mode, the trigger inputs switch between segments. To use this mode, define up to as many segments as trigger inputs are available and fill each segment with a stimulus pattern.

Now a trigger on trigger input 1 switches the STG to the first segment and starts all triggers in this segment. Likewise, a trigger on trigger input 2, 3 and 4 selects the respective segment and start all triggers in this segment So the Multi-File Mode can be used to predefine up to four different stimuli which can be selected without the need for a computer connection.

Parameters

submode

The submode. Submode 0 is regular Multi-File mode as described above, submode 1 is extended Multi-File mode, where the segment is selected based on the digital pattern on the digital inputs. In this mode, 256 different segments can be defined and used.

```
11.111.3.4 GetModuleCurrent() uint32_t GetModuleCurrent (
             unsigned int channel )
```

```
11.111.3.5 GetModuleTemp() int32_t GetModuleTemp (
            unsigned int channel )
```

```
11.111.3.6 PrepareAndAppendData() void PrepareAndAppendData (
```

```
uint32_t channel,
array< int32_t >^{\wedge} amplitude,
array< uint64_t >^{\wedge} duration,
STG_DestinationEnumNet destType )
```

Prepare and append data to a given channel on the STG.

Each datapoint is represented by an signed 32bit integer value. When using voltage stimulation, the values are in multiple of 1 uV, thus the possible range is += 2000 V. When using current stimulation, the values are in multiple of 1 nA, this the possible range is += 2000 mA.

The duration is given as a list of 64 bit integers. Durations are given in units of µs. The STG has a resolution of 20 μs.

Blocks of data which should repeat can be defined by prepending such a block with an entry in the arrays where both amplitude and duration is zero. The end of such an block is marked by an entry where the duration is set to zero and the amplitude beeing set to the number of times the block should run. Blocks can be nested.

channel	The channel number to send data to.
---------	-------------------------------------

Parameters

itude \mid A list of amplitudes in units of μV and nA in voltage and current mode, respectively.	ctively.
---	----------

Parameters

duration A list of durations in units of μs.		A list of durations in units of μs.
	destType	specifies wheather the data is for syncout, current or voltage stimulation.

11.111.3.7 PrepareAndSendData() void PrepareAndSendData (

```
uint32_t channel,
array< int32_t >^ amplitude,
array< uint64_t >^ duration,
STG_DestinationEnumNet destType )
```

Prepare and send data to a given channel on the STG. Previous data sent to that channel is erased first.

Each datapoint is represented by an signed 32bit integer value. When using voltage stimulation, the values are in multiple of 1 uV, thus the possible range is += 2000 V. When using current stimulation, the values are in multiple of 1 nA, this the possible range is += 2000 mA.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Blocks of data which should repeat can be defined by prepending such a block with an entry in the arrays where both amplitude and duration is zero. The end of such an block is marked by an entry where the duration is set to zero and the amplitude beeing set to the number of times the block should run. Blocks can be nested.

Parameters

channel The channel number to send data to	ata to.
--	---------

amplitude	A list of amplitudes in units of μV and nA in voltage and current mode, respectively.

duration	A list of durations in units of μs.
destType	specifies wheather the data is for syncout, current or voltage stimulation.

11.111.3.8 QueryTriggerstatus() StgStatusNet ^ QueryTriggerstatus ()

```
11.111.3.9 SendSegmentDefine() void SendSegmentDefine ( array < uint32_t >^{\land} segment\_list )
```

Defines the segment memory layout of the STG.

On reset, the STG has one segment containing all available memory.

With this command, the STG memory can be devided into several segments. Each segment can be filled with stimulus data.

Parameters

	segment_list	The List of memory sizes (one per segment).
--	--------------	---

Switchs segment.

Parameters

segment	The number of the segment to select.

Switchs segment and starts trigger.

Parameters

triggermap A bitmap of triggers that will be sta
--

Parameters

Parameters

segmentflags | A bitmap of flags, bit 1: assign all channels to the trigger number equal to the segment.

```
11.111.3.12 SetOutputMap() void SetOutputMap (

array< uint32_t >^ ChannelLayout )
```


11.111.4 Event Documentation

```
11.111.4.1 MwPollStatusEvent OnMwPollStatus^ MwPollStatusEvent [add], [remove], [raise]
```

```
11.111.4.2 Stg200xPollStatusEvent OnStgPollStatus^ Stg200xPollStatusEvent [add], [remove], [raise]
```

11.112 CStimulusFunctionNet Class Reference

Inheritance diagram for CStimulusFunctionNet:



Classes

- class SidebandData
- class StimulusDeviceDataAndUnrolledData

Public Member Functions

- CStimulusFunctionNet (CMcsUsbNet[∧] mcsusb, CMcsUsbFunctionPointerContainer[∧] stimulusFunction
 —
 PointerContainer)
- CStimulusFunctionNet (CMcsUsbNet[^] mcsusb)
- void StartPoll ()

Starts the interrupt fetching thread and delivers events

void StopPoll ()

Stops the interrupt fetching thread and delivers events

void ForceStatusEvent ()

Force a status event. Force the DLL to create a PollMessage event and to call the pPollCallback function, even if no new status information is available.

void SendStart (uint32_t triggermap)

Start (Trigger) the STG. The startup delay is in the range of a few ms.

void SendStop (uint32 t triggermap)

Stop some or all triggers of the STG.

void SendStop (uint32_t triggermap, int options)

Stop some or all triggers of the STG.

· void ClearChannelData (int channel)

Delete a Stimulus Pattern from STG memory

• void ClearSyncData (int channel)

Delete a Syncout Pattern from STG memory

 void PrepareAndSendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType)

Prepare and send data to a given channel on the STG. Previous data sent to that channel is erased first.

 void PrepareAndAppendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType)

Prepare and append data to a given channel on the STG.

- void ClearChannel_PrepareAndSendData (uint32_t channel, array< int32_t >^ amplitude, array< uint64_t >^ duration, STG_DestinationEnumNet destType, bool doClear)
- StimulusDeviceDataAndUnrolledData ^ PrepareData (int channel, array< int32_t > ^ amplitude, array< uint64 t > ^ duration, STG DestinationEnumNet destType)
- void SendPreparedData (int channel, StimulusDeviceDataAndUnrolledData[^] device_data_and_unrolled, S
 — TG_DestinationEnumNet destType)
- SidebandData ^ CreateSideband (array< int32_t >^ StimulusActive, array< int32_t >^ Syncout, array< uint64_t >^ Duration, uint32_t Bit0Time, uint32_t Bit3Time, uint32_t Bit4Time)

Creates the Sideband Channel for the MEA2100 device.

void ClearMultiplexedData ()

Clears the Stimulation Memory in the STG device.

void SendMultiplexedData (array< uint16 t >^ data)

Sends stimulus data in multiplexed form. All 16 bits words for the enabled DAC and digital channels are muxed together per time slice.

int GetMultiplexedDataChannelsInBlock ()

Gets the number of stimulus data channels to send per time slice. Might be greater than the number of configured channels. Fill unused channels with dummy data in SendMultiplexedData

• int GetDACResolution ()

Gets number of bits of the DAC resolution.

int GetVoltageRangeInMicroVolt (uint32_t channel)

Gets the Voltage Range of the specified channel in Microvolts.

int GetVoltageResolutionInMicroVolt (uint32 t channel)

Gets the Voltage Resolution of the specified channel in Microvolts.

• int GetCurrentRangeInNanoAmp (uint32_t channel)

Gets the Current Range of the specified channel in Nanoamps.

• int GetCurrentResolutionInNanoAmp (uint32_t channel)

Gets the Current Resolution of the specified channel in Nanoamps.

void SetupTrigger (uint32_t first_trigger, array< uint32_t >^ channelmap, array< uint32_t >^ syncoutmap, array< uint32_t >^ repeat)

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

void SetupTriggerSingle (uint32_t trigger, uint32_t channelmap, uint32_t syncoutmap, uint32_t repeat)

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

uint32_t GetTotalMemory ()

Get the total amount of memory available on the STG (all segments).

uint32 t GetAvailableMemory ()

Get the amount of memory available in the currently selected segment of the STG.

int GetNumberOfAnalogChannels ()

Get the number of STG channels.

Events

OnStgPollStatus^ PollStatusEvent

Additional Inherited Members

11.112.1 Constructor & Destructor Documentation

```
11.112.1.1 CStimulusFunctionNet() [1/2] CStimulusFunctionNet (

CMcsUsbNet^ mcsusb,

CMcsUsbFunctionPointerContainer^ stimulusFunctionPointerContainer)
```

```
11.112.1.2 CStimulusFunctionNet() [2/2] CStimulusFunctionNet (
CMcsUsbNet^ mcsusb )
```

11.112.2 Member Function Documentation

Delete a Stimulus Pattern from STG memory

Parameters

channel specifies the channel to clear.

11.112.2.3 ClearMultiplexedData() void ClearMultiplexedData ()

Clears the Stimulation Memory in the STG device.

Delete a Syncout Pattern from STG memory

Parameters

channel specifies the channel to clear.

Creates the Sideband Channel for the MEA2100 device.

Each datapoint is represented by an signed 32bit integer value. A value 0 means that the stimulation is active during that time. A value 1 means that the stimulation is not active during that time.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Parameters

StimulusActive	A list of datapoints which define weather the Stimulus is active or idle at that time as int32.
----------------	---

Parameters

Duration	A list of durations as uint64. The time is given in units of μs.
Bit0Time Time in μs for which Bit 0 (Blanking) is to be extended.	

Parameters

Bit3Time	Time in µs for which Bit 3 (Stimulus Enable) is to be extended.
----------	---

Parameters

Bit4Time	Time in μs for which Bit 4 (Stimulus Selector) is to be extended.
----------	---

Returns

Error Status. 0 on success.

11.112.2.6 ForceStatusEvent() void ForceStatusEvent ()

Force a status event. Force the DLL to create a PollMessage event and to call the pPollCallback function, even if no new status information is available.

11.112.2.7 GetAvailableMemory() uint32_t GetAvailableMemory ()

Get the amount of memory available in the currently selected segment of the STG.

Returns

The total memory available on the STG in bytes.

Gets the Current Range of the specified channel in Nanoamps.

Parameters

channel	Channel which is queried.
---------	---------------------------

Returns

The Current Range of the specified channel in Nanoamps.

11.112.2.9 GetCurrentResolutionInNanoAmp() int GetCurrentResolutionInNanoAmp (uint32_t channel)

Gets the Current Resolution of the specified channel in Nanoamps.

Parameters

channel	Channel which is queried.
---------	---------------------------

Returns

The Current Resolution of the specified channel in Nanoamps.

11.112.2.10 GetDACResolution() int GetDACResolution ()

Gets number of bits of the DAC resolution.

Returns

The DAC resolution in bits.

11.112.2.11 GetMultiplexedDataChannelsInBlock() int GetMultiplexedDataChannelsInBlock ()

Gets the number of stimulus data channels to send per time slice. Might be greater than the number of configured channels. Fill unused channels with dummy data in SendMultiplexedData

11.112.2.12 GetNumberOfAnalogChannels() int GetNumberOfAnalogChannels ()

Get the number of STG channels.

Returns

The number of STG channels.

11.112.2.13 GetTotalMemory() uint32_t GetTotalMemory ()

Get the total amount of memory available on the STG (all segments).

Returns

The total memory available on the STG in bytes.

11.112.2.14 GetVoltageRangeInMicroVolt() int GetVoltageRangeInMicroVolt (uint32_t channel)

Gets the Voltage Range of the specified channel in Microvolts.

Parameters

channel	Channel which is queried.
---------	---------------------------

Returns

The Voltage Range of the specified channel in Microvolts.

11.112.2.15 **GetVoltageResolutionInMicroVolt()** int GetVoltageResolutionInMicroVolt (uint32_t channel)

Gets the Voltage Resolution of the specified channel in Microvolts.

channel	Channel which is queried.
channel	Channel which is queried.

Returns

The Voltage Resolution of the specified channel in Microvolts.

11.112.2.16 PrepareAndAppendData() void PrepareAndAppendData (

```
uint32_t channel,
array< int32_t >^ amplitude,
array< uint64_t >^ duration,
STG_DestinationEnumNet destType )
```

Prepare and append data to a given channel on the STG.

Each datapoint is represented by an signed 32bit integer value. When using voltage stimulation, the values are in multiple of 1 uV, thus the possible range is += 2000 V. When using current stimulation, the values are in multiple of 1 nA, this the possible range is += 2000 mA.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Blocks of data which should repeat can be defined by prepending such a block with an entry in the arrays where both amplitude and duration is zero. The end of such an block is marked by an entry where the duration is set to zero and the amplitude beeing set to the number of times the block should run. Blocks can be nested.

Parameters

channel	The channel number to send data to.

Parameters

amplitude	A list of amplitudes in units of μV and nA in voltage and current mode, respectively.
-----------	--

Parameters

duration	A list of durations in units of μs.
destType	specifies wheather the data is for syncout, current or voltage stimulation.

Returns

Error Status. 0 on success.

11.112.2.17 PrepareAndSendData() void PrepareAndSendData (

```
uint32_t channel,
array< int32_t >^ amplitude,
array< uint64_t >^ duration,
STG_DestinationEnumNet destType )
```

Prepare and send data to a given channel on the STG. Previous data sent to that channel is erased first.

Each datapoint is represented by an signed 32bit integer value. When using voltage stimulation, the values are in multiple of 1 uV, thus the possible range is += 2000 V. When using current stimulation, the values are in multiple of 1 nA, this the possible range is += 2000 mA.

The duration is given as a list of 64 bit integers. Durations are given in units of μ s. The STG has a resolution of 20 μ s.

Blocks of data which should repeat can be defined by prepending such a block with an entry in the arrays where both amplitude and duration is zero. The end of such an block is marked by an entry where the duration is set to zero and the amplitude beeing set to the number of times the block should run. Blocks can be nested.

Parameters

nel The channel number to send data to.

Parameters

	amplitude	A list of amplitudes in units of μV and nA in voltage and current mode, respectively.
--	-----------	--

Parameters

duration	A list of durations in units of μs.
destType	specifies wheather the data is for syncout, current or voltage stimulation.

Returns

Error Status. 0 on success.

```
11.112.2.18 PrepareData() StimulusDeviceDataAndUnrolledData ^ PrepareData (
    int channel,
    array< int32_t >^ amplitude,
    array< uint64_t >^ duration,
    STG_DestinationEnumNet destType )
```

```
11.112.2.19 SendMultiplexedData() void SendMultiplexedData ( array < uint16_t >^{\wedge} data )
```

Sends stimulus data in multiplexed form. All 16 bits words for the enabled DAC and digital channels are muxed together per time slice.

Parameters

```
data Array of data to be sent.
```

```
11.112.2.21 SendStart() void SendStart ( uint32_t triggermap )
```

Start (Trigger) the STG. The startup delay is in the range of a few ms.

Parameters

triggermap A bitmap of triggers which will be started.

```
11.112.2.22 SendStop() [1/2] void SendStop ( uint32_t triggermap )
```

Stop some or all triggers of the STG.

Parameters

triggermap	A bitmap of triggers which will be stopped.

```
11.112.2.23 SendStop() [2/2] void SendStop (
          uint32_t triggermap,
          int options )
```

Stop some or all triggers of the STG.

Parameters

triggermap	A bitmap of triggers which will be stopped.
options	bitmap of options, currently only STOP_OPTION_SAVESTOP (0x80) is defined, which bypasses
	the stop commands when a syncout assossiated with a given sync-out has bit 1 (0x02) set. Can
	be used e.g. to prevent a stop while a biphasic stimulation pulse is active

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.

Parameters

irst_trigger The number of the first trigger to chan	e.
--	----

Parameters

channelmap	For each trigger, a bitmap of channels that belong to this trigger.

Parameters

syncoutmap	For each trigger, a bitmap of syncouts that belong to this trigger.
repeat	For each trigger, define the number of times this trigger should be repeated.

Configures the trigger settings for the STG. Note that all memory segments have their own trigger setting.	

trigger	The trigger to change.
---------	------------------------

Parameters

channelmap	A bitmap of channels that belong to this trigger.
------------	---

Parameters

syncoutmap	A bitmap of syncouts that belong to this trigger.
repeat	The number of times this trigger should be repeated.

11.112.2.26 StartPoll() void StartPoll ()

Starts the interrupt fetching thread and delivers events

$\textbf{11.112.2.27} \quad \textbf{StopPoll()} \quad \texttt{void StopPoll ()}$

Stops the interrupt fetching thread and delivers events


11.112.3 Event Documentation

11.112.3.1 PollStatusEvent OnStgPollStatus^ PollStatusEvent

11.113 CSw2to64DeviceNet Class Reference

The class to control the MCS-USB-Sw2to64 device.

Inheritance diagram for CSw2to64DeviceNet:



Public Member Functions

- CSw2to64DeviceNet ()
- ∼CSw2to64DeviceNet ()
- unsigned short GetNumber ()

Gets the number of channels that can be switched in this box.

array< unsigned char > ^ GetChannels ()

Gets the current switch positions as char array.

void SetChannels (array< unsigned char >^ pattern)

Sets the switch positions from a char array.

unsigned char GetChannel (unsigned short index)

Gets one current switch position.

void SetChannel (unsigned short index, unsigned char pattern)

Sets one switch position.

Additional Inherited Members

11.113.1 Detailed Description

The class to control the MCS-USB-Sw2to64 device.

This class controls the settings of the MCS-USB-Sw2to64. The box has two inputs for signals. Each of the 64 outputs can be connected to one of the input signals, could be held open or connected ground. Valid switch states are 0, 1, 2 or 3 for each of the settings.

11.113.2 Constructor & Destructor Documentation

```
11.113.2.1 CSw2to64DeviceNet() CSw2to64DeviceNet ()
```

```
11.113.2.2 ~CSw2to64DeviceNet() ~CSw2to64DeviceNet ()
```

11.113.3 Member Function Documentation

```
11.113.3.1 GetChannel() unsigned char GetChannel ( unsigned short index )
```

Gets one current switch position.

in	index	number of channel to read the switch position from	
----	-------	--	--

Returns

switch position of desired channel

11.113.3.2 GetChannels() array<unsigned char> ^ GetChannels ()

Gets the current switch positions as char array.

Returns

array of char with the size of the number of channels, each char has the setting of a channel

11.113.3.3 GetNumber() unsigned short GetNumber ()

Gets the number of channels that can be switched in this box.

The box can have a different number of channels it can switch. Up to now usually 64 channels are returned

```
11.113.3.4 SetChannel() void SetChannel (
unsigned short index,
unsigned char pattern )
```

Sets one switch position.

Parameters

i	n	index	number of channel to write the switch position to
i	n	pattern	switch position of the channel

11.113.3.5 SetChannels() void SetChannels ($array < unsigned char >^{\land} pattern$)

Sets the switch positions from a char array.


Parameters

in	pattern	array of char with the size of the number of channels, each char has the setting of a channel

11.114 CTcxDeviceNet Class Reference

Class to control a Temperature Controller (TCX)

Inheritance diagram for CTcxDeviceNet:



Public Member Functions

CTcxDeviceNet ()

Initializes a new instance of CTcxDeviceNet class.

- ∼CTcxDeviceNet ()
- unsigned int GetNumControlChannels ()

Gets the number of channels the device can control/regulate.

unsigned int GetNumMeasureChannels ()

Gets the number of channels the device can measure.

int GetValue (unsigned int channel)

Gets the temperate of the specified channel in units of 0.1 ℃.

• int GetValueHires (unsigned int channel)

Gets the temperate of the specified channel in units of 0.01 ℃.

int GetHeaterTemp (unsigned int channel)

Gets the temperate of the specified heater in units of 0.1 °C.

• int GetHeaterLimit (unsigned int device)

Gets the temperate limit of the specified heater in units of 0.1 ℃.

double GetMaxHeaterPowerMultiwell ()

queries the max. heater power that the Multiwell temperature controller will apply; unit: W; useful range: 5.2W..7.6W

void SetMaxHeaterPowerMultiwell (double MaxPowerWatt)

sets the max. heater power that the Multiwell temperature controller will apply; unit: W; useful range: 5.2W..7.6W

bool GetHasThermocouple ()

Gets weather the device supports a thermocouple.

- bool GetEnableHeaterLimit (unsigned int device)
- bool GetEnableThermocouple (unsigned int device)
- TcxSensorTypeEnumNet GetSensorType (unsigned int device)
- String ^ GetUnit (unsigned int channel)
- unsigned int GetBoardTemp ()

Gets the temperate of the mainboard in units of 0.1 °C.

- · unsigned int GetVolti (unsigned int channel)
- unsigned int GetNumDevices ()
- void SetSetpoint (unsigned int channel, int sp)

Sets the target temperate of specified channel in units of 0.1 °C.

- void SetDevice (unsigned int channel, int device)
- void SetOnOff (unsigned int channel, bool on)

Switches the specified channel on or off.

- · void SetCalibration (unsigned int channel, int calib)
- void SetP (unsigned int device, int p_coeff)

Sets the P-coefficient of the specified device.

void SetI (unsigned int device, int i coeff)

Sets the I-coefficient of the specified device.

• void SetD (unsigned int device, int d coeff)

Sets the D-coefficient of the specified device.

void SetMaxP (unsigned int device, int maxp)

Sets the maximum heater power of the specified device.

- void SetHeaterLimit (unsigned int device, int heater_limit)
- void SetEnableHeaterLimit (unsigned int device, bool enable)
- void SetEnableThermocouple (unsigned int device, bool enable)
- void SetSensorType (unsigned int device, TcxSensorTypeEnumNet type)
- void SetDevname (unsigned int device, String[^] Devicename)
- int GetSetpoint (unsigned int channel)

Gets the target temperate of specified channel in units of 0.1 ℃.

- int GetDevice (unsigned int channel)
- int GetOnOff (unsigned int channel)

Gets if the specified channel is on or off.

- int GetCalibration (unsigned int channel)
- int GetP (unsigned int device)

Gets the P-coefficient of the specified device.

· int GetI (unsigned int device)

Gets the I-coefficient of the specified device.

int GetD (unsigned int device)

Gets the D-coefficient of the specified device.

• int GetMaxP (unsigned int device)

Gets the maximum heater power of the specified device.

- String \(^\) GetDevname (unsigned int device)
- TcxDeviceTypeEnumNet GetDeviceType ()
- int GetSetpointMin (unsigned int channel)
- int GetCalibrationMin (unsigned int channel)
- int GetPMin (unsigned int device)
- int GetIMin (unsigned int device)
- int GetDMin (unsigned int device)
- int GetMaxpMin (unsigned int device)
- · int GetSetpointMax (unsigned int channel)
- · int GetCalibrationMax (unsigned int channel)
- int GetPMax (unsigned int device)
- int GetIMax (unsigned int device)
- int GetDMax (unsigned int device)
- int GetMaxpMax (unsigned int device)
- int GetSetpointDecp (unsigned int channel)
- int GetCalibrationDecp (unsigned int channel)
- int GetPDecp (unsigned int device)
- int GetIDecp (unsigned int device)
- int GetDDecp (unsigned int device)
- int GetMaxpDecp (unsigned int device)
- · int GetResX (unsigned int channel)
- int GetResS (unsigned int channel)
- int GetRes1 (unsigned int channel)
- int GetRes2 (unsigned int channel)
- int GetPwrSet (unsigned int channel)
- int GetPwrOut (unsigned int channel)
- int GetDuty (unsigned int channel)

Gets the duty cycle of the heating element.

• int GetUOut (unsigned int channel)

Gets the voltage on the heating element.

• int GetlOut (unsigned int channel)

Gets the current through the heating element.

int GetROut (unsigned int channel)

Gets the resistance of the heating element.

int GetPOut (unsigned int channel)

Gets the output power of the heating element.

- int GetCurrent (unsigned int channel)
- int GetThermocoupleTemp (unsigned int channel)
- int GetThermocoupleTempAbs (unsigned int channel)
- int GetThermocoupleReferenceTemp (unsigned int channel)
- unsigned int GetThermocoupleNanovoltPerKelvin (unsigned int channel)

Gets the proportional constant for the thermocouple.

void SetThermocoupleNanovoltPerKelvin (unsigned int channel, unsigned int value)

Sets the proportional constant for the thermocouple.

- int GetThermocoupleCalibration (unsigned int channel)
- void CalibrateThermocouple (unsigned int channel)
- void SetDeviceType (TcxDeviceTypeEnumNet devicetype)
- void FactoryReset ()

Additional Inherited Members

11.114.1 Detailed Description

Class to control a Temperature Controller (TCX)

11.114.2 Constructor & Destructor Documentation

```
11.114.2.1 CTcxDeviceNet() CTcxDeviceNet ()
```

Initializes a new instance of CTcxDeviceNet class.

```
11.114.2.2 ~CTcxDeviceNet() ~CTcxDeviceNet ()
```

11.114.3 Member Function Documentation

```
11.114.3.2 FactoryReset() void FactoryReset ( )
\textbf{11.114.3.3} \quad \textbf{GetBoardTemp()} \quad \texttt{unsigned int GetBoardTemp ()} \\
Gets the temperate of the mainboard in units of 0.1 °C.
\textbf{11.114.3.4} \quad \textbf{GetCalibration()} \quad \texttt{int GetCalibration (}
                unsigned int channel )
11.114.3.5 GetCalibrationDecp() int GetCalibrationDecp (
                unsigned int channel )
\textbf{11.114.3.6} \quad \textbf{GetCalibrationMax()} \quad \texttt{int GetCalibrationMax} \quad \textbf{(}
                unsigned int channel )
11.114.3.7 GetCalibrationMin() int GetCalibrationMin (
                unsigned int channel )
11.114.3.8 GetCurrent() int GetCurrent (
                unsigned int channel )
11.114.3.9 GetD() int GetD (
                unsigned int device )
Gets the D-coefficient of the specified device.
11.114.3.10 GetDDecp() int GetDDecp (
                unsigned int device )
```

```
11.114.3.11 GetDevice() int GetDevice (
             unsigned int channel )
11.114.3.12 GetDeviceType() TcxDeviceTypeEnumNet GetDeviceType ( )
11.114.3.13 GetDevname() String ^ GetDevname (
             unsigned int device )
11.114.3.14 GetDMax() int GetDMax (
             unsigned int device )
11.114.3.15 GetDMin() int GetDMin (
             unsigned int device )
11.114.3.16 GetDuty() int GetDuty (
             unsigned int channel )
Gets the duty cycle of the heating element.
Parameters
 channel The channel number.
Returns
     The duty cycle in percent, the value of 320 \ast 64 corresponds to 100 %.
11.114.3.17 GetEnableHeaterLimit() bool GetEnableHeaterLimit (
             unsigned int device )
```

11.114.3.18 GetEnableThermocouple() bool GetEnableThermocouple (

unsigned int device)

```
\textbf{11.114.3.19} \quad \textbf{GetHasThermocouple()} \quad \texttt{bool GetHasThermocouple ()} \quad \texttt{bool GetHasThermocouple ()}
```

Gets weather the device supports a thermocouple.

```
11.114.3.20 GetHeaterLimit() int GetHeaterLimit ( unsigned int device )
```

Gets the temperate limit of the specified heater in units of 0.1 °C.

```
11.114.3.21 GetHeaterTemp() int GetHeaterTemp (
unsigned int channel)
```

Gets the temperate of the specified heater in units of 0.1 °C.

```
11.114.3.22 Getl() int GetI ( unsigned int device )
```

Gets the I-coefficient of the specified device.

```
11.114.3.23 GetIDecp() int GetIDecp (
unsigned int device)
```

```
11.114.3.24 GetlMax() int GetlMax (
unsigned int device)
```

```
11.114.3.25 GetlMin() int GetIMin (
unsigned int device)
```

```
11.114.3.26 GetlOut() int GetlOut ( unsigned int channel)
```

Gets the current through the heating element.

channel The channel number.

Returns

The current in units of mA.

```
11.114.3.27 GetMaxHeaterPowerMultiwell() double GetMaxHeaterPowerMultiwell ( )
```

queries the max. heater power that the Multiwell temperature controller will apply; unit: W; useful range: 5.2W..7.6W

```
11.114.3.28 GetMaxP() int GetMaxP (
unsigned int device)
```

Gets the maximum heater power of the specified device.

```
11.114.3.29 GetMaxpDecp() int GetMaxpDecp ( unsigned int device )
```

```
11.114.3.30 GetMaxpMax() int GetMaxpMax ( unsigned int device )
```

```
11.114.3.31 GetMaxpMin() int GetMaxpMin ( unsigned int device )
```

11.114.3.32 GetNumControlChannels() unsigned int GetNumControlChannels ()

Gets the number of channels the device can control/regulate.

11.114.3.33 GetNumDevices() unsigned int GetNumDevices ()

```
\textbf{11.114.3.34} \quad \textbf{GetNumMeasureChannels()} \quad \texttt{unsigned int GetNumMeasureChannels ()}
```

Gets the number of channels the device can measure.

```
11.114.3.35 GetOnOff() int GetOnOff (
          unsigned int channel )
```

Gets if the specified channel is on or off.

```
11.114.3.36 GetP() int GetP ( unsigned int device )
```

Gets the P-coefficient of the specified device.

```
11.114.3.37 GetPDecp() int GetPDecp ( unsigned int device )
```

```
11.114.3.38 GetPMax() int GetPMax (
unsigned int device)
```

```
11.114.3.39 GetPMin() int GetPMin (
unsigned int device)
```

```
11.114.3.40 GetPOut() int GetPOut (
unsigned int channel)
```

Gets the output power of the heating element.

Parameters

```
channel The channel number.
```

Returns

The resistance in units of mW.

```
11.114.3.41 GetPwrOut() int GetPwrOut (
             unsigned int channel )
11.114.3.42 GetPwrSet() int GetPwrSet (
             unsigned int channel )
11.114.3.43 GetRes1() int GetRes1 (
             unsigned int channel )
11.114.3.44 GetRes2() int GetRes2 (
             unsigned int channel )
11.114.3.45 GetResS() int GetResS (
             unsigned int channel )
11.114.3.46 GetResX() int GetResX (
             unsigned int channel )
11.114.3.47 GetROut() int GetROut (
             unsigned int channel )
Gets the resistance of the heating element.
Parameters
 channel The channel number.
```

Returns

The resistance in units of 0.1 Ohm.

```
11.114.3.48 GetSensorType() TcxSensorTypeEnumNet GetSensorType ( unsigned int device )
```

```
11.114.3.49 GetSetpoint() int GetSetpoint ( unsigned int channel )
```

Gets the target temperate of specified channel in units of $0.1\,^{\circ}$ C.

```
11.114.3.53 GetThermocoupleCalibration() int GetThermocoupleCalibration ( unsigned int channel)
```

```
11.114.3.54 GetThermocoupleNanovoltPerKelvin() unsigned int GetThermocoupleNanovoltPerKelvin ( unsigned int channel)
```

Gets the proportional constant for the thermocouple.

Parameters

channel	Thermocouple channel number.
---------	------------------------------

Returns

The proportional constant in Nanovolt per Kelvin.

```
11.114.3.55 GetThermocoupleReferenceTemp() int GetThermocoupleReferenceTemp ( unsigned int channel)
```

Gets the voltage on the heating element.

Parameters

channel	The channel number.
---------	---------------------

Returns

The voltage in units of mV.

```
11.114.3.60 GetValue() int GetValue ( unsigned int channel )
```

Gets the temperate of the specified channel in units of 0.1 °C.

```
11.114.3.61 GetValueHires() int GetValueHires ( unsigned int channel)
```

Gets the temperate of the specified channel in units of 0.01 ℃.

```
11.114.3.62 GetVolti() unsigned int GetVolti ( unsigned int channel )
```

```
11.114.3.63 SetCalibration() void SetCalibration (
             unsigned int channel,
             int calib )
11.114.3.64 SetD() void SetD (
             unsigned int device,
             int d\_coeff )
Sets the D-coefficient of the specified device.
11.114.3.65 SetDevice() void SetDevice (
             unsigned int channel,
             int device )
11.114.3.66 SetDeviceType() void SetDeviceType (
             TcxDeviceTypeEnumNet devicetype )
11.114.3.67 SetDevname() void SetDevname (
             unsigned int device,
             String^ Devicename )
11.114.3.68 SetEnableHeaterLimit() void SetEnableHeaterLimit (
             unsigned int device,
             bool enable )
11.114.3.69 SetEnableThermocouple() void SetEnableThermocouple (
             unsigned int device,
             bool enable )
11.114.3.70 SetHeaterLimit() void SetHeaterLimit (
             unsigned int device,
             int heater_limit )
```

```
11.114.3.71 Setl() void SetI (
         unsigned int device,
         int i_coeff )
```

Sets the I-coefficient of the specified device.

```
11.114.3.72 SetMaxHeaterPowerMultiwell() void SetMaxHeaterPowerMultiwell (

double MaxPowerWatt )
```

sets the max. heater power that the Multiwell temperature controller will apply; unit: W; useful range: 5.2W..7.6W

```
11.114.3.73 SetMaxP() void SetMaxP (
unsigned int device,
int maxp )
```

Sets the maximum heater power of the specified device.

```
11.114.3.74 SetOnOff() void SetOnOff (
     unsigned int channel,
     bool on )
```

Switches the specified channel on or off.

Parameters

```
channel The channel number.
```

```
11.114.3.75 SetP() void SetP (
    unsigned int device,
    int p\_coeff)
```

Sets the P-coefficient of the specified device.

```
11.114.3.77 SetSetpoint() void SetSetpoint (
     unsigned int channel,
     int sp )
```

Sets the target temperate of specified channel in units of 0.1 ℃.

```
11.114.3.78 SetThermocoupleNanovoltPerKelvin() void SetThermocoupleNanovoltPerKelvin (
unsigned int channel,
unsigned int value)
```

Sets the proportional constant for the thermocouple.


Parameters

channel	Thermocouple channel number.
value	Proportinal constant in Nanovolt per Kelvin.

11.115 CTEERFunctionNet Class Reference

CTEERFunctionNet is the class to control the TEER device

Inheritance diagram for CTEERFunctionNet:



Public Member Functions

CTEERFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pTEERFunctionPointer
 — Container)

Initializes a new instance of the CTEERFunctionNet class.

- CTEERFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ∼CTEERFunctionNet ()
- !CTEERFunctionNet ()
- uint32_t GetPeriod_us ()

gets the period of TEER stimulation in us

void SetPeriod_us (uint32_t period_us)

sets the period of TEER stimulation in us

• uint32 t GetAmplitude nA ()

gets TEER stimulation amplitude in nA

void SetAmplitude_nA (uint32_t Amplitude_nA)

sets TEER stimulation amplitude in nA

• TeerWaveformEnumNet GetWaveform ()

gets TEER stimulation waveform (sine/rect)

```
    void SetWaveform (TeerWaveformEnumNet Waveform)

      sets TEER stimulation waveform (sine/rect)

    TeerClampModeEnumNet GetClampMode ()

      gets TEER clamp mode (voltage/current)

    void SetClampMode (TeerClampModeEnumNet ClampMode)

      sets TEER clamp mode (voltage/current)

    void StartSampling (uint32 t NumberOfCycles)

      starts TEER stimulation (duration: n cycles) and samples during last cycle

    void StopSampling ()

      stops TEER stimulation and sampling

    uint32_t IsSamplingFinished ()

      returns false iff stimulation/sampling is going on, otherwise true

    void SetControllerParams (uint32_t P, uint32_t I, uint32_t D)

      sets PID controller parameters for voltage clamp mode

    void GetControllerParams ([System::Runtime::InteropServices::Out]uint32 t% P, [System::Runtime::←

  InteropServices::Out]uint32_t% I, [System::Runtime::InteropServices::Out]uint32_t% D)
      gets PID controller parameters for voltage clamp mode

    array< int32 t > ^ GetSampleBufferChunk (int Buffer Length)

     private function to query max. 100 bytes of sample buffer; called internally

    array< int32 t > ^ GetSampleVoltageBuffer uV (int Buffer Length)

     returns voltage sample buffer (max. 500 values); unit: uV

    uint32_t GetMaxChunkSize_Byte ()

     private function to be called internally only

    uint32_t GetBytesPerSample ()

     private function to be called internally only
• uint32_t GetNumberOfAvailableSamples ()
     private function to be called internally only

    void SetBufferIndex (uint32 t NewBufferIndex)

      pre-selects sample buffer to be tranferred by GetSampleVoltageBuffer_uV()

    uint32_t GetAdapterCode ()

     gets the adapter code

    uint32 t GetRotaryPositionCode ()

      gets the rotary position code

    void SetExternalLED (uint32_t NewState)

     sets the external LED

    void SetCurrentEnable (bool NewCurrentEnable)

      when disabled, no current will flow through chamber

    bool GetCurrentEnable ()

      when disabled, no current will flow through chamber

    int32 t GetUptimeSeconds ()

      returns time in seconds since device was powered up
· void StartInternalCalibration ()
      starts determination of internal DAC-offset; result is used internally; NON-BLOCKING call

    bool IsInternalCalibrationFinished ()

      queries whether internal calibration has finished

    int GetDacZero ()

     returns DAC-offset (result of internal calibration); use to check for plausibility only

    void CancelInternalCalibration ()

      in case the internal calibration "hangs", this will cancel it

    void SetLiquidResistance (int32_t NewLiquidResistance_Ohm)

      sets the resistance of the liquid in ohms
```

```
    int32_t GetLiquidResistance ()
        gets the resitance of the liquid in ohms
    int GetScaleFactorU1 ()
        returns U1 scale factor times 10^6 (result of internal calibration)
    int GetScaleFactorU2 ()
        returns U2 scale factor times 10^6 (result of internal calibration)
    int GetAdcOffsetU1 ()
        returns ADC offset of U1 channel (result of internal calibration)
    int GetAdcOffsetU2 ()
        returns ADC offset of U2 channel (result of internal calibration)
    int GetFrameErrorCounter ()
        returns number of times (since bootup) sample memory got overwritten
    int GetSampleRate ()
        returns sample rate in Hz
```

Additional Inherited Members

11.115.1 Detailed Description

CTEERFunctionNet is the class to control the TEER device

11.115.2 Constructor & Destructor Documentation

Initializes a new instance of the CTEERFunctionNet class.

```
11.115.2.2 CTEERFunctionNet() [2/2] CTEERFunctionNet (
CMcsUsbNet^ mcsusb )

11.115.2.3 ~CTEERFunctionNet() virtual ~CTEERFunctionNet ( ) [virtual]

11.115.2.4 "!CTEERFunctionNet() !CTEERFunctionNet ( )
```

11.115.3 Member Function Documentation

```
11.115.3.1 CancelInternalCalibration() void CancelInternalCalibration ( )
in case the internal calibration "hangs", this will cancel it
11.115.3.2 GetAdapterCode() uint32_t GetAdapterCode ( )
gets the adapter code
Returns
     the adapter code
11.115.3.3 GetAdcOffsetU1() int GetAdcOffsetU1 ( )
returns ADC offset of U1 channel (result of internal calibration)
Returns
     the ADC offset for U1
11.115.3.4 GetAdcOffsetU2() int GetAdcOffsetU2 ( )
returns ADC offset of U2 channel (result of internal calibration)
Returns
     the ADC offset for U2
11.115.3.5 GetAmplitude_nA() uint32_t GetAmplitude_nA ( )
gets TEER stimulation amplitude in nA
Returns
```

current stimulation amplitude in nA

```
11.115.3.6 GetBytesPerSample() uint32_t GetBytesPerSample ( )
```

private function to be called internally only

Returns

not documented

```
11.115.3.7 GetClampMode() TeerClampModeEnumNet GetClampMode ( )
```

gets TEER clamp mode (voltage/current)

Returns

current TEER clamp mode

11.115.3.8 GetControllerParams() void GetControllerParams (

```
[System::Runtime::InteropServices::Out] uint32_t% P,
[System::Runtime::InteropServices::Out] uint32_t% I,
[System::Runtime::InteropServices::Out] uint32_t% D)
```

gets PID controller parameters for voltage clamp mode

Parameters

Р	the P value
1	the I value
D	the D value

11.115.3.9 GetCurrentEnable() bool GetCurrentEnable ()

when disabled, no current will flow through chamber

Returns

false when disabled, true when enabled

11.115.3.10 GetDacZero() int GetDacZero ()

returns DAC-offset (result of internal calibration); use to check for plausibility only

Returns

the DAC offset

```
11.115.3.11 GetFrameErrorCounter() int GetFrameErrorCounter ( )
returns number of times (since bootup) sample memory got overwritten
Returns
     the number of errors
11.115.3.12 GetLiquidResistance() int32_t GetLiquidResistance ( )
gets the resitance of the liquid in ohms
Returns
     the resistance in ohms
11.115.3.13 GetMaxChunkSize_Byte() uint32_t GetMaxChunkSize_Byte ()
private function to be called internally only
Returns
     not documented
11.115.3.14 GetNumberOfAvailableSamples() uint32_t GetNumberOfAvailableSamples ( )
private function to be called internally only
Returns
     not documented
11.115.3.15 GetPeriod_us() uint32_t GetPeriod_us ( )
gets the period of TEER stimulation in us
Returns
     the period in us
11.115.3.16 GetRotaryPositionCode() uint32_t GetRotaryPositionCode ( )
gets the rotary position code
Returns
     the rotary position code
11.115.3.17 GetSampleBufferChunk() array<int32_t> ^ GetSampleBufferChunk (
              int Buffer_Length )
```

private function to query max. 100 bytes of sample buffer; called internally

Buffer_Length	The maximal length of Buffer.
---------------	-------------------------------

Returns

not documented

11.115.3.18 GetSampleRate() int GetSampleRate ()

returns sample rate in Hz

Returns

the sample rate in Hz

11.115.3.19 GetSampleVoltageBuffer_uV() array<int32_t> $^{\land}$ GetSampleVoltageBuffer_uV (int Buffer_Length)

returns voltage sample buffer (max. 500 values); unit: uV

Parameters

Buffer Lenath Th	ne maximal length of Buffer.

Returns

the voltage sample buffer

11.115.3.20 GetScaleFactorU1() int GetScaleFactorU1 ()

returns U1 scale factor times 10^6 (result of internal calibration)

Returns

the U1 scale factor

```
11.115.3.21 GetScaleFactorU2() int GetScaleFactorU2 ( )
returns U2 scale factor times 10<sup>6</sup> (result of internal calibration)
Returns
     the U2 scale factor
11.115.3.22 GetUptimeSeconds() int32_t GetUptimeSeconds ()
returns time in seconds since device was powered up
Returns
     seconds since power-on
11.115.3.23 GetWaveform() TeerWaveformEnumNet GetWaveform ( )
gets TEER stimulation waveform (sine/rect)
Returns
     waveform enum
11.115.3.24 IsInternalCalibrationFinished() bool IsInternalCalibrationFinished ( )
queries whether internal calibration has finished
Returns
     true if calibration has finished
11.115.3.25 IsSamplingFinished() uint32_t IsSamplingFinished ( )
returns false iff stimulation/sampling is going on, otherwise true
Returns
     true if sampling is finished
11.115.3.26 SetAmplitude_nA() void SetAmplitude_nA (
              uint32_t Amplitude_nA )
```

sets TEER stimulation amplitude in nA

Amplitude_nA	new stimulation amplitude in nA
--------------	---------------------------------

```
11.115.3.27 SetBufferIndex() void SetBufferIndex ( uint32_t NewBufferIndex )
```

pre-selects sample buffer to be tranferred by GetSampleVoltageBuffer_uV()

Parameters

NewBufferIndex	0 - chamber voltage; 1 - compliance voltage
----------------	---

```
11.115.3.28 SetClampMode() void SetClampMode (
TeerClampModeEnumNet ClampMode)
```

sets TEER clamp mode (voltage/current)

Parameters

ClampMode	new TEER clamp mode
-----------	---------------------

sets PID controller parameters for voltage clamp mode

Parameters

Р	the P value
1	the I value
D	the D value

```
11.115.3.30 SetCurrentEnable() void SetCurrentEnable ( bool NewCurrentEnable )
```

when disabled, no current will flow through chamber

NewCurrentEnable disabled when false, enabled when true

```
11.115.3.31 SetExternalLED() void SetExternalLED (
uint32_t NewState)
```

sets the external LED

Parameters

NewState state

11.115.3.32 SetLiquidResistance() void SetLiquidResistance (int32_t NewLiquidResistance_Ohm)

sets the resistance of the liquid in ohms

Parameters

```
11.115.3.33 SetPeriod_us() void SetPeriod_us ( uint32_t period_us )
```

sets the period of TEER stimulation in us

Parameters

period_us the period in us

```
11.115.3.34 SetWaveform() void SetWaveform (
TeerWaveformEnumNet Waveform)
```

sets TEER stimulation waveform (sine/rect)

Parameters

Waveform waveform enum

11.115.3.35 StartInternalCalibration() void StartInternalCalibration ()

starts determination of internal DAC-offset; result is used internally; NON-BLOCKING call

```
11.115.3.36 StartSampling() void StartSampling (
uint32_t NumberOfCycles)
```

starts TEER stimulation (duration: n cycles) and samples during last cycle

Parameters


NumberOfCycles | number of cycles (sine or rect) to output (0 - loop forever)

11.115.3.37 StopSampling() void StopSampling ()

stops TEER stimulation and sampling

11.116 CTEERMachineDeviceNet Class Reference

Inheritance diagram for CTEERMachineDeviceNet:



Public Member Functions

- CTEERMachineDeviceNet ()
- ∼CTEERMachineDeviceNet ()

Properties

• CTEERFunctionNet [get]

Additional Inherited Members

11.116.1 Constructor & Destructor Documentation

```
11.116.1.1 CTEERMachineDeviceNet() CTEERMachineDeviceNet ( )
```

```
11.116.1.2 ~CTEERMachineDeviceNet() ~CTEERMachineDeviceNet ()
```


11.116.2 Property Documentation

```
11.116.2.1 TEERFunctionNet CTEERFunctionNet^ TEERFunctionNet [get]
```

11.117 CUsbDeviceConfigurationFunctionNet Class Reference

CUsbDeviceConfigurationFunctionNet is the class to configure the USB firmware

Inheritance diagram for CUsbDeviceConfigurationFunctionNet:



Public Member Functions

• CUsbDeviceConfigurationFunctionNet (CMcsUsbNet[∧] mcsusb, CMcsUsbFunctionPointerContainer[∧] p↔ UsbDeviceConfigurationFunctionPointerContainer)

Initializes a new instance of the CUsbDeviceConfigurationFunctionNet class.

- CUsbDeviceConfigurationFunctionNet (CMcsUsbNet[∧] mcsusb)
- virtual ~CUsbDeviceConfigurationFunctionNet ()
- !CUsbDeviceConfigurationFunctionNet ()
- void SetDeviceName (String^ name)

sets the USB device name for configurable devices

void SetDeviceId (ProductIdEnumNet id)

sets the USB device name for configurable devices

Additional Inherited Members

11.117.1 Detailed Description

CUsbDeviceConfigurationFunctionNet is the class to configure the USB firmware

11.117.2 Constructor & Destructor Documentation

```
11.117.2.1 CUsbDeviceConfigurationFunctionNet() [1/2] CUsbDeviceConfigurationFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pUsbDeviceConfigurationFunctionPointerContainer)
```

 $Initializes\ a\ new\ instance\ of\ the\ CUsbDeviceConfigurationFunctionNet\ class.$

```
11.117.2.2 CUsbDeviceConfigurationFunctionNet() [2/2] CUsbDeviceConfigurationFunctionNet (
CMcsUsbNet^ mcsusb )
```

```
11.117.2.3 \sim CUsbDeviceConfigurationFunctionNet() virtual \sim CUsbDeviceConfigurationFunctionNet () [virtual]
```

11.117.2.4 "!CUsbDeviceConfigurationFunctionNet() !CUsbDeviceConfigurationFunctionNet ()

11.117.3 Member Function Documentation

```
11.117.3.1 SetDeviceId() void SetDeviceId ( ProductIdEnumNet id )
```

sets the USB device name for configurable devices

Parameters

id

```
11.117.3.2 SetDeviceName() void SetDeviceName ( String^{\wedge} name )
```

sets the USB device name for configurable devices


Parameters

name

11.118 CUsbExceptionNet Class Reference

Exception class that is thrown in case of an USB error.

Inheritance diagram for CUsbExceptionNet:



Public Member Functions

- CUsbExceptionNet (uint32_t status)
 Constructor of a CUsbException.
- CUsbExceptionNet (uint32_t status, String^ message)

Properties

• uint32_t Status [get]

11.118.1 Detailed Description

Exception class that is thrown in case of an USB error.

11.118.2 Constructor & Destructor Documentation

```
11.118.2.1 CUsbExceptionNet() [1/2] CUsbExceptionNet ( uint32_t status )
```

Constructor of a CUsbException.

status the status number

11.118.3 Property Documentation

```
11.118.3.1 Status uint32_t Status [get]
```

11.119 CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet Class Reference

Public Member Functions

• CVoltageRangeInfoNet (int vr, String^ vrString)

Public Attributes

- int VoltageRangeInMicroVolt
- $\bullet \ \, \textbf{String} \ ^{\wedge} \ \, \textbf{VoltageRangeDisplayStringMilliVolt}$

11.119.1 Constructor & Destructor Documentation


11.119.2 Member Data Documentation

$\textbf{11.119.2.1} \quad \textbf{VoltageRangeDisplayStringMilliVolt} \quad \texttt{String} \; \wedge \; \texttt{VoltageRangeDisplayStringMilliVolt}$

11.119.2.2 VoltageRangeInMicroVolt int VoltageRangeInMicroVolt

11.120 CW2100_FunctionNet Class Reference

Inheritance diagram for CW2100_FunctionNet:



Classes

struct AudioChannelsNet

Public Member Functions

- CW2100_FunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] w2100_Function
 — PointerContainer)
- CW2100_FunctionNet (CMcsUsbNet[^] mcsusb)
- array< HeadStageIDType[^]> [^] GetAvailableHeadstages (unsigned int max_length)
- void SelectHeadstage (unsigned int IDorEntry, int TimeSlotNr)
- void DeselectHeadstage (int TimeSlotNr)
- void DeselectAllHeadstages ()
- HeadStageIDTypeState ^ GetSelectedHeadstageState (int TimeSlotNr)
- BatteryState ^ GetBatteryState (int TimeSlotNr)
- System::String \(^\) GetUserDefinedName (unsigned short ID)
- System::String ^ GetUserDefinedNameFromSelectedHS (int TimeSlotNr)
- System::String \(^\) GetUserDefinedNameCache (unsigned short ID)
- W2100_StimulusParametersNet ^ GetStiumlusParameters (unsigned short ID)
- W2100 StimulusParametersNet ^ GetStimulusParametersFromSelectedHS (int TimeSlotNr)
- W2100 StimulusParametersNet \(^\) GetStimulusParametersCache (unsigned int typeValue)
- uint32_t GetStimulusParametersCache (unsigned int typeValue, [System::Runtime::InteropServices::
 — Out]W2100_StimulusParametersNet^% StimulusParameters)
- void SetSelectedChannels (array< BYTE >^ channels, int TimeSlotNr)
- array< BYTE > ^ GetSelectedChannels (int TimeSlotNr)
- void SetMultiHeadstageMode (bool Mode)
- bool GetMultiHeadstageMode ()
- void SetHeadstageSamplingActive (bool Active, int TimeSlotNr)
- bool GetHeadstageSamplingActive (int TimeSlotNr)
- void SetHeadstageToSleep (unsigned int Sleep16ms, int TimeSlotNr)
- void SetHeadstageOnOff (unsigned short On, int TimeSlotNr)
- unsigned short GetHeadstageOnOff (int TimeSlotNr)
- unsigned int GetAnalogOutChannel ([System::Runtime::InteropServices::Out]int % automatic, unsigned short index)
- void SetAnalogOutChannel (int automatic, unsigned short index, unsigned int Channel)

- array< unsigned int > ^ GetAnalogOutFilter ([System::Runtime::InteropServices::Out]int % automatic)
- void SetAnalogOutFilter (int automatic, array< unsigned int >^ Coeffs)
- AnalogOut_DAC_Range_EnumNet GetDacRange ()
- void SetDacRange (AnalogOut_DAC_Range_EnumNet range)
- CFilterPropertyNet ^ GetFilterProperty (W2100DacqGroupChannelEnumNet GroupID, unsigned int index)
- array < CFilterPropertyNet^> ^ GetFilterProperties (W2100DacqGroupChannelEnumNet GroupID)
- void SetAccelGyroEnabled (W2100_Accel_Gyro_Select_EnumNet enable, int TimeSlotNr)
- W2100 Accel Gyro Select EnumNet GetAccelGyroEnabled (int TimeSlotNr)
- void SetAccelGyroDesiredRate (int rate, int TimeSlotNr)
- int GetAccelGyroDesiredRate (int TimeSlotNr)
- int GetAccelGyroCurrentRate (int TimeSlotNr)
- void SetAccelRange (int range, int TimeSlotNr)
- int GetAccelRange (int TimeSlotNr)
- void SetGyroRange (int range, int TimeSlotNr)
- int GetGyroRange (int TimeSlotNr)
- void SetAudioChannels (array< AudioChannelsNet[^]>[^] channels)
- array< AudioChannelsNet[^]> [^] GetAudioChannels ()
- unsigned int GetPicFirmwareType (int TimeSlotNr)
- unsigned int GetFPGAFirmwareType (int TimeSlotNr)

Static Public Member Functions

- static void ClearUserDefinedNameCache ()
- static void ClearUserDefinedNameCache (unsigned short ID)
- static void ClearStimulusParametersCache ()
- static void ClearStimulusParametersCache (unsigned short ID)

Properties

- CW2100_StimulatorFunctionNet^ Stimulator [get]
- CPulseGeneratorFunctionNet[^] PulseGenerator [get]

Additional Inherited Members

11.120.1 Constructor & Destructor Documentation

11.120.2 Member Function Documentation

```
11.120.2.1 ClearStimulusParametersCache() [1/2] static void ClearStimulusParametersCache ( )
[static]
11.120.2.2 ClearStimulusParametersCache() [2/2] static void ClearStimulusParametersCache (
              unsigned short ID ) [static]
11.120.2.3 ClearUserDefinedNameCache() [1/2] static void ClearUserDefinedNameCache ( ) [static]
\textbf{11.120.2.4} \quad \textbf{ClearUserDefinedNameCache() [2/2]} \quad \texttt{static void ClearUserDefinedNameCache ()} \\
              unsigned short ID ) [static]
11.120.2.5 DeselectAllHeadstages() void DeselectAllHeadstages ()
11.120.2.6 DeselectHeadstage() void DeselectHeadstage (
              int TimeSlotNr )
\textbf{11.120.2.7} \quad \textbf{GetAccelGyroCurrentRate()} \quad \texttt{int GetAccelGyroCurrentRate ()}
              int TimeSlotNr )
11.120.2.8 GetAccelGyroDesiredRate() int GetAccelGyroDesiredRate (
              int TimeSlotNr )
11.120.2.9 GetAccelGyroEnabled() W2100_Accel_Gyro_Select_EnumNet GetAccelGyroEnabled (
              int TimeSlotNr )
```

```
11.120.2.10 GetAccelRange() int GetAccelRange (
              int TimeSlotNr )
11.120.2.11 GetAnalogOutChannel() unsigned int GetAnalogOutChannel (
              [System::Runtime::InteropServices::Out] int % automatic,
              unsigned short index )
\textbf{11.120.2.12} \quad \textbf{GetAnalogOutFilter()} \quad \texttt{array} < \texttt{unsigned int} > \\ ^ \land \quad \texttt{GetAnalogOutFilter} \ \ (
              [System::Runtime::InteropServices::Out] int % automatic )
11.120.2.13 GetAudioChannels() array<AudioChannelsNet^> ^ GetAudioChannels ( )
11.120.2.14 GetAvailableHeadstages() array<HeadStageIDType^> ^ GetAvailableHeadstages (
              unsigned int max_length )
11.120.2.15 GetBatteryState() BatteryState ^ GetBatteryState (
             int TimeSlotNr )
11.120.2.16 GetDacRange() AnalogOut_DAC_Range_EnumNet GetDacRange ( )
11.120.2.17 GetFilterProperties() array<CFilterPropertyNet^> ^ GetFilterProperties (
              W2100DacqGroupChannelEnumNet GroupID )
11.120.2.18 GetFilterProperty() CFilterPropertyNet ^ GetFilterProperty (
              W2100DacqGroupChannelEnumNet GroupID,
              unsigned int index )
11.120.2.19 GetFPGAFirmwareType() unsigned int GetFPGAFirmwareType (
              int TimeSlotNr )
```


```
\textbf{11.120.2.20} \quad \textbf{GetGyroRange()} \quad \texttt{int GetGyroRange ()}
                                                                    int TimeSlotNr )
11.120.2.21 GetHeadstageOnOff() unsigned short GetHeadstageOnOff (
                                                                   int TimeSlotNr )
11.120.2.22 GetHeadstageSamplingActive() bool GetHeadstageSamplingActive (
                                                                    int TimeSlotNr )
11.120.2.23 GetMultiHeadstageMode() bool GetMultiHeadstageMode ( )
11.120.2.24 GetPicFirmwareType() unsigned int GetPicFirmwareType (
                                                                   int TimeSlotNr )
11.120.2.25 GetSelectedChannels() array<br/>
BYTE> ^{\wedge} GetSelectedChannels (
                                                                    int TimeSlotNr )
\textbf{11.120.2.26} \quad \textbf{GetSelectedHeadstageState()} \quad \texttt{HeadStageIDTypeState} \quad ^{\land} \quad \texttt{GetSelectedHeadstageState} \quad \texttt{(}
                                                                    int TimeSlotNr )
11.120.2.27 GetStimulusParametersCache() [1/2] W2100_StimulusParametersNet ^ GetStimulus↔
ParametersCache (
                                                                   unsigned int typeValue )
\textbf{11.120.2.28} \quad \textbf{GetStimulusParametersCache() [2/2]} \quad \texttt{uint32\_t GetStimulusParametersCache ()} \quad \textbf{(2/2)} \quad \textbf{(2/2)
                                                                  unsigned int typeValue,
                                                                    [System::Runtime::InteropServices::Out] W2100_StimulusParametersNet^% Stimulus↔
Parameters )
```

```
11.120.2.29 GetStimulusParametersFromSelectedHS() w2100_StimulusParametersNet ^ GetStimulus↔
ParametersFromSelectedHS (
             int TimeSlotNr )
11.120.2.30 GetStiumlusParameters() W2100_StimulusParametersNet ^ GetStiumlusParameters (
             unsigned short ID )
\textbf{11.120.2.31} \quad \textbf{GetUserDefinedName()} \quad \texttt{System::String} \; \wedge \; \texttt{GetUserDefinedName} \; \; \textbf{(}
             unsigned short ID )
11.120.2.32 GetUserDefinedNameCache() [1/2] System::String ^ GetUserDefinedNameCache (
             unsigned short ID )
11.120.2.33 GetUserDefinedNameCache() [2/2] uint32_t GetUserDefinedNameCache (
             unsigned short ID,
             [System::Runtime::InteropServices::Out] System::String^% Name )
11.120.2.34 GetUserDefinedNameFromSelectedHS() System::String ^ GetUserDefinedNameFrom←
SelectedHS (
             int TimeSlotNr )
11.120.2.35 SelectHeadstage() void SelectHeadstage (
             unsigned int IDorEntry,
             int TimeSlotNr )
11.120.2.36 SetAccelGyroDesiredRate() void SetAccelGyroDesiredRate (
             int rate,
             int TimeSlotNr )
11.120.2.37 SetAccelGyroEnabled() void SetAccelGyroEnabled (
             W2100_Accel_Gyro_Select_EnumNet enable,
             int TimeSlotNr )
```

```
11.120.2.38 SetAccelRange() void SetAccelRange (
             int range,
             int TimeSlotNr )
11.120.2.39 SetAnalogOutChannel() void SetAnalogOutChannel (
            int automatic,
             unsigned short index,
             unsigned int Channel )
11.120.2.40 SetAnalogOutFilter() void SetAnalogOutFilter (
            int automatic,
             array< unsigned int >^{\land} Coeffs )
11.120.2.41 SetAudioChannels() void SetAudioChannels (
             array< AudioChannelsNet^>^ channels )
11.120.2.42 SetDacRange() void SetDacRange (
             AnalogOut_DAC_Range_EnumNet range )
11.120.2.43 SetGyroRange() void SetGyroRange (
             int range,
             int TimeSlotNr )
11.120.2.44 SetHeadstageOnOff() void SetHeadstageOnOff (
             unsigned short On,
             int TimeSlotNr )
11.120.2.45 SetHeadstageSamplingActive() void SetHeadstageSamplingActive (
            bool Active,
             int TimeSlotNr )
```

11.121 CW2100_StimulatorFunctionNet Class Reference

Inheritance diagram for CW2100_StimulatorFunctionNet:



Public Member Functions

- CW2100 StimulatorFunctionNet (CMcsUsbNet[^] mcsusb)
- void SendStart (uint32_t triggermap)

Start (Trigger) the STG. The startup delay is in the range of a few ms.

void SendStop (uint32 t triggermap)

Stop some or all triggers of the STG.

- CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData ^ PrepareData (int channel, array< int32_t > ^ amplitude, array< uint64_t > ^ duration, STG_DestinationEnumNet destType, uint32_t repeat)
- CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData
 PrepareDataSync (int channel, array< int32_t >^ amplitude, array< uint32_t >^ Sync, array< uint64_t >^ duration, STG_DestinationEnumNet destType, uint32_t repeat)
- void SendPreparedData (int channel, CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData[^] device_data_and_unrolled, STG_DestinationEnumNet destType)
- · void ClearChannelData (int channel)

Delete a Stimulus Pattern from STG memory

int GetDACResolution ()

Gets number of bits of the DAC resolution.

int GetTimeResolutionInNanoSeconds ()

Gets number of bits of the DAC resolution.

• int GetVoltageRangeInMicroVolt (uint32_t channel)

Gets the Voltage Range of the specified channel in Microvolts.

int GetVoltageResolutionInMicroVolt (uint32 t channel)

Gets the Voltage Resolution of the specified channel in Microvolts.

int GetCurrentRangeInNanoAmp (uint32_t channel)

Gets the Current Range of the specified channel in Nanoamps.

int GetCurrentResolutionInNanoAmp (uint32_t channel)

Gets the Current Resolution of the specified channel in Nanoamps.

- uint32 t GetNumberOfAnalogChannels ()
- uint32 t GetNumberOfSyncoutChannels ()
- uint32 t GetNumberOfTriggerInputs ()
- void SelectTimeSlot (int TimeSlotNr)
- int GetTimeSlot ()
- uint32_t GetStimulationPatternMemory ()
- uint32 t GetBoostPreTime ()
- uint32_t GetBoostAlwaysOnMode ()
- void SetDigitalStimulatorTrigger (int TimeSlotNr, DigitalStimulatorTriggerEventEnumNet trigger_event, int trigger_number, W2100DigitalSourceEnumNet digstream_source, int bitnumber_offset)
- void GetDigitalStimulatorTrigger (int TimeSlotNr, DigitalStimulatorTriggerEventEnumNet trigger_event, int trigger_number, [System::Runtime::InteropServices::Out]W2100DigitalSourceEnumNet% digstream_source, [System::Runtime::InteropServices::Out]int% bitnumber_offset)
- void SetDigitalStimulatorTriggerSlope (int TimeSlotNr, DigitalStimulatorTriggerEventEnumNet trigger_event, int trigger_number, DigitalStimulatorTriggerSlopeEnumNet slope)
- DigitalStimulatorTriggerSlopeEnumNet GetDigitalStimulatorTriggerSlope (int TimeSlotNr, DigitalStimulator
 — TriggerEventEnumNet trigger_event, int trigger_number)
- void StartPoll ()
- void StopPoll ()

Static Public Attributes

- static const uint32_t BOOST_BIT = (1 << 0)
- static const uint32_t GND_SWITCH_BIT = (1 << 1)
- static const uint32_t SYNC_BIT0 = (1 << 2)
- static const uint32_t SYNC_BIT1 = (1 << 3)

Events

OnStgPollStatus[^] PollStatusEvent

Additional Inherited Members

11.121.1 Constructor & Destructor Documentation

```
11.121.1.1 CW2100_StimulatorFunctionNet() CW2100_StimulatorFunctionNet (
CMcsUsbNet^ mcsusb )
```

11.121.2 Member Function Documentation

Delete a Stimulus Pattern from STG memory

Parameters

channel specifies the channel to clear.

```
11.121.2.2 GetBoostAlwaysOnMode() uint32_t GetBoostAlwaysOnMode ( )
```

```
11.121.2.3 GetBoostPreTime() uint32_t GetBoostPreTime ( )
```

Gets the Current Range of the specified channel in Nanoamps.

Parameters

channel Channel which is queried.

Returns

The Current Range of the specified channel in Nanoamps.

11.121.2.5 GetCurrentResolutionInNanoAmp() int GetCurrentResolutionInNanoAmp (uint32_t channel)

Gets the Current Resolution of the specified channel in Nanoamps.

Parameters

channel	Channel which is queried.
---------	---------------------------

Returns

The Current Resolution of the specified channel in Nanoamps.

11.121.2.6 GetDACResolution() int GetDACResolution ()

Gets number of bits of the DAC resolution.

Returns

The DAC resolution in bits.

```
11.121.2.7 GetDigitalStimulatorTrigger() void GetDigitalStimulatorTrigger (
```

```
int TimeSlotNr,
DigitalStimulatorTriggerEventEnumNet trigger_event,
int trigger_number,
[System::Runtime::InteropServices::Out] W2100DigitalSourceEnumNet% digstream_
source,
[System::Runtime::InteropServices::Out] int% bitnumber_offset )
```

```
\textbf{11.121.2.8} \quad \textbf{GetDigitalStimulatorTriggerSlope()} \quad \texttt{DigitalStimulatorTriggerSlopeEnumNet} \quad \texttt{GetDigital} \leftarrow \\
```

```
\textbf{11.121.2.9} \quad \textbf{GetNumberOfAnalogChannels()} \quad \texttt{uint32\_t GetNumberOfAnalogChannels ()} \\
\textbf{11.121.2.10} \quad \textbf{GetNumberOfSyncoutChannels()} \quad \texttt{uint32\_t GetNumberOfSyncoutChannels ()}
11.121.2.11 GetNumberOfTriggerInputs() uint32_t GetNumberOfTriggerInputs ( )
11.121.2.12 GetStimulationPatternMemory() uint32_t GetStimulationPatternMemory ( )
11.121.2.13 GetTimeResolutionInNanoSeconds() int GetTimeResolutionInNanoSeconds ( )
Gets number of bits of the DAC resolution.
Returns
     The time resolution in ns.
11.121.2.14 GetTimeSlot() int GetTimeSlot ( )
11.121.2.15 GetVoltageRangeInMicroVolt() int GetVoltageRangeInMicroVolt (
               uint32_t channel )
Gets the Voltage Range of the specified channel in Microvolts.
Parameters
 channel
            Channel which is queried.
```

Generated by Doxygen

The Voltage Range of the specified channel in Microvolts.

Returns

11.121.2.16 GetVoltageResolutionInMicroVolt() int GetVoltageResolutionInMicroVolt (uint32_t *channel*)

Gets the Voltage Resolution of the specified channel in Microvolts.

channel	Channel which is queried.
---------	---------------------------

Returns

The Voltage Resolution of the specified channel in Microvolts.

```
11.121.2.17 PrepareData() CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData ^ Prepare↔
Data (
              int channel,
              array< int32_t >^{\land} amplitude,
              array< uint64_t >^{\land} duration,
              STG_DestinationEnumNet destType,
              uint32_t repeat )
11.121.2.18 PrepareDataSync() CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData ^ Prepare←
DataSync (
             int channel,
             array< int32_t >^{\wedge} amplitude,
              array< uint32_t >^{\wedge} Sync,
              array< uint64_t >^{\wedge} duration,
              STG_DestinationEnumNet destType,
              uint32_t repeat )
11.121.2.19 SelectTimeSlot() void SelectTimeSlot (
              int TimeSlotNr )
11.121.2.20 SendPreparedData() void SendPreparedData (
             int channel,
              CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData^ device_data_and_unrolled,
              {\tt STG\_DestinationEnumNet}\ \textit{destType}\ )
11.121.2.21 SendStart() void SendStart (
             uint32_t triggermap )
```

Generated by Doxygen

Start (Trigger) the STG. The startup delay is in the range of a few ms.

triggermap	A bitmap of triggers which will be started.

```
11.121.2.22 SendStop() void SendStop ( uint32_t triggermap)
```

Stop some or all triggers of the STG.

Parameters

tr	iggermap	A bitmap of triggers which will be stopped.
----	----------	---

```
11.121.2.23 SetDigitalStimulatorTrigger() void SetDigitalStimulatorTrigger (
```

```
int TimeSlotNr,
DigitalStimulatorTriggerEventEnumNet trigger_event,
int trigger_number,
W2100DigitalSourceEnumNet digstream_source,
int bitnumber_offset )
```

11.121.2.24 SetDigitalStimulatorTriggerSlope() void SetDigitalStimulatorTriggerSlope (

```
int TimeSlotNr,
DigitalStimulatorTriggerEventEnumNet trigger_event,
int trigger_number,
DigitalStimulatorTriggerSlopeEnumNet slope )
```

```
11.121.2.25 StartPoll() void StartPoll ()
```

```
11.121.2.26 StopPoll() void StopPoll ()
```

11.121.3 Member Data Documentation

```
11.121.3.1 BOOST_BIT const uint32_t BOOST_BIT = (1 << 0) [static]
```

```
\textbf{11.121.3.2} \quad \textbf{GND\_SWITCH\_BIT} \quad \texttt{const uint32\_t GND\_SWITCH\_BIT} = (1 << 1) \quad \texttt{[static]}
```

```
11.121.3.3 SYNC_BIT0 const uint32_t SYNC_BIT0 = (1 << 2) [static]
```

```
11.121.3.4 SYNC_BIT1 const uint32_t SYNC_BIT1 = (1 << 3) [static]
```

11.121.4 Event Documentation

11.121.4.1 PollStatusEvent OnStgPollStatus^ PollStatusEvent

11.122 CW2100DacqGroupChannelSelectionNet Class Reference

Inheritance diagram for CW2100DacqGroupChannelSelectionNet:



Public Member Functions

CW2100DacqGroupChannelSelectionNet (CMcsUsbNet[^] mcsusb)

Additional Inherited Members


11.122.1 Constructor & Destructor Documentation

11.122.1.1 CW2100DacqGroupChannelSelectionNet() CW2100DacqGroupChannelSelectionNet (CMcsUsbNet^ mcsusb)

11.123 CWarnerUssingDeviceNet Class Reference

CWarnerUssingDeviceNet is the class to control the Ussing device

Inheritance diagram for CWarnerUssingDeviceNet:



Public Member Functions

- CWarnerUssingDeviceNet ()
 - Initializes a new instance of the CWarnerUssingDeviceNet class.
- virtual ~CWarnerUssingDeviceNet ()
- !CWarnerUssingDeviceNet ()

Properties

CWarnerUssingFunctionNet[^] WarnerUssingFunction [get]

Additional Inherited Members

11.123.1 Detailed Description

CWarnerUssingDeviceNet is the class to control the Ussing device

11.123.2 Constructor & Destructor Documentation

11.123.2.1 CWarnerUssingDeviceNet() CWarnerUssingDeviceNet ()

Initializes a new instance of the CWarnerUssingDeviceNet class.

$\textbf{11.123.2.2} \quad \sim \textbf{CWarnerUssingDeviceNet()} \quad \text{virtual} \quad \sim \texttt{CWarnerUssingDeviceNet} \quad \textbf{()} \quad \text{[virtual]}$

11.123.2.3 "!CWarnerUssingDeviceNet() !CWarnerUssingDeviceNet ()


11.123.3 Property Documentation

11.123.3.1 WarnerUssingFunction CWarnerUssingFunctionNet^ WarnerUssingFunction [get]

11.124 CWarnerUssingFunctionNet Class Reference

CWarnerUssingFunctionNet is the class to control the Ussing device

Inheritance diagram for CWarnerUssingFunctionNet:



Public Member Functions

CWarnerUssingFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pWarner
 — UssingFunctionPointerContainer)

Initializes a new instance of the CWarnerUssingFunctionNet class.

- CWarnerUssingFunctionNet (CMcsUsbNet[∧] mcsusb)
- virtual ~CWarnerUssingFunctionNet ()
- !CWarnerUssingFunctionNet ()
- int32_t GetChannelsCountOfChamber (int32_t ChamberId)

gets number of channels in datastream from chamber amp with given index

int32_t GetNumberOfHardwareSlotsForChambers ()

gets number of physical hardware slots for chambers amps

int32 t GetNumberOfAvailableChambers ()

gets number of actually connected chamber amps

bool IsChamberAvailable (int32_t ChamberId)

checks whether chamber amp is connected to slot

void SetPulse (int32_t ChamberId, UssingClampModeEnumNet StgMode, int32_t NumberOfRepetitions, array< int >^ Amplitudes, array< int >^ Durations, array< int >^ PulseMarker)

defines stimulation pulse pattern for voltage or current stimulation; CAUTION: zero-length amplitude will be briefly applied -> choose matching to neighbour to avoid spikes

void SetVoltageClampControllerParam_P (int32_t ChamberId, uint32_t P)

sets P value of PID controller;

• void SetVoltageClampControllerParam_I (int32_t ChamberId, uint32_t I)

sets I value of PID controller;

void SetVoltageClampControllerParam_D (int32_t ChamberId, uint32_t D)

sets D value of PID controller;

• uint32_t GetVoltageClampControllerParam_P (int32_t ChamberId)

gets P value of PID controller;

```
    uint32_t GetVoltageClampControllerParam_I (int32_t ChamberId)

     gets I value of PID controller;

    uint32 t GetVoltageClampControllerParam D (int32 t Chamberld)

     gets D value of PID controller;

    void SetClampMode (int32 t Chamberld, UssingClampModeEnumNet NewClampMode)

     sets clamp mode (voltage, current or open clamp)

    UssingClampModeEnumNet GetClampMode (int32 t ChamberId)

     gets clamp mode (voltage, current or open clamp; do not use when device is in internal calibration mode)

    bool IsInternalCalibrationFinished (int32 t ChamberId)

      when internal calibration is finished, values for U1,2_offset and U1,2_reference and DAC_offset are available

    int32 t GetU1Offset (int32 t ChamberId)

    int32_t GetU2Offset (int32_t ChamberId)

    int32 t GetU1Reference (int32 t ChamberId)

    int32 t GetU2Reference (int32 t ChamberId)

    int32_t GetDacZero (int32_t ChamberId)

    void SetHighCurrentMode (int32 t ChamberId)

     switch to high-current mode

    void SetLowCurrentMode (int32_t ChamberId)

     switch to low-current mode

    bool IsHighCurrentMode (int32 t ChamberId)

    uint32 t GetLowCurrentRange (int32 t Chamberld)

      query the range of the low current mode

    uint32_t GetHighCurrentRange (int32_t ChamberId)

     query the range of the high current mode

    uint32 t GetDacPampsPerDigitLowCurrentRange (int32 t ChamberId)

     get the resolution of the low current mode

    uint32 t GetDacPampsPerDigitHighCurrentRange (int32 t ChamberId)

     get the resolution of the high current mode

    uint32 t GetUnitsPerDigit (int32 t Chamberld, int32 t Channelld)

     gets amps/volts per digit for specified chamber and channel

    int32 t GetUnitExponent (int32 t Chamberld, int32 t Channelld)

     gets the unit exponent for specified chamber and channel

    UssingUnitEnumNet GetUnitName (int32_t Chamberld, int32_t Channelld)

     gets the channel's unit name

    String \(^\) GetUnitDescription (int32 t Chamberld, int32 t Channelld)

     gets the description for the unit

    array< int > ^ GetAvailableChambers ()

     returns array with (zero-based) Chamberlds of all available chambers

    int32 t GetUptimeSeconds (int32 t ChamberId)

     gets the uptime in seconds

    void SetIdleModeOffset (int32_t ChamberId, UssingClampModeEnumNet ClampMode, int32_t NewIdle

  Offset)
     sets the offset (voltage or current) that will be applied when clamping is DISABLED

    int32 t GetIdleModeOffset (int32 t ChamberId, UssingClampModeEnumNet ClampMode)

     gets the offset (voltage or current) that will be applied when clamping is DISABLED

    void SetEnablePulse (int32 t ChamberId, UssingClampModeEnumNet ClampMode, bool Enable)

     enable pulse of given chamber and mode (voltage/current clamp) of this chamber
• bool IsPulseEnabled (int32_t Chamberld, UssingClampModeEnumNet ClampMode)
     returns true when pulse of given chamber and current mode (voltage/current clamp) of this chamber is ENABLED

    void SetLiquidResistance (int32_t ChamberId, int32_t NewLiquidResistance_Ohm)
```

sets the resistance of the liquid

• int32_t GetLiquidResistance (int32_t ChamberId)

gets the resistance of the liquid

• int32_t GetComplianceVoltageThreshold (int32_t ChamberId)

returns compliance voltage threshold in uV; when Uc is above, current source is overloaded

bool CompensateElectrodeOffset (int32_t Chamberld)

blocking call to compensate electrode offset of one chamber; returns true when successful

bool WaitForChamber (int32_t ChamberId)

blocking call that waits for chamber boot-up calibration to complete

bool WaitForAllChambers ()

blocking call that waits for ALL chambers' boot-up calibration to complete

Additional Inherited Members

11.124.1 Detailed Description

CWarnerUssingFunctionNet is the class to control the Ussing device

11.124.2 Constructor & Destructor Documentation

```
11.124.2.1 CWarnerUssingFunctionNet() [1/2] CWarnerUssingFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ pWarnerUssingFunctionPointerContainer)
```

Initializes a new instance of the CWarnerUssingFunctionNet class.

```
11.124.2.2 CWarnerUssingFunctionNet() [2/2] CWarnerUssingFunctionNet (
CMcsUsbNet^ mcsusb )
```

```
11.124.2.3 ~CWarnerUssingFunctionNet() virtual ~CWarnerUssingFunctionNet () [virtual]
```

```
11.124.2.4 "!CWarnerUssingFunctionNet() !CWarnerUssingFunctionNet ( )
```

11.124.3 Member Function Documentation

```
11.124.3.1 CompensateElectrodeOffset() bool CompensateElectrodeOffset ( int32_t ChamberId )
```

blocking call to compensate electrode offset of one chamber; returns true when successful

Chamber⊷	index of hardware chamber slot (zero-based)	
ld		

Returns

true if compensation succeeded

11.124.3.2 GetAvailableChambers() array<int> $^{\land}$ GetAvailableChambers ()

returns array with (zero-based) Chamberlds of all available chambers

gets number of channels in datastream from chamber amp with given index

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

return value of zero means that amp is not placed

11.124.3.4 GetClampMode() UssingClampModeEnumNet GetClampMode (int32_t ChamberId)

gets clamp mode (voltage, current or open clamp; do not use when device is in internal calibration mode)

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the current clamp mode

11.124.3.5 GetComplianceVoltageThreshold() int32_t GetComplianceVoltageThreshold (int32_t ChamberId)

returns compliance voltage threshold in uV; when Uc is above, current source is overloaded

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the compliance voltage threshold in uV

get the resolution of the high current mode

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)	
ld		

Returns

unit: pA/digit in high current mode

11.124.3.7 GetDacPampsPerDigitLowCurrentRange() uint32_t GetDacPampsPerDigitLowCurrentRange (int32_t ChamberId)

get the resolution of the low current mode

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

pA/digit in low current mode

 diagnostic function only -; gets real zero value of DAC in digits (0 -> neg. current; 32767 -> near zero; 65535 -> pos. current)

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the zero value of the DAC

query the range of the high current mode

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

low current range in nA

gets the offset (voltage or current) that will be applied when clamping is DISABLED

Parameters

Chamberld	index of hardware chamber slot (zero-based)
ClampMode	voltage or current clamp stimulation

Returns

unit: nA or uV

$\textbf{11.124.3.11} \quad \textbf{GetLiquidResistance()} \quad \texttt{int32_t} \;\; \texttt{GetLiquidResistance} \;\; \textbf{(}$

int32_t ChamberId)

gets the resistance of the liquid

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the liquid restistance in ohm

query the range of the low current mode

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

low current range in nA

11.124.3.13 GetNumberOfAvailableChambers() int32_t GetNumberOfAvailableChambers ()

gets number of actually connected chamber amps

Returns

the number of actually connected chambers

11.124.3.14 GetNumberOfHardwareSlotsForChambers() int32_t GetNumberOfHardwareSlotsFor← Chambers ()

gets number of physical hardware slots for chambers amps

Returns

the number of hardware chamber slots on the backplane

· diagnostic function only -

Chamber←	index of hardware chamber slot (zero-based)
ld	

Returns

U1 offset

· diagnostic function only -

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

U1 reference

· diagnostic function only -

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

U2 offset

· diagnostic function only -

Chamber←	index of hardware chamber slot (zero-based)
ld	

Returns

U2 reference

gets the description for the unit

Parameters

Chamber← Id	index of hardware chamber slot (zero-based)
Channelld	index of channel (zero-based)

Returns

the description of the unix

gets the unit exponent for specified chamber and channel

Parameters

Chamber← Id	index of hardware chamber slot (zero-based)
Channelld	index of channel (zero-based)

Returns

example: return value -9 means that amps/volts per digit is in nano

gets the channel's unit name

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	
Channelld	index of channel (zero-based)

Returns

the name of the unit

gets amps/volts per digit for specified chamber and channel

Parameters

Chamber← Id	index of hardware chamber slot (zero-based)
Channelld	index of channel (zero-based)

Returns

amps/volts per digit

```
11.124.3.23 GetUptimeSeconds() int32_t GetUptimeSeconds ( int32_t ChamberId )
```

gets the uptime in seconds

Parameters

Chamber←	index of hardware chamber slot (zero-based)
ld	

Returns

seconds since power-on

11.124.3.24 GetVoltageClampControllerParam_D() uint32_t GetVoltageClampControllerParam_D (int32_t ChamberId)

gets D value of PID controller;

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the D value

11.124.3.25 GetVoltageClampControllerParam_I() uint32_t GetVoltageClampControllerParam_I (int32_t ChamberId)

gets I value of PID controller;

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the I value

11.124.3.26 GetVoltageClampControllerParam_P() uint32_t GetVoltageClampControllerParam_P (int32_t ChamberId)

gets P value of PID controller;

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

the P value

11.124.3.27 IsChamberAvailable() bool IsChamberAvailable (int32_t *ChamberId*)

checks whether chamber amp is connected to slot

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)]
ld		

Returns

true if the chamber is available

11.124.3.28 IsHighCurrentMode() bool IsHighCurrentMode (int32_t *ChamberId*)

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

true if in hight current mode

11.124.3.29 IsInternalCalibrationFinished() bool IsInternalCalibrationFinished (int32_t ChamberId)

when internal calibration is finished, values for U1,2_offset and U1,2_reference and DAC_offset are available

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

Returns

true if finished

returns true when pulse of given chamber and current mode (voltage/current clamp) of this chamber is ENABLED

Chamberld	index of hardware chamber slot (zero-based)
ClampMode	voltage or current clamp stimulation

Returns

when ENABLED, previously defined pulse pattern will be applied, otherwise the chamber current/voltage will be kept at specified offset level

sets clamp mode (voltage, current or open clamp)

Parameters

Chamberld	index of hardware chamber slot (zero-based)
NewClampMode	the clamp mode to use

enable pulse of given chamber and mode (voltage/current clamp) of this chamber

Parameters

Chamberld	index of hardware chamber slot (zero-based)
ClampMode	voltage or current clamp stimulation
Enable	when ENABLED, previously defined pulse pattern will be applied, otherwise the chamber current/voltage will be kept at specified offset level

```
11.124.3.33 SetHighCurrentMode() void SetHighCurrentMode ( int32_t ChamberId )
```

switch to high-current mode

Chamber⊷	index of hardware chamber slot (zero-based)]
ld		

sets the offset (voltage or current) that will be applied when clamping is DISABLED

Parameters

Chamberld	index of hardware chamber slot (zero-based)	
ClampMode	voltage or current clamp stimulation	
NewIdleOffset	unit: nA or uV	

sets the resistance of the liquid

Parameters

Chamberld	index of hardware chamber slot (zero-based)
NewLiquidResistance_Ohm	the liquid resistiance in ohm

```
11.124.3.36 SetLowCurrentMode() void SetLowCurrentMode (
    int32_t ChamberId )
```

switch to low-current mode

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	

```
11.124.3.37 SetPulse() void SetPulse ( int32_t ChamberId,
```

```
UssingClampModeEnumNet StgMode, int32_t NumberOfRepetitions, array< int >^ Amplitudes, array< int >^ Durations, array< int >^ PulseMarker)
```

defines stimulation pulse pattern for voltage or current stimulation; CAUTION: zero-length amplitude will be briefly applied -> choose matching to neighbour to avoid spikes

Parameters

Chamberld	index of hardware chamber slot (zero-based); send pattern to connected amp
StgMode	voltage or current clamp stimulation
NumberOfRepetitions	number of repetitions for pulse pattern (-1 for infinite; n means pattern is applied n+1 times)
Amplitudes	amplitude; unit in voltage clamp: uV; unit in current clamp: nA
Durations	duration in 100us; CAUTION: first element is applied only one; auto-loop back to second element after last one
PulseMarker	defines values on digital channel for each step (positive: digital channel "01", neg: "10", zero: "00")

11.124.3.38 SetVoltageClampControllerParam_D() void SetVoltageClampControllerParam_D (int32_t ChamberId, uint32_t D)

sets D value of PID controller;

Parameters

Chamber⊷ Id	index of hardware chamber slot (zero-based)
D	useful range: 0700

```
11.124.3.39 SetVoltageClampControllerParam_I() void SetVoltageClampControllerParam_I ( int32_t ChamberId, uint32_t I)
```

sets I value of PID controller;

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	
1	useful range: 80000120000

$\textbf{11.124.3.40} \quad \textbf{SetVoltageClampControllerParam_P()} \quad \texttt{void SetVoltageClampControllerParam_P} \quad \textbf{(} \quad \textbf{(} \quad \textbf{)} \quad \textbf{(} \quad \textbf{)} \quad \textbf{(} \quad \textbf{(} \quad \textbf{)} \quad \textbf{(} \quad \textbf{)} \quad \textbf{(} \quad \textbf{(}$

```
int32_t ChamberId,
uint32_t P )
```

sets P value of PID controller;

Parameters

Chamber← Id	index of hardware chamber slot (zero-based)
Р	useful value: 130000

11.124.3.41 WaitForAllChambers() bool WaitForAllChambers ()

blocking call that waits for ALL chambers' boot-up calibration to complete

Returns

returns false when at least one chamber's calibration fails (e.g. timeout...)

blocking call that waits for chamber boot-up calibration to complete

Parameters

Chamber⊷	index of hardware chamber slot (zero-based)
ld	


Returns

returns false when calibration fails (e.g. timeout...)

11.125 CWarnerValveControllerDeviceNet Class Reference

CWarnerValveControllerDeviceNet is the class to access the Warner Valve Controller

Inheritance diagram for CWarnerValveControllerDeviceNet:



Public Member Functions

- delegate void OnGetValveActive (uint16_t valve, int valveActive)
- delegate void OnGetValveManualState (uint16_t valve, int32_t valveManualState)
- delegate void OnGetValveManualGroup (uint16_t valve, int32_t valveManualGroup)
- delegate void OnGetValveMode (uint16_t valve, WvcValveModeEnumNet ValveMode)
- delegate void OnGetAnalogThresholdLow (uint16_t valve, int32_t threshold)
- delegate void OnGetAnalogThresholdHigh (uint16 t valve, int32 t threshold)
- delegate void OnGetDigitalPortDirection (uint16 t port, PortDirectionEnumNet direction)
- delegate void OnlsValveDigitalInInverted (uint16_t valve, bool isInverted)
- delegate void OnGetValveDigitalInPort (uint16 t valve, uint32 t digitalInPort)
- delegate void OnlsDigitalOutPortInverted (uint16_t digitalOutPort, bool isInverted)
- delegate void OnGetDigitalOutPortValve (uint16_t digitalOutPort, uint32_t valve)
- delegate void OnlsValveOpen (uint16_t valve, bool valveOpen)
- delegate void OnlsValveOpenInDigitalMode (uint16_t valve, bool valveOpen)
- delegate void OnlsValveOpenInAnalogMode (uint16_t valve, bool valveOpen)
- delegate void OnGetAnalogVoltage (int32_t voltage)
- delegate void OnTableEntryChanged (uint16_t tableNumber)
- delegate void OnGetTableNamebyIndex (uint16 t tableNumber, String[^] tableName)
- delegate void OnGetActiveRunningTableNumber (uint32_t tableNumber)
- delegate void OnGetCurrentNumberOfValves (int32 t numberOfValves)
- delegate void OnGetValveBoardRevision (uint32_t revision)
- delegate void OnGetValveLedOn (bool ledon)
- delegate void OnGetDisplayMode (WvcDisplayModeEnumNet DisplayMode)
- CWarnerValveControllerDeviceNet ()

Initializes a new instance of the CWarnerValveControllerDeviceNet class.

- virtual ~CWarnerValveControllerDeviceNet ()
- !CWarnerValveControllerDeviceNet ()
- int GetValveActive (uint16_t valve)

Gets the valve active/inactive state

void SetValveActive (uint16 t valve, int valveActive)

Sets the valve active/inactive state

uint32_t GetValvesActiveMap ()

Gets the valves active/inactive states

void SetValvesActiveMap (uint32 t valvesActive)

Sets the valve active/inactive state

• int32 t GetValveManualState (uint16 t valve)

Gets the valve manual on/off state

• void SetValveManualState (uint16_t valve, int32_t valveManualState)

Sets the valve manual on/off state

uint32 t GetValvesManualStateMap ()

Gets the valves manual on/off states

void SetValvesManualStateMap (uint32_t valveaManualState)

Sets the valve manual on/off state

int32_t GetValveManualGroup (uint16_t valve)

Gets the valve manual group

void SetValveManualGroup (uint16 t valve, int32 t valveManualGroup)

Sets the valve manual group

WvcValveModeEnumNet GetValveMode (uint16_t valve)

Reads the valve mode

void SetValveMode (uint16_t valve, WvcValveModeEnumNet ValveMode)

Writes the valve mode

int32_t GetAnalogThresholdLow (uint16_t valve)

Gets the lower threshold for the analog in port per valve

void SetAnalogThresholdLow (uint16 t valve, int32 t threshold)

Sets the lower threshold for the analog in port per valve

int32 t GetAnalogThresholdHigh (uint16 t valve)

Gets the upper threshold for the analog in port per valve

void SetAnalogThresholdHigh (uint16_t valve, int32_t threshold)

Sets the upper threshold for the analog in port per valve

PortDirectionEnumNet GetDigitalPortDirection (uint16 t port)

Gets the direction of a digital port

• void SetDigitalPortDirection (uint16_t port, PortDirectionEnumNet direction)

Sets the direction of a digital port

bool IsValveDigitalInInverted (uint16 t valve)

Is digital in inverted

void SetValveDigitalInInvert (uint16_t valve, bool isInverted)

Invert digital in

uint32 t GetValveDigitalInPort (uint16 t valve)

Gets the number of the digital in port which is mapped to a valve

void SetValveDigitalInPort (uint16 t valve, uint32 t digitalInPort)

Map a digital in port to a valve

bool IsDigitalOutPortInverted (uint16 t digitalOutPort)

Gets the number of the valve which is mapped to a digital out port

void SetDigitalOutPortInvert (uint16 t digitalOutPort, bool isInverted)

Map a valve to a digital out port

uint32 t GetDigitalOutPortValve (uint16 t digitalOutPort)

Gets the number of the valve which is mapped to a digital out port

void SetDigitalOutPortValve (uint16_t digitalOutPort, uint32_t valve)

Map a valve to a digital out port

void SetDefault ()

Sets the settings of the valve controller to default

bool IsValveOpen (uint16 t valve)

Is valve open

bool IsValveOpenInDigitalMode (uint16_t valve)

True, if the valve would be open when the device is in digital mode

bool IsValveOpenInAnalogMode (uint16_t valve)

True, if the valve would be open when the device is in analog mode

int32_t GetAnalogVoltage ()

Reads the voltage on the analog in port

void GetValveTableEntry (uint16_t valve, uint16_t index, [System::Runtime::InteropServices::Out]uint32_t% duration, [System::Runtime::InteropServices::Out]bool% state)

Read an entry from the valve protocol table

• void SetValveTableEntry (uint16_t valve, uint16_t index, uint32_t duration, bool state)

Write an entry to the valve protocol table

void ClearValveTable (uint16_t valve)

Clear the valve protocol table

• void LoadValveTable ()

Load the current table from permanent memory

void StoreValveTable ()

Store the current table in permanent memory

String \(^\) GetTableNamebyIndex (uint16_t tableNumber)

Get the name of a protocol table

String ^ GetTableName ()

Get the name of the current protocol table

void SetTableName (String[^] tableName)

Set the name of the current protocol table

uint32 t GetActiveRunningTableNumber ()

Gets the number of the table that is active for running

void SetActiveRunningTableNumber (uint32_t tableNumber)

Sets the number of the tanle that is active for running

uint32 t GetCurrentEditTableNumber ()

Gets the number of the table that is current for editing

void SetCurrentEditTableNumber (uint32_t tableNumber)

Sets the number of the table that is current for editing

void ClearTableName ()

Clear the name of current protocol table

void SetTableStep (uint16_t valve, int32_t steps)

Skips the table protocol for a valve by steps

void SetTableStepAll (int32 t steps)

Skips the table protocol for all valves by steps

• int32 t GetTotalNumberOfValves ()

Get the total number of valves in the system

int32 t GetTotalNumberOfDigitalPorts ()

Get the total number of digital ports in the system

int32 t GetTotalTableSize ()

Get the total table size in the system

int32_t GetTotalNumberOfTables ()

Get the total number of tables in the system

int32_t GetCurrentNumberOfValves ()

Get the current number of valves connected to the system

• uint32 t GetValveBoardRevision ()

Gets the revision code of the valve board

• bool GetValveLedOn ()

Gets the LED state of the valve board

void SetValveLedOn (bool ledon)

Gets the LED state of the valve board

WvcDisplayModeEnumNet GetDisplayMode ()

Reads the display mode

void SetDisplayMode (WvcDisplayModeEnumNet DisplayMode, int32 t lockTimeMs)

Writes the display mode

String \(^\) GetValveBoardRevisionString ()

Gets the revision name of the valve board

Events

• OnGetValveActive^ GetValveActiveEvent [add, remove, raise]

Event fires when the valve state for the valve number has changed

• OnGetValveManualState GetValveManualStateEvent [add, remove, raise]

Event fires when the manual valve state for the valve number has changed

• OnGetValveManualGroup GetValveManualGroupEvent [add, remove, raise]

Event fires when the manual valve group for the valve number has changed

• OnGetValveMode GetValveModeEvent [add, remove, raise]

Event fires when the valve mode for the valve number has changed

- OnGetAnalogThresholdLow^ GetAnalogThresholdLowEvent [add, remove, raise]

 Event fires when the threshold in mV for the valve number has changed
- $\bullet \ \ On Get Analog Threshold High ^ \land \ Get Analog Threshold High Event \ \ [\texttt{add, remove, raise}]$

Event fires when the threshold in mV for the valve number has changed

OnGetDigitalPortDirection^ GetDigitalPortDirectionEvent [add, remove, raise]

Event fires when the direction for the port number has changed

• OnlsValveDigitalInInverted^ IsValveDigitalInInvertedEvent [add, remove, raise]

Event fires when is inverted for the valve number has changed

• OnGetValveDigitalInPort^ GetValveDigitalInPortEvent [add, remove, raise]

Event fires when the digital in port for the valve number has changed

• OnlsDigitalOutPortInverted^ IsDigitalOutPortInvertedEvent [add, remove, raise]

Event fires when is inverted for the digital out port has changed

OnGetDigitalOutPortValve^ GetDigitalOutPortValveEvent [add, remove, raise]

Event fires when the valve number for the digital out port has changed

• OnlsValveOpen^ IsValveOpenEvent [add, remove, raise]

Event fires when is open for the valve number has changed

OnlsValveOpenInDigitalMode^\lambda IsValveOpenInDigitalModeEvent [add, remove, raise]

Event fires when is open for the valve number has changed

OnlsValveOpenInAnalogMode^ IsValveOpenInAnalogModeEvent [add, remove, raise]

Event fires when is open for the valve number has changed

OnGetAnalogVoltage[^] GetAnalogVoltageEvent [add, remove, raise]

Event fires when the voltage in mV has changed

OnTableEntryChanged^ TableEntryChangedEvent [add, remove, raise]

Event fires when an entry of a table changed

• OnGetTableNamebyIndex^ GetTableNamebyIndexEvent [add, remove, raise]

Event fires when the name of the table for the table number has changed

OnGetActiveRunningTableNumber[^] GetActiveRunningTableNumberEvent [add, remove, raise]

Event fires when the table number has changed

• OnGetCurrentNumberOfValves^ GetCurrentNumberOfValvesEvent [add, remove, raise]

Event fires when the number of valves has changed

• OnGetValveBoardRevision^ GetValveBoardRevisionEvent [add, remove, raise]

Event fires when the revision code has changed

• OnGetValveLedOn^ GetValveLedOnEvent [add, remove, raise]

Event fires when the LED state has changed

OnGetDisplayMode^ GetDisplayModeEvent [add, remove, raise]

Event fires when the display mode has changed

Additional Inherited Members

11.125.1 Detailed Description

CWarnerValveControllerDeviceNet is the class to access the Warner Valve Controller

11.125.2 Constructor & Destructor Documentation

```
11.125.2.1 CWarnerValveControllerDeviceNet() CWarnerValveControllerDeviceNet ()
```

Initializes a new instance of the CWarnerValveControllerDeviceNet class.

```
11.125.2.2 ~CWarnerValveControllerDeviceNet() virtual ~CWarnerValveControllerDeviceNet () [virtual]
```

```
11.125.2.3 "!CWarnerValveControllerDeviceNet() !CWarnerValveControllerDeviceNet ( )
```

11.125.3 Member Function Documentation

```
11.125.3.1 ClearTableName() void ClearTableName ( )
```

Clear the name of current protocol table

```
11.125.3.2 ClearValveTable() void ClearValveTable ( uint16_t valve )
```

Clear the valve protocol table

Parameters

valve The valve number

11.125.3.3 GetActiveRunningTableNumber() uint32_t GetActiveRunningTableNumber ()

Gets the number of the table that is active for running

Returns

The table number

```
11.125.3.4 GetAnalogThresholdHigh() int32_t GetAnalogThresholdHigh ( uint16_t valve)
```

Gets the upper threshold for the analog in port per valve

Returns

The threshold in mV

11.125.3.5 GetAnalogThresholdLow() int32_t GetAnalogThresholdLow (uint16_t valve)

Gets the lower threshold for the analog in port per valve

Parameters

valve The valve number

Returns

The threshold in mV

11.125.3.6 GetAnalogVoltage() int32_t GetAnalogVoltage ()

Reads the voltage on the analog in port

Returns

The voltage in mV

11.125.3.7 GetCurrentEditTableNumber() uint32_t GetCurrentEditTableNumber ()

Gets the number of the table that is current for editing

Returns

The table number

11.125.3.8 GetCurrentNumberOfValves() int32_t GetCurrentNumberOfValves ()

Get the current number of valves connected to the system

Returns

The number of valves

Gets the number of the valve which is mapped to a digital out port

digitalOutPort T	he digital out port
------------------	---------------------

Returns

The valve number

11.125.3.10 **GetDigitalPortDirection()** PortDirectionEnumNet GetDigitalPortDirection (uint16_t port)

Gets the direction of a digital port

Parameters

```
port The port number
```

Returns

the direction

11.125.3.11 GetDisplayMode() WvcDisplayModeEnumNet GetDisplayMode ()

Reads the display mode

Returns

The display mode

11.125.3.12 GetTableName() String $^{\wedge}$ GetTableName ()

Get the name of the current protocol table

Returns

The name of the table

11.125.3.13 GetTableNamebyIndex() String $^{\land}$ GetTableNamebyIndex (uint16_t tableNumber)

Get the name of a protocol table

tableNumber T	he table number
---------------	-----------------

Returns

The name of the table

11.125.3.14 GetTotalNumberOfDigitalPorts() int32_t GetTotalNumberOfDigitalPorts ()

Get the total number of digital ports in the system

Returns

The number of digital ports

11.125.3.15 GetTotalNumberOfTables() int32_t GetTotalNumberOfTables ()

Get the total number of tables in the system

Returns

The number of tables

11.125.3.16 GetTotalNumberOfValves() int32_t GetTotalNumberOfValves ()

Get the total number of valves in the system

Returns

The number of valves

11.125.3.17 GetTotalTableSize() int32_t GetTotalTableSize ()

Get the total table size in the system

Returns

The table size

```
11.125.3.18 GetValveActive() int GetValveActive ( uint16\_t \ valve )
```

Gets the valve active/inactive state

valve	The valve number
-------	------------------

Returns

The valve state

11.125.3.19 GetValveBoardRevision() uint32_t GetValveBoardRevision ()

Gets the revision code of the valve board

Returns

The revision code

$\textbf{11.125.3.20} \quad \textbf{GetValveBoardRevisionString()} \quad \texttt{String} \; \wedge \; \texttt{GetValveBoardRevisionString} \; (\;)$

Gets the revision name of the valve board

Returns

The revision name

11.125.3.21 GetValveDigitalInPort() uint32_t GetValveDigitalInPort (uint16_t valve)

Gets the number of the digital in port which is mapped to a valve

Parameters

valve	The valve number

Returns

The digital in port

11.125.3.22 GetValveLedOn() bool GetValveLedOn ()

Gets the LED state of the valve board

Returns

The LED state

$\textbf{11.125.3.23} \quad \textbf{GetValveManualGroup()} \quad \texttt{int32_t GetValveManualGroup} \quad \textbf{(}$

 $uint16_t\ valve$)

Gets the valve manual group

Parameters

valve The valve no	umber
--------------------	-------

Returns

The manual valve group

11.125.3.24 GetValveManualState() int32_t GetValveManualState (uint16_t valve)

Gets the valve manual on/off state

Parameters

valve	The valve number

Returns

The manual valve state

11.125.3.25 GetValveMode() WvcValveModeEnumNet GetValveMode (uint16_t valve)

Reads the valve mode

Parameters

valve	The valve number

Returns

The valve mode

11.125.3.26 GetValvesActiveMap() uint32_t GetValvesActiveMap ()

Gets the valves active/inactive states

Returns

The valves states

11.125.3.27 GetValvesManualStateMap() uint32_t GetValvesManualStateMap ()

Gets the valves manual on/off states

Returns

The manual valves states

11.125.3.28 GetValveTableEntry() void GetValveTableEntry (

```
uint16_t valve,
uint16_t index,
[System::Runtime::InteropServices::Out] uint32_t% duration,
[System::Runtime::InteropServices::Out] bool% state )
```

Read an entry from the valve protocol table

Parameters

valve	The valve number
index	The index in the table
duration	the duration in ms
state	the state

11.125.3.29 IsDigitalOutPortInverted() bool IsDigitalOutPortInverted (uint16_t digitalOutPort)

Gets the number of the valve which is mapped to a digital out port

Parameters

digitalOutPort	The digital out port

Returns

is inverted

```
11.125.3.30 IsValveDigitalInInverted() bool IsValveDigitalInInverted ( uint16_t valve)
```

Is digital in inverted

Parameters

```
valve The valve number
```

Returns

is inverted

```
11.125.3.31 IsValveOpen() bool IsValveOpen ( uint16_t valve )
```

Is valve open

Parameters

valve The valve number

Returns

is open

11.125.3.32 IsValveOpenInAnalogMode() bool IsValveOpenInAnalogMode (uint16_t *valve*)

True, if the valve would be open when the device is in analog mode

Parameters

```
valve The valve number
```

Returns

is open

```
11.125.3.33 IsValveOpenInDigitalMode() bool IsValveOpenInDigitalMode ( uint16_t valve)
```

True, if the valve would be open when the device is in digital mode

Pa	ra	m	ρi	ŀρ	re

valve The valve number	r
------------------------	---

Returns

is open

```
11.125.3.34 LoadValveTable() void LoadValveTable ( )
```

Load the current table from permanent memory

```
11.125.3.35 OnGetActiveRunningTableNumber() delegate void OnGetActiveRunningTableNumber ( uint32_t tableNumber)
```

```
11.125.3.36 OnGetAnalogThresholdHigh() delegate void OnGetAnalogThresholdHigh ( uint16_t valve, int32_t threshold)
```

```
11.125.3.37 OnGetAnalogThresholdLow() delegate void OnGetAnalogThresholdLow ( uint16_t valve, int32_t threshold)
```

```
11.125.3.38 OnGetAnalogVoltage() delegate void OnGetAnalogVoltage ( int32_t voltage )
```

```
11.125.3.39 OnGetCurrentNumberOfValves() delegate void OnGetCurrentNumberOfValves ( int32_t numberOfValves )
```

```
11.125.3.40 OnGetDigitalOutPortValve() delegate void OnGetDigitalOutPortValve ( uint16_t digitalOutPort, uint32_t valve )
```

```
11.125.3.41 OnGetDigitalPortDirection() delegate void OnGetDigitalPortDirection (
             uint16_t port,
             PortDirectionEnumNet direction )
11.125.3.42 OnGetDisplayMode() delegate void OnGetDisplayMode (
             WvcDisplayModeEnumNet DisplayMode )
11.125.3.43 OnGetTableNamebyIndex() delegate void OnGetTableNamebyIndex (
             uint16_t tableNumber,
             String^{\wedge} tableName)
11.125.3.44 OnGetValveActive() delegate void OnGetValveActive (
             uint16_t valve,
             int valveActive )
11.125.3.45 OnGetValveBoardRevision() delegate void OnGetValveBoardRevision (
             uint32_t revision )
11.125.3.46 OnGetValveDigitalInPort() delegate void OnGetValveDigitalInPort (
             uint16_t valve,
             uint32_t digitalInPort )
11.125.3.47 OnGetValveLedOn() delegate void OnGetValveLedOn (
             bool ledon )
11.125.3.48 OnGetValveManualGroup() delegate void OnGetValveManualGroup (
             uint16_t valve,
             int32_t valveManualGroup )
11.125.3.49 OnGetValveManualState() delegate void OnGetValveManualState (
             uint16_t valve,
             int32_t valveManualState )
```

```
11.125.3.50 OnGetValveMode() delegate void OnGetValveMode (
              uint16_t valve,
              WvcValveModeEnumNet ValveMode )
11.125.3.51 OnlsDigitalOutPortInverted() delegate void OnIsDigitalOutPortInverted (
              uint16_t digitalOutPort,
              bool isInverted )
\textbf{11.125.3.52} \quad \textbf{OnlsValveDigitalInInverted()} \quad \texttt{delegate void OnlsValveDigitalInInverted ()}
              uint16_t valve,
              bool is Inverted )
11.125.3.53 OnlsValveOpen() delegate void OnIsValveOpen (
              uint16_t valve,
              bool valveOpen )
11.125.3.54 OnlsValveOpenInAnalogMode() delegate void OnIsValveOpenInAnalogMode (
              uint16_t valve,
              bool valveOpen )
\textbf{11.125.3.55} \quad \textbf{OnlsValveOpenInDigitalMode()} \quad \texttt{delegate void OnlsValveOpenInDigitalMode} \quad \texttt{(}
              uint16_t valve,
              bool valveOpen )
11.125.3.56 OnTableEntryChanged() delegate void OnTableEntryChanged (
              uint16_t tableNumber )
11.125.3.57 SetActiveRunningTableNumber() void SetActiveRunningTableNumber (
              uint32_t tableNumber )
Sets the number of the tanle that is active for running
Parameters
```

tableNumber

The table number

11.125.3.58 SetAnalogThresholdHigh() void SetAnalogThresholdHigh (uint16_t valve, int32_t threshold)

Sets the upper threshold for the analog in port per valve

Parameters

valve	The valve number
threshold	The threshold in mV

Sets the lower threshold for the analog in port per valve

Parameters

valve	The valve number
threshold	The threshold in mV

11.125.3.60 SetCurrentEditTableNumber() void SetCurrentEditTableNumber (uint32_t tableNumber)

Sets the number of the table that is current for editing

Parameters

tableNumber	The table number
-------------	------------------

11.125.3.61 SetDefault() void SetDefault ()

Sets the settings of the valve controller to default

Map a valve to a digital out port

digitalOutPort	The digital out port
isInverted	True if digital out is to be inverted

11.125.3.63 SetDigitalOutPortValve() void SetDigitalOutPortValve (uint16_t digitalOutPort, uint32_t valve)

Map a valve to a digital out port

Parameters

digitalOutPort	The digital out port
valve	The valve number

Sets the direction of a digital port

Parameters

port	The port number
direction	the direction

Writes the display mode

Parameters

DisplayMode	The display mode
lockTimeMs	Locks the display for ms

```
11.125.3.66 SetTableName() void SetTableName ( String^{ } tableName )
```

Set the name of the current protocol table

Parameters

tableName	The name of the table
-----------	-----------------------

Skips the table protocol for a valve by steps

Parameters

valve	The valve number
steps	Number of steps

Skips the table protocol for all valves by steps

Parameters

steps	Number of steps

Sets the valve active/inactive state

Parameters

valve	The valve number
valveActive	The valve state

Invert digital in

Parameters

valve	The valve number
isInverted	True if digital in is to be inverted

Map a digital in port to a valve

Parameters

valve	The valve number
digitalInPort	The digital in port

Gets the LED state of the valve board

Parameters

ledon The LED state

Sets the valve manual group

Parameters

valve	The valve number
valveManualGroup	The manual valve group

Sets the valve manual on/off state

Parameters

valve	The valve number
valveManualState	The manual valve state

Writes the valve mode

Parameters

valve	The valve number
ValveMode	The valve mode

11.125.3.76 SetValvesActiveMap() void SetValvesActiveMap (uint32_t valvesActive)

Sets the valve active/inactive state

Parameters

valvesActive	The valves states

11.125.3.77 SetValvesManualStateMap() void SetValvesManualStateMap (uint32_t valveaManualState)

Sets the valve manual on/off state

Parameters

valveaManualState The manual valves states

Write an entry to the valve protocol table

Parameters

valve	The valve number
index	The index in the table
duration	the duration in ms
state	the state

11.125.3.79 StoreValveTable() void StoreValveTable ()

Store the current table in permanent memory

11.125.4 Event Documentation

11.125.4.1 GetActiveRunningTableNumberEvent OnGetActiveRunningTableNumber^ GetActiveRunning← TableNumberEvent [add], [remove], [raise]

Event fires when the table number has changed

11.125.4.2 GetAnalogThresholdHighEvent OnGetAnalogThresholdHigh^ GetAnalogThresholdHighEvent [add], [remove], [raise]

Event fires when the threshold in mV for the valve number has changed

11.125.4.3 GetAnalogThresholdLowEvent OnGetAnalogThresholdLow^ GetAnalogThresholdLowEvent [add], [remove], [raise]

Event fires when the threshold in mV for the valve number has changed

11.125.4.4 GetAnalogVoltageEvent OnGetAnalogVoltage^ GetAnalogVoltageEvent [add], [remove], [raise]

Event fires when the voltage in mV has changed

11.125.4.5 GetCurrentNumberOfValvesEvent OnGetCurrentNumberOfValves^ GetCurrentNumberOf← ValvesEvent [add], [remove], [raise]

Event fires when the number of valves has changed

11.125.4.6 GetDigitalOutPortValveEvent OnGetDigitalOutPortValve^ GetDigitalOutPortValveEvent [add], [remove], [raise]

Event fires when the valve number for the digital out port has changed

11.125.4.7 GetDigitalPortDirectionEvent OnGetDigitalPortDirection^ GetDigitalPortDirectionEvent [add], [remove], [raise]

Event fires when the direction for the port number has changed

11.125.4.8 GetDisplayModeEvent OnGetDisplayMode^ GetDisplayModeEvent [add], [remove], [raise]

Event fires when the display mode has changed

11.125.4.9 GetTableNamebyIndexEvent OnGetTableNamebyIndex^ GetTableNamebyIndexEvent [add], [remove], [raise]

Event fires when the name of the table for the table number has changed

11.125.4.10 GetValveActiveEvent OnGetValveActive^ GetValveActiveEvent [add], [remove], [raise]

Event fires when the valve state for the valve number has changed

11.125.4.11 GetValveBoardRevisionEvent OnGetValveBoardRevision^ GetValveBoardRevisionEvent [add], [remove], [raise]

Event fires when the revision code has changed

11.125.4.12 GetValveDigitalInPortEvent OnGetValveDigitalInPort^ GetValveDigitalInPortEvent [add], [remove], [raise]

Event fires when the digital in port for the valve number has changed

11.125.4.13 GetValveLedOnEvent OnGetValveLedOn^ GetValveLedOnEvent [add], [remove], [raise]

Event fires when the LED state has changed

11.125.4.14 GetValveManualGroupEvent OnGetValveManualGroup^ GetValveManualGroupEvent [add], [remove], [raise]

Event fires when the manual valve group for the valve number has changed

11.125.4.15 GetValveManualStateEvent OnGetValveManualState^ GetValveManualStateEvent [add], [remove], [raise]

Event fires when the manual valve state for the valve number has changed

11.125.4.16 GetValveModeEvent onGetValveMode^ GetValveModeEvent [add], [remove], [raise]

Event fires when the valve mode for the valve number has changed

11.125.4.17 IsDigitalOutPortInvertedEvent OnIsDigitalOutPortInverted^ IsDigitalOutPortInverted← Event [add], [remove], [raise]

Event fires when is inverted for the digital out port has changed

11.125.4.18 IsValveDigitalInInvertedEvent OnIsValveDigitalInInverted^ IsValveDigitalInInverted← Event [add], [remove], [raise]

Event fires when is inverted for the valve number has changed

11.125.4.19 IsValveOpenEvent OnIsValveOpen^ IsValveOpenEvent [add], [remove], [raise]

Event fires when is open for the valve number has changed

11.125.4.20 IsValveOpenInAnalogModeEvent OnIsValveOpenInAnalogMode^ IsValveOpenInAnalogMode← Event [add], [remove], [raise]

Event fires when is open for the valve number has changed

11.125.4.21 IsValveOpenInDigitalModeEvent OnIsValveOpenInDigitalMode^ IsValveOpenInDigital← ModeEvent [add], [remove], [raise]

Event fires when is open for the valve number has changed


11.125.4.22 TableEntryChangedEvent OnTableEntryChanged^ TableEntryChangedEvent [add], [remove], [raise]

Event fires when an entry of a table changed

11.126 CWarnerValveControllerDeviceTesterFunctionNet Class Reference

CWarnerValveControllerDeviceTesterFunctionNet is the class to access the functions for the Warner Valve Controller Device Tester

Inheritance diagram for CWarnerValveControllerDeviceTesterFunctionNet:



Public Member Functions

CWarnerValveControllerDeviceTesterFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] pWarnerValveControllerDeviceTesterFunctionPointerContainer)

Initializes a new instance of the CWarnerValveControllerDeviceTesterFunctionNet class.

- CWarnerValveControllerDeviceTesterFunctionNet (CMcsUsbNet[^] mcsusb)
- virtual ~CWarnerValveControllerDeviceTesterFunctionNet ()
- !CWarnerValveControllerDeviceTesterFunctionNet ()
- void SetADC (uint32_t onoff)

Sets the ADC port of the tester

• uint32_t GetSync ()

Gets the output from the sync port

void SetTrigger (uint32_t trigger)

Sets the input to the trigger port

• void SetTriggerSyncDirection (uint32_t direction)

Sets the direction of the trigger/sync test port

· uint32_t GetIO ()

Gets the output from the io ports

void SetIO (uint32_t io)

Sets the input to the io ports

· void SetIODirection (int32_t direction)

Sets the direction of the IO test ports

Additional Inherited Members

11.126.1 Detailed Description

CWarnerValveControllerDeviceTesterFunctionNet is the class to access the functions for the Warner Valve Controller Device Tester

11.126.2 Constructor & Destructor Documentation

Initializes a new instance of the CWarnerValveControllerDeviceTesterFunctionNet class.

```
11.126.2.2 CWarnerValveControllerDeviceTesterFunctionNet() [2/2] CWarnerValveControllerDeviceTesterFunctionNet (

CMcsUsbNet^ mcsusb )
```

```
11.126.2.3 ~CWarnerValveControllerDeviceTesterFunctionNet() virtual ~CWarnerValveControllerDeviceTesterFunctio
( ) [virtual]
```

```
11.126.2.4 "!CWarnerValveControllerDeviceTesterFunctionNet() !CWarnerValveControllerDeviceTesterFunctionNet ( )
```

11.126.3 Member Function Documentation

```
11.126.3.1 GetIO() uint32_t GetIO ()
```

Gets the output from the io ports

Returns

The manual valves states

```
11.126.3.2 GetSync() uint32_t GetSync ()
```

Gets the output from the sync port

Returns

The sync state

```
11.126.3.3 SetADC() void SetADC ( uint32_t onoff )
```

Sets the ADC port of the tester

Parameters

```
onoff The port state
```

```
11.126.3.4 SetIO() void SetIO ( uint32_t io )
```

Sets the input to the io ports

Parameters

io The manual valves states

```
11.126.3.5 SetIODirection() void SetIODirection ( int32_t direction )
```

Sets the direction of the IO test ports

Parameters

```
direction The 16bit direction map: 1=IN 0=OUT
```

Sets the input to the trigger port

trigger	The trigger state
---------	-------------------

11.126.3.7 SetTriggerSyncDirection() void SetTriggerSyncDirection (uint32_t direction)


Sets the direction of the trigger/sync test port

Parameters

direction The direction: 1=IN 0=OUT

11.127 CWClassicFunctionNet Class Reference

Inheritance diagram for CWClassicFunctionNet:



Public Member Functions

- CWClassicFunctionNet (CMcsUsbNet[^] mcsusb, CMcsUsbFunctionPointerContainer[^] wClassicFuntion←
 PointerContainer)
- CWClassicFunctionNet (CMcsUsbNet[^] mcsusb)
- uint32 t ResetChannelmap (unsigned int virtualDevice)
- uint32 t SetChannelmap (unsigned char position, unsigned char channel, unsigned int Device)
- void SetHWSelectedChannels (array< bool >^ channels, unsigned int Device)
- void SetRFLostBehaviour (uint8 t stoponfailure, unsigned int Device)
- void SetHeadstageOnOff (uint16_t onoff)
- USHORT GetHeadstageOnOff ()
- void SetRFFrequencyHeadstage (uint8_t receiver_nb, unsigned short frequency)
- unsigned short GetRFFrequencyHeadstage (uint8_t receiver_nb)
- · void SetRFFrequencyReceiver (uint8 t receiver nb, uint8 t configuration, unsigned short frequency)
- void SetRFFrequencyReceiverEeprom (uint8_t receiver_nb, uint8_t configuration, unsigned short frequency)
- unsigned short GetRFFrequencyReceiver (uint8 t receiver nb, uint8 t configuration)
- void SetSerialNumberHeadstage (unsigned short number)
- unsigned short GetSerialNumberHeadstage ()
- void SetSelectedHeadstage (uint8_t number)
- uint8_t GetSelectedHeadstage ()
- void ScanForHeadstages ()
- uint8_t GetScanHeadstagesResult (int max_wait_for_ms)

- void SetFilterParametersHeadstage (unsigned short index, array< int >^ buffer)
- array< int > ^ GetFilterParametersHeadstage (unsigned short index)
- bool GetHasRedLedHeadstage ()
- · void SetHasChecksum (unsigned int has, unsigned int Device)
- unsigned int GetHasChecksum (unsigned int Device)
- void SetResetFilter (unsigned int reset, unsigned int Device)
- unsigned int GetResetFilter (unsigned int Device)
- void SetWPAType (unsigned short type, unsigned int Device)
- · unsigned short GetWPAType (unsigned int Device)
- void SetWPADebugMode (unsigned int mode, unsigned int Device)
- unsigned int GetWPADebugMode (unsigned int Device)
- void SetRFPower (unsigned short power)
- unsigned short GetRFPower ()
- unsigned int GetRFConnectionStatus ()

Additional Inherited Members

11.127.1 Constructor & Destructor Documentation

```
11.127.2.1 GetFilterParametersHeadstage() array<int> ^{\land} GetFilterParametersHeadstage ( unsigned short index )
```

```
11.127.2.2 GetHasChecksum() unsigned int GetHasChecksum ( unsigned int Device )
```

11.127.2.3 GetHasRedLedHeadstage() bool GetHasRedLedHeadstage ()


```
11.127.2.4 GetHeadstageOnOff() USHORT GetHeadstageOnOff ()
11.127.2.5 GetResetFilter() unsigned int GetResetFilter (
             unsigned int Device )
\textbf{11.127.2.6} \quad \textbf{GetRFConnectionStatus()} \quad \textbf{unsigned int GetRFConnectionStatus ()}
11.127.2.7 GetRFFrequencyHeadstage() unsigned short GetRFFrequencyHeadstage (
             uint8_t receiver_nb )
11.127.2.8 GetRFFrequencyReceiver() unsigned short GetRFFrequencyReceiver (
             uint8_t receiver_nb,
             uint8_t configuration )
11.127.2.9 GetRFPower() unsigned short GetRFPower ()
11.127.2.10 GetScanHeadstagesResult() uint8_t GetScanHeadstagesResult (
             int max\_wait\_for\_ms )
11.127.2.11 GetSelectedHeadstage() uint8_t GetSelectedHeadstage ( )
11.127.2.12 GetSerialNumberHeadstage() unsigned short GetSerialNumberHeadstage ( )
11.127.2.13 GetWPADebugMode() unsigned int GetWPADebugMode (
             unsigned int Device )
```

```
\textbf{11.127.2.14} \quad \textbf{GetWPAType()} \quad \texttt{unsigned short GetWPAType (}
             unsigned int Device )
11.127.2.15 ResetChannelmap() uint32_t ResetChannelmap (
             unsigned int virtualDevice )
11.127.2.16 ScanForHeadstages() void ScanForHeadstages ()
11.127.2.17 SetChannelmap() uint32_t SetChannelmap (
             unsigned char position,
             unsigned char channel,
             unsigned int Device )
11.127.2.18 SetFilterParametersHeadstage() void SetFilterParametersHeadstage (
             unsigned short index,
             array< int >^{\wedge} buffer )
11.127.2.19 SetHasChecksum() void SetHasChecksum (
             unsigned int has,
             unsigned int Device )
11.127.2.20 SetHeadstageOnOff() void SetHeadstageOnOff (
             uint16_t onoff )
11.127.2.21 SetHWSelectedChannels() void SetHWSelectedChannels (
             array< bool >^{\wedge} channels,
             unsigned int Device )
11.127.2.22 SetResetFilter() void SetResetFilter (
             unsigned int reset,
             unsigned int Device )
```

```
\textbf{11.127.2.23} \quad \textbf{SetRFF} \textbf{requencyHeadstage()} \quad \textbf{void SetRFF} \textbf{requencyHeadstage ()}
              uint8_t receiver_nb,
              unsigned short frequency )
11.127.2.24 SetRFFrequencyReceiver() void SetRFFrequencyReceiver (
              uint8_t receiver_nb,
              uint8_t configuration,
              unsigned short frequency )
11.127.2.25 SetRFFrequencyReceiverEeprom() void SetRFFrequencyReceiverEeprom (
              uint8_t receiver_nb,
              uint8_t configuration,
              unsigned short frequency )
11.127.2.26 SetRFLostBehaviour() void SetRFLostBehaviour (
              uint8_t stoponfailure,
              unsigned int Device )
11.127.2.27 SetRFPower() void SetRFPower (
              unsigned short power )
11.127.2.28 SetSelectedHeadstage() void SetSelectedHeadstage (
              uint8_t number )
11.127.2.29 SetSerialNumberHeadstage() void SetSerialNumberHeadstage (
              unsigned short number )
{\bf 11.127.2.30} \quad {\bf SetWPADebugMode()} \quad {\tt void SetWPADebugMode \ (}
              unsigned int mode,
              unsigned int Device )
11.127.2.31 SetWPAType() void SetWPAType (
              unsigned short type,
              unsigned int Device )
```

11.128 CWirelessBaseFunctionNet Class Reference

Inheritance diagram for CWirelessBaseFunctionNet:



Public Member Functions

• CWirelessBaseFunctionNet (CMcsUsbNet^ mcsusb, CMcsUsbFunctionPointerContainer^ mcsusbfunction)

Static Public Member Functions

• static String \(^\) CreateWirelessHeadstageSerialNumberString (unsigned short ID)

Additional Inherited Members

11.128.1 Constructor & Destructor Documentation

```
11.128.1.1 CWirelessBaseFunctionNet() CWirelessBaseFunctionNet (
CMcsUsbNet^ mcsusb,
CMcsUsbFunctionPointerContainer^ mcsusbfunction)
```

11.128.2 Member Function Documentation

```
11.128.2.1 CreateWirelessHeadstageSerialNumberString() static String ^ CreateWirelessHeadstage←
SerialNumberString (
    unsigned short ID ) [static]
```

11.129 DeviceIdNet Struct Reference

Device Id.

Public Member Functions

- DeviceIdNet ()
- DeviceIdNet (VendorIdEnumNet vendor, ProductIdEnumNet product, int bcd, McsBusTypeEnumNet bustype)
- DeviceIdNet (DeviceIdNet% deviceId)
- DeviceIdNet operator= (DeviceIdNet% deviceId)

Public Attributes

- VendorldEnumNet IdVendor
- ProductIdEnumNet IdProduct
- int BcdDevice
- McsBusTypeEnumNet BusType

11.129.1 Detailed Description

Device Id.

11.129.2 Constructor & Destructor Documentation

```
11.129.2.1 DeviceIdNet() [1/3] DeviceIdNet ( )
```

```
\textbf{11.129.2.2} \quad \textbf{DeviceIdNet()} \; \texttt{[2/3]} \quad \texttt{DeviceIdNet} \; \; \texttt{(}
```

```
VendorIdEnumNet vendor,
ProductIdEnumNet product,
int bcd,
McsBusTypeEnumNet bustype )
```

```
11.129.2.3 DeviceIdNet() [3/3] DeviceIdNet (

DeviceIdNet% deviceId)
```

11.129.3 Member Function Documentation

```
11.129.3.1 operator=() DeviceIdNet operator= (
DeviceIdNet% deviceId)
```

11.129.4 Member Data Documentation

11.129.4.1 BcdDevice int BcdDevice

11.129.4.2 BusType McsBusTypeEnumNet BusType

11.129.4.3 IdProduct ProductIdEnumNet IdProduct

11.129.4.4 IdVendor VendorIdEnumNet IdVendor

11.130 DigitalSource< digitalsourceenum > Class Template Reference

Public Member Functions

- DigitalSource ()
- DigitalSource (digitalsourceenum source)
- int MaxBitNumber ()
- int MaxBitNumber (digitalsourceenum Source)

Static Public Member Functions

- static int MaxBitNumberStatic (digitalsourceenum Source)
- static int size ()

Properties

• digitalsourceenum Source [get, set]

11.130.1 Constructor & Destructor Documentation

11.130.1.1 DigitalSource() [1/2] DigitalSource ()

```
11.130.1.2 DigitalSource() [2/2] DigitalSource (
             digitalsourceenum source )
11.130.2 Member Function Documentation
11.130.2.1 MaxBitNumber() [1/2] int MaxBitNumber ( )
11.130.2.2 MaxBitNumber() [2/2] int MaxBitNumber (
             digitalsourceenum Source )
11.130.2.3 MaxBitNumberStatic() static int MaxBitNumberStatic (
             digitalsourceenum Source ) [static]
11.130.2.4 size() static int size ( ) [static]
11.130.3 Property Documentation
11.130.3.1 Source digitalsourceenum Source [get], [set]
```

11.131 DigitalSourceGeneral Class Reference

Public Member Functions

- DigitalSourceGeneral (Type^ type)
- DigitalSourceGeneral (Type[^] type, int Source)
- int MaxBitNumber ()
- int MaxBitNumber (int Source)

Static Public Member Functions

- static int MaxBitNumber (Type[∧] type, int Source)
- static int size (Type[^] type)

Properties

```
• int Source [get, set]
```

11.131.1 Constructor & Destructor Documentation

11.131.2 Member Function Documentation

int Source)

```
11.131.2.1 MaxBitNumber() [1/3] int MaxBitNumber ( )
```

```
11.131.2.3 MaxBitNumber() [3/3] static int MaxBitNumber (

Type^ type,

int Source ) [static]
```

11.131.3 Property Documentation

```
\textbf{11.131.3.1} \quad \textbf{Source} \quad \texttt{int Source} \quad \texttt{[get], [set]}
```

11.132 DriverVersionNet Class Reference

Class gives firmware versions of the device's firmware destinations.

Public Member Functions

DriverVersionNet ()

Contructor.

∼DriverVersionNet ()

Destructor.

• unsigned int GetStatus (CFirmwareDestinationNet dest)

Get status of firmware destination.

unsigned int GetStatus (unsigned int index)

Get status of firmware destination.

• unsigned int GetVersionInt (CFirmwareDestinationNet dest)

Get the version number of firmware destination (major in high word, minor in low word)

unsigned int GetVersionInt (unsigned int index)

Get the version number of firmware destination (major in high word, minor in low word)

unsigned int GetMajor (CFirmwareDestinationNet dest)

Get the major version number of firmware destination.

• unsigned int GetMajor (unsigned int index)

Get the major version number of firmware destination.

unsigned int GetMinor (CFirmwareDestinationNet dest)

Get the minor version number of firmware destination.

• unsigned int GetMinor (unsigned int index)

Get the minor version number of firmware destination.

• unsigned int GetNumEntries ()

Get the number of available firmware destinations.

String \(^\) GetVersionString (CFirmwareDestinationNet dest)

Get the version as a string in the format Major. Minor.

String \(^\) GetVersionString (unsigned int index)

Get the version as a string in the format Major. Minor.

CFirmwareDestinationNet GetDestinationCode (unsigned int index)

Get CFirmwareDestinationNet.

String \(^\) GetDestinationName (CFirmwareDestinationNet dest)

Get firmware destination name.

String \(^\) GetDestinationName (unsigned int index)

Get firmware destination name.

String \(^\) GetSerialNumber (CFirmwareDestinationNet dest)

Get the serial number of the destination, when no serial number if found, return an empty string.

String \(^\) GetSerialNumber (unsigned int index)

Get the serial number of the destination, when no serial number if found, return an empty string.

Static Public Member Functions

static String \(^\) DriverVersionNet::FormatVersion (unsigned int v)

11.132.1 Detailed Description

Class gives firmware versions of the device's firmware destinations.

11.132.2 Constructor & Destructor Documentation

```
11.132.2.1 DriverVersionNet() DriverVersionNet ( )
```

Contructor.

```
11.132.2.2 \sim DriverVersionNet() \sim DriverVersionNet ()
```

Destructor.

11.132.3 Member Function Documentation

```
11.132.3.1 DriverVersionNet::FormatVersion() static String ^{\land} DriverVersionNet::FormatVersion ( unsigned int v ) [static]
```

```
11.132.3.2 GetDestinationCode() CFirmwareDestinationNet GetDestinationCode ( unsigned int index )
```

Get CFirmwareDestinationNet.

Parameters

index by index of firmware destination

```
11.132.3.3 GetDestinationName() [1/2] String ^{\land} GetDestinationName ( CFirmwareDestinationNet dest )
```

Get firmware destination name.

Parameters

dest	by CFirmwareDestionationNet
acci	by or initivate become tation in ter

11.132.3.4 GetDestinationName() [2/2] String $^{\land}$ GetDestinationName (unsigned int index)

Get firmware destination name.

Parameters

index by index of firmware destination	
--	--

11.132.3.5 GetMajor() [1/2] unsigned int GetMajor (CFirmwareDestinationNet dest)

Get the major version number of firmware destination.

Parameters

11.132.3.6 GetMajor() [2/2] unsigned int GetMajor (unsigned int index)

Get the major version number of firmware destination.

Parameters

Get the minor version number of firmware destination.

Parameters

```
dest by CFirmwareDestionationNet
```

11.132.3.8 **GetMinor()** [2/2] unsigned int GetMinor (unsigned int *index*)

Get the minor version number of firmware destination.

Parameters

index by index of firmware destination	
--	--

11.132.3.9 GetNumEntries() unsigned int GetNumEntries ()

Get the number of available firmware destinations.

```
11.132.3.10 GetSerialNumber() [1/2] String ^{\land} GetSerialNumber ( CFirmwareDestinationNet dest )
```

Get the serial number of the destination, when no serial number if found, return an empty string.

Parameters

```
dest by CFirmwareDestionationNet
```

```
11.132.3.11 GetSerialNumber() [2/2] String ^ GetSerialNumber ( unsigned int index )
```

Get the serial number of the destination, when no serial number if found, return an empty string.

Parameters

index by index of firmware destination

```
11.132.3.12 GetStatus() [1/2] unsigned int GetStatus (
CFirmwareDestinationNet dest)
```

Get status of firmware destination.

Parameters

dest	by CFirmwareDestionationNet
------	-----------------------------

```
11.132.3.13 GetStatus() [2/2] unsigned int GetStatus ( unsigned int index )
```

Get status of firmware destination.

Parameters

index by index of firmware destination

11.132.3.14 GetVersionInt() [1/2] unsigned int GetVersionInt (CFirmwareDestinationNet dest)

Get the version number of firmware destination (major in high word, minor in low word)

Parameters

dest by CFirmwareDestionationNet

11.132.3.15 GetVersionInt() [2/2] unsigned int GetVersionInt (unsigned int index)

Get the version number of firmware destination (major in high word, minor in low word)

Parameters

index by index of firmware destination

```
11.132.3.16 GetVersionString() [1/2] String ^{\land} GetVersionString ( CFirmwareDestinationNet dest )
```

Get the version as a string in the format Major.Minor.

Parameters

dest by CFirmwareDestionationNet

```
11.132.3.17 GetVersionString() [2/2] String ^{\land} GetVersionString ( unsigned int index )
```

Get the version as a string in the format Major.Minor.

Parameters

index by index of firmware

11.133 FirmwareDestinationNames Class Reference

Static Public Attributes

```
    static String \(^\text{DSP} = \text{gcnew String("DSP")}\)

• static String ^ USB = gcnew String( "USB" )

    static String \(^{\text{MCU1}} = \text{gcnew String("MCU1")}\)

    static String \(^\) Bootstrap = gcnew String( "Bootstrap" )

    static String \(^{\text{MCSBUS1}} = \text{gcnew String( "McsBus1" )}\)

    static String \(^{\text{MCSBUS2}} = \text{gcnew String( "McsBus2" )}\)

    static String \(^{\text{MCSBUS3}} = \text{gcnew String( "McsBus3" )}\)

    static String \(^\text{MCSBUS4} = \text{gcnew String( "McsBus4" )}\)

    static String \(^\text{MCSBUS5} = \text{gcnew String( "McsBus5" )}\)

    static String \(^{\text{MCSBUS6}} = \text{gcnew String( "McsBus6" )}\)

    static String \(^{\text{MCSBUS7}} = \text{gcnew String( "McsBus7" )}\)

    static String \(^{\text{MCSBUS8}} = \text{gcnew String( "McsBus8" )}\)

    static String \(^\text{MCSBUS9} = \text{gcnew String( "McsBus9" )}\)

    static String \(^{\text{MCSBUS10}} = \text{gcnew String("McsBus10")}\)

    static String \(^{\text{MCSBUS11}} = \text{gcnew String( "McsBus11" )}\)

    static String \(^{\text{MCSBUS12}} = \text{gcnew String("McsBus12")}\)

    static String \(^\text{MCSBUS13} = \text{gcnew String( "McsBus13" )}\)

    static String \(^\text{BUS1}_MCSBUS1\) = gcnew String("Bus1McsBus1")

    static String \(^\text{BUS1_MCSBUS2} = \text{gcnew String("Bus1McsBus2")}\)

    static String ^ PIC = gcnew String("PIC")

    static String \(^\text{PIC2} = \text{gcnew String("PIC2")}\)

    static String \(^\text{PIC3} = \text{gcnew String("PIC3")}\)

    static String ^ PIC4 = gcnew String( "PIC4" )

    static String ^ Altera = gcnew String( "Altera" )

    static String \(^\text{FPGA2} = \text{gcnew String("FPGA2")}\)

    static String \(^\text{FPGA3} = \text{gcnew String("FPGA3")}\)

    static String \(^\) FPGA4 = gcnew String( "FPGA4" )

    static String \(^{\text{FPGA5}} = \text{gcnew String( "FPGA5" )}\)

    static String \(^\text{FPGA6} = \text{gcnew String("FPGA6")}\)
```

11.133.1 Member Data Documentation

```
11.133.1.1 Altera String ^ Altera = gcnew String( "Altera" ) [static]
```

```
11.133.1.2 Bootstrap String ^ Bootstrap = gcnew String( "Bootstrap" ) [static]
```

```
11.133.1.3 BUS1_MCSBUS1 String ^ BUS1_MCSBUS1 = gcnew String( "Bus1McsBus1" ) [static]
11.133.1.4 BUS1_MCSBUS2 String ^ BUS1_MCSBUS2 = gcnew String( "Bus1McsBus2") [static]
11.133.1.5 DSP String ^{\wedge} DSP = gcnew String("DSP") [static]
11.133.1.6 FPGA2 String ^{\wedge} FPGA2 = gcnew String( "FPGA2" ) [static]
11.133.1.7 FPGA3 String ^ FPGA3 = gcnew String( "FPGA3" ) [static]
11.133.1.8 FPGA4 String ^ FPGA4 = gcnew String( "FPGA4" ) [static]
11.133.1.9 FPGA5 String ^ FPGA5 = gcnew String( "FPGA5" ) [static]
11.133.1.10 FPGA6 String ^ FPGA6 = gcnew String( "FPGA6" ) [static]
\textbf{11.133.1.11} \quad \textbf{MCSBUS1} \quad \texttt{String} \; \land \; \texttt{MCSBUS1} \; = \; \texttt{gcnew} \; \texttt{String} ( \; \texttt{"McsBus1"} \; ) \quad \texttt{[static]}
11.133.1.12 MCSBUS10 String ^ MCSBUS10 = gcnew String( "McsBus10" ) [static]
11.133.1.13 MCSBUS11 String ^{\wedge} MCSBUS11 = gcnew String( "McsBus11" ) [static]
```

```
11.133.1.14 MCSBUS12 String ^{\wedge} MCSBUS12 = gcnew String( "McsBus12" ) [static]
11.133.1.15 MCSBUS13 String ^{\wedge} MCSBUS13 = gcnew String( "McsBus13" ) [static]
11.133.1.16 MCSBUS2 String ^{\land} MCSBUS2 = gcnew String( "McsBus2" ) [static]
11.133.1.17 MCSBUS3 String ^{\wedge} MCSBUS3 = gcnew String( "McsBus3" ) [static]
11.133.1.18 MCSBUS4 String ^ MCSBUS4 = gcnew String( "McsBus4" ) [static]
11.133.1.19 MCSBUS5 String ^ MCSBUS5 = gcnew String( "McsBus5" ) [static]
11.133.1.20 MCSBUS6 String ^ MCSBUS6 = gcnew String( "McsBus6" ) [static]
11.133.1.21 MCSBUS7 String ^ MCSBUS7 = gcnew String( "McsBus7" ) [static]
\textbf{11.133.1.22} \quad \textbf{MCSBUS8} \quad \texttt{String} \; \land \; \texttt{MCSBUS8} \; = \; \texttt{gcnew} \; \texttt{String} ( \; \texttt{"McsBus8"} \; ) \quad \texttt{[static]}
11.133.1.23 MCSBUS9 String ^ MCSBUS9 = gcnew String( "McsBus9" ) [static]
11.133.1.24 MCU1 String ^{\wedge} MCU1 = gcnew String( "MCU1" ) [static]
```

```
11.133.1.25 PIC String ^ PIC = genew String( "PIC" ) [static]

11.133.1.26 PIC2 String ^ PIC2 = genew String( "PIC2" ) [static]


11.133.1.27 PIC3 String ^ PIC3 = genew String( "PIC3" ) [static]

11.133.1.28 PIC4 String ^ PIC4 = genew String( "PIC4" ) [static]

11.133.1.29 USB String ^ USB = genew String( "USB" ) [static]
```

11.134 HeadStageIDType Class Reference

Inheritance diagram for HeadStageIDType:



Public Types

enum HeadstageTypeEnum {
 Unknown,
 MeasuringOnly,
 OpticalStimulation,
 ElectricalStimulation }

Public Member Functions

- HeadStageIDType (unsigned int entry, CW2100_FunctionNet[^] device)
- virtual System::String ^ ToString () override
- virtual bool Equals (Object[^] obj) override
- virtual Int32 CompareTo (Object[^] obj)

Properties

```
bool Valid [get]
unsigned int Entry [get]
unsigned short ID [get]
System::String^ SN [get]
unsigned int TypeValue [get]
System::String^ Type [get]
HeadstageTypeEnum HeadstageType [get]
System::String^ UserDefinedName [get]
int NumberOfAnalogChannels [get]
int NumberOfStimulationChannels [get]
W2100_StimulusParametersNet^ StimulusParameters [get]
bool HasIMU [get]
bool W16IsW14 [get]
bool HasOptoCurrentMessurement [get]
```

11.134.1 Member Enumeration Documentation

11.134.1.1 HeadstageTypeEnum enum HeadstageTypeEnum [strong]

Enumerator

Unknown	
MeasuringOnly	
OpticalStimulation	
ElectricalStimulation	

11.134.2 Constructor & Destructor Documentation

11.134.3 Member Function Documentation

```
11.134.3.1 CompareTo() virtual Int32 CompareTo (
Object^ obj ) [virtual]
```

```
11.134.3.2 Equals() virtual bool Equals (
            Object^ obj ) [override], [virtual]
11.134.3.3 ToString() virtual System::String ^ ToString ( ) [override], [virtual]
11.134.4 Property Documentation
11.134.4.1 Entry unsigned int Entry [get]
11.134.4.2 HasIMU bool HasIMU [get]
11.134.4.3 HasOptoCurrentMessurement bool HasOptoCurrentMessurement [get]
11.134.4.4 HeadstageType HeadstageTypeEnum HeadstageType [get]
11.134.4.5 ID unsigned short ID [get]
11.134.4.6 NumberOfAnalogChannels int NumberOfAnalogChannels [get]
11.134.4.7 NumberOfStimulationChannels int NumberOfStimulationChannels [get]
11.134.4.8 SN System:: String^{\wedge} SN [get]
```

```
11.134.4.10 Type System:: String^ Type [get]

11.134.4.11 TypeValue unsigned int TypeValue [get]

11.134.4.12 UserDefinedName System:: String^ UserDefinedName [get]

11.134.4.13 Valid bool Valid [get]
```

11.135 HeadstageIDTypeObject Class Reference

Public Member Functions

- HeadstageIDTypeObject (HeadStageIDType^ idType)
- virtual String ^ ToString () override
- virtual bool Equals (Object[^] obj) override
- virtual int GetHashCode () override

Public Attributes

- HeadStageIDType ^ _IdType
- String ^ _AdditionalText

Properties

- String^ AdditionalText [get, set]

11.135.1 Constructor & Destructor Documentation

```
11.135.1.1 HeadstageIDTypeObject() HeadstageIDTypeObject (
            HeadStageIDType^ idType )
11.135.2 Member Function Documentation
11.135.2.1 Equals() virtual bool Equals (
            Object^ obj ) [override], [virtual]
11.135.2.2 GetHashCode() virtual int GetHashCode ( ) [override], [virtual]
11.135.2.3 ToString() virtual String ^{\wedge} ToString ( ) [override], [virtual]
11.135.3 Member Data Documentation
11.135.3.1 _AdditionalText String ^ _AdditionalText
11.135.3.2 _ldType HeadStageIDType ^ _IdType
11.135.4 Property Documentation
11.135.4.1 AdditionalText String AdditionalText [get], [set]
11.135.4.2 IdType HeadStageIDType^ IdType [get]
11.136 HeadStageIDTypeState Class Reference
Properties
```

• unsigned int State [get]

bool ControlState [get]bool DataState [get]

HeadStageIDType[^] IdType [get]

11.136.1 Property Documentation

```
11.136.1.1 ControlState bool ControlState [get]

11.136.1.2 DataState bool DataState [get]

11.136.1.3 IdType HeadStageIDType^ IdType [get]
```

11.137 mkfilterNet Class Reference

11.136.1.4 State unsigned int State [get]

Static Public Member Functions

- static int mkfilter (String[^] filtertype, double value, String[^] passtype, int order, double alpha1, double alpha2, [System::Runtime::InteropServices::Out] array< double >^% xcoeffs, [System::Runtime::InteropServices::Out] array< double >^% ycoeffs)
- static int mkfilter_MCS (int SamplesPerSecond, double R1, double R2, double C, double Amplification, double Correction, [System::Runtime::InteropServices::Out] array< double >^% xcoeffs, [System::Runtime::
 InteropServices::Out] array< double >^% ycoeffs)
- static int mkfilter_MCS (int SamplesPerSecond, double R1, double R2, double C, double Correction, [System::Runtime::InteropServices::Out] array< double >^% xcoeffs, [System::Runtime::InteropServices::Out] array< double >^% ycoeffs)
- static int mkfilter_MCS_k (int SamplesPerSecond, double R1, double R2, double C, double Amplification, double Correction, [System::Runtime::InteropServices::Out] array< double >^% coeffs)
- static int mkfilter_MCS_k (int SamplesPerSecond, double R1, double R2, double C, double Correction, [System::Runtime::InteropServices::Out] array< double >^% coeffs)
- static void mkfilter_coef_in_one_set (int n, [System::Runtime::InteropServices::In] array< double >^ xcoeffs, [System::Runtime::InteropServices::Out] array< double >^% out_coeffs)
- static void mkfilter_scale_coef_in_one_set (int n, double scale, [System::Runtime::InteropServices::In] array< double >^ xcoeffs, [System::Runtime::InteropServices::In] array< double >^ ycoeffs, [System::

 Runtime::InteropServices::Out] array< double >^% out_coeffs)
- static void mkfilter_normalize_coeffs_short (short maxvalue, [System::Runtime::InteropServices::In] array
 double >^ coeffs, [System::Runtime::InteropServices::Out] array< short >^% out_coeffs)
- static void mkfilter_normalize_coeffs_int (int maxvalue, [System::Runtime::InteropServices::In] array< double
 coeffs, [System::Runtime::InteropServices::Out] array< int >^% out_coeffs)
- static void mkfilter_normalize_scale_coeffs_int (int maxvalue, [System::Runtime::InteropServices::In] array
 double >^ coeffs, [System::Runtime::InteropServices::Out] array< int >^% out_coeffs)
- static double mkfilter_highpass_coeff (int SamplesPerSecond, double Frequency)
- static double mkfilter_highpass_k (int SamplesPerSecond, double Frequency)
- static double mkfilter_highpass_frequency_from_coeff (int SamplesPerSecond, double coeff)
- static double mkfilter_highpass_frequency_from_k (int SamplesPerSecond, double k)

11.137.1 Member Function Documentation

```
11.137.1.1 mkfilter() static int mkfilter (
             String^{\wedge} filtertype,
             double value,
             String^ passtype,
             int order,
             double alpha1,
             double alpha2,
             [System::Runtime::InteropServices::Out] array< double >^{\%} xcoeffs,
             [System::Runtime::InteropServices::Out] array<br/> double >^{\%} ycoeffs ) [static]
11.137.1.2 mkfilter_coef_in_one_set() static void mkfilter_coef_in_one_set (
             int n_{i}
             [System::Runtime::InteropServices::In] array< double >^{\wedge} xcoeffs,
             [System::Runtime::InteropServices::In] array< double >^{\wedge} ycoeffs,
             [System::Runtime::InteropServices::Out] array< double >^{\%} out_coeffs ) [static]
11.137.1.3 mkfilter_highpass_coeff() static double mkfilter_highpass_coeff (
             int SamplesPerSecond,
             double Frequency ) [static]
11.137.1.4 mkfilter_highpass_frequency_from_coeff() static double mkfilter_highpass_frequency_←
from_coeff (
             int SamplesPerSecond,
             double coeff ) [static]
11.137.1.5 mkfilter_highpass_frequency_from_k() static double mkfilter_highpass_frequency_from←
_k (
             int SamplesPerSecond,
             double k ) [static]
11.137.1.6 mkfilter_highpass_k() static double mkfilter_highpass_k (
             int SamplesPerSecond,
             double Frequency ) [static]
```

```
11.137.1.7 mkfilter_MCS() [1/2] static int mkfilter_MCS (
              int SamplesPerSecond,
             double R1,
             double R2,
             double C,
             double Amplification,
             double Correction,
             [System::Runtime::InteropServices::Out] array< double >^% xcoeffs,
             [System::Runtime::InteropServices::Out] array<br/> double >^{\%} ycoeffs ) [static]
11.137.1.8 mkfilter\_MCS() [2/2] static int mkfilter\_MCS (
             int SamplesPerSecond,
             double R1.
             double R2,
             double C,
             double Correction,
             [System::Runtime::InteropServices::Out] array< double >^% xcoeffs,
             [System::Runtime::InteropServices::Out] array< double >^{\%} ycoeffs ) [static]
11.137.1.9 mkfilter\_MCS\_k() [1/2] static int mkfilter\_MCS\_k (
             int SamplesPerSecond,
             double R1,
             double R2,
             double C_{\prime}
             double Amplification,
             double Correction,
             [System::Runtime::InteropServices::Out] array< double >^{\%} coeffs ) [static]
11.137.1.10 mkfilter_MCS_k() [2/2] static int mkfilter_MCS_k (
             int SamplesPerSecond,
             double R1,
             double R2,
             double C,
             double Correction,
             [System::Runtime::InteropServices::Out] array< double >^{^{\land}} % coeffs ) [static]
11.137.1.11 mkfilter_normalize_coeffs_int() static void mkfilter_normalize_coeffs_int (
             int maxvalue,
              [System::Runtime::InteropServices::In] array< double >^{\land} coeffs,
              [System::Runtime::InteropServices::Out] array< int >^% out_coeffs ) [static]
```

11.138 CRoboDeviceNet::RoboMainLowLevelCommands Class Reference

Public Member Functions

- void SetParameter (unsigned short command, unsigned short index, unsigned int value)
- void SetParameter (unsigned short command, unsigned short index, unsigned int value1, unsigned int value2)
- void SetUserParameter (unsigned short index, unsigned int value)

Stores persistently 32 bit integer values on RoboMain

· void SetUserParameter (unsigned short index, int value)

Stores persistently 32 bit integer values on RoboMain

- void GetParameter (unsigned short command, unsigned short index, [System::Runtime::InteropServices::
 Out]unsigned int% value)
- void GetParameter (unsigned short command, unsigned short index, [System::Runtime::InteropServices::
 Out]unsigned int% value1, [System::Runtime::InteropServices::Out]unsigned int% value2)
- void GetUserParameter (unsigned short index, [System::Runtime::InteropServices::Out]unsigned int% value)

Reads 32 bit integer values stored persistently on RoboMain

void GetUserParameter (unsigned short index, [System::Runtime::InteropServices::Out]int% value)

Reads 32 bit integer values stored persistently on RoboMain

- void FindReferencePhase0 (unsigned char busaddress, char axes)
- void FindReferencePhase0 (unsigned char busaddress, char axes, int timeout)
- unsigned char HasRef (unsigned char busaddress, char axes)
- void SetHWRevision (unsigned int revision)
- unsigned int GetHWRevision ()
- · void SetHWConfig (unsigned int config)
- unsigned int GetHWConfig ()
- void SetMinPressureWaitTime (unsigned int t)
- unsigned int GetMinPressureWaitTime ()
- void SetMinPressure (unsigned int pressure)
- unsigned int GetMinPressure ()
- void SetMaxPressureWaitTime (unsigned int t)

- unsigned int GetMaxPressureWaitTime ()
- void SetMinNoPressureWaitTime (unsigned int t)
- unsigned int GetMinNoPressureWaitTime ()
- void SetMaxNoPressure (unsigned int pressure)
- unsigned int GetMaxNoPressure ()
- void SetMaxNoPressureWaitTime (unsigned int t)
- unsigned int GetMaxNoPressureWaitTime ()
- void SetSearchReferenceMethod (unsigned char busaddress, char axes, unsigned int method)
- unsigned int GetSearchReferenceMethod (unsigned char busaddress, char axes)
- void SetSearchReferenceOffsetPos (unsigned char busaddress, char axes, int offsetpos)
- int GetSearchReferenceOffsetPos (unsigned char busaddress, char axes)
- void SetSearchReferenceFastSpeed (unsigned char busaddress, char axes, unsigned short speed)
- unsigned short GetSearchReferenceFastSpeed (unsigned char busaddress, char axes)
- void SetSearchReferenceFastAccel (unsigned char busaddress, char axes, unsigned short accel)
- unsigned short GetSearchReferenceFastAccel (unsigned char busaddress, char axes)
- void SetSearchReferenceFineSpeed (unsigned char busaddress, char axes, unsigned short speed)
- unsigned short GetSearchReferenceFineSpeed (unsigned char busaddress, char axes)
- void SetSearchReferenceFineAccel (unsigned char busaddress, char axes, unsigned short accel)
- unsigned short GetSearchReferenceFineAccel (unsigned char busaddress, char axes)
- void SetSearchReferenceMoveOut (unsigned char busaddress, char axes, int move)
- int GetSearchReferenceMoveOut (unsigned char busaddress, char axes)
- void SetAxisConfig (unsigned char busaddress, char axes, unsigned int config)
- unsigned int GetAxisConfig (unsigned char busaddress, char axes)
- void GetPhases (unsigned char busaddress, char axes, [System::Runtime::InteropServices::Out] unsigned short% phase0, [System::Runtime::InteropServices::Out] unsigned short% lastphase)

11.138.1 Member Function Documentation

```
11.138.1.4 GetHWConfig() unsigned int GetHWConfig ( )
11.138.1.5 GetHWRevision() unsigned int GetHWRevision ()
11.138.1.6 GetMaxNoPressure() unsigned int GetMaxNoPressure ( )
11.138.1.7 GetMaxNoPressureWaitTime() unsigned int GetMaxNoPressureWaitTime ( )
11.138.1.8 GetMaxPressureWaitTime() unsigned int GetMaxPressureWaitTime ( )
11.138.1.9 GetMinNoPressureWaitTime() unsigned int GetMinNoPressureWaitTime ()
11.138.1.10 GetMinPressure() unsigned int GetMinPressure ( )
11.138.1.11 GetMinPressureWaitTime() unsigned int GetMinPressureWaitTime ( )
11.138.1.12 GetParameter() [1/2] void GetParameter (
             unsigned short command,
             unsigned short index,
             [System::Runtime::InteropServices::Out] unsigned int% value )
11.138.1.13 GetParameter() [2/2] void GetParameter (
             unsigned short command,
             unsigned short index,
             [System::Runtime::InteropServices::Out] unsigned int% value1,
             [System::Runtime::InteropServices::Out] unsigned int% value2 )
```

```
11.138.1.14 GetPhases() void GetPhases (
              unsigned char busaddress,
              char axes,
              [System::Runtime::InteropServices::Out] unsigned short% phase0,
              [{\tt System::Runtime::InteropServices::Out}] \ {\tt unsigned \ short} \ {\tt lastphase} \ )
11.138.1.15 GetSearchReferenceFastAccel() unsigned short GetSearchReferenceFastAccel (
             unsigned char busaddress,
              char axes )
11.138.1.16 GetSearchReferenceFastSpeed() unsigned short GetSearchReferenceFastSpeed (
              unsigned char busaddress,
              char axes )
11.138.1.17 GetSearchReferenceFineAccel() unsigned short GetSearchReferenceFineAccel (
              unsigned char busaddress,
              char axes )
11.138.1.18 GetSearchReferenceFineSpeed() unsigned short GetSearchReferenceFineSpeed (
              unsigned char busaddress,
              char axes )
11.138.1.19 GetSearchReferenceMethod() unsigned int GetSearchReferenceMethod (
             unsigned char busaddress,
              char axes )
11.138.1.20 GetSearchReferenceMoveOut() int GetSearchReferenceMoveOut (
             unsigned char busaddress,
              char axes )
\textbf{11.138.1.21} \quad \textbf{GetSearchReferenceOffsetPos()} \quad \texttt{int GetSearchReferenceOffsetPos} \quad \textbf{(}
              unsigned char busaddress,
              char axes )
11.138.1.22 GetUserParameter() [1/2] void GetUserParameter (
              unsigned short index,
              [System::Runtime::InteropServices::Out] int% value )
Reads 32 bit integer values stored persistently on RoboMain
```

Generated by Doxygen

intention: provide free persistent user memory space on motor controller

Parameters

index	address offset of parameter; range: 015
value	data buffer

Reads 32 bit integer values stored persistently on RoboMain

intention: provide free persistent user memory space on motor controller

Parameters

index address of		address offset of parameter; range: 015
ĺ	value	data buffer

```
11.138.1.24 HasRef() unsigned char HasRef (
unsigned char busaddress,
char axes )
```

```
11.138.1.26 SetHWConfig() void SetHWConfig ( unsigned int config )
```

```
11.138.1.27 SetHWRevision() void SetHWRevision ( unsigned int revision)
```

```
11.138.1.28 SetMaxNoPressure() void SetMaxNoPressure ( unsigned int pressure )
```

```
11.138.1.29 SetMaxNoPressureWaitTime() void SetMaxNoPressureWaitTime (
             unsigned int t )
11.138.1.30 SetMaxPressureWaitTime() void SetMaxPressureWaitTime (
             unsigned int t )
11.138.1.31 SetMinNoPressureWaitTime() void SetMinNoPressureWaitTime (
             unsigned int t )
11.138.1.32 SetMinPressure() void SetMinPressure (
             unsigned int pressure )
11.138.1.33 SetMinPressureWaitTime() void SetMinPressureWaitTime (
             unsigned int t )
11.138.1.34 SetParameter() [1/2] void SetParameter (
             unsigned short command,
             unsigned short index,
             unsigned int value )
11.138.1.35 SetParameter() [2/2] void SetParameter (
             unsigned short command,
             unsigned short index,
             unsigned int value1,
             unsigned int value2 )
11.138.1.36 SetSearchReferenceFastAccel() void SetSearchReferenceFastAccel (
             unsigned char busaddress,
             char axes,
             unsigned short \mathit{accel} )
```

```
\textbf{11.138.1.37} \quad \textbf{SetSearchReferenceFastSpeed()} \quad \texttt{void SetSearchReferenceFastSpeed} \quad \textbf{(}
              unsigned char busaddress,
              char axes,
              unsigned short speed )
11.138.1.38 SetSearchReferenceFineAccel() void SetSearchReferenceFineAccel (
              unsigned char busaddress,
              char axes,
              unsigned short accel )
11.138.1.39 SetSearchReferenceFineSpeed() void SetSearchReferenceFineSpeed (
              unsigned char busaddress,
              char axes,
              unsigned short speed )
11.138.1.40 SetSearchReferenceMethod() void SetSearchReferenceMethod (
              unsigned char busaddress,
              char axes,
              unsigned int method )
\textbf{11.138.1.41} \quad \textbf{SetSearchReferenceMoveOut()} \quad \texttt{void SetSearchReferenceMoveOut} \quad (
              unsigned char busaddress,
              char axes,
              int move )
11.138.1.42 SetSearchReferenceOffsetPos() void SetSearchReferenceOffsetPos (
              unsigned char busaddress,
              char axes,
              int offsetpos )
11.138.1.43 SetUserParameter() [1/2] void SetUserParameter (
              unsigned short index,
              int value )
```

Stores persistently 32 bit integer values on RoboMain

intention: provide free persistent user memory space on RoboMain

Parameters

index	address offset of parameter; range: 015
value	data to be stored

```
11.138.1.44 SetUserParameter() [2/2] void SetUserParameter (
unsigned short index,
unsigned int value)
```

Stores persistently 32 bit integer values on RoboMain

intention: provide free persistent user memory space on RoboMain

Parameters

index	address offset of parameter; range: 015
value	data to be stored

11.139 CRoboStatorDeviceNet::RoboMainStatorLowLevelCommands Class Reference

Public Member Functions

- void FindReferencePhase0XY ()
- void FindReferencePhase0XY (int timeout)

11.139.1 Member Function Documentation

```
11.139.1.1 FindReferencePhaseOXY() [1/2] void FindReferencePhaseOXY ( )
```

11.140 CFilterCoefficientsNet::s_FilterAttributesNet Struct Reference

Public Member Functions

- s_FilterAttributesNet (s_FilterAttributes attrib)
- s_FilterAttributes ToCpp ()

Public Attributes

- uint32_t PreCommaB
- uint32_t PostCommaB
- uint32_t CommaPositionB
- uint32_t PreCommaA
- uint32_t PostCommaA
- uint32_t CommaPositionA

11.140.1 Constructor & Destructor Documentation

```
11.140.1.1 s_FilterAttributesNet() s_FilterAttributesNet ( s_FilterAttributes attrib )
```

11.140.2 Member Function Documentation

```
\textbf{11.140.2.1} \quad \textbf{ToCpp()} \quad \texttt{s\_FilterAttributes ToCpp ()}
```

11.140.3 Member Data Documentation

- 11.140.3.1 CommaPositionA uint32_t CommaPositionA
- 11.140.3.2 CommaPositionB uint32_t CommaPositionB
- 11.140.3.3 PostCommaA uint32_t PostCommaA
- 11.140.3.4 PostCommaB uint32_t PostCommaB

```
11.140.3.5 PreCommaA uint32_t PreCommaA
```

```
11.140.3.6 PreCommaB uint32_t PreCommaB
```

11.141 CMeaAudioFunctionNet::s_setaudionet Struct Reference

Public Attributes

- int channel
- · int amplification

11.141.1 Member Data Documentation

```
11.141.1.1 amplification int amplification
```

```
11.141.1.2 channel int channel
```

11.142 CStimulusFunctionNet::SidebandData Class Reference

Public Member Functions

- SidebandData ()
- ∼SidebandData ()

Destructor: called by Dispose()

• !SidebandData ()

Finalizer: called by GC before collecting

Properties

```
• array< int32_t >^ Sideband [get]
```

• array< uint64_t >^ Duration [get]

11.142.1 Constructor & Destructor Documentation

11.142.1.1 SidebandData() SidebandData ()

```
11.142.1.2 ~SidebandData() ~SidebandData ()
```

Destructor: called by Dispose()

11.142.1.3 "!SidebandData() !SidebandData ()

Finalizer: called by GC before collecting

11.142.2 Property Documentation

```
11.142.2.1 Duration array< uint64_t>^{\wedge} Duration [get]
```

```
11.142.2.2 Sideband array< int32_t>^{\land} Sideband [get]
```

11.143 StgStatusNet Class Reference

Static Public Member Functions

- static StgStatusNet ^ FromIntPtr (IntPtr stgstatus)
- static StgStatusNet ^ FromPtr (stgstatus t *stgstatus)

Public Attributes

- array< Stg200xTriggerStatusEnumNet > ^ TiggerStatus
- array< uint32_t > ^ ListOfChangedTriggers

11.143.1 Member Function Documentation

```
11.143.1.2 FromPtr() static StgStatusNet ^ FromPtr (
    stgstatus_t * stgstatus ) [static]
```

11.143.2 Member Data Documentation

 $\textbf{11.143.2.1} \quad \textbf{ListOfChangedTriggers} \quad \texttt{array} < \texttt{uint32_t} > \\ \land \quad \texttt{ListOfChangedTriggers}$

 $\textbf{11.143.2.2} \quad \textbf{TiggerStatus} \quad \texttt{array} < \texttt{Stg200xTriggerStatusEnumNet} > \\ \land \quad \texttt{TiggerStatus}$

11.144 CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData Class Reference

Public Member Functions

- StimulusDeviceDataAndUnrolledData ()
- ~StimulusDeviceDataAndUnrolledData ()

Destructor: called by Dispose()

• !StimulusDeviceDataAndUnrolledData ()

Finalizer: called by GC before collecting

Properties

```
    array< uint8_t >^ DeviceData [get]
```

- int DeviceDataLength [get]
- array< int32_t >^ UnrolledAmplitude [get]
- array< uint32_t >^ UnrolledSync [get]
- array< uint64_t >^ UnrolledDuration [get]

11.144.1 Constructor & Destructor Documentation

11.144.1.1 StimulusDeviceDataAndUnrolledData() StimulusDeviceDataAndUnrolledData ()

11.144.1.2 ~StimulusDeviceDataAndUnrolledData() ~StimulusDeviceDataAndUnrolledData ()

Destructor: called by Dispose()

11.144.1.3 "!StimulusDeviceDataAndUnrolledData() !StimulusDeviceDataAndUnrolledData () Finalizer: called by GC before collecting 11.144.2 Property Documentation 11.144.2.1 DeviceData array< uint8_t>^ DeviceData [get] 11.144.2.2 DeviceDataLength int DeviceDataLength [get] $\textbf{11.144.2.3} \quad \textbf{UnrolledAmplitude} \quad \texttt{array} < \text{int32_t} > ^{\wedge} \text{UnrolledAmplitude} \quad \texttt{[get]}$ $\textbf{11.144.2.4} \quad \textbf{UnrolledDuration} \quad \texttt{array} < \ \texttt{uint64_t} > ^{\wedge} \ \texttt{UnrolledDuration} \quad \texttt{[get]}$ 11.144.2.5 UnrolledSync array< uint32_t>^ UnrolledSync [get]

11.145 usbSetupPacket_t Class Reference

Public Attributes

- uint8_t bmRequestType
- uint8_t bRequest
- uint16 t wValue
- uint16_t wlndex
- uint16_t wLength

11.145.1 Member Data Documentation

11.145.1.1 bmRequestType uint8_t bmRequestType

- 11.145.1.2 **bRequest** uint8_t bRequest
- 11.145.1.3 windex uint16_t wIndex
- 11.145.1.4 wLength uint16_t wLength
- **11.145.1.5 wValue** uint16_t wValue

11.146 W2100 StimulusParametersNet Struct Reference

Public Attributes

- int DACResolution
- int TimeResolutionInNanoSeconds
- int VoltageRangeInMicroVolt
- int VoltageResolutionInMicroVolt
- int CurrentRangeInNanoAmp
- int CurrentResolutionInNanoAmp

11.146.1 Member Data Documentation

- 11.146.1.1 CurrentRangeInNanoAmp int CurrentRangeInNanoAmp
- 11.146.1.2 CurrentResolutionInNanoAmp int CurrentResolutionInNanoAmp
- 11.146.1.3 DACResolution int DACResolution
- $\textbf{11.146.1.4} \quad \textbf{TimeResolutionInNanoSeconds} \quad \texttt{int TimeResolutionInNanoSeconds}$
- 11.146.1.5 VoltageRangeInMicroVolt int VoltageRangeInMicroVolt
- 11.146.1.6 VoltageResolutionInMicroVolt int VoltageResolutionInMicroVolt

Index

!CDacCalibrationFunctionNet	CUsbDeviceConfigurationFunctionNet, 524
CDacCalibrationFunctionNet, 51	!CWarnerUssingDeviceNet
!CDigOutStimulatorFunctionNet	CWarnerUssingDeviceNet, 543
CDigOutStimulatorFunctionNet, 61	!CWarnerUssingFunctionNet
!CExternDTesterDeviceNet	CWarnerUssingFunctionNet, 546
CExternDTesterDeviceNet, 65	!CWarnerValveControllerDeviceNet
•	
!CGrapheneFunctionNet	CWarnerValveControllerDeviceNet, 565
CGrapheneFunctionNet, 102	!CWarnerValveControllerDeviceTesterFunctionNet
!CInterfaceboard2FunctionNet	CWarnerValveControllerDeviceTesterFunctionNet,
CInterfaceboard2FunctionNet, 118	586
!CInterfaceboardFunctionNet	!SidebandData
CInterfaceboardFunctionNet, 120	CStimulusFunctionNet::SidebandData, 624
!CLIH3DeviceNet	!StimulusDeviceDataAndUnrolledData
CLIH3DeviceNet, 123	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData,
!CMEA2100x256FunctionNet	625
CMEA2100x256FunctionNet, 265	_AdditionalText
!CMcsUsbFunctionNet	HeadstageIDTypeObject, 610
CMcsUsbFunctionNet, 233	_ldType
!CMcsUsbListNet	HeadstageIDTypeObject, 610
CMcsUsbListNet, 240	\sim CCMOSMeaDeviceNet
!CMcsUsbNet	CCMOSMeaDeviceNet, 46
CMcsUsbNet, 246	\sim CChannelTestDeviceNet
!CMeFunctionNet	CChannelTestDeviceNet, 33
CMeFunctionNet, 298	\sim CCreateFilterNet
!CMeaCleanDeviceNet	CCreateFilterNet, 48
CMeaCleanDeviceNet, 270	\sim CDacCalibrationFunctionNet
!CMeaCoatDeviceNet	CDacCalibrationFunctionNet, 51
CMeaCoatDeviceNet, 274	\sim CDigOutStimulatorFunctionNet
!CMultiBatteryChargerDeviceNet	CDigOutStimulatorFunctionNet, 61
CMultiBatteryChargerDeviceNet, 300	~CExternDTesterDeviceNet
!CMultiwellCallbackFunctionNet	CExternDTesterDeviceNet, 65
CMultiwellCallbackFunctionNet, 307	~CFilterCoefficientsNet
!CMultiwellDeviceNet	CFilterCoefficientsNet, 67
CMultiwellDeviceNet, 309	~CFilterPropertyNet
!CMultiwellOptoStimFunctionNet	CFilterPropertyNet, 73
CMultiwellOptoStimFunctionNet, 315	~CFluidControlDeviceNet
!CPPCFunctionNet	CFluidControlDeviceNet, 75
CPPCFunctionNet, 349	~CGenericDevelopDeviceNet
!CPedoterDeviceNet	GGenericDevelopDeviceNet, 88
CPedoterDeviceNet, 330	~CGilsonDeviceNet
!CPositionIIDeviceNet	CGilsonDeviceNet, 99
CPositionIIDeviceNet, 336	~CGrapheneFunctionNet
!CPositionImpDeviceNet	CGrapheneFunctionNet, 102
CPositionImpDeviceNet, 344	~CInterfaceboard2FunctionNet
!CProgramPressureCurveNet	CInterfaceboard2FunctionNet, 118
CProgramPressureCurveNet, 363	~CInterfaceboardFunctionNet
!CPulseGeneratorFunctionNet	CInterfaceboardFunctionNet, 120
CPulseGeneratorFunctionNet, 364	~CLIH3DeviceNet
!CRFFunctionNet	CLIH3DeviceNet, 123
CRFFunctionNet, 372	~CMEA2100x256FunctionNet
!CSCUFunctionNet	CMEA2100x256FunctionNet, 265
CSCUFunctionNet, 425	\sim CMcsBusNet
!CTEERFunctionNet	CMcsBusNet, 169
CTEERFunctionNet, 514	\sim CMcsBus_AxisParametersNet
!CUsbDeviceConfigurationFunctionNet	CMcsBus_AxisParametersNet, 131

\sim CMcsBus_ExtensionNet	\sim CPgaDeviceNet
CMcsBus_ExtensionNet, 132	CPgaDeviceNet, 332
~CMcsBus_FYIExtensionNet	\sim CPositionIIDeviceNet
CMcsBus_FYIExtensionNet, 133	CPositionIIDeviceNet, 336
~CMcsBus_MotorControlNet	\sim CPositionImpDeviceNet
CMcsBus_MotorControlNet, 138	CPositionImpDeviceNet, 344
~CMcsBus_SensorNet	\sim CProgramPressureCurveNet
CMcsBus_SensorNet, 154	CProgramPressureCurveNet, 362
\sim CMcsBus_TempSensorNet	\sim CPulseGeneratorFunctionNet
CMcsBus_TempSensorNet, 163	CPulseGeneratorFunctionNet, 364
\sim CMcsBus_VoltageModeNet	\sim CRFFunctionNet
CMcsBus_VoltageModeNet, 165	CRFFunctionNet, 372
\sim CMcsUsbDacqNet	\sim CRetinaLedDeviceNet
CMcsUsbDacqNet, 179	CRetinaLedDeviceNet, 370
\sim CMcsUsbFactoryNet	\sim CRoboDeviceNet
CMcsUsbFactoryNet, 225	CRoboDeviceNet, 396
\sim CMcsUsbFunctionNet	\sim CRoboFluidDeviceNet
CMcsUsbFunctionNet, 233	CRoboFluidDeviceNet, 408
\sim CMcsUsbListEntryNet	\sim CSCUFunctionNet
CMcsUsbListEntryNet, 234	CSCUFunctionNet, 425
\sim CMcsUsbListNet	\sim CSafeISDeviceNet
CMcsUsbListNet, 240	CSafeISDeviceNet, 420
\sim CMcsUsbNet	~CStg200xBasicNet
CMcsUsbNet, 246	CStg200xBasicNet, 443
\sim CMeFunctionNet	\sim CStg200xDownloadNet
CMeFunctionNet, 298	CStg200xDownloadNet, 479
\sim CMeaCleanDeviceNet	~CSw2to64DeviceNet
CMeaCleanDeviceNet, 269	CSw2to64DeviceNet, 497
\sim CMeaCoatDeviceNet	\sim CTEERFunctionNet
CMeaCoatDeviceNet, 274	CTEERFunctionNet, 514
\sim CMeaDeviceNet	\sim CTEERMachineDeviceNet
CMeaDeviceNet, 280	CTEERMachineDeviceNet, 523
\sim CMealmpedanceDeviceNet	\sim CTcxDeviceNet
CMealmpedanceDeviceNet, 292	CTcxDeviceNet, 501
\sim CMeaSwitchDeviceNet	\sim CUsbDeviceConfigurationFunctionNet
CMeaSwitchDeviceNet, 295	CUsbDeviceConfigurationFunctionNet, 524
\sim CMeaUSBDeviceNet	\sim CWarnerUssingDeviceNet
CMeaUSBDeviceNet, 297	CWarnerUssingDeviceNet, 543
~CMultiBatteryChargerDeviceNet	~CWarnerUssingFunctionNet
CMultiBatteryChargerDeviceNet, 300	CWarnerUssingFunctionNet, 546
\sim CMultiwellCallbackFunctionNet	\sim CWarnerValveControllerDeviceNet
CMultiwellCallbackFunctionNet, 307	CWarnerValveControllerDeviceNet, 565
\sim CMultiwellDeviceNet	\sim CWarnerValveControllerDeviceTesterFunctionNet
CMultiwellDeviceNet, 309	CWarnerValveControllerDeviceTesterFunctionNet,
\sim CMultiwellOptoStimFunctionNet	586
CMultiwellOptoStimFunctionNet, 315	\sim DriverVersionNet
~CNF_GenDeviceNet	DriverVersionNet, 599
CNF_GenDeviceNet, 319	\sim SidebandData
~COkuvisionStimulatorDeviceNet	CStimulusFunctionNet::SidebandData, 624
COkuvisionStimulatorDeviceNet, 325	\sim StimulusDeviceDataAndUnrolledData
~CPPCFunctionNet	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData,
CPPCFunctionNet, 349	625
\sim CPathIdentDeviceNet	
CPathIdentDeviceNet, 329	A
~CPedoterDeviceNet	CFilterCoefficientsNet, 68
CPedoterDeviceNet, 330	AdditionalText
~CPeristalticPumpDeviceNet	HeadstageIDTypeObject, 610
CPeristalticPumpDeviceNet, 331	AddLoopEntry
1	CRetinal edDeviceNet 370

AddSelectedChannelsQueue	CW2100_StimulatorFunctionNet, 541
CMcsUsbDacqNet, 179–181	Bootstrap
AddSoftwareKey	FirmwareDestinationNames, 603
CMcsUsbNet, 246	bRequest
AddTableEntry	usbSetupPacket_t, 626
CRetinaLedDeviceNet, 370	BurnAdcOffset
Altera	COctoPotDeviceNet, 321
FirmwareDestinationNames, 603	BurnDacOffset
amplification	CDacCalibrationFunctionNet, 51
CMeaAudioFunctionNet::s_setaudionet, 623	COctoPotDeviceNet, 321
CW2100_FunctionNet::AudioChannelsNet, 29	BUS1_MCSBUS1
AmplifierSettle	FirmwareDestinationNames, 603
CIntanMea_FunctionNet, 116	BUS1_MCSBUS2
AnalogGain	FirmwareDestinationNames, 604
CMeaDeviceNet, 285	BusType
ApplyGains	DeviceIdNet, 595
CPgaDeviceNet, 333	ButterworthFilterHighPassNet, 31
AreTransistorVoltagesSet	ButterworthFilterHighPassNet, 32
CCMOSMea_FunctionNet, 36	ButterworthFilterLowPassNet, 32
AssociateToThis	ButterworthFilterLowPassNet, 32
CMcsUsbNet, 246	
AutomaticAnalogOut	CalibrateThermocouple
CSCUFunctionNet, 425	CFluidControlDeviceNet, 75
Axes_I	CTcxDeviceNet, 501
CRoboDeviceNet, 403	CancelInternalCalibration
Axes_X	CTEERFunctionNet, 514
CRoboDeviceNet, 403	CancelPoolLoop
Axes Y	CRoboDeviceNet, 396
CRoboDeviceNet, 403	CancelPoolLoopAndStopMovement
Axes Z	CRoboDeviceNet, 396
CRoboDeviceNet, 403	CancelTableLoop
Axis_I	CRoboDacqNet, 382
CRoboDeviceNet, 403	CancelTableLoopAndStopTable
Axis X	CRoboDacqNet, 382
CRoboDeviceNet, 403	CapacityTest
Axis Y	CMultiBatteryChargerDeviceNet, 300
CRoboDeviceNet, 403	CatchAmpGetAdcMean
Axis_Z	CMcsBus_SensorNet, 154
CRoboDeviceNet, 403	CatchAmpGetAdcValue
	CMcsBus_SensorNet, 154
В	CatchAmpGetAdcValueH
CFilterCoefficientsNet, 68	CMcsBus_SensorNet, 154
BatteryState, 29	CatchAmpGetAdcValueL
Charge, 30	CMcsBus_SensorNet, 154
ChargeRegionString, 30	CatchAmpGetDacAmplitude
ChargeString, 30	CMcsBus_SensorNet, 155
Voltage, 30	CatchAmpGetDacEnable
VoltageString, 30	CMcsBus_SensorNet, 155
BcdDevice	CatchAmpGetDacOffset
DeviceIdNet, 595	CMcsBus_SensorNet, 155
BeginImpedanceCheck	CatchAmpGetPwmEnable
CIntanMea_FunctionNet, 116	CMcsBus_SensorNet, 155
BesselFilterHighPassNet, 30	CatchAmpSetDacAmplitude
BesselFilterHighPassNet, 30	CMcsBus_SensorNet, 155
BesselFilterLowPassNet, 31	CatchAmpSetDacEnable
BesselFilterLowPassNet, 31	CMcsBus_SensorNet, 155
bmRequestType	CatchAmpSetDacOffset
usbSetupPacket_t, 626	CMcsBus_SensorNet, 155
BOOST_BIT	CatchAmpSetPwmEnable

OM D O N 1 455	0.10
CMcsBus_SensorNet, 155	SetSourceBulk, 44
CChannelTestDeviceNet, 33	SetSourceDrain, 44
~CChannelTestDeviceNet, 33	SetSourceGate, 44
CChannelTestDeviceNet, 33	SetStimulusSites, 44
SetAmplitude, 33	UpdateTransistorVoltages, 44
SetAttenuation, 33	VOPSTimerSetResetTimes, 44, 45
SetFrequency, 33	CCMOSMeaDeviceNet, 45
SetWaveform, 34	\sim CCMOSMeaDeviceNet, 46
CCMOSMea_FunctionNet, 34	CCMOSMeaDeviceNet, 46
AreTransistorVoltagesSet, 36	CMosMea, 47
CCMOSMea_FunctionNet, 36	GetAvailableBaseSamplerates, 46
ClearSTGOutput, 36	GetBaseSamplerate, 46
DetectChipType, 36	GetChannelDatal16, 46
EnableChannelsInGroup, 36, 37	GetChannelDatal32, 46
GetADCInputOffset, 37	GetChannelDataUI16, 46
GetBath, 37	GetChannelDataUI32, 47
GetBathMode, 37	GetCMOSDataDictionary, 47
GetEnabledChannelsInGroup, 37	SetBaseSamplerate, 47
GetGate, 37	SetRegionOfInterests, 47
GetGNDI, 37	Stimulus, 47
GetGroupADCBits, 37	UpdateChannelBlock, 47
GetGroupChannelBitmaskBySelect, 38	CCMOSMeaDeviceNet::CRegionOfInterestRect, 368
GetGroupChannelBitmaskHS1NCBathCurrent, 38	CRegionOfInterestRect, 368
GetGroupChannelBitmaskHS1NCCol2Current, 38	DeepCopy, 368
GetGroupChannelBitmaskHS1NChipTemp, 38	m_Bottom, 369
GetGroupChannelBitmaskHS1Sidebands, 39	m_Left, 369
GetGroupChannelBitmaskHS1TriggerStatus, 39	m_Right, 369
GetGroupChannelBitmaskIFDigChannels, 39	m_Top, 369
GetGroupChannelBitmaskInterfaceADC, 39	CCreateFilterNet, 48
GetGroupChannelBitmaskPacketFrameContext, 40	~CCreateFilterNet, 48
GetGroupChannelBitmaskSTG1DACSignal, 40	CCreateFilterNet, 48, 49
GetGroupDCOffset, 40	CutoffFrequency, 49
GetGroupID, 40	FindFilter, 49
•	
GetGroupNumberOfChannels, 41 GetGroupResolutionPerDigit, 41	GetBiQuad, 49
•	GetBiQuads, 49
GetGroupSampleSize, 41	NumCoefSets, 49
GetGroupType, 41	Order, 49
GetGroupUnit, 41, 42	SampleRate, 50
GetMaxNumOfColumns, 42	Scale, 50
GetNeurochipMemoryData, 42	CDacCalibrationFunctionNet, 50
GetNeurochipMemorySize, 42	!CDacCalibrationFunctionNet, 51
GetNumberOfSupportedGroups, 42	~CDacCalibrationFunctionNet, 51
GetSourceBulk, 42	BurnDacOffset, 51
GetSourceDrain, 42	CDacCalibrationFunctionNet, 50, 51
GetSourceGate, 42	GetDacOffset, 51
GetStimulusSites, 43	SetDacOffset, 51
GetVDD3I, 43	CDacqGroupChannelGenericSelectionNet, 52
GetVDDI, 43	CDacqGroupChannelGenericSelectionNet, 52
IsChipPowered, 43	CDacqGroupChannelSelectionNet, 52
IsGateFloating, 43	CDacqGroupChannelSelectionNet, 53
PowerChip, 43	CDacqGroupChannelSelectionTemplateNet
SetADCInputOffset, 43	CDacqGroupChannelSelectionTemplateNet< Dac-
SetBath, 43	qGroupChannelEnumTemplateNet, Dac-
SetBathMode, 43	qGroupChannelEnumTemplate, CDevice-
SetGate, 43	GroupChannelInfoTemplateNet >, 54
SetGateFloating, 44	CDacqGroupChannelSelectionTemplateNet< Dac-
SetGateToVOP, 44	qGroupChannelEnumTemplateNet, Dac-
SetNeurochipMemoryData, 44	qGroupChannelEnumTemplate, CDevice-

GroupChannelInfoTemplateNet >, 53	CFilterCoefficientsNet, 67
CDacqGroupChannelSelectionTemplateNet, 54	GetUintA, 68
EnableChannelsInGroup, 54	GetUintB, 68
GetDeviceGroupChannelInfos, 54	IsEqual, 68
GetEnabledChannelsInGroup, 54	CFilterCoefficientsNet::s_FilterAttributesNet, 621
GetGroupID, 54, 55	CommaPositionA, 622
GetGroupNumberOfChannels, 55	CommaPositionB, 622
GetGroupSampleSize, 55	PostCommaA, 622
GetGroupType, 55	PostCommaB, 622
GetNumberOfSupportedGroups, 55	PreCommaA, 622
CDeviceGroupChannelInfoGenericNet, 56	PreCommaB, 623
CDeviceGroupChannelInfoGenericNet, 56	s_FilterAttributesNet, 622
CDeviceGroupChannelInfoMEA2100_256Net, 56	ToCpp, 622
CDeviceGroupChannelInfoMEA2100_256Net, 57	CFilterConfigurationNet, 68
CDeviceGroupChannelInfoNet, 57	CFilterConfigurationNet, 69
CDeviceGroupChannelInfoNet, 57	EraseFilterParameterPermanent, 69
CDeviceGroupChannelInfoSCUNet, 57	GetFilterAttributes, 69
CDeviceGroupChannelInfoSCUNet, 58	GetHighpassFilterEnable, 69
CDeviceGroupChannelInfoTemplateNet	ResetHighpassFilter, 70
CDeviceGroupChannelInfoTemplateNet< Dacq-	SetFilterParameter, 70
GroupChannelEnumTemplateNet >, 58	SetFilterParameterPermanent, 70
CDeviceGroupChannelInfoTemplateNet Dacq-	SetHighpassFilterEnable, 70
GroupChannelEnumTemplateNet >, 58	CFilterConfigurationRegisterNet, 70
CDeviceGroupChannelInfoTemplateNet, 58	CFilterConfigurationRegisterNet, 70
·	
GroupID, 58	EraseFilterParameterPermanent, 71
GroupType, 58	SetFilterParameter, 71, 72
NumberOfChannels, 59	SetFilterParameterPermanent, 72
CDeviceGroupChannelInfoW2100Net, 59	CFilterPropertyNet, 72
CDeviceGroupChannelInfoW2100Net, 59	~CFilterPropertyNet, 73
CDigOutStimulatorFunctionNet, 59	CFilterPropertyNet, 73
!CDigOutStimulatorFunctionNet, 61	CornerFrequency, 73
~CDigOutStimulatorFunctionNet, 61	CornerFrequencymHz, 73
CDigOutStimulatorFunctionNet, 60	FilterActive, 73
ClearChannel, 61	FilterBand, 73
GetGlobalRepeat, 61	FilterFamily, 73
GetNumberOfChannels, 61	FilterType, 73
GetStartTriggerSlope, 61	Order, 74
GetStopTriggerSlope, 62	ToString, 73
PrepareChannelData, 62	CFluidControlDeviceNet, 74
SendChannelData, 62	~CFluidControlDeviceNet, 75
SetGlobalRepeat, 63	CalibrateThermocouple, 75
SetStartTriggerSlope, 63	CFluidControlDeviceNet, 75
SetStopTriggerSlope, 63	GetAdc, 76
CEncapsulatorDeviceNet, 64	GetDigin, 76
CEncapsulatorDeviceNet, 64	GetDigout, 76
GetRoboFluidDevice, 64	GetPWM, 76
CExternDTesterDeviceNet, 64	GetReferenceTemperature, 76
!CExternDTesterDeviceNet, 65	GetSingleValve, 77
~CExternDTesterDeviceNet, 65	GetThermocoupleCalibration, 77
CExternDTesterDeviceNet, 65	GetThermocoupleNanovoltPerKelvin, 77
Read, 65	GetThermocoupleTemperature, 78
Read2, 66	GetValve, 78
Write, 66	McsBus_VoltageMode, 80
Write2, 66	SetDigout, 78
CFilterCoefficientsNet, 66	SetPWM, 78
~CFilterCoefficientsNet, 67	SetSingleValve, 79
A, 68	SetThermocoupleNanovoltPerKelvin, 79
B, 68	SetValve, 79

CFYIDeviceNet, 80	CW2100_FunctionNet::AudioChannelsNet, 29
CFYIDeviceNet, 80	ChannelBlock_AvailFrames
FYIProgram, 81	CMcsUsbDacqNet, 182
FYITemp, 81	ChannelBlock_ReadAsFrameArrayI16
Sensor, 81	CMcsUsbDacqNet, 182, 183
CGenericDevelopDeviceNet, 81	ChannelBlock_ReadAsFrameArrayl32
∼CGenericDevelopDeviceNet, 88	CMcsUsbDacqNet, 184
CGenericDevelopDeviceNet, 88	ChannelBlock_ReadAsFrameArrayUI16
ClosePipe, 88	CMcsUsbDacqNet, 185, 186
GetBuffer, 88	ChannelBlock_ReadAsFrameArrayUI32
GetByteBuffer, 89	CMcsUsbDacqNet, 186, 187
GetIntBuffer, 89	ChannelBlock_ReadFramesDictI16
GetShortBuffer, 90	CMcsUsbDacqNet, 188
GetUByteBuffer, 91	ChannelBlock_ReadFramesDictl32
GetUIntBuffer, 91	CMcsUsbDacqNet, 188
GetUShortBuffer, 92	ChannelBlock_ReadFramesDictUI16
OpenPipe, 93	CMcsUsbDacqNet, 189
ReadPipe, 93	ChannelBlock_ReadFramesDictUI32
ResetPipe, 93	CMcsUsbDacqNet, 189
SetBuffer, 94	ChannelBlock_ReadFramesI16
SetByteBuffer, 94	CMcsUsbDacqNet, 190, 191
SetIntBuffer, 94	ChannelBlock_ReadFramesl32
SetShortBuffer, 95	CMcsUsbDacqNet, 191, 192
SetUByteBuffer, 96	ChannelBlock_ReadFramesUI16
SetUIntBuffer, 96	CMcsUsbDacqNet, 193
SetUShortBuffer, 97	ChannelBlock_ReadFramesUI32
SetValue, 97	CMcsUsbDacqNet, 194, 195
WritePipe, 98	ChannelDataEvent
CGilsonDeviceNet, 98	CMcsUsbDacqNet, 221
\sim CGilsonDeviceNet, 99	ChannelReset
CGilsonDeviceNet, 99	CMultiBatteryChargerDeviceNet, 301
ConnectSlave, 99	Charge
GetLastAnswer, 100	BatteryState, 30
m_pGilsonDevice, 100	ChargeRegionString
SendBuffered, 100	BatteryState, 30
SendImmediate, 100	ChargeString
SendImmediateGetResponse, 100	BatteryState, 30
CGrapheneFunctionNet, 100	CHiClampDeviceNet, 109
!CGrapheneFunctionNet, 102	CHiClampDeviceNet, 110
\sim CGrapheneFunctionNet, 102	RoboDacq, 110
CGrapheneFunctionNet, 102	CHLADacqNet, 110
GetDACOffset, 102, 103	CHLADacqNet, 111
GetVdVs, 103	CHLADeviceNet, 111
GetVdVsDAC, 103, 104	CHLADeviceNet, 111
GetVoltageRange, 104	HLADacq, 112
GetVoltageReached, 105	SerialPort, 112
GetVoltageResolution, 105	CHWInfo
SetDACOffset, 106	CMcsUsbDacqNet::CHWInfo, 113
SetVds, 106	CIntanMea_FunctionNet, 115
SetVdVs, 107	AmplifierSettle, 116
SetVdVsDAC, 107	BeginImpedanceCheck, 116
SetVgs, 108	CIntanMea_FunctionNet, 115
SetVoltageRange, 108	GetDSPHighPassByIndex, 116
SetVoltageResolution, 109	GetImpedanceResult, 116
ChangeSerialNumber	GetIntanRegister, 116
CMcsUsbFactoryNet, 225	GetLowerFrequencyByIndex, 116
channel	GetUpperFrequencyByIndex, 116
CMeaAudioFunctionNet::s_setaudionet, 623	SetBandwidthByIndex, 116

CatDiagnosticMade 116	CatDagOffact 10E
SetDiagnosticMode, 116	GetDacOffset, 125
SetDSPHighPassByIndex, 117	GetDacqRunStatus, 125
SetIntanRegister, 117	GetDacUseIdleValue, 125
CInterfaceboard2FunctionNet, 117	GetDigInState, 126
!CInterfaceboard2FunctionNet, 118	GetEEpromPage, 126
\sim CInterfaceboard2FunctionNet, 118	GetSampleInterval, 126
CInterfaceboard2FunctionNet, 118	IsUserTriggerEnabled, 126
GetloVoltage, 118	ReadClipping, 127
SetloVoltage, 118	ReadUARTData, 127
CInterfaceboardFunctionNet, 119	SendCommand, 127
!CInterfaceboardFunctionNet, 120	SetAdcOffset, 127
~CInterfaceboardFunctionNet, 120	SetAdcOffsetPermanent, 128
CInterfaceboardFunctionNet, 119	SetAudioOutDacParameter, 128
GetCardinalDacqSamplerate, 120	SetDacIdleValue, 128
GetCardinalStgOutputrate, 120	SetDacoffset, 128
- · · · · · · · · · · · · · · · · · · ·	
SetCardinalDacqSamplerate, 120	SetDacOffsetPermanent, 129
SetCardinalStgOutputrate, 120	SetDacUseIdleValue, 129
ClampAmpRestart	SetDigOutState, 129
CRoboDacqNet, 382	SetEEpromPage, 129
ClearBuffers	SetSampleInterval, 130
CMcsUsbDacqNet, 195	StimulusFunction, 130
ClearChannel	WriteUARTData, 130
CDigOutStimulatorFunctionNet, 61	CloseAllValves
ClearChannel_PrepareAndSendData	CRoboFluidDeviceNet, 408
CStg200xDownloadNet, 480	ClosePipe
CStimulusFunctionNet, 487	CGenericDevelopDeviceNet, 88
ClearChannelData	ClosePlateClamp
CStg200xDownloadBasicNet, 471	CMultiwellDeviceNet, 310
CStimulusFunctionNet, 487	CMcsBus_AxisParametersNet, 130
CW2100_StimulatorFunctionNet, 536	~CMcsBus_AxisParametersNet, 131
ClearMultiplexedData	CMcsBus_AxisParametersNet, 131
CStimulusFunctionNet, 487	GetAxisParametersSignedEeprom, 131
ClearSTGOutput	GetAxisParametersUnsignedEeprom, 131
CCMOSMea FunctionNet, 36	SetAxisParametersEeprom, 131, 132
ClearStimulusParametersCache	•
	CMcsBus_ExtensionNet, 132
CW2100_FunctionNet, 529	~CMcsBus_ExtensionNet, 132
ClearSyncData	CMcsBus_ExtensionNet, 132
CStg200xDownloadBasicNet, 471	GetLEDSwitch, 133
CStimulusFunctionNet, 487	SetLEDSwitch, 133
ClearTable	CMcsBus_FYIExtensionNet, 133
CRetinaLedDeviceNet, 370	∼CMcsBus_FYIExtensionNet, 133
ClearTableName	CMcsBus_FYIExtensionNet, 133
CWarnerValveControllerDeviceNet, 565	GetDIO, 134
ClearUserDefinedNameCache	GetSingleHeater, 134
CW2100_FunctionNet, 529	GetValves, 134
ClearValveTable	SetDIO, 134
CWarnerValveControllerDeviceNet, 565	SetSingleHeater, 134
CLIH3DeviceNet, 121	SetValves, 134
!CLIH3DeviceNet, 123	CMcsBus MotorControlNet, 135
~CLIH3DeviceNet, 123	~CMcsBus_MotorControlNet, 138
CLIH3DeviceNet, 123	CMcsBus_MotorControlNet, 138
DummyCommand, 123	GetMCAcceleration, 138
-	
EnableUserTrigger, 123	GetMCAccelerationEeprom, 138
ErasePermanentAdcOffset, 124	GetMCAccelerationShortCommand, 138
ErasePermanentDacOffset, 124	GetMCAxisRevisionEeprom, 138
GetAdcOffset, 124	GetMCBreakCurrent, 138
GetAudioOutDacParameter, 124	GetMCBreakCurrentEeprom, 139
GetDacIdleValue, 125	GetMCConfig, 139

CatMCCanfigEngram 100	CatMCMayCacad 140
GetMCConfigEeprom, 139	SetMCMaxSpeed, 148
GetMCCurrent, 139	SetMCMaxSpeedEeprom, 148
GetMCCurrentEeprom, 139	SetMCMaxTravel, 148
GetMCCurrentMode, 139	SetMCMaxTravelEeprom, 148
GetMCCurrentModeEeprom, 139	SetMCMaxTravelShortCommand, 148
GetMCCurrentModeShortCommand, 140	SetMCNewPosition, 148
GetMCCurrentPosition, 140	SetMCOutputOnOff, 149
GetMCCurrentShortCommand, 140	SetMCReference, 149
GetMCCurrentSpeed, 140	SetMCReferenceCurrent, 149
GetMCMaxAcceleration, 140	SetMCReferenceCurrentEeprom, 149
GetMCMaxAccelerationEeprom, 140	SetMCRegulatorGain, 149
GetMCMaxCurrent, 140	SetMCRegulatorGainEeprom, 149
GetMCMaxCurrentEeprom, 141	SetMCRotation, 150
GetMCMaxSpeed, 141	SetMCScalingFactor, 150
GetMCMaxSpeedEeprom, 141	SetMCScalingFactorEeprom, 150
GetMCMaxTravel, 141	SetMCSpeed, 150
GetMCMaxTravelEeprom, 141	SetMCSpeedEeprom, 150
GetMCMaxTravelShortCommand, 141	SetMCSpeedShortCommand, 150
GetMCMovement, 141	SetMCSpeedUnitEeprom, 151
GetMCNewPosition, 142	· · · · · · · · · · · · · · · · · · ·
	SetMCStandbyCurrent, 151
GetMCOutputOnOff, 142	SetMCStandbyCurrentEeprom, 151
GetMCPhase, 142	SetMCStandbyTime, 151
GetMCPhaseOffset, 142	SetMCStandbyTimeEeprom, 151
GetMCReference, 142	SetSubChannel, 151
GetMCReferenceCurrent, 142	StartMCMovement, 152
GetMCReferenceCurrentEeprom, 142	StopMCMovement, 152
GetMCRegulatorGain, 143	CMcsBus_SensorNet, 152
GetMCRegulatorGainEeprom, 143	\sim CMcsBus_SensorNet, 154
GetMCScalingFactor, 143	CatchAmpGetAdcMean, 154
GetMCScalingFactorEeprom, 143	CatchAmpGetAdcValue, 154
GetMCSpeed, 143	CatchAmpGetAdcValueH, 154
GetMCSpeedEeprom, 143	CatchAmpGetAdcValueL, 154
GetMCSpeedShortCommand, 143	CatchAmpGetDacAmplitude, 155
GetMCSpeedUnitEeprom, 144	CatchAmpGetDacEnable, 155
GetMCStandbyCurrent, 144	CatchAmpGetDacOffset, 155
GetMCStandbyCurrentEeprom, 144	CatchAmpGetPwmEnable, 155
GetMCStandbyTime, 144	CatchAmpSetDacAmplitude, 155
GetMCStandbyTimeEeprom, 144	CatchAmpSetDacEnable, 155
GetSubChannel, 144	CatchAmpSetDacOffset, 155
SetMCAcceleration, 144	CatchAmpSetPwmEnable, 155
SetMCAccelerationEeprom, 145	CMcsBus SensorNet, 154
SetMCAccelerationShortCommand, 145	Get2AnalogInput, 156
SetMCAxisRevisionEeprom, 145	Get2DigitalInput, 156
SetMCBreakCurrent, 145	Get4ADC, 156
SetMCBreakCurrentEeprom, 145	Get4ADCAverage, 156
SetMCConfig, 145	Get4ADCCatchampAverageShift, 156
SetMCConfigEeprom, 146	Get4ADCMode, 156
SetMCCurrent, 146	Get4DAC, 156
SetMCCurrentEeprom, 146	GetADCs, 156
SetMCCurrentMode, 146	GetADCsLoop, 157
SetMCCurrentModeEeprom, 146	GetBubbleStatus, 157
SetMCCurrentModeShortCommand, 146	GetDACs, 157
SetMCCurrentPosition, 147	GetDetectionThreshold, 157
SetMCCurrentShortCommand, 147	GetDetectorValue, 157
SetMCMaxAcceleration, 147	GetLatency, 157
SetMCMaxAccelerationEeprom, 147	GetLatencyCounter, 157
SetMCMaxCurrent, 147	GetMinimalThreshold, 157
SetMCMaxCurrentEeprom, 147	GetMovePump, 158

GetPiezoState, 158	SetVMMaxPositiveCurrentEeprom, 168
GetPressure, 158	SetVMMaxPositiveVoltage, 168
GetPressureOffset, 158	SetVMMaxPositiveVoltageEeprom, 168
GetRegulationTimeouts, 158	SetVMOutputOnOff, 168
GetRegulatorFactor, 159	SetVMVoltage, 168
GetRegulatorOnOff, 159	CMcsBusNet, 169
GetRegulatorStatus, 159	~CMcsBusNet, 169
GetRotatePump, 159	CMcsBusNet, 169
GetSamplePeriode, 159	CMcsBusNet::GetMode, 170
GetSollPressure, 159	CMcsBusNet::GetModeEeprom, 170
GetSyncState, 159	CMcsBusNet::SetMode, 170
Set4ADCCatchampAverageShift, 160	CMcsBusNet::SetModeEeprom, 170
Set4ADCMode, 160	GetBusAddress, 170
Set4DAC, 160	GetBusAddressEeprom, 170
SetDACs, 160	GetCommand, 170, 171
SetDetectionThreshold, 160	GetHWRevisionEeprom, 171
SetLatency, 160	SetBusAddress, 171
SetMinimalThreshold, 160	SetBusAddressEeprom, 171
SetMovePump, 161	SetCommand, 172
SetPiezoState, 161	SetHWRevisionEeprom, 172
SetPressureOffset, 161	CMcsBusNet::GetMode
SetRegulationTimeouts, 161	CMcsBusNet, 170
SetRegulatorFactor, 161	CMcsBusNet::GetModeEeprom
SetRegulatorOnOff, 161	CMcsBusNet, 170
SetRotatePump, 161	CMcsBusNet::SetMode
SetSamplePeriode, 162	CMcsBusNet, 170
SetSollPressure, 162	CMcsBusNet::SetModeEeprom
StartSync, 162	CMcsBusNet, 170
TactSwitchGetState, 162	CMcsUsbDacqNet, 173
TactSwitchSetDisplay, 162	~CMcsUsbDacqNet, 179
CMcsBus_TempSensorNet, 162	AddSelectedChannelsQueue, 179–181
~CMcsBus_TempSensorNet, 163	ChannelBlock_AvailFrames, 182
CMcsBus TempSensorNet, 163	ChannelBlock ReadAsFrameArrayI16, 182, 183
GetNanoVoltsPerKelvin, 163	ChannelBlock ReadAsFrameArray132, 184
GetTemperatur, 163	ChannelBlock_ReadAsFrameArrayUI16, 185, 186
GetThermoOffset, 163	ChannelBlock ReadAsFrameArrayUI32, 186, 187
GetThermoTemp, 164	ChannelBlock ReadFramesDict116, 188
GetThermoVoltage, 164	ChannelBlock_ReadFramesDict132, 188
SetThormo Officet 164	ChannelBlock_ReadFramesDictUI16, 189
SetThermoOffset, 164	ChannelBlock_ReadFramesDictUI32, 189
CMcsBus_VoltageModeNet, 164	ChannelBlock_ReadFramesI16, 190, 191
~CMcsBus_VoltageModeNet, 165	ChannelBlock_ReadFramesI32, 191, 192 ChannelBlock_ReadFramesUI16, 193
CMcsBus_VoltageModeNet, 165 GetVMMaxNegativeCurrent, 166	ChannelBlock_ReadFramesUl32, 194, 195
•	
GetVMMaxNegativeVoltage 166	ChannelDataEvent, 221 ClearBuffers, 195
GetVMMaxNegativeVoltage, 166	
GetVMMaxNegativeVoltageEeprom, 166	CMcsUsbDacqNet, 179
GetVMMaxPositiveCurrent, 166	CMcsUsbDacqNet::GetFilterProperties, 195
GetVMMaxPositiveCurrentEeprom, 166	Error_Callback_Aquisition_Stopped, 220
GetVMMaxPositiveVoltage, 166	Error_Callback_Data_lost, 220
GetVMOxtoutOnOff 167	Error_Callback_Frames_Lost, 220
GetVMOutputOnOff, 167	Error_Callback_Packet_Error, 220
GetVMVoltage, 167	Error_Callback_Queue_Full, 220
SetVMMaxNegativeCurrent, 167	Error_Callback_RingQueue_Full, 220
SetVMMaxNegativeCurrentEeprom, 167	ErrorEvent, 221
SetVMMaxNegativeVoltage, 167	GetAdapterType, 196
SetVMMaxNegativeVoltageEeprom, 167	GetAdcDataFormat, 196
SetVMMaxPositiveCurrent, 168	GetAdcZero, 196

GetAnalogValueUnit, 196	VoltageRangeInMicroVolt, 526
GetChannelDataFillSize, 196	CMcsUsbDacqNet::GetFilterProperties
GetChannelLayout, 196	CMcsUsbDacqNet, 195
GetChannelsInBlock, 197	CMcsUsbDeviceStatePushFunctionNet, 221
GetDataFormat, 197	CMcsUsbDeviceStatePushFunctionNet, 222
	McsUsbDeviceStateEvent, 222
GetDataMode, 197	•
GetDigitalSource, 197–199	TriggerStatus, 222
GetFilterProperty, 199	CMcsUsbDeviceStatePushNet, 222
GetGroupChannelDatal16, 199	CMcsUsbDeviceStatePushNet, 223
GetGroupChannelDatal32, 200	McsUsbDeviceStateEvent, 223
GetGroupChannelDataUI16, 200	TriggerStatus, 223
GetGroupChannelDataUI32, 201	CMcsUsbFactoryNet, 223
GetHardwareMaxRange, 202	~CMcsUsbFactoryNet, 225
GetHardwareMinRange, 202	ChangeSerialNumber, 225
GetMaxSamplingFrequency, 202	CMcsUsbFactoryNet, 225
GetMeaLayout, 202	Coldstart, 225
GetMinSamplingFrequencyStepsize, 202	CompareFirmware, 226
GetNumberOfDataBits, 203	DownloadFirmware, 226
GetPoti, 203	FindFirmwareVersionMagicInBuffer, 226
GetResolutionPerDigit, 203	FX3MCSDataAddress, 231
GetSamplerate, 203	FX3MCSDataDeviceIdOffset, 231
GetVoltageRangeIndex, 203	FX3MCSDatalFB1ImageOffset, 231
GetVoltageRangeInMicroVolt, 203	FX3MCSDatalFB2ImageOffset, 231
GetVoltageRangeInMilliVolt, 204	FX3MCSDataVersionOffset, 231
HWInfo, 204	GetChecksumFromFX3Image, 226
Samplerate, 221	GetDestination, 226
SendStartDacq, 204	GetDestinationDisplayLabel, 226
SendStartStgAndDacq, 204	GetDestinationName, 226, 227
SendStopDacq, 205	GetDestinationSerialNumber, 227
SendStopStgAndDacq, 205	GetDestinationTargetAddress, 227
SendStopStgAndDacqWithOptions, 205	GetFirmwareVersionFromFile, 227
SetDataMode, 206	GetFirmwareVersionFromHexFile, 227
SetDigitalSource, 206–208	GetNumDestinations, 227
SetPoti, 208	GetUSBDeviceIDFromFX3Image, 228
SetSamplerate, 208	GetUsercodeFromBitFile, 228
SetSelectedChannels, 208–210	GetUsercodeFromFlash, 228
SetSelectedChannelsQueue, 211–213	GetXilinxFlashOffset, 228
SetSelectedData, 213–215	GetXilinxFlashReadCommand, 228
SetupGroupDacqQueue, 215	LoadUserFirmware, 228, 229
SetVoltageRangeByIndex, 215	ReadBlockFromFlash, 229
SetVoltageRangeInMicroVolt, 215	ReadBlockFromIFBGlobalEEprom, 229
StartDacq, 216, 217	ReadBlockFromNVMEM, 229
StartLoop, 218, 219	SetDestinationSerialNumber, 229
StopDacq, 219, 220	UpdateFirmware, 229–231
StopLoop, 220	CMcsUsbFunctionNet, 232
CMcsUsbDacqNet::CHWInfo, 112	!CMcsUsbFunctionNet, 233
CHWInfo, 113	~CMcsUsbFunctionNet, 233
GetAvailableSampleRates, 113	CMcsUsbFunctionNet, 232, 233
	• • •
GetAvailableVoltageRangesInMicroVolt, 113	m_pMcsUsb, 233
GetAvailableVoltageRangesInMicroVoltAnd-	m_pMcsUsbFunction, 233
StringsInMilliVolt, 113	ThrowCUsbExceptionNetOnError, 233
GetNumberOfHWADCChannels, 114	CMcsUsbFunctionPointerContainer, 233
GetNumberOfHWDigitalChannels, 114	CMcsUsbListEntryNet, 233
IsDigitalChannelDedicated, 114	~CMcsUsbListEntryNet, 234
CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet,	Deviceld, 238
526	DeviceName, 238
CVoltageRangeInfoNet, 526	Equals, 235
VoltageRangeDisplayStringMilliVolt, 526	GetEntry, 236

	GetEntryCount, 237	IsConnected, 253
	HwVersion, 238	IsDeviceHighSpeed, 253
	Manufacturer, 238	IsDeviceHighSpeedCapable, 254
	Product, 238	IsExceptionsEnabled, 254
	SerialNumber, 238	MultibootGetCypressImageId, 254
	SetStringFormat, 237	MultibootGetImageId, 254
	ToString, 238	MultibootGetSelectedImage, 254
CMcs	sUsbListNet, 239	MultibootSelectImage, 254
	!CMcsUsbListNet, 240	ReadEepromRegisterPreconfig, 255
	\sim CMcsUsbListNet, 240	ReadRegister, 255
	CMcsUsbListNet, 239	ReadRegister32, 255
	Count, 241	ReadRegisterTimeSlot, 255
	DeviceArrival, 241	RemoveSoftwareKey, 256
	DeviceRemoval, 241	RescanHeadstage, 256
	GetNumberOfDevices, 240	SerialNumber, 263
	GetUsbListEntries, 240	SetConfiguration, 256
	GetUsbListEntry, 240	SetSoftwareKey, 256
	IsDeviceTypeOf, 241	Status_AlreadyConfigured, 258
	SetStringFormat, 241	Status_BadStartFrame, 259
CMcs	sUsbNet, 242	Status_Btstuff, 259
	!CMcsUsbNet, 246	Status_BufferOverrun, 259
	∼CMcsUsbNet, 246	Status_BufferUnderrun, 259
	AddSoftwareKey, 246	Status_Canceled, 259
	AssociateToThis, 246	Status_Canceling, 259
	CMcsUsbNet, 246	Status_ConnectedPipes, 259
	Connect, 246, 247	Status_ControlNotOwned, 259
	CyclePort, 248	Status_Crc, 259
	Disconnect, 248	Status_DataOverrun, 259
	EmptyKey, 248	Status_DataToggleMismatch, 259
	EnableExceptions, 248	Status_DataUnderrun, 260
	EraseEepromRegisterPreconfig, 248	Status_DeviceLocked, 260
	GetConfiguration, 249	Status_DeviceNotFound, 260
	GetDeviceCannotStallOutRequests, 249	Status DeviceRemoved, 260
	GetDeviceCapableSpeed, 249	Status_DevNotResponding, 260
	GetDeviceEnum, 249	Status_EndpointHalted, 260
	GetDeviceId, 249	Status_ErrorBusy, 260
	GetDeviceRootHubVendorEnum, 249	Status ErrorShortTransfer, 260
	GetDeviceRootHubVendorID, 249	Status_Fifo, 260
	GetDeviceRootHubVendorName, 249	Status FrameControlOwned, 260
	GetDeviceSpeed, 250	Status InternalHcError, 260
	GetErrorText, 250	Status InvalidDeviceHandle, 261
	GetFirmwareVersion, 250	Status_InvalidHandle, 261
	GetHardwareRevision, 250	Status Invalid Parameter, 261
	GetHeadstageActive, 251	Status_InvalidPipeHandle, 261
	GetHeadstageID, 251	Status InvalidUrbFunction, 261
	-	_
	GetHeadstagePresent, 251	Status_loPending, 261
	GetIdent, 251	Status_loads Degree Failed 001
	GetLastUSBError, 252	Status_IsochRequestFailed, 261
	GetMea21UsbPort, 252	Status_LastUsbErrorMismatch, 261
	GetNumConfigurations, 252	Status_NoBandwidth, 261
	GetSerialNumber, 252	Status_NoMemory, 262
	GetSoftwareKey, 252	Status_NoSuchDevice, 262
	GetSoftwareKeyString, 252	Status_NotAccessed, 262
	GetStatus, 252	Status_NotSupported, 262
	GetStatusOfLastCommand, 253	Status_PidCheckFailure, 262
	GetUsbListEntry, 253	Status_PipeNotLinked, 262
	GetVersion, 253	Status_RequestFailed, 262
	HasSoftwareKey, 253	Status_RequestMutexFailed, 262

0 D	0 10 11 074
Status_RequestMutexTimeout, 262	GetDuration, 274
Status_Stall, 262	GetMaxCurrent, 274
Status_Unconfigured, 262	GetOffsetCurrent, 275
Status_UnexpectedPid, 263	GetOutputCurrent, 275
ThrowCUsbExceptionNetOnError, 256	GetPauseDuration, 275
TxnGetSerialNumber, 256	GetSlope, 275
TxnSetSerialNumber, 256	GetTimeInPause, 275
TxnTestMemoryReadAndCheck, 256	GetTimeInPlateau, 276
TxnTestMemoryWrite, 257	IsRunning, 276
ValidKey, 257	SetCycles, 276
WPAError_ScanningIsPending, 263	SetDuration, 276
WriteEepromRegisterPreconfig, 257	SetMaxCurrent, 277
WriteRegister, 257, 258	SetOffsetCurrent, 277
WriteRegister32, 258	SetPauseDuration, 277
WriteRegisterArray, 258	SetSlope, 277
WriteRegisterTimeSlot, 258	Start, 277
WriteRegisterValue, 258	Stop, 278
CMcsUsbPointerContainer, 263	CMeaDeviceNet, 278
CMEA2100_256DacqGroupChannelSelectionNet, 263	~CMeaDeviceNet, 280
CMEA2100_256DacqGroupChannelSelectionNet,	AnalogGain, 285
263	CMeaDeviceNet, 279
CMEA2100x256FunctionNet, 264	EnableChecksum, 280
!CMEA2100x256FunctionNet, 265	EnableDigitalIn, 280, 281
~CMEA2100x256FunctionNet, 265	EnableTimestamp, 281
CMEA2100x256FunctionNet, 264	Gain, 285
GetLayoutConfiguration, 265	GetAnalogGain, 282
SetLayoutConfiguration, 265	GetEnumerationSpeed, 282
CMeaAudioFunctionNet, 265	GetGain, 282
CMeaAudioFunctionNet, 266	
	MeaAudioFunctionNet, 285
GetAudioChannels, 266, 267	MeaDigitalDataFunctionNet, 285
GetNumberOfAudioChannels, 267	MeaFeedbackFunctionNet, 285
SetAudioChannels, 267, 268	MeFunctionNet, 286
CMeaAudioFunctionNet::s_setaudionet, 623	SetDigitalOut, 282
amplification, 623	SetNumberOfAnalogChannels, 282
channel, 623	SetNumberOfChannels, 283, 284
CMeaCleanDeviceNet, 268	SetTriggerMaskValue, 284
!CMeaCleanDeviceNet, 270	SetTriggerPeriod, 285
~CMeaCleanDeviceNet, 269	W2100_FunctionNet, 286
CMeaCleanDeviceNet, 269	WClassicFunctionNet, 286
GetCycle, 270	CMeaDigitalDataFunctionNet, 286
GetCycles, 270	CMeaDigitalDataFunctionNet, 286
GetMaxVoltage, 270	GetDigitalData, 287
GetMinVoltage, 270	SetDigitalData, 287
GetOutputVoltage, 270	CMeaFeedbackFunctionNet, 288
GetSlope, 271	CMeaFeedbackFunctionNet, 289
IsRunning, 271	FeedbackGetSampleTimerCount, 289
SetCycles, 271	FeedbackSetAnalogSource, 289
SetMaxVoltage, 271	FeedbackSetChannelFilter, 289
SetMinVoltage, 272	FeedbackSetDigitalMapping, 289
SetSlope, 272	FeedbackSetFeedback, 289
Start, 272	FeedbackSetFilterOff, 290
Stop, 272	FeedbackSetFilterParameter, 290
CMeaCoatDeviceNet, 272	FeedbackSetFilterParameter32, 290
!CMeaCoatDeviceNet, 274	FeedbackSetGlobalChannelFilter, 290
~CMeaCoatDeviceNet, 274	FeedbackSetIIRFilterParameter, 290
CMeaCoatDeviceNet, 274	FeedbackSetLogic, 290
GetCurrentCycle, 274	FeedbackSetMkFilter, 290
GetCycles, 274	FeedbackSetNumberOfLogics, 291

FeedbackSetNumberOfRateCounter, 291	SetChargingPCoefficient, 304
FeedbackSetNumberOfRateDetectors, 291	SetDischargeCurrentSetPoint, 305
FeedbackSetNumberOfSpikeDetectors, 291	SetFinalDischargeVoltage, 305
FeedbackSetNumberOfTriggers, 291	SetRatedCapacity, 305
FeedbackSetRateCounter, 291	SetRatedCapacityVolatile, 305
FeedbackSetRateDetector, 291	CMultiwellCallbackFunctionNet, 306
FeedbackSetSpikeDetectorThreshold, 291	!CMultiwellCallbackFunctionNet, 307
FeedbackSetTrigger, 292	~CMultiwellCallbackFunctionNet, 307
CMealmpedanceDeviceNet, 292	CMultiwellCallbackFunctionNet, 307
~CMealmpedanceDeviceNet, 292	GetPlateClampStateByHeadstage, 307
CMealmpedanceDeviceNet, 292	GetPlateClampStateByHeadstageEvent, 308
GetAdapterCode, 293	OnGetPlateClampStateByHeadstage, 307
GetArraySize, 293	CMultiwellDeviceNet, 308
GetImpedanceTestFrequency, 293	!CMultiwellDeviceNet, 309
GetReady, 293	∼CMultiwellDeviceNet, 309
GetResult, 293	ClosePlateClamp, 310
SetImpedanceTestFrequency, 293	CMultiwellDeviceNet, 309
StartMeasurement, 293	GetPlateClampLockState, 310
CMeasureTableDeviceNet, 293	GetPlateClampState, 310
CMeasureTableDeviceNet, 294	GetPlateMux, 310, 311
Sensor, 294	GetPlateType, 311
CMeaSwitchDeviceNet, 294	GetPowerMuxPlate, 311
~CMeaSwitchDeviceNet, 295	IsPlateTypeValid, 312
CMeaSwitchDeviceNet, 295	LockPlateClamp, 312
GetNumber, 295	OpenPlateClamp, 312
GetPattern, 295	SetPlateMux, 312, 313
GetPatternBool, 296	SetPlateType, 313
SetPattern, 296	SetPowerMuxPlate, 313
SetPatternBool, 296	StopPlateClamp, 314
CMeaUSBDeviceNet, 296	UnlockPlateClamp, 314
CMeaUSBDeviceNet, 296 ∼CMeaUSBDeviceNet, 297	UnlockPlateClamp, 314 CMultiwellOptoStimFunctionNet, 314
	•
~CMeaUSBDeviceNet, 297	CMultiwellOptoStimFunctionNet, 314
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 299 !CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 CMeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 CMeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDurationHighCurrentInMicroSec, 318
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 ~CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 299 !CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300 ChannelReset, 301 CMultiBatteryChargerDeviceNet, 300 GetBatteryVoltage, 301 GetChannelState, 301 GetChannelState, 301	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 ~CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 299 !CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300 ChannelReset, 301 CMultiBatteryChargerDeviceNet, 300 GetBatteryVoltage, 301 GetChannels, 301 GetChannelState, 301 GetChargeCapacity, 302	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319
~CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDuryCycleHighCurrent, 318 SetMaxDuryCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ○MeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDuryCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDurationHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 CNF_GenDeviceNet, 319
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 CMeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 Set_Values, 319
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 Set_Values, 319 COctoPotDeviceNet, 320
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 ~CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300 ChannelReset, 301 CMultiBatteryChargerDeviceNet, 300 GetBatteryVoltage, 301 GetChannelState, 301 GetChargeCapacity, 302 GetChargingMode, 302 GetChargingPCoefficient, 302 GetDischargeCurrent, 303 GetDischargeCurrent, 303 GetDischargeCurrent, 303	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 CNF_GenDeviceNet, 319 Set_Values, 319 COctoPotDeviceNet, 320 BurnAdcOffset, 321
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 ~CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 299 !CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300 ChannelReset, 301 CMultiBatteryChargerDeviceNet, 300 GetBatteryVoltage, 301 GetChannels, 301 GetChannelState, 301 GetChargeCapacity, 302 GetChargingMode, 302 GetChargingPCoefficient, 302 GetDischargeCurrent, 303 GetDischargeCurrent, 303 GetDischargeCurrentSetPoint, 303 GetDischargeCurrentSetPoint, 303	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 CotPermanentCurrentInMicroAmp, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 CNF_GenDeviceNet, 319 Set_Values, 319 COctoPotDeviceNet, 320 BurnAdcOffset, 321 BurnDacOffset, 321
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 ∼CMeFunctionNet, 298 SetTransformer, 298 CMosMea	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetWaveLengthInNanometer, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 CNF_GenDeviceNet, 319 Set_Values, 319 COctoPotDeviceNet, 320 BurnAdcOffset, 321 BurnDacOffset, 321 COctoPotDeviceNet, 320, 321
∼CMeaUSBDeviceNet, 297 CMeaUSBDeviceNet, 297 CMeFunctionNet, 297 !CMeFunctionNet, 298 ~CMeFunctionNet, 298 ~CMeFunctionNet, 298 SetTransformer, 298 CMosMea CCMOSMeaDeviceNet, 47 CMultiBatteryChargerDeviceNet, 299 !CMultiBatteryChargerDeviceNet, 300 ~CMultiBatteryChargerDeviceNet, 300 CapacityTest, 300 ChannelReset, 301 CMultiBatteryChargerDeviceNet, 300 GetBatteryVoltage, 301 GetChannels, 301 GetChannelState, 301 GetChargeCapacity, 302 GetChargingMode, 302 GetChargingPCoefficient, 302 GetDischargeCurrent, 303 GetDischargeCurrent, 303 GetDischargeCurrentSetPoint, 303 GetDischargeCurrentSetPoint, 303	CMultiwellOptoStimFunctionNet, 314 !CMultiwellOptoStimFunctionNet, 315 ~CMultiwellOptoStimFunctionNet, 315 CMultiwellOptoStimFunctionNet, 315 GetAbsMaxCurrentInMicroAmp, 315 GetColorRgb, 315 GetColorStr, 316 GetMaxDurationHighCurrentInMicroSec, 316 GetMaxDutyCycleHighCurrent, 316 GetPermanentCurrentInMicroAmp, 317 GetWaveLengthInNanometer, 317 SetAbsMaxCurrentInMicroAmp, 317 SetColorRgb, 317 SetColorRgb, 317 SetColorStr, 318 SetMaxDurationHighCurrentInMicroSec, 318 SetMaxDutyCycleHighCurrent, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 SetPermanentCurrentInMicroAmp, 318 CotPermanentCurrentInMicroAmp, 319 CNF_GenDeviceNet, 319 ~CNF_GenDeviceNet, 319 CNF_GenDeviceNet, 319 Set_Values, 319 COctoPotDeviceNet, 320 BurnAdcOffset, 321 BurnDacOffset, 321

EnableTimestamp, 321	CPositionImpDeviceNet, 344
GetAdcOffset, 321	ConnectImp
GetDacOffset, 321	CPositionImpDeviceNet, 344
PatternListStart, 321	ConnectSlave
RampStart, 322	CGilsonDeviceNet, 99
ResetAdcOffset, 322	ControlState
ResetDacOffset, 322	HeadStageIDTypeState, 611
SetAdcOffset, 322	CornerFrequency
SetAmplificationSwitch, 322	CFilterPropertyNet, 73
SetBathclamp, 322	CornerFrequencymHz
SetChannelSwitch, 322	CFilterPropertyNet, 73
SetDacAutoControl, 322	Count
SetDacOffset, 322	CMcsUsbListNet, 241
SetDacValue, 323	CPatchServerDeviceNet, 327
SetNumberOfChannels, 323	CPatchServerDeviceNet, 328
SetOutputRate, 323	Sensor, 328
SetPatternListEntry, 323	CPathIdentDeviceNet, 328
SetPidParameter, 323	
	~CPathIdentDeviceNet, 329
SetRampParameter, 323	CPathIdentDeviceNet, 329
SetSineParameter, 323	Measure, 329
SineStart, 324	Set_Values, 329
COkuvisionStimulatorDeviceNet, 324	CPedoterDeviceNet, 329
~COkuvisionStimulatorDeviceNet, 325	!CPedoterDeviceNet, 330
COkuvisionStimulatorDeviceNet, 325	~CPedoterDeviceNet, 330
GetCheckVoltage, 325	CPedoterDeviceNet, 330
GetCurrentFactor, 325	GetCommand, 330
GetDACOffset, 325	SetCommand, 330
GetMaxPower, 325	CPeristalticPumpDeviceNet, 331
GetMaxVoltage, 325	\sim CPeristalticPumpDeviceNet, 331
GetPulseform, 325	CPeristalticPumpDeviceNet, 331
GetRTC, 326	McsBus_MotorControl, 332
GetStimulatorStatus, 326	CPgaDeviceNet, 332
GetVoltage, 326	\sim CPgaDeviceNet, 332
SetCheckVoltage, 326	ApplyGains, 333
SetCurrentFactor, 326	CPgaDeviceNet, 332
SetDACOffset, 326	DefineAmplification, 333
SetMaxPower, 327	DefineFrequencyRange, 333
SetMaxVoltage, 327	DefineNumAmplifications, 333
SetPulseform, 327	DefineNumFrequencyRanges, 333
SetRTC, 327	GetAmplification, 333
Coldstart	GetFrequencyRange, 333
CMcsUsbFactoryNet, 225	GetGain, 334
CommaPositionA	GetNumAmplifications, 334
CFilterCoefficientsNet::s_FilterAttributesNet, 622	GetNumFrequencyRanges, 334
CommaPositionB	SetGain, 334
CFilterCoefficientsNet::s_FilterAttributesNet, 622	CPositionIIDeviceNet, 334
CompareFirmware	!CPositionIIDeviceNet, 336
CMcsUsbFactoryNet, 226	\sim CPositionIIDeviceNet, 336
CompareTo	CPositionIIDeviceNet, 336
HeadStageIDType, 607	GetCoilCommunication, 336
CompensateElectrodeOffset	GetDebugData, 336
CWarnerUssingFunctionNet, 546	GetEventData, 337
Connect	GetImplantCurrentSetpoint, 337
CMcsUsbNet, 246, 247	GetImplantResult, 338
CRFFunctionNet, 372	GetImplantState, 338
ConnectDevice	GetOnOff, 338
CRadioControledDevicesNet, 367	GetPowerStrength, 339
ConnectedImp	GetRTC, 339
•	

GetStateDebugData, 339	GetAnalogVoltages, 358
GetStateEventData, 340	GetBubbleState, 358
RFFunction, 343	GetDigitalIn, 358
SetImplantCurrentSetpoint, 340	GetPumpCouple, 358
SetPowerStrength, 340	GetPumpEnableSpeedRatio, 358
SetRTC, 340	GetPumpFastOnOff, 358
SetStateDebugData, 342	GetPumpFastSpeed, 358
SetStateEventData, 342	GetPumpFunctionSpeeds, 359
SwitchOnOff, 342	GetPumpManualOnOff, 359
CPositionImpDeviceNet, 343	· · · · · · · · · · · · · · · · · · ·
•	GetPumpMaxSpeed, 359
!CPositionImpDeviceNet, 344	GetPumpModeType, 359
~CPositionImpDeviceNet, 344	GetPumpSpeedRatio, 359
ConnectedImp, 344	GetPumpSpeedUnit, 359
ConnectImp, 344	GetSupplyVoltage, 359
CPositionImpDeviceNet, 344	GetUseBubble, 359
GetDeviceList, 344	SetAnalogVoltages, 359
GetImpId, 345	SetPumpCouple, 360
GetRFFrequency, 345	SetPumpEnableSpeedRatio, 360
SetDeviceList, 345	SetPumpFastOnOff, 360
SetImpld, 345	SetPumpFastSpeed, 360
SetRFFrequency, 346	SetPumpFunctionSpeeds, 360
CPPCDeviceNet, 346	SetPumpManualOnOff, 360
CPPCDeviceNet, 346	SetPumpMaxSpeed, 360
McsBus, 347	SetPumpModeType, 360
McsBus_MotorControl, 347	SetPumpSpeedRatio, 361
McsBus_Sensor, 347	SetPumpSpeedUnit, 361
PPCFunction, 347	SetUseBubble, 361
CPPCFunctionNet, 347	CPPSDeviceNet, 361
!CPPCFunctionNet, 349	CPPSDeviceNet, 362
~CPPCFunctionNet, 349	CProgramPressureCurveNet, 362
CPPCFunctionNet, 348	!CProgramPressureCurveNet, 363
FirePressurePulse, 349	~CProgramPressureCurveNet, 362
GetAnalogVoltage, 349	CProgramPressureCurveNet, 362
GetAnalogVoltageRange, 349	GetRepeats, 363
GetDigitalIn, 351	Program, 363
GetPressureRange, 351	SetRepeats, 363
GetPumpModeType, 351	CPulseGeneratorFunctionNet, 363
GetPumpSpeedUnit, 352	!CPulseGeneratorFunctionNet, 364
GetSupplyVoltage, 352	\sim CPulseGeneratorFunctionNet, 364
GetValveActive, 352	CPulseGeneratorFunctionNet, 364
IsBusy, 352	GetModeSelect, 365
LoadPressure, 354	GetPeriod, 365
MeasureReservoir, 354	GetPulseLength, 365
SetAnalogVoltageRange, 354	SetModeSelect, 366
SetPressureOffset, 354	SetPeriod, 366
SetPressureRange, 354	SetPulseLength, 366
_	CRadioControledDevicesNet, 366
SetPumpModeType, 355	
SetPumpSpeedUnit, 355	ConnectDevice, 367
SetValveActive, 355	CRadioControledDevicesNet, 367
CPPS_DeviceNet, 356	DisConnectDevice, 367
CPPS_DeviceNet, 356	GetDeviceNames, 367
McsBus, 356	GetFrequency, 368
McsBus_MotorControl, 356	HasRadioControl, 368
McsBus_Sensor, 356	SetFrequency, 368
PPS_Function, 356	StillConnected, 368
CPPS_FunctionNet, 357	CreateSideband
CPPS_FunctionNet, 357, 358	CStimulusFunctionNet, 487
GetAnalogVoltage, 358	CreateWirelessHeadstageSerialNumberString

OM/mala as Dana Franchis as Nat. 500	OIT- -
CWirelessBaseFunctionNet, 593	CancelTableLoopAndStopTable, 382
CRegionOfInterestRect	ClampAmpRestart, 382
CCMOSMeaDeviceNet::CRegionOfInterestRect,	CRoboDacqNet, 382
368	DoRamp, 382
CRetinaLedDeviceNet, 369	Emu_GetCellCapacity, 382
\sim CRetinaLedDeviceNet, 370	Emu_GetCellPotential, 382
AddLoopEntry, 370	Emu_GetCellResists, 382
AddTableEntry, 370	Emu_GetElectrodeResists, 383
ClearTable, 370	Emu_GetNoise, 383
CRetinaLedDeviceNet, 370	Emu_SetCellCapacity, 383
GetTablepointer, 370	Emu_SetCellPotential, 383
SetLED, 370	Emu_SetCellResists, 383
SetLumi, 370	Emu_SetElectrodeResists, 383
SetPersistency, 370	Emu_SetNoise, 383
SetRepeat, 371	GetAllDigout, 383
SetTablepointer, 371	GetCapacityC, 383
SetTrigger, 371	GetCapacityV, 383
CRFFunctionNet, 371	GetCapacityX, 384
!CRFFunctionNet, 372	GetClampAmpSerialNumber, 384
~CRFFunctionNet, 372	GetCommand, 384
Connect, 372	GetConfigurationBit, 384
CRFFunctionNet, 372	GetConfigurationBit, 364 GetConfigurationBitAxc, 384
GetAvailableDeviceList, 373	GetConfigurationBitBlu_Led, 384
•	GetConfigurationBitBlu LedToggleFast, 384
GetAvailableDeviceListEx, 373	
GetAvailableStateList, 373	GetConfigurationBitBlu_LedToggleSlow, 384
GetAvailableStateListEx, 373	GetConfigurationBitCC_Gen, 384
GetBaseFrequency, 374	GetConfigurationBitCV_Gen, 384
GetConnectedDevice, 374	GetConfigurationBitRC_Gen, 384
GetState, 374	GetConfigurationBitRed_Led, 385
GetTestMode, 374	GetConfigurationBitRed_LedSaturation, 385
GetWorkingFrequency, 375	GetConfigurationBitRed_LedToggleFast, 385
SetBaseFrequency, 375	GetConfigurationBitRed_LedToggleSlow, 385
SetTestMode, 375	GetConfigurationBitRelais, 385
SetWorkingFrequency, 375	GetConfigurationBitRV_Gen, 385
CRobo_FYIProgram_FunctionNet, 376	GetConfigurationBits, 385
CRobo_FYIProgram_FunctionNet, 376	GetConfigurationBitStream, 385
GetLength, 376	GetConfigurationBitSupply, 385
GetState, 377	GetCrossTalkOffset, 385
GetValve1, 377	GetCrossTalkOptimum, 385
GetValve2, 377	GetDigout, 386
SetLength, 377	GetDisplayText, 386
SetValve1, 377	GetDownsampleFactor, 386
SetValve2, 377	GetFilter, 386
Start, 377	GetFilterCoeffs, 386
CRobo_FYITemp_FunctionNet, 377	GetIC, 386
CRobo_FYITemp_FunctionNet, 378	GetlClamp, 386
GetlCoeff, 378	GetICOffset, 386
GetMaxPower, 378	GetlGain, 386
GetPCoeff, 378	GetNIC MS, 386
GetRegulatorOnOff, 378	GetNUC_MS, 386
GetSollTemp, 378	GetNUV_MS, 387
SetICoeff, 379	GetPGain, 387
SetMaxPower, 379	GetRecordingNumber, 387
SetPCoeff, 379	GetResistanceC, 387
SetRegulatorOnOff, 379	GetResistanceV, 387
SetSollTemp, 379	GetScreen, 387
CRoboDacqNet, 379	GetSimulation, 387
CancelTableLoop, 382	GetUC, 387
oancenableLoop, soz	UGIUU, 301

GetUClamp, 387	Axis_Z, 403
GetUCOffset, 387	CancelPoolLoop, 396
GetUpdateDisplay, 387	CancelPoolLoopAndStopMovement, 396
GetUV, 388	CRoboDeviceNet, 396
GetUVOffset, 388	FindReference, 396
GetXGain, 388	GetAirpressure, 397
RunTable, 388	GetAirpressureLimit, 397
SetAllDigout, 388	GetAirValve, 397
SetCommand, 388	GetCurrentAirvalve, 397
SetConfigurationBit, 388	GetCurrentAirvalveLimit, 397
SetConfigurationBitAxc, 388	GetCurrentPosition, 397
SetConfigurationBitBlu_Led, 388	GetErrorAirpressure, 398
SetConfigurationBitBlu_LedToggleFast, 389	GetErrorCurrentAirvalve, 398
SetConfigurationBitBlu_LedToggleSlow, 389	GetErrorVoltage12V, 398
SetConfigurationBitCC_Gen, 389	GetErrorVoltage5V, 398
SetConfigurationBitCV_Gen, 389	GetErrorVoltageAirvalve, 398
SetConfigurationBitRC_Gen, 389	GetErrorVoltageRs485A, 398
SetConfigurationBitRed_Led, 389	GetErrorVoltageRs485B, 398
SetConfigurationBitRed_LedSaturation, 389	GetErrorVoltageValves, 398
SetConfigurationBitRed_LedToggleFast, 389	GetInMovement, 398
SetConfigurationBitRed_LedToggleSlow, 389	GetMinPressure, 398
SetConfigurationBitRelais, 389	GetMovementError, 398
SetConfigurationBitRV_Gen, 390	GetVoltage12V, 399
SetConfigurationBitStream, 390	GetVoltage12VLimit, 399
SetConfigurationBitSupply, 390	GetVoltage5V, 399
SetCrossTalkOffset, 390	GetVoltage5VLimit, 399
SetCrossTalkOptimum, 390	GetVoltageAirvalve, 399
SetDigout, 390	GetVoltageAirvalveLimit, 399
SetDisplayText, 390	GetVoltageRs485A, 399
SetDownsampleFactor, 390	GetVoltageRs485ALimit, 399
SetFilter, 390	GetVoltageRs485B, 399
SetFilterCoeffs, 391	GetVoltageRs485BLimit, 400
SetlClamp, 391	GetVoltageValves, 400
SetICOffset, 391	GetVoltageValvesLimit, 400
SetlGain, 391	McsBus, 407
SetNoFilterCoeffs, 391	McsBus_MotorControl, 407
SetPGain, 391	McsBus_XY, 404
SetRecordingNumber, 391	McsBus_ZI, 404
SetScreen, 391	MoveAbs, 400
SetSimulation, 391	RoboError AnotherMaster, 404
SetUClamp, 392	RoboError_Base, 404
SetUCOffset, 392	RoboError_CannotEscapeEndSwitch, 404
SetUVOffset, 392	RoboError_CommandAlreadyInProgress, 404
SetXGain, 392	RoboError CommandNotPossible, 404
StopTable, 392	RoboError_CommunicationTimeout, 404
Table Wait, 392	RoboError DacqNotReady, 404
TableDefBegin, 392	RoboError DLLMovementTimeout, 405
TableDefEnd, 392	RoboError FindReferenceMethod, 405
UpdateDisplay, 392	RoboError GilsonCommandPending, 405
CRoboDeviceNet, 393	RoboError_GilsonTimeout, 405
∼CRoboDeviceNet, 396	RoboError_GilsonWrondID, 405
Axes_I, 403	RoboError_McsBus_UnknownCommand, 405
Axes_X, 403	RoboError NoEndSwitch, 405
Axes Y, 403	RoboError NoMoreData, 405
Axes_Z, 403	RoboError NoReference, 405
Axis_I, 403	RoboError_NoSpeedOrAcceleration, 405
Axis X, 403	RoboError_OverPressure, 406
Axis_Y, 403	RoboError_ParameterNotAllowed, 406
_	

RoboError_PeristalticTimeout, 406 RoboError_Phase0OutOfRange, 406 RoboError_PollLoopCanceled, 406 RoboError_PollLoopCanceledAndStopMovement,	SetSearchReferenceFastSpeed, 619 SetSearchReferenceFineAccel, 620 SetSearchReferenceFineSpeed, 620 SetSearchReferenceMethod, 620
406	SetSearchReferenceMoveOut, 620
RoboError_Pressure, 406	SetSearchReferenceOffsetPos, 620
RoboError_RangeExceeded, 406	SetUserParameter, 620, 621
RoboError_StateChangeNotPossible, 406	CRoboFluidDeviceNet, 407
RoboError_Timeout, 406	\sim CRoboFluidDeviceNet, 408
RoboError_UnknownCommand, 407	CloseAllValves, 408
RoboMainLowLevelCommand, 407	CRoboFluidDeviceNet, 408
RoboStatusEvent, 407	GetPumpSpeed, 409
SetAirpressureLimit, 400	GetSingleValve, 409
SetAirValve, 401	GetValve, 409
SetCurrentAirvalveLimit, 401	IsPumpMotorOn, 409
SetCurrentAndAir, 401	m_pMcsBus_MotorControlNet, 410
SetInMovement, 401	m_pRoboFluidDevice, 410
SetMinPressure, 401	McsBus_MotorControl, 410
SetVoltage12VLimit, 401	PumpOff, 409
SetVoltage5VLimit, 401	PumpOn, 409
SetVoltageAirvalveLimit, 402	SetPumpSpeed, 409
SetVoltageRs485ALimit, 402	SetSingleValve, 410
SetVoltageRs485BLimit, 402	SetValve, 410
SetVoltageValvesLimit, 402	CRobolnjectDeviceNet, 411
StopMovement, 402	CRobolnjectDeviceNet, 411
CRoboDeviceNet::RoboMainLowLevelCommands, 614	CRoboocyte2DeviceNet, 411
FindReferencePhase0, 615	CRoboocyte2DeviceNet, 412
GetAxisConfig, 615	GetAxisLED, 412
GetHWConfig, 615	GetGilsonDevice, 412
GetHWRevision, 616	GetMcsBus_Extension, 412
GetMaxNoPressure, 616	GetRoboDacq, 412
GetMaxNoPressureWaitTime, 616	GetRoboFluidDevice, 412
GetMaxPressureWaitTime, 616	SetAxisLED, 413
GetMinNoPressureWaitTime, 616	CRoboStatorDeviceNet, 413
GetMinPressure, 616	CRoboStatorDeviceNet, 414
GetMinPressureWaitTime, 616	FindReferencel, 414
GetParameter, 616	FindReferenceXY, 414, 415
GetPhases, 616	FindReferenceZ, 415
GetSearchReferenceFastAccel, 617	GetCurrentPositionI, 415
GetSearchReferenceFastSpeed, 617	GetCurrentPositionXY, 415
GetSearchReferenceFineAccel, 617	GetCurrentPositionZ, 415
GetSearchReferenceFineSpeed, 617	HasRefl, 415
GetSearchReferenceMethod, 617	HasRefXY, 415
GetSearchReferenceMoveOut, 617	HasRefZ, 415
GetSearchReferenceOffsetPos, 617	MoveAbsl, 415, 416
GetUserParameter, 617, 618	MoveAbsXY, 416
HasRef, 618	MoveAbsZ, 416
SetAxisConfig, 618	RoboMainStatorLowLevelCommand, 419
SetHWConfig, 618	SetAccelerationI, 416
SetHWRevision, 618	SetAccelerationNativel, 416
SetMaxNoPressure, 618	SetAccelerationNativeXY, 416
SetMaxNoPressureWaitTime, 618	SetAccelerationNativeZ, 416
SetMaxPressureWaitTime, 619	SetAccelerationXY, 417
SetMinNoPressureWaitTime, 619	SetAccelerationZ, 417
SetMinPressure, 619	SetCurrentAndAirXY, 417
SetMinPressureWaitTime, 619	SetSpeedl, 417
SetParameter, 619	SetSpeedNativel, 417
SetSearchReferenceFastAccel, 619	SetSpeedNativeXY, 417

SetSpeedNativeZ, 417	GetReferenceElectrodeMode, 432
SetSpeedXY, 417	GetReferenceElectrodeSwitchState, 433
SetSpeedZ, 418	HasAnalogOut, 433
SetVelocityI, 418	HasGalvanicIsolation, 433
SetVelocityXY, 418	HasHSPowerSwitch, 433
SetVelocityZ, 418	IsAnalogOutEnabled, 434
StopMovementI, 418	IsAutomaticAnalogOut, 434
StopMovementXY, 418	IsHeadstageAvailable, 434
StopMovementZ, 418	IsHeadstageAvailableEvent, 437
CRoboStatorDeviceNet::RoboMainStatorLowLevelComm	
621	IsInDacqLegacyMode, 435
FindReferencePhase0XY, 621	OnGetAvailableHeadstages, 435
CSafeISDeviceNet, 419	OnIsHeadstageAvailable, 435
~CSafeISDeviceNet, 420	PowerHS, 435
CSafeISDeviceNet, 420	SetAnalogOutADCRange, 435
DacqDevice, 421	SetAnalogOutChannels, 436
FluidControlDevice, 422	SetAnalogOutDACRange, 436
RoboDevice, 422	SetDacqLegacyMode, 436
SetAdcChannels, 420	SetHeadstagePowerStateAtStart, 436
SetAdcSamplePos, 420	SetReferenceElectrodeMode, 437
SetDacMode, 420	SetReferenceElectrodeSwitchState, 437
SetDacPeriode, 421	CSerialPortNet, 438
SetDacPulseform, 421	CSerialPortNet, 438
SetSwitches, 421	GetBytesAvailable, 438
CSCUDacqGroupChannelSelectionNet, 422	Receive, 438
CSCUDacqGroupChannelSelectionNet, 422	ReceiveString, 438, 439
CSCUFunctionNet, 423	Send, 439
!CSCUFunctionNet, 425	CStg200xBasicNet, 439
~CSCUFunctionNet, 425	~CStg200xBasicNet, 443
AutomaticAnalogOut, 425	GetAnalogRanges, 443
CSCUFunctionNet, 425	GetAnalogResolution, 444
EnableAnalogOut, 426	GetAutocalibrationDisabled, 444
GetAnalogOutADCRange, 426	GetAvailableMemory, 444
GetAnalogOutChannels, 426	GetBlankingEnable, 444, 445
GetAnalogOutDACRange, 426	GetCurrentRangeInNanoAmp, 445
GetAvailableHeadstages, 426	GetCurrentResolutionInNanoAmp, 445
GetAvailableHeadstagesEvent, 437	GetDacAmplificationFactor, 446
GetFilterProperties, 427	GetDACResolution, 446
GetFilterProperty, 427	GetDACHesolution, 440 GetDiginValue, 446
GetHeadstageAdcBits, 427	GetDigitValue, 446 GetDigoutMode, 446
GetHeadstageAdcBits, 427 GetHeadstageAdcRangeInMicroVolt, 428	GetDigoutWalue, 447
GetHeadstageDacBits, 428	GetElectrodeDacMux, 447
GetHeadstageDacCurrentRangeInMicroAmpere,	GetElectrodeEnable, 448
428	GetElectrodeEnable, 449
GetHeadstageDacCurrentResolutionInNanoAm-	GetEnableAmplifierProtectionSwitch, 449, 450
pere, 429	GetExternalElectrodeEnable, 450
GetHeadstageDacVoltageRangeInMilliVolt, 429	GetFAAmplification, 451
GetHeadstageDacVoltageResolutionInMicroVolt,	GetHeadstage, 451
429	GetListmodeIndexRange, 451
GetHeadstageGainInPermille, 430	GetListmodeTriggerSource, 451
GetHeadstageID, 430	GetNumberOfINADACRethe 451
GetHeadstageNumberOfAnalogChannels, 430	GetNumberOfHWDACPaths, 451
GetHeadstageNumberOfStimulationChannels, 431	GetNumberOfStimulationElectrodes, 451
GetHeadstagePowerStateAtStart, 431	GetNumberOfStimulationSourcesPerElectrode,
GetHeadstageSamplerate, 431	452
GetHeadstageSerialNumber, 432	GetNumberOfSyncoutChannels, 452
GetMaxNumberOfHeadstages, 432	GetNumberOfTriggerInputs, 452
GetMaxStimulusChannelsPerHeadstage, 432	GetOutputRate, 452

	GetStgProgramInfo, 452, 453	PrepareAndAppendData, 481
	GetStgVersionInfo, 453	PrepareAndSendData, 482
	GetSyncoutMap, 453	QueryTriggerstatus, 483
	GetTotalMemory, 454	SendSegmentDefine, 483
	GetTriggerSource, 454	SendSegmentSelect, 483
	GetVoltageRangeInMicroVolt, 454	SendSegmentStart, 484
	GetVoltageResolutionInMicroVolt, 454	SetOutputMap, 484
	ListModeSendStart, 455	Stg200xPollStatusEvent, 484
	ListModeSendStop, 455	CStimulusFunctionNet, 485
	SendStart, 455	ClearChannel_PrepareAndSendData, 487
	•	ClearChannelData, 487
	SendStop, 455	
	SetAutocalibrationDisabled, 455	ClearMultiplexedData, 487
	SetBlankingEnable, 456, 457	ClearSyncData, 487
	SetCurrentMode, 457	CreateSideband, 487
	SetDacAmplificationFactor, 457	CStimulusFunctionNet, 486
	SetDigoutMode, 458	ForceStatusEvent, 488
	SetDigoutValue, 458	GetAvailableMemory, 488
	SetElectrodeDacMux, 458–460	GetCurrentRangeInNanoAmp, 489
	SetElectrodeEnable, 460–462	GetCurrentResolutionInNanoAmp, 489
	SetElectrodeMode, 463, 464	GetDACResolution, 489
	SetEnableAmplifierProtectionSwitch, 464, 465	GetMultiplexedDataChannelsInBlock, 489
	SetExternalElectrodeEnable, 466	GetNumberOfAnalogChannels, 490
	SetFAAmplification, 467	GetTotalMemory, 490
	SetHeadstage, 467	GetVoltageRangeInMicroVolt, 490
	SetListmodeIndexRange, 467	GetVoltageResolutionInMicroVolt, 490
	SetListmodeTriggerSource, 467	PollStatusEvent, 496
	SetMeasurementMode, 467	PrepareAndAppendData, 491
	SetOutputRate, 468	PrepareAndSendData, 492
	SetStgProgramInfo, 468	PrepareData, 492
	SetSyncoutMap, 468	SendMultiplexedData, 493
	SetTriggerSource, 468, 469	SendPreparedData, 493
	SetVoltageMode, 469	SendStart, 493
CS	tg200xDownloadBasicNet, 469	SendStop, 493
	ClearChannelData, 471	SetupTrigger, 494
	ClearSyncData, 471	SetupTriggerSingle, 494
	DisableAutoReset, 471	StartPoll, 496
	EnableAutoReset, 472	StopPoll, 496
	ForceStatusEvent, 472	CStimulusFunctionNet::SidebandData, 623
	GetMemoryUsageDAC, 472	!SidebandData, 624
	GetMemoryUsageSyncout, 472	∼SidebandData, 624
	GetSweepCount, 473	Duration, 624
	GetTrigger, 473	Sideband, 624
	ResetStatus, 473	SidebandData, 623
	SendChannelData, 475	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData
	SendSyncData, 475	625
	SetupRetriggerMode, 476	!StimulusDeviceDataAndUnrolledData, 625
	SetupTrigger, 476	~StimulusDeviceDataAndUnrolledData, 625
	SetupTriggerSingle, 477	DeviceData, 626
	Stimulus, 478	DeviceDataLength, 626
CS	tg200xDownloadNet, 478	StimulusDeviceDataAndUnrolledData, 625
	~CStg200xDownloadNet, 479	UnrolledAmplitude, 626
	ClearChannel_PrepareAndSendData, 480	UnrolledDuration, 626
	CStg200xDownloadNet, 479	UnrolledSync, 626
	DisableMultiFileMode, 480	CSw2to64DeviceNet, 496
	EnableMultiFileMode, 480	~CSw2to64DeviceNet, 497
	GetModuleCurrent, 481	CSw2to64DeviceNet, 497
	GetModuleTemp, 481	GetChannel, 497
	MwPollStatusEvent. 484	GetChannels, 498
	IVIVI OHOLULUOLIVOIILI TUT	MOLOTIGITION IVV

GetNumber, 498	GetThermocoupleCalibration, 508
SetChannel, 498	GetThermocoupleNanovoltPerKelvin, 508
SetChannels, 498	GetThermocoupleReferenceTemp, 508
CTcxDeviceNet, 499	GetThermocoupleTemp, 508
∼CTcxDeviceNet, 501	GetThermocoupleTempAbs, 509
CalibrateThermocouple, 501	GetUnit, 509
CTcxDeviceNet, 501	GetUOut, 509
FactoryReset, 501	GetValue, 509
GetBoardTemp, 502	GetValueHires, 509
GetCalibration, 502	GetVolti, 509
GetCalibrationDecp, 502	SetCalibration, 509
GetCalibrationMax, 502	SetD, 510
GetCalibrationMin, 502	SetDevice, 510
GetCurrent, 502	SetDeviceType, 510
GetD, 502	SetDevname, 510
GetDDecp, 502	SetEnableHeaterLimit, 510
GetDevice, 502	SetEnableThermocouple, 510
GetDeviceType, 503	SetHeaterLimit, 510
GetDevname, 503	Setl, 510
GetDMax, 503	SetMaxHeaterPowerMultiwell, 511
	•
GetDMin, 503 GetDuty, 503	SetMaxP, 511 SetOnOff, 511
• .	
GetEnableHeaterLimit, 503	SetP, 511
Get Las Thormosouple, 503	SetSetneint 511
GetHasThermocouple, 503	SetSetpoint, 511
GetHeaterLimit, 504	SetThermocoupleNanovoltPerKelvin, 512
GetHeaterTemp, 504	CTEERFunctionNet, 512
Getl, 504	!CTEERFunctionNet, 514
GetIDecp, 504	~CTEERFunctionNet, 514
GetIMax, 504	CancelInternalCalibration, 514
GetIMin, 504	CTEERFunctionNet, 514
GetlOut, 504	GetAdapterCode, 515
GetMaxHeaterPowerMultiwell, 505	GetAdcOffsetU1, 515
GetMaxP, 505	GetAdcOffsetU2, 515
GetMaxpDecp, 505	GetAmplitude_nA, 515
GetMaxpMax, 505	GetBytesPerSample, 515
GetMaxpMin, 505	GetClampMode, 516
GetNumControlChannels, 505	GetControllerParams, 516
GetNumDevices, 505	GetCurrentEnable, 516
GetNumMeasureChannels, 505	GetDacZero, 516
GetOnOff, 506	GetFrameErrorCounter, 516
GetP, 506	GetLiquidResistance, 517
GetPDecp, 506	GetMaxChunkSize_Byte, 517
GetPMax, 506	GetNumberOfAvailableSamples, 517
GetPMin, 506	GetPeriod_us, 517
GetPOut, 506	GetRotaryPositionCode, 517
GetPwrOut, 506	GetSampleBufferChunk, 517
GetPwrSet, 507	GetSampleRate, 518
GetRes1, 507	GetSampleVoltageBuffer uV, 518
GetRes2, 507	GetScaleFactorU1, 518
GetResS, 507	GetScaleFactorU2, 518
GetResX, 507	GetUptimeSeconds, 519
GetROut, 507	GetWaveform, 519
GetSensorType, 507	IsInternalCalibrationFinished, 519
GetSetpoint, 507	IsSamplingFinished, 519
GetSetpointDecp, 508	SetAmplitude nA, 519
GetSetpointMax, 508	SetBufferIndex, 520
•	
GetSetpointMin, 508	SetClampMode, 520

SetControllerParams, 520	GetStimulusParametersFromSelectedHS, 531
SetCurrentEnable, 520	GetStiumlusParameters, 532
SetExternalLED, 521	GetUserDefinedName, 532
SetLiquidResistance, 521	GetUserDefinedNameCache, 532
SetPeriod_us, 521	GetUserDefinedNameFromSelectedHS, 532
SetWaveform, 521	PulseGenerator, 534
StartInternalCalibration, 522	SelectHeadstage, 532
StartSampling, 522	SetAccelGyroDesiredRate, 532
StopSampling, 522	SetAccelGyroEnabled, 532
CTEERMachineDeviceNet, 522	SetAccelRange, 532
\sim CTEERMachineDeviceNet, 523	SetAnalogOutChannel, 533
CTEERMachineDeviceNet, 523	SetAnalogOutFilter, 533
TEERFunctionNet, 523	SetAudioChannels, 533
CurrentRangeInNanoAmp	SetDacRange, 533
W2100_StimulusParametersNet, 627	SetGyroRange, 533
CurrentResolutionInNanoAmp	SetHeadstageOnOff, 533
W2100_StimulusParametersNet, 627	SetHeadstageSamplingActive, 533
CUsbDeviceConfigurationFunctionNet, 523	SetHeadstageToSleep, 533
!CUsbDeviceConfigurationFunctionNet, 524	SetMultiHeadstageMode, 534
~CUsbDeviceConfigurationFunctionNet, 524	SetSelectedChannels, 534
CUsbDeviceConfigurationFunctionNet, 524	Stimulator, 534
SetDeviceId, 524	CW2100_FunctionNet::AudioChannelsNet, 29
SetDeviceName, 524	amplification, 29
CUsbExceptionNet, 525	channel, 29
CUsbExceptionNet, 525, 526	dacqgroup, 29
Status, 526	CW2100_StimulatorFunctionNet, 534
CutoffFrequency	BOOST_BIT, 541
CCreateFilterNet, 49	ClearChannelData, 536
CVoltageRangeInfoNet	CW2100_StimulatorFunctionNet, 536
CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNe	
526	GetBoostPreTime, 536
CW2100_FunctionNet, 527	GetCurrentRangeInNanoAmp, 536
ClearStimulusParametersCache, 529	GetCurrentResolutionInNanoAmp, 537
ClearUserDefinedNameCache, 529	GetDACResolution, 537
CW2100_FunctionNet, 528	GetDigitalStimulatorTrigger, 537
DeselectAllHeadstages, 529	GetDigitalStimulatorTriggerSlope, 537
DeselectHeadstage, 529	GetNumberOfAnalogChannels, 537
GetAccelGyroCurrentRate, 529	GetNumberOfSyncoutChannels, 538
GetAccelGyroDesiredRate, 529	GetNumberOfTriggerInputs, 538
GetAccelGyroEnabled, 529	GetStimulationPatternMemory, 538
GetAccelRange, 529	GetTimeResolutionInNanoSeconds, 538
GetAnalogOutChannel, 530	GetTimeSlot, 538
GetAnalogOutFilter, 530	GetVoltageRangeInMicroVolt, 538
GetAudioChannels, 530	GetVoltageResolutionInMicroVolt, 538
GetAvailableHeadstages, 530	GND_SWITCH_BIT, 541
GetBatteryState, 530	PollStatusEvent, 542
GetDacRange, 530	PrepareData, 540
GetFilterProperties, 530	PrepareDataSync, 540
GetFilterProperty, 530	SelectTimeSlot, 540
GetFPGAFirmwareType, 530	SendPreparedData, 540
GetGyroRange, 530	SendStart, 540
GetHeadstageOnOff, 531	SendStop, 541
GetHeadstageSamplingActive, 531	SetDigitalStimulatorTrigger, 541
GetMultiHeadstageMode, 531	SetDigitalStimulatorTriggerSlope, 541
GetPicFirmwareType, 531	StartPoll, 541
GetSelectedChannels, 531	StopPoll, 541
GetSelectedHeadstageState, 531	SYNC_BIT0, 542
GetStimulusParametersCache, 531	SYNC BIT1, 542

CW2100DacqGroupChannelSelectionNet, 542	CWarnerValveControllerDeviceNet, 564
CW2100DacqGroupChannelSelectionNet, 542	GetActiveRunningTableNumber, 565
CWarnerUssingDeviceNet, 543	GetActiveRunningTableNumberEvent, 582
!CWarnerUssingDeviceNet, 543	GetAnalogThresholdHigh, 565
~CWarnerUssingDeviceNet, 543	GetAnalogThresholdHighEvent, 582
CWarnerUssingDeviceNet, 543	GetAnalogThresholdLow, 566
WarnerUssingFunction, 544	GetAnalogThresholdLowEvent, 582
CWarnerUssingFunctionNet, 544	GetAnalogVoltage, 566
!CWarnerUssingFunctionNet, 546	GetAnalogVoltageEvent, 582
~CWarnerUssingFunctionNet, 546	GetCurrentEditTableNumber, 566
CompensateElectrodeOffset, 546	GetCurrentNumberOfValves, 566
CWarnerUssingFunctionNet, 546	GetCurrentNumberOfValvesEvent, 582
GetAvailableChambers, 547	GetDigitalOutPortValve, 567
GetChannelsCountOfChamber, 547	GetDigitalOutPortValveEvent, 583
GetClampMode, 547	GetDigitalPortDirection, 568
GetComplianceVoltageThreshold, 547	GetDigitalPortDirectionEvent, 583
GetDacPampsPerDigitHighCurrentRange, 548	GetDisplayMode, 568
GetDacPampsPerDigitLowCurrentRange, 548	GetDisplayModeEvent, 583
GetDacZero, 548	GetTableName, 568
GetHighCurrentRange, 549	
	GetTableNamebyIndex, 568
GetIdleModeOffset, 549	GetTableNamebyIndexEvent, 583
GetLiquidResistance, 549	GetTotalNumberOfDigitalPorts, 569
GetLowCurrentRange, 550	GetTotalNumberOfTables, 569
GetNumberOfAvailableChambers, 550	GetTotalNumberOfValves, 569
GetNumberOfHardwareSlotsForChambers, 550	GetTotalTableSize, 569
GetU1Offset, 550	GetValveActive, 569
GetU1Reference, 552	GetValveActiveEvent, 583
GetU2Offset, 552	GetValveBoardRevision, 570
GetU2Reference, 552	GetValveBoardRevisionEvent, 583
GetUnitDescription, 553	GetValveBoardRevisionString, 570
GetUnitExponent, 553	GetValveDigitalInPort, 570
GetUnitName, 553	GetValveDigitalInPortEvent, 583
GetUnitsPerDigit, 554	GetValveLedOn, 570
GetUptimeSeconds, 554	GetValveLedOnEvent, 584
GetVoltageClampControllerParam_D, 554	GetValveManualGroup, 571
GetVoltageClampControllerParam_I, 555	GetValveManualGroupEvent, 584
GetVoltageClampControllerParam_P, 555	GetValveManualState, 571
IsChamberAvailable, 555	GetValveManualStateEvent, 584
IsHighCurrentMode, 556	GetValveMode, 571
IsInternalCalibrationFinished, 556	GetValveModeEvent, 584
IsPulseEnabled, 556	GetValvesActiveMap, 571
SetClampMode, 557	GetValvesManualStateMap, 572
SetEnablePulse, 557	GetValveTableEntry, 572
SetHighCurrentMode, 557	IsDigitalOutPortInverted, 572
SetIdleModeOffset, 558	IsDigitalOutPortInvertedEvent, 584
SetLiquidResistance, 558	IsValveDigitalInInverted, 572
SetLowCurrentMode, 558	IsValveDigitalInInvertedEvent, 584
SetPulse, 558	IsValveOpen, 573
SetVoltageClampControllerParam_D, 559	IsValveOpenEvent, 584
SetVoltageClampControllerParam_I, 559	IsValveOpenInAnalogMode, 573
SetVoltageClampControllerParam_P, 559	IsValveOpenInAnalogModeEvent, 584
WaitForAllChambers, 560	IsValveOpenInDigitalMode, 573
WaitForChamber, 560	IsValveOpenInDigitalModeEvent, 585
CWarnerValveControllerDeviceNet, 560	LoadValveTable, 574
!CWarnerValveControllerDeviceNet, 565	OnGetApplogThrosholdHigh 574
~CWarnerValveControllerDeviceNet, 565	OnGetAnalogThresholdHigh, 574
Clear Table Name, 565	OnGetAnalogThresholdLow, 574
ClearValveTable, 565	OnGetAnalogVoltage, 574

OnGetCurrentNumberOfValves, 574	CatEiltar Paramatara Haadataga 500
,	GetFilterParametersHeadstage, 589
OnGetDigitalOutPortValve, 574	GetHasChecksum, 589
OnGetDigitalPortDirection, 574	GetHasRedLedHeadstage, 589
OnGetDisplayMode, 575	GetHeadstageOnOff, 589
OnGetTableNamebyIndex, 575	GetResetFilter, 590
OnGetValveActive, 575	GetRFConnectionStatus, 590
OnGetValveBoardRevision, 575	GetRFFrequencyHeadstage, 590
OnGetValveDigitalInPort, 575	GetRFFrequencyReceiver, 590
OnGetValveLedOn, 575	GetRFPower, 590
OnGetValveManualGroup, 575	GetScanHeadstagesResult, 590
OnGetValveManualState, 575	GetSelectedHeadstage, 590
OnGetValveMode, 575	GetSerialNumberHeadstage, 590
OnIsDigitalOutPortInverted, 576	GetWPADebugMode, 590
OnIsValveDigitalInInverted, 576	GetWPAType, 590
OnIsValveOpen, 576	ResetChannelmap, 591
OnlsValveOpenInAnalogMode, 576	ScanForHeadstages, 591
OnlsValveOpenInDigitalMode, 576	SetChannelmap, 591
OnTableEntryChanged, 576	SetFilterParametersHeadstage, 591
SetActiveRunningTableNumber, 576	SetHasChecksum, 591
SetAnalogThresholdHigh, 577	SetHeadstageOnOff, 591
SetAnalogThresholdLow, 577	_
· ·	SetHWSelectedChannels, 591
SetCurrentEditTableNumber, 577	SetResetFilter, 591
SetDefault, 577	SetRFFrequencyHeadstage, 591
SetDigitalOutPortInvert, 577	SetRFFrequencyReceiver, 592
SetDigitalOutPortValve, 578	SetRFFrequencyReceiverEeprom, 592
SetDigitalPortDirection, 578	SetRFLostBehaviour, 592
SetDisplayMode, 578	SetRFPower, 592
SetTableName, 578	SetSelectedHeadstage, 592
SetTableStep, 579	SetSerialNumberHeadstage, 592
SetTableStepAll, 579	SetWPADebugMode, 592
SetValveActive, 579	SetWPAType, 592
SetValveDigitaIInInvert, 579	CWirelessBaseFunctionNet, 593
SetValveDigitaIInPort, 580	CreateWirelessHeadstageSerialNumberString,
SetValveLedOn, 580	593
SetValveManualGroup, 580	CWirelessBaseFunctionNet, 593
SetValveManualState, 580	CyclePort
SetValveMode, 581	CMcsUsbNet, 248
SetValvesActiveMap, 581	
SetValvesManualStateMap, 581	DacqDevice
SetValveTableEntry, 581	CSafeISDeviceNet, 421
Store Valve Table 1.111 y, 361	dacqgroup
TableEntryChangedEvent, 585	CW2100_FunctionNet::AudioChannelsNet, 29
	DACResolution
CWarnerValveControllerDeviceTesterFunctionNet, 585	W2100_StimulusParametersNet, 627
!CWarnerValveControllerDeviceTesterFunctionNet,	DataState
586	HeadStageIDTypeState, 611
$\sim\!\!CWarnerValveControllerDeviceTesterFunctionNet,$	DeepCopy
586	• • • •
CWarner Valve Controller Device Tester Function Net,	CCMOSMeaDeviceNet::CRegionOfInterestRect
586	368
GetIO, 586	DefineAmplification
GetSync, 586	CPgaDeviceNet, 333
SetADC, 587	DefineFrequencyRange
SetIO, 587	CPgaDeviceNet, 333
SetIODirection, 587	DefineNumAmplifications
SetTrigger, 587	CPgaDeviceNet, 333
SetTriggerSyncDirection, 588	DefineNumFrequencyRanges
CWClassicFunctionNet, 588	CPgaDeviceNet, 333
CWClassicFunctionNet, 589	DeselectAllHeadstages
2.1.0.000.0. 0.1.00000000	CW2100_FunctionNet, 529

DeselectHeadstage	GetMajor, 600
CW2100_FunctionNet, 529	GetMinor, 600
DetectChipType	GetNumEntries, 601
CCMOSMea_FunctionNet, 36	GetSerialNumber, 601
DEVICE_NOT_FOUND	GetStatus, 601
Mcs::Usb, 27	GetVersionInt, 602
DeviceArrival	GetVersionString, 602
CMcsUsbListNet, 241	DriverVersionNet::FormatVersion
DeviceData	DriverVersionNet, 599
CStimulusFunctionNet::StimulusDeviceDataAndUnro	ol ଢ6 ₽ata,
626	FirmwareDestinationNames, 604
DeviceDataLength	DummyCommand
CStimulusFunctionNet::StimulusDeviceDataAndUnro	olledD &a JH3DeviceNet, 123
626	Duration
DeviceId	CStimulusFunctionNet::SidebandData, 624
CMcsUsbListEntryNet, 238	
DeviceIdNet, 593	ElectricalStimulation
BcdDevice, 595	HeadStageIDType, 607
BusType, 595	EmptyKey
DeviceIdNet, 594	CMcsUsbNet, 248
IdProduct, 595	Emu_GetCellCapacity
IdVendor, 595	CRoboDacqNet, 382
operator=, 594	Emu_GetCellPotential
DeviceName	CRoboDacqNet, 382
CMcsUsbListEntryNet, 238	Emu_GetCellResists
DeviceRemoval	CRoboDacqNet, 382
CMcsUsbListNet, 241	Emu_GetElectrodeResists
DigitalSource	CRoboDacqNet, 383
DigitalSource< digitalsourceenum >, 595	Emu_GetNoise
DigitalSource< digitalsourceenum >, 595	CRoboDacqNet, 383
DigitalSource, 595	Emu_SetCellCapacity
MaxBitNumber, 596	CRoboDacqNet, 383
MaxBitNumberStatic, 596	Emu_SetCellPotential
size, 596	CRoboDacqNet, 383
Source, 596	Emu_SetCellResists
DigitalSourceGeneral, 596	CRoboDacqNet, 383
DigitalSourceGeneral, 597	Emu_SetElectrodeResists
MaxBitNumber, 597	CRoboDacqNet, 383
size, 597	Emu_SetNoise
Source, 597	CRoboDacqNet, 383
DisableAutoReset	EnableAnalogOut
CStg200xDownloadBasicNet, 471	CSCUFunctionNet, 426
DisableMultiFileMode	EnableAutoReset
CStg200xDownloadNet, 480	CStg200xDownloadBasicNet, 472
Disconnect	EnableChannelsInGroup
CMcsUsbNet, 248	CCMOSMea_FunctionNet, 36, 37
DisConnectDevice	CDacqGroupChannelSelectionTemplateNet< Dac
CRadioControledDevicesNet, 367	qGroupChannelEnumTemplateNet, Dac
DoRamp	qGroupChannelEnumTemplate, CDevice
CRoboDacqNet, 382	GroupChannelInfoTemplateNet >, 54
DownloadFirmware	EnableChecksum
CMcsUsbFactoryNet, 226	CMeaDeviceNet, 280
DriverVersionNet, 598	COctoPotDeviceNet, 321
~DriverVersionNet, 599	EnableDigitalIn
Driver VersionNet, 599	CMeaDeviceNet, 280, 281
DriverVersionNet::FormatVersion, 599	COctoPotDeviceNet, 321
GetDestinationCode, 599	EnableExceptions
GetDestinationName, 599, 600	CMcsUsbNet, 248
Goldosination value, 555, 666	EnableMultiFileMode

CStg200xDownloadNet, 480	CMeaFeedbackFunctionNet, 290
EnableTimestamp	FeedbackSetIIRFilterParameter
CMeaDeviceNet, 281	CMeaFeedbackFunctionNet, 290
COctoPotDeviceNet, 321	FeedbackSetLogic
EnableUserTrigger	CMeaFeedbackFunctionNet, 290
CLIH3DeviceNet, 123	FeedbackSetMkFilter
enCMosMeaChipType	CMeaFeedbackFunctionNet, 290
Mcs::Usb, 27	FeedbackSetNumberOfLogics
EnSTG200x_STATUS	CMeaFeedbackFunctionNet, 291
Mcs::Usb, 27	FeedbackSetNumberOfRateCounter
Entry	CMeaFeedbackFunctionNet, 291
HeadStageIDType, 608	FeedbackSetNumberOfRateDetectors
Equals	CMeaFeedbackFunctionNet, 291
CMcsUsbListEntryNet, 235	FeedbackSetNumberOfSpikeDetectors
HeadStageIDType, 607	CMeaFeedbackFunctionNet, 291
HeadstageIDTypeObject, 610	FeedbackSetNumberOfTriggers
EraseEepromRegisterPreconfig	CMeaFeedbackFunctionNet, 291
CMcsUsbNet, 248	FeedbackSetRateCounter
EraseFilterParameterPermanent	CMeaFeedbackFunctionNet, 291
CFilterConfigurationNet, 69	FeedbackSetRateDetector
CFilterConfigurationRegisterNet, 71	CMeaFeedbackFunctionNet, 291
ErasePermanentAdcOffset	FeedbackSetSpikeDetectorThreshold
CLIH3DeviceNet, 124	CMeaFeedbackFunctionNet, 291
ErasePermanentDacOffset	FeedbackSetTrigger
CLIH3DeviceNet, 124 From Callback, Agricultura, Standard	CMeaFeedbackFunctionNet, 292 FilterActive
Error_Callback_Aquisition_Stopped	
CMcsUsbDacqNet, 220	CFilterPropertyNet, 73 FilterBand
Error_Callback_Data_lost	
CMcsUsbDacqNet, 220	CFilterPropertyNet, 73
Error_Callback_Frames_Lost	FilterFamily CFilterPreparts Not. 72
CMcsUsbDacqNet, 220	CFilterPropertyNet, 73
Error_Callback_Packet_Error	FilterType
CMcsUsbDacqNet, 220	CFilterPropertyNet, 73
Error_Callback_Queue_Full	FindFilter
CMcsUsbDacqNet, 220	CCreateFilterNet, 49
Error_Callback_RingQueue_Full	FindFirmwareVersionMagicInBuffer
CMcsUsbDacqNet, 220	CMcsUsbFactoryNet, 226
ErrorEvent	FindReference
CMcsUsbDacqNet, 221	CRoboDeviceNet, 396
FactoryReset	FindReferencel
CTcxDeviceNet, 501	CRoboStatorDeviceNet, 414
FeedbackGetSampleTimerCount	FindReferencePhase0
CMeaFeedbackFunctionNet, 289	CRoboDeviceNet::RoboMainLowLevelCommands,
FeedbackSetAnalogSource	615
CMeaFeedbackFunctionNet, 289	FindReferencePhase0XY
FeedbackSetChannelFilter	CRoboStatorDeviceNet::RoboMainStatorLowLevelCommands,
CMeaFeedbackFunctionNet, 289	621
FeedbackSetDigitalMapping	FindReferenceXY
CMeaFeedbackFunctionNet, 289	CRoboStatorDeviceNet, 414, 415
FeedbackSetFeedback	FindReferenceZ
	CRoboStatorDeviceNet, 415
CMeaFeedbackFunctionNet, 289 FeedbackSetFilterOff	FirePressurePulse
CMeaFeedbackFunctionNet, 290	CPPCFunctionNet, 349
FeedbackSetFilterParameter	FirmwareDestinationNames, 603
CMeaFeedbackFunctionNet, 290	Altera, 603
FeedbackSetFilterParameter32	Bootstrap, 603
CMeaFeedbackFunctionNet, 290	BUS1_MCSBUS1, 603
FeedbackSetGlobalChannelFilter	BUS1_MCSBUS2, 604
i eeubauksetsiubaioilailileifiitei	

P.O.P. 404	
DSP, 604	Gain
FPGA2, 604	CMeaDeviceNet, 285
FPGA3, 604	Get2AnalogInput
FPGA4, 604	CMcsBus_SensorNet, 156
FPGA5, 604	Get2DigitalInput
FPGA6, 604	CMcsBus_SensorNet, 156
MCSBUS1, 604	Get4ADC
MCSBUS10, 604	CMcsBus_SensorNet, 156
MCSBUS11, 604	Get4ADCAverage
MCSBUS12, 604	CMcsBus_SensorNet, 156
MCSBUS13, 605	Get4ADCCatchampAverageShift
MCSBUS2, 605	CMcsBus_SensorNet, 156
MCSBUS3, 605	Get4ADCMode
MCSBUS4, 605	CMcsBus_SensorNet, 156
MCSBUSS, 605	Get4DAC
MCSBUS6, 605	CMcsBus_SensorNet, 156
MCSBUS7, 605	GetAbsMaxCurrentInMicroAmp
MCSBUS8, 605	CMultiwellOptoStimFunctionNet, 315
MCSBUS9, 605	GetAccelGyroCurrentRate
MCU1, 605	CW2100_FunctionNet, 529
PIC, 605	GetAccelGyroDesiredRate
PIC2, 606	CW2100_FunctionNet, 529
PIC3, 606	GetAccelGyroEnabled
PIC4, 606	CW2100_FunctionNet, 529
USB, 606	GetAccelRange
FluidControlDevice	CW2100_FunctionNet, 529
CSafeISDeviceNet, 422	GetActiveRunningTableNumber
ForceStatusEvent	CWarnerValveControllerDeviceNet, 565
CStg200xDownloadBasicNet, 472	GetActiveRunningTableNumberEvent
CStimulusFunctionNet, 488 FPGA2	CotAdenterCode
FirmwareDestinationNames, 604	GetAdapterCode CMeaImpedanceDeviceNet, 293
FPGA3	CTEERFunctionNet, 515
FirmwareDestinationNames, 604	GetAdapterType
FPGA4	CMcsUsbDacqNet, 196
FirmwareDestinationNames, 604	GetAdc
FPGA5	CFluidControlDeviceNet, 76
FirmwareDestinationNames, 604	GetAdcDataFormat
FPGA6	CMcsUsbDacqNet, 196
FirmwareDestinationNames, 604	GetADCInputOffset
FromIntPtr	CCMOSMea_FunctionNet, 37
StgStatusNet, 624	GetAdcOffset
FromPtr	CLIH3DeviceNet, 124
StgStatusNet, 624	COctoPotDeviceNet, 321
FX3MCSDataAddress	GetAdcOffsetU1
CMcsUsbFactoryNet, 231	CTEERFunctionNet, 515
FX3MCSDataDeviceIdOffset	GetAdcOffsetU2
CMcsUsbFactoryNet, 231	CTEERFunctionNet, 515
FX3MCSDataIFB1ImageOffset	GetADCs
CMcsUsbFactoryNet, 231	CMcsBus_SensorNet, 156
FX3MCSDataIFB2ImageOffset	GetADCsLoop
CMcsUsbFactoryNet, 231	CMcsBus_SensorNet, 157
FX3MCSDataVersionOffset	GetAdcZero
CMcsUsbFactoryNet, 231	CMcsUsbDacqNet, 196
FYIProgram	GetAirpressure
CFYIDeviceNet, 81	CRoboDeviceNet, 397
FYITemp	GetAirpressureLimit
CFYIDeviceNet, 81	CRoboDeviceNet, 397
51 115011001101, 01	5.13555 577661761, 50 7

0-14:4/-1	ODEE wasting Net 070
GetAirValve	CRFFunctionNet, 373
CRoboDeviceNet, 397	GetAvailableDeviceListEx
GetAllDigout	CRFFunctionNet, 373
CRoboDacqNet, 383	GetAvailableHeadstages
GetAmplification	CSCUFunctionNet, 426
CPgaDeviceNet, 333	CW2100_FunctionNet, 530
GetAmplitude_nA	GetAvailableHeadstagesEvent
CTEERFunctionNet, 515	CSCUFunctionNet, 437
GetAnalogGain	GetAvailableMemory
CMeaDeviceNet, 282	CStg200xBasicNet, 444
GetAnalogOutADCRange	CStimulusFunctionNet, 488
CSCUFunctionNet, 426	GetAvailableSampleRates
GetAnalogOutChannel	CMcsUsbDacqNet::CHWInfo, 113
CW2100_FunctionNet, 530	GetAvailableStateList
GetAnalogOutChannels	CRFFunctionNet, 373
CSCUFunctionNet, 426	GetAvailableStateListEx
GetAnalogOutDACRange	CRFFunctionNet, 373
CSCUFunctionNet, 426	GetAvailableVoltageRangesInMicroVolt
GetAnalogOutFilter	CMcsUsbDacqNet::CHWInfo, 113
CW2100 FunctionNet, 530	GetAvailableVoltageRangesInMicroVoltAndStringsInMilliVolt
GetAnalogRanges	CMcsUsbDacqNet::CHWInfo, 113
CStg200xBasicNet, 443	GetAxisConfig
GetAnalogResolution	CRoboDeviceNet::RoboMainLowLevelCommands.
CStg200xBasicNet, 444	615
GetAnalogThresholdHigh	GetAxisLED
CWarnerValveControllerDeviceNet, 565	CRoboocyte2DeviceNet, 412
GetAnalogThresholdHighEvent	GetAxisParametersSignedEeprom
CWarnerValveControllerDeviceNet, 582	CMcsBus_AxisParametersNet, 131
GetAnalogThresholdLow	GetAxisParametersUnsignedEeprom
	- ,
CWarnerValveControllerDeviceNet, 566	CMcsBus_AxisParametersNet, 131
GetAnalogThresholdLowEvent	GetBaseFrequency
CWarnerValveControllerDeviceNet, 582	CRFFunctionNet, 374
GetAnalogValueUnit	GetBaseSamplerate
CMcsUsbDacqNet, 196	CCMOSMeaDeviceNet, 46
GetAnalogVoltage	GetBath
CPPCFunctionNet, 349	CCMOSMea_FunctionNet, 37
CPPS_FunctionNet, 358	GetBathMode
CWarnerValveControllerDeviceNet, 566	CCMOSMea_FunctionNet, 37
GetAnalogVoltageEvent	GetBatteryState
CWarnerValveControllerDeviceNet, 582	CW2100_FunctionNet, 530
GetAnalogVoltageRange	GetBatteryVoltage
CPPCFunctionNet, 349	CMultiBatteryChargerDeviceNet, 301
GetAnalogVoltages	GetBiQuad
CPPS_FunctionNet, 358	CCreateFilterNet, 49
GetArraySize	GetBiQuads
CMealmpedanceDeviceNet, 293	CCreateFilterNet, 49
GetAudioChannels	GetBlankingEnable
CMeaAudioFunctionNet, 266, 267	CStg200xBasicNet, 444, 445
CW2100_FunctionNet, 530	GetBoardTemp
GetAudioOutDacParameter	CTcxDeviceNet, 502
CLIH3DeviceNet, 124	GetBoostAlwaysOnMode
GetAutocalibrationDisabled	CW2100_StimulatorFunctionNet, 536
CStg200xBasicNet, 444	GetBoostPreTime
GetAvailableBaseSamplerates	CW2100_StimulatorFunctionNet, 536
CCMOSMeaDeviceNet, 46	GetBubbleState
GetAvailableChambers	
	CPPS_FunctionNet, 358
CWarnerUssingFunctionNet, 547	GetBubbleStatus
GetAvailableDeviceList	CMcsBus SensorNet, 157

GetBuffer	CMultiBatteryChargerDeviceNet, 302
CGenericDevelopDeviceNet, 88	GetChargingPCoefficient
GetBusAddress	CMultiBatteryChargerDeviceNet, 302
CMcsBusNet, 170	GetChecksumFromFX3Image
GetBusAddressEeprom	CMcsUsbFactoryNet, 226
CMcsBusNet, 170	GetCheckVoltage
GetByteBuffer	COkuvisionStimulatorDeviceNet, 325
CGenericDevelopDeviceNet, 89	GetClampAmpSerialNumber
GetBytesAvailable	CRoboDacqNet, 384
CSerialPortNet, 438	GetClampMode
GetBytesPerSample	CTEERFunctionNet, 516
CTEERFunctionNet, 515	CWarnerUssingFunctionNet, 547
GetCalibration	GetCMOSDataDictionary
CTcxDeviceNet, 502	CCMOSMeaDeviceNet, 47
GetCalibrationDecp	GetCoilCommunication
CTcxDeviceNet, 502	CPositionIIDeviceNet, 336
GetCalibrationMax	GetColorRgb
CTcxDeviceNet, 502	CMultiwellOptoStimFunctionNet, 315
GetCalibrationMin	GetColorStr
CTcxDeviceNet, 502	CMultiwellOptoStimFunctionNet, 316
GetCapacityC	GetCommand
CRoboDacqNet, 383	CMcsBusNet, 170, 171
GetCapacityV	CPedoterDeviceNet, 330
CRoboDacqNet, 383	CRoboDacqNet, 384
GetCapacityX	GetComplianceVoltageThreshold
CRoboDacqNet, 384	CWarnerUssingFunctionNet, 547
GetCardinalDacqSamplerate	GetConfiguration
CInterfaceboardFunctionNet, 120	CMcsUsbNet, 249
GetCardinalStgOutputrate	GetConfigurationBit
CInterfaceboardFunctionNet, 120	CRoboDacqNet, 384
GetChannel	GetConfigurationBitAxc
CSw2to64DeviceNet, 497	CRoboDacqNet, 384
GetChannelDataFillSize	GetConfigurationBitBlu_Led
CMcsUsbDacqNet, 196	CRoboDacqNet, 384
GetChannelDatal16	GetConfigurationBitBlu_LedToggleFast
CCMOSMeaDeviceNet, 46	CRoboDacqNet, 384
GetChannelDatal32	GetConfigurationBitBlu_LedToggleSlow
CCMOSMeaDeviceNet, 46	CRoboDacqNet, 384
GetChannelDataUI16	GetConfigurationBitCC Gen
CCMOSMeaDeviceNet, 46	CRoboDacqNet, 384
GetChannelDataUI32	GetConfigurationBitCV_Gen
CCMOSMeaDeviceNet, 47	CRoboDacqNet, 384
GetChannelLayout	GetConfigurationBitRC_Gen
CMcsUsbDacqNet, 196	CRoboDacqNet, 384
GetChannels	GetConfigurationBitRed Led
CMultiBatteryChargerDeviceNet, 301	CRoboDacqNet, 385
CSw2to64DeviceNet, 498	GetConfigurationBitRed_LedSaturation
GetChannelsCountOfChamber	CRoboDacqNet, 385
CWarnerUssingFunctionNet, 547	GetConfigurationBitRed_LedToggleFast
GetChannelsInBlock	CRoboDacqNet, 385
CMcsUsbDacqNet, 197	GetConfigurationBitRed_LedToggleSlow
GetChannelState	CRoboDacqNet, 385
CMultiBatteryChargerDeviceNet, 301	GetConfigurationBitRelais
GetChargeCapacity	CRoboDacqNet, 385
CMultiBatteryChargerDeviceNet, 302	GetConfigurationBitRV_Gen
GetChargeCurrent	CRoboDacqNet, 385
CMultiBatteryChargerDeviceNet, 302	GetConfigurationBits
GetChargingMode	CRoboDacqNet, 385
~ ~	• • • •

GetConfigurationBitStream	CGrapheneFunctionNet, 102, 103
CRoboDacqNet, 385	COkuvisionStimulatorDeviceNet, 325
GetConfigurationBitSupply	GetDacOffset
CRoboDacqNet, 385	CDacCalibrationFunctionNet, 51
GetConnectedDevice	CLIH3DeviceNet, 125
CRFFunctionNet, 374	COctoPotDeviceNet, 321
GetControllerParams	GetDacPampsPerDigitHighCurrentRange
CTEERFunctionNet, 516	CWarnerUssingFunctionNet, 548
GetCrossTalkOffset	GetDacPampsPerDigitLowCurrentRange
CRoboDacqNet, 385	CWarnerUssingFunctionNet, 548
GetCrossTalkOptimum	GetDacqRunStatus
CRoboDacqNet, 385	CLIH3DeviceNet, 125
GetCurrent	GetDacRange
CTcxDeviceNet, 502	CW2100_FunctionNet, 530
GetCurrentAirvalve	GetDACResolution
CRoboDeviceNet, 397	CStg200xBasicNet, 446
GetCurrentAirvalveLimit	CStimulusFunctionNet, 489
CRoboDeviceNet, 397	CW2100_StimulatorFunctionNet, 537
GetCurrentCycle	GetDACs
CMeaCoatDeviceNet, 274	CMcsBus_SensorNet, 157
GetCurrentEditTableNumber	GetDacUseIdleValue
CWarnerValveControllerDeviceNet, 566	CLIH3DeviceNet, 125
GetCurrentEnable	GetDacZero
CTEERFunctionNet, 516	CTEERFunctionNet, 516
GetCurrentFactor	CWarnerUssingFunctionNet, 548
COkuvisionStimulatorDeviceNet, 325	GetDataFormat
GetCurrentNumberOfValves	CMcsUsbDacqNet, 197
CWarnerValveControllerDeviceNet, 566	GetDataMode
GetCurrentNumberOfValvesEvent	CMcsUsbDacqNet, 197
CWarnerValveControllerDeviceNet, 582	GetDDecp
GetCurrentPosition	CTcxDeviceNet, 502
CRoboDeviceNet, 397	GetDebugData
GetCurrentPositionI	CPositionIIDeviceNet, 336
CRoboStatorDeviceNet, 415	GetDestination
GetCurrentPositionXY	CMcsUsbFactoryNet, 226
CRoboStatorDeviceNet, 415	GetDestinationCode
GetCurrentPositionZ	DriverVersionNet, 599
CRoboStatorDeviceNet, 415	GetDestinationDisplayLabel
GetCurrentRangeInNanoAmp	CMcsUsbFactoryNet, 226
CStg200xBasicNet, 445	GetDestinationName
CStimulusFunctionNet, 489	CMcsUsbFactoryNet, 226, 227
CW2100_StimulatorFunctionNet, 536	DriverVersionNet, 599, 600
GetCurrentResolutionInNanoAmp	GetDestinationSerialNumber
CStg200xBasicNet, 445	CMcsUsbFactoryNet, 227
CStimulusFunctionNet, 489	GetDestinationTargetAddress
CW2100_StimulatorFunctionNet, 537	CMcsUsbFactoryNet, 227
	GetDetectionThreshold
GetCycle CMacClean Davisa Nat. 270	
CMeaCleanDeviceNet, 270	CMcsBus_SensorNet, 157
GetCycles CMagClean Davisa Not. 270	GetDetectorValue
CMeaClearDeviceNet, 270	CMcsBus_SensorNet, 157 GetDevice
CotD	
GetD	CTcxDeviceNet, 502
CTcxDeviceNet, 502	GetDeviceCannotStallOutRequests
GetDacAmplificationFactor	CMcsUsbNet, 249
CStg200xBasicNet, 446	GetDeviceCapableSpeed
GetDacIdleValue	CMcsUsbNet, 249
CLIH3DeviceNet, 125	GetDeviceEnum
GetDACOffset	CMcsUsbNet, 249

GetDeviceGroupChannelInfos	CMultiBatteryChargerDeviceNet, 303
CDacqGroupChannelSelectionTemplateNet< Dac-	GetDischargeCurrent
qGroupChannelEnumTemplateNet, Dac-	CMultiBatteryChargerDeviceNet, 303
qGroupChannelEnumTemplate, CDevice-	GetDischargeCurrentSetPoint
GroupChannelInfoTemplateNet >, 54	CMultiBatteryChargerDeviceNet, 303
GetDeviceId	GetDisplayMode
CMcsUsbNet, 249	CWarnerValveControllerDeviceNet, 568
GetDeviceList	GetDisplayModeEvent
CPositionImpDeviceNet, 344	CWarnerValveControllerDeviceNet, 583
GetDeviceNames	GetDisplayText
CRadioControledDevicesNet, 367	CRoboDacqNet, 386
GetDeviceRootHubVendorEnum	GetDMax
CMcsUsbNet, 249	CTcxDeviceNet, 503
GetDeviceRootHubVendorID	GetDMin
CMcsUsbNet, 249	CTcxDeviceNet, 503
GetDeviceRootHubVendorName	GetDownsampleFactor
CMcsUsbNet, 249	CRoboDacqNet, 386
GetDeviceSpeed	GetDSPHighPassByIndex
CMcsUsbNet, 250	CIntanMea_FunctionNet, 116
GetDeviceType	GetDuration
CTcxDeviceNet, 503	CMeaCoatDeviceNet, 274
GetDevname	GetDuty
CTcxDeviceNet, 503	CTcxDeviceNet, 503
GetDigin	GetEEpromPage
CFluidControlDeviceNet, 76	CLIH3DeviceNet, 126
GetDigInState	GetElectrodeDacMux
CLIH3DeviceNet, 126	CStg200xBasicNet, 447
GetDiginValue	GetElectrodeEnable
CStg200xBasicNet, 446	CStg200xBasicNet, 448
GetDigitalData	GetElectrodeMode
CMeaDigitalDataFunctionNet, 287	CStg200xBasicNet, 449
GetDigitalIn	GetEnableAmplifierProtectionSwitch
CPPCFunctionNet, 351	CStg200xBasicNet, 449, 450
CPPS_FunctionNet, 358	GetEnabledChannelsInGroup
GetDigitalOutPortValve	CCMOSMea_FunctionNet, 37
CWarnerValveControllerDeviceNet, 567	CDacqGroupChannelSelectionTemplateNet< Dac
GetDigitalOutPortValveEvent	qGroupChannelEnumTemplateNet, Dac
CWarnerValveControllerDeviceNet, 583	qGroupChannelEnumTemplate, CDevice
GetDigitalPortDirection	GroupChannelInfoTemplateNet >, 54
CWarnerValveControllerDeviceNet, 568	GetEnableHeaterLimit
GetDigitalPortDirectionEvent	CTcxDeviceNet, 503
CWarnerValveControllerDeviceNet, 583	GetEnableThermocouple
GetDigitalSource	CTcxDeviceNet, 503
CMcsUsbDacqNet, 197-199	GetEntry
GetDigitalStimulatorTrigger	CMcsUsbListEntryNet, 236
CW2100_StimulatorFunctionNet, 537	GetEntryCount
GetDigitalStimulatorTriggerSlope	CMcsUsbListEntryNet, 237
CW2100_StimulatorFunctionNet, 537	GetEnumerationSpeed
GetDigout	CMeaDeviceNet, 282
CFluidControlDeviceNet, 76	GetErrorAirpressure
CRoboDacqNet, 386	CRoboDeviceNet, 398
GetDigoutMode	GetErrorCurrentAirvalve
CStg200xBasicNet, 446	CRoboDeviceNet, 398
GetDigoutValue	GetErrorText
CStg200xBasicNet, 447	CMcsUsbNet, 250
GetDIO	GetErrorVoltage12V
CMcsBus_FYIExtensionNet, 134	CRoboDeviceNet, 398
GetDischargeCapacity	GetErrorVoltage5V

CRoboDeviceNet, 398	CCMOSMea FunctionNet, 37
GetErrorVoltageAirvalve	GetGroupChannelBitmaskBySelect
CRoboDeviceNet, 398	CCMOSMea_FunctionNet, 38
GetErrorVoltageRs485A	GetGroupChannelBitmaskHS1NCBathCurrent
CRoboDeviceNet, 398	CCMOSMea FunctionNet, 38
GetErrorVoltageRs485B	GetGroupChannelBitmaskHS1NCCol2Current
CRoboDeviceNet, 398	CCMOSMea_FunctionNet, 38
GetErrorVoltageValves	GetGroupChannelBitmaskHS1NChipTemp
CRoboDeviceNet, 398	CCMOSMea FunctionNet, 38
GetEventData	GetGroupChannelBitmaskHS1Sidebands
CPositionIIDeviceNet, 337	CCMOSMea_FunctionNet, 39
GetExternalElectrodeEnable	GetGroupChannelBitmaskHS1TriggerStatus
CStg200xBasicNet, 450	CCMOSMea_FunctionNet, 39
GetFAAmplification	GetGroupChannelBitmaskIFDigChannels
CStg200xBasicNet, 451	CCMOSMea_FunctionNet, 39
GetFilter	GetGroupChannelBitmaskInterfaceADC
CRoboDacqNet, 386	CCMOSMea FunctionNet, 39
GetFilterAttributes	GetGroupChannelBitmaskPacketFrameContext
CFilterConfigurationNet, 69	CCMOSMea FunctionNet, 40
GetFilterCoeffs	GetGroupChannelBitmaskSTG1DACSignal
CRoboDacqNet, 386	CCMOSMea_FunctionNet, 40
GetFilterParametersHeadstage	GetGroupChannelDatal16
CWClassicFunctionNet, 589	CMcsUsbDacqNet, 199
GetFilterProperties	GetGroupChannelDatal32
CSCUFunctionNet, 427	CMcsUsbDacqNet, 200
CW2100_FunctionNet, 530	GetGroupChannelDataUI16
GetFilterProperty	CMcsUsbDacqNet, 200
CMcsUsbDacqNet, 199	GetGroupChannelDataUI32
CSCUFunctionNet, 427	CMcsUsbDacqNet, 201
CW2100_FunctionNet, 530	GetGroupDCOffset
GetFinalDischargeVoltage	CCMOSMea_FunctionNet, 40
CMultiBatteryChargerDeviceNet, 304	GetGroupID
GetFirmwareVersion	CCMOSMea_FunctionNet, 40
CMcsUsbNet, 250	CDacqGroupChannelSelectionTemplateNet< Dac-
GetFirmwareVersionFromFile	qGroupChannelEnumTemplateNet, Dac-
CMcsUsbFactoryNet, 227	qGroupChannelEnumTemplate, CDevice-
GetFirmwareVersionFromHexFile	GroupChannelInfoTemplateNet >, 54, 55
CMcsUsbFactoryNet, 227	GetGroupNumberOfChannels
GetFPGAFirmwareType	CCMOSMea_FunctionNet, 41
CW2100_FunctionNet, 530	CDacqGroupChannelSelectionTemplateNet< Dac-
GetFrameErrorCounter	qGroupChannelEnumTemplateNet, Dac-
CTEERFunctionNet, 516	qGroupChannelEnumTemplate, CDevice-
GetFrequency	GroupChannelInfoTemplateNet >, 55
CRadioControledDevicesNet, 368	GetGroupResolutionPerDigit
GetFrequencyRange	CCMOSMea_FunctionNet, 41
CPgaDeviceNet, 333	GetGroupSampleSize
GetGain	CCMOSMea_FunctionNet, 41
CMeaDeviceNet, 282	CDacqGroupChannelSelectionTemplateNet< Dac-
CPgaDeviceNet, 334	qGroupChannelEnumTemplateNet, Dac-
GetGate	qGroupChannelEnumTemplate, CDevice-
CCMOSMea_FunctionNet, 37	GroupChannelInfoTemplateNet >, 55
GetGilsonDevice	GetGroupType
CRoboocyte2DeviceNet, 412	CCMOSMea_FunctionNet, 41
GetGlobalRepeat	CDacqGroupChannelSelectionTemplateNet< Dac-
CDigOutStimulatorFunctionNet, 61	qGroupChannelEnumTemplateNet, Dac-
GetGNDI	qGroupChannelEnumTemplate, CDevice-
CCMOSMea_FunctionNet, 37	GroupChannelInfoTemplateNet >, 55
GetGroupADCBits	GetGroupUnit
-	ereeer eegs e

CCMOSMea_FunctionNet, 41, 42	CTcxDeviceNet, 504
GetGyroRange	GetHeaterTemp
CW2100_FunctionNet, 530	CTcxDeviceNet, 504
GetHardwareMaxRange	GetHighCurrentRange
CMcsUsbDacqNet, 202	CWarnerUssingFunctionNet, 549
GetHardwareMinRange	GetHighpassFilterEnable
CMcsUsbDacqNet, 202	CFilterConfigurationNet, 69
GetHardwareRevision	GetHWConfig
CMcsUsbNet, 250	CRoboDeviceNet::RoboMainLowLevelCommands,
GetHasChecksum	615
CWClassicFunctionNet, 589	GetHWRevision
GetHashCode	CRoboDeviceNet::RoboMainLowLevelCommands,
HeadstageIDTypeObject, 610	616
GetHasRedLedHeadstage	GetHWRevisionEeprom
CWClassicFunctionNet, 589	CMcsBusNet, 171
GetHasThermocouple	Getl
CTcxDeviceNet, 503	CTcxDeviceNet, 504
GetHeadstage	GetIC
CStg200xBasicNet, 451	CRoboDacqNet, 386
GetHeadstageActive	GetIClamp
CMcsUsbNet, 251	CRoboDacqNet, 386
GetHeadstageAdcBits	GetlCoeff
CSCUFunctionNet, 427	CRobo_FYITemp_FunctionNet, 378
GetHeadstageAdcRangeInMicroVolt	GetICOffset
CSCUFunctionNet, 428	CRoboDacqNet, 386
GetHeadstageDacBits	GetIDecp
CSCUFunctionNet, 428	CTcxDeviceNet, 504
GetHeadstageDacCurrentRangeInMicroAmpere	GetIdent
CSCUFunctionNet, 428	CMcsUsbNet, 251
GetHeadstageDacCurrentResolutionInNanoAmpere	GetIdleModeOffset
CSCUFunctionNet, 429	CWarnerUssingFunctionNet, 549
GetHeadstageDacVoltageRangeInMilliVolt	GetlGain
CSCUFunctionNet, 429	CRoboDacqNet, 386
GetHeadstageDacVoltageResolutionInMicroVolt	GetIMax
CSCUFunctionNet, 429	CTcxDeviceNet, 504
GetHeadstageGainInPermille	GetIMin
CSCUFunctionNet, 430	CTcxDeviceNet, 504
GetHeadstageID	GetImpedanceResult
CMcsUsbNet, 251	CIntanMea_FunctionNet, 116
CSCUFunctionNet, 430	GetImpedanceTestFrequency
GetHeadstageNumberOfAnalogChannels	CMealmpedanceDeviceNet, 293
CSCUFunctionNet, 430	GetImpld
GetHeadstageNumberOfStimulationChannels	CPositionImpDeviceNet, 345
CSCUFunctionNet, 431	GetImplantCurrentSetpoint
GetHeadstageOnOff	CPositionIIDeviceNet, 337
•	
CW2100_FunctionNet, 531	GetImplantResult
Cott Loadstage Payer State At Start	CPositionIIDeviceNet, 338
GetHeadstagePowerStateAtStart	GetImplantState
CSCUFunctionNet, 431	CPositionIIDeviceNet, 338 GetInMovement
GetHeadstagePresent CMcsUsbNet, 251	CRoboDeviceNet, 398
GetHeadstageSamplerate	GetIntanRegister
CSCUFunctionNet, 431	CIntanMea_FunctionNet, 116
GetHeadstageSamplingActive	GetIntBuffer Connection Devices Not 80
CW2100_FunctionNet, 531	CGenericDevelopDeviceNet, 89
GetHeadstageSerialNumber	GetIO
CSCUFunctionNet, 432	CWarnerValveControllerDeviceTesterFunctionNet,
GetHeaterLimit	586

GetlOut	CTcxDeviceNet, 505
CTcxDeviceNet, 504	GetMaxPower
GetloVoltage	COkuvisionStimulatorDeviceNet, 325
CInterfaceboard2FunctionNet, 118	CRobo_FYITemp_FunctionNet, 378
GetLastAnswer	GetMaxPressureWaitTime
CGilsonDeviceNet, 100	CRoboDeviceNet::RoboMainLowLevelCommands
GetLastUSBError	616
CMcsUsbNet, 252	GetMaxSamplingFrequency
GetLatency	CMcsUsbDacqNet, 202
CMcsBus_SensorNet, 157	GetMaxStimulusChannelsPerHeadstage
GetLatencyCounter	CSCUFunctionNet, 432
CMcsBus_SensorNet, 157	GetMaxVoltage
GetLayoutConfiguration	CMeaCleanDeviceNet, 270
CMEA2100x256FunctionNet, 265	COkuvisionStimulatorDeviceNet, 325
GetLEDSwitch	GetMCAcceleration
CMcsBus_ExtensionNet, 133	CMcsBus_MotorControlNet, 138
GetLength	GetMCAccelerationEeprom
CRobo_FYIProgram_FunctionNet, 376	CMcsBus_MotorControlNet, 138
GetLiquidResistance	GetMCAccelerationShortCommand
CTEERFunctionNet, 517	CMcsBus_MotorControlNet, 138
CWarnerUssingFunctionNet, 549	GetMCAxisRevisionEeprom
GetListmodeIndexRange	CMcsBus_MotorControlNet, 138
CStg200xBasicNet, 451	GetMCBreakCurrent
GetListmodeTriggerSource	CMcsBus_MotorControlNet, 138
CStg200xBasicNet, 451	GetMCBreakCurrentEeprom
GetLowCurrentRange	CMcsBus_MotorControlNet, 139
CWarnerUssingFunctionNet, 550	GetMCConfig
GetLowerFrequencyByIndex	CMcsBus_MotorControlNet, 139
CIntanMea_FunctionNet, 116	GetMCConfigEeprom
GetMajor	CMcsBus_MotorControlNet, 139
DriverVersionNet, 600	GetMCCurrent
GetMaxChunkSize_Byte	CMcsBus_MotorControlNet, 139
CTEERFunctionNet, 517	GetMCCurrentEeprom
GetMaxCurrent	CMcsBus_MotorControlNet, 139
CMeaCoatDeviceNet, 274	GetMCCurrentMode
GetMaxDurationHighCurrentInMicroSec	CMcsBus_MotorControlNet, 139
CMultiwellOptoStimFunctionNet, 316	GetMCCurrentModeEeprom
GetMaxDutyCycleHighCurrent	CMcsBus_MotorControlNet, 139
CMultiwellOptoStimFunctionNet, 316	GetMCCurrentModeShortCommand
GetMaxHeaterPowerMultiwell	CMcsBus MotorControlNet, 140
CTcxDeviceNet, 505	GetMCCurrentPosition
GetMaxNoPressure	CMcsBus_MotorControlNet, 140
CRoboDeviceNet::RoboMainLowLevelCommands,	GetMCCurrentShortCommand
616	CMcsBus_MotorControlNet, 140
GetMaxNoPressureWaitTime	GetMCCurrentSpeed
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBus_MotorControlNet, 140
616	GetMCMaxAcceleration
GetMaxNumberOfHeadstages	CMcsBus_MotorControlNet, 140
CSCUFunctionNet, 432	GetMCMaxAccelerationEeprom
GetMaxNumOfColumns	CMcsBus_MotorControlNet, 140
CCMOSMea_FunctionNet, 42	GetMCMaxCurrent
GetMaxP	CMcsBus_MotorControlNet, 140
CTcxDeviceNet, 505	GetMCMaxCurrentEeprom
GetMaxpDecp	CMcsBus_MotorControlNet, 141
CTcxDeviceNet, 505	GetMCMaxSpeed
GetMaxpMax	CMcsBus_MotorControlNet, 141
CTcxDeviceNet, 505	GetMCMaxSpeedEeprom
GetMaxpMin	CMcsBus MotorControlNet. 141

GetMCMaxTravel	GetMinNoPressureWaitTime
CMcsBus_MotorControlNet, 141	CRoboDeviceNet::RoboMainLowLevelCommands,
GetMCMaxTravelEeprom	616
CMcsBus_MotorControlNet, 141	GetMinor
GetMCMaxTravelShortCommand	DriverVersionNet, 600
CMcsBus MotorControlNet, 141	GetMinPressure
GetMCMovement	CRoboDeviceNet, 398
CMcsBus_MotorControlNet, 141	CRoboDeviceNet::RoboMainLowLevelCommands,
GetMCNewPosition	616
CMcsBus_MotorControlNet, 142	GetMinPressureWaitTime
GetMCOutputOnOff	CRoboDeviceNet::RoboMainLowLevelCommands,
CMcsBus MotorControlNet, 142	616
GetMCPhase	GetMinSamplingFrequencyStepsize
CMcsBus_MotorControlNet, 142	CMcsUsbDacqNet, 202
GetMCPhaseOffset	GetMinVoltage
CMcsBus_MotorControlNet, 142	CMeaCleanDeviceNet, 270
GetMCReference	GetModeSelect
CMcsBus MotorControlNet, 142	CPulseGeneratorFunctionNet, 365
GetMCReferenceCurrent	GetModuleCurrent
CMcsBus MotorControlNet, 142	
——————————————————————————————————————	CStg200xDownloadNet, 481
GetMCReferenceCurrentEeprom CMcsBus MotorControlNet, 142	GetModuleTemp
——————————————————————————————————————	CStg200xDownloadNet, 481
GetMCRegulatorGain	GetMovementError
CMcsBus_MotorControlNet, 143	CRoboDeviceNet, 398
GetMCRegulatorGainEeprom	GetMovePump
CMcsBus_MotorControlNet, 143	CMcsBus_SensorNet, 158
GetMcsBus_Extension	GetMultiHeadstageMode
CRoboocyte2DeviceNet, 412	CW2100_FunctionNet, 531
GetMCScalingFactor	GetMultiplexedDataChannelsInBlock
CMcsBus_MotorControlNet, 143	CStimulusFunctionNet, 489
GetMCScalingFactorEeprom	GetNanoVoltsPerKelvin
CMcsBus_MotorControlNet, 143	CMcsBus_TempSensorNet, 163
GetMCSpeed	GetNeurochipMemoryData
CMcsBus_MotorControlNet, 143	CCMOSMea_FunctionNet, 42
GetMCSpeedEeprom	GetNeurochipMemorySize
CMcsBus_MotorControlNet, 143	CCMOSMea_FunctionNet, 42
GetMCSpeedShortCommand	GetNIC_MS
CMcsBus_MotorControlNet, 143	CRoboDacqNet, 386
GetMCSpeedUnitEeprom	GetNUC_MS
CMcsBus_MotorControlNet, 144	CRoboDacqNet, 386
GetMCStandbyCurrent	GetNumAmplifications
CMcsBus_MotorControlNet, 144	CPgaDeviceNet, 334
GetMCStandbyCurrentEeprom	GetNumber
CMcsBus_MotorControlNet, 144	CMeaSwitchDeviceNet, 295
GetMCStandbyTime	CSw2to64DeviceNet, 498
CMcsBus_MotorControlNet, 144	GetNumberOfAnalogChannels
GetMCStandbyTimeEeprom	CStg200xBasicNet, 451
CMcsBus_MotorControlNet, 144	CStimulusFunctionNet, 490
GetMea21UsbPort	CW2100_StimulatorFunctionNet, 537
CMcsUsbNet, 252	GetNumberOfAudioChannels
GetMeaLayout	CMeaAudioFunctionNet, 267
CMcsUsbDacqNet, 202	GetNumberOfAvailableChambers
GetMemoryUsageDAC	CWarnerUssingFunctionNet, 550
CStg200xDownloadBasicNet, 472	GetNumberOfAvailableSamples
GetMemoryUsageSyncout	CTEERFunctionNet, 517
CStg200xDownloadBasicNet, 472	GetNumberOfChannels
GetMinimalThreshold	CDigOutStimulatorFunctionNet, 61
CMcsRus SensorNet 157	GetNumberOfDataRits

014 II I B	
CMcsUsbDacqNet, 203	CRoboDeviceNet::RoboMainLowLevelCommands,
GetNumberOfDevices	616
CMcsUsbListNet, 240	GetPattern
GetNumberOfHardwareSlotsForChambers	CMeaSwitchDeviceNet, 295
CWarnerUssingFunctionNet, 550	GetPatternBool
GetNumberOfHWADCChannels	CMeaSwitchDeviceNet, 296
CMcsUsbDacqNet::CHWInfo, 114	GetPauseDuration
GetNumberOfHWDACPaths	CMeaCoatDeviceNet, 275
CStg200xBasicNet, 451	GetPCoeff
GetNumberOfHWDigitalChannels	CRobo_FYITemp_FunctionNet, 378
CMcsUsbDacqNet::CHWInfo, 114	GetPDecp
GetNumberOfStimulationElectrodes	CTcxDeviceNet, 506
CStg200xBasicNet, 451	GetPeriod
GetNumberOfStimulationSourcesPerElectrode	CPulseGeneratorFunctionNet, 365
CStg200xBasicNet, 452	GetPeriod_us
GetNumberOfSupportedGroups	CTEERFunctionNet, 517
CCMOSMea FunctionNet, 42	GetPermanentCurrentInMicroAmp
CDacqGroupChannelSelectionTemplateNet< Dac-	CMultiwellOptoStimFunctionNet, 317
qGroupChannelEnumTemplateNet, Dac-	GetPGain
qGroupChannelEnumTemplate, CDevice-	CRoboDacqNet, 387
GroupChannelInfoTemplateNet >, 55	GetPhases
GetNumberOfSyncoutChannels	CRoboDeviceNet::RoboMainLowLevelCommands,
CStg200xBasicNet, 452	616
CW2100_StimulatorFunctionNet, 538	GetPicFirmwareType
GetNumberOfTriggerInputs	CW2100_FunctionNet, 531
	Gw2100_FunctionNet, 531 GetPiezoState
CStg200xBasicNet, 452	
CW2100_StimulatorFunctionNet, 538	CMcsBus_SensorNet, 158
GetNumConfigurations	GetPlateClampLockState
CMcsUsbNet, 252	CMultiwellDeviceNet, 310
GetNumControlChannels	GetPlateClampState
CTcxDeviceNet, 505	CMultiwellDeviceNet, 310
GetNumDestinations	GetPlateClampStateByHeadstage
CMcsUsbFactoryNet, 227	CMultiwellCallbackFunctionNet, 307
GetNumDevices	GetPlateClampStateByHeadstageEvent
CTcxDeviceNet, 505	CMultiwellCallbackFunctionNet, 308
GetNumEntries	GetPlateMux
DriverVersionNet, 601	CMultiwellDeviceNet, 310, 311
GetNumFrequencyRanges	GetPlateType
CPgaDeviceNet, 334	CMultiwellDeviceNet, 311
GetNumMeasureChannels	GetPMax
CTcxDeviceNet, 505	CTcxDeviceNet, 506
GetNUV_MS	GetPMin
CRoboDacqNet, 387	CTcxDeviceNet, 506
GetOffsetCurrent	GetPoti
CMeaCoatDeviceNet, 275	CMcsUsbDacqNet, 203
GetOnOff	GetPOut
CPositionIIDeviceNet, 338	CTcxDeviceNet, 506
CTcxDeviceNet, 506	GetPowerMuxPlate
GetOutputCurrent	CMultiwellDeviceNet, 311
CMeaCoatDeviceNet, 275	GetPowerStrength
GetOutputRate	CPositionIIDeviceNet, 339
CStg200xBasicNet, 452	GetPressure
GetOutputVoltage	CMcsBus_SensorNet, 158
	GetPressureOffset
CMeaCleanDeviceNet, 270	
GetP CTovPoviceNet 506	CotProcure Pongo
CTcxDeviceNet, 506	GetPressureRange
GetParameter	CPPCFunctionNet, 351
	GetPulseform

COkuvisionStimulatorDeviceNet, 325	GetRes2
GetPulseLength	CTcxDeviceNet, 507
CPulseGeneratorFunctionNet, 365	GetResetFilter
GetPumpCouple	CWClassicFunctionNet, 590
CPPS_FunctionNet, 358	GetResistanceC
GetPumpEnableSpeedRatio	CRoboDacqNet, 387
CPPS_FunctionNet, 358	GetResistanceV
GetPumpFastOnOff	CRoboDacqNet, 387
CPPS FunctionNet, 358	GetResolutionPerDigit
GetPumpFastSpeed	CMcsUsbDacqNet, 203
CPPS FunctionNet, 358	GetResS
GetPumpFunctionSpeeds	CTcxDeviceNet, 507
CPPS_FunctionNet, 359	GetResult
GetPumpManualOnOff	CMealmpedanceDeviceNet, 293
CPPS_FunctionNet, 359	GetResX
GetPumpMaxSpeed	CTcxDeviceNet, 507
CPPS_FunctionNet, 359	GetRFConnectionStatus
GetPumpModeType	CWClassicFunctionNet, 590
CPPCFunctionNet, 351	GetRFFrequency
CPPS_FunctionNet, 359	CPositionImpDeviceNet, 345
GetPumpSpeed	GetRFFrequencyHeadstage
CRoboFluidDeviceNet, 409	CWClassicFunctionNet, 590
GetPumpSpeedRatio	GetRFFrequencyReceiver
CPPS_FunctionNet, 359	CWClassicFunctionNet, 590
GetPumpSpeedUnit	GetRFPower
CPPCFunctionNet, 352	CWClassicFunctionNet, 590
CPPS_FunctionNet, 359	GetRoboDacq
GetPWM	CRoboocyte2DeviceNet, 412
	GetRoboFluidDevice
CFluidControlDeviceNet, 76 GetPwrOut	
	CEncapsulatorDeviceNet, 64
CTcxDeviceNet, 506 GetPwrSet	CRoboocyte2DeviceNet, 412
CTcxDeviceNet, 507	GetRotaryPositionCode
	CTEERFunctionNet, 517
GetRatedCapacity CMultiPatters Charges Payloa Not 2004	GetRotatePump
CMultiBatteryChargerDeviceNet, 304	CMcsBus_SensorNet, 159
GetReady CMaplemedanas Pavijas Nat. 200	GetROut 507
CMealmpedanceDeviceNet, 293	CTcxDeviceNet, 507
GetRecordingNumber	GetRTC
CRoboDacqNet, 387	COkuvisionStimulatorDeviceNet, 326
GetReferenceElectrodeMode	CPositionIIDeviceNet, 339
CSCUFunctionNet, 432	GetSampleBufferChunk
GetReferenceElectrodeSwitchState	CTEERFunctionNet, 517
CSCUFunctionNet, 433	GetSampleInterval
GetReferenceTemperature	CLIH3DeviceNet, 126
CFluidControlDeviceNet, 76	GetSamplePeriode
GetRegulationTimeouts	CMcsBus_SensorNet, 159
CMcsBus_SensorNet, 158	GetSampleRate
GetRegulatorFactor	CTEERFunctionNet, 518
CMcsBus_SensorNet, 159	GetSamplerate
GetRegulatorOnOff	CMcsUsbDacqNet, 203
CMcsBus_SensorNet, 159	GetSampleVoltageBuffer_uV
CRobo_FYITemp_FunctionNet, 378	CTEERFunctionNet, 518
GetRegulatorStatus	GetScaleFactorU1
CMcsBus_SensorNet, 159	CTEERFunctionNet, 518
GetRepeats	GetScaleFactorU2
CProgramPressureCurveNet, 363	CTEERFunctionNet, 518
GetRes1	GetScanHeadstagesResult
CTcxDeviceNet, 507	CWClassicFunctionNet, 590

0.10	0 10 % 1/ 01 1
GetScreen	GetSoftwareKeyString
CRoboDacqNet, 387	CMcsUsbNet, 252
GetSearchReferenceFastAccel	GetSollPressure
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBus_SensorNet, 159
617	GetSollTemp
GetSearchReferenceFastSpeed	CRobo_FYITemp_FunctionNet, 378
CRoboDeviceNet::RoboMainLowLevelCommands,	GetSourceBulk
617	CCMOSMea_FunctionNet, 42
GetSearchReferenceFineAccel	GetSourceDrain
CRoboDeviceNet::RoboMainLowLevelCommands,	CCMOSMea_FunctionNet, 42
617	GetSourceGate
GetSearchReferenceFineSpeed	CCMOSMea_FunctionNet, 42
CRoboDeviceNet::RoboMainLowLevelCommands,	GetStartTriggerSlope
617	CDigOutStimulatorFunctionNet, 61
GetSearchReferenceMethod	GetState
CRoboDeviceNet::RoboMainLowLevelCommands,	CRFFunctionNet, 374
617	CRobo_FYIProgram_FunctionNet, 377
GetSearchReferenceMoveOut	GetStateDebugData
CRoboDeviceNet::RoboMainLowLevelCommands,	CPositionIIDeviceNet, 339
617	GetStateEventData
GetSearchReferenceOffsetPos	CPositionIIDeviceNet, 340
CRoboDeviceNet::RoboMainLowLevelCommands,	GetStatus
617	CMcsUsbNet, 252
GetSelectedChannels	DriverVersionNet, 601
CW2100_FunctionNet, 531	GetStatusOfLastCommand
GetSelectedHeadstage	CMcsUsbNet, 253
CWClassicFunctionNet, 590	GetStgProgramInfo
GetSelectedHeadstageState	CStg200xBasicNet, 452, 453
CW2100_FunctionNet, 531	GetStgVersionInfo
GetSensorType	CStg200xBasicNet, 453
CTcxDeviceNet, 507	GetStimulationPatternMemory
GetSerialNumber	CW2100_StimulatorFunctionNet, 538
CMcsUsbNet, 252	GetStimulatorStatus
DriverVersionNet, 601	COkuvisionStimulatorDeviceNet, 326
GetSerialNumberHeadstage	GetStimulusParametersCache
CWClassicFunctionNet, 590	CW2100_FunctionNet, 531
GetSetpoint	GetStimulusParametersFromSelectedHS
CTcxDeviceNet, 507	CW2100 FunctionNet, 531
GetSetpointDecp	GetStimulusSites
CTcxDeviceNet, 508	CCMOSMea_FunctionNet, 43
	GetStiumlusParameters
GetSetpointMax	
CTcxDeviceNet, 508	CW2100_FunctionNet, 532
GetSetpointMin	GetStopTriggerSlope
CTcxDeviceNet, 508	CDigOutStimulatorFunctionNet, 62
GetShortBuffer	GetSubChannel
CGenericDevelopDeviceNet, 90	CMcsBus_MotorControlNet, 144
GetSimulation	GetSupplyVoltage
CRoboDacqNet, 387	CPPCFunctionNet, 352
GetSingleHeater	CPPS_FunctionNet, 359
CMcsBus_FYIExtensionNet, 134	GetSweepCount
GetSingleValve	CStg200xDownloadBasicNet, 473
CFluidControlDeviceNet, 77	GetSync
CRoboFluidDeviceNet, 409	CWarnerValveControllerDeviceTesterFunctionNet,
GetSlope	586
CMeaCleanDeviceNet, 271	GetSyncoutMap
	•
CMeaCoatDeviceNet, 275	CStg200xBasicNet, 453
GetSoftwareKey	GetSyncState
CMcsl lshNet 252	CMcsRus SensorNet 159

GetTableName	CWarnerUssingFunctionNet, 552
CWarnerValveControllerDeviceNet, 568	GetU2Offset
GetTableNamebyIndex	CWarnerUssingFunctionNet, 552
CWarnerValveControllerDeviceNet, 568	GetU2Reference
GetTableNamebyIndexEvent	CWarnerUssingFunctionNet, 552
CWarnerValveControllerDeviceNet, 583	GetUByteBuffer
GetTablepointer	CGenericDevelopDeviceNet, 91
CRetinaLedDeviceNet, 370	GetUC
GetTemperatur	CRoboDacqNet, 387
CMcsBus_TempSensorNet, 163	GetUClamp
GetTestMode	CRoboDacqNet, 387
CRFFunctionNet, 374	GetUCOffset
GetThermocoupleCalibration	CRoboDacqNet, 387
CFluidControlDeviceNet, 77	GetUintA
CTcxDeviceNet, 508	CFilterCoefficientsNet, 68
GetThermocoupleNanovoltPerKelvin	GetUintB
CFluidControlDeviceNet, 77	CFilterCoefficientsNet, 68
CTcxDeviceNet, 508	GetUIntBuffer
GetThermocoupleReferenceTemp	CGenericDevelopDeviceNet, 91
CTcxDeviceNet, 508	GetUnit
GetThermocoupleTemp	CTcxDeviceNet, 509
CTcxDeviceNet, 508	GetUnitDescription
GetThermocoupleTempAbs	CWarnerUssingFunctionNet, 553
CTcxDeviceNet, 509	GetUnitExponent
GetThermocoupleTemperature	CWarnerUssingFunctionNet, 553
CFIuidControlDeviceNet, 78	GetUnitName
GetThermoOffset	CWarnerUssingFunctionNet, 553
CMcsBus_TempSensorNet, 163	GetUnitsPerDigit
GetThermoTemp	CWarnerUssingFunctionNet, 554
CMcsBus_TempSensorNet, 164	GetUOut
GetThermoVoltage CMcsBus_TempSensorNet, 164	CTcxDeviceNet, 509 GetUpdateDisplay
GetTimeInPause	CRoboDacqNet, 387
CMeaCoatDeviceNet, 275	GetUpperFrequencyByIndex
GetTimeInPlateau	CIntanMea_FunctionNet, 116
CMeaCoatDeviceNet, 276	GetUptimeSeconds
GetTimeResolutionInNanoSeconds	CTEERFunctionNet, 519
CW2100_StimulatorFunctionNet, 538	CWarnerUssingFunctionNet, 554
GetTimeSlot	GetUSBDeviceIDFromFX3Image
CW2100_StimulatorFunctionNet, 538	CMcsUsbFactoryNet, 228
GetTotalMemory	GetUsbListEntries
CStg200xBasicNet, 454	CMcsUsbListNet, 240
CStimulusFunctionNet, 490	GetUsbListEntry
GetTotalNumberOfDigitalPorts	CMcsUsbListNet, 240
CWarnerValveControllerDeviceNet, 569	CMcsUsbNet, 253
GetTotalNumberOfTables	GetUseBubble
CWarnerValveControllerDeviceNet, 569	CPPS_FunctionNet, 359
GetTotalNumberOfValves	GetUsercodeFromBitFile
CWarnerValveControllerDeviceNet, 569	CMcsUsbFactoryNet, 228
GetTotalTableSize	GetUsercodeFromFlash
CWarnerValveControllerDeviceNet, 569	CMcsUsbFactoryNet, 228
GetTrigger	GetUserDefinedName
CStg200xDownloadBasicNet, 473	CW2100_FunctionNet, 532
GetTriggerSource	GetUserDefinedNameCache
CStg200xBasicNet, 454	CW2100_FunctionNet, 532
GetU1Offset	GetUserDefinedNameFromSelectedHS
CWarnerUssingFunctionNet 550	CW2100 FunctionNet 532
CWarnerUssingFunctionNet, 550 GetU1Reference	CW2100_FunctionNet, 532 GetUserParameter

CPahaDaviaaNat::PahaMaiaLawLavalCammanda	GetVDD3I
CRoboDeviceNet::RoboMainLowLevelCommands,	
617, 618 Catt ShortPuffor	CCMOSMea_FunctionNet, 43
GetUShortBuffer	GetVDDI
CGenericDevelopDeviceNet, 92	CCMOSMea_FunctionNet, 43
GetUV	GetVdVs
CRoboDacqNet, 388	CGrapheneFunctionNet, 103
GetUVOffset	GetVdVsDAC
CRoboDacqNet, 388	CGrapheneFunctionNet, 103, 104
GetValue	GetVersion
CTcxDeviceNet, 509	CMcsUsbNet, 253
GetValueHires	GetVersionInt
CTcxDeviceNet, 509	DriverVersionNet, 602
GetValve	GetVersionString
CFluidControlDeviceNet, 78	DriverVersionNet, 602
CRoboFluidDeviceNet, 409	GetVMMaxNegativeCurrent
GetValve1	CMcsBus_VoltageModeNet, 166
CRobo_FYIProgram_FunctionNet, 377	GetVMMaxNegativeCurrentEeprom
GetValve2	CMcsBus_VoltageModeNet, 166
CRobo_FYIProgram_FunctionNet, 377	GetVMMaxNegativeVoltage
GetValveActive	CMcsBus VoltageModeNet, 166
CPPCFunctionNet, 352	GetVMMaxNegativeVoltageEeprom
CWarnerValveControllerDeviceNet, 569	CMcsBus_VoltageModeNet, 166
GetValveActiveEvent	GetVMMaxPositiveCurrent
CWarnerValveControllerDeviceNet, 583	CMcsBus VoltageModeNet, 166
GetValveBoardRevision	GetVMMaxPositiveCurrentEeprom
CWarnerValveControllerDeviceNet, 570	CMcsBus_VoltageModeNet, 166
GetValveBoardRevisionEvent	GetVMMaxPositiveVoltage
CWarnerValveControllerDeviceNet, 583	CMcsBus_VoltageModeNet, 166
GetValveBoardRevisionString	GetVMMaxPositiveVoltageEeprom
CWarnerValveControllerDeviceNet, 570	CMcsBus_VoltageModeNet, 167
GetValveDigitalInPort	GetVMOutputOnOff
CWarnerValveControllerDeviceNet, 570	CMcsBus_VoltageModeNet, 167
	— -
GetValveDigitalInPortEvent	GetVMVoltage CMcsBus VoltageModeNet, 167
CWarnerValveControllerDeviceNet, 583	<u> </u>
GetValveLedOn	GetVoltage
CWarnerValveControllerDeviceNet, 570	COkuvisionStimulatorDeviceNet, 326
GetValveLedOnEvent	GetVoltage12V
CWarnerValveControllerDeviceNet, 584	CRoboDeviceNet, 399
GetValveManualGroup	GetVoltage12VLimit
CWarnerValveControllerDeviceNet, 571	CRoboDeviceNet, 399
GetValveManualGroupEvent	GetVoltage5V
CWarnerValveControllerDeviceNet, 584	CRoboDeviceNet, 399
GetValveManualState	GetVoltage5VLimit
CWarnerValveControllerDeviceNet, 571	CRoboDeviceNet, 399
GetValveManualStateEvent	GetVoltageAirvalve
CWarnerValveControllerDeviceNet, 584	CRoboDeviceNet, 399
GetValveMode	GetVoltageAirvalveLimit
CWarnerValveControllerDeviceNet, 571	CRoboDeviceNet, 399
GetValveModeEvent	GetVoltageClampControllerParam_D
CWarnerValveControllerDeviceNet, 584	CWarnerUssingFunctionNet, 554
GetValves	GetVoltageClampControllerParam_I
CMcsBus_FYIExtensionNet, 134	CWarnerUssingFunctionNet, 555
GetValvesActiveMap	GetVoltageClampControllerParam_P
CWarnerValveControllerDeviceNet, 571	CWarnerUssingFunctionNet, 555
GetValvesManualStateMap	GetVoltageRange
CWarnerValveControllerDeviceNet, 572	CGrapheneFunctionNet, 104
GetValveTableEntry	GetVoltageRangeIndex
CWarnerValveControllerDeviceNet, 572	CMcsUsbDacqNet, 203

GetVoltageRangeInMicroVolt	CSCUFunctionNet, 433
CMcsUsbDacqNet, 203	HasIMU
CStg200xBasicNet, 454	HeadStageIDType, 608
CStimulusFunctionNet, 490	HasOptoCurrentMessurement
CW2100_StimulatorFunctionNet, 538	HeadStageIDType, 608
GetVoltageRangeInMilliVolt	HasRadioControl
CMcsUsbDacqNet, 204	CRadioControledDevicesNet, 368
GetVoltageReached	HasRef
CGrapheneFunctionNet, 105	CRoboDeviceNet::RoboMainLowLevelCommands
GetVoltageResolution	618
CGrapheneFunctionNet, 105	HasRefl
GetVoltageResolutionInMicroVolt	CRoboStatorDeviceNet, 415
CStg200xBasicNet, 454	HasRefXY
CStimulusFunctionNet, 490	CRoboStatorDeviceNet, 415
CW2100_StimulatorFunctionNet, 538	HasRefZ
GetVoltageRs485A	CRoboStatorDeviceNet, 415
CRoboDeviceNet, 399	HasSoftwareKey
GetVoltageRs485ALimit	CMcsUsbNet, 253
CRoboDeviceNet, 399	HeadStageIDType, 606
GetVoltageRs485B	CompareTo, 607
CRoboDeviceNet, 399	ElectricalStimulation, 607
GetVoltageRs485BLimit	Entry, 608
CRoboDeviceNet, 400	Equals, 607
GetVoltageValves	HasIMU, 608
CRoboDeviceNet, 400	HasOptoCurrentMessurement, 608
GetVoltageValvesLimit	HeadStageIDType, 607
CRoboDeviceNet, 400	HeadstageType, 608
GetVolti	HeadstageTypeEnum, 607
CTcxDeviceNet, 509	ID, 608
GetWaveform	MeasuringOnly, 607
CTEERFunctionNet, 519	NumberOfAnalogChannels, 608
GetWaveLengthInNanometer	NumberOfStimulationChannels, 608
CMultiwellOptoStimFunctionNet, 317	OpticalStimulation, 607
GetWorkingFrequency	SN, 608
CRFFunctionNet, 375	StimulusParameters, 608
GetWPADebugMode	ToString, 608
CWClassicFunctionNet, 590	Type, 609
GetWPAType	TypeValue, 609
CWClassicFunctionNet, 590	Unknown, 607
GetXGain	UserDefinedName, 609
CRoboDacqNet, 388	Valid, 609
GetXilinxFlashOffset	W16IsW14, 609
CMcsUsbFactoryNet, 228	HeadstageIDTypeObject, 609
GetXilinxFlashReadCommand	_AdditionalText, 610
CMcsUsbFactoryNet, 228	IdType, 610
GND_SWITCH_BIT	AdditionalText, 610
CW2100_StimulatorFunctionNet, 541	Equals, 610
GroupID	GetHashCode, 610
CDeviceGroupChannelInfoTemplateNet< Dacq-	HeadstageIDTypeObject, 609
GroupChannelEnumTemplateNet >, 58	IdType, 610
GroupType	ToString, 610
CDeviceGroupChannelInfoTemplateNet< Dacq-	HeadStageIDTypeState, 610
GroupChannelEnumTemplateNet >, 58	ControlState, 611
Grouponamentminiemplateriet /, 30	DataState, 611
HasAnalogOut	IdType, 611
CSCUFunctionNet, 433	State, 611
HasGalvanicIsolation	HeadstageType
CSCUFunctionNet, 433	HeadStageIDType, 608
HasHSPowerSwitch	rieaustageib type, 000

HeadstageTypeEnum	CTEERFunctionNet, 519
HeadStageIDType, 607	CWarnerUssingFunctionNet, 556
HLADacq	IsPlateTypeValid
CHLADeviceNet, 112	CMultiwellDeviceNet, 312
HWInfo	IsPulseEnabled
CMcsUsbDacqNet, 204	CWarnerUssingFunctionNet, 556
HwVersion	IsPumpMotorOn
CMcsUsbListEntryNet, 238	CRoboFluidDeviceNet, 409
SMSSSSEISTERN J. TOO, 200	IsRunning
ID	CMeaCleanDeviceNet, 271
HeadStageIDType, 608	CMeaCoatDeviceNet, 276
IdProduct	IsSamplingFinished
DeviceIdNet, 595	• •
IdType	CTEERFunctionNet, 519
HeadstageIDTypeObject, 610	IsUserTriggerEnabled
HeadStageIDTypeState, 611	CLIH3DeviceNet, 126
IdVendor	IsValveDigitalInInverted
DeviceIdNet, 595	CWarnerValveControllerDeviceNet, 572
•	IsValveDigitaIInInvertedEvent
IsAnalogOutEnabled	CWarnerValveControllerDeviceNet, 584
CSCUFunctionNet, 434	IsValveOpen
IsAutomaticAnalogOut	CWarnerValveControllerDeviceNet, 573
CSCUFunctionNet, 434	IsValveOpenEvent
IsBusy	CWarnerValveControllerDeviceNet, 584
CPPCFunctionNet, 352	IsValveOpenInAnalogMode
IsChamberAvailable	CWarnerValveControllerDeviceNet, 573
CWarnerUssingFunctionNet, 555	IsValveOpenInAnalogModeEvent
IsChipPowered	CWarnerValveControllerDeviceNet, 584
CCMOSMea_FunctionNet, 43	IsValveOpenInDigitalMode
IsConnected	CWarnerValveControllerDeviceNet, 573
CMcsUsbNet, 253	IsValveOpenInDigitalModeEvent
IsDeviceHighSpeed	CWarnerValveControllerDeviceNet, 585
CMcsUsbNet, 253	GWaitier valveControllerDevicemen, 363
IsDeviceHighSpeedCapable	ListModeSendStart
CMcsUsbNet, 254	CStg200xBasicNet, 455
IsDeviceTypeOf	ListModeSendStop
CMcsUsbListNet, 241	CStg200xBasicNet, 455
IsDigitalChannelDedicated	
-	ListOfChangedTriggers
CMcsUsbDacqNet::CHWInfo, 114	StgStatusNet, 625
IsDigitalOutPortInverted	LoadPressure
CWarnerValveControllerDeviceNet, 572	CPPCFunctionNet, 354
IsDigitalOutPortInvertedEvent	LoadUserFirmware
CWarnerValveControllerDeviceNet, 584	CMcsUsbFactoryNet, 228, 229
IsEqual	LoadValveTable
CFilterCoefficientsNet, 68	CWarnerValveControllerDeviceNet, 574
IsExceptionsEnabled	LockPlateClamp
CMcsUsbNet, 254	CMultiwellDeviceNet, 312
IsGateFloating	
CCMOSMea_FunctionNet, 43	m_Bottom
IsHeadstageAvailable	CCMOSMeaDeviceNet::CRegionOfInterestRect,
CSCUFunctionNet, 434	369
IsHeadstageAvailableEvent	m_Left
CSCUFunctionNet, 437	CCMOSMeaDeviceNet::CRegionOfInterestRect,
IsHighCurrentMode	369
CWarnerUssingFunctionNet, 556	m_pGilsonDevice
IsHSPowered	CGilsonDeviceNet, 100
CSCUFunctionNet, 434	m_pMcsBus_MotorControlNet
IsInDacqLegacyMode	CRoboFluidDeviceNet, 410
CSCUFunctionNet, 435	m_pMcsUsb
IsInternalCalibrationFinished	CMcsUsbFunctionNet, 233
isinteniaivaliviativiii iliisileu	ONICSOSDI UNCUONINGL, 200

m_pMcsUsbFunction	FirmwareDestinationNames, 605
CMcsUsbFunctionNet, 233	MCSBUS4
m_pRoboFluidDevice	FirmwareDestinationNames, 605
CRoboFluidDeviceNet, 410	MCSBUS5
m_Right	FirmwareDestinationNames, 605
CCMOSMeaDeviceNet::CRegionOfInterestRect,	MCSBUS6
369	FirmwareDestinationNames, 605
m_Top	MCSBUS7
CCMOSMeaDeviceNet::CRegionOfInterestRect,	FirmwareDestinationNames, 605
369	MCSBUS8
Manufacturer	FirmwareDestinationNames, 605
CMcsUsbListEntryNet, 238	MCSBUS9
MaxBitNumber	FirmwareDestinationNames, 605
DigitalSource< digitalsourceenum >, 596	McsBus_MotorControl
DigitalSourceGeneral, 597	CPeristalticPumpDeviceNet, 332
MaxBitNumberStatic	CPPCDeviceNet, 347
DigitalSource< digitalsourceenum >, 596	CPPS_DeviceNet, 356
Mcs, 22	CRoboDeviceNet, 407
Mcs::Usb, 22	CRoboFluidDeviceNet, 410
DEVICE_NOT_FOUND, 27	McsBus_Sensor
enCMosMeaChipType, 27	CPPCDeviceNet, 347
EnSTG200x STATUS, 27	CPPS_DeviceNet, 356
nMos16LV, 27	McsBus VoltageMode
nMos32LV, 27	CFluidControlDeviceNet, 80
nMos36LN, 27	McsBus XY
nMos64LN, 27	CRoboDeviceNet, 404
NOT_CONNECTED, 27	McsBus ZI
OK, 27	CRoboDeviceNet, 404
OnChannelData, 27	McsUsbDeviceStateEvent
OnDeviceArrivalRemoval, 27	CMcsUsbDeviceStatePushFunctionNet, 222
OnError, 28	CMcsUsbDeviceStatePushNet, 223
OnMcsUsbDeviceState, 28	MCU1
OnMcsUsbDeviceStateCallback, 28	FirmwareDestinationNames, 605
OnMwPollStatus, 28	MeaAudioFunctionNet
OnStg200xDataHandler, 28	CMeaDeviceNet, 285
OnStg200xErrorHandler, 28	MeaDigitalDataFunctionNet
OnStgPollStatus, 28	CMeaDeviceNet, 285
OnUpdateFirmwareProgress, 28	MeaFeedbackFunctionNet
OnUpdateFirmwareStatusChange, 29	CMeaDeviceNet, 285
•	,
RoboStatusEventDelegate, 29	Measure CPathIdentDoviceNet 320
unknown, 27 McsBus	CPathIdentDeviceNet, 329 MeasureReservoir
	CPPCFunctionNet, 354
CPPCDeviceNet, 347	•
CPPS_DeviceNet, 356	MeasuringOnly
CRoboDeviceNet, 407	HeadStageIDType, 607
MCSBUS1	MeFunctionNet
FirmwareDestinationNames, 604	CMeaDeviceNet, 286
MCSBUS10	mkfilter
FirmwareDestinationNames, 604	mkfilterNet, 612
MCSBUS11	mkfilter_coef_in_one_set
FirmwareDestinationNames, 604	mkfilterNet, 612
MCSBUS12	mkfilter_highpass_coeff
FirmwareDestinationNames, 604	mkfilterNet, 612
MCSBUS13	mkfilter_highpass_frequency_from_coeff
FirmwareDestinationNames, 605	mkfilterNet, 612
MCSBUS2	mkfilter_highpass_frequency_from_k
FirmwareDestinationNames, 605	mkfilterNet, 612
MCSBUS3	mkfilter_highpass_k

mkfilterNet, 612	CDeviceGroupChannelInfoTemplateNet< Dacq-
mkfilter_MCS	GroupChannelEnumTemplateNet >, 59
mkfilterNet, 612, 613	NumberOfStimulationChannels
mkfilter_MCS_k	HeadStageIDType, 608
mkfilterNet, 613	NumCoefSets
mkfilter_normalize_coeffs_int	CCreateFilterNet, 49
mkfilterNet, 613	
mkfilter_normalize_coeffs_short	OK
mkfilterNet, 613	Mcs::Usb, 27
mkfilter_normalize_scale_coeffs_int	OnChannelData
mkfilterNet, 614	Mcs::Usb, 27
mkfilter_scale_coef_in_one_set	OnDeviceArrivalRemoval
mkfilterNet, 614	Mcs::Usb, 27
mkfilterNet, 611	OnError
mkfilter, 612	Mcs::Usb, 28
mkfilter_coef_in_one_set, 612	OnGetActiveRunningTableNumber
mkfilter_highpass_coeff, 612	CWarnerValveControllerDeviceNet, 574
mkfilter_highpass_frequency_from_coeff, 612	OnGetAnalogThresholdHigh
mkfilter_highpass_frequency_from_k, 612	CWarnerValveControllerDeviceNet, 574
mkfilter highpass k, 612	OnGetAnalogThresholdLow
mkfilter_MCS, 612, 613	CWarnerValveControllerDeviceNet, 574
mkfilter_MCS_k, 613	OnGetAnalogVoltage
mkfilter_normalize_coeffs_int, 613	CWarnerValveControllerDeviceNet, 574
mkfilter_normalize_coeffs_short, 613	OnGetAvailableHeadstages
mkfilter_normalize_scale_coeffs_int, 614	CSCUFunctionNet, 435
mkfilter_scale_coef_in_one_set, 614	OnGetCurrentNumberOfValves
MoveAbs	CWarnerValveControllerDeviceNet, 574
CRoboDeviceNet, 400	OnGetDigitalOutPortValve
MoveAbsI	CWarnerValveControllerDeviceNet, 574
CRoboStatorDeviceNet, 415, 416	OnGetDigitalPortDirection
MoveAbsXY	CWarnerValveControllerDeviceNet, 574
CRoboStatorDeviceNet, 416	OnGetDisplayMode
MoveAbsZ	CWarnerValveControllerDeviceNet, 575
CRoboStatorDeviceNet, 416	OnGetPlateClampStateByHeadstage
	CMultiwellCallbackFunctionNet, 307
MultibootGetCypressImageId	OnGetTableNamebyIndex
CMcsUsbNet, 254	CWarnerValveControllerDeviceNet, 575
MultibootGetImageId	OnGetValveActive
CMcsUsbNet, 254	CWarnerValveControllerDeviceNet, 575
MultibootGetSelectedImage	OnGetValveBoardRevision
CMcsUsbNet, 254	CWarnerValveControllerDeviceNet, 575
MultibootSelectImage	OnGetValveDigitalInPort
CMcsUsbNet, 254	CWarnerValveControllerDeviceNet, 575
MwPollStatusEvent	OnGetValveLedOn
CStg200xDownloadNet, 484	CWarnerValveControllerDeviceNet, 575
nMos16LV	OnGetValveManualGroup
Mcs::Usb, 27	CWarnerValveControllerDeviceNet, 575
nMos32LV	OnGetValveManualState
	CWarnerValveControllerDeviceNet, 575
Mcs::Usb, 27	OnGetValveMode
nMos36LN	
Mcs::Usb, 27	CWarnerValveControllerDeviceNet, 575
nMos64LN	OnlsDigitalOutPortInverted
Mcs::Usb, 27	CWarnerValveControllerDeviceNet, 576
NOT_CONNECTED	OnlsHeadstageAvailable
Mcs::Usb, 27	CSCUFunctionNet, 435
NumberOfAnalogChannels	OnlsValveDigitalInInverted
HeadStageIDType, 608	CWarnerValveControllerDeviceNet, 576
NumberOfChannels	OnlsValveOpen
	CWarnerValveControllerDeviceNet, 576

OnIsValveOpenInAnalogMode	PreCommaA
CWarnerValveControllerDeviceNet, 576	CFilterCoefficientsNet::s_FilterAttributesNet, 622
OnIsValveOpenInDigitalMode	PreCommaB
CWarnerValveControllerDeviceNet, 576	CFilterCoefficientsNet::s_FilterAttributesNet, 623
OnMcsUsbDeviceState	PrepareAndAppendData
Mcs::Usb, 28	CStg200xDownloadNet, 481
OnMcsUsbDeviceStateCallback	CStimulusFunctionNet, 491
Mcs::Usb, 28	PrepareAndSendData
OnMwPollStatus	CStg200xDownloadNet, 482
Mcs::Usb, 28	CStimulusFunctionNet, 492
OnStg200xDataHandler	PrepareChannelData
<u> </u>	CDigOutStimulatorFunctionNet, 62
Mcs::Usb, 28	•
OnStg200xErrorHandler	PrepareData
Mcs::Usb, 28	CStimulusFunctionNet, 492
OnStgPollStatus	CW2100_StimulatorFunctionNet, 540
Mcs::Usb, 28	PrepareDataSync
OnTableEntryChanged	CW2100_StimulatorFunctionNet, 540
CWarnerValveControllerDeviceNet, 576	Product
OnUpdateFirmwareProgress	CMcsUsbListEntryNet, 238
Mcs::Usb, 28	Program
OnUpdateFirmwareStatusChange	CProgramPressureCurveNet, 363
Mcs::Usb, 29	PulseGenerator
OpenPipe	CW2100 FunctionNet, 534
CGenericDevelopDeviceNet, 93	PumpOff
OpenPlateClamp	CRoboFluidDeviceNet, 409
CMultiwellDeviceNet, 312	PumpOn
operator=	CRoboFluidDeviceNet, 409
DeviceIdNet, 594	Criobor laidbeviocriet, 400
	QueryTriggerstatus
OpticalStimulation	CStg200xDownloadNet, 483
HeadStageIDType, 607	Cotg200xDownloadivet, 400
Order	RampStart
CCreateFilterNet, 49	COctoPotDeviceNet, 322
CFilterPropertyNet, 74	Read
PatternListStart	CExternDTesterDeviceNet, 65
COctoPotDeviceNet, 321	Read2
	CExternDTesterDeviceNet, 66
PIC	•
FirmwareDestinationNames, 605	ReadBlockFromFlash
PIC2	CMcsUsbFactoryNet, 229
FirmwareDestinationNames, 606	ReadBlockFromIFBGlobalEEprom
PIC3	CMcsUsbFactoryNet, 229
FirmwareDestinationNames, 606	ReadBlockFromNVMEM
PIC4	CMcsUsbFactoryNet, 229
FirmwareDestinationNames, 606	ReadClipping
PollStatusEvent	CLIH3DeviceNet, 127
CStimulusFunctionNet, 496	ReadEepromRegisterPreconfig
CW2100_StimulatorFunctionNet, 542	CMcsUsbNet, 255
PostCommaA	ReadPipe
CFilterCoefficientsNet::s_FilterAttributesNet, 622	CGenericDevelopDeviceNet, 93
PostCommaB	ReadRegister
CFilterCoefficientsNet::s_FilterAttributesNet, 622	CMcsUsbNet, 255
PowerChip	ReadRegister32
CCMOSMea_FunctionNet, 43	CMcsUsbNet, 255
PowerHS	ReadRegisterTimeSlot
	CMcsUsbNet, 255
CSCUFunctionNet, 435	ReadUARTData
PPCFunction CRRCRoviceNet 247	
CPPCDeviceNet, 347	CLIH3DeviceNet, 127
PPS_Function	Receive
CPPS_DeviceNet, 356	CSerialPortNet, 438

ReceiveString	RoboError_OverPressure
CSerialPortNet, 438, 439	CRoboDeviceNet, 406
RemoveSoftwareKey	RoboError_ParameterNotAllowed
CMcsUsbNet, 256	CRoboDeviceNet, 406
RescanHeadstage	RoboError_PeristalticTimeout
CMcsUsbNet, 256	CRoboDeviceNet, 406
ResetAdcOffset	RoboError_Phase0OutOfRange
COctoPotDeviceNet, 322	CRoboDeviceNet, 406
ResetChannelmap	RoboError_PollLoopCanceled
CWClassicFunctionNet, 591	CRoboDeviceNet, 406
ResetDacOffset	RoboError_PollLoopCanceledAndStopMovement
COctoPotDeviceNet, 322	CRoboDeviceNet, 406
ResetHighpassFilter	RoboError_Pressure
CFilterConfigurationNet, 70	CRoboDeviceNet, 406
ResetPipe	RoboError_RangeExceeded
CGenericDevelopDeviceNet, 93	CRoboDeviceNet, 406
ResetStatus	RoboError_StateChangeNotPossible
CStg200xDownloadBasicNet, 473	CRoboDeviceNet, 406
RFFunction	RoboError Timeout
CPositionIIDeviceNet, 343	CRoboDeviceNet, 406
RoboDacq	RoboError_UnknownCommand
CHiClampDeviceNet, 110	CRoboDeviceNet, 407
RoboDevice	RoboMainLowLevelCommand
CSafeISDeviceNet, 422	CRoboDeviceNet, 407
RoboError_AnotherMaster	RoboMainStatorLowLevelCommand
CRoboDeviceNet, 404	CRoboStatorDeviceNet, 419
RoboError_Base	RoboStatusEvent
CRoboDeviceNet, 404	CRoboDeviceNet, 407
RoboError_CannotEscapeEndSwitch	RoboStatusEventDelegate
CRoboDeviceNet, 404	Mcs::Usb, 29
RoboError_CommandAlreadyInProgress	RunTable
CRoboDeviceNet, 404	CRoboDacqNet, 388
RoboError_CommandNotPossible	On lobobacquet, ooo
CRoboDeviceNet, 404	s FilterAttributesNet
RoboError_CommunicationTimeout	CFilterCoefficientsNet::s_FilterAttributesNet, 622
CRoboDeviceNet, 404	SampleRate
RoboError_DacqNotReady	CCreateFilterNet, 50
CRoboDeviceNet, 404	Samplerate
RoboError_DLLMovementTimeout	CMcsUsbDacqNet, 221
CRoboDeviceNet, 405	Scale
RoboError_FindReferenceMethod	CCreateFilterNet, 50
CRoboDeviceNet, 405	ScanForHeadstages
RoboError_GilsonCommandPending	CWClassicFunctionNet, 591
	SelectHeadstage
CRoboDeviceNet, 405	CW2100 FunctionNet, 532
RoboError_GilsonTimeout	SelectTimeSlot
CRoboDeviceNet, 405	CW2100_StimulatorFunctionNet, 540
RoboError_GilsonWrondID	Send
CRoboDeviceNet, 405	CSerialPortNet, 439
RoboError_McsBus_UnknownCommand	SendBuffered
CRoboDeviceNet, 405	CGilsonDeviceNet, 100
RoboError_NoEndSwitch	SendChannelData
CRoboDeviceNet, 405	CDigOutStimulatorFunctionNet, 62
RoboError_NoMoreData	CStg200xDownloadBasicNet, 475
CRoboDeviceNet, 405	SendCommand
RoboError_NoReference	CLIH3DeviceNet, 127
CRoboDeviceNet, 405	SendImmediate
RoboError_NoSpeedOrAcceleration	CGilsonDeviceNet, 100
CRoboDeviceNet, 405	
	SendImmediateGetResponse

CGilsonDeviceNet, 100	SetAccelerationNativeZ
SendMultiplexedData	CRoboStatorDeviceNet, 416
CStimulusFunctionNet, 493	SetAccelerationXY
SendPreparedData	CRoboStatorDeviceNet, 417
CStimulusFunctionNet, 493	SetAccelerationZ
CW2100_StimulatorFunctionNet, 540	CRoboStatorDeviceNet, 417
SendSegmentDefine	SetAccelGyroDesiredRate
CStg200xDownloadNet, 483	CW2100_FunctionNet, 532
SendSegmentSelect	SetAccelGyroEnabled
CStg200xDownloadNet, 483	CW2100_FunctionNet, 532
SendSegmentStart	SetAccelRange
CStg200xDownloadNet, 484	CW2100_FunctionNet, 532
SendStart	SetActiveRunningTableNumber
CStg200xBasicNet, 455	CWarnerValveControllerDeviceNet, 576
CStimulusFunctionNet, 493	SetADC
CW2100_StimulatorFunctionNet, 540	CWarnerValveControllerDeviceTesterFunctionNet
SendStartDacq	587
CMcsUsbDacqNet, 204	SetAdcChannels
SendStartStgAndDacq	CSafeISDeviceNet, 420
CMcsUsbDacqNet, 204	SetADCInputOffset
SendStop	CCMOSMea_FunctionNet, 43
CStg200xBasicNet, 455	SetAdcOffset
CStimulusFunctionNet, 493	CLIH3DeviceNet, 127
CW2100_StimulatorFunctionNet, 541	COctoPotDeviceNet, 322
SendStopDacq	SetAdcOffsetPermanent
CMcsUsbDacqNet, 205	CLIH3DeviceNet, 128
SendStopStgAndDacq	SetAdcSamplePos
CMcsUsbDacqNet, 205	CSafeISDeviceNet, 420
SendStopStgAndDacqWithOptions	SetAirpressureLimit
CMcsUsbDacqNet, 205	CRoboDeviceNet, 400
SendSyncData	SetAirValve
CStg200xDownloadBasicNet, 475	CRoboDeviceNet, 401
Sensor	SetAllDigout
CFYIDeviceNet, 81	CRoboDacqNet, 388
CMeasureTableDeviceNet, 294	SetAmplificationSwitch
CPatchServerDeviceNet, 328	COctoPotDeviceNet, 322
SerialNumber	SetAmplitude
CMcsUsbListEntryNet, 238	CChannelTestDeviceNet, 33
CMcsUsbNet, 263	SetAmplitude_nA
SerialPort	CTEERFunctionNet, 519
CHLADeviceNet, 112	SetAnalogOutADCRange
Set4ADCCatchampAverageShift	CSCUFunctionNet, 435
CMcsBus_SensorNet, 160	SetAnalogOutChannel
Set4ADCMode	CW2100_FunctionNet, 533
CMcsBus SensorNet, 160	SetAnalogOutChannels
Set4DAC	CSCUFunctionNet, 436
CMcsBus SensorNet, 160	SetAnalogOutDACRange
Set Values	CSCUFunctionNet, 436
_	SetAnalogOutFilter
CNF_GenDeviceNet, 319 CPathIdentDeviceNet, 329	CW2100_FunctionNet, 533
SetAbsMaxCurrentInMicroAmp	SetAnalogThresholdHigh
•	CWarnerValveControllerDeviceNet, 577
CMultiwellOptoStimFunctionNet, 317	
SetAccelerationI	SetAnalogThresholdLow
CRoboStatorDeviceNet, 416	CWarnerValveControllerDeviceNet, 577
SetAccelerationNativel	SetAnalogVoltageRange
CRoboStatorDeviceNet, 416	CPPCFunctionNet, 354
SetAccelerationNativeXY	SetAnalogVoltages
CRoboStatorDeviceNet, 416	CPPS FunctionNet, 359

SetAttenuation	SetCheckVoltage
CChannelTestDeviceNet, 33	COkuvisionStimulatorDeviceNet, 326
SetAudioChannels	SetClampMode
CMeaAudioFunctionNet, 267, 268	CTEERFunctionNet, 520
CW2100_FunctionNet, 533	CWarnerUssingFunctionNet, 557
SetAudioOutDacParameter	SetColorRgb
CLIH3DeviceNet, 128	CMultiwellOptoStimFunctionNet, 317
SetAutocalibrationDisabled	SetColorStr
CStg200xBasicNet, 455	CMultiwellOptoStimFunctionNet, 318
SetAxisConfig	SetCommand
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBusNet, 172
618	CPedoterDeviceNet, 330
SetAxisLED	CRoboDacqNet, 388
CRoboocyte2DeviceNet, 413	SetConfiguration
SetAxisParametersEeprom	CMcsUsbNet, 256
CMcsBus_AxisParametersNet, 131, 132	SetConfigurationBit
SetBandwidthByIndex	CRoboDacqNet, 388
CIntanMea_FunctionNet, 116	SetConfigurationBitAxc
SetBaseFrequency	CRoboDacqNet, 388
CRFFunctionNet, 375	SetConfigurationBitBlu_Led
SetBaseSamplerate	CRoboDacqNet, 388
CCMOSMeaDeviceNet, 47	SetConfigurationBitBlu_LedToggleFast
SetBath	
	CRoboDacqNet, 389
CCMOSMea_FunctionNet, 43	SetConfigurationBitBlu_LedToggleSlow
SetBathclamp	CRoboDacqNet, 389
COctoPotDeviceNet, 322	SetConfigurationBitCC_Gen
SetBathMode	CRoboDacqNet, 389
CCMOSMea_FunctionNet, 43	SetConfigurationBitCV_Gen
SetBlankingEnable	CRoboDacqNet, 389
CStg200xBasicNet, 456, 457	SetConfigurationBitRC_Gen
SetBuffer	CRoboDacqNet, 389
CGenericDevelopDeviceNet, 94	SetConfigurationBitRed_Led
SetBufferIndex	CRoboDacqNet, 389
CTEERFunctionNet, 520	SetConfigurationBitRed_LedSaturation
SetBusAddress	CRoboDacqNet, 389
CMcsBusNet, 171	SetConfigurationBitRed_LedToggleFast
SetBusAddressEeprom	CRoboDacqNet, 389
CMcsBusNet, 171	SetConfigurationBitRed_LedToggleSlow
SetByteBuffer	CRoboDacqNet, 389
CGenericDevelopDeviceNet, 94	SetConfigurationBitRelais
SetCalibration	CRoboDacqNet, 389
CTcxDeviceNet, 509	SetConfigurationBitRV_Gen
SetCardinalDacqSamplerate	CRoboDacqNet, 390
CInterfaceboardFunctionNet, 120	SetConfigurationBitStream
SetCardinalStgOutputrate	CRoboDacqNet, 390
CInterfaceboardFunctionNet, 120	SetConfigurationBitSupply
SetChannel	CRoboDacqNet, 390
CSw2to64DeviceNet, 498	SetControllerParams
SetChannelmap	CTEERFunctionNet, 520
CWClassicFunctionNet, 591	SetCrossTalkOffset
SetChannels	CRoboDacqNet, 390
CSw2to64DeviceNet, 498	SetCrossTalkOptimum
SetChannelSwitch	CRoboDacqNet, 390
COctoPotDeviceNet, 322	SetCurrentAirvalveLimit
SetChargingMode	CRoboDeviceNet, 401
CMultiBatteryChargerDeviceNet, 304	SetCurrentAndAir
SetChargingPCoefficient CMultiPetters/ChargerPessionNet 204	CRoboDeviceNet, 401
CMultiBatteryChargerDeviceNet, 304	SetCurrentAndAirXY

CRoboStatorDeviceNet, 417	CPositionImpDeviceNet, 345
SetCurrentEditTableNumber	SetDeviceName
CWarnerValveControllerDeviceNet, 577	CUsbDeviceConfigurationFunctionNet, 524
SetCurrentEnable	SetDeviceType
CTEERFunctionNet, 520	CTcxDeviceNet, 510
SetCurrentFactor	SetDevname
COkuvisionStimulatorDeviceNet, 326	CTcxDeviceNet, 510
SetCurrentMode	SetDiagnosticMode
CStg200xBasicNet, 457	CIntanMea_FunctionNet, 116
SetCycles	SetDigitalData
CMeaCleanDeviceNet, 271	CMeaDigitalDataFunctionNet, 287
CMeaCoatDeviceNet, 276	SetDigitalOut
SetD	CMeaDeviceNet, 282
CTcxDeviceNet, 510	SetDigitalOutPortInvert
SetDacAmplificationFactor	CWarnerValveControllerDeviceNet, 577
CStg200xBasicNet, 457	SetDigitalOutPortValve
SetDacAutoControl	CWarnerValveControllerDeviceNet, 578
	SetDigitalPortDirection
COctoPotDeviceNet, 322 SetDacIdleValue	•
	CWarnerValveControllerDeviceNet, 578
CLIH3DeviceNet, 128	SetDigitalSource
SetDacMode CSetalSDavissNet 400	CMcsUsbDacqNet, 206–208
CSafeISDeviceNet, 420	SetDigitalStimulatorTrigger
SetDACOffset	CW2100_StimulatorFunctionNet, 541
CGrapheneFunctionNet, 106	SetDigitalStimulatorTriggerSlope
COkuvisionStimulatorDeviceNet, 326	CW2100_StimulatorFunctionNet, 541
SetDacOffset	SetDigout
CDacCalibrationFunctionNet, 51	CFluidControlDeviceNet, 78
CLIH3DeviceNet, 128	CRoboDacqNet, 390
COctoPotDeviceNet, 322	SetDigoutMode
SetDacOffsetPermanent	CStg200xBasicNet, 458
CLIH3DeviceNet, 129	SetDigOutState
SetDacPeriode	CLIH3DeviceNet, 129
CSafeISDeviceNet, 421	SetDigoutValue
SetDacPulseform	CStg200xBasicNet, 458
CSafeISDeviceNet, 421	SetDIO
SetDacqLegacyMode	CMcsBus_FYIExtensionNet, 134
CSCUFunctionNet, 436	SetDischargeCurrentSetPoint
SetDacRange	CMultiBatteryChargerDeviceNet, 305
CW2100_FunctionNet, 533	SetDisplayMode
SetDACs	CWarnerValveControllerDeviceNet, 578
CMcsBus_SensorNet, 160	SetDisplayText
SetDacUseIdleValue	CRoboDacqNet, 390
CLIH3DeviceNet, 129	SetDownsampleFactor
SetDacValue	CRoboDacqNet, 390
COctoPotDeviceNet, 323	SetDSPHighPassByIndex
SetDataMode	CIntanMea_FunctionNet, 117
CMcsUsbDacqNet, 206	SetDuration
SetDefault	CMeaCoatDeviceNet, 276
CWarnerValveControllerDeviceNet, 577	SetEEpromPage
SetDestinationSerialNumber	CLIH3DeviceNet, 129
CMcsUsbFactoryNet, 229	SetElectrodeDacMux
SetDetectionThreshold	CStg200xBasicNet, 458-460
CMcsBus_SensorNet, 160	SetElectrodeEnable
SetDevice	CStg200xBasicNet, 460–462
CTcxDeviceNet, 510	SetElectrodeMode
SetDeviceId	CStg200xBasicNet, 463, 464
CUsbDeviceConfigurationFunctionNet, 524	SetEnableAmplifierProtectionSwitch
SetDeviceList	CStg200xBasicNet, 464, 465

SetEnableHeaterLimit	SetHighpassFilterEnable
CTcxDeviceNet, 510	CFilterConfigurationNet, 70
SetEnablePulse	SetHWConfig
CWarnerUssingFunctionNet, 557	CRoboDeviceNet::RoboMainLowLevelCommands,
SetEnableThermocouple	618
CTcxDeviceNet, 510	SetHWRevision
SetExternalElectrodeEnable	CRoboDeviceNet::RoboMainLowLevelCommands,
CStg200xBasicNet, 466	618
SetExternalLED	SetHWRevisionEeprom
CTEERFunctionNet, 521	CMcsBusNet, 172
SetFAAmplification	SetHWSelectedChannels
CStg200xBasicNet, 467	CWClassicFunctionNet, 591
SetFilter	SetI
CRoboDacqNet, 390	CTcxDeviceNet, 510
SetFilterCoeffs	SetIClamp
CRoboDacqNet, 391	CRoboDacqNet, 391
SetFilterParameter	SetICoeff
CFilterConfigurationNet, 70	CRobo_FYITemp_FunctionNet, 379
CFilterConfigurationRegisterNet, 71, 72	SetICOffset
SetFilterParameterPermanent	CRoboDacqNet, 391
CFilterConfigurationNet, 70	SetIdleModeOffset
CFilterConfigurationRegisterNet, 72	CWarnerUssingFunctionNet, 558
SetFilterParametersHeadstage	SetIGain
CWClassicFunctionNet, 591	CRoboDacqNet, 391
SetFinalDischargeVoltage	SetImpedanceTestFrequency
CMultiBatteryChargerDeviceNet, 305	CMealmpedanceDeviceNet, 293
SetFrequency	SetImpId
CChannelTestDeviceNet, 33	CPositionImpDeviceNet, 345
CRadioControledDevicesNet, 368	SetImplantCurrentSetpoint
SetGain	CPositionIIDeviceNet, 340
CPgaDeviceNet, 334	SetInMovement
SetGate	CRoboDeviceNet, 401
CCMOSMea_FunctionNet, 43	SetIntanRegister
	-
SetGateFloating	CIntanMea_FunctionNet, 117 SetIntBuffer
CCMOSMea_FunctionNet, 44 SetGateToVOP	
	CGenericDevelopDeviceNet, 94
CCMOSMea_FunctionNet, 44	SetIO
SetGlobalRepeat	CWarnerValveControllerDeviceTesterFunctionNet,
CDigOutStimulatorFunctionNet, 63	587
SetGyroRange	SetIODirection
CW2100_FunctionNet, 533	CWarnerValveControllerDeviceTesterFunctionNet,
SetHasChecksum	587
CWClassicFunctionNet, 591	Setlo Voltage
SetHeadstage	CInterfaceboard2FunctionNet, 118
CStg200xBasicNet, 467	SetLatency
SetHeadstageOnOff	CMcsBus_SensorNet, 160
CW2100_FunctionNet, 533	SetLayoutConfiguration
CWClassicFunctionNet, 591	CMEA2100x256FunctionNet, 265
SetHeadstagePowerStateAtStart	SetLED
CSCUFunctionNet, 436	CRetinaLedDeviceNet, 370
SetHeadstageSamplingActive	SetLEDSwitch
CW2100_FunctionNet, 533	CMcsBus_ExtensionNet, 133
SetHeadstageToSleep	SetLength
CW2100_FunctionNet, 533	CRobo_FYIProgram_FunctionNet, 377
SetHeaterLimit	SetLiquidResistance
CTcxDeviceNet, 510	CTEERFunctionNet, 521
SetHighCurrentMode	CWarnerUssingFunctionNet, 558
CWarnerUssingFunctionNet, 557	SetListmodeIndexRange

00: 000 B : N : 407	0 11400 10 11
CStg200xBasicNet, 467	SetMCCurrentPosition
SetListmodeTriggerSource	CMcsBus_MotorControlNet, 147
CStg200xBasicNet, 467	SetMCCurrentShortCommand
SetLowCurrentMode	CMcsBus_MotorControlNet, 147
CWarnerUssingFunctionNet, 558	SetMCMaxAcceleration
SetLumi	CMcsBus_MotorControlNet, 147
CRetinaLedDeviceNet, 370	SetMCMaxAccelerationEeprom
SetMaxCurrent	CMcsBus_MotorControlNet, 147
CMeaCoatDeviceNet, 277	SetMCMaxCurrent
SetMaxDurationHighCurrentInMicroSec	CMcsBus_MotorControlNet, 147
CMultiwellOptoStimFunctionNet, 318	SetMCMaxCurrentEeprom
SetMaxDutyCycleHighCurrent	CMcsBus_MotorControlNet, 147
CMultiwellOptoStimFunctionNet, 318	SetMCMaxSpeed
•	•
SetMaxHeaterPowerMultiwell	CMcsBus_MotorControlNet, 148
CTcxDeviceNet, 511	SetMCMaxSpeedEeprom
SetMaxNoPressure	CMcsBus_MotorControlNet, 148
CRoboDeviceNet::RoboMainLowLevelCommands,	SetMCMaxTravel
618	CMcsBus_MotorControlNet, 148
SetMaxNoPressureWaitTime	SetMCMaxTravelEeprom
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBus_MotorControlNet, 148
618	SetMCMaxTravelShortCommand
SetMaxP	CMcsBus_MotorControlNet, 148
CTcxDeviceNet, 511	SetMCNewPosition
SetMaxPower	CMcsBus_MotorControlNet, 148
COkuvisionStimulatorDeviceNet, 327	SetMCOutputOnOff
CRobo_FYITemp_FunctionNet, 379	CMcsBus_MotorControlNet, 149
SetMaxPressureWaitTime	SetMCReference
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBus_MotorControlNet, 149
619	SetMCReferenceCurrent
SetMaxVoltage	CMcsBus_MotorControlNet, 149
CMeaCleanDeviceNet, 271	SetMCReferenceCurrentEeprom
	•
COkuvisionStimulatorDeviceNet, 327	CMcsBus_MotorControlNet, 149
SetMCAcceleration	SetMCRegulatorGain
CMcsBus_MotorControlNet, 144	CMcsBus_MotorControlNet, 149
SetMCAccelerationEeprom	SetMCRegulatorGainEeprom
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 149
SetMCAccelerationShortCommand	SetMCRotation
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 150
SetMCAxisRevisionEeprom	SetMCScalingFactor
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 150
SetMCBreakCurrent	SetMCScalingFactorEeprom
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 150
SetMCBreakCurrentEeprom	SetMCSpeed
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 150
SetMCConfig	SetMCSpeedEeprom
CMcsBus_MotorControlNet, 145	CMcsBus_MotorControlNet, 150
SetMCConfigEeprom	SetMCSpeedShortCommand
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 150
SetMCCurrent	SetMCSpeedUnitEeprom
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 151
	SetMCStandbyCurrent
SetMCCurrentEeprom CMacRus MaterControlNet 146	
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 151
SetMCCurrentMode	SetMCStandbyCurrentEeprom
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 151
SetMCCurrentModeEeprom	SetMCStandbyTime
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 151
SetMCCurrentModeShortCommand	SetMCStandbyTimeEeprom
CMcsBus_MotorControlNet, 146	CMcsBus_MotorControlNet, 151

SetMeasurementMode	CPulseGeneratorFunctionNet, 366
CStg200xBasicNet, 467	SetPeriod_us
SetMinimalThreshold	CTEERFunctionNet, 521
CMcsBus_SensorNet, 160	SetPermanentCurrentInMicroAmp
SetMinNoPressureWaitTime	CMultiwellOptoStimFunctionNet, 318
CRoboDeviceNet::RoboMainLowLevelCommands,	SetPersistency
619	CRetinaLedDeviceNet, 370
SetMinPressure	SetPGain
CRoboDeviceNet, 401	CRoboDacqNet, 391
CRoboDeviceNet::RoboMainLowLevelCommands,	SetPidParameter
619	COctoPotDeviceNet, 323
SetMinPressureWaitTime	SetPiezoState
CRoboDeviceNet::RoboMainLowLevelCommands,	CMcsBus_SensorNet, 161
619	SetPlateMux
SetMinVoltage	CMultiwellDeviceNet, 312, 313
CMeaCleanDeviceNet, 272	SetPlateType
SetModeSelect	CMultiwellDeviceNet, 313
CPulseGeneratorFunctionNet, 366	SetPoti
SetMovePump	CMcsUsbDacqNet, 208
CMcsBus_SensorNet, 161	SetPowerMuxPlate
SetMultiHeadstageMode	CMultiwellDeviceNet, 313
CW2100_FunctionNet, 534	SetPowerStrength
SetNanoVoltsPerKelvin	CPositionIIDeviceNet, 340
CMcsBus_TempSensorNet, 164	SetPressureOffset
SetNeurochipMemoryData	CMcsBus_SensorNet, 161
CCMOSMea_FunctionNet, 44	CPPCFunctionNet, 354
SetNoFilterCoeffs	SetPressureRange
CRoboDacqNet, 391	CPPCFunctionNet, 354
SetNumberOfAnalogChannels	SetPulse
CMeaDeviceNet, 282	CWarnerUssingFunctionNet, 558
SetNumberOfChannels	SetPulseform
CMeaDeviceNet, 283, 284	COkuvisionStimulatorDeviceNet, 327
COctoPotDeviceNet, 323	SetPulseLength
SetOffsetCurrent	CPulseGeneratorFunctionNet, 366
CMeaCoatDeviceNet, 277	SetPumpCouple
SetOnOff	CPPS_FunctionNet, 360
CTcxDeviceNet, 511	SetPumpEnableSpeedRatio
SetOutputMap	CPPS_FunctionNet, 360
CStg200xDownloadNet, 484	SetPumpFastOnOff
SetOutputRate	CPPS FunctionNet, 360
COctoPotDeviceNet, 323	SetPumpFastSpeed
CStg200xBasicNet, 468	CPPS_FunctionNet, 360
SetP	SetPumpFunctionSpeeds
CTcxDeviceNet, 511	CPPS_FunctionNet, 360
SetParameter	SetPumpManualOnOff
CRoboDeviceNet::RoboMainLowLevelCommands,	CPPS_FunctionNet, 360
619	SetPumpMaxSpeed
SetPattern	CPPS_FunctionNet, 360
CMeaSwitchDeviceNet, 296	SetPumpModeType
SetPatternBool	CPPCFunctionNet, 355
CMeaSwitchDeviceNet, 296	CPPS_FunctionNet, 360
SetPatternListEntry	SetPumpSpeed
COctoPotDeviceNet, 323	CRoboFluidDeviceNet, 409
SetPauseDuration	SetPumpSpeedRatio
CMeaCoatDeviceNet, 277	CPPS_FunctionNet, 361
SetPCoeff	SetPumpSpeedUnit
CRobo_FYITemp_FunctionNet, 379	CPPCFunctionNet, 355
SetPeriod	CPPS_FunctionNet, 361

SetPWM	CRoboDeviceNet::RoboMainLowLevelCommands
CFluidControlDeviceNet, 78	619
SetRampParameter	SetSearchReferenceFineAccel
COctoPotDeviceNet, 323	CRoboDeviceNet::RoboMainLowLevelCommands
SetRatedCapacity	620
CMultiBatteryChargerDeviceNet, 305	SetSearchReferenceFineSpeed
SetRatedCapacityVolatile	CRoboDeviceNet::RoboMainLowLevelCommands
CMultiBatteryChargerDeviceNet, 305	620
SetRecordingNumber	SetSearchReferenceMethod
CRoboDacqNet, 391	CRoboDeviceNet::RoboMainLowLevelCommands
SetReferenceElectrodeMode	620
CSCUFunctionNet, 437	SetSearchReferenceMoveOut
SetReferenceElectrodeSwitchState	CRoboDeviceNet::RoboMainLowLevelCommands
CSCUFunctionNet, 437	620
	SetSearchReferenceOffsetPos
SetRegionOfInterests	
CCMOSMeaDeviceNet, 47	CRoboDeviceNet::RoboMainLowLevelCommands
SetRegulationTimeouts	620
CMcsBus_SensorNet, 161	SetSelectedChannels
SetRegulatorFactor	CMcsUsbDacqNet, 208–210
CMcsBus_SensorNet, 161	CW2100_FunctionNet, 534
SetRegulatorOnOff	SetSelectedChannelsQueue
CMcsBus_SensorNet, 161	CMcsUsbDacqNet, 211–213
CRobo_FYITemp_FunctionNet, 379	SetSelectedData
SetRepeat	CMcsUsbDacqNet, 213–215
CRetinaLedDeviceNet, 371	SetSelectedHeadstage
SetRepeats	CWClassicFunctionNet, 592
CProgramPressureCurveNet, 363	SetSensorType
SetResetFilter	CTcxDeviceNet, 511
CWClassicFunctionNet, 591	SetSerialNumberHeadstage
SetRFFrequency	CWClassicFunctionNet, 592
CPositionImpDeviceNet, 346	SetSetpoint
SetRFFrequencyHeadstage	CTcxDeviceNet, 511
CWClassicFunctionNet, 591	SetShortBuffer
SetRFFrequencyReceiver	CGenericDevelopDeviceNet, 95
CWClassicFunctionNet, 592	SetSimulation
SetRFFrequencyReceiverEeprom	CRoboDacqNet, 391
CWClassicFunctionNet, 592	SetSineParameter
SetRFLostBehaviour	COctoPotDeviceNet, 323
CWClassicFunctionNet, 592	SetSingleHeater
SetRFPower	CMcsBus FYIExtensionNet, 134
CWClassicFunctionNet, 592	SetSingleValve
SetRotatePump	CFluidControlDeviceNet, 79
•	
CMcsBus_SensorNet, 161 SetRTC	CRoboFluidDeviceNet, 410
	SetSlope
COkuvisionStimulatorDeviceNet, 327	CMeaCleanDeviceNet, 272
CPositionIIDeviceNet, 340	CMeaCoatDeviceNet, 277
SetSampleInterval	SetSoftwareKey
CLIH3DeviceNet, 130	CMcsUsbNet, 256
SetSamplePeriode	SetSollPressure
CMcsBus_SensorNet, 162	CMcsBus_SensorNet, 162
SetSamplerate	SetSollTemp
CMcsUsbDacqNet, 208	CRobo_FYITemp_FunctionNet, 379
SetScreen	SetSourceBulk
CRoboDacqNet, 391	CCMOSMea_FunctionNet, 44
SetSearchReferenceFastAccel	SetSourceDrain
CRoboDeviceNet::RoboMainLowLevelCommands,	CCMOSMea_FunctionNet, 44
619	SetSourceGate
SetSearchBeferenceFastSneed	CCMOSMea FunctionNet 44

SetSpeedI	SetTriggerSource
CRoboStatorDeviceNet, 417	CStg200xBasicNet, 468, 469
SetSpeedNativel	SetTriggerSyncDirection
CRoboStatorDeviceNet, 417	CWarnerValveControllerDeviceTesterFunctionNet,
SetSpeedNativeXY	588
CRoboStatorDeviceNet, 417	SetUByteBuffer
SetSpeedNativeZ	CGenericDevelopDeviceNet, 96
CRoboStatorDeviceNet, 417	SetUClamp
SetSpeedXY	CRoboDacqNet, 392
CRoboStatorDeviceNet, 417	SetUCOffset
,	
SetSpeedZ	CRoboDacqNet, 392
CRoboStatorDeviceNet, 418	SetUIntBuffer
SetStartTriggerSlope	CGenericDevelopDeviceNet, 96
CDigOutStimulatorFunctionNet, 63	SetupGroupDacqQueue
SetStateDebugData	CMcsUsbDacqNet, 215
CPositionIIDeviceNet, 342	SetupRetriggerMode
SetStateEventData	CStg200xDownloadBasicNet, 476
CPositionIIDeviceNet, 342	SetupTrigger
SetStgProgramInfo	CStg200xDownloadBasicNet, 476
CStg200xBasicNet, 468	CStimulusFunctionNet, 494
SetStimulusSites	SetupTriggerSingle
CCMOSMea_FunctionNet, 44	CStg200xDownloadBasicNet, 477
SetStopTriggerSlope	CStimulusFunctionNet, 494
CDigOutStimulatorFunctionNet, 63	SetUseBubble
SetStringFormat	CPPS_FunctionNet, 361
CMcsUsbListEntryNet, 237	SetUserParameter
CMcsUsbListNet, 241	CRoboDeviceNet::RoboMainLowLevelCommands,
SetSubChannel	620, 621
CMcsBus_MotorControlNet, 151	SetUShortBuffer
SetSwitches	CGenericDevelopDeviceNet, 97
CSafeISDeviceNet, 421	SetUVOffset
SetSyncoutMap	CRoboDacqNet, 392
CStg200xBasicNet, 468	SetValue
SetTableName	CGenericDevelopDeviceNet, 97
CWarnerValveControllerDeviceNet, 578	SetValve
SetTablepointer	CFluidControlDeviceNet, 79
CRetinaLedDeviceNet, 371	CRoboFluidDeviceNet, 410
SetTableStep	SetValve1
CWarnerValveControllerDeviceNet, 579	CRobo_FYIProgram_FunctionNet, 377
SetTableStepAll	SetValve2
CWarnerValveControllerDeviceNet, 579	CRobo_FYIProgram_FunctionNet, 377
SetTestMode	SetValveActive
CRFFunctionNet, 375	CPPCFunctionNet, 355
SetThermocoupleNanovoltPerKelvin	CWarnerValveControllerDeviceNet, 579
CFluidControlDeviceNet, 79	SetValveDigitalInInvert
CTcxDeviceNet, 512	CWarnerValveControllerDeviceNet, 579
SetThermoOffset	SetValveDigitalInPort
CMcsBus_TempSensorNet, 164	CWarnerValveControllerDeviceNet, 580
SetTransformer	SetValveLedOn
CMeFunctionNet, 298	CWarnerValveControllerDeviceNet, 580
SetTrigger	SetValveManualGroup
	CWarnerValveControllerDeviceNet, 580
CRetinaLedDeviceNet, 371	SetValveManualState
CWarnerValveControllerDeviceTesterFunctionNet,	
587	CWarnerValveControllerDeviceNet, 580
SetTriggerMaskValue	SetValveMode
CMeaDeviceNet, 284	CWarnerValveControllerDeviceNet, 581
SetTriggerPeriod	SetValves
CMeaDeviceNet, 285	CMcsBus FYIExtensionNet, 134

SetValvesActiveMap	SetVoltageRangeInMicroVolt
CWarnerValveControllerDeviceNet, 581	CMcsUsbDacqNet, 215
SetValvesManualStateMap	SetVoltageResolution
CWarnerValveControllerDeviceNet, 581	CGrapheneFunctionNet, 109
SetValveTableEntry	SetVoltageRs485ALimit
CWarnerValveControllerDeviceNet, 581	CRoboDeviceNet, 402
SetVds	SetVoltageRs485BLimit
CGrapheneFunctionNet, 106	CRoboDeviceNet, 402
SetVdVs	SetVoltageValvesLimit
CGrapheneFunctionNet, 107	CRoboDeviceNet, 402
SetVdVsDAC	SetWaveform
CGrapheneFunctionNet, 107	CChannelTestDeviceNet, 34
SetVelocityI	CTEERFunctionNet, 521
CRoboStatorDeviceNet, 418	SetWaveLengthInNanometer
SetVelocityXY	CMultiwellOptoStimFunctionNet, 319
CRoboStatorDeviceNet, 418	SetWorkingFrequency
	- · · · · · · · · · · · · · · · · · · ·
SetVelocityZ	CRFFunctionNet, 375
CRoboStatorDeviceNet, 418	SetWPADebugMode
SetVgs	CWClassicFunctionNet, 592
CGrapheneFunctionNet, 108	SetWPAType
SetVMMaxNegativeCurrent	CWClassicFunctionNet, 592
CMcsBus_VoltageModeNet, 167	SetXGain
SetVMMaxNegativeCurrentEeprom	CRoboDacqNet, 392
CMcsBus_VoltageModeNet, 167	Sideband
SetVMMaxNegativeVoltage	CStimulusFunctionNet::SidebandData, 624
CMcsBus_VoltageModeNet, 167	SidebandData
SetVMMaxNegativeVoltageEeprom	CStimulusFunctionNet::SidebandData, 623
CMcsBus_VoltageModeNet, 167	SineStart
SetVMMaxPositiveCurrent	COctoPotDeviceNet, 324
CMcsBus_VoltageModeNet, 168	size
SetVMMaxPositiveCurrentEeprom	DigitalSource< digitalsourceenum >, 596
CMcsBus_VoltageModeNet, 168	DigitalSourceGeneral, 597
SetVMMaxPositiveVoltage	SN
CMcsBus_VoltageModeNet, 168	HeadStageIDType, 608
SetVMMaxPositiveVoltageEeprom	Source
CMcsBus_VoltageModeNet, 168	DigitalSource< digitalsourceenum >, 596
SetVMOutputOnOff	DigitalSourceGeneral, 597
CMcsBus_VoltageModeNet, 168	Start
SetVMVoltage	CMeaCleanDeviceNet, 272
CMcsBus_VoltageModeNet, 168	CMeaCoatDeviceNet, 277
SetVoltage12VLimit	CRobo_FYIProgram_FunctionNet, 377
CRoboDeviceNet, 401	StartDacq
SetVoltage5VLimit	CMcsUsbDacqNet, 216, 217
-	StartInternalCalibration
CRoboDeviceNet, 401	
SetVoltageAirvalveLimit	CTEERFunctionNet, 522
CRoboDeviceNet, 402	StartLoop
SetVoltageClampControllerParam_D	CMcsUsbDacqNet, 218, 219
CWarnerUssingFunctionNet, 559	StartMCMovement
SetVoltageClampControllerParam_I	CMcsBus_MotorControlNet, 152
CWarnerUssingFunctionNet, 559	StartMeasurement
SetVoltageClampControllerParam_P	CMealmpedanceDeviceNet, 293
CWarnerUssingFunctionNet, 559	StartPoll
SetVoltageMode	CStimulusFunctionNet, 496
CStg200xBasicNet, 469	CW2100_StimulatorFunctionNet, 541
SetVoltageRange	StartSampling
CGrapheneFunctionNet, 108	CTEERFunctionNet, 522
SetVoltageRangeByIndex	StartSync
CMcsUsbDacqNet, 215	CMcsBus_SensorNet, 162

State	Status_InvalidUrbFunction
HeadStageIDTypeState, 611	CMcsUsbNet, 261
Status	Status_IoPending
CUsbExceptionNet, 526	CMcsUsbNet, 261
Status_AlreadyConfigured	Status_loTimeout
CMcsUsbNet, 258	CMcsUsbNet, 261
Status_BadStartFrame	Status_IsochRequestFailed
CMcsUsbNet, 259	CMcsUsbNet, 261
Status_Btstuff	Status_LastUsbErrorMismatch
CMcsUsbNet, 259	CMcsUsbNet, 261
Status_BufferOverrun	Status_NoBandwidth
CMcsUsbNet, 259	CMcsUsbNet, 261
Status_BufferUnderrun	Status_NoMemory
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_Canceled	Status_NoSuchDevice
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_Canceling	Status_NotAccessed
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_ConnectedPipes	Status_NotSupported
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_ControlNotOwned	Status_PidCheckFailure
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_Crc	Status_PipeNotLinked
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_DataOverrun	Status_RequestFailed
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_DataToggleMismatch	Status_RequestMutexFailed
CMcsUsbNet, 259	CMcsUsbNet, 262
Status_DataUnderrun	Status_RequestMutexTimeout
CMcsUsbNet, 260	CMcsUsbNet, 262
Status_DeviceLocked	Status_Stall
CMcsUsbNet, 260	CMcsUsbNet, 262
Status_DeviceNotFound	Status_Unconfigured
CMcsUsbNet, 260	CMcsUsbNet, 262
Status_DeviceRemoved	Status_UnexpectedPid
CMcsUsbNet, 260	CMcsUsbNet, 263
Status_DevNotResponding	Stg200xPollStatusEvent
CMcsUsbNet, 260	CStg200xDownloadNet, 484
Status_EndpointHalted	StgStatusNet, 624
CMcsUsbNet, 260	FromIntPtr, 624 FromPtr, 624
Status_ErrorBusy	
CMcsUsbNet, 260 Status_ErrorShortTransfer	ListOfChangedTriggers, 625 TiggerStatus, 625
CMcsUsbNet, 260	StillConnected
Status_Fifo	CRadioControledDevicesNet, 368
CMcsUsbNet, 260	Stimulator
Status_FrameControlOwned	CW2100_FunctionNet, 534
CMcsUsbNet, 260	Stimulus
Status_InternalHcError	CCMOSMeaDeviceNet, 47
CMcsUsbNet, 260	CStg200xDownloadBasicNet, 478
Status_InvalidDeviceHandle	StimulusDeviceDataAndUnrolledData
CMcsUsbNet, 261	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData,
Status_InvalidHandle	625
CMcsUsbNet, 261	StimulusFunction
Status_InvalidParameter	CLIH3DeviceNet, 130
CMcsUsbNet, 261	StimulusParameters
Status_InvalidPipeHandle	HeadStageIDType, 608
CMcsUsbNet, 261	Stop
511100000110t, <u>201</u>	p

CMeaCleanDeviceNet, 272	CMcsUsbListEntryNet, 238
CMeaCoatDeviceNet, 278	HeadStageIDType, 608
StopDacq	HeadstageIDTypeObject, 610
CMcsUsbDacqNet, 219, 220	TriggerStatus
StopLoop	CMcsUsbDeviceStatePushFunctionNet, 222
CMcsUsbDacqNet, 220	CMcsUsbDeviceStatePushNet, 223
StopMCMovement	TxnGetSerialNumber
CMcsBus_MotorControlNet, 152	CMcsUsbNet, 256
StopMovement	TxnSetSerialNumber
CRoboDeviceNet, 402	CMcsUsbNet, 256
StopMovementI	TxnTestMemoryReadAndCheck
CRoboStatorDeviceNet, 418	CMcsUsbNet, 256
StopMovementXY	TxnTestMemoryWrite
CRoboStatorDeviceNet, 418	CMcsUsbNet, 257
StopMovementZ	Type
CRoboStatorDeviceNet, 418	HeadStageIDType, 609
StopPlateClamp	TypeValue
CMultiwellDeviceNet, 314	HeadStageIDType, 609
StopPoll	rieadotageib type, 000
CStimulusFunctionNet, 496	Unknown
	HeadStageIDType, 607
CW2100_StimulatorFunctionNet, 541	unknown
StopSampling CTETPE wasting Net 500	Mcs::Usb, 27
CTEERFunctionNet, 522	UnlockPlateClamp
StopTable	
CRoboDacqNet, 392	CMultiwellDeviceNet, 314 UnrolledAmplitude
StoreValveTable	
CWarnerValveControllerDeviceNet, 582	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData
SwitchOnOff	626
CPositionIIDeviceNet, 342	UnrolledDuration
SYNC_BIT0	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData
CW2100_StimulatorFunctionNet, 542	626
SYNC_BIT1	UnrolledSync
CW2100_StimulatorFunctionNet, 542	CStimulusFunctionNet::StimulusDeviceDataAndUnrolledData 626
Table_Wait	UpdateChannelBlock
CRoboDacqNet, 392	CCMOSMeaDeviceNet, 47
TableDefBegin	UpdateDisplay
CRoboDacqNet, 392	CRoboDacqNet, 392
TableDefEnd	UpdateFirmware
CRoboDacqNet, 392	CMcsUsbFactoryNet, 229-231
TableEntryChangedEvent	UpdateTransistorVoltages
CWarnerValveControllerDeviceNet, 585	CCMOSMea_FunctionNet, 44
TactSwitchGetState	USB
CMcsBus_SensorNet, 162	FirmwareDestinationNames, 606
TactSwitchSetDisplay	usbSetupPacket_t, 626
CMcsBus_SensorNet, 162	bmRequestType, 626
TEERFunctionNet	bRequest, 626
CTEERMachineDeviceNet, 523	wIndex, 627
ThrowCUsbExceptionNetOnError	wLength, 627
CMcsUsbFunctionNet, 233	wValue, 627
CMcsUsbNet, 256	UserDefinedName
TiggerStatus	HeadStageIDType, 609
StgStatusNet, 625	J 71 /
TimeResolutionInNanoSeconds	Valid
W2100_StimulusParametersNet, 627	HeadStageIDType, 609
ToCpp	ValidKey
CFilterCoefficientsNet::s_FilterAttributesNet, 622	CMcsUsbNet, 257
ToString	Voltage
CFilterPropertyNet, 73	BatteryState, 30

```
VoltageRangeDisplayStringMilliVolt
                                                    wValue
    CMcsUsbDacqNet::CHWInfo::CVoltageRangeInfoNet,
                                                         usbSetupPacket_t, 627
         526
VoltageRangeInMicroVolt
    CMcsUsbDacqNet:: CHWInfo:: CVoltageRangeInfoNet,\\
    W2100 StimulusParametersNet, 627
VoltageResolutionInMicroVolt
    W2100_StimulusParametersNet, 627
VoltageString
    BatteryState, 30
VOPSTimerSetResetTimes
    CCMOSMea_FunctionNet, 44, 45
W16lsW14
    HeadStageIDType, 609
W2100_FunctionNet
    CMeaDeviceNet, 286
W2100_StimulusParametersNet, 627
    CurrentRangeInNanoAmp, 627
    CurrentResolutionInNanoAmp, 627
    DACResolution, 627
    TimeResolutionInNanoSeconds, 627
    VoltageRangeInMicroVolt, 627
    VoltageResolutionInMicroVolt, 627
WaitForAllChambers
    CWarnerUssingFunctionNet, 560
WaitForChamber
    CWarnerUssingFunctionNet, 560
WarnerUssingFunction
    CWarnerUssingDeviceNet, 544
WClassicFunctionNet
    CMeaDeviceNet, 286
wIndex
    usbSetupPacket t, 627
wLength
    usbSetupPacket_t, 627
WPAError_ScanningIsPending
    CMcsUsbNet, 263
Write
    CExternDTesterDeviceNet, 66
Write2
    CExternDTesterDeviceNet, 66
WriteEepromRegisterPreconfig
    CMcsUsbNet, 257
WritePipe
    CGenericDevelopDeviceNet, 98
WriteRegister
    CMcsUsbNet, 257, 258
WriteRegister32
    CMcsUsbNet, 258
WriteRegisterArray
    CMcsUsbNet, 258
WriteRegisterTimeSlot
    CMcsUsbNet, 258
WriteRegisterValue
    CMcsUsbNet, 258
WriteUARTData
```

CLIH3DeviceNet, 130